CSR WRAP

EN CSR Wrap Single Ply

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EN PACKAGING INSTRUCTIONS

Wrap your medical device set in double layers using a sequential folding (Figure 1) according to your standard procedure.

Double sequential envelope folding is recommended to ensure optimal bacterial barrier efficiency and aseptic opening at the point of use.

STERILIZATION INSTRUCTIONS

When loading the autoclaves, comply with the autoclaves manufacturer's instructions.

It is critical that process parameters are validated for each individual type of sterilization equipment and load configuration.

1. Steam sterilization (moist heat)

In compliance with EN 285 and EN ISO 17665.

Sterilization temperature 269.6°F (132°C) may be used with plateau duration of 6 min.

2. EO Sterilization

In compliance with EN 1422 and EN ISO 11135.

All of our wraps have been validated under the following typical sterilization cycle conditions:

- EO Concentration: 725 mg/L
- Temperature: 130°F/54°C
- Exposure Time: 90 minutes
- Humidity: 60%
- Aeration: 12 hours

For a different sterilization process please refer to your local regulation and/or contact your distributor. It remains the responsibility of the processor to ensure that the processing is actually performed using equipment, materials and trained personnel in the processing facility to achieve the desired result.

Please keep these instructions with product for reference.

Figure 1 - Double sequential envelope folding

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