

CRITERION™ NITRILE EXAM N250 GLOVES

NEXT-LEVEL COMFORT AND PERFORMANCE

ULTRATHIN FORMULATION IMPROVES
TACTILE SENSITIVITY AND DEXTERITY



SAME GREAT PRODUCTS, NEW PACKAGING



For more information about the full Criterion line of personal protection equipment scan here or visit HENRYSCHEINDENTAL.COM/CRITERION

CRITERION[™] NITRILE EXAM N250 GLOVES

- Ultrathin formulation enhances tactile sensitivity and dexterity
- Chlorinated for reduced friction and easy donning
- Textured fingers provide enhanced grip
- Beaded cuff. 2.0-mil/0.05mm thickness at palm



LILAC COLOR

Small	250/box	570-1947
Medium	250/box	570-1948
Large	250/box	570-1949
X-Large	230/box	570-1950


CHEMOTHERAPY TESTED TO ASTM D6978 STANDARD

Caution: gloves used for protection against chemotherapy drug exposure must be selected specifically for the type of chemicals used. This glove has been tested for resistance to permeation of various chemotherapy drugs per ASTM D6978, "standard practice for assessment of resistance of medical gloves to permeation by chemotherapy drugs."

CHEMOTHERAPY DRUGS & CONCENTRATION	MINIMUM BREAKTHROUGH DETECTION TIME
Carmustine (BCNU) (3.3mg/mL)	14.7min.
Cisplatin (1.0mg/mL)	>240min
Cyclophosphamide (Cytoxan) (20mg/mL)	>240min
Cytarabine (100mg/mL)	>240min
Dacarbazine (DTIC) (10.0mg/mL)	>240min
Doxorubicin Hydrochloride (2.0mg/mL)	>240min
Etoposide (Toposar) (20.0mg/mL)	>240min
Fluorouracil (50.0mg/mL)	>240min
Ifosfamide (50.0mg/mL)	>240min
Methotrexate (25mg/mL)	>240min
Mitomycin C (0.5mg/mL)	>240min
Mitoxantrone (2.0mg/mL)	>240min
Paclitaxel (Taxol) (6.0mg/mL)	>240min
Thiotepa (10.0mg/mL)	58.8min
Vincristine Sulfate (1.0mg/mL)	>240min

Warning: Do not use Carmustine (BCNU) (3.3mg/mL) and Thiotepa (10mg/mL)
When chemotherapy drugs are present, gloves selection should be based on the specific types of drugs used. Users should review drug labeling or material safety data sheets for the drugs being used to determine an adequate level of protection.

SAME GREAT PRODUCTS, NEW PACKAGING




1. UPDATED PACKAGING

Makes it easy to know what's inside the box with attractive imagery and bold design elements.



2. INFO ZONE

Appears in the bottom right of all packaging. The icons indicate key product information: PPE standards, size, color and quantity.



3. REORDER PANEL

Helps to simplify the reordering process by clearly stating the full product name and reorder item number.