



Flex^{VeinViewer} & AESTHETICS

GREATER SUCCESS WITH THE FIRST PROCEDURE.

HOW IT WORKS

VeinViewer[®] uses near-infrared (NIR) light and other technologies to detect subcutaneous blood and create a digital image of the patient's superficial vein pattern projected directly onto the skin's surface in real time.

The benefits of VeinViewer could improve the success of your practice.

- See veins up to 10mm deep, visualizing the procedure as it happens.
- Fully assess the patient and identify "invisible" deeper feeder veins that are the root cause of the problem.
- Visualize the procedure as it happens; watching fluids clear the blood column.
- Potentially avoid an infiltration.
- Potentially increase the chance of success of the procedure and avoid additional treatments.
- Potentially reduce costs and improve patient throughput.
- Expand marketing opportunities with increased patient satisfaction and word of mouth referrals.
- Use during consultation to help educate patients on the procedure.

VeinViewer **Flex** & AESTHETICS



EVIDENCE: AESTHETICS MARKET

Feeder vein targets marked for phlebectomy using naked-eye visualization and palpation versus VeinViewer® only.

Total Naked Eye Marks	Total VeinViewer® Marks	Total Combined Marks	% of Marks Made With VeinViewer® Only
103	211	314	67%

Peer Review Journal, Aug. 2006, Dermatologic Surgery, Vein Imaging: A new method of near-infrared imaging, when a processed image is projected onto the surface of the skin for the enhancement of vein treatment.

SPECIFICATIONS

Depth of Visualization	Up to 10mm
Brightness	6 lumens
Resolution	Under testing; designed with same accuracy specifications as VeinViewer Vision, ~97% accurate vein width against ultrasound
Power Source	Battery or AC (outlet)
Battery Life	Up to 2 hours continuous run time (4 hours with additional battery)
Dimensions	11.5 x 1.5 x 3.8" (29.21 x 3.81 x 9.65cm)
Weight	1.6lbs (0.73kg)
Durability	Drop tested per IEC/UL 60601-1

TESTIMONIALS

"We have been using our VeinViewer for about 18 months and love it. We use it predominantly for reticular veins and identifying the veins which feed the large clusters of spider telangiectasias. It is also very helpful for the lateral subdermal venous plexus. Most importantly, the patients absolutely love the technology and I feel that it legitimizes us as vein specialists."

—William Schuh, MD
Vein Care Institute
Denver, CO

"Dr. Miyake says his VeinViewer boosted referrals by 20%, cut costs by 30%, and reduced the need for anesthesia."

—Economist.com
Augmented Reality: Reality, Only Better
December 2007

See the VeinViewer difference for yourself!
Call today for a free demonstration.

CHRISTIE+ www.VeinViewer.com

CHRISTIE MEDICAL HOLDINGS, INC. | 1256 UNION AVE. 3rd FL. MEMPHIS, TN 38104 | 877-SEE-VEIN (877-733-8346) TOLL FREE