

P R O D U C T   C A T A L O G



**busse**  
Hospital Disposables

1 . 8 0 0 . 6 4 5 . 6 5 2 6   W W W . B U S S E I N C . C O M

# THE BUSSE COMMITMENT

Busse Hospital Disposables has a proven track record for over 40 years of providing the health care industry with superb products and services:

- High-quality products precisely manufactured
- Strict adherence to industry standards and guidelines
- Flexibility with both standard and custom solutions to meet your needs
- Responsiveness with prompt price quotes, delivery of samples, and processing and delivering of orders
- Cost-effective solutions with fair and competitive pricing
- ISO certified and CE compliant, where necessary
- And, Busse's manufacturing facility is located in the U.S.A. to assure our consumers the quality and performance they've come to expect from us.

**Visit our website at [www.busseinc.com](http://www.busseinc.com) for the most up-to-date product information.**

If you don't see a product you require, please contact your local Busse representative, or call our main office at 1.800.645.6526.

Unless stated as reusable, all Busse products are for single use only.



ISO 13485-2003  
Certified  
FM 77400

**busse**  
Hospital Disposables

# CONTENTS

Pain Management Trays	
Single Shot Epidural Tray.....	4
Continuous Epidural Tray.....	4
Basic Nerve Block Tray.....	4
Pain Management Components for Custom Trays.....	5
Titan Spinal Needles and LOR Syringes.....	6
Tumescent Delivery and Vein Closure Trays.....	7
Diagnostic Procedure Trays	
Lumbar Puncture Trays.....	8, 9
Manometer.....	9
Basic Soft Tissue Biopsy Trays.....	10
Amniocentesis Trays.....	10
Myelogram and Arthrogram Trays.....	11
Bone Marrow Biopsy/Aspiration Trays.....	12
Bone Marrow Biopsy/Aspiration Needle.....	13
Thoracentesis and Paracentesis Trays.....	14
Central Line Dressing Change Kits and Trays.....	15, 16
I.V. Start Kits.....	17
Incision and Drainage, Wound Closure, General Purpose Trays.....	18
Debridement and Nosebleed Trays.....	18
Minor Laceration and General Purpose Instrument Trays, Minor Procedure Pack.....	19
E.R. Laceration and Suturing Trays, Irrigation Splash Guard.....	20
Skin Staple Remover Kits, Nurses' Scissors, Forceps, Hemostats and Sterile Floor-Grade, Single-Use Instrument.....	21
Suture Removal Kits.....	22
Tracheostomy Care Sets.....	23, 24
Suction Products.....	25, 26
Tubing.....	26, 28
Tubing Connectors and Catheter Plug.....	27
Curettes and Laminaria Cervical Dilator.....	28
Mucus Specimen Traps.....	29
Urine Collection and Calculi Strainer.....	30
Ear/Ulcer Syringes.....	31
Labor and Delivery.....	31, 32
Sterile Fields, Towels and Drapes.....	33, 34, 35
Epidural/Pain Management and Ophthalmic Drapes.....	36
Ice Bags.....	37
O.R. Head and Shoe Coverings.....	38
Isolation/Procedural Masks and Staff Protection Kit.....	39
Staff Protection Gowns.....	40
Isolation Gowns.....	41
Protective Lab Coats and Coveralls.....	42
Post Mortem Bags and Kits.....	43, 44, 45
Disaster Bag, Shroud Sheets and Kits.....	46
Bowls and Containers.....	47
Patient Care.....	48
Accessories.....	49

# STERILE PAIN MANAGEMENT TRAYS

CUSTOMIZED SOLUTIONS AVAILABLE



At Busse, we recognize that pain management is the fastest growing sector of the medical field today. Hospitals, clinics and surgi-centers are increasingly choosing pain management techniques over more costly and more invasive procedures. Increasing awareness of the benefits of regional anesthesia vs. general anesthesia has created a demand for a variety of pain management products. To meet this demand, Busse offers a line of Pain Management Trays.

## BENEFITS OF CHOOSING BUSSE:

Premium products that fit your requirements

Flexibility – standard and custom trays

Personalized products and service

Engineered sharps protection components available

Quick turnaround time allows you to keep inventory low

Prompt price quote and samples



## SINGLE SHOT EPIDURAL TRAY

### PREP COMPONENTS

- One 1 oz. PVP solution
- Three sponge sticks
- Four 3" x 3" gauze dressings
- One absorbent towel
- One integral prep well
- One clear, fenestrated drape with adhesive strip

### CATALOG NO. 670

Packed 10/carton

### PROCEDURE COMPONENTS

- Two multi-compartment trays
- One 18g x 3-1/2" Tuohy epidural needle with wings
- One 5cc LOR glass syringe with metal Luer slip
- One 3cc plastic Luer lock syringe
- One 20cc plastic Luer lock syringe
- One 18g x 1-1/2" needle
- One 25g x 1-1/2" needle
- One filter needle
- One foam needle stick block
- One 5mL ampule 1% lidocaine hydrochloride, USP
- One 10mL 0.9% sodium chloride
- One integral mixing well
- One adhesive bandage
- One CSR wrap

## CONTINUOUS EPIDURAL TRAY

### PREP COMPONENTS

- One 1 oz. PVP solution
- Three sponge sticks
- Four 3" x 3" gauze dressings
- One absorbent towel
- One integral prep well
- One clear, fenestrated drape with adhesive strip

### CATALOG NO. 671

Packed 10/carton

### PROCEDURE COMPONENTS

- Two multi-compartment trays
- One 20g closed-end epidural catheter
- One thread assist guide
- One Tuohy Borst adapter
- One catheter pad
- One epidural filter
- One 18g x 3-1/2" Tuohy epidural needle with wings
- One 5cc LOR glass syringe with metal Luer slip
- One 3cc plastic Luer lock syringe
- One 20cc plastic Luer lock syringe
- One 18g x 1-1/2" needle
- One 25g x 1-1/2" needle
- One filter needle
- One foam needle stick block
- One 1mL ampule 0.1% epinephrine
- One 5mL ampule 1% lidocaine hydrochloride, USP
- One 5mL ampule 1.5% lidocaine with epinephrine
- One 10mL 0.9% sodium chloride
- One integral mixing well
- One adhesive bandage
- One CSR wrap

## BASIC NERVE BLOCK TRAY

### PREP COMPONENTS

- One 1 oz. PVP solution
- Three sponge sticks
- Four 3" x 3" gauze dressings
- Two absorbent towels
- One skin marker
- One 6" flexible ruler
- One integral prep well
- One clear, fenestrated drape with adhesive strip

### CATALOG NO. 672

Packed 10/carton

### PROCEDURE COMPONENTS

- Two multi-compartment trays
- One 1cc plastic Luer lock syringe
- One 3cc plastic Luer lock syringe
- One 5cc plastic Luer lock syringe
- One 20cc plastic Luer lock control syringe
- One 18g x 1-1/2" needle
- One 25g x 1-1/2" needle
- One filter needle
- One foam needle stick block
- One 5mL ampule 1% lidocaine hydrochloride, USP
- One integral mixing well
- One adhesive bandage
- One CSR wrap

# CUSTOMIZED SOLUTIONS FOR PAIN MANAGEMENT

The Busse Pain Management Tray program enables you to build the exact tray for your clinical needs.

## STEP 1 – DETERMINE YOUR NEEDS

What procedure(s) will you be performing?

- Single Shot Epidural
- Continuous Epidural
- Combined Spinal/Epidural (CSE)
- Nerve Block

What components do you need?

- Compare your component requirements to our standard tray listings on the previous page.

What if you need components not included in our standard trays?

- Simple... build your own custom tray.
- Using one of our standard trays as your foundation, you may add or remove any components to customize a tray to your facility's specifications.

## STEP 2 – SELECT A STANDARD TRAY

Choose one of our three standard trays (on the previous page) packed with the preferred prep and procedure components. The layouts of these trays are specially designed to facilitate quick identification and easy access. Our standard trays serve as the building blocks to unlimited custom trays.

## FOR ASSISTANCE

If you do not see the components you require, contact your local Busse representative, or one of our Pain Management Product Specialists at

**1-800-645-6526, Ext 217**

## STEP 3 – BUILD YOUR OWN CUSTOM STERILE PAIN MANAGEMENT TRAY

Add optional components to our standard trays or create a total custom solution. **The following list provides a sample of components to build your custom tray. If you do not see the component that you require on this list, please call us.**

### EPIDURAL NEEDLES

Tuohy, Weiss (fixed wing, detachable wing)

Available in various lengths, *such as*

18g x 3-1/2"  
20g x 3-1/2"  
22g x 3-1/2"

### SPINAL NEEDLES

(Assorted styles including Whitacre, Quincke)

Available in various lengths, *such as*

18g x 3-1/2"    22g x 5"  
20g x 3-1/2"    22g x 7"  
22g x 3-1/2"

### NEEDLES

(Standard, Safety), *such as*

Needle	18g x 1-1/2"
Needle	20g x 1-1/2"
Needle	21g x 1-1/2"
Needle	22g x 1-1/2"
Needle	25g x 5/8"
Needle	25g x 1-1/2"
Needle	27g x 1-1/2"
Filter Needle	19g x 1-1/2"

### LOSS OF RESISTANCE (LOR) SYRINGES

(Luer lock, Luer slip), *such as*

Glass syringe	5cc
(metal/glass tip)	
Glass syringe	10cc
(metal/glass tip)	
Plastic syringe	10cc

### CONTROL SYRINGES

*such as*

Control syringe	10cc
Control syringe	20cc

### SYRINGES

(Luer lock, Luer slip), *such as*

Syringe	1cc
Syringe	3cc
Syringe	5cc
Syringe	10cc
Syringe	20cc
Syringe	60cc

### EPIDURAL CATHETERS

(Open or closed-end, with or without stylet)  
Available in various lengths, *such as*

19g Epidural catheter  
20g Epidural catheter

### MEDICATIONS

*such as*

Lidocaine 1%  
Sodium chloride 0.9%  
Epinephrine  
Epinephrine and lidocaine

### DRAPES & TOWELS

Available in a variety of materials, fenestrations, adhesive applications and sizes, *such as*

Blue, fenestrated drape  
Blue, fenestrated drape with adhesive strip  
Clear, fenestrated drape  
Clear, fenestrated drape with adhesive strip  
Blue, absorbent towel

### EXTENSION SETS

Assorted lengths available, *such as*

Extension set 6"  
Extension set 12"  
Extension set 20"

### SAFETY COMPONENTS


*such as*

Saf-T Cube  
Safety Scalpel  
Safety Needles  
Filter Needles





# STERILE TITAN SPINAL NEEDLES & POSI-SPACE LOR SYRINGES



## LOSS OF RESISTANCE (LOR) PLASTIC SYRINGES

Our LOR syringes have been designed with smooth plunger movement to improve user sensitivity and provide accurate detection of the epidural space.

Sterile plastic 10mL LOR syringe with Luer lock tip.

### CATALOG NO. 595

Packed 5 dispenser boxes of 10 each, 50/carton


Sterile plastic 10mL LOR syringe with Luer slip tip.

### CATALOG NO. 596

Packed 5 dispenser boxes of 10 each, 50/carton



Titan Spinal Needles are designed with the quality and precision necessary for superlative performance. To ensure patient safety and comfort, the needles maintain stability with minimal deflection during insertion. Featuring translucent hubs for unobstructed viewing during injection. Trust Titan Spinal Needles to fully support your procedures and pre-ops with the quality they demand.



## QUINCKE-STYLE SPINAL NEEDLES

### STANDARD LENGTH

#### CATALOG NO. 560

18g x 3-1/2"  
Packed 4 dispenser boxes of 25 each, 100/carton

#### CATALOG NO. 561

20g x 3-1/2"  
Packed 4 dispenser boxes of 25 each, 100/carton

#### CATALOG NO. 562

22g x 3-1/2"  
Packed 4 dispenser boxes of 25 each, 100/carton

#### CATALOG NO. 563

25g x 3-1/2"  
Packed 4 dispenser boxes of 25 each, 100/carton

### LONG LENGTH

#### CATALOG NO. 565

22g x 5"  
Packed 10 dispenser boxes of 5 each, 50/carton


#### CATALOG NO. 566

22g x 7"  
Packed 10 dispenser boxes of 5 each, 50/carton

### SHORT LENGTH

#### CATALOG NO. 564

22g x 1-1/2"  
Packed 4 dispenser boxes of 25 each, 100/carton



# STERILE TUMESCENT DELIVERY & VEIN CLOSURE

CUSTOMIZED  
SOLUTIONS AVAILABLE



Allows multiple doses of fluid to be delivered.  
No manual reloading of syringe is necessary.

Polycarbonate syringe assists in the controlled  
delivery of tumescent anesthesia with less  
hand fatigue.

Dual-check valve allows repetitive injection of  
anesthesia solution.

## TUMESCENT DELIVERY SYSTEM

One 10cc syringe with dual-check valve  
One 10' I.V. tubing

**CATALOG NO. 620**

Packed 10/carton



## ALL-IN-ONE CONVENIENCE



Busse's Vein Closure Trays are  
designed to save time and money.

Our trays include the key  
components most commonly used  
in a vein closure procedure.

Busse's trays are manufactured  
under strict guidelines for  
uncompromising quality.

Custom trays available to meet your  
requirements upon request.  
(Minimum quantities apply)

## DELUXE VEIN CLOSURE TRAY

One tray  
One Tumescent Delivery System  
(10cc tumescent syringe with  
dual-check valve and 10' tubing)  
One 44" x 88" table cover  
One 77" x 110" split drape  
One 55" x 44" half drape  
One 16 oz. (500cc) plastic bowl  
One 5mL ampule 1% lidocaine  
hydrochloride, USP  
Twenty 4"x 4" gauze sponges  
One 3cc Luer lock syringe  
One 10cc Luer lock syringe  
One 20cc Luer lock syringe  
One 25g x 5/8" needle

One 25g x 1-1/2" needle  
One 22g x 3-1/2" spinal needle for  
tumescent injection  
One #11 scalpel blade  
One Saf-T Cube  
Two ChloraPrep® 10.5mL applicators  
Three ultrasound gel, 20g packets  
One probe cover  
One card of tape strips, 3 per card  
One surgical skin marker  
Four blue O.R. towels  
One surgical gown, XL, with towel  
One surgical gown, LG, with towel

**CATALOG NO. 622**

Packed 5/carton



## BASIC VEIN CLOSURE TRAY

One tray  
One 44" x 88" table cover  
One 77" x 110" split drape  
One 55" x 44" half drape  
One 16 oz. (500cc) plastic bowl  
One 5mL ampule 1% lidocaine  
hydrochloride, USP  
Twenty 4"x 4" gauze sponges  
One 3cc Luer lock syringe  
One 10cc Luer lock syringe  
One 25g x 1-1/2" needle  
One 22g x 3-1/2" spinal needle for  
tumescent injection

One #11 scalpel blade  
One Saf-T Cube  
Two ChloraPrep® 10.5mL applicators  
Three ultrasound gel, 20g packets  
One probe cover  
One card of tape strips, 3 per card  
One surgical skin marker  
Four blue O.R. towels  
One surgical gown, XL, with towel  
One surgical gown, LG, with towel

**CATALOG NO. 621**

Packed 5/carton







**STERILE  
DIAGNOSTIC  
PROCEDURE TRAYS**  
LUMBAR PUNCTURE TRAYS

CUSTOMIZED  
SOLUTIONS AVAILABLE

**ADULT TRAYS**



One multi-compartment tray with prep well  
One 2mL ampule 1% lidocaine hydrochloride, USP  
One 3cc syringe with 25g x 5/8" needle  
Three sponge sticks  
Three gauze pads  
One towel  
One fenestrated drape with adhesive strip  
One bandage  
One 20g x 3-1/2" spinal needle with stylet  
One 3-way stopcock  
One 5" extension tube  
One 2-piece 550mm manometer  
One 22g x 1-1/2" infiltration needle  
Four 10mL specimen vials, prenumbered, with caps  
Four patient I.D. labels  
One CSR wrap  
One set of 10 Universal Medication labels

**CATALOG NO. 648**

Packed 10/carton

Same components as Catalog No. 648, but with  
18g x 3-1/2" spinal needle.

**CATALOG NO. 647**

Packed 10/carton

Same components as Catalog No. 648, but with  
22g x 3-1/2" spinal needle

**CATALOG NO. 649**

Packed 10/carton



**SAFETY-DELUXE ADULT LUMBAR  
PUNCTURE TRAYS**

One multi-compartment tray with prep well  
One 20g x 3-1/2" spinal needle  
One 2-piece 550mm manometer  
One 3-way stopcock  
One 5" extension set  
Four 10mL specimen vials, prenumbered, with caps  
and patient I.D. labels  
One 5mL ampule 1% lidocaine hydrochloride, USP  
One 25g x 1" BD SafetyGlide™ needle  
One 22g x 1-1/2" BD SafetyGlide™ needle  
One 5cc Luer lock syringe  
One ear loop mask  
One Filter Straw®  
One fenestrated drape with adhesive strip  
One towel  
One fine tip marker  
Three gauze pads  
One needle foam stop  
Two Point-Lok® needle retention safety devices  
One PVP triple swabstick  
Three sponge sticks  
One time-out label  
One bandage  
One CSR wrap  
One set of 10 Universal Medication labels


**CATALOG NO. 974**

Packed 10/carton

Same components as Catalog No. 974, but with a  
22g x 3-1/2" spinal needle

**CATALOG NO. 975**

Packed 10/carton



**STERILE MANOMETER**

One 1-piece 550mm sterile manometer with  
3-way stopcock

**CATALOG NO. 644**

Individually packaged

Packed 5 dispenser boxes of  
10 each, 50/carton



BD SafetyGlide™ is a trademark of Becton, Dickinson and Company.  
Filter Straw® is a registered trademark of B.Braun.  
Point-Lok® is a registered trademark of Smiths Medical

# STERILE BASIC SOFT TISSUE BIOPSY & AMNIOCENTESIS TRAYS

CUSTOMIZED  
SOLUTIONS AVAILABLE



## BASIC SOFT TISSUE BIOPSY TRAY

This tray is complete with all of the common components required for skin preparation and handling of the aspirated sample. Trays are provided without a needle to allow the practitioner the freedom to use the needle of their choice. This tray may be used during CT-guided biopsies.

- One multi-compartment tray with prep well
- One 3cc syringe with 25g x 5/8" needle
- One 5cc Luer lock syringe
- One 20cc Luer lock syringe
- One 22g x 1-1/2" needle
- One 20g x 1-1/2" needle
- One 5mL ampule 1% lidocaine hydrochloride, USP
- One scalpel
- Two 10mL specimen vials, prenumbered, with caps
- Six frosted microscope slides
- One 6" ruler
- Two sponge sticks
- Four gauze pads
- Two towels
- One fenestrated drape with adhesive strip
- One bandage
- Two patient I.D. labels
- One CSR wrap
- One set of 10 Universal Medication labels

### CATALOG NO. 651

Packed 10/carton



Busse Safety-Deluxe Biopsy and Amniocentesis Trays are designed to help you comply with today's safety requirements.

## SAFETY-DELUXE BASIC SOFT TISSUE BIOPSY TRAY

- One multi-compartment tray with prep well
- One 6" ruler
- Two 10mL specimen vials, prenumbered, with caps and patient I.D. labels
- Two 5mL ampules 1% lidocaine hydrochloride, USP
- One 25g x 1" BD SafetyGlide™ needle
- One 22g x 1-1/2" BD SafetyGlide™ needle
- Two 10cc Luer lock syringes
- One ear loop mask
- One Filter Straw®
- One fenestrated drape with adhesive strip
- Two towels
- One fine tip marker
- Ten gauze pads
- One needle foam stop
- One Point-Lok® needle retention safety device
- One NeedleVISE® needle retention safety device
- One 3mL ChloraPrep® Applicator Hi Lite Orange™
- One time-out label
- One bandage
- One alcohol prep pad
- Ten frosted microscope slides
- One safety scalpel
- One CSR wrap
- One set of 10 Universal Medication labels

### CATALOG NO. 977

Packed 10/carton



## AMNIOCENTESIS TRAY

- One multi-compartment tray with prep well
- One 20g x 3-1/2" spinal needle with stylet
- One 6cc Luer lock syringe
- One 20cc Luer lock syringe
- One 25g x 5/8" needle
- One 22g x 1-1/2" needle
- One 5mL ampule 1% lidocaine hydrochloride, USP
- Three 10mL amber specimen vials, prenumbered, with caps
- Three sponge sticks
- One 6" ruler
- Three gauze pads
- One towel
- One fenestrated drape with adhesive strip
- One bandage
- Three patient I.D. labels
- One CSR wrap
- One set of 10 Universal Medication labels

### CATALOG NO. 650

Packed 10/carton

## SAFETY-DELUXE AMNIOCENTESIS TRAY

- One multi-compartment tray with prep well
- One 20g x 3-1/2" spinal needle
- One 6" ruler
- Three 10mL amber specimen vials, prenumbered, with caps and patient I.D. labels
- One 5mL ampule 1% lidocaine hydrochloride, USP
- One 25g x 1" BD SafetyGlide™ needle
- One 22g x 1-1/2" BD SafetyGlide™ needle
- One 5cc Luer lock syringe
- One 20cc Luer lock syringe
- One ear loop mask
- One Filter Straw®
- One fenestrated drape with adhesive strip
- One towel
- One fine tip marker
- Three gauze pads
- One needle foam stop
- Two Point-Lok® needle retention safety devices
- One 3mL ChloraPrep® Applicator Hi Lite Orange™
- Three sponge sticks
- One time-out label
- One bandage
- One CSR wrap
- One set of 10 Universal Medication labels

### CATALOG NO. 976

Packed 10/carton



BD SafetyGlide™ is a trademark of Becton, Dickinson and Company. Filter Straw® is a registered trademark of B. Braun. Point-Lok® is a registered trademark of Smiths Medical. NeedleVISE® is a registered trademark of Atrion Medical. ChloraPrep® and ChloraPrep® Applicator Hi Lite Orange™ are trademarks of Carefusion.

# STERILE DIAGNOSIS MYELOGRAM & ARTHROGRAM TRAYS



## BUSSE'S MYELOGRAM TRAYS

Busse's Myelogram Trays contain the components necessary to safely sample cerebral spinal fluid and the injection of water-based contrast media (not for use with oil-based contrast media) with engineered sharps protection components.



## BUSSE'S ARTHROGRAM TRAYS

Busse's Arthrogram Trays provide all the components for successful double and single contrast arthrography procedures.

**Busse Safety-Deluxe Myelogram and Arthrogram Trays are designed to help comply with today's safety requirements.**

## MYELOGRAM TRAY

One multi-compartment tray with prep well  
 One 22g x 3-1/2" spinal needle with stylet  
 One 5" extension tube  
 One 20" extension tube  
 One 5cc Luer lock syringe  
 One 20cc Luer lock syringe  
 One 18g x 1-1/2" needle  
 One 20g x 1-1/2" needle  
 One 25g x 5/8" needle  
 One 5mL ampule 1% lidocaine hydrochloride, USP  
 Three 10mL specimen vials, prenumbered, with caps  
 Three sponge sticks  
 Three gauze pads  
 One towel  
 One fenestrated drape with adhesive strip  
 One bandage  
 Three patient I.D. labels  
 One CSR wrap  
 One set of 10 Universal Medication labels

### CATALOG NO. 656

Packed 10/carton

## SAFETY MYELOGRAM TRAY

Contains all of the components of Catalog No. 656, but with sharps protection components.

### CATALOG NO. 956

Packed 10/carton

## SAFETY ARTHROGRAM TRAY

One multi-compartment tray with prep well  
 One 20" extension tube  
 Two 10cc Luer lock syringe  
 One 60cc Luer lock syringe  
 One 18g x 1-1/2" needle  
 One 20g x 1-1/2" needle  
 One 21g x 1-1/2" needle  
 One 25g x 5/8" needle  
 One Saf-T Cube  
 One 5mL ampule 1% lidocaine hydrochloride, USP  
 Three sponge sticks  
 Three gauze pads  
 One absorbent towel  
 One fenestrated drape with adhesive strip  
 One bandage  
 One CSR wrap  
 One set of 10 Universal Medication labels

### CATALOG NO. 957

Packed 10/carton

## SAFETY-DELUXE MYELOGRAM TRAY

One multi-compartment tray with prep well  
 One 22g x 3-1/2" spinal needle  
 One 5" extension set  
 One 20" extension set  
 Three 10mL specimen vials, prenumbered, with caps and patient I.D. labels  
 One 5mL ampule 1% lidocaine hydrochloride, USP  
 One 25g x 1" BD SafetyGlide™ needle  
 One 22g x 1-1/2" BD SafetyGlide™ needle  
 One 18g x 1-1/2" BD SafetyGlide™ needle  
 One 5cc Luer lock syringe  
 One 20cc Luer lock syringe  
 One ear loop mask  
 One Filter Straw®  
 One fenestrated drape with adhesive strip  
 One towel  
 One fine tip marker  
 Three gauze pads  
 One needle foam stop  
 Two Point-Lok® needle retention safety devices  
 One PVP triple swabstick  
 Three sponge sticks  
 One time-out label  
 One bandage  
 One CSR wrap  
 One set of 10 Universal Medication labels

### CATALOG NO. 978

Packed 10/carton



## SAFETY-DELUXE ARTHROGRAM TRAY

One multi-compartment tray with prep well  
 One 20" extension set  
 One 5mL ampule 1% lidocaine hydrochloride, USP  
 One 25g x 1" BD SafetyGlide™ needle  
 One 22g x 1-1/2" BD SafetyGlide™ needle  
 One 18g x 1-1/2" BD SafetyGlide™ needle  
 Two 10cc Luer lock syringes  
 One 60cc Luer lock syringe  
 One ear loop mask  
 One Filter Straw®  
 One fenestrated drape with adhesive strip  
 One towel  
 One fine tip marker  
 Three gauze pads  
 One needle foam stop  
 One 3mL ChloroPrep® Applicator Hi Lite Orange™  
 Three sponge sticks  
 One time-out label  
 One bandage  
 One CSR wrap  
 One set of 10 Universal Medication labels

### CATALOG NO. 979

Packed 10/carton



1-800-645-6526

11

BD SafetyGlide™ is a trademark of Becton, Dickinson and Company. Filter Straw® is a registered trademark of B. Braun. Point-Lok® is a registered trademark of Smiths Medical. NeedleWISE® is a registered trademark of Atrion Medical. ChloroPrep® and ChloroPrep® Applicator Hi Lite Orange™ are trademarks of Carefusion.

# STERILE BONE MARROW BIOPSY/ASPIRATION TRAYS

CUSTOMIZED SOLUTIONS AVAILABLE



The bone marrow trays offer both Illinois-style Sternal-Iliac needles and our exclusive, enhanced J-style bone marrow biopsy/aspiration needles.



- Two multi-compartment trays
- One PVP triple swabstick
- One 5mL ampule 1% lidocaine hydrochloride, USP
- One scalpel
- One 5cc Luer lock syringe with 20g x 1-1/2" needle
- One 20cc Luer lock syringe
- One 25g x 5/8" needle
- One 21g x 1-1/2" needle
- One towel
- One fenestrated drape with adhesive strip
- One bandage
- One 11g x 4" J-style bone marrow biopsy needle with cap, extractor and guide
- One 15g I-style aspiration needle, adjustable from 3/8" to 1-7/8"
- One Luer lock adapter
- Two specimen tubes containing edetate (EDTA) potassium
- Ten frosted microscope slides
- One amber specimen vial with cap
- One stirring rod
- Five gauze pads
- One patient I.D. label
- One CSR wrap
- One set of 10 Universal Medication labels

## CATALOG NO. 652

Packed 10/carton

## CATALOG NO. 653

Packed 10/carton

- Two multi-compartment trays
- One PVP triple swabstick
- One 5mL ampule 1% lidocaine hydrochloride, USP
- One scalpel
- One 5cc Luer lock syringe with 20g x 1-1/2" needle
- One 20cc Luer lock syringe
- One 25g x 5/8" needle
- One 21g x 1-1/2" needle
- One towel
- One fenestrated drape with adhesive strip
- One bandage

- One 11g x 4" J-style bone marrow biopsy needle, with cap, extractor and guide
- One Luer lock adapter
- Two specimen tubes containing edetate (EDTA) potassium
- Ten frosted microscope slides
- One amber specimen vial with cap
- One stirring rod
- Five gauze pads
- One patient I.D. label
- One CSR wrap
- One set of 10 Universal Medication labels

## CATALOG NO. 654

Packed 10/carton

- Two multi-compartment trays
- One PVP triple swabstick
- One 5mL ampule 1% lidocaine hydrochloride, USP
- One scalpel
- One 5cc Luer lock syringe with 20g x 1-1/2" needle
- One 20cc Luer lock syringe
- One 25g x 5/8" needle
- One 21g x 1-1/2" needle
- One towel
- One fenestrated drape with adhesive strip

- One bandage
- One 11g x 4" J-style bone marrow biopsy needle, with cap, extractor and guide
- One 15g, I-style aspiration needle, adjustable from 3/8" to 1-7/8"
- One Luer lock adapter
- Eight gauze pads
- One CSR wrap
- One set of 10 Universal Medication labels

## CATALOG NO. 655

Packed 10/carton

- Two multi-compartment trays
- One PVP triple swabstick
- One 5mL ampule 1% lidocaine hydrochloride, USP
- One scalpel
- One 5cc Luer lock syringe with 20g x 1-1/2" needle
- One 20cc Luer lock syringe
- One 25g x 5/8" needle
- One 21g x 1-1/2" needle
- One towel
- One fenestrated drape with adhesive strip
- One bandage

- One 15g I-style aspiration needle, adjustable from 3/8" to 1-7/8"
- One Luer lock adapter
- Two specimen tubes containing edetate (EDTA) potassium
- Ten frosted microscope slides
- One amber specimen vial with cap
- One stirring rod
- Five gauze pads
- One patient I.D. label
- One CSR wrap
- One set of 10 Universal Medication labels



## SAFETY DELUXE BONE MARROW BIOPSY/ASPIRATION TRAYS

CUSTOMIZED SOLUTIONS AVAILABLE

Two multi-compartment trays  
 One 11g x 4" J-Style bone marrow biopsy needle, with cap, extractor and guide  
 One 15g I-Style aspiration needle, adjustable from 3/8" to 1-7/8"  
 One 5cc Luer lock syringe  
 One 20cc Luer lock syringe  
 One 25g x 1" BD SafetyGlide™ needle  
 One 22g x 1-1/2" BD SafetyGlide™ needle  
 One safety scalpel  
 One Filter Straw®  
 One time out label  
 Two Saf-T Cubes  
 One skin marker  
 One ear loop mask  
 One 5mL ampule 1% lidocaine hydrochloride, USP  
 Ten frosted microscope slides  
 Five gauze pads  
 One towel  
 One fenestrated drape with adhesive strip  
 Two specimen tubes containing edetate (EDTA) potassium  
 One stirring rod  
 One amber specimen vial with cap  
 One patient I.D. label  
 One Luer lock adapter  
 One PVP triple swabstick  
 One bandage  
 One CSR Wrap  
 One set of 10 Universal Medication labels

### CATALOG NO. 985

Packed 10/carton

Same components as Catalog No. 985, but with only one Saf-T Cube and without I-style aspiration needle.

### CATALOG NO. 986

Packed 10/carton


Two multi-compartment trays  
 One 11g x 4" J-Style bone marrow biopsy needle, with cap, extractor and guide  
 One 15g I-Style aspiration needle, adjustable from 3/8" to 1-7/8"  
 One 5cc Luer lock syringe  
 One 20cc Luer lock syringe  
 One 25g x 1" BD SafetyGlide™ needle  
 One 22g x 1-1/2" BD SafetyGlide™ needle  
 One safety scalpel  
 One Filter Straw®  
 One time out label  
 Two Saf-T-Cubes  
 One skin marker  
 One ear loop mask  
 One 5mL ampule 1% lidocaine hydrochloride, USP  
 Eight gauze pads  
 One towel  
 One fenestrated drape with adhesive strip  
 One Luer lock adapter  
 One PVP triple swabstick  
 One bandage  
 One CSR Wrap  
 One set of 10 Universal Medication labels

### CATALOG NO. 987

Packed 10/carton


BD SafetyGlide™ is a trademark of Becton, Dickinson and Company  
 Filter Straw® is a registered trademark of B.Braun.

## J-STYLE BONE MARROW BIOPSY/ASPIRATION NEEDLE



Ultra-sharp and manufactured to the highest tolerances, the Busse bone marrow biopsy needle is designed to provide outstanding and consistent performance.

With its substantial, easily gripped handle, our 11g x 4" J-style needle provides the physician with greater torque than similar instruments. Combined with the ultra-sharp trocar tip on the stylet, Busse's exclusive design offers much easier penetration of bone tissue. Once the stylet is withdrawn, the cannula tip presents four sharp, beveled planes with pointed serrations at their junctures, and produces tightly sheared marrow samples as the needle is advanced.




J-style bone marrow/ aspiration needle, 11g x 4"

### CATALOG NO. 630

Packed 10/carton



## I-STYLE STERNAL-ILIAIC BONE MARROW BIOPSY/ASPIRATION NEEDLE



Manufactured in the popular Illinois Sternal-Iliac design, Busse's 15g I-style needle allows variable depth adjustment from 3/8" to 1-7/8". Luer lock allows for easy connection to syringe.

I-style Sternal-Iliac aspiration needle, 15g x 3/8" to 1-7/8" adjustable

### CATALOG NO. 631

Packed 10/carton





# STERILE THORACENTESIS/ PARACENTESIS TRAYS

CUSTOMIZED  
SOLUTIONS AVAILABLE

Designed for image-guided procedures, these three trays are designed specifically for safe and cost-effective drainage and collection of fluid in pleural space and peritoneal cavity.

## THORACENTESIS TRAY

For use with drainage bag or vacuum bottle

- One tray
- One 14g catheter over 17g x 6" needle
- One drainage bag
- One tubing set for drainage bag
- One tubing set with needle for vacuum bottle
- One 3-way stopcock
- One 25g x 5/8" safety needle
- One 21g x 1-1/2" safety needle
- One 5cc Luer lock syringe
- One 60cc Luer lock syringe
- One 5mL ampule 1% lidocaine hydrochloride, USP
- Three specimen vials
- One fenestrated drape with adhesive strip
- One bandage
- One underpad/towel
- One 3mL ChloraPrep® One-Step Applicator
- Six gauze dressings
- One safety scalpel
- One Saf-T-cube
- One CSR wrap
- One set of 10 Universal Medication labels

### CATALOG NO. 960

Packed 10/carton



ChloraPrep® is a registered trademark of Carefusion.

## THORACENTESIS TRAY



### CATALOG NO. 961

Packed 10/carton



Same components as Catalog No. 960, without drainage bag and tubing set for drainage bag. **For use with vacuum bottle.**

## PARACENTESIS TRAY



### CATALOG NO. 962

Packed 10/carton



Same components as Catalog No. 960, but with 14g catheter over 17g x 3-1/4" (5.7 cm) needle. **For use with drainage bag or vacuum bottle.**

# STERILE CENTRAL LINE DRESSING CHANGE KITS & TRAYS WITH CHLORAPREP®

FOR TPN/CVP AND OTHER CENTRAL LINE MAINTENANCE

CUSTOMIZED SOLUTIONS AVAILABLE



## CENTRAL LINE KIT WITH TEGADERM™ DRESSING AND CHLORAPREP® ONE-STEP APPLICATOR

- One soft pack
- One polylined, sterile field drape
- One mask
- One pair gloves
- One ChloraPrep® One-Step applicator, 3mL
- One skin protectant prep pad
- One isopropyl alcohol swabstick
- One 4" x 4-3/4" Tegaderm™ transparent dressing
- One 2" x 2" non-adhering dressing
- One 2" x 2" non-adhering split dressing
- One 2" x 2" gauze dressing
- One 4" x 4" gauze dressing
- One roll medical tape
- One dressing change label

### CATALOG NO. 836

Packed 50 kits/carton



## CENTRAL LINE TRAY WITH TEGADERM™ DRESSING AND CHLORAPREP® ONE-STEP APPLICATOR

- One tray
- One mask
- One pair gloves
- One ChloraPrep® One-Step applicator, 3mL
- One isopropyl alcohol swabstick
- One 4" x 4-3/4" Tegaderm™ transparent dressing
- One 2" x 2" non-adhering dressing
- One 2" x 2" non-adhering split dressing
- One 2" x 2" gauze dressing
- One 4" x 4" gauze dressing
- One roll medical tape
- One dressing change label
- One CSR wrap

### CATALOG NO. 835

Packed 30 trays/carton



## CENTRAL LINE KIT WITH TEGADERM™ DRESSING WITH TRANSPARENT BORDERS AND CHLORAPREP® ONE-STEP APPLICATOR



- One soft pack
- One polylined, sterile field drape
- One mask
- One pair gloves
- One ChloraPrep® One-Step applicator, 3mL
- One skin protectant swabstick
- One isopropyl alcohol swabstick
- One 4" x 4-3/4" Tegaderm™ transparent dressing with absorbent borders
- One 2" x 2" non-adhering dressing
- One 2" x 2" non-adhering split dressing
- One 2" x 2" gauze dressing
- One 4" x 4" gauze dressing
- One roll medical tape
- One dressing change label

### CATALOG NO. 837

Packed 50 kits/carton



ChlorPrep® is a registered trademark of Carefusion.  
Tegaderm™ is a trademark of 3M, St. Paul, Minnesota.

# STERILE CENTRAL LINE DRESSING CHANGE KITS & TRAYS

FOR TPN/CVP AND OTHER CENTRAL LINE MAINTENANCE

CUSTOMIZED SOLUTIONS AVAILABLE



## CENTRAL VENOUS CATHETER DRESSING CHANGE TRAY WITH TEGADERM™ DRESSING AND ISOPROPYL ALCOHOL SWABSTICKS

- One tray
- One mask
- One pair gloves
- One povidone iodine ointment
- Three iodophor PVP swabsticks
- Three isopropyl alcohol swabsticks
- One 2" x 2" non-adhering dressing
- One 2" x 2" gauze dressing
- One 4" x 4" gauze dressing
- One 2" x 2" non-adhering split dressing
- One 4" x 4-3/4" Tegaderm™ transparent dressing
- One roll medical tape
- One dressing change label
- One CSR wrap

### CATALOG NO. 829

Packed 30 trays/carton



Tegaderm™ is a trademark of 3M, St. Paul, Minnesota.

## CENTRAL VENOUS CATHETER DRESSING CHANGE KIT WITH TEGADERM™ DRESSING AND SKIN PROTECTANT SWABSTICKS



### CATALOG NO. 832

Packed 50 kits/carton



- One soft pack
- One polylined, sterile field drape
- One pair gloves
- One iodophor PVP ointment
- Three isopropyl alcohol swabsticks
- Three iodophor PVP swabsticks
- One skin protectant swabstick
- One 2" x 2" non-woven dressing
- One 2" x 2" non-woven split dressing
- One 3" x 3" gauze dressing
- One 4" x 4-3/4" Tegaderm™ transparent dressing
- One roll medical tape
- One dressing change label

## CENTRAL LINE DRESSING CHANGE KIT WITH TEGADERM™ DRESSING AND SKIN PROTECTANT SWABSTICKS



### CATALOG NO. 834

Packed 50 kits/carton



- One soft pack
- One polylined, sterile field drape
- One mask
- One pair gloves
- Three isopropyl alcohol swabsticks
- Three iodophor PVP swabsticks
- One skin protectant swabstick
- One 2" x 2" non-woven dressing
- One 2" x 2" gauze dressing
- One 3" x 3" gauze dressing
- One 4" x 5" Tegaderm™ transparent dressing
- One roll medical tape

## PRIMARY DRESSING CHANGE TRAY WITH SALINE SOLUTION AND STRETCH GAUZE DRESSING



### CATALOG NO. 1443

Packed 20 trays/carton



- One three-compartment tray
- One pair gloves
- One polylined, sterile field drape
- One abdominal pad
- Eight 4" x 4" gauze sponges
- One 100mL container sterile saline solution
- One polybag with twist tie
- One roll medical tape
- One medical ruler
- One sharp/blunt scissors
- One stretch gauze dressing
- One dressing change label

## DRESSING CHANGE TRAY WITH INSTRUMENTS AND ABDOMINAL PAD



### CATALOG NO. 759

Packed 20 trays/carton



- One tray
- One polylined, sterile field drape
- One abdominal pad
- One safety pin
- One iodophor PVP swabstick (triples)
- Four 4" x 4" gauze sponges
- One sharp/blunt scissors
- One polybag with twist tie
- One roll medical tape
- One transfer forceps
- One Posi-Grip™ plastic forceps



# STERILE I.V. START KITS

CUSTOMIZED SOLUTIONS AVAILABLE



An economical I.V. start kit that provides the infection control needed.

## I.V. START KIT WITH CATHETER SECUREMENT DEVICE AND CHLORASCRUB® SWAB

- One soft pack
- One Chlorascrub® swab
- One Posi-Guard catheter securement device
- One transparent window dressing
- One latex-free tourniquet
- Two 2" x 2" gauze dressings
- One roll of medical tape
- One I.V. site label

### CATALOG NO. 822

Packed 50 kits/carton



Chlorascrub® swab is a registered trademark of Professional Disposables International Inc.  
Chloraprep® is a registered trademark of Carefusion.  
Tegaderm™ is a trademark of 3M, St. Paul, Minnesota.

## STANDARD I.V. START KIT WITH CHLORAPREP® SEPP® APPLICATOR



- One soft pack
- One Chloraprep® Sepp® applicator, .67mL
- One latex-free tourniquet
- One 2-3/8" x 2-3/4" Tegaderm™ transparent dressing
- Two 2" x 2" gauze dressings
- One roll medical tape
- One I.V. site label



### CATALOG NO. 818

Packed 50 kits/carton

## I.V. START KIT WITH TEGADERM™ DRESSING



- One soft pack
- One 2-3/8" x 2-3/4" Tegaderm™ dressing
- One latex-free tourniquet
- One alcohol prep pad
- One PVP prep pad
- Two 2" x 2" gauze sponges
- One roll medical tape

### CATALOG NO. 820

Packed 50 kits/carton



# STERILE INCISION & DRAINAGE TRAY, WOUND CLOSURE TRAYS, SHARP DEBRIDEMENT TRAY & NOSEBLEED TRAY

CUSTOMIZED SOLUTIONS AVAILABLE



## INCISION AND DRAINAGE TRAY

- One tray
- One 5-1/2" floor-grade, straight Kelly hemostat
- One scalpel, #11 blade
- One Posi-Grip™ plastic forceps
- One polylined, sterile field
- One polylined, fenestrated drape
- One polybag with twist tie
- Two large iodophor PVP prep pads
- Four cotton-filled sponges

### CATALOG NO. 758

Packed 20 trays/carton

## SAFETY INCISION AND DRAINAGE TRAY

Contains all the components of Catalog No. 758, but with sharps protection components.

### CATALOG NO. 558

Packed 20 trays/carton



## DELUXE WOUND CLOSURE TRAY WITH INSTRUMENTS



### CATALOG NO. 756

Packed 20 trays/carton



- One tray
- One 5" floor-grade, Baumgartner (serrated jaw) needle holder
- One 5-1/2" floor-grade, curved Kelly hemostat
- One 5-1/4" sharp/blunt scissors
- One 5" floor-grade tissue forceps
- One polybag with twist tie
- One polylined, sterile field
- One polylined, fenestrated drape
- One large iodophor PVP prep pad
- One large alcohol prep pad
- Four gauze sponges

## DELUXE FACIAL WOUND CLOSURE TRAY WITH INSTRUMENTS



### CATALOG NO. 757

Packed 20 trays/carton



- One tray
- One 5" floor-grade, Webster (smooth jaw) needle holder
- One 5" floor-grade, curved Mosquito hemostat
- One 4-1/2" sharp/sharp scissors
- One 5" metal tissue forceps
- One polybag with twist tie
- One polylined, sterile field
- One polylined, fenestrated drape
- One large iodophor PVP prep swabstick
- One large alcohol prep pad
- Four gauze sponges

## SHARP DEBRIDEMENT TRAY




### CATALOG NO. 744

Packed 50 trays/carton



- One tray
- One pair Nitrile gloves
- One alcohol prep pad
- One PVP prep pad
- One gauze sponge
- One Iris scissors
- One scalpel #15 blade
- One metal insert forceps

## NOSEBLEED TRAY



### CATALOG NO. 640

Packed 10 trays/carton



- One tray
- Two medicine cups
- Three cotton balls
- Five gauze dressings
- One 10 Fr Fraiser instrument
- One Hartman nasal speculum
- One Gruenwald nasal forceps
- One Iris scissors
- One instrument holder
- One CSR wrap



# STERILE MINOR LACERATION TRAYS AND GENERAL PURPOSE INSTRUMENT TRAYS

CUSTOMIZED SOLUTIONS AVAILABLE



## MINOR LACERATION TRAY WITH INSTRUMENTS

- One tray
- One 5" floor-grade, Baumgartner (serrated jaw) needle holder
- One 4-1/2" sharp/sharp scissors
- One 5" metal tissue forceps
- Two absorbent towels
- Five gauze sponges
- One polylined, sterile field
- One polylined, fenestrated drape

### CATALOG NO. 751

Packed 20 trays/carton



## MINOR LACERATION TRAY WITHOUT INSTRUMENTS

Contains all the components of Catalog No. 751, excluding the metal instruments.

### CATALOG NO. 750

Packed 20 trays/carton



## GENERAL PURPOSE INSTRUMENT TRAY



- One tray
- One 5" floor-grade, straight Mosquito hemostat
- One 4-1/2" sharp/sharp scissors
- One 5" metal thumb forceps
- One safety pin
- Two cotton-filled sponges
- Two cotton-tipped applicators
- One iodophor PVP prep pad
- One large alcohol prep pad

### CATALOG NO. 753

Packed 20 trays/carton



## GENERAL PURPOSE INSTRUMENT TRAY

Contains all the components of Catalog No. 753, but uses a floor-grade, curved Mosquito hemostat in place of the straight Mosquito hemostat.

### CATALOG NO. 752

Packed 20 trays/carton



## GENERAL PURPOSE INSTRUMENT TRAY



- One tray
- One 5-1/2" floor-grade, curved Kelly hemostat
- One 4-1/2" sharp/blunt scissors
- One 5" metal thumb forceps
- One safety pin
- Two cotton-filled sponges
- Two cotton-tipped applicators
- One iodophor PVP prep pad
- One large alcohol prep pad

### CATALOG NO. 754

Packed 20 trays/carton



## GENERAL PURPOSE INSTRUMENT TRAY

Contains all the components of Catalog No. 754, but uses a floor-grade, straight Kelly hemostat in place of the curved Kelly hemostat.

### CATALOG NO. 755

Packed 20 trays/carton



## STER-L-UNIT™ MINOR PROCEDURE PACK WITHOUT INSTRUMENTS – CSR WRAPPED



- One tray
- One fenestrated drape, 18" x 26"
- One polylined, sterile field drape
- Five gauze sponges
- Two absorbent towels



### CATALOG NO. 3220

Packed 8 dispensers of 10 each, 80/carton

1-800-645-6526

# STERILE E.R. LACERATION TRAYS AND IRRIGATION SPLASH GUARD

CUSTOMIZED  
SOLUTIONS AVAILABLE



## E.R. LACERATION TRAY WITH FLOOR-GRADE POLISHED INSTRUMENTS

One tray  
One 5" floor-grade, Webster (smooth) needle holder  
One 4-1/2" floor-grade, straight Iris scissors  
One 5" curved, floor-grade, Mosquito hemostat  
One 1" x 2" floor-grade, Adson tissue forceps  
Two 2 oz. medicine cups  
One 10cc Luer lock syringe  
One 25g x 5/8" needle  
One 27g x 1-1/2" needle  
One 18g x 1-1/2" needle  
Five 2" x 2" gauze sponges  
Ten 4" x 4" gauze sponges  
One polylined, fenestrated drape  
Two absorbent towels  
One CSR wrap  
One set of 10 Universal Medication labels

### CATALOG NO. 749

Packed 20 trays/carton



Contains all the components of Catalog No. 749, but with sharps protection components.

### CATALOG NO. 549

Packed 20 trays/carton

## SAF-SHIELD™ IRRIGATION SPLASH GUARD



Fits all syringes – Luer lock or slip-fit – and the protective dome acts as a shield against splash back.

### CATALOG NO. 766 Sterile

Packed 50/carton

### CATALOG NO. 767 Non-Sterile

Bulk packed 200/carton

## E.R. LACERATION TRAY WITH FLOOR-GRADE POLISHED INSTRUMENTS AND CLOTH TOWELS



### CATALOG NO. 746

Packed 20 trays/carton

### CATALOG NO. 546

Packed 20 trays/carton



Contains all the components of Catalog No. 746, but with sharps protection components.

One tray  
One 5" floor-grade Webster (smooth) needle holder  
One 4-1/2" floor-grade, straight Iris scissors  
One 5" curved, floor-grade, Mosquito hemostat  
One 1" x 2" floor-grade, Adson tissue forceps  
Two 2 oz medicine cups  
One 10cc Luer lock syringe  
One 25g x 5/8" needle  
One 18g x 1-1/2" needle  
One 27g x 1-1/2" needle  
Five 2" x 2" gauze sponges  
Ten 4" x 4" gauze sponges  
One blue cotton O.R. fenestrated drape  
Two blue cotton O.R. towels  
One CSR wrap  
One set of 10 Universal Medication labels

## SUTURING TRAY WITH FLOOR-GRADE SATIN FINISHED INSTRUMENTS



### CATALOG NO. 747

Packed 20 trays/carton



One tray  
One 5" floor-grade, Halsey needle holder, fine point, satin finish  
One 4-1/2" floor-grade, straight Iris scissors  
One 1" x 2" floor-grade, Adson tissue forceps, fine point, satin finish  
Two 2 oz. medicine cups  
Five 2" x 2" gauze sponges  
Eight 4" x 4" gauze sponges  
One polylined, fenestrated drape  
One absorbent towel  
One CSR wrap  
One set of 10 Universal Medication labels

## E.R. LACERATION KIT WITH SAF-SHIELD™ IRRIGATION SPLASH GUARD



### CATALOG NO. 748

Packed 20 trays/carton



One tray  
One 5" floor-grade, Webster (smooth) needle holder  
One 4-1/2" floor-grade, straight Iris scissors  
One 5" curved, floor-grade, Mosquito hemostat  
One 1" x 2" floor-grade, Adson tissue forceps  
Two 2 oz. medicine cups  
One Saf-Shield™ Irrigation Splash Guard  
One 10cc Luer lock syringe  
One 18g x 1-1/2" needle  
One 27g x 1-1/2" needle  
Five 2" x 2" gauze sponges  
Ten 4" x 4" gauze sponges  
One polylined, fenestrated drape  
Two absorbent towels  
One CSR wrap  
One set of 10 Universal Medication labels

# STERILE SKIN STAPLE REMOVAL KITS AND INDIVIDUALLY PACKED INSTRUMENTS



## STERILE SKIN STAPLE REMOVER KIT

Made of stainless steel, our skin staple remover allows fast, easy removal with minimal discomfort to the patient. Each sterile, single-use instrument is packed along with a 3" x 3" gauze sponge, in a sealed plastic tray with an easy-to-open, peel-back lid.

One tray  
One skin staple remover  
One 3" x 3" gauze sponge

**CATALOG NO. 716**  
Packed 48 kits/carton



## STERILE SKIN STAPLE REMOVER KIT – DISPENSER PACKED

Same as Catalog No. 716, but packed as four dispenser packs of 12 each, 48 total.

**CATALOG NO. 770**



## STERILE SKIN STAPLE REMOVER WITH PLASTIC HANDLE



Made of stainless steel, this skin staple remover provides fast, easy removal of staples with minimal discomfort to patient. The scissors-style handles allow for a more controlled grip. Each sterile, single-use instrument is in a sealed plastic tray with an easy-to-open, peel-back lid.

One tray  
One stainless steel skin staple remover with plastic handles

**CATALOG NO. 772**  
Individually packed 30/carton



## STERILE GENERAL PURPOSE NURSES' SCISSORS



Sterile, 5-1/4" sharp/blunt scissors, convenient for use by entire staff. Cuts gauze or tape easily. Each scissor is individually packed in a sterile tray.

**CATALOG NO. 7246**  
Individually packed 100 trays/carton



## 4" POSI-GRIP™ PLASTIC FORCEPS

**CATALOG NO. 7190**  
Sterile, packed 100/carton



## 4" POSI-GRIP™ PLASTIC FORCEPS

Same as Cat. No. 7190, but non-sterile.

**CATALOG NO. 7195**  
Bulk packed 200/carton



## STERILE FLOOR-GRADE, SINGLE-USE INSTRUMENTS



### SHARP/BLUNT NURSES' SCISSORS 5-1/2"

**CATALOG NO. 773**  
Individually packed 50/carton

### CURVED KELLY HEMOSTAT 5-1/2"

**CATALOG NO. 774**  
Individually packed 50/carton

### STRAIGHT KELLY HEMOSTAT 5-1/2"

**CATALOG NO. 775**  
Individually packed 50/carton





# STERILE SUTURE REMOVAL KITS

CUSTOMIZED SOLUTIONS AVAILABLE



## STERILE SUTURE REMOVAL KITS FEATURING FLOOR-GRADE INSTRUMENTS

- One tray
- One floor-grade, Iris suture scissors
- One floor-grade, Adson serrated forceps
- One 3" x 3" gauze sponge
- One large alcohol prep pad
- One iodophor PVP prep pad

### CATALOG NO. 723

Packed 50 kits/carton



## STERILE SUTURE REMOVAL KIT

Same as Catalog No. 723, with floor-grade, Littauer (hook-tipped) suture scissors in place of the Iris.

### CATALOG NO. 726

Packed 50 kits/carton



## STERILE SUTURE REMOVAL KIT

Same as Catalog No. 723, with metal forceps in place of the Adson serrated forceps.

### CATALOG NO. 729

Packed 50 kits/carton



All Suture Removal Kits are packed in a rigid tray.

## STERILE "CLASSIC" SUTURE REMOVAL KIT



- One tray
- One plastic handle Littauer (hook-tipped) scissors
- One Posi-Grip™ plastic forceps
- One 3" x 3" gauze dressing

### CATALOG NO. 737

Packed 48 kits/carton



Catalog No. 737 features the Littauer hook-tipped blade.

## STERILE SUTURE REMOVAL KIT



- One tray
- One metal Littauer (hook-tipped) scissors
- One Posi-Grip™ plastic forceps
- One 3" x 3" gauze sponge
- One large alcohol prep pad
- One iodophor PVP prep pad

### CATALOG NO. 724

Packed 50 kits/carton



## STERILE SUTURE REMOVAL KIT WITH METAL FORCEPS



- One tray
- One metal Littauer (hook-tipped) scissors
- One 4" metal forceps
- One 3" x 3" gauze sponge

### CATALOG NO. 718

Packed 50 kits/carton



## STERILE SUTURE REMOVAL KIT WITH PLASTIC FORCEPS



- One tray
- One metal Littauer (hook-tipped) scissors
- One Posi-Grip™ plastic forceps
- One 3" x 3" gauze sponge

### CATALOG NO. 717

Packed 50 kits/carton



## STERILE SUTURE REMOVAL KIT WITH STITCH CUTTER



- One tray
- One stitch cutter
- One Posi-Grip™ plastic forceps
- One 3" x 3" gauze sponge

### CATALOG NO. 732

Packed 50 kits/carton



# STERILE TRACHEOSTOMY CARE SETS

## THREE-COMPARTMENT TRAYS

CUSTOMIZED SOLUTIONS AVAILABLE



### TRACHEOSTOMY CARE SET – STANDARD

- One three-compartment tray
- One pair gloves
- One pre-cut gauze tracheostomy dressing
- One polylined towel/drape
- Four 4" x 4" gauze sponges
- Four pipe cleaners
- Two cotton-tipped applicators
- One flexible nylon bristle brush
- One 36" twill tape

#### CATALOG NO. 715

Packed 20 sets/carton



### TRACHEOSTOMY CARE SET WITH HYDROGEN PEROXIDE AND SALINE SOLUTION



#### CATALOG NO. 800

Packed 20 sets/carton



- One three-compartment tray
- One pair gloves
- One pre-cut gauze tracheostomy dressing
- One polylined towel/drape
- Four 4" x 4" gauze sponges
- Two pipe cleaners
- Three cotton-tipped applicators
- One flexible nylon bristle brush
- One 30" twill tape
- One foam neck pad
- One 100mL container sterile saline solution
- Two 1 oz. packets hydrogen peroxide

### TRACHEOSTOMY CARE SET WITH HYDROGEN PEROXIDE



#### CATALOG NO. 714

Packed 20 sets/carton



- One three-compartment tray
- One pair gloves
- One pre-cut gauze tracheostomy dressing
- One polylined towel/drape
- Four 4" x 4" gauze sponges
- Two pipe cleaners
- Two cotton-tipped applicators
- One flexible nylon bristle brush
- One 30" twill tape
- Two 1 oz. packets hydrogen peroxide

### TRACHEOSTOMY CARE SET WITH FORCEPS



#### CATALOG NO. 703

Packed 20 sets/carton



- One three-compartment tray
- One pair gloves
- One pre-cut gauze tracheostomy dressing
- Two polylined towel/drapes
- One Posi-Grip™ plastic forceps
- Three 4" x 4" gauze sponges
- Three pipe cleaners
- Three cotton-tipped applicators
- One flexible nylon bristle brush
- One 30" twill tape

### TRACHEOSTOMY CARE SET WITH FORCEPS WITHOUT GLOVES

Same as Catalog No. 703, without gloves.



#### CATALOG NO. 799

Packed 20 sets/carton

1-800-645-6526



# STERILE TRACHEOSTOMY CARE SETS

## THREE-COMPARTMENT TRAYS

CUSTOMIZED SOLUTIONS AVAILABLE



### TRACHEOSTOMY CARE SET WITH EXTRA DRAPE

- One three-compartment tray
- One pair gloves
- One pre-cut gauze tracheostomy dressing
- Two polylined towel/drapes
- Four 4" x 4" gauze sponges
- Six pipe cleaners
- Two cotton-tipped applicators
- One flexible nylon bristle brush
- One 30" twill tape

#### CATALOG NO. 711

Packed 24 sets/carton



### TRACHEOSTOMY CARE SET WITH EXTRA DRAPE WITHOUT GLOVES

Same as Catalog No. 711, except without gloves.

#### CATALOG NO. 710

Packed 24 sets/carton



## TWO-COMPARTMENT TRAYS

### TRACHEOSTOMY CARE SET WITH REMOVABLE BASIN



- One two-compartment tray with removable basin
- One pair gloves
- One pre-cut gauze tracheostomy dressing
- One polylined towel/drape
- Four 4" x 4" gauze sponges
- Four pipe cleaners
- Two cotton-tipped applicators
- One flexible nylon bristle brush
- One 36" twill tape

#### CATALOG NO. 713

Packed 20 sets/carton



### TRACHEOSTOMY CARE SET WITH REMOVABLE BASIN AND SUCTION CATHETER



- One two-compartment tray with removable basin
- One 14 Fr catheter
- One pair gloves
- Three pre-cut gauze tracheostomy dressings
- Two polylined towel/drapes
- Four 4" x 4" gauze sponges
- Four pipe cleaners
- Two cotton-tipped applicators
- One flexible nylon bristle brush
- One 30" twill tape

#### CATALOG NO. 706

Packed 20 sets/carton



### STANDARD TRACHEOSTOMY DRESSING



Busse's individually wrapped, pre-cut gauze dressings complement every tracheostomy care set.

Highly absorbent quality gauze provides patient comfort and protection.

Sterile single dressing offers the ultimate in economy, yet assures maximum patient protection for each tracheostomy dressing replacement.

Contents:


- One 4" x 4" tracheostomy dressing

#### CATALOG NO. 707

Individually packed sterile 300/carton (6 trays of 50 each)



# STERILE SUCTION PRODUCTS



## STERILE POOLE SUCTION TIPS

The Busse Poole is made of tough, shatter-resistant plastic that features a retractable screen. It is a versatile, three-position unit which features a one-hand thumb "click and slide" operation for the outer screen to extend or retract the pinpoint suction tip. With the screen in place, the Poole will help eliminate tissue grab and deep abdominal trauma. With the screen retracted, the Poole presents a fine suction tip. With the screen removed, the Poole can reach into confined areas.

### CATALOG NO. 290

Instrument, 50/carton

### CATALOG NO. 291

With 6' x 1/4" ID connective tubing, 20/carton



## STERILE YANKAUER SUCTION TIPS



Open Tip



Bulb Tip

The Busse Yankauer is a rigid, curved, one-piece design, made from a blend of clear, blue-tinted, shatterproof plastic. Edges are smooth and clean, with no flash or parting lines. Clear tips permit visual monitoring of suction fluids. A ribbed, universal connector fits securely on all types of tubing. Busse Yankauers are available with or without a suction control vent, in two tip styles.

## STERILE YANKAUER SUCTION TIPS

### CATALOG NO. 296

Open Tip, w/o vent  
Packed 50/carton

### CATALOG NO. 297

Open Tip, w/ vent  
Packed 50/carton

### CATALOG NO. 298

Bulb Tip, w/o vent  
Packed 50/carton

### CATALOG NO. 299

Bulb Tip, w/ vent  
Packed 50/carton



## STERILE YANKAUER SUCTION TIPS WITH TUBING

### CATALOG NO. 302

Bulb Tip, w/o vent,  
6' x 1/4" ID connective tubing  
Packed 20/carton

### CATALOG NO. 303

Bulb Tip, w/ vent,  
6' x 1/4" ID connective tubing  
Packed 20/carton

### CATALOG NO. 305

Bulb Tip, w/o vent,  
10' x 1/4" ID connective tubing  
Packed 20/carton



## STERILE SUCTION PRODUCTS



### STERILE FRAZIER SUCTION TIPS

Comes complete with a stylette that allows the malleable shaft to be bent to the desired position without occlusion, and a control vent permits continuous or intermittent suctioning. Available in 8, 10, 12 French sizes.

#### CATALOG NO. 162

8 Fr, 5 dispenser boxes of 10 each, 50/carton

#### CATALOG NO. 163

10 Fr, 5 dispenser boxes of 10 each, 50/carton

#### CATALOG NO. 164

12 Fr, 5 dispenser boxes of 10 each, 50/carton



### STERILE SIGMOIDOSCOPIC SUCTION TIPS



Clear plastic shaft allows visual monitoring of suctioning, while a sliding control vent permits continuous or intermittent suction. Universal type connector facilitates easy connection of all tubing. Packed sterile with or without connective tubing.

#### CATALOG NO. 165

18 Fr 30 cm, 50/carton

#### CATALOG NO. 167

18 Fr 30 cm, with 6' x 3/16" ID connective tubing, 20/carton



### STERILE SUCTION TUBING FOR COSMETIC SURGERY



High volume suction tubing. Thick enough to eliminate collapse, but still flexible and easy to use. Able to withstand high negative pressure and is transparent for ease of visualization. The tubing features a large bore to maximize volume flow and drainage.

#### CATALOG NO. 1539

8' x 3/8" ID connective tubing, 10/carton

#### CATALOG NO. 1540

10' x 3/8" ID connective tubing, 10/carton



### STERILE CONNECTIVE TUBING



#### CATALOG NO. 154

6' x 1/4" ID connective tubing, 50/carton

#### CATALOG NO. 155

10' x 1/4" ID connective tubing, 50/carton

#### CATALOG NO. 156

12' x 1/4" ID connective tubing, 20/carton


#### CATALOG NO. 157

6' x 3/16" ID connective tubing, 50/carton



Soft, low durometer tubing is memory-free, and made of clear plastic to permit monitoring of suction fluids. A ribbed outer wall gives added strength and slip resistance. Soft and durable connectors are permanently affixed to tubing.

# TUBING CONNECTORS AND CATHETER PLUG




Tough, shatterproof polypropylene connectors adapt to almost any hospital use, and permit visual monitoring of fluids. Sterile connectors are packed in individual plasti-packs; clean connectors are packed in convenient dispenser cartons. Select from nine different types and sizes.



## 5-IN-1 TUBING CONNECTOR

Fits five sizes of tubing. Connect same size tubing or any combination of 5 to 11mm. (3/16", 1/4", 5/16", 3/8", 7/16")




**CATALOG NO. 501** Sterile  
Packed 25/box; 8 boxes/carton

**CATALOG NO. 511** Clean  
Packed 50/box; 10 boxes/carton

## CATHETER PLUG WITH DRAINAGE TUBE COVER

Tapered plug with easy-grip handle offers fast, easy plugging of any size catheter



**CATALOG NO. 510** Sterile  
Packed 200/carton

## SIMS TUBING CONNECTOR/ IRRIGATION NOZZLE

Nozzle fits any size tubing from 3 to 6mm (1/8" to 1/4") in seven size increments. Connect 1/4" to 1/4" or use as a reducer from 1/4" to 1/8" tubing.




**CATALOG NO. 502** Sterile  
Packed 50/box; 8 boxes/carton

**CATALOG NO. 512** Clean  
Packed 50/box; 10 boxes/carton

## "T" TUBING CONNECTOR

Joins three sections of tubing. Five steps on each leg provide connections of any size tubing from 5 to 11mm (3/16" to 7/16").



**CATALOG NO. 503** Sterile  
Packed 15/box; 8 boxes/carton

**CATALOG NO. 513** Clean  
Packed 15/box; 10 boxes/carton

## STRAIGHT TUBING CONNECTORS


Joins two sections of same size tubing.



### 6MM (1/4")

**CATALOG NO. 504** Sterile


**CATALOG NO. 514** Clean



### 8MM (5/16")

**CATALOG NO. 505** Sterile

**CATALOG NO. 515** Clean



### 10MM (3/8")


**CATALOG NO. 506** Sterile

**CATALOG NO. 516** Clean

Sterile: Packed 50/box; 8 boxes/carton  
Clean: Packed 50/box; 10 boxes/carton

## "Y" TUBING CONNECTOR


Joins three sections of tubing. Steps on each leg permit joining of multiple sizes for each connector.



Fits tubing from 5 to 8mm (3/16" to 5/16").

**CATALOG NO. 507** Sterile  
Packed 50/box; 8 boxes/carton

**CATALOG NO. 517** Clean  
Packed 50/box; 10 boxes/carton




Fits tubing from 8 to 11mm (5/16" to 7/16").

**CATALOG NO. 508** Sterile  
Packed 25/box; 8 boxes/carton

**CATALOG NO. 518** Clean  
Packed 50/box; 10 boxes/carton

## 6-IN-1 "Y" TUBING CONNECTOR

Tapered legs fit six sizes of tubing, from 6mm to 15mm (1/4" to 9/16").




**CATALOG NO. 509** Sterile  
Packed 15/box; 8 boxes/carton

**CATALOG NO. 519** Clean  
Packed 15/box; 10 boxes/carton




# STERILE UTERINE ASPIRATION SYSTEM

## STERILE VACUUM CURETTES




Disposable vacuum curettes are made of clear, rigid plastic, with molded markings for size and tip orientation, and a visible, frosted sounding line to gauge proper insertion depth. Our curettes are strong, firm, and finely finished; wall thickness is uniform throughout; and there are no sharp edges anywhere. Select from seven sizes for minimal uterine dilation of each patient.

Busse vacuum curettes are available with open end tips, in both curved or straight designs, that fit all popular styles of permanent or reusable handles.



**CURVED  
OPEN TIP**



**STRAIGHT  
OPEN TIP**

CATALOG		CATALOG	
NO.	SIZE	NO.	SIZE
<b>279</b>	7mm	<b>133</b>	7mm
<b>280</b>	8mm	<b>134</b>	8mm
<b>281</b>	9mm	<b>135</b>	9mm
<b>282</b>	10mm	<b>136</b>	10mm
<b>283</b>	11mm	<b>137</b>	11mm
<b>284</b>	12mm	<b>138</b>	12mm
<b>285</b>	14mm	<b>139</b>	14mm

Each single-use, disposable curette is packed in an individual, sterile, see-through peel pack, properly sized for easy aseptic removal.

Packed 25/carton.



## SUCTION HANDLES AND TUBING



Swivel handle with 6" x 1/2" ID tubing, with male adaptor to fit Berkeley pump.

### CATALOG NO. 150

Packed 25/carton

Swivel handle with 6" x 1/2" ID tubing without adaptor

### CATALOG NO. 151

Packed 25/carton

Swivel handle with 6" x 3/8" ID tubing, with male adaptor to fit Berkeley pump.

### CATALOG NO. 153

Packed 25/carton



## STERILE SUCTION TUBING FOR COSMETIC SURGERY



High-volume suction tubing. Thick enough to eliminate collapse, but still flexible and is easy to use. Able to withstand high negative pressure and transparent for ease of visualization. The tubing features a large bore to maximize volume flow and drainage.

8' x 3/8" ID tubing  
**CATALOG NO. 1539**  
10/carton

10' x 3/8" ID tubing  
**CATALOG NO. 1540**  
10/carton



## LAMINARIA CERVICAL DILATOR



The gradual and uniform swelling of this natural substance, a byproduct of seaweed, is far less traumatic than rapid dilation, making Laminaria the preferred method of dilation or softening of the cervix. Normal dilation is two-to-three times the diameter of the selected size, with most swelling occurring within the first six hours after insertion. Ethylene Oxide processed. All approximately 60mm long.

### CATALOG NO.

- 206**
- 207**
- 208**
- 209**

### SIZE

- 3±mm/small
- 4±mm/medium
- 5±mm/large
- 6±mm/extra large

Packed 4 dispenser boxes of 12 each,  
48/carton





# STERILE MUCUS SPECIMEN TRAPS

## STERILE MUCUS SPECIMEN TRAPS

A simple yet efficient way to collect specimens using any suction pump or device. Two see-through plastic container sizes, 40cc and 80cc, with graduations clearly marked. Disposable adaptor prevents contamination from non-sterile suction tube. Screw-on cap protects lips and helps maintain sterile pouring surface. Collapse-proof, 5" long heavy wall, non-latex tubing.

Each packed sterile, in an easy-to-peel envelope that assures aseptic handling. Write-on bottle.



### STERILE MUCUS SPECIMEN TRAP WITH TRANSPORT CAP

Eliminate the possibility of accidental tube separation during pneumatic tube transport with this sterile mucus specimen trap. This trap is packed with an extra cap – designed to seal the trap for transport. After collection, simply remove the standard tubed cap and replace it with the transport cap.

#### CATALOG NO. 405

40cc  
Packed 24/carton



#### CATALOG NO. 406

40cc  
Packed 24/carton



### STERILE MUCUS SPECIMEN TRAP, WHEN AN 80CC CAPACITY IS DESIRED

#### CATALOG NO. 407

80cc  
Packed 4 boxes of 25 each, 100/carton



# URINE COLLECTION AND CALCULI STRAINER



## CALCULI STRAINER

Mesh screen is optimal size for filtering the smallest stones while allowing unobstructed flow-through

Built-in tab for easier handling and sanitary conditions

Durable, white plastic construction

Disposable

Rinsable plastic allows one patient to use strainer for several days

Fits all 24-hour specimen collection containers

For use in hospitals and homes

### CATALOG NO. 460

Packed 100/carton



## STERILE URINE COLLECTION



Smooth, tapered catheter tip and eyes for low-trauma insertions

Sheath protects catheter from bending during shipping

Unique vented lid promotes urine flow without the need to loosen lid

Closed system reduces the risk of contamination and false readings

Printed graduations on tube makes it easier for technician to see levels of urine

15cc centrifuge-ready tubes

Kits are packed in a pouch organizer, procedure-ready



## FEMALE CATHETER KIT

- One 6" rigid catheter, 8 Fr
- One pair vinyl exam gloves
- Three iodophor PVP swabsticks
- One 3 gm. lubricant packet
- One towel
- One patient I.D. label
- One blue, plastic pouch

### CATALOG NO. 410

Packed 4 inner cartons of 25 each, 100/carton

## ADULT CATHETER KIT (FEMALE OR MALE USE)

- One 9" soft catheter, 8 Fr
- One pair vinyl exam gloves
- Three iodophor PVP swabsticks
- One 3 gm. lubricant packet
- One towel
- One patient I.D. label
- One blue, plastic pouch

### CATALOG NO. 411

Packed 50/carton

## PEDIATRIC CATHETER KIT

- One 9" soft catheter, 8 Fr
- One pair vinyl exam gloves
- Three iodophor PVP swabsticks
- One 3 gm. lubricant packet
- One towel
- One patient I.D. label
- One blue, plastic pouch
- Instruction sheet

### CATALOG NO. 412

Packed 50/carton

## INFANT CATHETER KIT

- One 9" soft catheter, 5 Fr
- One pair vinyl exam gloves
- Three iodophor PVP swabsticks
- One 3 gm. lubricant packet
- One towel
- One patient I.D. label
- One blue, plastic pouch
- Instruction sheet

### CATALOG NO. 413

Packed 50/carton



## LABOR AND DELIVERY AND GENERAL IRRIGATION



### REUSABLE ALL-PURPOSE EAR/ULCER SYRINGE

Sturdy, flexible, latex-free Kraton™ allows increased suction, better control and a better grip when wet. Packed sterile in a peel-back blister pack.

**CATALOG NO. 148**

**2 OZ.** Sterile

Packed 24/carton

**CATALOG NO. 149**

**3 OZ.** Sterile

Packed 12/carton

**CATALOG NO. 140**

**2 OZ.** Non-sterile

Bulk packed, 200/carton

**CATALOG NO. 147**

**3 OZ.** Non-sterile

Bulk packed, 100/carton



### ALL-PURPOSE DISPOSABLE EAR/ULCER SYRINGE Long tapered tip and increased suction power



Can be used for any aspiration or irrigation procedure. Pronounced ribs provide a secure grip for more accurate volume control. Flat bottom allows syringe to stand on its own. Sterile blister pack assures aseptic technique.

**CATALOG NO. 141 1 OZ.** Green, Sterile

Packed 50/carton

**CATALOG NO. 142 2 OZ.** Blue, Sterile

Packed 50/carton

**CATALOG NO. 143 3 OZ.** Green, Sterile

Packed 50/carton

**CATALOG NO. 144 2 OZ.** Blue, Non-Sterile

Bulk packed, 200/carton

**CATALOG NO. 145 3 OZ.** Green, Non-Sterile

Bulk packed, 100/carton

**CATALOG NO. 146 1 OZ.** Green, Non-Sterile

Bulk packed, 200/carton



# LABOR AND DELIVERY



## POSI-GRIP™ UMBILICAL CORD CLAMP

The single-use clamp is designed to seal the umbilical cord, preventing loss of blood as the cord dries beneath it. It can be easily placed with one hand, even if gloved and slippery. The interlocking teeth hold the cord securely and keep it stationary. Once the protected closure is in place and closed, it cannot open accidentally. Included in each box of 64 is one cord clamp cutter, 2 cord clamp cutters per case of 128. Works with any style cord cutter. Available sterile and non-sterile.

### CATALOG NO. 384 Sterile

Individual packets, 128/carton  
(2 boxes of 64 each)

### CATALOG NO. 385 Non-Sterile

Bulk packed, 250/carton



## TAPE MEASURES

Newborn 24" – shows inches and millimeters/centimeters

### CATALOG NO. 791

Packed 1,000/carton

Pediatric 36" – shows inches and millimeters/centimeters

### CATALOG NO. 795

Packed 1,000/carton



## OBSTETRICAL KITS



### CATALOG NO. 388 Pouch

Packed 10 kits/carton

### CATALOG NO. 389 Box

Packed 10 kits/carton



A complete disposable obstetrical kit available in cardboard box or soft pouch packaging. Contains all necessary equipment for emergency delivery situations.

- One pair sterile exam gloves
- One disposable plastic apron
- One 17" x 24" underpad
- One 23" x 36" receiving blanket
- One 48" x 36" drape sheet
- Four absorbent towels
- Two OB towelettes
- Four sterile 4" x 3" cover sponges
- One sterile #10 scalpel
- One sterile OB pad
- One sterile 2 oz. bulb syringe
- Two alcohol prep pads
- Two sterile Posi-Grip™ umbilical cord clamps
- One clear polybag with tie
- One white polybag with tie



## STERILE CIRCUMCISION TRAY



### CATALOG NO. 386 Pouch

Packed 10 trays/carton

- One 16" x 18" blue/white drape with 5/8" fenestration
- Two 13" x 19" white absorbent towels
- Six 4" x 4" 12 ply gauze dressings
- One 1" x 8" petrolatum gauze dressing
- Two 5" straight, Mosquito hemostats with large finger loops
- One 5" curved, Mosquito hemostat with large finger loops
- One 5-1/2" probe with eye
- One 4-1/2" S/B straight operating dressing scissors
- One Adson tissue forceps
- One scalpel with #10 blade
- One safety pin
- One foam instrument holder/10 slit die cut / 5-instrument
- One 2 oz. medicine cup
- One CSR wrap



## AMNIOPERF™ — AMNIOTIC MEMBRANE PERFORATOR



The Busse Sterile AmnioPerf™ is designed to rupture the amniotic membrane without any injury to the fetus or mother. This device allows medical personnel to carefully snag the membrane, without scratching or causing any damage.

### CATALOG NO. 387

Packed 100/carton





# STERILE FIELDS, TOWELS AND DRAPES



Multi-layered, all-purpose surgical drapes bring convenience, performance and added protection. These drapes are absorbent yet impervious to strikethrough. Polylining, affixed by a heat-sealing process, prevents passage of liquids, air and bacteria, preserves integrity of sterile field.

## STERILE FIELD

18" x 26"  
Sterile  
Tissue/poly/tissue  
White/non-glare blue tissue

### CATALOG NO. 696

Packed 6 dispensers of 50 each, 300/carton



## STERILE FENESTRATED DRAPE

18" x 26"  
Sterile  
Tissue/poly/tissue  
White/non-glare blue tissue

### CATALOG NO. 697

Packed 6 dispensers of 50 each, 300/carton



## SCRIM-REINFORCED TOWELS



### Blue, general purpose, scrim-reinforced towels

19-1/2" x 23"  
Non-sterile  
Absorbent

### CATALOG NO. 6710

Bulk packed, 200/carton

### CATALOG NO. 6711

Same as Catalog No. 6710,  
but bulk packed 500/carton

19-1/2" x 28"  
Non-sterile  
Absorbent

### CATALOG NO. 6935

Bulk packed, 500/carton



## AUTOCLAVABLE TOWELS/DRAPES

White tissue with blue, non-glare  
polypropylene film center  
**Autoclavable, for use with individual  
treatment trays**

18" x 26"  
Non-sterile  
Tissue/poly/tissue

### CATALOG NO. 695

Bulk packed, 500/carton

### CATALOG NO. 3270

Same as Catalog No. 695,  
but bulk packed 300/carton


Same as Catalog No. 695,  
but with 3" round fenestration.

### CATALOG NO. 694

Bulk packed, 300/carton



# GENERAL SURGICAL DRAPES



## STERILE TISSUE/POLY/TISSUE DRAPE WITH 3" FENESTRATION AND ADHESIVE PATCH

Polyethylene film is sandwiched between two absorbent layers of white tissue to both absorb fluids and prevent strikethrough.

Adhesive around fenestration ensures secure positioning on surgical site.

**CATALOG NO. 698** 18" x 26"  
Packed 4 dispensers of 50 each, 200/carton



## STERILE SMS GENERAL PURPOSE UTILITY DRAPE WITH ADHESIVE STRIP




Blue multi-purpose drapes are perfect add-ons to other drapes, helping to complete the procedural draping field.

Adhesive strip is edge-aligned to minimize drape lift and to prevent any fluid accumulation under the drape.

Ideal for use as a squaring drape to define a sterile field.

General surgical and minor surgical procedures.

**CATALOG NO. 680** 15-1/2" x 25"  
Packed 2 drapes/package, 50 packages/carton



## STERILE NON-WOVEN LAMINATED GENERAL PURPOSE DRAPES

The blue, laminated fabric absorbs fluids on the operative side while providing an impervious barrier to the patient.

Adhesive around fenestration on Catalog No. 687 ensures secure positioning on surgical site.

General surgical and minor surgical procedures.

With 2-1/2" fenestration and adhesive patch  
**CATALOG NO. 687** 22" x 25"  
Packed 4 dispenser boxes of 25 each, 100/carton

Without fenestration and patch  
**CATALOG NO. 686** 22" x 25"  
Packed 2 drapes/package, 80 packages/carton



## GENERAL SURGICAL DRAPES




### STERILE PLASTIC U-DRAPE

U-shaped aperture allows for a snug fit around extremity or head.

Polyethylene film creates an impervious barrier.

Matte finish reduces glare from drape surfaces.

**CATALOG NO. 688** 47" x 51"  
Packed 4 dispensers of 10 each, 40/carton



### STERILE PLASTIC GENERAL PURPOSE UTILITY DRAPES WITH ADHESIVE STRIP

Polyethylene film creates an impervious barrier.

Adhesive strip is edge-aligned to minimize drape lift and to prevent any fluid accumulation under the drape.

Matte finish reduces glare from drape surfaces.

Tinted adhesive liner allows for easy identification and removal.

**CATALOG NO. 683** 17" x 11"  
Packed 4 dispensers of 10 each, 40/carton

**CATALOG NO. 684** 17" x 23"  
Packed 4 dispensers of 10 each, 40/carton



# EPIDURAL/PAIN MANAGEMENT & OPHTHALMIC DRAPES



## STERILE, CLEAR PLASTIC EPIDURAL/PAIN MANAGEMENT DRAPE WITH 4" X 8" RECTANGULAR FENESTRATION AND ADHESIVE STRIP


Large fenestration allows greater access without the need to reposition the drape.

Clear polyethylene film creates an impervious barrier while allowing for visual access to the operative site.

Adhesive strip ensures proper drape placement and prevents drape from moving.

**CATALOG NO. 690** 24" x 30"

Packed 100/carton



## STERILE, CLEAR PLASTIC EPIDURAL/PAIN MANAGEMENT DRAPE WITH 4" FENESTRATION AND ADHESIVE PATCH WITH ADHESIVE STRIPS


Clear polyethylene film creates an impervious barrier while allowing for visual access to the operative site.

Adhesive around fenestration ensures secure positioning on surgical site.

Separate adhesive strips allow for easier handling and placement on patient.

**CATALOG NO. 691** 24" x 30"

Packed 100/carton



## STERILE PLASTIC OPHTHALMIC DRAPE WITH 2-1/2" FENESTRATION AND ADHESIVE PATCH

Adhesive around fenestration ensures secure positioning on surgical site.

Tinted adhesive liner allows for easy identification and removal.

Ultra comfortable film allows drape to adhere to operative site.

Matte finish reduces glare from drape surface

**CATALOG NO. 685** – 15" x 15"

Packed 4 dispensers of 10 each,  
40/carton



## STERILE PLASTIC OPHTHALMIC DRAPE WITH INCISE FILM (3-7/8" X 4-7/8" INCISE AREA) AND FLUID COLLECTION POUCH

Adhesive field isolates the surgical site to reduce the risk of contamination.

Fluid collection pouch minimizes fluid run-off.

Design enables drape to be used for a variety of surgical procedures.

Matte finish reduces glare from drape surface.

Protects patient clothing from contamination during procedures.

**CATALOG NO. 689** – 51" x 51"

Packed 4 dispensers of 10 each,  
40/carton





# ICE BAGS



## DUPONT™ SOFTESSE® ICE BAG

Busse ice bags are offered in both a comfortable, non-woven "Softesse®" synthetic fabric, and a soft, woven "Stockinet" knit fabric. Both styles, available in large and extra large sizes, with and without ties, allow for different size, shape and location applications.

## DUPONT™ SOFTESSE® SOFT-LOC ICE BAG WITH TIES – LARGE

**CATALOG NO. 614** 11" x 8"  
25 bags/box, 2 boxes/carton

## DUPONT™ SOFTESSE® SOFT-LOC ICE BAG WITH TIES – X-LARGE

**CATALOG NO. 615** 11" x 16"  
20 bags/box, 2 boxes/carton

## DUPONT™ SOFTESSE® SOFT-LOC ICE BAG WITHOUT TIES – LARGE

**CATALOG NO. 608** 11" x 8"  
80 bags/carton

## DUPONT™ SOFTESSE® SOFT-LOC ICE BAG WITHOUT TIES – X-LARGE

**CATALOG NO. 609** 11" x 16"  
80 bags/carton

## STOCKINET ICE BAG WITH TIES – LARGE

**CATALOG NO. 618** 9" x 10"  
20 bags/box, 2 boxes/carton

## STOCKINET ICE BAG WITH TIES – X-LARGE

**CATALOG NO. 619** 11" x 13"  
20 bags/box, 2 boxes/carton



The new and improved leakproof ice bags from Busse have been designed to make the caregiver's job faster and easier while providing the patient with cold, yet dry therapy. This innovative ice bag has been thoroughly researched to meet all the clinician's and patient's needs.

## OPENS WIDE

Wider opening allows for easier filling with no spills

Opening fits right under standard ice chute

No need to scoop ice with a funnel or cup

One person can fill the bag without assistance



## LEAKPROOF/COMFORTABLE FIT



### LEAKPROOF

Busse's press, slide and seal closure virtually eliminates leaks. These inner bags are heat-sealed for added strength and assurance.

### COMFORTABLE FIT

Both styles of ice bags are available with soft, durable tie strings designed for secure attachment, as well as maximum comfort. And, there's no plastic clip closures that can irritate skin.

## SIZING



## HOLDS MORE/ SAVES MONEY

One extra large bag does the work of two traditional bags and comfortably conforms to body contours. No discomfort and awkwardness of clip closures.

# O.R. HEAD COVERINGS AND SHOE COVERS



## SPUNBONDED BOUFFANT KOOL-CAPS™

Made of lightweight, but durable spunbonded fabric that permits increased airflow, keeping the wearer cool and comfortable during long surgical procedures.

### 21" Large Bouffant

**CATALOG NO. 327** White

**CATALOG NO. 328** Blue

**CATALOG NO. 329** Green

Packed 5 dispensers of 100 each, 500/carton



### 24" Extra Large Bouffant

for fuller hair styles

**CATALOG NO. 330** White

**CATALOG NO. 331** Blue

**CATALOG NO. 332** Green

Packed 5 dispensers of 100 each, 500/carton



## RAYON BOUFFANT CAPS

Made of lightweight, non-woven rayon material that is comfortable and cool to wear and meets FTC requirements.

### 21" Large Bouffant

**CATALOG NO. 268** Blue

Packed 5 dispensers of 100 each, 500/carton



## SURGEON'S CAP



Made of soft and absorbent material, with wide panel sides, ventilated crown, and adjustable ties.

### CATALOG NO. 335

Packed 5 dispensers of 100 each, 500 per/carton



## SUR-STEP™ ANTI-SKID STYLE SHOE COVERS



Designed with exclusive anti-skid strips applied to the soles that cut down on accidental slips and falls.

Made of durable spunbonded polypropylene

Fits all sizes – small to large

Elastic top ensures snug fit

Liquid repellent and abrasion resistant

Cool and comfortable to wear

### BULK

#### CATALOG NO. 346

Bulk packed, 150 pairs/carton



### DISPENSER PACKED

#### CATALOG NO. 348

Packed 3 dispensers of 50 each, 150 pairs/carton



## EXTRA LARGE SUR-STEP™ ANTI-SKID STYLE SHOE COVER

Fits sizes 13 and above

### CATALOG NO. 347

Bulk packed, 100 pairs/carton



## STANDARD STYLE SHOE COVER



Features quality and durability at an economical price.

### BULK

#### CATALOG NO. 338

Bulk packed, 150 pairs/carton



### DISPENSER PACKED

#### CATALOG NO. 340

Packed 3 dispensers of 50 each, 150 pairs/carton



# HIGH PERFORMANCE ISOLATION / PROCEDURAL MASKS & STAFF PROTECTION KIT



## HIGH FILTRATION, FLUID RESISTANT BLUE MASK WITH FACE SHIELD

Exceeds 99% filtration efficiency at .1 microns. Combines the performance of the high filtration fluid resistant mask with an optically clear plastic wraparound face shield. Provides peripheral protection for the face and eyes against body fluids. Lightly tinted to help eliminate glare. Anti-fog and anti-static treatment. Earloop design

### CATALOG NO. 373

Packed 4 dispensers of 25 each, 100/carton



## HIGH FILTRATION, FLUID RESISTANT BLUE MASK

Fluid resistant mask affords the highest level of fluid and respiratory protection.

Exceeds 99% filtration efficiency at .1 microns. Earloop design. Enclosed nosepiece.

### CATALOG NO. 372

Packed 10 dispensers of 40 each, 400/carton



## HIGH PERFORMANCE YELLOW ISOLATION MASK



Exceeds 99% filtration efficiency at .1 microns.

Ultra soft earloop design. Enclosed nosepiece.

Each layer is ultrasonically bonded to eliminate the sewing pinholes that permit bacteria to migrate from one side to the other.

### CATALOG NO. 370

Packed 10 dispensers of 50 each, 500/carton



## HIGH PERFORMANCE BLUE PROCEDURAL MASK



Exceeds 99% filtration efficiency at .1 microns.

Ultra soft earloop design is comfortable to wear, easy to use.

Enclosed nosepiece for extra softness, extra comfort.

### CATALOG NO. 371

Packed 10 dispensers of 50 each, 500/carton



## DELUXE STAFF PROTECTION KIT



One white, fluid-repellent polyethylene-coated, spunbonded polypropylene gown, with full back, elastic cuffs, and ties at neck and waist

One fluid resistant earloop mask with clear plastic wraparound face shield

One pair of examination gloves

One pair of Sur-Step™ non-skid shoe covers

One 24" blue spunbonded cap

One red waste bag with tie

### CATALOG NO. 198

Packed 20/carton





# STAFF PROTECTION GOWNS



## PERSONAL PROTECTION GOWN

Impervious – Level 4

Thumbhooks prevent gown from riding up and exposing wrists

Pullover design eliminates ties at neck

Perforations in back for fast removal

Universal size

Latex-free

Blue, embossed polyethylene

Meets ASTM F1670-97 requirements for synthetic blood test penetration

Meets ASTM F1671-97B requirements for viral penetration

Ideal for the E.R. or patient rooms

### CATALOG NO. 235

Individually packed  
(15/box, 5 boxes/carton)  
75/carton



## FEATURES IN OUR PERSONAL PROTECTION GOWNS



**THUMBHOOK STIRRUPS**  
Makes glove donning easier and protects against wrist exposure.



**COMFORTABLE & CONVENIENT**  
Open back keeps wearer cool and comfortable.

Perforations allow for fast, easy removal.



**EASY-TO-USE PACKAGING**  
Individually packed 15/box

Hooks onto mounting bracket on most personal protection stations.

## FLUID REPELLENT STAFF PROTECTION GOWNS



Busse staff protection gowns are made of polyethylene-coated, spunbonded polypropylene that is fluid repellent. All gowns have long sleeves, with either elastic or knit cuffs that fit snugly under gloves. They are available in either full back or open back styles.

Gowns are unisex and unisize to fit both men and women, with neck and waist ties for individual size adjustment.

Fluid repellent poly-coated gown, with open back and elastic cuffs

**CATALOG NO. 199** White  
Packed 50/carton (10/bag, 5 bags/carton)

Fluid repellent poly-coated gown, with full back and elastic cuffs.

**CATALOG NO. 200** White  
Packed 50/carton (10/bag, 5 bags/carton)

Fluid repellent, poly-coated gown, with full back and knit cuffs.

**CATALOG NO. 205** White  
Packed 50/carton (10/bag, 5 bags/carton)





# ISOLATION GOWNS



Made with fluid resistant, spunbond polypropylene. Has elastic cuffs and attached neck and waist ties to assure snug fit. Full back. Ideal for visitors in isolation areas.

## STANDARD ISOLATION GOWN

**CATALOG NO. 202** Yellow  
Packed 50/carton (10/bag, 5 bags/carton)



## FLUID REPELLENT ISOLATION GOWN

For maximum protection, fluid repellent spunbond polypropylene poly-coated with full back and knit cuffs.

**CATALOG NO. 203** Yellow  
Packed 50/carton (10/bag, 5 bags/carton)



## X-LARGE ISOLATION GOWN



**Extra room in chest and shoulders provides a more generous fit**

Longer length provides more coverage protection

Tape back placed at top of gown makes it easier and quicker to put on and remove

Elastic cuffs and full back provides complete protection of arms and back

Made from fluid-resistant, spunbonded, polypropylene gown

Tailored fit for a more professional look

**CATALOG NO. 211** Yellow  
Packed 50/carton (10/bag, 5 bags/carton)



## SMS ISOLATION GOWN




SMS is made of meltblown polypropylene sandwiched between two layers of spunbonded polypropylene. **More fluid resistant than common spunbonded material**, the multi-layered SMS apparel offers better protection in a variety of health care situations. Sized generously and softer than regular spunbonded materials, latex-free SMS apparel drapes easily.

Exposure Anticipated: Potential for light fluid contact and contact isolation.

- Elastic cuffs
- Tape back neck closure
- Full back
- Waist ties
- One size fits all

**CATALOG NO. 242** Yellow  
Packed 100/carton (10/bag, 10 bags/carton)



Cross section of SMS apparel showing the meltblown polypropylene sandwiched between two layers of spunbonded polypropylene.

# PROTECTIVE LAB COATS & COVERALLS



## SMS PROTECTIVE LAB COATS

SMS is made of meltblown polypropylene sandwiched between two layers of spunbonded polypropylene. **More fluid resistant than common spunbonded material**, the multi-layered SMS apparel offers better protection for any level of exposure in a variety of health care situations. Sized generously and softer than regular spunbonded materials, latex-free SMS apparel drapes easily.

Exposure Anticipated: Potential for moderate fluid contact and extended use.

- Knit cuffs
- Traditional collar
- Six snap-front closures
- Three reinforced pockets
- Back vent
- Neck loop for hanging
- White

**CATALOG NO. 222**  
Small/Medium

**CATALOG NO. 225**  
Large/X-Large

**CATALOG NO. 229**  
XX-Large

**CATALOG NO. 231**  
XXX-Large

**CATALOG NO. 236**  
XXXX-Large

Packed 25/carton



## LAB COATS




Busse disposable lab coats are made of fluid resistant, spunbonded polypropylene and feature ribbed knit cuffs and collar. Lab coats have snap-front closure and pockets. Available in blue or white.

**CATALOG NO. 224**  
Small/Medium  
Blue  
Packed 30/carton

**CATALOG NO. 226**  
Small/Medium  
White  
Packed 30/carton

**CATALOG NO. 227**  
Large/X-Large  
Blue  
Packed 30/carton

**CATALOG NO. 228**  
Large/X-Large  
White  
Packed 30/carton



Cross section of SMS apparel showing the meltblown polypropylene sandwiched between two layers of spunbonded polypropylene.



## COVERALLS



Busse coveralls are made of white, spunbonded polypropylene. Select catalog numbers available with polyethylene coating to increase fluid repellency.

Generously cut and styled for maximum comfort and mobility.

Nylon zipper front. Available with open or elastic cuffs at both wrist and ankle.

Fluid resistant, open cuffs.  
**CATALOG NO. 212** Large  
**CATALOG NO. 213** X-Large  
Packed 25/carton


Fluid resistant, elastic cuffs.  
**CATALOG NO. 214** Large  
**CATALOG NO. 215** X-Large  
Packed 25/carton

Fluid repellent, poly-coated, elastic cuffs.  
**CATALOG NO. 219** X-Large  
Packed 25/carton



# POST MORTEM BAGS & KITS

## STRAIGHT ZIPPER STYLE



The straight zipper bag is designed with no side seams and top and bottom seams are heat sealed to prevent leakage. The bags are made of strong grade plastic and heavy-duty, rust-proof, nylon zippers. They are all latex-free. Our designs conform to OSHA's universal precautions regarding containment of body fluids and protection from bloodborne pathogens.

### POST MORTEM BAG WITH TAGS

#### CATALOG NO. 901

Packed 10 bags/carton, each in individual poly bag

36" x 90" straight zipper bag  
White  
Three white I.D. tags

### POST MORTEM KIT

#### CATALOG NO. 904

Packed 10 bags/carton, each in individual poly bag

36" x 90" straight zipper bag  
White  
Chin strap  
One underpad  
Two 60" ties  
Three 36" ties  
Three white I.D. tags  
Individual polyethylene bag to hold personal effects



**POST MORTEM BAG KIT COMPONENTS**

### POST MORTEM KIT WITH SHROUD

#### CATALOG NO. 913

Packed 10 bags/carton, each in individual poly bag

Same components of Catalog No. 904, but with 72" x 108" waterproof, opaque linen finish shroud sheet.



### POST MORTEM BAG WITH TAGS

#### CATALOG NO. 902

Packed 10 bags/carton, each in individual poly bag

36" x 90" straight zipper bag, yellow  
Three white I.D. tags  
Three yellow I.D. tag

### POST MORTEM KIT

#### CATALOG NO. 903

Packed 10 bags/carton, each in individual poly bag

36" x 90" straight zipper bag, yellow  
Same components as Catalog No. 904, plus three yellow tags

### INFANT/PEDIATRIC POST MORTEM BAGS

#### CATALOG NO. 905

Packed 10 bags/carton, each in individual poly bag

18-1/2" x 30-1/2" infant straight zipper bag, white  
Three white I.D. tags

#### CATALOG NO. 906

Packed 10 bags/carton, each in individual poly bag

30" x 48" pediatric straight zipper bag, white  
Three white I.D. tags



### PVC-FREE ENVIRO-SAF POST MORTEM KITS

Our Enviro-Saf post mortem bags comply with EPA chlorine-free burn requirements which allow for safe incineration. They also conform to OSHA standards.

#### CATALOG NO. 914

Packed 10 bags/carton, each in individual poly bag

Same components as Catalog No. 901, but with PVC-free components.

#### CATALOG NO. 915

Packed 10 bags/carton, each in individual poly bag

Same components as Catalog No. 904, but with PVC-free components.



# POST MORTEM BAGS & KITS

CURVED ZIPPER STYLE



The curved zipper style bag opens fully for ease of access. All four sides are heat sealed to prevent leakage. The bags are made of strong, high-grade plastic with bound edges and heavy-duty, rust-proof nylon zippers. They are all latex-free. Our designs conform to OSHA's universal precautions regarding containment of body fluids and protection from bloodborne pathogens.

## POST MORTEM BAG WITH TAGS

### CATALOG NO. 907

Packed 10 bags/carton,  
each in individual poly bag

36" x 90" curved zipper bag, white  
Three white I.D. tags

## POST MORTEM KIT

### CATALOG NO. 908

Packed 10 bags/carton,  
each in individual poly bag

36" x 90" curved zipper bag, white  
Chin strap  
One underpad  
Two 60" ties  
Three 36" ties  
Three white I.D. tags  
Individual polyethylene bag to hold personal effects



## POST MORTEM BAG KIT COMPONENTS

## PVC-FREE ENVIRO-SAF POST MORTEM BAG

Our Enviro-Saf post mortem bag complies with EPA chlorine-free burn requirements which allow for safe incineration. It also conforms to OSHA standards.

### CATALOG NO. 916

Packed 10 bags/carton,  
each in individual poly bag

Same components as Catalog No. 907,  
but with PVC-free curved zipper bag.





# POST MORTEM BAGS & KITS

## HEAVY-DUTY WRAP-AROUND ZIPPER STYLE



This white, 36" x 90," heavy-duty post mortem bag features a wrap-around side zipper along the edge of the bag. The bag is made of strong, high-grade latex-free plastic with bound edges and a heavy-duty, rust-proof nylon zipper. Fluid flap is designed to catch excess body fluid run off.

### HEAVY-DUTY POST MORTEM KIT

#### CATALOG NO. 911

10 bags/carton,  
each in individual poly bag

36" x 90" heavy-duty, wrap-around side zipper bag, white

50% heavier than other style bags

Chin strap

Two underpads

Two 60" ties

Three 36" ties

Three white I.D. tags

One set of precaution labels (four per set)

Individual polyethylene bag to hold personal effects



### POST MORTEM BAG KIT COMPONENTS

### X-LARGE, HEAVY DUTY POST MORTEM KIT

#### CATALOG NO. 909

6 bags/carton  
each in individual poly bag

#### 72" x 90" Extra Large heavy-duty

wrap-around side zipper bag

Used for the extremely obese - holds up  
to 550 lbs.

White

Same components as Catalog No. 911



### PVC-FREE ENVIRO-SAF HEAVY DUTY POST MORTEM KIT

Our ENVIRO-SAF heavy-duty post mortem bag complies with EPA chlorine-free burn requirements which allow for safe incineration. It also conforms to OSHA standards.

#### CATALOG NO. 918

6 bags/carton,  
each in individual poly bag

Same components as Catalog No. 911, but  
with PVC Free **heavy-duty** wrap-around  
zipper bag.



# DISASTER BAG AND I.D. KIT/ SHROUD KITS AND SHEETS



## DISASTER BAG

The Busse Disaster Bag is made to meet government specifications as well as OSHA safety standards.

**34" x 96" Black Disaster Bag with I.D. Kit Holds up to 450 lbs.**

Meets US Government Specifications and Department of Defense (DOD) Specifications NSN 9930-01-331-6244

Complies with OSHA Regulation 3130, Occupational Exposure To Blood borne Pathogens – Precautions for Emergency Responders  
Heat sealed seams prevent leakage of fluid from the sides during transport  
Envelope style easy access zipper opening  
20 mil black scrim reinforced PVC for strength  
6 handles for ease of carrying  
Personal effects bag

### CATALOG NO. 920

Packed 5 each/carton



## DISASTER BAG IDENTIFICATION KIT

Used to maintain identification and chain of custody

### CATALOG NO. 921

Packed 20 each/carton

One red biohazard tag  
Three yellow containment tags  
Three identification tags  
One self-adhesive document pouch with biohazard symbol



Our durable shroud sheets and kits are made from PVC-free PE film. All come individually packaged in an outer PE bag that can be used for holding personal effects.

## ADULT KITS AND SHEETS

### CATALOG NO. 725

Packed 50/carton

54" x 108"  
White, opaque linen plastic sheet  
Waterproof  
44" chin strap  
One underpad  
Two 60" ties  
Three 36" ties  
Three white I.D. tags with string ties  
Plastic bag to hold personal effects

### CATALOG NO. 700

Same components as Catalog No. 725, but packed 10/carton

### CATALOG NO. 727

Packed 24/carton

72" x 108" EXTRA LARGE size  
Same components as Catalog No. 725

## CHILD KITS

### CATALOG NO. 701

Packed 10/carton

54" x 72"  
Same components as Catalog No. 725

## SHROUD SHEETS

### CATALOG NO. 702

Packed 50/carton

54" x 108"  
Plastic shroud sheet

### CATALOG NO. 728

Packed 24/carton

72" x 108" Extra Large Size  
Plastic shroud sheet

# BOWLS AND CONTAINERS



## SPONGE BOWLS

Our multi-purpose sponge/solution/utility bowls are made from durable, heavy-weight polypropylene and are extremely versatile. A broad lip at the top of each bowl offers an excellent surface on which to rest instruments, while molded ridges that circle the inside allow instrument tips to rest in solution without immersing the entire instrument. A broad base provides excellent stability.

Busse disposable sponge bowls are an excellent, economical alternative to stainless steel, eliminating all washing, sterilizing, storing, and inventory control problems.

Available in two sizes of opaque, blue polypropylene, 16 oz. and 35 oz.

**CATALOG NO. 400** 16 oz.  
Sterile, individually packed, 75/carton

**CATALOG NO. 401** 35 oz.  
Sterile, individually packed, 50/carton

**CATALOG NO. 402** 16 oz.  
Non-sterile, bulk packed, 500/carton

**CATALOG NO. 403** 35 oz.  
Non-sterile, bulk packed, 200/carton



Each sterile Busse bowl is packed in a see-through plastic bag and sealed with a peel-back strip that opens easily for a struggle-free, aseptic delivery.



## 8 OZ. SPECIMEN CONTAINERS



Plastic container  
8 oz. capacity  
Seal-tite lid to contain odors  
Red container with white lid  
Engraved "Name" and "Room"

**CATALOG NO. 420**

Packed 250/carton



## PLASTIC DENTURE CUP



Lasts for the entire stay of the patient  
Turquoise cup with clear seal-tite lid  
Engraved "Name" and "Room"

**CATALOG NO. 490**

Packed 250/carton



# PATIENT CARE



## PLASTIC APRON

28" x 46"  
White polyethylene with opaque linen finish  
Water repellent  
Oil and chemical resistant  
Helps protect wearer from contact with body fluids  
Individually folded and packed in clear polybag

**CATALOG NO. 250**

Packed 144/carton



## GERIATRIC PLASTIC BIB WITH POCKET



16" x 22"  
White polyethylene with opaque linen finish  
Waterproof  
Attached neckties and dribble pocket  
across bottom

**CATALOG NO. 810**

Packed 500/carton



## PLASTIC PILLOW CASE PROTECTOR



20" x 28"  
Soft, clothlike  
Opaque white linen embossed polyethylene  
Tuck-in end  
Use under cotton pillow case to protect  
pillow from stains

**CATALOG NO. 815**

Packed 200/carton





# CSR WRAPS & ACCESSORIES



## CSR WRAPS

A comprehensive, single-use wrapping system – unparalleled for quality, versatility, and cost-effectiveness. Busse CSR wrap offers superior barrier properties against both air and water-borne bacteria, combined with excellent water repellency and high wet or dry tensile strength.

## STANDARD PACKS

An economical wrap ideally suited for general purpose use in hospitals, physician's offices, clinics, etc.

CAT NO.	SIZE	CARTON PACK
850	12" x 12"	1000
851	15" x 15"	500
852	20" x 20"	500
853	24" x 24"	500
854	30" x 30"	250
855	36" x 36"	250
856	40" x 40"	250
857	45" x 45"	250
858	54" x 54"	100

## PHYSICIAN CONVENIENCE PACKS

Each pack contains 100 sheets in the six most commonly used sizes.

CAT NO.	SIZE	CARTON PACK
845	12" x 12"	100
843	15" x 15"	100
846	18" x 18"	100
847	24" x 24"	100
848	30" x 30"	100
849	36" x 36"	100



## TAPE MEASURES



Newborn 24" – shows inches and millimeters/centimeters

### CATALOG NO. 791

Packed 1,000/carton



Pediatric 36" – shows inches and millimeters/centimeters

### CATALOG NO. 795

Packed 1,000/carton



## 1/2 OZ. FECAL CONTAINER



Opaque black plastic  
Tight fitting, hinged lid  
Disposable spatula and envelope  
Envelopes imprinted: Number, Date, Address, Age, Sex, and Result

### CATALOG NO. 450

Packed 1,000/carton



## PROFESSIONAL TOWELS – ABSORBENT



13" x 19" folded to 5" x 7".

3-ply wet strength facial quality tissue


Suitable for a multitude of uses, including burping cloths, pillow covers, wipers, bibs, bassinet sheets, nursery scale liners, bedpan drapes, napkins, sterile wraps, tray covers

### CATALOG NO. 802

Packed 500/carton



## PLASTIC MICROSCOPE COVER SLIPS



Clear, non-curling, unbreakable plastic cover slips. For use in immediate wet preparations such as urine, direct feces smears, fats, parasites, etc.

### CATALOG NO. 736

Size 22mm x 22mm  
Packed 10 boxes of 100 each per dispenser  
15 dispensers per carton  
15,000 per carton



# PRODUCT INDEX

CAT. NO. DESCRIPTION	PAGE
133 Curette, Straight	28
134 Curette, Straight	28
135 Curette, Straight	28
136 Curette, Straight	28
137 Curette, Straight	28
138 Curette, Straight	28
139 Curette, Straight	28
140 Ear/Ulcer Syringe	31
141 Ear/Ulcer Syringe	31
142 Ear/Ulcer Syringe	31
143 Ear/Ulcer Syringe	31
144 Ear/Ulcer Syringe	31
145 Ear/Ulcer Syringe	31
146 Ear/Ulcer Syringe	31
147 Ear/Ulcer Syringe	31
148 Ear/Ulcer Syringe	31
149 Ear/Ulcer Syringe	31
150 Tubing, Suction	28
151 Tubing, Suction	28
153 Tubing, Suction	28
154 Tubing, Connective	26
155 Tubing, Connective	26
156 Tubing, Connective	26
157 Tubing, Connective	26
162 Suction Tip, Frazier	26
163 Suction Tip, Frazier	26
164 Suction Tip, Frazier	26
165 Suction Tip, Sigmoidoscopic	26
167 Suction Tip, Sigmoidoscopic	26
198 Staff Protection Kit	39
199 Staff Protection Gown	40
200 Staff Protection Gown	40
202 Isolation Gown	41
203 Isolation Gown	41
205 Staff Protection Gown	40
206 Laminaria Cervical Dilator	28
207 Laminaria Cervical Dilator	28
208 Laminaria Cervical Dilator	28
209 Laminaria Cervical Dilator	28
211 Isolation Gown	41
212 Coveralls	42
213 Coveralls	42
214 Coveralls	42
215 Coveralls	42
219 Coveralls	42
222 SMS Protective Lab Coat	42
224 Lab Coat	42

CAT. NO. DESCRIPTION	PAGE
225 SMS Protective Lab Coat	42
226 Lab Coat	42
227 Lab Coat	42
228 Lab Coat	42
229 SMS Protective Lab Coat	42
231 SMS Protective Lab Coat	42
235 Personal Protection Gown	40
236 SMS Protective Lab Coat	42
242 SMS Isolation Gown	41
250 Apron	48
268 Rayon Bouffant Cap	38
279 Curette, Curved	28
280 Curette, Curved	28
281 Curette, Curved	28
282 Curette, Curved	28
283 Curette, Curved	28
284 Curette, Curved	28
285 Curette, Curved	28
290 Suction Tip, Poole	25
291 Suction Tip, Poole	25
296 Suction Tip, Yankauer	25
297 Suction Tip, Yankauer	25
298 Suction Tip, Yankauer	25
299 Suction Tip, Yankauer	25
302 Suction Tip with Tubing, Yankauer	25
303 Suction Tip with Tubing, Yankauer	25
305 Suction Tip with Tubing, Yankauer	25
327 Spunbonded Bouffant Kool-Cap	38
328 Spunbonded Bouffant Kool-Cap	38
329 Spunbonded Bouffant Kool-Cap	38
330 Spunbonded Bouffant Kool-Cap	38
331 Spunbonded Bouffant Kool-Cap	38
332 Spunbonded Bouffant Kool-Cap	38
335 Surgeon's Cap	38
338 Shoe Cover	38
340 Shoe Cover	38
346 Shoe Cover	38
347 Shoe Cover	38
348 Shoe Cover	38
370 Isolation Mask	39
371 Procedural Mask	39
372 Fluid-Resistant Mask	39
373 Fluid-Resistant Mask w/ Face Shield	39
384 Umbilical Cord Clamp	32
385 Umbilical Cord Clamp	32
386 Circumcision Tray	32
387 Amniotic Membrane Perforator	32
388 Obstetrical Kit	32
389 Obstetrical Kit	32
400 Sponge Bowls	47
401 Sponge Bowls	47
402 Sponge Bowls	47
403 Sponge Bowls	47

CAT. NO. DESCRIPTION	PAGE
405 Mucus Specimen Traps	29
406 Mucus Specimen Traps	29
407 Mucus Specimen Traps	29
410 Catheter Kit	30
411 Catheter Kit	30
412 Catheter Kit	30
413 Catheter Kit	30
420 Specimen Container	47
450 Fecal Container	49
460 Calculi Strainer	30
490 Denture Cup	47
501 Tubing Connector, 5-in-1	27
502 Tubing Connector, Sims	27
503 Tubing Connector, "T"	27
504 Tubing Connector, Straight	27
505 Tubing Connector, Straight	27
506 Tubing Connector, Straight	27
507 Tubing Connector, "Y"	27
508 Tubing Connector, "Y"	27
509 Tubing Connector, 6-in-1 "Y"	27
510 Catheter Plug	27
511 Tubing Connector, 5-in-1	27
512 Tubing Connector, Sims	27
513 Tubing Connector, "T"	27
514 Tubing Connector, Straight	27
515 Tubing Connector, Straight	27
516 Tubing Connector, Straight	27
517 Tubing Connector, "Y"	27
518 Tubing Connector, "Y"	27
519 Tubing Connector, 6-in-1 "Y"	27
546 E.R. Laceration Tray	20
549 E.R. Laceration Tray	20
558 Incision and Drainage Tray	18
560 Spinal Needle, Quincke-Style	6
561 Spinal Needle, Quincke-Style	6
562 Spinal Needle, Quincke-Style	6
563 Spinal Needle, Quincke-Style	6
564 Spinal Needle, Quincke-Style	6
565 Spinal Needle, Quincke-Style	6
566 Spinal Needle, Quincke-Style	6
595 Loss of Resistance Syringe	6
596 Loss of Resistance Syringe	6
608 Ice Bag, Softesse®	37
609 Ice Bag, Softesse®	37
614 Ice Bag, Softesse®	37
615 Ice Bag, Softesse®	37
618 Ice Bag, Stockinet	37
619 Ice Bag, Stockinet	37
620 Tumescant Delivery System	7
621 Vein Closure Tray, Basic	7
622 Vein Closure Tray, Deluxe	7
630 Bone Marrow Biopsy/Aspiration Ndl	13
631 Bone Marrow Biopsy/Aspiration Ndl	13

CAT. NO. DESCRIPTION	PAGE	CAT. NO. DESCRIPTION	PAGE	CAT. NO. DESCRIPTION	PAGE
640 Nosebleed Tray	18	732 Suture Removal Kit	22	854 CSR Wrap	49
644 Manometer	9	736 Microscope Cover Slip	49	855 CSR Wrap	49
645 Lumbar Puncture Tray	8	737 Suture Removal Kit	22	856 CSR Wrap	49
646 Lumbar Puncture Tray	8	744 Sharp Debridement Tray	18	857 CSR Wrap	49
647 Lumbar Puncture Tray	9	746 E.R. Laceration Tray	20	858 CSR Wrap	49
648 Lumbar Puncture Tray	9	747 Suturing Tray	20	901 Post Mortem Kit	43
649 Lumbar Puncture Tray	9	748 E.R. Laceration Kit	20	902 Post Mortem Bag	43
650 Amniocentesis Tray	10	749 E.R. Laceration Tray	20	903 Post Mortem Kit	43
651 Basic Soft Tissue Biopsy Tray	10	750 Minor Laceration Tray	19	904 Post Mortem Kit	43
652 Bone Marrow Biopsy/Aspiration Tray	12	751 Minor Laceration Tray	19	905 Post Mortem Bag	43
653 Bone Marrow Biopsy/Aspiration Tray	12	752 General Purpose Instrument Tray	19	906 Post Mortem Bag	43
654 Bone Marrow Biopsy/Aspiration Tray	12	753 General Purpose Instrument Tray	19	907 Post Mortem Bag	44
655 Bone Marrow Biopsy/Aspiration Tray	12	754 General Purpose Instrument Tray	19	908 Post Mortem Kit	44
656 Myelogram Tray	11	755 General Purpose Instrument Tray	19	909 Post Mortem Kit	45
670 Epidural Tray, Single Shot	4	756 Deluxe Wound Closure Tray	18	911 Post Mortem Kit	45
671 Epidural Tray, Continuous	4	757 Deluxe Wound Closure Tray, Facial	18	913 Post Mortem Kit	43
672 Basic Nerve Block Tray	4	758 Incision and Drainage Tray	18	914 Post Mortem Kit	43
680 Drape, General Purpose, Sterile	34	759 Dressing Change Tray	16	915 Post Mortem Kit	43
683 Drape, General Purpose, Sterile	35	766 Irrigation Splash Guard	20	916 Post Mortem Bag	44
684 Drape, General Purpose, Sterile	35	767 Irrigation Splash Guard	20	918 Post Mortem Bag	45
685 Drape with Fenestration, Sterile	36	770 Skin Staple Remover Kit	21	920 Disaster Bag	46
686 Drape, General Purpose, Sterile	34	772 Skin Staple Remover	21	921 Disaster Bag I.D. Kit	46
687 Drape, General Purpose, Sterile	34	773 Nurses' Scissors	21	946 Lumbar Puncture Tray	8
688 U-Drape, Sterile	35	774 Hemostat	21	956 Myelgram	11
689 Drape, Ophthalmic, Sterile	36	775 Hemostat	21	957 Arthrogram Tray	11
690 Drape with Fenestration, Sterile	36	791 Tape Measure	32, 49	960 Thoracentesis Tray	14
691 Drape with Fenestration, Sterile	36	795 Tape Measure	32, 49	961 Thoracentesis Tray	14
694 Autoclavable Towels/Drapes	33	799 Tracheostomy Care Set	23	962 Paracentesis Tray	14
695 Autoclavable Towels/Drapes	33	800 Tracheostomy Care Set	23	974 Lumbar Puncture Tray	9
696 Sterile Field	33	802 Towel	49	975 Lumbar Puncture Tray	9
697 Drape, General Purpose	33	810 Bib	48	976 Amniocentesis Tray	10
698 Drape with Fenestration, Sterile	34	815 Pillow Case Protector	48	977 Basic Soft Tissue Biopsy Tray	10
700 Shroud Kit	46	818 I.V. Start Kit	17	978 Myelogram Tray	11
701 Shroud Kit	46	820 I.V. Start Kit	17	979 Arthrogram Tray	11
702 Shroud Sheet	46	822 I.V. Start Kit	17	985 Bone Marrow Biopsy/Aspiration Tray	13
703 Tracheostomy Care Set	23	829 Central Venous Catheter Dressing Change Tray	16	986 Bone Marrow Biopsy/Aspiration Tray	13
706 Tracheostomy Care Set	24	832 Central Venous Catheter Dressing Change Tray	16	987 Bone Marrow Biopsy/Aspiration Tray	13
707 Tracheostomy Dressing	24	834 Central Line Dressing Change Tray	16	1443 Primary Dressing Change Tray	16
710 Tracheostomy Care Set	24	835 Central Line Kit	15	1539 Tubing, Cosmetic Surgery	26, 28
711 Tracheostomy Care Set	24	836 Central Line Kit	15	1540 Tubing, Cosmetic Surgery	26, 28
713 Tracheostomy Care Set	24	837 Central Line Kit	15	3220 Minor Procedure Pack	19
714 Tracheostomy Care Set	23	843 CSR Wrap	49	3270 Autoclavable Towels/Drape	33
715 Tracheostomy Care Set	23	845 CSR Wrap	49	6710 Scrim-Reinforced Towel	33
716 Skin Staple Remover Kit	21	846 CSR Wrap	49	6711 Scrim-Reinforced Towel	33
717 Suture Removal Kit	22	847 CSR Wrap	49	6935 Scrim-Reinforced Towel	33
718 Suture Removal Kit	22	848 CSR Wrap	49	7190 Forceps	21
723 Suture Removal Kit	22	849 CSR Wrap	49	7195 Forceps	21
724 Suture Removal Kit	22	850 CSR Wrap	49	7246 Nurses' Scissors	21
725 Shroud Kit	46	851 CSR Wrap	49		
726 Suture Removal Kit	22	852 CSR Wrap	49		
727 Shroud Kit	46	853 CSR Wrap	49		
728 Shroud Sheet	46				
729 Suture Removal Kit	22				

**busse**  
**Hospital Disposables**

P.O. Box 11067, Hauppauge, NY 11788-0920  
800.645.6526 • Fax: 631.435.4721  
mail@busseinc.com  
[www.busseinc.com](http://www.busseinc.com)