

# ORTHOPEDIC CATALOG



# TABLE OF CONTENTS

Foam Collars.....	1-4
Philadelphia Collars.....	5
Extrication Collars .....	6-7
Post-Op Collars.....	8
Cervical Traction .....	9-12
Clavicle Straps .....	13-14
Shoulder Immobilizers .....	15-18
Shoulder Abduction Pillows.....	19-20
Arm Slings/Elevators.....	21-25
Elbow .....	26-28
Post-Op Elbow .....	29
Wrist & Wrist/Forearm .....	30-36
Thumb .....	37-38
Aluminum Splints.....	39-47
Abdominal Binders.....	48-49
Rib Belts.....	50
Lumbo Sacral Supports .....	51-53
Pelvic Traction.....	54
Spine.....	55-57
Knee Immobilizers .....	58-64
Knee Supports .....	65-68
Post-Op Knee.....	69-70
Bucks Traction Boots .....	71-73
Ankle .....	74-78
Walkers .....	79-82
Post-Op Shoes .....	83-86
Cast Shoes .....	87-89
Foot & Heel.....	90
Pressure Reduction.....	91-94
Foam Positioners .....	95-98
Hot & Cold Therapy .....	99-102
Specialty Nursing .....	103-106
Bandages & Wraps .....	107
Ambulatory Aids .....	108
Pediatric.....	109

Dear Valued Customer:

For over thirty-five years, we have been privileged to stand quietly behind you in your emergency departments and orthopedic hospital floors, your offices and clinics, your rehab centers and physical therapy labs.

We are proud to be the largest supplier of orthopedic softgoods to hospitals and clinics in the nation. The diversity of our product line, the quality inherent in our manufacturing process and our committed staff enables us to offer our customers exceptional products at reasonable prices.

Our commitment continues - a commitment to providing patients with the best products available today. Using our own raw materials and many modern manufacturing techniques, we make comfortable products that fit, products patients can afford and products which help improve their quality of life.

We hope this catalog provides you with all of the information you need to order and use our products. For more information, please contact your local DeRoyal Sales Representative. To order, either fax the order to 800.543.2182, e-mail [customerservice@deroyal.com](mailto:customerservice@deroyal.com) or call Customer Service directly at 800.251.9864.

And, as in any successful partnership, please let us know how we can work better with you. We strive for continuous improvement based on customer input, and we very much look forward to hearing from you.



Pete DeBusk  
Chairman and CEO, DeRoyal

# FOAM COLLARS

## Contoured Cervical Collar

- Medium density foam construction
- Breathable cotton stockinette liner
- Removable cover can be washed for better patient hygiene
- Deep contour holds the head in a neutral position

**Indications:** Cervical facet syndrome; Cervical sprains and strains; Discogenic pain; Degenerative disc disease with and without radiculopathy

Size	Length	2½"	3"	4"
Universal	22"	1061-00	1057-00	1056-00
X-Small	13"	1061-01	1057-01	1056-01
Small	15"	1061-02	1057-02	1056-02
Medium	17"	1061-03	1057-03	1056-03
Large	18"	1061-04	1057-04	1056-04
X-Large	20"	1061-05	1057-05	1056-05



## Medium Density Cervical Collar

- Medium density foam construction
- Breathable cotton stockinette liner
- Removable cover can be washed for better patient hygiene

**Indications:** Cervical facet syndrome; Cervical sprains and strains; Discogenic pain; Degenerative disc disease with and without radiculopathy

Size	Length	2½"	3"	3½"	4"
Universal	22"	1000251	1000301	1000351	1000401
Small	16"	1010251	1010301	1010351	1010401
Medium	18"	1020251	1020301	1020351	1020401
Large	20"	1030251	1030301	1030351	1030401
X-Large	22"	1040251	1040301	1040351	1040401



## Medium/Firm Density Cervical Collar

- Medium/firm density foam construction
- Breathable cotton stockinette liner
- Removable cover can be washed for better patient hygiene
- Contour holds the head in a neutral position

**Indications:** Cervical facet syndrome; Cervical sprains and strains; Discogenic pain; Degenerative disc disease with and without radiculopathy

Size	Length	2½"	3"	3½"	4"
Universal	22"	1050-00	1050-10	1050-20	1050-30
Small	16"	1050-01	1050-11	1050-21	1050-31
Medium	18"	1050-02	1050-12	1050-22	1050-32
Large	20"	1050-03	1050-13	1050-23	1050-33
X-Large	22"	1050-04	1050-14	1050-24	1050-34



# FOAM COLLARS



## Comfo-Eze Cervical Collar

- Highest quality medium density foam construction
- Breathable synthetic stockinette liner
- Removable cover can be washed for better patient hygiene
- Available in wide and narrow styles

**Indications:** Cervical facet syndrome; Cervical sprains and strains; Discogenic pain; Degenerative disc disease with and without radiculopathy

Size	Dimensions	Wide
X-Small	3½" x 19½"	1001-01
Small	4" x 20½"	1001-02
Medium	4" x 22½"	1001-03
Large	4" x 23½"	1001-04
X-Large	4½" x 25½"	1001-05
Universal	4¾" x 22"	1010-01

Size	Dimensions	Regular
X-Small	3" x 19½"	1002-01
Small	3½" x 20½"	1002-02
Medium	3½" x 22½"	1002-03
Large	3½" x 23½"	1002-04
X-Large	3½" x 25½"	1002-05
Universal	3½" x 22"	1010-02

Size	Dimensions	Narrow
X-Small	2" x 19½"	1003-01
Small	2" x 20½"	1003-02
Medium	2" x 22½"	1003-03
Large	2" x 23½"	1003-04
X-Large	2" x 25½"	1003-05

Size	Dimensions	Product #
Pediatric Small	2½" x 14"	1004-01
Pediatric Large	3½" x 17"	1004-02



## Medium Density Cervical Collar

- Medium density foam construction
- Breathable synthetic stockinette liner
- Removable cover can be washed for better patient hygiene

**Indications:** Cervical facet syndrome; Cervical sprains and strains; Discogenic pain; Degenerative disc disease with and without radiculopathy

Size	Dimensions	Product #
Small	3" x 16½"	1005-02
Small, X-Long	3" x 19"	1005-06
Small, X-Long, X-Wide	3½" x 19½"	1005-01
Medium	4" x 20½"	1005-03
Medium, X-Long	4" x 22"	1005-05
Large	5" x 21½"	1005-04

# FOAM COLLARS

## Serpentine Cervical Collar

- Firm density foam construction
- Breathable synthetic stockinette liner
- Removable cover can be washed for better patient hygiene
- Contoured serpentine design for patient comfort

**Indications:** Cervical facet syndrome; Cervical sprains and strains; Discogenic pain; Degenerative disc disease with and without radiculopathy

Size	Dimensions	Product #
Small	3½" x 17"	1006-02
Small, X-Wide	3¾" x 18"	1006-01
Small, X-Long	3¼" x 24½"	1006-05
Medium	4" x 20"	1006-03
Large	4" x 24"	1006-04



## Cerv-o-Collar

- Firm density foam construction
- Breathable synthetic stockinette liner
- Removable cover can be washed for better patient hygiene
- Designed to hyperextend chin

**Indications:** Cervical facet syndrome; Cervical sprains and strains; Discogenic pain; Degenerative disc disease with and without radiculopathy

Size	Dimensions	Product #
Small	3" x 15½"	1007-02
Medium	3½" x 21½"	1007-03
Medium, Short	3½" x 19"	1007-04
Medium, X-Long	3½" x 23½"	1007-06
Large	4½" x 21½"	1007-05



## Medium Density Cervical Collar

- Medium density foam construction
- Breathable synthetic stockinette liner
- Removable cover can be washed for better patient hygiene
- Contoured design holds the head in a neutral position

**Indications:** Cervical facet syndrome; Cervical sprains and strains; Discogenic pain; Degenerative disc disease with and without radiculopathy

Size	Dimensions	Product #
Small, Narrow	2" x 22"	1063-06
Small	3½" x 19"	1063-01
Small, X-Long	3½" x 23"	1063-02
Medium	4" x 21"	1063-03
Medium, X-Long	4" x 24"	1063-04
Large	4½" x 24"	1063-05



# FOAM COLLARS



## Soft Foam Cervical Collar

- Soft density foam construction
- Breathable cotton stockinette liner
- Packaged 6/Case

**Indications:** Cervical facet syndrome; Cervical sprains and strains; Discogenic pain; Degenerative disc disease with and without radiculopathy

Size	Dimensions	Product #
X-Small	3" x 19½"	1018-01
Small	3½" x 20½"	1018-02
Medium	3½" x 22½"	1018-03
Large	3½" x 23½"	1018-04
X-Large	3½" x 25½"	1018-05
Universal	3" x 22"	1018-10
Universal Narrow	3½" x 22"	1018-11
Universal Wide	4¼" x 22"	1018-12



# PHILADELPHIA COLLARS

## Philadelphia® Cervical Collar

- 2-piece design restricts flexion, extension and rotation while supporting the occipital and mandibular areas
- Non-toxic, hypo-allergenic and latex-safe Plastazote® foam reduces skin irritation
- Large trachea opening allows quick access for carotid pulse monitoring, emergency airway procedures and visualization
- Posterior ventilation holes reduce heat and moisture
- X-ray, CT and MRI lucent

**Indications:** Provides superior immobilization for both extrication and rehabilitation; Post-operative discectomy and cervical fusion



Size	Circumference	1 1/4"	1 3/4"	2 1/4"	3 1/4"	4 1/4"	5 1/4"
Infant	6" - 8"	1023-01	—	—	—	—	—
Pediatric	8" - 11"	—	1023-00	—	—	—	—
Small	10" - 13"	—	—	1024-11	1034-11	1044-11	1054-11
Medium	13" - 16"	—	—	1024-12	1034-12	1044-12	1054-12
Large	16" - 19"	—	—	1024-13	1034-13	1044-13	1054-13
X-Large	19" - 22"	—	—	1024-14	1034-14	1044-14	1054-14

## Definitive Collar

- 2-piece design restricts flexion, extension and rotation while supporting the occipital and mandibular areas
- Large trachea opening allows quick access for carotid pulse monitoring, emergency airway procedures and visualization
- Pre-formed molded foam is non-toxic and hypo-allergenic
- Contoured edges to increase patient comfort
- X-ray, CT and MRI lucent

**Indications:** Provides superior immobilization for both extrication and rehabilitation; Post-operative discectomy and cervical fusion



Size	Circumference	2"	2 1/4"	3 1/4"	4 1/4"	5 1/4"
Pediatric	8" - 11"	1023-50	—	—	—	—
Small	10" - 13"	—	1024-51	1034-51	1044-51	1054-51
Medium	13" - 16"	—	1024-52	1034-52	1044-52	1054-52
Large	16" - 19"	—	1024-53	1034-53	1044-53	1054-53



# EXTRICATION COLLARS



## Philly EMT Select Extrication Collar

- 1-piece flat design allows efficient storage where space is limited
- Universal/adjustable to four adult and three pediatric sizes
- Locking tabs, slide adjustment and reflective sizing indicators for maximum hold
- Closed cell foam liner to help prevent skin breakdown
- Large trachea opening allows easy access for carotid pulse monitoring, emergency airway procedures and visualization
- X-ray, CT and MRI lucent

**Indications:** Provides superior immobilization for both extrication and rehabilitation; Cervical trauma

Size	Front	Circumference	Qty	Product #
Adult	3½"	11" - 23"	1 Each	1044-31
			20/Case	1044-30
			100/Case	1044-90
Pediatric	1½"	8" - 18"	1 Each	1044-41
			20/Case	1044-40
			100/Case	1044-91



## Philly One-Piece Extrication Collar

- 1-piece flat design allows efficient storage where space is limited
- Four adult and two pediatric sizes ensure accurate fit
- Color-coded for quick size reference
- Large trachea opening allows easy access for carotid pulse monitoring, emergency airway procedures and visualization
- X-ray, CT and MRI lucent

**Indications:** Provides superior immobilization for both extrication and rehabilitation; Cervical trauma

Size	Front	Circumference	Product #
Infant	1½"	8" - 18"	1044-80
Pediatric	1¾"	8" - 18"	1044-81
Low	1¾"	11" - 23"	1044-82
Short	2½"	11" - 23"	1044-83
Regular	3"	11" - 23"	1044-84
Tall	3½"	11" - 23"	1044-85

# EXTRICATION COLLARS

## STAT® Neck Extrication Collar

- 2-piece design and construction
- Maneuverable chin piece allows for a precise fit
- Non-absorbent padding for patient comfort
- Color-coded for quick size reference
- X-ray, CT and MRI lucent

**Indications:** Provides superior immobilization for both extrication and rehabilitation; Cervical trauma

Size	Front	Circumference	Product #
Pediatric	2½"	10"- 15"	1043-01
X-Small	3½"	10"- 16"	1042-00
Small	3½"	12"- 17"	1043-02
Small, X-Long	3½"	14"- 19"	1043-00
Medium	4½"	13"- 18"	1043-03
Large	5½"	14"- 19"	1043-04



# POST-OP COLLARS



## XTW Extended Wear Collar

- Designed to help restrict flexion, extension, lateral bending and axial rotation
- High density polyethylene frame for maximum rigidity and support
- Shoulder and mandible contours for increased stability
- Anterior and posterior flex tabs for individual comfort
- Latex-safe, non-toxic, open-cell liner helps prevent skin breakdown
- Five adult and five child sizes to ensure best patient fit
- X-ray, CT and MRI lucent
- Replacement pad set available upon request

**Indications:** Cervical sprains and strains; Cervical spinal stenosis; Cervical disc displacement; Cervical arthritis, Whiplash injury; Cervical spondylosis and spondylolisthesis; Degenerative disc disease; Cervicalgia

Size	Collar	Collar & Pad Set
Adult X-Short	11201000	11202000
Adult Short	11211000	11212000
Adult Regular	11221000	11222000
Adult Tall	11231000	11232000
Adult X-Tall	11241000	11242000
Infant Short	11251000	11252000
Infant Regular	11261000	11262000
Child Tot	11271000	11272000
Child Short	11281000	11282000
Child Regular	11291000	11292000

# CERVICAL TRACTION

## Deluxe Head Halter

- Constructed of perforated foam material with tricot lining
- Designed to provide anterior tilt cervical traction
- Halter is easily adjustable with hook and loop closures

**Indications:** Degenerative disc disease; Cervical facet syndrome; Herniated nucleus pulposus with or without radiculopathy; Chronic cervical strains

Size	Product #
Small	6101-01
Medium	6101-02
Large	6101-03



## Repro Type Universal Head Halter

- Made of durable canvas material with tricot lining
- Pulls from the occipital area to eliminate chin pressure
- Halter tilts head slightly during traction
- Chin straps have buckle closures for easy adjustment

**Indications:** Degenerative disc disease; Cervical facet syndrome; Herniated nucleus pulposus with or without radiculopathy; Chronic cervical strains

Size	w/o Spreader Bar	w/ 12" Spreader Bar
Universal	6102-00	6114-00



## Upper-7 Type Universal Head Halter

- Made of durable canvas material with tricot lining
- Unique design does not cover ears
- Minimizes amount of pressure to the patient's head
- Chin straps have buckle closures for easy adjustment

**Indications:** Degenerative disc disease; Cervical facet syndrome; Herniated nucleus pulposus with or without radiculopathy; Chronic cervical strains

Size	Product #
Universal	6103-00



# CERVICAL TRACTION



## Discard Head Halter

- Available in cotton or flannel material
- Product can be applied for hyperextension or flexion
- Spreader bar sold separately
- Packaged 12/Case

**Indications:** Degenerative disc disease; Cervical facet syndrome; Herniated nucleus pulposus with or without radiculopathy; Chronic cervical strains

Size	Cotton	Flannel
Universal	6104-00	6105-00



## Standard Over-Door Traction Unit

- Unit comes with standard discard flannel head halter
- 12" spreader bar, traction cord and vinyl water bag included
- Fits easily on any door

**Indications:** Degenerative disc disease; Cervical facet syndrome; Herniated nucleus pulposus with or without radiculopathy; Chronic cervical strains

Size	w/ 10 lb. Bag	w/ 20 lb. Bag
Universal	7102-00	7102-22



## Exo-Static Over-Door/Exo-Bed Traction Unit

- Unit comes with standard discard flannel head halter
- Includes 12" spreader bar, traction cord and vinyl water bag
- Extra support on over-door bar and bracket for increased security
- Exo-bed unit designed for bed-ridden patients

**Indications:** Degenerative disc disease; Cervical facet syndrome; Herniated nucleus pulposus with or without radiculopathy; Chronic cervical strains

<b>Over Door Unit</b>	<b>w/ 10 lb. Bag</b>	<b>w/ 20 lb. Bag</b>
Universal	7101-00	7101-20
<b>Bed Unit</b>	<b>w/ 10 lb. Bag</b>	
Universal	7103-00	



# CERVICAL TRACTION

## Over-Door Traction Kit w/ Deluxe Head Halter

- Unit comes with deluxe head halter
- 12" spreader bar, traction cord and vinyl water bag included
- Extra support on over-door bar and bracket for increased security

**Indications:** Degenerative disc disease; Cervical facet syndrome; Herniated nucleus pulposus with or without radiculopathy; Chronic cervical strains

Size	w/ 10 lb. Bag	w/ 20 lb. Bag
Medium	7104-02	7104-22
Large	—	7104-23



## Over-Door Traction Kit w/ Repro Head Halter

- Unit comes with universal repro head halter
- 12" spreader bar, traction cord and vinyl water bag included
- Extra support on over-door bar and bracket for increased security

**Indications:** Degenerative disc disease; Cervical facet syndrome; Herniated nucleus pulposus with or without radiculopathy; Chronic cervical strains

Size	w/ 10 lb. Bag	w/ 20 lb. Bag
Universal	7115-00	7115-20



## Traction Cord

- Inner core constructed of polyester
- Polypropylene outer weave
- 600 pound tensile strength
- 1/4" diameter

**Indications:** Traction

Length	Product #
100'	7116-01



# CERVICAL TRACTION



## Weight Bag

- Durable vinyl bag
- Graduated marks for weight adjustment

**Indications:** Traction

Size	Product #
10 lb.	7117-01
20 lb.	7117-02



## Spreader Bar

- Extra strong and stable
- Includes 78" of traction cord

**Indications:** Traction

Length	Product #
12"	6117-00



# CLAVICLE STRAPS

## Premium Clavicle Strap

- Synthetic stockinette over foam padding does not absorb water and can be worn while bathing
- Serpentine stitch design to prevent migration
- Available in buckle or hook and loop closure

**Indications:** Clavicle fractures; Grade II and III acromioclavicular separations; Osteolysis of the clavicle; AC degenerative joint pain; Mechanical cervical disease; Postural dysfunction; Sternoclavicular dislocations

Size	Buckle Closure	Hook & Loop Closure
Infant	3001-01	3002-01
Child	3001-02	3002-02
Small	3001-06	3002-06
Medium	3001-03	3002-03
Large	3001-04	3002-04
X-Large	3001-05	3002-05



## Heavy-Duty Clavicle Strap

- Constructed of foam padding with stockinette covering
- Extra long straps for easy adjustment
- Available in buckle or hook and loop closure

**Indications:** Clavicle fractures; Grade II and III acromioclavicular separations; Osteolysis of the clavicle; AC degenerative joint pain; Mechanical cervical disease; Postural dysfunction; Sternoclavicular dislocations

Size	Buckle Closure	Hook & Loop Closure
Infant	3004-01	3004-11
Child	3004-02	3004-12
Small	3004-06	3004-16
Medium	3004-03	3004-13
Large	3004-04	3004-14
X-Large	3004-05	3004-15



## 3-Way Clavicle Strap

- Foam construction with cotton stockinette cover
- 3-way posterior vector helps prevent migration
- Available in buckle or hook and loop closure

**Indications:** Clavicle fractures; Grade II and III acromioclavicular separations; Osteolysis of the clavicle; AC degenerative joint pain; Mechanical cervical disease; Postural dysfunction; Sternoclavicular dislocations

Size	Buckle Closure	Hook & Loop Closure
Child	3015-01	3015-11
Small	3015-02	3015-12
Medium	3015-03	3015-13
Large	3015-04	3015-14
X-Large	3015-05	3015-15



# CLAVICLE STRAPS



## Acromioclavicular Splint

- Constructed of closed cell material that does not absorb moisture and comfortably pads shoulder area
- Provides upward pressure on forearm and elbow for better alignment of the G-H joint
- D-ring/clip closures
- Universal right/left and size to fit majority of patients

**Indications:** S/P glenohumeral acromioplasty; Surgical repairs of the rotator cuff; Anterior glenohumeral capsular reconstruction; Anterior dislocation of the G-H joint

**Size**  
Universal

**Product #**  
3007-00

# SHOULDER IMMOBILIZERS

## Deluxe Sling & Swathe

- Lightweight perforated foam material provides patient comfort
- Adjustable shoulder straps attach to hook receivable sling to ensure accurate fit
- Open or closed elbow designs
- Swathe attaches at any point on sling
- Available in pediatric and extra long versions
- Universal right/left

**Indications:** S/P glenohumeral acromioplasty; Surgical repairs of the rotator cuff; Anterior glenohumeral capsular reconstruction; Anterior dislocation of the G-H joint

Size	Swathe Size	Pouch Size	Product #
Universal	5" x 54"	13" x 19"	9006-00
Universal Enclosed Elbow	5" x 54"	17" x 19"	9014-00
Universal w/ X-Long Swathe	5" x 58"	13" x 19"	9006-50
Universal Pediatric	4½" x 44"	9" x 16½"	9007-00
Swathe Only	5" x 54"	—	9006-10



## Foam Sling & Swathe

- Lightweight perforated foam material provides patient comfort
- Adjustable shoulder straps attach to hook receivable sling to ensure accurate fit
- Swathe attaches at any point on sling
- Available in universal and sized options
- Universal right/left

**Indications:** S/P glenohumeral acromioplasty; Surgical repairs of the rotator cuff; Anterior glenohumeral capsular reconstruction; Anterior dislocation of the G-H joint

Size	Swathe Size	Pouch Size	Product #
Universal	5" x 54"	13" x 19"	9013-00
Small	5" x 46"	10" x 16½"	9013-02
Medium	5" x 48"	11" x 18½"	9013-03
Large	5" x 51"	12" x 19"	9013-04



## Universal Foam Sling & Swathe

- Constructed of perforated foam material for patient comfort
- Adjustable shoulder straps attach to hook receivable sling to ensure accurate fit
- Swathe attaches at any point on sling
- Extra-long swathe size available
- Universal right/left

**Indications:** S/P glenohumeral acromioplasty; Surgical repairs of the rotator cuff; Anterior glenohumeral capsular reconstruction; Anterior dislocation of the G-H joint

Size	Swathe Size	Pouch Size	Product #
Universal	5" x 54"	13" x 19"	4442-00
Universal w/ X-Long Swathe	5" x 58"	13" x 19"	4442-10



# SHOULDER IMMOBILIZERS



## Foam Sling & Swathe

- Constructed of perforated foam material
- Adjustable shoulder straps attach to hook receivable sling to ensure accurate fit
- Swathe attaches at any point on sling
- Available in universal and sized options
- Universal right/left

**Indications:** S/P glenohumeral acromioplasty; Surgical repairs of the rotator cuff; Anterior glenohumeral capsular reconstruction; Anterior dislocation of the G-H joint

Size	Swathe Size	Pouch Size	Product #
Universal	5" x 54"	13" x 19"	TX4442-00
Small	5" x 46"	10" x 16½"	TX4442-02
Medium	5" x 48"	11" x 18½"	TX4442-03
Large	5" x 51"	12" x 19"	TX4442-04



## Velpeau Shoulder Immobilizer

- Constructed of heavy duty elastic material
- Durable foam padded wrist and humeral cuffs
- Contoured design with hook and loop closure for patient comfort
- Universal right/left

**Indications:** S/P glenohumeral acromioplasty; Surgical repair of the rotator cuff; Anterior glenohumeral capsular reconstruction; Anterior dislocation of the G-H joint

Size	Belt Size	Product #
X-Small	6" x 29"	9008-00
Small	6" x 32"	9008-01
Medium	6" x 35"	9008-02
Large	6" x 38"	9008-03
X-Large	6" x 41"	9008-04
XX-Large	6" x 44"	9008-05
XXX-Large	6" x 47"	9008-06

Best Choice



## Elastic Shoulder Immobilizer

- Constructed of lightweight elastic material
- Foam wrist and humeral cuffs
- Low profile design is contoured for patient comfort
- Hook and loop closure
- Universal right/left

**Indications:** S/P glenohumeral acromioplasty; Surgical repair of the rotator cuff; Anterior glenohumeral capsular reconstruction; Anterior dislocation of the G-H joint

Size	Belt Size	Product #
X-Small	6" x 29"	1000159
Small	6" x 32"	1010159
Medium	6" x 35"	1020159
Large	6" x 38"	1030159
X-Large	6" x 41"	1040159
XX-Large	6" x 44"	1050159

# SHOULDER IMMOBILIZERS

## Shoulder Immobilizer w/ Canvas Swathe

- Sling with attached swathe constructed of navy canvas material
- Shoulder strap includes o-ring/hook and loop closure
- Shoulder pad for added comfort
- Thumb loop helps prevent wrist drop
- Universal right/left

**Indications:** S/P glenohumeral acromioplasty; Surgical repairs of the rotator cuff; Anterior glenohumeral capsular reconstruction; Anterior dislocation of the G-H joint

Size	Pouch Size	Product #
X-Small	6" x 12½"	9001-01
Small	7" x 13½"	9001-02
Medium	8" x 16½"	9001-03
Large	9" x 19½"	9001-04
X-Large	9" x 22½"	9001-05



## Canvas Shoulder Immobilizer w/ Waist Strap

- Constructed of navy canvas sling with web shoulder and waist straps
- Shoulder strap includes o-ring/hook and loop closure
- Adjustable buckle on waist strap
- Thumb loop helps prevent wrist drop
- Universal right/left

**Indications:** S/P glenohumeral acromioplasty; Surgical repairs of the rotator cuff; Anterior glenohumeral capsular reconstruction; Anterior dislocation of the G-H joint

Size	Pouch Size	Product #
X-Small	6" x 12½"	9002-01
Small	7" x 13½"	9002-02
Medium	8" x 16½"	9002-03
Large	9" x 19½"	9002-04
X-Large	9" x 22½"	9002-05



## Muslin Shoulder Immobilizer w/ Waist Strap

- Constructed of lightweight muslin sling with web shoulder and waist straps
- Adjustable buckle closures
- Closed pouch design
- Thumb loop helps prevent wrist drop
- Universal right/left

**Indications:** S/P glenohumeral acromioplasty; Surgical repairs of the rotator cuff; Anterior glenohumeral capsular reconstruction; Anterior dislocation of the G-H joint

Size	Pouch Size	Product #
X-Small	6" x 12½"	9003-01
Small	7" x 13½"	9003-02
Medium	8" x 16½"	9003-03
Large	9" x 19½"	9003-04
X-Large	9" x 22½"	9003-05





# SHOULDER IMMOBILIZERS



## Shoulder Immobilizer w/ Foam Straps

- Constructed of canvas or lightweight tietex material with foam shoulder and waist straps
- Adjustable straps feature o-ring/hook and loop closure
- Thumb loop helps prevent wrist drop
- Universal right/left

**Indications:** S/P glenohumeral acromioplasty; Surgical repairs of the rotator cuff; Anterior glenohumeral capsular reconstruction; Anterior dislocation of the G-H joint

Size	Pouch Size	Canvas	Tietex
X-Small	6" x 12½"	9023-00	TX9023-00
Small	7" x 13½"	9023-01	TX9023-01
Medium	8" x 16½"	9023-02	TX9023-02
Large	9" x 19½"	9023-03	TX9023-03
X-Large	9" x 22½"	9023-04	TX9023-04



## Universal Foam Shoulder Immobilizer

- Contoured perforated foam construction ensures patient comfort
- Wrist cuff, humeral cuff and shoulder strap are adjustable and can be placed anywhere on belt to achieve desired fit
- Belt available in light or dark foam\*
- Provides shoulder and upper arm immobilization following dislocation or post-surgical application
- Universal right/left

**Indications:** S/P glenohumeral acromioplasty; Surgical repairs of the rotator cuff; Anterior glenohumeral capsular reconstruction; Anterior dislocation of the G-H joint

Size	Belt Size	Wrist Cuff Size	Product #
Universal	7½" x 48"	5" x 12"	9009-00
Universal	8" x 54"	3" x 9"	9009-50*



## Universal Sling Style Shoulder Immobilizer

- Perforated foam construction ensures patient comfort
- Sling style wrist cuff for optimal comfort and support
- Wrist cuff, humeral cuff and padded shoulder strap are adjustable and can be placed anywhere on belt to achieve desired fit
- Provides shoulder and upper arm immobilization following dislocation or post-surgical application
- Universal right/left

**Indications:** S/P glenohumeral acromioplasty; Surgical repairs of the rotator cuff; Anterior glenohumeral capsular reconstruction; Anterior dislocation of the G-H joint

Size	Belt Size	Wrist Cuff Size	Product #
Adult	5" x 51"	5½" x 15"	9012-00
Pediatric	4¾" x 34½"	2" x 9¾"	9011-00

# SHOULDER ABDUCTION PILLOWS

## Shoulder P.A.D. II

- Constructed of hook receivable material for infinite adjustability
- Soft breathable liner wicks away moisture
- Foam pillow provides 35° of abduction
- Shoulder and waist straps can be trimmed to desired length
- Adjustable for open or closed pouch design
- Quick-release buckle on waist strap for ease of application and removal
- Adjustable thumb loop to prevent wrist drop
- Arm can be positioned more anteriorly or posteriorly depending on the post-operative requirements
- Exercise ball packaged with product
- Universal right/left

**Indications:** Surgical rotator cuff repairs; Anterior glenohumeral capsular reconstruction; Dislocation/subluxation of G-H joint; S/P glenohumeral acromioplasty; Bankart repairs; Soft tissue repairs; Shoulder sprains and strains

Size	Pouch Size	Product #
Small	7" x 16"	11630005
Medium	7" x 17"	11630006
Large	7" x 18"	11630007
X-Large	7" x 19"	11630008



## Shoulder P.A.D.

- Constructed of lightweight canvas material
- Foam pillow provides 35° of abduction
- Shoulder and waist straps can be trimmed to desired length
- Thumb loop helps prevent wrist drop
- Allows arm to be placed in an optimal neutral position for effective healing
- Exercise ball packaged with product
- Universal right/left

**Indications:** Surgical rotator cuff repairs; Anterior glenohumeral capsular reconstruction; Dislocation/subluxation of G-H joint; S/P glenohumeral acromioplasty; Bankart repairs; Soft tissue repairs; Shoulder sprains and strains

Size	Pouch Size	Product #
Small	5½" x 14"	11640005
Medium	6½" x 17"	11640006
Large	7" x 19"	11640007
X-Large	7½" x 20½"	11640008





# SHOULDER ABDUCTION PILLOWS



## Shoulder Abduction Pillow

- Constructed of hook receivable foam
- Foam pillow provides 55° of abduction
- Adjustable wrist and humeral cuffs for proper positioning
- Padded strap for patient comfort
- Designed for abduction and immobilization of the shoulder
- Universal right/left

**Indications:** Acute shoulder instability; Surgical rotator cuff repair; Anterior glenohumeral capsular reconstruction; Acute rotator cuff tendinopathy; Dislocation/subluxation of G-H joint; Acromioclavicular joint dislocations; S/P glenohumeral acromioplasty; Bankart repairs; Soft tissue repairs; Shoulder sprains and strains

Size	Pillow Dimensions	Product #
Small	25 $\frac{1}{4}$ " x 23 $\frac{3}{4}$ "	11650005
Medium	29 $\frac{1}{2}$ " x 26 $\frac{3}{4}$ "	11650006
Large	30 $\frac{1}{2}$ " x 36"	11650007

# ARM SLINGS/ELEVATORS

## Specialty Arm Sling

- Durable navy canvas construction with deep pocket design
- Infant and child sizes feature cartoon print
- Choice of foam or vinyl shoulder pad for patient comfort
- Thumb loop helps prevent wrist drop
- Adjustable slide buckle on strap for best fit
- Hook and loop closure with o-rings

**Indications:** Cast support; Injuries to the arm, wrist or hand; Acromioclavicular sprains; Anterior glenohumeral dislocations

Size	Pouch Size	Strap Length	Foam Pad	Vinyl Pad
Infant	5" x 9"	1" x 21"	8004-01	8004-11
Child	6" x 11"	1" x 27"	8004-02	8004-12
Small	7" x 13"	2" x 34"	8004-03	8004-13
Medium	8" x 17"	2" x 37"	8004-04	8004-14
Large	9" x 20"	2" x 37"	8004-05	8004-15
X-Large	9" x 23"	2" x 42"	8004-06	8004-16



## Premium Arm Sling

- Deep pocket pouch design constructed of black cotton material
- Infant and child sizes feature cartoon print
- Foam shoulder pad provides additional patient comfort
- Thumb loop helps prevent wrist drop
- Adjustable slide buckle on strap for best fit
- Hook and loop closure with o-rings

**Indications:** Cast support; Injuries to the arm, wrist or hand; Acromioclavicular sprains; Anterior glenohumeral dislocations

Size	Pouch Size	Strap Length	Product #
Infant	5" x 9"	1" x 21"	11690001
Child	6" x 11"	1" x 27"	11690002
Small	7" x 13"	2" x 34"	11690005
Medium	8" x 17"	2" x 37"	11690006
Large	9" x 20"	2" x 37"	11690007
X-Large	9" x 23"	2" x 42"	11690008



## Universal Arm Sling

- Constructed of heavy-duty canvas material
- Features easy fold-back universal design
- Foam shoulder pad provides additional comfort
- Hook and loop closure with o-rings

**Indications:** Cast support; Injuries to the arm, wrist or hand; Acromioclavicular sprains; Anterior glenohumeral dislocations

Size	Pouch Size	Strap Length	Product #
Adult	8" x 19"	2" x 37"	8020-00
Child	5" x 9"	1" x 31"	8020-01



# ARM SLINGS/ELEVATORS



## Envelope Arm Sling

- Constructed of lightweight navy canvas material
- Deep pocket pouch design
- Foam shoulder pad for added comfort
- Thumb loop helps prevent wrist drop
- Buckle closure

**Indications:** Cast support; Injuries to the arm, wrist or hand; Acromioclavicular sprains; Anterior glenohumeral dislocations

Size	Pouch Size	Strap Length	Product #
Child	6" x 11"	1" x 27"	8001-01
Small	7" x 13"	1" x 34"	8001-02
Medium	8" x 17"	1" x 37"	8001-03
Large	9" x 20"	1" x 42"	8001-04
X-Large	9" x 23"	1" x 45"	8001-05



## Light Blue Arm Sling

- Constructed of light blue cotton material with deep pocket design
- Foam shoulder pad provides additional comfort
- Thumb loop helps prevent wrist drop
- Hook and loop closure

**Indications:** Cast support; Injuries to the arm, wrist or hand; Acromioclavicular sprains; Anterior glenohumeral dislocations

Size	Pouch Size	Strap Length	Product #
Infant	5" x 9"	1" x 21"	10101
Child	6" x 11"	1" x 27"	10102
Small	7" x 13"	2" x 34"	10103
Medium	8" x 17"	2" x 37"	10104
Large	9" x 20"	2" x 37"	10105
X-Large	9" x 23"	2" x 42"	10106



## Narrow Pouch Arm Sling

- Lightweight cotton material with narrow, closed pouch design
- Buckle closure
- Available in each or 12/Case quantity

**Indications:** Cast support; Injuries to the arm, wrist or hand; Acromioclavicular sprains; Anterior glenohumeral dislocations

Size	Pouch Size	Strap Length	Each	Case
Child	5" x 10"	1" x 28"	—	8017-00
Small	5" x 13"	1" x 34"	8017-51	8017-01
Medium	6" x 15"	1" x 38"	8017-52	8017-02
Large	6" x 18"	1" x 42"	8017-53	8017-03
X-Large	7" x 20"	1" x 46"	8017-54	8017-04

# ARM SLINGS/ELEVATORS

## Arm Sling w/ Foam Strap

- Constructed of canvas or lightweight tietex material
- Deep pocket, closed pouch design
- Foam shoulder strap provides maximum comfort
- Hook closure fastens anywhere on foam strap for superior adjustment
- Available in sized or easy fold back universal style

**Indications:** Cast support; Injuries to the arm, wrist or hand; Acromioclavicular sprains; Anterior glenohumeral dislocations



Size	Pouch Size	Strap Length	Canvas	Tietex
Universal	8" x 19"	2" x 37"	—	TX8020-00
Small	7" x 13"	2" x 37"	8023-01	TX8023-21
Medium	8" x 17"	2" x 37"	8023-02	TX8023-22
Large	9" x 20"	2" x 37"	8023-03	TX8023-23
X-Large	9" x 22"	2" x 40"	8023-04	TX8023-24
XX-Large	10" x 23"	2" x 40"	8023-05	—

## Deep Pocket Arm Sling

- Deep pocket, open pouch design constructed of breathable tietex material
- Foam shoulder pad provides additional comfort
- Thumb loop helps prevent wrist drop
- Hook and loop closure with o-rings and slide buckle adjustment
- Available in each or 100/Case quantity

**Indications:** Cast support; Injuries to the arm, wrist or hand; Acromioclavicular sprains; Anterior glenohumeral dislocations



Size	Pouch Size	Strap Length	Each	Case
Infant	5" x 9"	1" x 21"	TX9902-01	—
Child	6" x 11"	1" x 27"	TX9902-02	TX9902-02/100
Small	7" x 13"	2" x 37"	TX9902-03	TX9902-03/100
Medium	8" x 17"	2" x 37"	TX9902-04	TX9902-04/100
Large	9" x 20"	2" x 37"	TX9902-05	TX9902-05/100
X-Large	9" x 22"	2" x 42"	TX9902-06	TX9902-06/100

## Strap Arm Sling

- 2" strap construction
- Vinyl neck pad for patient comfort
- Adjustable buckles

**Indications:** Cast support; Injuries to the arm, wrist or hand

Size	Strap Length	Product #
Child	1" x 52"	8006-01
Small	2" x 62"	8006-02
Medium	2" x 68"	8006-03
Large	2" x 75"	8006-04
X-Large	2" x 86"	8006-05



# ARM SLINGS/ELEVATORS



## CVA Sling

- Constructed of durable canvas material
- Tubular foam neck pad for patient comfort
- Strap features adjustable slide buckle closure

**Indications:** Stroke or paralysis affecting arm or shoulder

Size	Pouch Size	Strap Length	Product #
Small/Medium	5½" x 7½"	2" x 30"	6157D02
Medium/Large	5½" x 7½"	2" x 38½"	6157D03



## Collar and Cuff

- Foam-lined leatherette cuff is perforated for breathability
- Double d-ring hook closure provides open cuff access
- Vinyl neck pad for patient comfort
- Slide buckle adjustment

**Indications:** Cast support; Wrist/forearm elevation

Size	Cuff Size	Strap Length	Product #
Small/Medium	5" x 11"	2" x 44"	8013-01
Medium/Large	8" x 11"	2" x 52"	8013-02



## Triangle Bandage Sling

- Constructed of muslin material
- Includes safety pin for fastening
- Available in each or 12/Case quantity

**Indications:** Cast support; Injuries to the arm, wrist or hand

Size	Each	Case
Universal	8875-00	8875-100

# ARM SLINGS/ELEVATORS

## Perforated Arm Elevator/Sling

- Constructed of hook receivable perforated foam
- Converts easily to sling with shoulder pad included
- Can be trimmed for best patient fit
- Hook and loop closures for easy application

**Indications:** Post upper extremity surgery; Upper extremity edema

<b>Size</b>	<b>Pouch Size</b>	<b>Product #</b>
Universal	20" x 20½"	8788-50



## Convuluted Foam Arm Elevator/Sling

- Constructed of convuluted foam inner with hook receivable outer lining
- Converts easily to sling with shoulder pad included
- Slide buckle for adjustment of shoulder strap
- Hook and loop closures for easy application

**Indications:** Post upper extremity surgery; Upper extremity edema

<b>Size</b>	<b>Pouch Size</b>	<b>Product #</b>
Universal	10" x 16"	8788-03



## Arm Elevator/Sling w/ Ice Pockets

- Heavy duty canvas construction
- Includes two ice bags for cold therapy
- Converts easily to sling

**Indications:** Post upper extremity surgery; Upper extremity edema

<b>Size</b>	<b>Pouch Size</b>	<b>Product #</b>
Small	10" x 14"	8011-01
Medium	10" x 16"	8011-02
Large	12" x 17"	8011-03





# ELBOW



## Premium Tennis Elbow Support

- Made of lightweight, ventilated material
- Support is low profile increasing patient compliance
- Inner pad can be strategically placed for maximum pain relief
- D-ring closure for easy application
- Universal one size fits all

**Indications:** Medial and lateral epicondylitis

Size	Measurement	Product #
Universal	2" x 16"	12870000



## PneuGel™ Tennis Elbow Support

- Constructed of soft foam material for patient comfort
- Concentration of pressure allows compression to be applied where needed while not constricting uninvolved muscles or blood flow
- Air/gel removable bladder can be placed in freezer to provide cold therapy
- Inflation bulb included to adjust bladder
- Universal one size fits all

**Indications:** Medial and lateral epicondylitis

Size	Measurement	Product #
Universal	2" x 20"	12840000



## Universal Tennis Elbow Support

- Perforated foam material is breathable for patient comfort
- Counterforce strap with loop lock contact closure for adjustable compression
- Universal one size fits all

**Indications:** Medial and lateral epicondylitis

Size	Measurement	Product #
Universal	3" x 13"	6005-00



## Sized Tennis Elbow Support

- Constructed of surgical webbing with foam lining
- Counterforce strap with loop lock contact closure for adjustable compression
- Sized to ensure a precise fit

**Indications:** Medial and lateral epicondylitis

Size	Measurement	Product #
Small	3" x 9"	6004-01
Medium	3" x 10"	6004-02
Large	3" x 11"	6004-03
X-Large	3" x 12"	6004-04



## Universal Neoprene Tennis Elbow Support

- 1/8" neoprene construction provides compression, warmth and support
- Adjustable felt buttress can be placed just below the epicondyle for pain relief
- Loop lock contact closure for adjustable compression and application
- Universal one size fits all

**Indications:** Medial and lateral epicondylitis

Size	Measurement	Product #
Universal	3" x 16"	NE7730-70



## Neoprene Tennis Elbow Support w/ Pad

- 1/8" neoprene construction provides compression, warmth and support
- Incorporates 1/4" pad to alleviate stress over the lateral or medial epicondyle
- Elastic counterforce strap and loop lock contact closure allow for maximum compression
- Sized to ensure a precise fit

**Indications:** Medial and lateral epicondylitis

Size	Forearm Circumference	Product #
X-Small	8" - 9"	NE7731-71
Small	9" - 10"	NE7731-72
Medium	10" - 11"	NE7731-73
Large	11" - 12"	NE7731-74
X-Large	12" - 14"	NE7731-75
XX-Large	14" - 16"	NE7731-76



# ELBOW



## Neoprene Elbow Sleeve

- 1/8" neoprene construction provides compression, warmth and support
- Anatomically contours to the elbow
- Available with or without elbow strap

**Indications:** Tendonitis; Osteoarthritis; Arthritis; Olecranon bursitis; (Medial and lateral epicondylitis w/ strap)

Size	Forearm Circumference	w/o Strap	w/ Strap
X-Small	8" - 9"	NE7727-71	NE7728-71
Small	9" - 10"	NE7727-72	NE7728-72
Medium	10" - 11"	NE7727-73	NE7728-73
Large	11" - 12"	NE7727-74	NE7728-74
X-Large	12" - 14"	NE7727-75	NE7728-75
XX-Large	14" - 16"	NE7727-76	NE7728-76



## Elastic Elbow Support

- Pull on cotton elastic design
- Provides snug fit to achieve maximum comfort and compression

**Indications:** Minor elbow sprains and strains; Edema; Osteoarthritis; Olecranon bursitis

Size	Forearm Circumference	Product #
Small	7" - 10"	6003-01
Medium	10" - 13"	6003-02
Large	13" - 15"	6003-03
X-Large	15" - 17"	6003-04


## Transition Elbow Brace

- Breathable foam softgoods are trimmable for a customized fit
- Hinge can be easily adjusted in 10-degree increments without tools (Flexion: 120° to 0° and Extension: 90° to 0°)
- Loop lock strapping is trimmable to allow for superior fit
- Length can be adjusted from 19" to 16"
- Universal right/left

**Indications:** Post-operative rehabilitation

Size	Measurement*	Product #
Small	6" - 8"	EB5000-01
Medium	8" - 10"	EB5000-02
Large	10" +	EB5000-03

\*Measure from axis of elbow to axis of wrist



## Slimline Elbow Brace

- Hinge features range of motion settings in 10-degree increments
- Contoured softgoods at biceps and forearm to provide precise fit
- Universal one size fits all

**Indications:** Post-operative rehabilitation

Size	Product #
Universal	EB2040-00



## Wrist & Elbow Splint

- Constructed of durable leatherette material
- Provides immobilization of elbow at 90°
- Aluminum frame provides protection to the elbow

**Indications:** Olecranon fractures; Radial head fractures; Hairline fracture of the radius, humerus and ulna; Ligamentous reconstruction procedures of the elbow

Size	Biceps Circumference	Forearm Circumference	Product #
Small	11½"	9¾"	6001-01
Medium	12½"	10"	6001-02
Large	13½"	11"	6001-03
X-Large	14½"	12"	6001-04



# WRIST & WRIST/FOREARM



## Suede Leatherette Wrist & Wrist/Forearm Splint

- Suede leatherette material reduces migration
- Unique wedge contour design ensures an intimate fit increasing patient comfort
- Extra-wide loop lock strap at wrist for compression and support
- Removable/malleable palmar and dorsal stays
- Does not extend past the distal palmar crease allowing full finger range of motion

**Indications:** Carpal tunnel syndrome; Sprains and strains of the wrist; Post-operative immobilization and rehabilitation; Stable fractures of the distal radius and carpal joints

Best Choice

Size	Wrist			Wrist/Forearm		
	Length	Right	Left	Length	Right	Left
Pediatric	4"	12981202	12981102	6"	12991202	12991102
Child	8"	—	—	9"	12992202	12992102
Small	8"	12982205	12982105	11½"	12993205	12993105
Medium	8"	12982206	12982106	11½"	12993206	12993106
Large	8"	12982207	12982107	11½"	12993207	12993107
X-Large	8"	12982208	12982108	11½"	12993208	12993108



## Leatherette Wrist & Wrist/Forearm Splint

- Constructed of durable leatherette material with foam tricot lining
- Unique wedge contour design ensures an intimate fit increasing patient comfort
- Extra-wide loop lock strap at wrist for compression and support
- Removable/malleable palmar and dorsal stays
- Does not extend past the distal palmar crease allowing full finger range of motion

**Indications:** Carpal tunnel syndrome; Sprains and strains of the wrist; Post-operative immobilization and rehabilitation; Stable fractures of the distal radius and carpal joints

Size	Wrist			Wrist/Forearm		
	Length	Right	Left	Length	Right	Left
Pediatric	4"	8745-01	8745-06	6"	8824-01	8824-06
Child	8"	—	—	9"	8824-00	8824-50
Small	8"	8745-02	8745-07	11"	8824-02	8824-07
Medium	8"	8745-03	8745-08	11"	8824-03	8824-08
Large	8"	8745-04	8745-09	11"	8824-04	8824-09
X-Large	8"	8745-05	8745-10	11"	8824-05	8824-10



## Canvas Wrist & Wrist/Forearm Splint

- Constructed of sturdy blue canvas material with foam tricot lining
- Unique wedge contour design ensures an intimate fit increasing patient comfort
- Extra-wide loop lock strap at wrist for compression and support
- Removable/malleable palmar and dorsal stays
- Does not extend past the distal palmar crease allowing full finger range of motion

**Indications:** Carpal tunnel syndrome; Sprains and strains of the wrist; Post-operative immobilization and rehabilitation; Stable fractures of the distal radius and carpal joints

Size	Wrist			Wrist/Forearm		
	Length	Right	Left	Length	Right	Left
Pediatric	4"	1010115	1060115	6"	1010025	1060025
Small	8"	1020115	1070115	11"	1020025	1070025
Medium	8"	1030115	1080115	11"	1030025	1080025
Large	8"	1040115	1090115	11"	1040025	1090025
X-Large	8"	1050115	1100115	11"	1050025	1100025

# WRIST & WRIST/FOREARM

## Black Foam Wrist & Wrist/Forearm Splint

- Comfortable black foam material with tricot lining
- Extra-wide loop lock strap at wrist for compression and support
- Removable/malleable palmar and dorsal stays
- Does not extend past the distal palmar crease allowing full finger range of motion

**Indications:** Carpal tunnel syndrome; Sprains and strains of the wrist; Post-operative immobilization and rehabilitation; Stable fractures of the distal radius and carpal joints



Size	Length	Wrist		Length	Wrist/Forearm	
		Right	Left		Right	Left
Pediatric	4"	BF5011-01	BF5011-06	6"	BF5002-01	BF5002-06
Child	—	—	—	9"	BF5030-01	BF5030-02
Small	8"	BF5011-02	BF5011-07	11"	BF5002-02	BF5002-07
Medium	8"	BF5011-03	BF5011-08	11"	BF5002-03	BF5002-08
Large	8"	BF5011-04	BF5011-09	11"	BF5002-04	BF5002-09
X-Large	8"	BF5011-05	BF5011-10	11"	BF5002-05	BF5002-10

## Tietex Wrist & Wrist/Forearm Splint

- Tietex material construction is lightweight and comfortable for patient
- Unique wedge contour design ensures an intimate fit increasing patient comfort
- Extra-wide loop lock strap at wrist for compression and support
- Removable/malleable palmar and dorsal stays
- Does not extend past the distal palmar crease allowing full finger range of motion

**Indications:** Carpal tunnel syndrome; Sprains and strains of the wrist; Post-operative immobilization and rehabilitation; Stable fractures of the distal radius and carpal joints



Size	Length	Wrist		Length	Wrist/Forearm	
		Right	Left		Right	Left
Pediatric	4"	TX9903-01	TX9903-06	6"	TX9904-01	TX9904-06
Small	8"	TX9903-02	TX9903-07	11"	TX9904-02	TX9904-07
Medium	8"	TX9903-03	TX9903-08	11"	TX9904-03	TX9904-08
Large	8"	TX9903-04	TX9903-09	11"	TX9904-04	TX9904-09
X-Large	8"	TX9903-05	TX9903-10	11"	TX9904-05	TX9904-10

## Universal Leatherette Wrist/Forearm Splint

- Constructed of durable white leatherette material with foam tricot lining
- Malleable palmar and dorsal stays provide anatomical fit and support
- Loop lock closures add compression and support
- Universal right/left and size to fit majority of patients

**Indications:** Carpal tunnel syndrome; Sprains and strains of the wrist; Post-operative immobilization and rehabilitation; Stable fractures of the distal radius and carpal joints



Size	Length	Product #
Universal	11"	5007-00

# WRIST & WRIST/FOREARM



## Universal Intera® Wrist & Wrist/Forearm Splint

- Constructed of a unique odor resistant Intera® material with moisture transport properties to keep patient's skin cool and dry
- Malleable palmar and dorsal stays provide anatomical fit and proper support
- Movable dorsal stay makes product universal in size
- D-ring closures for easy application and removal

**Indications:** Sprains and strains of the wrist; Tendonitis; Tenosynovitis of the wrist and distal forearm

Size	Length	Right	Left
Universal Wrist	8"	5066-92	5066-93
Universal Wrist/Forearm	10"	5066-90	5066-91



## Universal Foam Wrist & Wrist/Forearm Splint

- Perforated foam material allows skin to breathe
- Malleable palmar and dorsal stays provide anatomical fit and proper support
- Moveable dorsal stay makes product universal in size
- D-ring closures for easy application and removal
- Available without binding around edges for trimming

**Indications:** Sprains and strains of the wrist; Tendonitis; Tenosynovitis of the wrist and distal forearm

Size	Length	w/ Binding Right	w/ Binding Left
Universal	8"	5066-82	5066-83
Universal	10"	5066-80	5066-81

Size	Length	w/o Binding Right	w/o Binding Left
Universal	8½"	5005-03	5005-04
Universal	11½"	5005-01	5005-02



## Universal Black Foam Wrist & Wrist/Forearm Splint

- Comfortable black foam material with tricot lining
- Malleable palmar and dorsal stays provide anatomical fit and proper support
- Moveable dorsal stay makes product universal in size
- D-ring closures for easy application and removal

**Indications:** Sprains and strains of the wrist; Tendonitis; Tenosynovitis of the wrist and distal forearm

Size	Length	Right	Left
Universal	8"	BF5066-82	BF5066-83
Universal	10"	BF5066-80	BF5066-81

# WRIST & WRIST/FOREARM

## Universal Cutaway Wrist & Wrist/Forearm Splint

- Perforated foam construction allows skin to breathe
- Malleable palmar and dorsal stays provide anatomical fit and proper support
- Loop lock closure at wrist provides compression and support
- Foam can be trimmed for desired fit

**Indications:** Sprains and strains of the wrist; Tendonitis; Tenosynovitis of the wrist and distal forearm

Size	Length	Right	Left
Universal	8"	5040-03	5040-04
Universal	11"	5040-01	5040-02



## Universal Two-Panel Wrist & Wrist/Forearm Splint

- Blue canvas construction with foam tricot lining
- Malleable palmar and dorsal stays limit dorsiflexion motion
- Loop lock closure at wrist provides optimal immobilization and compression
- Universal right/left and size to fit majority of patients

**Indications:** Sprains or strains of the wrist; Tendonitis; Tenosynovitis of the wrist and distal forearm

Size	Length	Product #
Universal	8"	1000275
Universal	11"	1000285



## Black Wrist Splint

- Constructed of comfortable black foam material
- Malleable palmar stay provides anatomical fit and rigid support
- Low profile and lightweight design
- D-ring closures for easy application and removal
- Does not extend past distal palmar crease allowing full finger range of motion

**Indications:** Mild sprains and strains of the wrist; Tendonitis, Carpal tunnel syndrome; Overuse injuries; Tenosynovitis

Size	8" Right	8" Left	10" Right	10" Left
Pediatric	5073-11	5073-12	—	—
X-Small	5073-01	5073-06	5074-01	5074-06
Small	5073-02	5073-07	5074-02	5074-07
Medium	5073-03	5073-08	5074-03	5074-08
Large	5073-04	5073-09	5074-04	5074-09
X-Large	5073-05	5073-10	5074-05	5074-10





# WRIST & WRIST/FOREARM



## Premium Wrist Splint

- Breathable material with four way stretch provides patient comfort
- Removable/malleable palmar stay
- Low profile and lightweight design
- D-ring closures for easy application and removal
- Does not extend past distal palmar crease allowing full finger range of motion

**Indications:** Mild sprains and strains of the wrist; Tendonitis, Carpal tunnel syndrome; Overuse injuries; Tenosynovitis

Size	Length	Right	Left
X-Small	8"	351XSR	351XSL
Small	8"	351SR	351SL
Medium	8"	351MR	351ML
Large	8"	351LR	351LL
X-Large	8"	351XLR	351XLL
XX-Large	8"	351XXLR	351XXLL

Best Choice



## Lace Up Suede Leatherette Wrist Splint

- Soft perforated suede leatherette material with nylon tricot lining
- Removable/malleable palmar stay
- Single pull lace closure for compression and easy adjustment
- Does not extend past distal palmar crease allowing full finger range of motion

**Indications:** Mild sprains and strains of the wrist; Tendonitis; Carpal tunnel syndrome; Overuse injuries; Tenosynovitis

Size	Length	Right	Left
X-Small	8"	5023-01	5023-06
Small	8"	5023-02	5023-07
Medium	8"	5023-03	5023-08
Large	8"	5023-04	5023-09
X-Large	8"	5023-05	5023-10



## Lace Up Canvas Wrist Splint

- Constructed of durable canvas material
- Removable/malleable palmar spoon stay
- Laces are sewn onto hook and loop closure system to allow for compression adjustment and easy application
- Does not extend past distal palmar crease allowing full finger range of motion

**Indications:** Mild sprains and strains of the wrist; Tendonitis; Carpal tunnel syndrome; Overuse injuries; Tenosynovitis

Size	6" Right	6" Left	8" Right	8" Left
XX-Small	5017-00	5017-50	5022-00	5022-50
X-Small	5017-01	5017-05	5022-01	5022-05
Small	5017-02	5017-06	5022-02	5022-06
Medium	5017-03	5017-07	5022-03	5022-07
Large	5017-04	5017-08	5022-04	5022-08
X-Large	5017-09	5017-10	5022-09	5022-10

# WRIST & WRIST/FOREARM

## Elastic Wrist Splint

- Cotton elastic construction provides compression
- Removable/malleable palmar stay
- Hook and loop closure for easy application and removal
- Does not extend past distal palmar crease allowing full finger range of motion

**Indications:** Mild sprains and strains of the wrist; Tendonitis, Carpal tunnel syndrome; Overuse injuries; Tenosynovitis

Size	6" Right	6" Left	8" Right	8" Left
Pediatric	5015-01	5015-06	5016-01	5016-06
Small	5015-02	5015-07	5016-02	5016-07
Medium	5015-03	5015-08	5016-03	5016-08
Large	5015-04	5015-09	5016-04	5016-09
X-Large	5015-05	5015-10	5016-05	5016-10
XX-Large	—	—	5016-15	5016-20



## Neoprene Wrist Support

- 1/8" neoprene material provides compression, warmth and support
- Removable palmar stay can be molded for a custom fit
- Hook and loop closure for easy application
- Low profile at the palmar crease for easy thumb-finger opposition

**Indications:** Mild sprains and strains of the wrist; Tendonitis; Carpal tunnel syndrome; Overuse injuries; Tenosynovitis

Size	Length	Right	Left
Small	8"	NE7739-71	NE7739-76
Medium	8"	NE7739-72	NE7739-77
Large	8"	NE7739-73	NE7739-78
X-Large	8"	NE7739-74	NE7739-79



## Universal Ambidextrous Wristlet

- 1/8" neoprene material provides compression, warmth and support
- Wraparound design for universal sizing
- Hook and loop closure for easy application
- Universal one size fits all

**Indications:** Mild sprains and strains of the wrist; Tendonitis; Carpal tunnel syndrome; Overuse injuries; Tenosynovitis

Size	Product #
Universal	NE7740-70



# WRIST & WRIST/FOREARM



## Universal Neoprene Wristlet

- Neoprene construction provides compression, warmth and support
- Lightweight and low profile design
- Hook and loop contact closure for a secure fit
- Universal one size fits all

**Indications:** Wrist sprains and strains

**Size**  
Universal

**Product #**  
5019-00



## Elastic Wrist Support

- Cotton elastic material with thumb opening for precise fit
- Provides snug fit for maximum comfort and compression
- Universal right/left

**Indications:** Wrist sprains and strains

**Size**  
Small  
Medium  
Large  
X-Large

**Product #**  
5019-01  
5019-02  
5019-03  
5019-04

## Premium Wrist & Thumb Splint

- Breathable material with four way stretch provides patient comfort
- Removable/adjustable thumb and palmar stay
- Low profile and lightweight design
- D-ring closures for easy application and removal
- Does not extend past distal palmar crease allowing full finger range of motion

**Indications:** De Quervain's syndrome; Scaphoid injuries; Carpometacarpal instabilities; Metacarpal fractures; Post cast healing; Gamekeeper's thumb; Cumulative trauma injuries

Size	Length	Right	Left
X-Small	8"	350XSR	350XSL
Small	8"	350SR	350SL
Medium	8"	350MR	350ML
Large	8"	350LR	350LL
X-Large	8"	350XLR	350XLL
XX-Large	8"	350XXLR	350XXLL



## Black Foam Wrist & Thumb

- Constructed of comfortable black foam material
- Removable/adjustable thumb and palmar stay
- D-ring closures for easy application and removal
- Does not extend past distal palmar crease allowing full finger range of motion

**Indications:** De Quervain's syndrome; Scaphoid injuries; Carpometacarpal instabilities; Metacarpal fractures; Post cast healing; Gamekeeper's thumb Cumulative trauma injuries

Size	Length	Right	Left
X-Small	8"	A125204	A125104
Small	8"	A125205	A125105
Medium	8"	A125206	A125106
Large	8"	A125207	A125107
X-Large	8"	A125208	A125108
XX-Large	8"	A125209	A124109



## Universal Ambidextrous Wristlet w/ Abducted Thumb

- 1/8" neoprene construction provides warmth, compression and support
- Wrap around design for universal sizing
- Designed to hold thumb in abducted position
- Hook and loop closure for easy application

**Indications:** Arthritis; Thumb sprains; Tendonitis

Size	Product #
Universal	NE7738-70



# THUMB



## Universal Leatherette Thumb Splint

- Constructed of durable perforated leatherette material
- Can be used alone or in combination with wrist and forearm splint for thumb support
- Universal one size fits all

**Indications:** Thumb intrinsic muscle injuries; Post-thumb fracture; Collateral ligament injuries of the thumb

**Size**  
Universal

**Product #**  
5026-00



## Thumb & Wrist/Thumb Spica

- Lightweight splint constructed of perforated foam
- Thumb spica supports thumb while allowing IP joint to remain free
- Wrist & thumb spica supports the wrist and IP joint of the thumb
- Soft inner lining wicks away perspiration for comfort and hygiene
- Stretch loop straps designed for a precise anatomical fit allowing hand and finger function
- May be trimmed for a custom fit

**Indications:** De Quervain's syndrome; CMC arthritis; Post-thumb fractures; Tendon injury/repair



### THUMB

Size	Measurement	Right	Left
Small/Medium	2½" - 3½"	340S/MR	340S/ML
Medium/Large	3½" - 4¼"	340M/LR	340M/LL

### WRIST/THUMB

Size	Measurement	Right	Left
Small/Medium	2½" - 3½"	345S/MR	345S/ML
Medium/Large	3½" - 4¼"	345M/LR	345M/LL

# ALUMINUM SPLINTS

## Baseball Finger Splints

- Lightweight aluminum can be easily adjusted for proper fit
- Foam padding for patient comfort
- Easily moldable to maintain desired position
- X-ray lucent
- Packaged 12/Case

**Indications:** Fractures of the distal phalanx

Size	Length	Product #
Small	4¼"	9109-01
Medium	4¾"	9109-02
Large	5¾"	9109-03



## Toad Finger Splints

- Lightweight aluminum can be easily adjusted for proper fit
- Foam padding for patient comfort
- Malleable arms fold to hold finger in desired alignment
- X-ray lucent
- Packaged 12/Case

**Indications:** Distal and proximal interphalangeal joint positioning

Size	Length	Product #
Small	2¼"	9110-01
Medium	2¾"	9110-02
Large	2⅞"	9110-03



## Fold Over Finger Splints

- Lightweight aluminum can be adjusted for proper fit
- Available with or without foam padding
- X-ray lucent
- Packaged 12/Case

**Indications:** Protects injured nails and fingertips; Baseball or transverse fractures of the finger

Size	Length	w/ Foam	w/o Foam
Small	1"	9111-01	9111-04
Medium	2"	9111-02	9111-05
Large	3"	9111-03	9111-06
X-Large	3½"	9111-07	9111-08



# ALUMINUM SPLINTS



## Four Prong Finger Splints

- Lightweight aluminum can be adjusted for proper fit
- Available with or without foam padding
- X-ray lucent
- Packaged 12/Case

**Indications:** Provides protection and immobilization to injured fingers

Size	Length	w/ Foam	w/o Foam
Small	2"	9112-01	9112-04
Medium	3"	9112-02	9112-05
Large	3½"	9112-03	9112-06



## Gutter Splints

- Lightweight aluminum can be easily adjusted for proper fit
- Available with or without foam padding
- Features slight angle at tip to address a variety of phalangeal injuries
- X-ray lucent
- Packaged 12/Case

**Indications:** Phalangeal and metacarpal fractures; Finger dislocations

Measurement	w/ Foam	w/o Foam
1½" x 7/8"	9113-01	9114-01
3" x 7/8"	9113-02	9114-02
4" x 7/8"	9113-05	—
5½" x 7/8"	9113-03	9114-03
7" x 7/8"	9113-04	9114-04



## Padded Finger Splint Assortment Kit

- Includes an assortment of 57 splints packaged in convenient storage container
- Constructed of lightweight aluminum with foam padding
- X-ray lucent

Description	Qty	Product #
Complete Kit	1 Kit	9100-00
<i>Kit Includes:</i>		
2¾" Toad Splints	12 Each	9110-02
2⅞" Toad Splints	12 Each	9110-03
1½" Gutter Splints	6 Each	9113-01
3" Gutter Splints	6 Each	9113-02
2" Fold Over Finger Splints	6 Each	9111-02
3" Fold Over Finger Splints	6 Each	9111-03
3½" Fold Over Finger Splints	9 Each	9111-07
<i>(All sizes sold individually)</i>		

# ALUMINUM SPLINTS

## Foam Padded Finger Strips

- Constructed of aluminum with foam padding
- Easy to form and can be cut to desired length with scissors
- X-ray lucent
- Packaged 12/Case

**Indications:** Phalangeal fractures

Measurement	Product #
1/2" x 9"	9115-01
1/2" x 18"	9115-02
3/4" x 18"	9115-03
1" x 18"	9115-04
1 1/2" x 18"	9115-05
2" x 18"	9115-06



## Foam Padded Splinting Material

- Constructed of aluminum with foam padding
- Can be trimmed with scissors to achieve desired length for comfortable immobilization
- X-ray lucent
- Packaged 3/Case

**Indications:** Fractures and dislocations

Measurement	Product #
2" x 16"	3302D39
3" x 16"	3302D59
4" x 16"	3302D79



## Finger Splints

- Unique aluminum foam construction can be easily adjusted for proper fit
- X-ray lucent
- Universal right/left
- Packaged 6/Case

**Indications:** Phalangeal fractures

Size	Product #
Small	9125-01
Medium	9125-02
Large	9125-03
X-Large	9125-04





# ALUMINUM SPLINTS



## Burnham Thumb Splint

- Pre-molded aluminum provides proper support for thumb
- Foam padding for patient comfort
- X-ray lucent
- Universal one size fits all

**Indications:** Immobilizes wrist to alleviate muscle-tendon tensions that tend to distort fracture segments

Size	Right	Left
Universal	9130-00	9131-00



## Burnham Finger Splint

- Pre-molded aluminum provides proper support for finger
- Foam padding for patient comfort
- X-ray lucent
- Universal right/left and size

**Indications:** Immobilizes wrist to alleviate muscle-tendon tensions that tend to distort fracture segments

Size	Product #
Universal	9132-00



## Finger Dressing/Finger Cage

- 25 yard roll of tubular bandage
- Perfect for all types of dressings

Size	Width	Product #
Small Fingers & Toes	5/8"	FD1
Fingers & Toes, Infants Hands & Feet	1"	FD2
Small Adult Hands, Arms & Feet, Infant & Child Arms & Feet	1 1/2"	FD3

- Cage constructed of durable aluminum
- Assists in application of tubular finger dressing
- Resuable

Size	Diameter	Length	Product #
Small	1"	6"	9150-01
Large	1 3/4"	6"	9150-02



# ALUMINUM SPLINTS

## Plastalume Splints

- Includes an assortment of 54 splints packaged in convenient storage container
- Constructed of blue plastic attached to white foam padding

Description	Qty	Product #
Complete Kit	1 Kit	9120-00
<i>Kit Includes:</i>		
1¼" Length	6 Each	9119-01
2¼" Length	6 Each	9119-02
3¼" Length	6 Each	9119-03
4¼" Length	6 Each	9119-04
3¼" Length, circular base	6 Each	9119-05
4¼" Length, circular base	6 Each	9119-06
5¼" Length, circular base	6 Each	9119-07
6¼" Length, circular base	6 Each	9119-08
7¼" Length, circular base	6 Each	9119-09
<i>(All sizes sold individually)</i>		



## STAX Finger Splint Kit

- Includes an assortment of 30 splints packaged in convenient storage container
- Protects distal phalanx and extends the DIP joint
- Open tip permits examination of fingertip and nailbed

**Indications:** Fingertip injuries; Nailbed injuries; Mallet finger

Description	Qty	Product #
Complete Kit	1 Kit	9122-00
<i>Kit Includes:</i>		
Size 1	5 Each	9121-01
Size 2	5 Each	9121-02
Size 3	5 Each	9121-03
Size 4	5 Each	9121-04
Size 5	5 Each	9121-05
Size 5½	2 Each	9121-07
Size 6	2 Each	9121-06
Size 7	1 Each	9121-08
<i>(All sizes sold individually)</i>		



## Velcro® Finger Locks

- Soft material for patient comfort
- Enables a stronger finger to assist a weaker one
- Packaged 12/Case

**Indications:** Fractures of phalanges; Ligament injuries; Injuries that require buddy taping

Size	Product #
Universal	9126-00



# ALUMINUM SPLINTS



## Wrist & Forearm Splint

- Pre-molded aluminum design for an anatomical fit and proper support
- Perforated design promotes air circulation
- Available with or without foam padding
- X-ray lucent

**Indications:** Post-operative rehabilitation; Severe sprains and strains of the wrist; Wrist and distal forearm immobilization

Size	w/ Foam Right	w/ Foam Left	w/o Foam Right	w/o Foam Left
Child	9101-01	9101-04	9102-01	9102-04
Youth	9101-02	9101-05	9102-02	9102-05
Adult	9101-03	9101-06	9102-03	9102-06



## Cock-Up Hand Splint

- Pre-molded aluminum design for an anatomical fit and proper support
- Perforated design promotes air circulation
- Available with or without foam padding
- X-ray lucent
- Universal right/left

**Indications:** Post-operative rehabilitation; Severe sprains and strains of the wrist; Wrist and distal forearm immobilization

Size	w/ Foam	w/o Foam
Child	9103-01	9104-01
Youth	9103-02	9104-02
Adult	9103-03	9104-03



## Colles' Splint

- Pre-molded aluminum design for an anatomical fit and proper support
- Perforated design for air circulation
- Available with or without foam padding
- X-ray lucent

**Indications:** Colles' injury; Wrist immobilization

Size	w/ Foam Right	w/ Foam Left	w/o Foam Right	w/o Foam Left
Infant	9105-07	9105-08	—	—
Child	9105-01	9105-04	9106-01	9106-04
Youth	9105-02	9105-05	9106-02	9106-05
Adult	9105-03	9105-06	9106-03	9106-06

# ALUMINUM SPLINTS

## Plastic Colles' Splint

- Pre-molded plastic design for an anatomical fit and proper support
- Perforated design for air circulation
- Foam padding for patient comfort
- X-ray lucent

**Indications:** Colles' injury; Wrist immobilization

Size	Right	Left
Child	9105-24	9105-20
Youth	9105-25	9105-21
Adult	9105-26	9105-22



## Hand Grip Colles' Splint

- Pre-molded aluminum design for an anatomical fit and proper support
- Perforated design promotes air circulation
- Available with or without foam padding
- X-ray lucent

**Indications:** Colles' injury; Wrist immobilization

Size	w/ Foam Right	w/ Foam Left	w/o Foam Right	w/o Foam Left
Child	9107-01	9107-04	9108-01	9108-04
Youth	9107-02	9107-05	9108-02	9108-05
Adult	9107-03	9107-06	9108-03	9108-06



## Short Metacarpal Splint

- Pre-molded aluminum design for an anatomical fit and proper support
- Contoured design to fit in palm of hand
- Perforated design for air circulation
- Foam padding for patient comfort
- X-ray lucent

**Indications:** Phalangeal and metacarpal fractures

Size	Right	Left
Child	9136-01	9136-02
Youth	9134-01	9134-02
Adult	9135-01	9135-02



# ALUMINUM SPLINTS



## Long Metacarpal Splint

- Pre-molded aluminum design for an anatomical fit and proper support
- Perforated design promotes air circulation
- Foam padding for patient comfort
- X-ray lucent

**Indications:** Post-operative rehabilitation; Severe sprains and strains of the wrist; Wrist and distal forearm immobilization

Size	Right	Left
Child	9148-01	9148-04
Youth	9148-02	9148-05
Adult	9148-03	9148-06



## Full Hand Splint

- Pre-molded aluminum design for an anatomical fit and proper support
- Perforated design for air circulation
- Foam padding for patient comfort
- X-ray lucent

**Indications:** Positioning of hand, fingers and thumb

Size	Right	Left
Child	9151-01	9151-04
Youth	9151-02	9151-05
Adult	9151-03	9151-06



## Mason Allen Hand Splint

- Pre-molded aluminum design provides an anatomical fit and proper support
- Available with or without foam padding
- Universal right/left
- X-ray lucent

**Indications:** Immobilization of entire hand and wrist

Size	w/ Foam	w/o Foam
Small	9133-01	9133-04
Universal	9133-03	9133-06

# ALUMINUM SPLINTS

## Elbow Splint

- Pre-molded aluminum design for an anatomical fit and proper support
- Aluminum bar can be bent to adjust angle
- Perforated design for air circulation
- Foam padding for patient comfort
- Universal right/left
- X-ray lucent

**Indications:** Elbow immobilization

Size	Product #
Child	9118-01
Youth	9118-02
Adult	9118-03



## Posterior Leg Splint

- Pre-molded aluminum design for an anatomical fit and proper support
- Open front for easy application
- Short or full leg length options
- Available with or without foam padding
- X-ray lucent

**Indications:** Short: Transverse fractures of the tibia/fibula; Ankle fractures; Lower leg immobilization; Full: Shaft fractures of the femur; Knee immobilization

Short Leg	Length	w/ Foam	w/o Foam
Child	12"	9127-01	9116-01
Youth	14"	9127-02	9116-02
Adult	16"	9127-03	9116-03
Full Leg	Length	w/ Foam	w/o Foam
Child	19½"	9128-01	9117-01
Youth	25"	9128-02	9117-02
Adult	31"	9128-03	9117-03



# ABDOMINAL BINDERS



## Premium Sized Abdominal Binder

- Constructed of latex-safe knit elastic with soft lining for extended wear and comfort
- Three (9") or four (12") individual panel design for best patient fit
- Unique stitching prevents individual panels from rolling
- Hook and loop closure

**Indications:** Abdominal muscle strain; Soreness in abdominal region; Post-abdominal surgery

Size	Measurement*	9"	12"
Small/Medium	30" - 45"	13991000	13981000
Medium/Large	46" - 62"	13992000	13982000
Large	63" - 74"	13993000	13983000
X-Large	75" - 84"	13994000	13984000
XX-Large	85" - 94"	13995000	13985000

\*Measure circumference of male waist or female hips



## Premium Universal Abdominal Binder

- Constructed of latex-safe knit elastic
- Choose between three (9") or four (12") individual or solid (10") panel design
- Unique stitching design prevents individual panels from rolling
- Hook and loop closure
- Universal sizing reduces inventory needs
- Five panel (15") ilio bypass version for ultimate support

**Indications:** Abdominal muscle strain; Soreness in abdominal region; Post-abdominal surgery

Size	Measurement	9"	10"	12"	15"
Universal	27" - 48"	13970000	13950000	13960000	—
X-Long	up to 54"	—	—	—	13871000
XX-long	up to 70"	—	—	—	13872900



## Sized Abdominal Binder

- Latex-safe elastic material
- Three (9") or four (12") panel design
- Hook and loop closure

**Indications:** Abdominal muscle strain; Soreness in abdominal region; Post-abdominal surgery

Size	Measurement*	9"	12"
Small/Medium	30" - 45"	13661056	13651056
Medium/Large	46" - 62"	13662067	13652067
Large	63" - 74"	13663008	13653008
X-Large	75" - 84"	13664009	13654009
XX-Large	85" - 94"	13665010	13655010

\*Measure circumference of male waist or female hips



# ABDOMINAL BINDERS

## Universal Abdominal Binder

- Latex-safe elastic material
- Choose between three (9") or four (12") individual or solid (10") panel design
- Hook and loop closure
- Universal sizing reduces inventory needs

**Indications:** Abdominal muscle strain; Soreness in abdominal region;  
Post-abdominal surgery



Size	Measurement	9"	10"	12"
Universal	27" - 48"	13660000	13640000	13650000
X-Long	up to 54"	—	13600000	—
XX-Long	up to 70"	—	13590000	13630000

# RIB BELTS



## Premium 6" Rib Belt

- Constructed of latex-safe knit elastic with soft inner liner
- Two panel design with unique stitching prevents material from rolling
- Hook and loop closure for ease of application
- Male and female (contoured) styles for best fit
- Available in universal and sized options

**Indications:** Fractured or bruised ribs; Intercostal strains; Costovertebral sprains

Size	Dimensions	Male	Female
Universal	6" x 40"	13780300	13780400
Small	6" x 26"	13790305	13790405
Medium	6" x 31"	13790306	13790406
Large	6" x 36"	13790307	13790407
X-Large	6" x 43"	13790308	13790408
XX-Large	6" x 50"	13790309	13790409



## Premium 8" Rib Belt

- Constructed of latex-safe knit elastic with soft inner liner
- Three panel design with unique stitching prevents material from rolling
- Hook and loop closure for ease of application
- Male and female (contoured) styles for best fit
- Available in universal and sized options

**Indications:** Fractured or bruised ribs; Intercostal strains; Costovertebral sprains

Size	Dimensions	Male	Female
Universal	8" x 38"	13760300	13760400
Small	8" x 26"	13760305	13760405
Medium	8" x 31"	13760306	13760406
Large	8" x 36"	13760307	13760407
X-Large	8" x 43"	13760308	13760408



## Universal Rib Belt

- Constructed of latex-safe knit elastic
- 6" features solid panel design; 8" features three panel design
- Hook and loop closure for ease of application
- Available in universal, gender specific styles for best fit

**Indications:** Fractured or bruised ribs; Intercostal strains; Costovertebral sprains

Size	Dimensions	Qty	Male	Female
Universal	6" x 40"	12/Case	13740300	13740400
Universal X-Long	6" x 43"	12/Case	13749300	13749400
Universal	8" x 38"	1 Each	13060300	13060400

# LUMBO SACRAL SUPPORTS

## Criss-Cross Lumbo Sacral Support

- Durable latex-safe elastic construction with hook and loop closure
- Contoured in waist area for maximum comfort
- Two elastic side pulls for optimal compression
- Four adjustable posterior stays for stability and support
- 12" posterior height

**Indications:** Mild low back pain; Lumbar muscle weakness; Lumbar sprains and strains; Mechanical or discogenic thoracolumbar pain

Size	Waist Circumference	Product #
Universal	26" - 46"	13850000
Small	24" - 28"	13850005
Small/Medium	28" - 32"	13850056
Medium/Large	32" - 36"	13850067
Large	36" - 40"	13850007
X-Large	40" - 44"	13850008
XX-Large	44" - 48"	13850009
XXX-Large	48" - 52"	13850010
XXXX-Large	52" - 56"	13850011



## Elastic Lumbo Sacral Support

- Constructed of durable latex-safe elastic with hook and loop closure
- Two elastic side pulls allow optimal compression
- Six flexible posterior stays for stability and support
- 10" posterior height

**Indications:** Mild low back pain; Lumbar muscle weakness; Lumbar sprains and strains; Mechanical or discogenic thoracolumbar pain

Size	Waist Circumference	Product #
Universal	26" - 46"	13810000
X-Small	28" - 32"	13810004
Small	32" - 36"	13810005
Medium	36" - 40"	13810006
Large	40" - 44"	13810007
Large, X-Long	40" - 56"	13810907
X-Large	44" - 48"	13810008



# LUMBO SACRAL SUPPORTS



## Heat & Mold Lumbo Sacral Support

- Constructed of soft latex-safe elastic material
- Posterior pocket accommodates moldable thermoplastic insert (sold separately)
- 8" posterior height

**Indications:** Moderate to severe low back pain; Lumbar muscle weakness; Lumbar sprains and strains; Lumbar instability

Size	Waist Circumference	Product #
Small	24" - 32"	13800005
Medium	32" - 38"	13800006
Large	38" - 46"	13800007
X-Large	46" - 54"	13800008
Moldable Insert		13860900



## Elastic/Neoprene Back Support

- Constructed of beige or black elastic material
- Contoured shape for increased comfort in seated position
- Two side pulls for optimal compression
- Neoprene posterior pocket accommodates moldable thermoplastic insert (sold separately)
- 9" posterior height

**Indications:** Moderate to severe low back pain; Lumbar muscle weakness; Lumbar sprains and strains; Lumbar instability

Size	Waist Circumference	Beige	Black
Small	26" - 30"	NE7741-02	NE7741-12
Medium	30" - 34"	NE7741-03	NE7741-13
Large	34" - 38"	NE7741-04	NE7741-14
X-Large	38" - 42"	NE7741-05	NE7741-15
XX-Large	42" - 46"	NE7741-06	NE7741-16
XXX-Large	46" - 50"	NE7741-07	NE7741-17
Moldable Insert		13860900	

# LUMBO SACRAL SUPPORTS

## S-I Support Belt

- Constructed of non-slip neoprene material to provide compression and allow stretch for movement
- Inner adjustable pads can be customized to provide anterior and posterior support to S-I joint
- Provides stability to S-I joint reducing strain on the S-I ligaments decreasing pain and promoting faster healing
- Thin material is low profile allowing belt to be worn under clothes
- Simple hook closure eliminates need for bulky buckles

**Indications:** Sacroiliac pain; Sacroiliac joint stability; Sacroiliac ligament strains

Size	Length	Product #
Small/Medium	42"	13480056
Large/X-Large	54"	13480078



## Workers Back Support

- Constructed of high strength, breathable nylon material
- Five self-contouring stays
- Rubber lining on stay pockets help prevent migration
- 1½" suspension straps adjust easily and lay flat
- Hook and loop closure with vinyl covered side pulls for maximum support
- Available in both regular (9.5") and concise (7") posterior heights
- Support available with or without suspenders

**Indications:** Mild low back pain; Lumbar muscle weakness; Lumbar sprains and strains

Size	Waist Circumference	w/ Suspenders	w/o Suspenders	Concise w/ Suspenders
X-Small	28" - 32"	BS1000-10	BS1000-50	—
Small	32" - 36"	BS1000-12	BS1000-52	BS5000-12
Medium	36" - 40"	BS1000-14	BS1000-54	BS5000-14
Large	40" - 44"	BS1000-16	BS1000-56	BS5000-16
X-Large	44" - 48"	BS1000-18	BS1000-58	BS5000-18
XX-Large	48" - 52"	BS1000-20	BS1000-59	BS5000-20
XXX-Large	52" - 56"	BS1000-99	BS1000-40	—



# PELVIC TRACTION



## Paneled Pelvic Traction Belt

- 9", three panel elastic construction prevents material from rolling
- Soft lining for extended wear comfort
- Double strap design
- Hook and loop closure for easy application

**Indications:** Lumbar spine traction; Degenerative disc disease; Herniated nucleus pulposus with and without radiculopathy; Spinal stenosis

Size	Hip Circumference	Product #
Small	20" - 40"	13731000
Medium	40" - 60"	13732000
Large	60" - 80"	13733000



## Specialty Pelvic Traction Belt

- Constructed of perforated foam material
- Double strap design with buckle adjustment
- Removable straps adhere to any position on foam belt
- Hook and loop closure for easy application

**Indications:** Lumbar spine traction; Degenerative disc disease; Herniated nucleus pulposus with and without radiculopathy; Spinal stenosis

Size	Measurement	Product #
Universal	9" x 56"	13691000

## Solace™ Low Profile Back Brace

- Soft, breathable material provides cool, comfortable fit increasing patient compliance
- Adjustable pulley style closure system is easy to secure and provides maximum compression
- Anterior and posterior ABS plastic inserts provide control and support
- Low profile, lightweight design fits easily under patient's clothing
- Plastic d-rings for ease of application and removal

**Indications:** Chronic low back pain; Lumbar muscle weakness; Lumbar sprains and strains; Mechanical or discogenic lumbar pain



Size	Waist Circumference	Product #
Small	26" - 34"	13400005
Medium	32" - 40"	13400006
Large	38" - 46"	13400007
X-Large	44" - 51"	13400008

## Prolign® Lumbar Orthosis

- Constructed of hook receivable material with soft breathable liner that wicks away moisture to help prevent skin breakdown
- Patented dual dynamic closure system provides instant full range adjustment and maximum compression
- Rigid anterior insert provides support to intra-abdominal and viscera region helping to unload the spine under gravitational load
- Available in 15° and 25° spot moldable posterior ABS plastic inserts
- Plastic d-rings for ease of application and removal

**Indications:** Post-operative lumbar immobilization; Chronic low back pain; Lumbar muscle weakness; Lumbar strains and sprains; Mechanical or discogenic disc disease; Postural support; Degenerative disc disease; IDET procedure; Spinal blocks



Size	Waist Circumference	15°	25°
X-Small	25" - 30"	13501004	13500004
Small	30" - 35"	13501005	13500005
Medium	35" - 40"	13501006	13500006
Large	40" - 45"	13501007	13500007
X-Large	45" - 50"	13501008	13500008
XX-Large	50" - 55"	13501009	13500009



## UltraAlign® + Lumbar Sacral Orthosis



Best Choice



- Constructed of hook receivable material with soft, breathable liner that provides a cooler, lighter fit
- Patented dual dynamic closure system allows infinite adjustability and maximum compression
- Combines heat-moldable plastic inserts (anterior, posterior, lateral) for rigid support
- Available in 0°, 15° or 25° lordosis
- Non-Tapered (male) and tapered (female) versions available for best patient fit
- Mid and low profile versions for short-waisted patients
- LSO meets requirements for most L1 to L5 pathologies

**Indications:** Post-operative lumbar fusion; Laminectomy or discectomy; Lumbar compression fractures; Degenerative disc disease; Osteoporosis; Spondylolisthesis; Spondylolysis; Spinal stenosis; Facet syndrome

### Sizing Guide

Sizing determined by three circumferential measurements:

1. 1" below tip of breast bone or xyphoid process
2. Around navel or waist
3. Around fullest area of trochanter or hips

Non-Tapered	Xyphoid	Waist	Trochanter
Small	31" - 35"	31" - 35"	31" - 35"
Medium	34" - 38"	32" - 38"	34" - 39"
Large	37" - 41"	35" - 41"	38" - 42"
X-Large	38" - 45"	38" - 45"	38" - 45"
XX-Large	45" - 53"	45" - 53"	45" - 53"

Tapered	Xyphoid	Waist	Trochanter
Small	25" - 31"	25" - 31"	27" - 33"
Medium	28" - 35"	29" - 36"	32" - 38"
Large	32" - 38"	32" - 38"	35" - 41"
X-Large	35" - 42"	36" - 42"	38" - 46"
XX-Large	40" - 46"	39" - 47"	43" - 49"

### Product #

Non-Tapered	0°	15°	25°
Small	UPL1330-0	UPL1330-15	UPL1330-25
Medium	UPL1430-0	UPL1430-15	UPL1430-25
Large	UPL1530-0	UPL1530-15	UPL1530-25
X-Large	UPL1630-0	UPL1630-15	UPL1630-25
XX-Large	UPL1730-0	UPL1730-15	UPL1730-25

Tapered	0°	15°	25°
Small	UPL2330-0	UPL2330-15	UPL2330-25
Medium	UPL2430-0	UPL2430-15	UPL2430-25
Large	UPL2530-0	UPL2530-15	UPL2530-25
X-Large	UPL2630-0	UPL2630-15	UPL2630-25
XX-Large	UPL2730-0	UPL2730-15	UPL2730-25

Low Profile Tapered	15°	25°	Mid Profile Non-Tapered	15°
Small	UPLP2330-15	UPLP2330-25	Small	UPMP1330-15
Medium	UPLP2430-15	UPLP2430-25	Medium	UPMP1430-15
Large	UPLP2530-15	UPLP2530-25	Large	UPMP1530-15
X-Large	UPLP2630-15	UPLP2630-25	X-Large	UPMP1630-15
XX-Large	UPLP2730-15	UPLP2730-25	XX-Large	UPMP1730-15

## Ultralign®+ Thoracic Lumbar Sacral Orthosis

- Constructed of hook receivable material with soft, breathable liner that provides a cooler, lighter fit
- Patented dual dynamic closure system allows infinite adjustability and maximum compression
- Combines heat-moldable plastic inserts (anterior, posterior, lateral) for rigid support
- Available in 15° or 25° lordosis
- Non-tapered (male) and tapered (female) versions available for best patient fit
- DLE (Dorsal Lumbar Extension) or SPK (Sternal Pad Kit) available for maximum control
- Low profile version for short waisted patients
- TLSO meets requirements for most T6 to L5 pathologies



SPK



**Indications:** Post-operative thoracic/lumbar fusion; Laminectomy or discectomy; Compression fractures; Degenerative disc disease; Osteoporosis; Single column instability immobilization; Facet syndrome

### Sizing Guide

Sizing determined by three circumferential measurements:

1. 1" below tip of breast bone or xyphoid process
2. Around navel or waist
3. Around fullest area of trochanter or hips

Non-Tapered	Xyphoid	Waist	Trochanter
Small	31" - 35"	31" - 35"	31" - 35"
Medium	34" - 38"	32" - 38"	34" - 39"
Large	37" - 41"	35" - 41"	38" - 42"
X-Large	38" - 45"	38" - 45"	38" - 45"
XX-Large	45" - 53"	45" - 53"	45" - 53"

Tapered	Xyphoid	Waist	Trochanter
Small	25" - 31"	25" - 31"	27" - 33"
Medium	28" - 35"	29" - 36"	32" - 38"
Large	32" - 38"	32" - 38"	35" - 41"
X-Large	35" - 42"	36" - 42"	38" - 46"
XX-Large	40" - 46"	39" - 47"	43" - 49"



DLE



Low Profile

### Product #

#### Non-Tapered

	15°	25°
Small	UPT1333-15	UPT1333-25
Medium	UPT1433-15	UPT1433-25
Large	UPT1533-15	UPT1533-25
X-Large	UPT1633-15	UPT1633-25
XX-Large	UPT1733-15	UPT1733-25

#### Tapered

	15°	25°
Small	UPT2333-15	UPT2333-25
Medium	UPT2433-15	UPT2433-25
Large	UPT2533-15	UPT2533-25
X-Large	UPT2633-15	UPT2633-25
XX-Large	UPT2733-15	UPT2733-25

#### Low Profile Tapered

	15°	25°
Small	UPLT2333-15	UPLT2333-25
Medium	UPLT2433-15	UPLT2433-25
Large	UPLT2533-15	UPLT2533-25
X-Large	UPLT2633-15	UPLT2633-25
XX-Large	UPLT2733-15	UPLT2733-25

### Accessories

- Dorsal Lumbar Extension
- Sternal Pad Kit
- 3 1/2" Extension Panels
- 5 1/2" Extension Panels

### Product #

- UDLE
- USPK
- UEP2031
- UEP2051

# KNEE IMMOBILIZERS



## Deluxe Tri-Panel Knee Immobilizer

- Constructed of foam body with durable blue canvas panels
- T-bar stays provide optimal strength
- Contoured posterior stays for correct anatomical fit
- Adjustable medial and lateral panels allow for universal sizing
- Patella strap assists with application
- Also available in black foam construction

**Indications:** Mild to severe knee injuries; Post-operative immobilization; Grade II and III collateral ligament sprain; Patellar dislocations; Ostrochondral fractures

Length	Canvas	Foam
12"	1120917	BF1120917
14"	1140917	—
16"	1160917	BF1160917
19"	1190917	BF1190917
22"	1220917	BF1220917
24"	1240917	BF1240917

\*Universal size fits thigh circumference 12" - 24"



## Superlite Tri-Panel Knee Immobilizer

- Constructed of foam body with white perforated canvas panels
- T-bars stays provide optimal strength
- Contoured posterior stays for correct anatomical fit
- Adjustable medial and lateral panels allow for universal sizing
- Patella strap assists with application

**Indications:** Mild to severe knee injuries; Post-operative immobilization; Grade II and III collateral ligament sprain; Patellar dislocations; Ostrochondral fractures

Length	Product #
12"	7091-50
16"	7091-56
19"	7091-51
22"	7091-55
24"	7091-52
26"	7091-57
28"	7091-58

\*Universal size fits thigh circumference 12" - 24"

# KNEE IMMOBILIZERS

## Foam Tri-Panel Knee Immobilizer

- Constructed of blue perforated foam for patient comfort
- T-bar stays provide optimal strength
- Contoured posterior stays for correct anatomical fit
- Adjustable medial and lateral panels allow for universal sizing
- Trimmable foam allows for a more customized fit

**Indications:** Mild to severe knee injuries; Post-operative immobilization; Grade II and III collateral ligament sprain; Patellar dislocations; Ostrochondral fractures

Length	Product #
12"	7091-03
16"	7091-04
19"	7091-05
22"	7091-06
24"	7091-07
26"	7091-08

\*Universal size fits thigh circumference 12" - 24"



## Tietex Tri-Panel Knee Immobilizer

- Constructed of black foam body with durable tietex panels
- T-bar stays provide optimal strength
- Contoured posterior stays for correct anatomical fit
- Adjustable medial and lateral panels allow for universal sizing
- Patella strap assists with application

**Indications:** Mild to severe knee injuries; Post-operative immobilization; Grade II and III collateral ligament sprain; Patellar dislocations; Ostrochondral fractures

Length	Product #
12"	A141000
16"	A142000
19"	A143000
22"	A144000
24"	A145000
26"	A146000

\*Universal size fits thigh circumference 12" - 24"



# KNEE IMMOBILIZERS



## Cutaway Knee Immobilizer

- Perforated foam construction for patient comfort and breathability
- T-bar stays provide optimal strength
- Available with straight or contoured posterior stays
- Adjustable medial and lateral stays allow for universal sizing
- Loop lock straps affix posterior to encircle entire leg
- Also available in black foam construction

**Indications:** Mild to severe knee injuries; Post-operative immobilization; Grade II and III collateral ligament sprain; Patellar dislocations; Ostrochondral fractures

Length	T-bars	Contoured	Straight	Black Foam
12"	—	—	—	BF4433-12
16"	4433-11	4433-10	—	BF4433-16
19"	—	4433-01	4433-02	BF4433-19
22"	4433-22	4433-05	—	BF4433-22
24"	4433-33	4433-03	4433-04	BF4433-24
26"	—	—	—	BF4433-26

\*Universal size fits thigh circumference 12" - 24"



## Knee Immobilizer w/ Elastic Straps

- Perforated foam construction for patient comfort and breathability
- T-bar stays provide optimal strength
- Adjustable medial and lateral stays allow for universal sizing
- Elastic straps provide compression helping to reduce edema

**Indications:** Mild to severe knee injuries; Post-operative immobilization; Grade II and III collateral ligament sprain; Patellar dislocations; Ostrochondral fractures

Length	Blue	Black
12"	—	BF7072-01
16"	7072-02	BF7072-02
20"	7072-03	BF7072-03
22"	7072-04	BF7072-04
24"	7072-05	BF7072-05
26"	—	BF7072-06

\*Universal size fits thigh circumference 12" - 24"



## Cutaway Knee Immobilizer

- Perforated foam construction for patient comfort and breathability
- T-bar stays provide optimal strength
- Available with straight or contoured posterior stays
- Adjustable medial and lateral stays allow for universal sizing
- Quick fit tabs assist with application

**Indications:** Mild to severe knee injuries; Post-operative immobilization; Grade II and III collateral ligament sprain; Patellar dislocations; Ostrochondral fractures

Length	Straight	Contoured
12"	7041-01	7041-11
16"	7041-02	7041-22
20"	7041-03	7041-33
24"	7041-04	7041-44
28"	7041-05	7041-55

\*Universal size fits thigh circumference 12" - 24"

# KNEE IMMOBILIZERS

## Cutaway Knee Immobilizer w/ Quick Fit Tabs

- Perforated foam construction for patient comfort and breathability
- T-bar stays provide optimal strength
- Contoured posterior stays for correct anatomical fit
- Adjustable medial and lateral stays allow for universal sizing
- Quick fit tabs assist with application

**Indications:** Mild to severe knee injuries; Post-operative immobilization; Grade II and III collateral ligament sprain; Patellar dislocations; Osteochondral fractures

Length	w/ Tabs	w/o Tabs
12"	—	7081-00
16"	7082-22	7082-00
19"	7083-33	7083-00
22"	7084-44	7084-00
24"	7085-55	7085-00
26"	7086-66	7086-00
28"	—	7087-00

\*Universal size fits thigh circumference 12" - 24"



## Knee Immobilizer w/ Elastic Straps

- Perforated foam construction for patient comfort and breathability
- Available with or without t-bar stays
- Adjustable medial and lateral stays allow for universal sizing
- Elastic straps provide compression helping to reduce edema

**Indications:** Mild to severe knee injuries; Post-operative immobilization; Grade II and III collateral ligament sprain; Patellar dislocations; Osteochondral fractures

Length	w/o T-bars	w/ T-bars	40" Circumference
13"	7071-01	7071-11	—
16"	7071-02	7071-12	—
18"	—	—	7071-55
19"	7071-03	7071-13	—
23"	7071-05	7071-15	7071-57
28"	7071-07	7071-17	7071-59

\*Universal size fits thigh circumference 12" - 24"



# KNEE IMMOBILIZERS



## Cutaway Knee Immobilizer w/ Patella Strap

- Perforated foam construction for patient comfort and breathability
- Available with or without t-bar stays
- Contoured posterior stays for correct anatomical fit
- Adjustable medial and lateral stays allow for universal sizing
- Elastic patella strap helps reduce knee edema
- Also available in black foam construction (w/o t-bars)

**Indications:** Mild to severe knee injuries; Post-operative immobilization; Grade II and III collateral ligament sprain; Patellar dislocations; Osteochondral fractures

Length	w/ T-bars	w/o T-bars
13"	—	9307-01
14"	7066-14	—
16"	7066-16	—
18"	7066-18	9307-02
20"	7066-20	9307-03
22"	7066-22	9307-06
24"	7066-24	9307-04
26"	7066-26	—
28"	7066-28	9307-05

\*Universal size fits thigh circumference 12" - 24"



## Sized Canvas Knee Immobilizer

- Constructed of durable blue canvas
- Available with or without t-bar stays
- Contoured posterior stays for correct anatomical fit
- Alternating web straps with loop lock closure

**Indications:** Mild to severe knee injuries; Post-operative immobilization; Grade II and III collateral ligament sprain; Patellar dislocations; Osteochondral fractures

w/o T-bar	Measurement*	12"	16"	19"	22"	24"
Small	14" - 16"	1010127	1010167	1010197	1010227	1010247
Medium	16" - 18"	1020127	1020167	1020197	1020227	1020247
Large	18" - 20"	1030127	1030167	1030197	1030227	1030247
X-Large	20" - 22"	1040127	1040167	1040197	1040227	1040247
XX-Large	22" - 24"	1050127	1050167	1050197	1050227	1050247

w/ T-bar	Measurement*	12"	16"	19"	22"
Small	14" - 16"	1110127	1110167	1110197	1110227
Medium	16" - 18"	1120127	1120167	1120197	1120227
Large	18" - 20"	1130127	1130167	1130197	1130227
X-Large	20" - 22"	1140127	1140167	1140197	1140227
XX-Large	22" - 24"	1150127	1150167	1150197	1150227

\*Measure thigh circumference 6" above mid-patella



# KNEE IMMOBILIZERS

## Sized Black Foam Knee Immobilizer

- Constructed of black foam with tricot liner
- T-bar stays provide optimal strength
- Contoured posterior stays for correct anatomical fit
- Alternating straps with loop lock closure provide maximum compression

**Indications:** Mild to severe knee injuries; Post-operative immobilization; Grade II and III collateral ligament sprain; Patellar dislocations; Ostrochondral fractures



Size	Measurement*	10"	12"	16"	20"	24"
Small	14" - 16"	—	BF7006-01	BF7007-01	BF7008-01	BF7009-01
Medium	16" - 18"	—	BF7006-02	BF7007-02	BF7008-02	BF7009-02
Large	18" - 20"	—	BF7006-03	BF7007-03	BF7008-03	BF7009-03
X-Large	20" - 22"	—	BF7006-04	BF7007-04	BF7008-04	BF7009-04
XX-Large	22" - 24"	—	BF7006-05	BF7007-05	BF7008-05	BF7009-05
XXX-Large	24" - 26"	—	BF7006-06	BF7007-06	BF7008-06	BF7009-06
Pediatric Small	8" - 10"	14970005				
Pediatric Medium	10" - 12"	14970006				
Pediatric Large	12" - 14"	14970007				

\*Measure thigh circumference 6" above mid-patella

## Sized Superlite Knee Immobilizer

- Constructed of lightweight perforated canvas material with tricot liner
- T-bar stays provide optimal strength
- Available with straight or contoured posterior stays
- Foam popliteal pad included with straight stay design
- Alternating web straps with loop lock closure

**Indications:** Mild to severe knee injuries; Post-operative immobilization; Grade II and III collateral ligament sprain; Patellar dislocations; Ostrochondral fractures



Straight Stays	Measurement*	12"	16"	20"	24"
Small	14" - 16"	7001-01	7002-01	7003-01	7004-01
Medium	16" - 18"	7001-02	7002-02	7003-02	7004-02
Large	18" - 20"	7001-03	7002-03	7003-03	7004-03
X-Large	20" - 22"	7001-04	7002-04	7003-04	7004-04
XX-Large	22" - 24"	7001-05	7002-05	7003-05	7004-05
Contoured Stays	Measurement*	12"	16"	20"	24"
Small	14" - 16"	7006-01	7007-01	7008-01	7009-01
Medium	16" - 18"	7006-02	7007-02	7008-02	7009-02
Large	18" - 20"	7006-03	7007-03	7008-03	7009-03
X-Large	20" - 22"	7006-04	7007-04	7008-04	7009-04
XX-Large	22" - 24"	7006-05	7007-05	7008-05	7009-05

\*Measure thigh circumference 6" above mid-patella

# KNEE IMMOBILIZERS



## Sized Foam Knee Immobilizer

- Perforated foam construction for patient comfort and breathability
- T-bar stays provide optimal strength
- Contoured posterior stays for correct anatomical fit
- Alternating straps with loop lock closure provide maximum compression

**Indications:** Mild to severe knee injuries; Post-operative immobilization; Grade II and III collateral ligament sprain; Patellar dislocations; Osteochondral fractures

Size	Measurement*	16"	20"	22"	24"
Small	14" - 16"	4461-02	4462-02	4464-02	4463-02
Medium	16" - 18"	4461-03	4462-03	4464-03	4465-03
Large	18" - 20"	4461-04	4462-04	4464-04	4463-04
X-Large	20" - 22"	4461-05	4462-05	4464-05	4463-05
XX-Large	22" - 24"	—	4462-06	—	4463-06

\*Measure thigh circumference 6" above mid-patella



## Sized Tietex Knee Immobilizer

- Constructed of durable tietex foam
- T-bar stays provide optimal strength
- Contoured posterior stays for correct anatomical fit
- Alternating straps with loop lock closure provide maximum compression

**Indications:** Mild to severe knee injuries; Post-operative immobilization; Grade II and III collateral ligament sprain; Patellar dislocations; Osteochondral fractures

Size	Measurement*	14"	16"	20"	22"	24"
Small	14" - 16"	A145005	A146005	A142005	A143005	A144005
Medium	16" - 18"	A145006	A146006	A142006	A143006	A144006
Large	18" - 20"	A145007	A146007	A142007	A143007	A144007
X-Large	20" - 22"	A145008	A146008	A142008	A143008	A144008
XX-Large	22" - 24"	A145009	A146009	A142009	A143009	A144009

\*Measure thigh circumference 6" above mid-patella

## Closed Patella Knee Support

- 1/8" neoprene construction provides warmth, compression and support
- Available with oval pad for additional patella protection

**Indications:** Grade I collateral ligament sprain; Knee sprains and strains; Osteoarthritis; Mild joint effusion; Baker's cyst; Bursitis; Pes anserinus

Size	Measurement*	w/o Pad	w/ Pad
X-Small	13½" - 15½"	NE7701-71	NE7701-01
Small	15½" - 18"	NE7701-72	NE7701-02
Medium	18" - 20½"	NE7701-73	NE7701-03
Large	20½" - 23"	NE7701-74	NE7701-04
X-Large	23" - 25½"	NE7701-75	NE7701-05
XX-Large	25½" - 28"	NE7701-76	NE7701-06
XXX-Large	28" - 30½"	NE7701-77	NE7701-07

\*Measure circumference 6" above mid-patella



## Open Patella Knee Support

- 1/8" neoprene construction provides warmth, compression and support
- Available with oval pad for additional patella protection
- Also available in breathable Tri-tex construction (w/o pad)
- Open patella for patient comfort

**Indications:** Grade I collateral ligament sprain; Knee sprains and strains; Osteoarthritis; Mild joint effusion; Baker's cyst; Bursitis; Pes anserinus

Size	Measurement*	w/o Pad	w/ Pad	Tri-tex
X-Small	13½" - 15½"	NE7702-71	NE7702-01	14811004
Small	15½" - 18"	NE7702-72	NE7702-02	14811005
Medium	18" - 20½"	NE7702-73	NE7702-03	14811006
Large	20½" - 23"	NE7702-74	NE7702-04	14811007
X-Large	23" - 25½"	NE7702-75	NE7702-05	14811008
XX-Large	25½" - 28"	NE7702-76	NE7702-06	14811009
XXX-Large	28" - 30½"	NE7702-77	NE7702-07	14811010
XXXX-Large	30½" - 33"	—	NE7702-08	—

\*Measure circumference 6" above mid-patella



# KNEE SUPPORTS



## Universal Knee Support

- 1/8" neoprene construction provides warmth, compression and support
- Anterior wrap around strapping above and below knee helps prevent migration
- Enclosed tubular buttress for desired patella position
- Universal adult and pediatric sizes available

**Indications:** Grade I collateral ligament sprain; Patellar subluxation; Chondromalacia; Patella tendonitis; Patellofemoral tracking dysfunction; Quadriceps tendonitis; Plica syndrome

Size	Measurement*	Product #
Universal	up to 25"	NE7710-70
Pediatric	up to 12"	NE7710-71

\*Measure circumference 6" above mid-patella



## Knee Support w/ Trimmable Buttress

- 1/8" neoprene construction provides warmth, compression and support
- Enclosed tubular buttress can be trimmed for desired patella position

**Indications:** Grade I collateral ligament sprain; Patellar subluxation; Chondromalacia; Patella tendonitis; Patellofemoral tracking dysfunction; Quadriceps tendonitis; Plica syndrome

Size	Measurement*	Product #
X-Small	13 1/2" - 15 1/2"	NE7717-71
Small	15 1/2" - 18"	NE7717-72
Medium	18" - 20 1/2"	NE7717-73
Large	20 1/2" - 23"	NE7717-74
X-Large	23" - 25 1/2"	NE7717-75
XX-Large	25 1/2" - 28"	NE7717-76
XXX-Large	28" - 30 1/2"	NE7717-77

\*Measure circumference 6" above mid-patella

## Deluxe Knee Support w/ Trimmable Buttress

- 1/8" neoprene construction provides warmth, compression and support
- Anterior strapping above and below knee helps prevent migration while medial/lateral stays provide additional stability
- Enclosed tubular buttress can be trimmed for desired patella position
- Open popliteal option for patient comfort

**Indications:** Grade I collateral ligament sprain; Patellar subluxation; Chondromalacia; Patella tendonitis; Patellofemoral tracking dysfunction; Patellar dislocation; Quadriceps tendonitis; Plica syndrome



Size	Measurement*	Closed Popliteal	Open Popliteal
X-Small	13½" - 15½"	NE7718-71	14882004
Small	15½" - 18"	NE7718-72	14882005
Medium	18" - 20½"	NE7718-73	14882006
Large	20½" - 23"	NE7718-74	14882007
X-Large	23" - 25½"	NE7718-75	14882008
XX-Large	25½" - 28"	NE7718-76	14882009
XXX-Large	28" - 30½"	NE7718-77	14882010

\*Measure circumference 6" above mid-patella

## Fan Knee Stabilizer

- 1/8" neoprene construction provides warmth, compression and support
- Features four criss-cross elastic straps for additional patella support
- Strapping above and below the knee helps prevent migration
- Medial/lateral stays provide stability
- Enclosed tubular buttress and adjustable felt buttress included

**Indications:** Grade I collateral ligament sprain; Patella subluxation; Chondromalacia; Patella tendonitis; Patellofemoral tracking dysfunction; Patellar dislocation; Quadriceps tendonitis; Plica syndrome



Size	Measurement*	Product #
Small	15½" - 18"	NE7707-72
Medium	18" - 20½"	NE7707-73
Large	20½" - 23"	NE7707-74
X-Large	23" - 25½"	NE7707-75
XX-Large	25½" - 28"	NE7707-76

\*Measure circumference 6" above mid-patella

# KNEE SUPPORTS



## Elastic Knee Support

- Pull up cotton elastic design
- Provides snug fit for maximum compression
- Available in open and closed patella styles

**Indications:** Grade I collateral ligament sprain; Knee sprains and strains; Osteoarthritis; Mild joint effusion; Baker's cyst; Bursitis; Pes anserinus

Size	Measurement*	Closed Patella	Open Patella
Small	12" - 15"	7042-01	7043-01
Medium	15" - 18"	7042-02	7043-02
Large	18" - 21"	7042-03	7043-03
X-Large	21" - 24"	7042-04	7043-04
XX-Large	24" - 27"	7042-05	7043-05

\*Measure circumference 6" above mid-patella



## Visco Elastic Knee Support

- Visco-elastic construction with four way stretch for maximum compression
- Silicone buttress option features spiral stays
- Closed patella design

**Indications:** Grade I collateral ligament sprain; Knee sprains and strains; Osteoarthritis; Mild joint effusion; Baker's cyst; Bursitis; Pes anserinus

Size	Measurement*	w/o Buttress	w/ Buttress
X-Small	13½" - 15½"	14761004	14801004
Small	15½" - 18"	14761005	14801005
Medium	18" - 20½"	14761006	14801006
Large	20½" - 23"	14761007	14801007
X-Large	23" - 25½"	14761008	14801008
XX-Large	25½" - 28"	14761009	14801009
XXX-Large	28" - 30½"	14761010	14801010

\*Measure circumference 6" above mid-patella



## Arthroscopic Foam Knee Wrap

- Close fitting foam wrap provides compression to help control edema
- Optional pressure-sensitive adhesive foam pads can be placed in desired position after lateral patellar retinacular release or as adjunct for patellar stabilization
- Hook tab closures attach anywhere on wrap

**Indications:** Post-arthroscopic or arthroplasty knee surgery; Moderate to severe joint effusion

Size	Product #
Universal	FF4061-00

## Transition Knee Brace

- Drop-lock design offers easy access for pin adjustment
- Range of motion can be set in 10° increments
- Superior and inferior score zones for length adjustment
- Dial can be locked to encourage patient compliance
- Circumferential loop lock straps allow for even compression
- Trimmable softgoods available in full or cool foam options

**Indications:** Post-operative stabilization and rehabilitation; Post total knee replacement; Post ORIF; Fractures of the femur, tibia, fibula or patella; Ligament reconstructive procedures, Meniscal repairs; Post-fracture management

Description	Length Adjustment	Product #
Full Foam	23" - 17"	KB5000-00
Cool Foam	23" - 17"	KB5000-01
Full Foam, X-Long	26" - 20"	KB5000-02



## Slimline II Knee Brace

- Snap-lock hinge cover design with simple pin adjustment
- Range of motion can be set at 5°, 15°, 30°, 60° or 90° of flexion and 0°, 15°, 30°, 60° or 90° of extension
- Superior and inferior score zones for length adjustment
- Circumferential loop lock straps allow for even compression
- Trimmable softgoods available in full or cool foam options

**Indications:** Post-operative stabilization and rehabilitation; Post total knee replacement; Post ORIF; Fractures of the femur, tibia, fibula or patella; Ligament reconstructive procedures, Meniscal repairs; Post-fracture management

Description	Length Adjustment	Product #
Full Foam	23" - 17"	KB6000-00
Cool Foam	23" - 17"	KB6000-01



## Slimline Knee Brace

- Posterior offset hinge for correct anatomical fit
- Hinge can be easily adjusted in 10° degree increments of flexion and extension without tools
- Paddles allow for easy length adjustment
- Unique calf strap helps prevent migration
- Dial can be locked to encourage patient compliance

**Indications:** Post-operative stabilization and rehabilitation; Post total knee replacement; Post ORIF; Fractures of the femur, tibia, fibula or patella; Ligament reconstructive procedures, Meniscal repairs; Post-fracture management

Description	Length Adjustment	Product #
Telescopic Paddles	27" - 22"	KB2040-00
Fixed Paddles	22"	KB2040-10
No Paddles	21"	KB2040-40





# POST-OP KNEE



## Key Knee Brace

- Posterior offset hinge for best anatomical fit
- Hinge can be adjusted in 10° degree increments of flexion and extension
- Length is adjustable with easy-to-use key
- Unique calf strap helps prevent migration

**Indications:** Post-operative stabilization and rehabilitation; Post total knee replacement; Post ORIF; Fractures of the femur, tibia, fibula or patella; Ligament reconstructive procedures, Meniscal repairs; Post-fracture management

Description	Length Adjustment	Product #
Telescopic Paddles	27" - 22"	KB3001-00
Fixed Paddles	22"	KB3001-10
No Paddles	21"	KB3001-40

# BUCKS TRACTION BOOTS

## Total Universal Bucks Traction Boot

- Constructed of convoluted foam and laminated hook receivable material
- Medial and lateral stays provide stability and are removable for ease of application
- Full-foot design helps prevent foot drop
- Trimmable foam body and loop lock closure for more customized fit
- Spreader bar attached

**Indications:** Hip fractures; Hip dislocations; Femur fractures

**Size**  
Universal

**Product #**  
8140-00




## Full Foot Bucks Traction Boot

- Convoluted foam construction helps prevent migration
- Medial and lateral stays provide additional stability
- Full foot design helps prevent foot drop
- 1" ventilation holes increase airflow
- Trimmable loop lock closures
- Spreader bar attached

**Indications:** Hip fractures; Hip dislocations; Femur fractures

**Size**  
Universal  
Pediatric  
X-Small  
Small  
Medium  
Large  
X-Large

**Product #**  
8130-00  
8106-01  
8106-02  
8106-03  
8106-04  
8106-05  
8106-06



## Universal Bucks Traction Boot

- Convoluted foam construction helps prevent migration
- Medial and lateral stays provide additional stability
- 1" ventilation holes increase airflow
- Trimmable foam body for more customized fit
- Universal size fits most patients
- Spreader bar attached

**Indications:** Hip fractures; Hip dislocations; Femur fractures

**Size**  
Universal

**Product #**  
8121-00



# BUCKS TRACTION BOOTS




## Sized Bucks Traction Boot

- Non-convoluted foam construction helps reduce pressure points
- Medial and lateral stays provide additional stability
- 1" ventilation holes increase airflow
- Trimmable loop lock closures
- Spreader bar attached

**Indications:** Hip fractures; Hip dislocations; Femur fractures

Size	Product #
Universal	8120-00
Pediatric	8105-01
X-Small	8105-02
Small	8105-03
Medium	8105-04
Large	8105-05
X-Large	8105-06



## Bucks Traction Boot w/ Triangle Buckle

- Convoluted foam construction helps prevent migration
- Medial and lateral stays provide additional stability
- 1" ventilation holes increase airflow
- Trimmable hook closures
- Traction cord attaches to triangle buckle

**Indications:** Hip fractures; Hip dislocations; Femur fractures

Size	Product #
X-Small	19402
Small	19403
Medium	19404
Large	19405
X-Large	19406



## Foam/Canvas Bucks Traction Boot

- Foam foot and canvas leg construction
- Medial and lateral stays provide additional stability
- 1" ventilation holes increase airflow
- Trimmable loop lock closures
- Spreader bar attached

**Indications:** Hip fractures; Hip dislocations; Femur fractures

Size	Product #
Pediatric	8103-01
X-Small	8103-02
Small	8103-03
Medium	8103-04
Large	8103-05
X-Large	8103-06

# BUCKS TRACTION BOOTS

## Speciality Bucks Traction Boot

- Convoluted foam construction helps prevent migration
- Medial and lateral stays provide additional stability
- Features full-foot design and dorsiflexion assist strap to help prevent foot drop
- 1" ventilation holes increase airflow
- Trimmable loop lock closures
- Spreader bar attached

**Indications:** Hip fractures; Hip dislocations; Femur fractures

Size	Product #
X-Small	8108-02
Small	8108-03
Medium	8108-04
Large	8108-05
X-Large	8108-06



# ANKLE



## Confor® Ankle Stirrup

- Confor® open cell polyurethane foam provides a superior fit and has proven to be one of the most conforming and shock absorbing materials on the market
- Rigid medial and lateral shells minimize ankle rotation and inversion/eversion
- Elasticized Vel-stretch® straps provide compression to help reduce edema
- Adjustable outside heel strap for maximum fit and comfort
- Low profile design fits easily in shoe
- Universal right/left
- Available in each or 6/Case quantity

**Indications:** Acute and chronic ankle sprains; Osteoarthritis; Syndesmosis sprain

Size	Height	Each	Case
Standard	10"	AB2340-00	AB2340-06
Small	8½"	AB2341-00	AB2341-06



## Air Ankle Stirrup

- Rigid medial and lateral shells minimize ankle rotation and inversion/eversion
- Air bladder allows for simple adjustment of air with inflation bulb
- Swivel straps allow even compression when applied in opposite directions
- Adjustable outside heel strap for maximum fit and comfort
- Low profile design fits easily in shoe
- Universal right/left
- Kit features: Standard air ankle stirrup, ice pack, elastic bandage and exercise band

**Indications:** Acute and chronic ankle sprains; Osteoarthritis; Syndesmosis sprain

Size	Height	Product #
Standard	10"	AB2360-00
Small	8½"	AB2361-00
Pediatric	6"	AB2362-00
Kit		SC1000-00



## Air/Gel Ankle Stirrup

- Rigid medial and lateral shells minimize ankle rotation and inversion/eversion
- Air/gel bladder can be placed in freezer for cold therapy
- Swivel straps allow even compression when applied in opposite directions
- Adjustable outside heel strap for maximum fit and comfort
- Low profile design fits easily in shoe
- Universal right/left

**Indications:** Acute and chronic ankle sprains; Osteoarthritis; Syndesmosis sprain

Size	Height	Product #
Standard	10"	AB2350-00
Small	8½"	AB2351-00

## DeRoyal® Functional Ankle Brace

- Anatomically designed uprights for a secure fit
- Hinged stirrup maintains support throughout range of motion
- Padded foam inlay for high energy absorption during activity
- Polyurethane foam padding for patient comfort and shock absorption
- Vel-stretch® straps for maximum support and edema control
- Universal right/left

**Indications:** Acute and chronic ankle sprains; Chronic instability; Prophylactic talocrural support

Size	Male Shoe Size	Female Shoe Size	Product #
X-Small	—	5 - 7	AB2000-10
Small	5 - 8	7½ - 10	AB2000-12
Medium	8½ - 12	10½ - 13	AB2000-14
Large	12½ - 16	13½ +	AB2000-16



## DeRoyal® Sports Ankle Brace 2

- Constructed of durable nylon for long lasting wear
- Universal stabilization strap forces ankle into inversion or eversion
- Cuff closure for added support and compression
- Lace-up design for secure fit and primary stabilization
- Low profile design allows use in athletic shoe
- Universal right/left

**Indications:** Chronic ankle instability; Prophylactic talocrural support; Prevention or treatment of acute ankle sprains; Subtalar joint sprain; Syndesmosis sprain

Size	Ankle Circumference*	Product #
X-Small	10" - 11"	15600004
Small	11" - 12"	15600005
Medium	12" - 13"	15600006
Large	13" - 14"	15600007
X-Large	14" - 15"	15600008
XX-Large	15" - 16"	15600009

\*Measure circumference from the base of heel around the ankle



## DeRoyal® Sports Ankle Brace

- Constructed of durable nylon for long lasting wear
- Figure 8 strapping configuration, in conjunction with medial to lateral stirrup straps, places foot in slightly everted position for maximum ankle stability
- Removable medial and lateral plastic stays with malleoli cut-out provides comfort and support
- Lace-up design for a secure fit
- Low profile design allows use in athletic shoe

**Indications:** Chronic ankle instability; Prophylactic talocrural support; Prevention or treatment of acute ankle sprains; Subtalar joint sprain; Syndesmosis sprain

Size	Male Shoe Size	Female Shoe Size	Right	Left
X-Small	4 - 6	5 - 7	AB2801-20	AB2801-10
Small	6½ - 8	7½ - 9	AB2801-22	AB2801-12
Medium	8½ - 10	9½ - 11	AB2801-24	AB2801-14
Large	10½ - 12	11½ - 13	AB2801-26	AB2801-16
X-Large	12½ +	13½ +	AB2801-28	AB2801-18



# ANKLE



## R.I.C.E. Ankle Splint

- Constructed of perforated vinyl material with foam lining
- Removable medial and lateral plastic stays with malleoli cut-out for comfort and support
- Lace-up design with figure 8 strap for a secure fit
- Universal right/left

**Indications:** Acute and chronic ankle sprains; Ankle rehabilitation

Size	Ankle Circumference*	Product #
Child	6" - 7"	4050-00
Small	7" - 8"	4050-01
Medium	8" - 10"	4050-02
Large	10" - 12"	4050-03
X-Large	12" - 14"	4050-04

\*Measure circumference from the base of heel around the ankle



## Canvas Lace-Up Ankle Brace

- Constructed of durable canvas with flannel lining
- Contoured design for best anatomical fit
- Removable medial and lateral metal stays provide stability
- Lace-up design for a secure fit
- Universal right/left

**Indications:** Acute and chronic ankle sprains; Ankle rehabilitation

Size	Ankle Circumference*	Product #
Small	7" - 8"	8761-01
Medium	8" - 10"	8761-02
Large	10" - 12"	8761-03
X-Large	12" - 14"	8761-04

\*Measure circumference from the base of heel around the ankle



## Leatherette Ankle Brace

- Constructed of durable leatherette material with flannel lining
- Contoured design for best anatomical fit
- Removable medial and lateral metal stays provide stability
- Loop lock closure for compression and support
- Universal right/left

**Indications:** Mild to moderate ankle sprains and strains; Ankle edema; Ankle rehabilitation

Size	Ankle Circumference*	Product #
Small	7" - 8"	4013-01
Medium	8" - 10"	4013-02
Large	10" - 12"	4013-03
X-Large	12" - 14"	4013-04

\*Measure circumference from the base of heel around the ankle



## Neoprene Ankle Support

- 1/8" neoprene construction provides warmth, compression and support
- Open heel, pull on design
- Universal right/left

**Indications:** Mild ankle sprains and strains; Chronic ankle instability; Ankle edema; Osteoarthritis; Joint effusion

Size	Male Shoe Size	Female Shoe Size	Product #
Small	4½ - 6	5½ - 7	NE7732-72
Medium	6½ - 8	7½ - 9	NE7732-73
Large	8½ - 10	9½ - 11	NE7732-74
X-Large	10½ - 12	11½ - 13	NE7732-75
XX-Large	12½ +	13½ +	NE7732-76



## Figure 8 Wrap Ankle Support

- Constructed of a durable, breathable mesh sock material
- Elastic strap provides compression and edema control
- Figure 8 wrap with hook and loop closure provides additional support
- Universal right/left

**Indications:** Mild ankle sprains and strains; Ankle edema; Peroneal and anterior tibialis tendonitis; Osteoarthritis of the talocrural and subtalar joints

Size	Ankle Circumference*	Product #
Small	7" - 8"	4016-01
Medium	8" - 10"	4016-02
Large	10" - 12"	4016-03
X-Large	12" - 14"	4016-04

\*Measure circumference from the base of heel around the ankle



## Elastic Ankle Sleeve

- Pull on cotton elastic design
- Provides snug fit for maximum comfort and compression
- Low profile design fits easily in shoe
- Universal right/left

**Indications:** Mild ankle sprains and strains; Ankle edema; Osteoarthritis of the talocrural and subtalar joints

Size	Ankle Circumference*	Product #
Small	7" - 8"	4004-01
Medium	8" - 10"	4004-02
Large	10" - 12"	4004-03
X-Large	12" - 14"	4004-04

\*Measure circumference from the base of heel around the ankle



# ANKLE



## Closed Heel Elastic Ankle Sleeve

- Closed heel, elastic construction
- Provides snug fit for maximum comfort and compression
- Pull-on design for ease of application
- Low profile fits easily in shoe
- Universal right/left

**Indications:** Mild ankle sprains and strains; Ankle edema; Osteoarthritis of the talocrural and subtalar joints

Size	Ankle Circumference*	Product #
Small	7" - 8"	4005-01
Medium	8" - 10"	4005-02
Large	10" - 12"	4005-03
X-Large	12" - 14"	4005-04

\*Measure circumference from the base of heel around the ankle



## Deluxe Ankle Foot Orthosis

- Closed heel support provides additional security
- Can be modified with heat gun and trimmed for customized fit
- Hook and loop closure for easy application
- Low profile design fits easily in shoe
- Available in right/left design for best fit

**Indications:** Patients suffering from neurological or orthopedic conditions; Lower leg instabilities

Size	Male Shoe Size	Female Shoe Size	Right	Left
Small	6 - 7½	7½ - 9	AF1001-01	AF1001-04
Medium	8 - 9 ½	9½ - 11	AF1001-02	AF1001-05
Large	10 - 11½	11½ - 13	AF1001-03	AF1001-06
X-Large	12 +	13½ +	AF1001-07	AF1001-08



## Ankle Foot Orthosis

- Open back and heel design
- Can be modified with heat gun and trimmed for customized fit
- Hook and loop closure for easy application
- Low profile design fits easily in shoe
- Universal in size

**Indications:** Patients suffering from neurological or orthopedic conditions; Lower leg instabilities

Description	Male	Female
Right	AF1000-01	AF1000-03
Left	AF1000-02	AF1000-04

## Inline Pneumatic Walker

- Features inflatable air bladder wrap for compression and stabilization
- Heavy duty flared uprights are color-coded for easy size identification
- Low profile rocker sole promotes natural gait
- Wide foot plate for patient comfort and stability
- Treaded sole for patient safety
- Available in open and closed heel design
- Offered in both standard and low top versions

**Indications:** Moderate to severe ankle sprains and strains; Stable foot and/or ankle fractures; Syndesmosis sprain; Post ORIF



Best Choice

Size	Male Shoe Size	Female Shoe Size	Standard Open Heel	Standard Closed Heel	Low Top Open Heel	Low Top Closed Heel
Small	3½ - 7	5 - 8½	15470005	15490005	15480005	15500005
Medium/Large	7½ - 10½	9 - 12	15470067	15490067	15480067	15500067
X-Large	11 +	12½ +	15470008	15490008	15480008	15500008
Pediatric	8 - 12 (Child)		15470002	15490002	15480002	15500002

## Pacesetter™ Air Walker

- Inflatable air bladder design for compression and stabilization
- Unique non-welded air bladder reduces chance of over-inflation to the point of bursting or tearing
- Wide foot plate for patient comfort and stability
- Rocker sole helps promote natural gait
- Heel strap helps prevent foot migration

**Indications:** Moderate to severe ankle sprains and strains; Stable foot and/or ankle fractures; Syndesmosis sprain; Post ORIF



Size	Male Shoe Size	Female Shoe Size	Product #
Small	3 - 6½	4 - 7½	15780005
Medium/Large	7 - 10½	8 - 12½	15780067
X-Large	11 - 16	13 +	15780008

# WALKERS



## Inline Walker

- Heavy duty flared uprights are color-coded for easy size identification
- Low profile rocker bottom promotes natural gait
- Flat, wide foot plate provides patient comfort and stability
- Treaded sole for patient safety
- Available in open or closed heel design
- Offered in both standard and low top versions

**Indications:** Moderate to severe ankle sprains and strains; Stable foot and/or ankle fractures; Syndesmosis sprain; Post ORIF

Best Choice

Size	Male Shoe Size	Female Shoe Size	Standard Open Heel	Standard Closed Heel	Low Top Open Heel	Low Top Closed Heel
Small	3½ - 7	5 - 8½	15520005	15530005	15540005	15550005
Medium/Large	7½ - 10½	9 - 12	15520067	15530067	15540067	15550067
X-Large	11 +	12½ +	15520008	15530008	15540008	15550008
Pediatric	8 - 12 (Child)		15520002	15530002	15540002	15550002



## Pacesetter™ II Walker

- Rocker sole helps promote rolling gait pattern
- Double laminated foam wrap is breathable to prevent skin breakdown
- Wide foot plate for patient comfort and stability
- Available in open or closed heel design
- Offered in both standard and low top versions

**Indications:** Moderate to severe ankle sprains and strains; Stable foot and/or ankle fractures; Syndesmosis sprain; Post ORIF

Size	Male Shoe Size	Female Shoe Size	Standard Open Heel	Standard Closed Heel	Low Top Open Heel
Small	3 - 6½	4 - 7½	CSG53510-A	CSG53550-A	CSG53520-A
Medium/Large	7 - 10½	8 - 12½	CSG53510-B	CSG53550-B	CSG53520-B
X-Large	11 - 16	13 +	CSG53510-C	CSG53550-C	CSG53520-C



## Tracker EX Walker

- Lightweight design for patient comfort
- Low profile rocker bottom sole helps promote natural gait
- Wide foot plate for patient comfort and stability
- Closed heel design helps prevent heel migration
- Extra long wrap provides full coverage
- Available in standard and low top versions

**Indications:** Moderate to severe ankle sprains and strains; Stable foot and/or ankle fractures; Syndesmosis sprain; Post ORIF

Size	Male Shoe Size	Female Shoe Size	Standard	Low Top
X-Small	up to 5½	up to 4½	15590002	15580002
Small	6 - 8½	5 - 8½	15590005	15580005
Medium	9 - 10½	9 - 11½	15590006	15580006
Large	11 - 13	12 - 13½	15590067	15580067
X-Large	13½ +	14 +	15590008	15580008

## Three-D® Lowest Profile Walker

- Extra low heel height helps prevent external hip rotation
- Tri-radiate sole encourages natural gait pattern
- Offset uprights at the malleoli accommodate swelling
- Padded foam inlay for increased energy absorption
- Flat foot bed keeps foot in neutral position
- Criss-cross straps help eliminate forward foot migration

**Indications:** Moderate to severe ankle sprains and strains; Stable foot and/or ankle fractures; Syndesmosis sprain; Post ORIF

Size	Male Shoe Size	Female Shoe Size	Product #
X-Small	up to 5	up to 6	SL5002-01
Small	5 - 7	6 - 8	SL5002-03
Medium	7 - 10	8 - 11	SL5002-05
Large	10 - 12	11 - 13	SL5002-07
X-Large	12 +	13 +	SL5002-09



Best Choice

## Three-D® Low Profile Walker

- Low heel height helps prevent external hip rotation
- Tri-radiate sole encourages natural gait pattern
- Offset uprights at the malleoli accommodate swelling
- Two layer memory foam inlay conforms to foot
- Flat foot bed keeps foot in neutral position
- Criss-cross straps help eliminate forward foot migration
- Available in boot or open plantar wrap design

**Indications:** Moderate to severe ankle sprains and strains; Stable foot and/or ankle fractures; Syndesmosis sprain; Post ORIF

Size	Male Shoe Size	Female Shoe Size	Boot Wrap	Open Wrap
X-Small	up to 5	up to 6	SL2001-01	SL2001-08
Small	5 - 7	6 - 8	SL2001-03	SL2001-10
Medium	7 - 10	8 - 11	SL2001-05	SL2001-12
Large	10 - 12	11 - 13	SL2001-07	SL2001-14
X-Large	12 +	13 +	SL2001-09	SL2001-16



## One-Piece Samson Walker

- Low profile anatomical design promotes natural gait pattern
- Two layer memory foam inlay conforms to foot
- Adjustable straps for even compression
- Wrap material allows skin to breathe and helps to prevent skin breakdown
- Available in standard, mid and low top designs

**Indications:** Moderate to severe ankle sprains and strains; Stable foot and/or ankle fractures; Syndesmosis sprain; Post ORIF

Size	Male Shoe Size	Female Shoe Size	Standard	Mid	Low Top
X-Small	up to 5	up to 6	SL8001-01	—	—
Small	5 - 7	6 - 8	SL8001-03	SL8001-43	SL8001-33
Medium	7 - 10	8 - 11	SL8001-05	SL8001-45	SL8001-35
Large	10 - 12	11 - 13	SL8001-07	SL8001-47	SL8001-37
X-Large	12 +	13 +	SL8001-09	—	—





## Splintsrite™

- Alternative to posterior splinting in the emergency room
- Medial and lateral stays provide stability
- Rocker bottom sole provides natural gait pattern
- Foam wrap provides compression to help reduce edema

**Indications:** Grade II and III ankle sprains; Syndesmosis sprains; Ankle ligament reconstructive procedures

Size	Male Shoe Size	Female Shoe Size	Product #
X-Small	up to 5	up to 6	AB2700-10
Small	5 - 7	6 - 8	AB2700-12
Medium	7 - 10	8 - 11	AB2700-14
Large	10 - 12	11 - 13	AB2700-16
X-Large	12 +	13 +	AB2700-18

# POST-OP SHOES

## Original Darby Shoe

- Semi-rigid metatarsal shank provides controlled forefoot motion
- Reinforced padded heel counter to help reduce slippage
- Padded upper for patient comfort
- Loop-lock closure for easy application
- Universal right/left

**Indications:** Post-operative foot protection and support; Minor foot fractures

Size	Male Shoe Size	Product #	Female Shoe Size	Product #
Small	6 - 8	DS1014-10	4 - 6	DS1013-10
Medium	8½ - 10	DS1014-12	6½ - 8	DS1013-12
Large	10½ - 12	DS1014-14	8½ - 10	DS1013-14
X-Large	12½ - 14	DS1014-16	—	—
Pediatric	12 - 1 (Child)	DS1012-00	—	—



## Aegis Post-Op Shoe

- Semi-rigid rubber sole for support
- Reinforced padded heel counter to help reduce slippage
- Padded upper for patient comfort
- Loop-lock closure for easy application
- Universal right/left

**Indications:** Post-operative foot protection and support; Minor foot fractures

Size	Male Shoe Size	Product #	Female Shoe Size	Product #
Small	6 - 8	2044-01	4 - 6	2044-07
Medium	8½ - 10	2044-02	6½ - 8	2044-08
Large	10½ - 12	2044-03	8½ - 10	2044-09
X-Large	12½ - 14	2044-04	—	—



## Canvas Post-Op Shoe

- Rigid rubber sole provides controlled forefoot motion
- Canvas upper for durability
- Hook and loop closure
- Universal right/left

**Indications:** Post-operative foot protection and support; Minor foot fractures

Size	Male Shoe Size	Product #	Female Shoe Size	Product #
Small	6 - 8	2041-01	4 - 6	2041-07
Medium	8½ - 10	2041-02	6½ - 8	2041-08
Large	10½ - 12	2041-03	8½ - 10	2041-09
Pediatric	6 - 8 (Child)	2041-00	—	—





# POST-OP SHOES



## Deluxe Wooden Post-Op Shoe

- Rigid wood sole with rubber bottom for support
- Perforated suede upper for comfort and durability
- Available in hook and loop or lace up closure design
- Universal right/left

**Indications:** Post-operative foot protection and support; Minor foot fractures

Size	Male Shoe Size	Hook/Loop Closure	Lace Up Closure
Child	2 - 4½	2007-01	2006-01
Small	5 - 6½	2007-02	2006-02
Medium	7 - 8½	2007-03	2006-03
Large	9 - 11	2007-04	2006-04
X-Large	11½ - 13	2007-05	2006-05

Size	Female Shoe Size	Hook/Loop Closure	Lace Up Closure
Child	2 - 3½	2007-06	2006-06
Small	4 - 6	2007-07	2006-07
Medium	6½ - 8	2007-08	2006-08
Large	8½ - 10	2007-09	2006-09



## Canvas Wooden Post-Op Shoe

- Rigid wood sole with rubber bottom for support
- Blue canvas upper
- Available in hook and loop or lace up closure design
- Universal right/left

**Indications:** Post-operative foot protection and support; Minor foot fractures

Size	Male Shoe Size	Hook/Loop Closure	Lace Up Closure
Small	5 - 6½	2020-11	2020-01
Medium	7 - 8½	2020-12	2020-02
Large	9 - 11	2020-13	2020-03
X-Large	11½ - 13	2020-14	2020-04

Size	Female Shoe Size	Hook/Loop Closure	Lace Up Closure
Small	4 - 6	2021-11	2021-01
Medium	6½ - 8	2021-12	2021-02
Large	8½ - 10	2021-13	2021-03
X-Large	10½ - 12	2021-14	2021-04

## Healing Shoe

- Non-slip rubber sole for patient safety
- Durable vinyl upper construction
- Padded liner and foam insole for patient comfort
- Features heel strap and lace-up design for a secure fit
- Universal right/left

**Indications:** Post-operative foot protection and support; Minor foot fractures

Size	Male Shoe Size	Female Shoe Size	Product #
X-Small	—	4 - 6½	HS1003-08
Small	5 - 6½	7 - 9	HS1003-10
Medium	7 - 8	9½ - 11½	HS1003-12
Large	8½ - 11	12 - 13	HS1003-14
X-Large	11½ +	—	HS1003-16



## Darco® All Purpose Boot

- Closed toe design provides all-season protection
- High top design and ankle strap help eliminate heel slippage
- Hook and loop closure for easy application
- Universal right/left

**Indications:** Post-operative foot protection and support; Minor foot fractures

Size	Male Shoe Size	Female Shoe Size	Product #
X-Small	3 - 4½	3 - 6½	2045-01
Small	5 - 6½	7 - 9½	2045-02
Medium	7 - 8½	10 - 12½	2045-03
Large	9 - 11	13 +	2045-04
X-Large	11½ - 14	—	2045-05



## Darco® Med-Surg™ Shoe

- MetaShank™ protection provides rigid control under metatarsal heads
- Strapless forefoot closure expands to accommodate bulky dressings
- Square toe design serves as a bumper for additional protection
- Rocker sole reduces plantar pressure on forefoot and heel
- Ankle strap helps eliminate heel migration
- Universal right/left

**Indications:** Post-operative foot protection and support; Minor foot fractures

Size	Male Shoe Size	Product #	Female Shoe Size	Product #
Small	6 - 8	DS1016-10	4 - 6	DS1015-10
Medium	8½ - 10	DS1016-12	6½ - 8	DS1015-12
Large	10½ - 12	DS1016-14	8½ - 10	DS1015-14
X-Large	12½ - 14	DS1016-16	—	—



# POST-OP SHOES



## Darco® Orthowedge™ Shoe

- Reduces weight bearing pressure from the forefoot with 15° dorsi-angle
- Removable 10 mm insole allows customization for offloading specific areas
- Square toe design serves as a bumper for additional protection
- Hook and loop closure for easy application
- Universal right/left

**Indications:** Post-operative foot protection and support; Minor foot fractures

Size	Male Shoe Size	Female Shoe Size	Product #
X-Small	—	4 - 7	HS1005-08
Small	6 - 8	7½ - 10	HS1005-10
Medium	8½ - 10	10½ - 13	HS1005-12
Large	10½ - 12	13½ +	HS1005-14
X-Large	12½ - 14	—	HS1005-16

## Comfo-Eze Rocker Cast Shoe

- Constructed of durable navy canvas material
- Non-slip rubber sole for patient safety
- Rocker bottom promotes rolling gait pattern
- Heel strap helps prevent foot migration
- Hook and loop closure for easy application
- Available in open and closed toe design

**Indications:** Lower extremity cast protection

Size	Length	Width	Open Toe	Closed Toe
Pediatric	5½"	2½"	2004-00	—
X-Small	7"	3½"	2004-01	2038-01
Small	8½"	4½"	2004-02	2038-02
Medium	10"	4¾"	2004-03	2038-03
Large	11½"	5¼"	2004-04	2038-04
X-Large	13"	5½"	2004-05	2038-05



## Canvas Cast Shoe

- Constructed of navy canvas material
- Padded canvas footbed for maximum comfort
- Non-slip rubber sole for patient safety
- Rocker bottom promotes rolling gait pattern
- Heel strap helps prevent foot migration
- Hook and loop closure for easy application

**Indications:** Lower extremity cast protection

Size	Length	Width	Product #
Child	5½"	2½"	2035-01
Small	8½"	4½"	2035-02
Medium	10"	4¾"	2035-03
Large	11½"	5¼"	2035-04
X-Large	13"	5½"	2035-05



## Black Canvas Cast Shoe

- Constructed of durable black canvas material
- Non-slip rubber sole for patient safety
- Rocker bottom promotes rolling gait pattern
- Heel strap helps prevent foot migration
- Hook and loop closure for easy application

**Indications:** Lower extremity cast protection

Size	Length	Width	Product #
X-Small	7"	3½"	2033-01
Small	8½"	4½"	2033-02
Medium	11½"	5¼"	2033-03
Large	13"	5½"	2033-04



# CAST SHOES



## Vinyl Cast Boot

- Constructed of durable vinyl material
- Rocker bottom promotes rolling gait pattern
- Closed heel design to prevent foot migration
- Open or closed toe design
- Hook and loop closure for easy application

**Indications:** Lower extremity cast protection

Size	Length	Width	Open Toe	Closed Toe
X-Small	7"	3½"	2027-01	2037-01
Small	8½"	4½"	2027-02	2037-02
Medium	10"	4¾"	2027-03	2037-03
Large	11½"	5¼"	2027-04	2037-04
X-Large	13"	5½"	2027-05	2037-05



## E-Z Walk Cast Boot

- Constructed of durable vinyl material
- Foam foot bed helps provide shock absorption during ambulation
- Closed heel design helps prevent foot migration
- Open or closed toe design
- Hook and loop closure for easy application

**Indications:** Lower extremity cast protection

Size	Length	Width	Open Toe	Closed Toe
Pediatric	5½"	2½"	2001-00	—
Child	7"	3½"	2001-01	2005-01
Small	8½"	4½"	2001-02	2005-02
Medium	10"	4¾"	2001-03	2005-03
Large	11½"	5¼"	2001-04	2005-04
X-Large	13"	5½"	2001-05	2005-05



## Rainbow Cast Sandal

- Features low profile sandal design for patient comfort
- Rocker sole promotes rolling gait pattern
- Heel strap helps prevent foot migration
- Elasticized straps with hook and loop closure for easy application

**Indications:** Lower extremity cast protection

Size	Length	Width	Product #
X-Small	7"	3½"	2032-00
Small	8½"	4½"	2032-01
Medium	10"	4¾"	2032-02
Large	13"	5½"	2032-03

## Darco® SlimLine Cast Boot

- Square toe design serves as a bumper for protection and provides room for larger casts
- Stretchable upper material helps boot conform to cast allowing precise fit
- Rocker sole promotes rolling gait pattern
- Padded insole to cushion patient's foot
- Strapless forefoot closure adjusts easily to cast size

**Indications:** Lower extremity cast protection



Size	Male Shoe Size	Female Shoe Size	Product #
X-Small	—	3 - 5	2039-01
Small	—	5½ - 8	2039-02
Medium	6 - 8	8½ - 10	2039-03
Large	8½ - 11	10½ - 13	2039-04
X-Large	11½ - 14	13+	2039-05

Size	Child Shoe Size	Product #
Child/Small	8 - 10½	2039-98
Child/Medium	11 - 1½	2039-99
Child/Large	2 - 4	2039-00

# FOOT & HEEL



## DeRoyal® Night Splint

- Dorsiflexion design applies stretch for effective pain relief
- Lightweight and low profile design for patient comfort
- Adjustable straps to vary amount of stretch
- Removable foam liner for laundering
- Optional foam wedges provided for additional stretch
- Kit features: Night splint, ice pack, exercise band, heel cups, foot control strap and self-management guide

**Indications:** Plantar fasciitis; Achilles tendonitis; Tendonosis

Size	Male Shoe Size	Female Shoe Size	Splint	Kit
Small	4 - 6	5 - 7	NS2000-10	SC2000-10
Medium	6½ - 9	7½ - 10	NS2000-12	SC2000-12
Large	9½ +	10½ +	NS2000-14	SC2000-14



## Silicone Heel Cups

- Constructed of medical grade silicone for long-lasting support and comfort
- Soft blue spot helps relieve heel strike fatigue and discomfort through shock absorption
- Packaged 1/Pair

**Indications:** Heel spurs; Bruises; Callosities; Plantar fasciitis and apophyseal injuries; Retrocalcaneal bursitis; Fat pad atrophy

Size	Product #
Small	HC2000-10
Medium	HC2000-12
Large	HC2000-14



## Grey Spot Heel Cups

- Constructed of Visco™ elastic polymers to dissipate heel strike energy
- "U-shaped" design helps divert weight bearing forces away from plantar fascia
- Packaged 1/Pair

**Indications:** Heel spurs; Bruises; Callosities; Plantar fasciitis and apophyseal injuries; Retrocalcaneal bursitis; Fat pad atrophy

Size	Product #
Regular	HC1000-01
Large	HC1000-02



# PRESSURE REDUCTION

## Ankle Contracture Boot

- Static ankle foot orthosis used for positioning and pressure reduction
- Available in fleece or vel-foam material with or without boot sole
- Softgoods completely cover splint leaving no exposed metal against skin
- "No strap" closures make it easy for patients to apply and remove the boot
- Anti-rotation bar, padded toe piece and mesh laundry bag included
- Accessories include hip abduction bar, boot sole and extra softgoods
- Universal right/left

**Indications:** Hip rotation; Foot drop; Not recommended for severe plantar flexion contractures, open wounds, severe spasticity or thrombosis

Description	Vel-Foam	Fleece
Ankle Contracture Boot	4302	4307
Ankle Contracture Boot w/ Boot Sole	4301	4306
Ankle Contracture Boot w/ Boot Sole & Extra Liner	4300	4305

Accessories	Product #
Extra Liner - Vel-Foam	4300K
Extra Liner - Fleece	4305K
Hip Abduction Bar	4300/05AB
Boot Sole	4300/05WK

### Sizing Guide

Add letter to end of product # to indicate size

Measure calf circumference 8" up from bottom of heel

- A\* = 8" - 12"
- B = 10" - 14"
- C = 14" - 16"
- D = 16" - 21"
- E = 16" - 21"

Measure foot circumference around the widest part of foot

- A\* = 6" - 8"
- B = 7" - 9"
- C = 9" - 11"
- D = 11" - 13"
- E = 13" +

\*Size A available in 4301 style only



Vel-foam

Best  
Choice



Fleece



# PRESSURE REDUCTION



## Bunny Boot

- Breathable polyfiber liner allows air circulation for increased comfort
- Adjustable foam wedge at base of frame for additional comfort
- Rigid plastic frame holds foot at a 90° angle and can act as a foot rotator with internal or external rotation of 10° or 20° arcs
- Hook and loop/d-ring closures

**Indications:** Toe and heel pressure reduction; Foot drop; Static stretching of foot

Size	Qty	Product #
Universal	1 Each	M0501D01
Universal	24/Case	M0501D05



## Heel Pillow

- Soft, cotton padding provides maximum heel cushioning for protection against skin breakdown
- Hook and loop closure for easy application and adjustment

**Indications:** Equalizes pressure on heel to help prevent decubitus formation; Prevent skin shearing; Prevent skin friction

Size	Qty	Product #
Universal	1 Pair	M3040



## DeRoyal® Heel Suspension Boot

- Smooth or convoluted foam construction
- Posterior panel helps reduce friction against bed surface
- Ventilation holes help promote air circulation and dissipate heat
- Sloped heel suspension wedge pre-attached or free for precise placement
- Extra pad included for additional elevation and rotation protection
- Hook/d-ring strap closures

**Indications:** Pressure reduction of heel; Pressure protection; Hip rotation; Foot drop; Hip fractures (pre and post operative); Heel decubitus treatment and prevention

Description	Qty	Product #
Convoluted	6/Case	M60-100
Smooth	6/Case	M60-200
Smooth w/Attached Wedge	6/Case	M60-300

Best Choice

# PRESSURE REDUCTION

## Foam Heel/Ankle Protector

- Breathable convoluted foam construction for patient comfort
- Available in high or low-density foam
- Protects heel and ankle area from decubitus formation
- Hook and loop closure for easy application and adjustment

**Indications:** Equalizes pressure on heel to help prevent decubitus formation; Prevent skin shearing; Prevent skin friction

Description	Qty	Product #
High-density	1 Pair	M190FF
High-density	24 Pair/Case	M190FF-24
Low-density	20 Pair/Case	M190FFE



## Heel & Elbow Protector - Full Foot Design

- Available in breathable, deluxe or standard polyfiber material
- Full foot design
- Ventilation holes help promote air circulation
- Hook and loop closure for easy application and adjustment

**Indications:** Equalizes pressure on heel and/or elbow to help prevent decubitus formation; Prevent skin shearing; Prevent skin friction

Deluxe	Size	Qty	Product #
Heel	Universal	1 Pair	M3014D
Heel	Universal	6 Pair/Case	M3014D-6

Standard	Size	Qty	Product #
Heel	Universal	1 Pair	M3015
Heel	Universal	6 Pair/Case	M3015B
Elbow	Universal	1 Pair	M3006



## Gelbo™ Heel & Elbow Protector

- Nylon/spandex sleeve conforms around the heel or elbow to hold product in place
- Lined with natural mineral oil-based gel formula for shock absorption and cushioning
- Moisturizes dry, cracked skin at heel and/or elbow
- Available in standard and large sizes

**Indications:** Equalizes pressure on heel and/or elbow to help prevent decubitus formation; Prevent skin shearing; Prevent skin friction; Shock absorption

Size	Qty	Product #
Standard (>110 lbs)	1 Pair	M3022
Large (<110 lbs)	1 Pair	M3022-L



# PRESSURE REDUCTION



## Padded Heel & Elbow Protector

- Knitted sleeve conforms around heel or elbow to help hold product in place
- Breathable foam or polyfiber pad sewn into cushioned sock for maximum pressure protection with no irritating seam

**Indications:** Equalizes pressure on heel and/or elbow to help prevent decubitus formation; Prevent skin shearing; Prevent skin friction; Shock absorption

Size	Qty	Foam Pad	Polyfiber Pad
Universal	12 Pair/Case	M3000U	M3000KU
Universal	36 Pair/Case	M3001U	M3001KU

# FOAM POSITIONERS


## Concave Hip Abduction Pillow - High Density

- Breathable, high-density foam construction
- Concave sides hold legs in secure position
- Long foam straps with hook closure for easy application and adjustment
- Available in compressed packaging for efficient storage or non-compressed packaging

**Indications:** Post surgical hip abduction

Size	Dimensions				Qty	Compressed
	A	B	C	D		
Infant	4½"	5"	2"	2"	10/Case	M60-076-I
X-Small	9"	10"	3"	4"	5/Case	M60-076-XS
Small	14"	18"	7"	5"	4/Case	M60-076-S
Medium	16½"	21"	7½"	6"	4/Case	M60-076-M
Large	18"	24"	9"	6"	4/Case	M60-076-L

Size	Dimensions				Qty	Non-Compressed
	A	B	C	D		
Small	14"	18"	7"	5"	2/Case	M60-078-S
Medium	16½"	21"	7½"	6"	2/Case	M60-078-M
Large	18"	24"	9"	6"	2/Case	M60-078-L




## Concave Hip Abduction Pillow - Wide


- Breathable, high-density foam construction
- Concave sides hold legs in secure position
- Wider angle for increased abduction
- Long foam straps with hook closure for easy application and adjustment
- Compressed packaging for efficient storage

**Indications:** Post surgical hip abduction

Size	Dimensions				Qty	Product #
	A	B	C	D		
Small	16"	16"	10"	7"	5/Case	M60-024-S
Medium	18"	20"	10"	9"	5/Case	M60-024-M
Large	20"	24"	10"	9"	5/Case	M60-024-L



# FOAM POSITIONERS




## Concave Hip Abduction Pillow - Wide w/ Liner

- Breathable, high-density foam construction
- Concave sides hold legs in secure position
- Wider angle for increased abduction
- Polyfiber liner provides additional patient comfort
- Long foam straps with hook closure for easy application and adjustment
- Compressed packaging for efficient storage

**Indications:** Post surgical hip abduction

Size	Dimensions				Qty	Product #
	A	B	C	D		
Small	16"	16"	10"	7"	5/Case	M60-025-S
Medium	18"	20"	10"	9"	5/Case	M60-025-M
Large	20"	24"	10"	9"	5/Case	M60-025-L




## Flat Hip Abduction Pillow - Narrow

- Breathable, high-density foam construction
- Long foam straps with hook closure for easy application and adjustment
- Compressed packaging for efficient storage

**Indications:** Post surgical hip abduction

Size	Dimensions				Qty	Product #
	A	B	C	D		
Small	10"	16"	6"	7"	5/Case	M60-031-S
Medium	12"	21"	6"	7"	5/Case	M60-031-M
Large	14"	23"	6"	7"	5/Case	M60-031-L



## Hip Abduction Pillow - Adjustable

- Breathable, high-density foam construction
- 4" tear-off blocks on both ends for length adjustment
- Long, perforated foam straps with hook closure for easy application and adjustment
- Compressed packaging for efficient storage

**Indications:** Post surgical hip abduction

Size	Dimensions				Qty	Product #
	A	B	C	D		
Universal	10"	18"-27"	6"	5"	4/Case	M60-224


# FOAM POSITIONERS

## Arm Cradle

- Breathable, high-density foam construction
- Ventilation holes promote air circulation and dissipate heat
- Compressed packaging for efficient storage

**Indications:** Pressure reduction; Arm elevation

Size	Dimensions			Qty	Product #
	A	B	C		
Universal	12½"	9"	15½"	5/Case	M60-034




## Foot Cradle

- Breathable, high-density foam construction
- Ventilation holes promote air circulation and dissipate heat
- Compressed packaging for efficient storage

**Indications:** Pressure reduction; Foot drop support; Heel elevation

Size	Dimensions		Qty	Product #
	A	B		
Universal	14"	9"	10/Case	M60-045




## Foot Drop Support

- Breathable, high-density foam construction
- Adjustable anti-rotation strap holds foot in neutral position
- Hook receivable foam strap
- Compressed packaging for efficient storage

**Indications:** Pressure reduction; Foot drop support; Heel elevation


Size	Dimensions		Qty	Product #
	A	B		
Universal	8"	10"	5/Case	M60-018





# FOAM POSITIONERS

## Leg Positioner




- Breathable, high-density standard or flat foam construction
- Available with or without 45° angle
- Stabilizes knee and lower leg while promoting better circulation
- Minimizes level of immobilization
- Foam straps with hook and loop closure for easy application and adjustment
- Compressed packaging for efficient storage

**Indications:** Pressure reduction of heel, ankle, mid-foot and mid-calf; Pressure protection; Static back pain relief

Size	Dimensions			Qty	Product #
	A	B	C		
Standard	11"	29"	11"	2/Case	M60-015A
Standard w/ 45° angle	11"	28"	10"	1 Each	M60-015B
Flat w/ 45° angle	8	31"	15"	2/Case	M60-015CD

## Leg Positioner w/ Foot Drop Strap




- Breathable, high-density foam construction
- Adjustable foot drop strap holds foot in neutral position
- Stabilizes knee and lower leg while promoting better circulation
- Minimizes level of immobilization
- Foam straps with hook and loop closure for easy application and adjustment
- Compressed packaging for efficient storage

**Indications:** Pressure reduction of heel, ankle, mid-foot and mid-calf; Pressure protection; Static back pain relief

Size	Dimensions			Qty	Product #
	A	B	C		
Universal	15"	32"	11"	4/Case	M60-015

## Body Aligner



- Breathable, high-density foam construction
- Non-skid surface maintains proper positioning while reducing pressure
- Available in compressed packaging for efficient storage or non-compressed packaging

**Indications:** Pressure reduction; Positioning

Size	Dimensions			Qty	Non Compressed	
	A	B	C		Compressed	Non Compressed
Universal	24"	8"	11"	5/Case	M60-032	—
Universal	24"	8"	11"	10/Case	M60-032B	—
Universal	26"	27"	12"	2/Case	M60-033	—
Universal	24"	8"	11"	2/Case	—	M60-032NC

# HOT & COLD THERAPY

## DeRoyal® Hot/Cold Therapy Unit

- Automatic temperature controlled unit provides hours of hot or cold therapy
- Digital temperature read-out
- Temperature control reduces risk of nerve damage or frostbite
- Measures, adjusts and regulates motor speed for precision temperature control without manual adjustment
- Portable system suitable for any medical setting and continuous treatment at home
- Lightweight, low-profile design is patient friendly
- Automatic shut-off feature lets user know when to add ice
- Hot therapy range (90° - 120°F); cold therapy range (45° - 60°F)
- Hot/cold therapy blankets available in both sterile and non-sterile versions

**Indications:** Pain and soft tissue swelling as a result of surgery or trauma

**Description**  
Hot/Cold Therapy Unit

**Product #**  
T600



## DeRoyal® Cold Therapy Unit

- Automatic temperature controlled unit provides hours of cold therapy
- Temperature control reduces risk of nerve damage or frostbite
- Measures, adjusts and regulates motor speed for precision temperature control without manual adjustment
- Portable system suitable for any medical setting and continuous treatment at home
- Lightweight, low-profile design is patient friendly
- Automatic shut-off feature lets user know when to add ice
- Cold therapy range (45° - 60°F)
- Cold therapy blankets available in sterile or non-sterile versions
- Unit available in combo set with most popular cold therapy blankets

**Indications:** Pain and soft tissue swelling as a result of surgery or trauma

**Description**  
Cold Therapy Unit

**Product #**  
T505



## DeRoyal® Cold Therapy Unit Combos

- T505 cold therapy unit packaged with hot/cold therapy blanket
- Available in non-sterile or sterile versions

**Description**

Chest/Lumbar  
Knee/Shoulder  
Knee/Shoulder (w/ straps)  
Knee (w/ attached straps)  
Shoulder (w/ attached straps)  
Foot/Ankle

**Non-Sterile**

T520NS  
T530NS  
T531NS  
T532NS  
T533NS  
T580NS

**Sterile**

T520ST  
T530ST  
T531ST  
T532ST  
T533ST  
T580ST



# HOT & COLD THERAPY



## Hot/Cold Therapy Blankets

- Blankets feature unique channel construction to move water directly to the injury site
- Available in non-sterile or sterile versions
- Universal right/left

**Indications:** Pain and soft tissue swelling as a result of surgery or trauma



**Description**  
Arm/Leg

**Measurement**  
10" x 5"

**Qty**  
1 Each

**Non-Sterile**  
T5010NS

**Sterile**  
T5010ST

Lumbar

10½" x 11½"

1 Each

T5020NS

T5020ST

Knee/Shoulder

10½" x 10½"

1 Each  
10/Case

T5030NS  
T5030NS-10

T5030ST  
T5030ST-10

Knee/Shoulder - Full

10½" x 10½"

1 Each  
10/Case

T5031NS  
T5031NS-10

T5031ST  
T5031ST-10

Knee - Full, Attached Straps

10½" x 10½"

1 Each  
10/Case

T5032NS  
T5032NS-10

T5032ST  
T5032ST-10

Shoulder - Full, Attached Straps

10½" x 10½"

1 Each  
10/Case

T5033NS  
T5033NS-10

T5033ST  
T5033ST-10

Hand/Wrist - Double

3½" x 5"

2 Each

T5040NS

T5040ST

Hand/Wrist

3½" x 5"

1 Each

T5050NS

T5050ST

Hip

15½" x 18"

1 Each

T5060NS

T5060ST

Cervical

10" x 21½"

1 Each

T5070NS

T5070ST

Foot/Ankle

10" x 11½"

1 Each

T5080NS

T5080ST

Facial

5½" x 3¾"

1 Each

T5090NS

T5090ST



## Hot/Cold Therapy Wraps

- Wrap keeps blanket well insulated and provides a simple way to ensure precise blanket placement

**Description**

**Measurement**

**Product #**

Lumbar

12" x 58"

T6100

Knee

15" x 27"

T6110

Shoulder

10½" x 10½"

T6120

Hand

6" x 7¼"

T6130

Facial

5" x 30"

T6140

Hip

17" x 58"

T6160

# HOT & COLD THERAPY

## Foam Hot/Cold Therapy Wraps

- Laminated foam construction keeps hot/cold pack insulated
- Anatomically designed to hold hot/cold pack in place
- Foam is compressive and can be tightened to apply even compression and help reduce swelling
- Available with or without hot/cold packs
- Universal right/left

**Indications:** Pain and soft tissue swelling as a result of surgery or trauma

Description	Measurement	Packs	Product #
Shoulder	47" x 13"	0	9381-50
		2	9381-00
		4	9381-10
Knee	25½" x 14"	0	9383-50
		2	9383-00
		4	9383-10
Knee w/Straps	30" x 16"	0	9387-50
		4	9387-00
Back	58" x 10"	0	9384-50
		1	9384-00
Ankle	14½" x 5"	1	9385-00
Wrist	14" x 5¼"	1	9386-00



# HOT & COLD THERAPY



## Hot/Cold Pack

- Reusable for both hot or cold therapy
- Packaged each or 20/Case

**Indications:** Pain and soft tissue swelling as a result of surgery or trauma

Measurement	Each	Case
4½" x 10"	9380-00	9380-20



## Chill-n-Wrap Ice Bag

- Constructed of soft outer cover with moisture absorbing middle layer
- Leak-proof plaster bladder prevents patient exposure
- Plastic closure clip is completely enclosed
- Conformable opening holds shape while ice is being added
- Packaged 10/Case

**Indications:** Pain and soft tissue swelling as a result of surgery or trauma

Size	Pouch Dimensions	Strap Length	Product #
Small	6" x 13"	16½"	8726-40
Large	6" x 16"	16½"	8726-50

## Surgical Bra

- Constructed of Lycra® material
- Front hook and loop closure for easy access to dressings and application
- Wide adjustable shoulder straps for maximum patient comfort

**Indications:** Post operative compression; Dressing holder

Size	Chest Measurement	Product #
X-Small	28" - 31"	M5001-XS
Small	32" - 35"	M5001-S
Medium	36" - 39"	M5001-M
Large	40" - 43"	M5001-L
X-Large	44" - 47"	M5001-XL
XX-Large	48" - 51"	M5001-XXL
XXX-Large	52" - 55"	M5001-XXXL



## Trach Tube Holder

- Soft, foam one-piece holder secures most tracheostomy tubes
- Available with wide or narrow fasteners for secure fit (adult only)
- Adjustable to accommodate different neck circumferences
- Packaged 12/Case

**Indications:** Patients with tracheostomy

Size	Neck Circumference	Product #
Adult	up to 20"	M1151
Pediatric	up to 10"	M1152
Bariatric	up to 30"	M1153
Adult, Wide	up to 23"	M1150



Narrow



Wide

## Catheter Strap

- Soft elastic or foam construction for patient comfort
- Excellent substitute to messy, skin-irritating tape
- Hook and loop fasteners for easy application and placement
- Available in standard or extra long length
- Ideal for any area where foley catheters are applied

**Indications:** Securing foley catheters; Catheter positioning

Size	Dimensions	Qty	Elastic	Foam
Universal	22"L x 2"W	36/Case	M1133	—
Universal, X-Long	38"L x 2"W	24/Case	M1133-38	—
Universal	28"L x 1"W	50/Case	—	M1128



Elastic



Foam





## Flexible I.V. Arm Board

- Malleable board for custom positioning with soft foam cover
- Adjustable foam straps eliminate the need for skin-irritating tape
- Straps attach to board to keep arm immobilized during I.V. therapy
- Available in adult, child and infant sizes

**Indications:** Positioning aid for I.V. therapy; Arterial support

Size	Dimensions	Qty	Product #
Adult	9"L x 4"W	25/Case	M8125-A
Adult, X-Long	18"L x 4"W	25/Case	M8125-18
Preemie	2"L x 1"W	30/Case	M8125-P
Infant	4"L x 1"W	30/Case	M8125-I
Child	6"L x 2"W	30/Case	M8125-C



## Elbow Immobilizer

- Soft, perforated foam construction for patient comfort and breathability
- Metal t-bar stays for optimal immobilization and strength
- Immobilizes the elbow joint to prevent pulling NG tubes, tracheostomy tubes or I.V. lines
- Hook and loop and elastic closures for fast and accurate placement

**Indications:** Limits arm movement; Limb positioning

Size	Dimensions	Qty	Product #
Adult	14"L x 17"W	1 Pair	M2035
Adult	14"L x 17"W	50 Eaches	M2035B
Infant	7½"L x 12½"W	1 Pair	M7035-S
Pediatric	7½"L x 17½"W	1 Pair	M7035-L



## Cardiac Cushion

- Heart diagram printed on pillow with permanent marker included for bedside education purposes
- Cotton cover with polyester fill provides superior cushioning
- Packaged 10/Case

**Indications:** Post-cardiac surgery; Support during expectoration; Sternal support to aid in healing; Helps reduce pain when coughing and deep breathing; Marketing tool\*

Size	Product #
Adult	M13-100
Pediatric	M13-500

\*Contact your local DeRoyal sales representative for facility's custom logo information.



## Disposable Telemetry Pouch

- Durable, non-woven material construction
- Neck and torso tie straps for easy application and adjustment
- Open or closed pouch design
- Packaged 50/Case

**Indications:** Patients with telemetry units

Dimensions	Closed Pouch	Open Pouch
7½"L x 6½"W	M4210D06	M4209D06



## Polyester Telemetry Pouch

- Durable polyester pouch with foam neck strap for patient comfort
- Hook receivable material for easy application and adjustability
- Packaged 50/Case

**Indications:** Patients with telemetry units

Dimensions	Product #
6½"L x 5"W	M13-055



## Bed Rail Covers

- Foam padded, vinyl covered protectors help prevent patient injury
- Vinyl covering wipes clean for easy disinfecting to help prevent cross contamination issues
- Hook and loop closure for easy application and removal

**Indications:** Helps prevent side rail injury and bed exit

Description	Dimensions	Product #
2 Piece Set	69"L x 12"H	M3000D08
4 Piece Set	69"L x 12"H & 32"L x 12"H	M3002D08





## Cervical Tension Pillow

- Constructed of soft foam with cotton cover
- Designed to help relieve pressure and tension on cervical spine
- Anatomically shaped to fit back of neck and head area

**Indications:** Cervical spine positioning and support; Neck tension relief

**Size**  
Universal

**Product #**  
9222-00



## Cerv-O-Rest Pillow

- Constructed of soft collapsible foam with cotton cover
- Compresses to fill in space behind neck when in a supine position

**Indications:** Cervical spine positioning and support

**Size**  
Universal

**Product #**  
9221-00



## Orthopedic Felt

- Ideal for extra padding

**Measurement**  
36" x 21" x 1/4"


**Product #**  
8829-01

## Invalid Ring

- Soft foam construction provides pressure-free comfort while sitting
- Available with or without cloth cover

**Indications:** Post-natal discomfort; Post-operative pain; Hemorrhoids; Bruised coccyx

Size	Dimensions			w/ Cover	w/o Cover
	A	B	C		
Universal	14½"	12"	2½"	M10-070-C	M10-070



# BANDAGES & WRAPS



## Elastic Bandage

- Durable elastic construction
- Available in standard length (5 yds) and double length (10 yds)
- Packaged individually
- Non-sterile
- Available in Velcro® or clip closure
- Premium option features hook closure on both ends of roll

Size	Quantity	Premium Clip Closure	Premium Velcro® Closure	Econo Clip Closure	Econo Velcro® Closure
2" x 5 yds	10 rolls/case	EB2	EBV2	DEB2	DEBV2
3" x 5 yds	10 rolls/case	EB3	EBV3	DEB3	DEBV3
4" x 5 yds	10 rolls/case	EB4	EBV4	DEB4	DEBV4
6" x 5 yds	10 rolls/case	EB6	EBV6	DEB6	DEBV6
4" x 10 yds	6 rolls/case	HEB4	—	—	—
6" x 10 yds	6 rolls/case	HEB6	—	—	—

# AMBULATORY AIDS

## Crutches

- Lightweight, sturdy aluminum construction
- Push-button height adjustment
- Rubber tips for maximum surface traction
- Padded underarm rest and hand grips for patient comfort
- Packaged 10 Pair/Case

Size	Height	Product #
Short	4'8" - 5'2"	CR1000
Regular	5'2" - 5'10"	CR2000
Tall	5'10" - 6'6"	CR3000



## Folding Walkers

- Lightweight, sturdy aluminum construction
- Push-button height adjustment
- Push-button folding for easy storage and transport
- Rubber tips for maximum surface traction
- Padded hand grips for patient comfort
- Available with or without 3" wheels
- Packaged 2/Case

Size	w/o Wheels	w/ Wheels
Universal	WLK02	WLK02-W



## Walking Canes

- Lightweight and sturdy aluminum construction
- Push-button height adjustable from 30" - 40"
- Rubber tip for maximum surface traction
- Padded hand grip for patient comfort
- Packaged 20/Case

Size	Product #
Universal	CA298



## Pediatric Size Products

Product	Page #
Comfo-Eze Cervical Collar.....	2
Philadelphia® Cervical Collar.....	5
Definitive Collar.....	5
Philly EMT Select Extrication Collar.....	6
Philly One-Piece Extrication Collar.....	6
STAT® Neck Extrication Collar.....	7
XTW Extended Wear Collar.....	8
Premium Clavicle Strap.....	13
Heavy-Duty Clavicle Strap.....	13
3-Way Clavicle Strap.....	13
Deluxe Sling & Swathe.....	15
Universal Sling Style Shoulder Immobilizer.....	18
Premium Arm Sling.....	21
Specialty Arm Sling.....	21
Universal Arm Sling.....	21
Envelope Arm Sling.....	22
Light Blue Arm Sling.....	22
Narrow Pouch Arm Sling.....	22
Deep Pocket Arm Sling.....	23
Strap Arm Sling.....	23
Suede Leatherette Wrist & Wrist/Forearm Splint.....	30
Leatherette Wrist & Wrist/Forearm Splint.....	30
Canvas Wrist & Wrist/Forearm Splint.....	30
Black Foam Wrist & Wrist/Forearm Splint.....	31
Tietex Wrist & Wrist/Forearm Splint.....	31
Black Wrist Splint.....	33
Elastic Wrist Splint.....	35
Finger Dressing/Finger Cage.....	42
Wrist & Forearm Splint.....	44
Cock Up Hand Splint.....	44
Colles' Splint.....	44
Plastic Colles' Splint.....	45
Hand Grip Colles' Splint.....	45
Short Metacarpal Splint.....	45
Long Metacarpal Splint.....	46
Full Hand Splint.....	46
Mason Allen Hand Splint.....	46
Elbow Splint.....	47
Posterior Knee Splint.....	47
Posterior Leg Splint.....	47
Sized Black Foam Knee Immobilizer.....	63
Universal Knee Support.....	66
Full Foot Bucks Traction Boot.....	71
Sized Bucks Traction Boot.....	72
Foam/Canvas Bucks Traction Boot.....	72
Air Ankle Stirrup.....	74
R.I.C.E. Ankle Splint.....	76
Inline Pneumatic Walker.....	79
Inline Walker.....	80
Original Darby Shoe.....	83
Canvas Post-Op Shoe.....	83
Deluxe Wooden Post-Op Shoe.....	84
Comfo-Eze Rocker Shoe.....	87
Canvas Cast Shoe.....	87
E-Z Walk Cast Boot.....	88
Darco® SlimLine Cast Boot.....	89
Ankle Contracture Boot.....	91
Concave Hip Abduction Pillow.....	95
Trach Tube Holder.....	103
Flexible I.V. Arm Board.....	104
Elbow Immobilizer.....	104
Cardiac Cushion.....	104

## ™ & ® Products

Confor® is a registered trademark of DeRoyal, Inc.  
Darco® is a registered trademark of Darco International, Inc.  
DeRoyal® is a registered trademark of DeRoyal, Inc.  
Gelbo™ is a registered trademark of Knit-Rite, Inc.  
Lycra® is a registered trademark of E.I. du Pont de Nemours  
Med-Surg™ is a trademark of Darco International, Inc.  
Metashank™ is a trademark of Darco International, Inc.  
Pacesetter™ is a trademark of DeRoyal, Inc.  
Philadelphia® is a registered trademark of Philadelphia Cervical Collar Company  
Plastazote® is a registered trademark of Zotefoams PLC  
PneuGel™ is a trademark of Sroufe Healthcare Products, Inc.  
Prolign® is a registered trademark of DeRoyal, Inc.  
Orthowedge™ is a trademark of Darco International, Inc.  
Solace™ is a trademark of DeRoyal, Inc.  
Splintsrite™ is a trademark of DeRoyal, Inc.  
STAT® is a registered trademark of DeRoyal, Inc.  
Three-D® is a registered trademark of DeRoyal, Inc.  
Ultralign® is a registered trademark of DeRoyal, Inc.  
Velcro® is a registered trademark of Velcro Industries B.V.  
Versa-Pac™ is a trademark of Tyco Healthcare Group  
Vel-stretch® is a registered trademark of Velcro Industries B.V.

## Patents

Ultralign® – #5,853,378; #5,967,998 and #6,478,759  
Prolign® – #6,666,838  
Ankle Contracture Boot - #6,428,493 and #D435294



# INDEX BY PRODUCT NUMBER

1000159.....	16	1018-03 .....	4	1043-02 .....	7	1056-05 .....	1
1000251.....	1	1018-04 .....	4	1043-03 .....	7	1057-00 .....	1
1000275.....	33	1018-05 .....	4	1043-04 .....	7	1057-01 .....	1
1000285.....	33	1018-10 .....	4	1044-11 .....	5	1057-02 .....	1
1000301.....	1	1018-11 .....	4	1044-12 .....	5	1057-03 .....	1
1000351.....	1	1018-12 .....	4	1044-13 .....	5	1057-04 .....	1
1000401.....	1	1020025.....	30	1044-14 .....	5	1057-05 .....	1
1001-01 .....	2	1020115.....	30	1044-30 .....	6	1060025.....	30
1001-02 .....	2	1020127.....	62	1044-31 .....	6	1060115.....	30
1001-03 .....	2	1020159.....	16	1044-40 .....	6	1061-00 .....	1
1001-04 .....	2	1020167.....	62	1044-41 .....	6	1061-01 .....	1
1001-05 .....	2	1020197.....	62	1044-51 .....	5	1061-02 .....	1
1002-01 .....	2	1020227.....	62	1044-52 .....	5	1061-03 .....	1
1002-02 .....	2	1020247.....	62	1044-53 .....	5	1061-04 .....	1
1002-03 .....	2	1020251.....	1	1044-80 .....	6	1061-05 .....	1
1002-04 .....	2	1020301.....	1	1044-81 .....	6	1063-01 .....	3
1002-05 .....	2	1020351.....	1	1044-82 .....	6	1063-02 .....	3
1003-01 .....	2	1020401.....	1	1044-83 .....	6	1063-03 .....	3
1003-02 .....	2	1023-00 .....	5	1044-84 .....	6	1063-04 .....	3
1003-03 .....	2	1023-01 .....	5	1044-85 .....	6	1063-05 .....	3
1003-04 .....	2	1023-50 .....	5	1044-90 .....	6	1063-06 .....	3
1003-05 .....	2	1024-11 .....	5	1044-91 .....	6	1070025.....	30
1004-01 .....	2	1024-12 .....	5	1050-00 .....	1	1070115.....	30
1004-02 .....	2	1024-13 .....	5	1050-01 .....	1	1080025.....	30
1005-01 .....	2	1024-14 .....	5	1050-02 .....	1	1080115.....	30
1005-02 .....	2	1024-51 .....	5	1050025.....	30	1090025.....	30
1005-03 .....	2	1024-52 .....	5	1050-03 .....	1	1090115.....	30
1005-04 .....	2	1024-53 .....	5	1050-04 .....	1	1100025.....	30
1005-05 .....	2	1030025.....	30	1050-10 .....	1	1100115.....	30
1005-06 .....	2	1030115.....	30	1050-11 .....	1	1110127.....	62
1006-01 .....	3	1030127.....	62	1050115.....	30	1110167.....	62
1006-02 .....	3	1030159.....	16	1050-12 .....	1	1110197.....	62
1006-03 .....	3	1030167.....	62	1050127.....	62	1110227.....	62
1006-04 .....	3	1030197.....	62	1050-13 .....	1	11201000.....	8
1006-05 .....	3	1030227.....	62	1050-14 .....	1	1120127.....	62
1007-02 .....	3	1030247.....	62	1050159.....	16	1120167.....	62
1007-03 .....	3	1030251.....	1	1050167.....	62	1120197.....	62
1007-04 .....	3	1030301.....	1	1050197.....	62	11202000.....	8
1007-05 .....	3	1030351.....	1	1050-20 .....	1	1120227.....	62
1007-06 .....	3	1030401.....	1	1050-21 .....	1	1120917.....	58
1010-01 .....	2	1034-11 .....	5	1050-22 .....	1	11211000.....	8
1010-02 .....	2	1034-12 .....	5	1050227.....	62	11212000.....	8
1010025.....	30	1034-13 .....	5	1050-23 .....	1	11221000.....	8
10101 .....	22	1034-14 .....	5	1050-24 .....	1	11222000.....	8
1010115.....	30	1034-51 .....	5	1050247.....	62	11231000.....	8
1010127.....	62	1034-52 .....	5	1050-30 .....	1	11232000.....	8
1010159.....	16	1034-53 .....	5	1050-31 .....	1	11241000.....	8
1010167.....	62	1040025.....	30	1050-32 .....	1	11242000.....	8
1010197.....	62	1040115.....	30	1050-33 .....	1	11251000.....	8
10102 .....	22	1040127.....	62	1050-34 .....	1	11252000.....	8
1010227.....	62	1040159.....	16	1054-11 .....	5	11261000.....	8
1010247.....	62	1040167.....	62	1054-12 .....	5	11262000.....	8
1010251.....	1	1040197.....	62	1054-13 .....	5	11271000.....	8
10103.....	22	1040227.....	62	1054-14 .....	5	11272000.....	8
1010301.....	1	1040247.....	62	1054-51 .....	5	11281000.....	8
1010351.....	1	1040251.....	1	1054-52 .....	5	11282000.....	8
10104.....	22	1040301.....	1	1054-53 .....	5	11291000.....	8
1010401.....	1	1040351.....	1	1056-00 .....	1	11292000.....	8
10105.....	22	1040401.....	1	1056-01 .....	1	1130127.....	62
10106.....	22	1042-00 .....	7	1056-02 .....	1	1130167.....	62
1018-01 .....	4	1043-00 .....	7	1056-03 .....	1	1130197.....	62
1018-02 .....	4	1043-01 .....	7	1056-04 .....	1	1130227.....	62

# INDEX BY PRODUCT NUMBER

1140127.....	62	13500004.....	55	13800006.....	52	14882009.....	67
1140167.....	62	13500005.....	55	13800007.....	52	14882010.....	67
1140197.....	62	13500006.....	55	13800008.....	52	14970005.....	63
1140227.....	62	13500007.....	55	13810000.....	51	14970006.....	63
1140917.....	58	13500008.....	55	13810004.....	51	14970007.....	63
1150127.....	62	13500009.....	55	13810005.....	51	15470002.....	79
1150167.....	62	13501004.....	55	13810006.....	51	15470005.....	79
1150197.....	62	13501005.....	55	13810007.....	51	15470008.....	79
1150227.....	62	13501006.....	55	13810008.....	51	15470067.....	79
1160917.....	58	13501007.....	55	13810907.....	51	15480002.....	79
11630005.....	19	13501008.....	55	13850000.....	51	15480005.....	79
11630006.....	19	13501009.....	55	13850005.....	51	15480008.....	79
11630007.....	19	13590000.....	49	13850007.....	51	15480067.....	79
11630008.....	19	13600000.....	49	13850008.....	51	15490002.....	79
11640005.....	19	13630000.....	49	13850009.....	51	15490005.....	79
11640006.....	19	13640000.....	49	13850010.....	51	15490008.....	79
11640007.....	19	13650000.....	49	13850011.....	51	15490067.....	79
11640008.....	19	13651056.....	48	13850056.....	51	15500002.....	79
11650005.....	20	13652067.....	48	13850067.....	51	15500005.....	79
11650006.....	20	13653008.....	48	13860900.....	52	15500008.....	79
11650007.....	20	13654009.....	48	13860900.....	52	15500067.....	79
11690001.....	21	13655010.....	48	13871000.....	48	15520002.....	80
11690002.....	21	13660000.....	49	13872900.....	48	15520005.....	80
11690005.....	21	13661056.....	48	13950000.....	48	15520008.....	80
11690006.....	21	13662067.....	48	13960000.....	48	15520067.....	80
11690007.....	21	13663008.....	48	13970000.....	48	15530002.....	80
11690008.....	21	13664009.....	48	13981000.....	48	15530005.....	80
1190917.....	58	13665010.....	48	13982000.....	48	15530008.....	80
1220917.....	58	13691000.....	54	13983000.....	48	15530067.....	80
1240917.....	58	13721000.....	57	13984000.....	48	15540002.....	80
12840000.....	26	13722000.....	57	13985000.....	48	15540005.....	80
12870000.....	26	13723000.....	57	13991000.....	48	15540008.....	80
12981102.....	30	13731000.....	54	13992000.....	48	15540067.....	80
12981202.....	30	13732000.....	54	13993000.....	48	15550002.....	80
12982105.....	30	13733000.....	54	13994000.....	48	15550005.....	80
12982106.....	30	13740300.....	50	13995000.....	48	15550008.....	80
12982107.....	30	13740400.....	50	14761004.....	68	15550067.....	80
12982108.....	30	13749300.....	50	14761005.....	68	15580067.....	80
12982205.....	30	13749400.....	50	14761006.....	68	15580002.....	80
12982206.....	30	13760300.....	50	14761007.....	68	15580005.....	80
12982207.....	30	13760305.....	50	14761008.....	68	15580006.....	80
12982208.....	30	13760306.....	50	14761009.....	68	15580008.....	80
12991102.....	30	13760307.....	50	14761010.....	68	15590002.....	80
12991202.....	30	13760308.....	50	14801004.....	68	15590005.....	80
12992102.....	30	13760400.....	50	14801005.....	68	15590006.....	80
12992202.....	30	13760405.....	50	14801006.....	68	15590008.....	80
12993105.....	30	13760406.....	50	14801007.....	68	15590067.....	80
12993106.....	30	13760407.....	50	14801008.....	68	15600004.....	75
12993107.....	30	13760408.....	50	14801009.....	68	15600005.....	75
12993108.....	30	13780300.....	50	14801010.....	68	15600006.....	75
12993205.....	30	13780400.....	50	14811004.....	65	15600007.....	75
12993206.....	30	13790305.....	50	14811005.....	65	15600008.....	75
12993207.....	30	13790306.....	50	14811006.....	65	15600009.....	75
12993208.....	30	13790307.....	50	14811007.....	65	15780005.....	79
13060300.....	50	13790308.....	50	14811008.....	65	15780008.....	79
13060400.....	50	13790309.....	50	14811009.....	65	15780067.....	79
13400005.....	55	13790405.....	50	14811010.....	65	19402.....	72
13400006.....	55	13790406.....	50	14882004.....	67	19403.....	72
13400007.....	55	13790407.....	50	14882005.....	67	19404.....	72
13400008.....	55	13790408.....	50	14882006.....	67	19405.....	72
13480056.....	53	13790409.....	50	14882007.....	67	19406.....	72
13480078.....	53	13800005.....	52	14882008.....	67	2001-00.....	88

# INDEX BY PRODUCT NUMBER

2001-01	88	2033-03	88	3004-11	13	4013-03	76
2001-02	88	2033-04	88	3004-12	13	4013-04	76
2001-03	88	2035-01	87	3004-13	13	4016-01	77
2001-04	88	2035-02	87	3004-14	13	4016-02	77
2001-05	88	2035-03	87	3004-15	13	4016-03	77
2004-00	87	2035-04	87	3004-16	13	4016-04	77
2004-01	87	2035-05	87	3007-00	14	4050-00	76
2004-02	87	2037-01	88	3015-01	13	4050-01	76
2004-03	87	2037-02	88	3015-02	13	4050-02	76
2004-04	87	2037-03	88	3015-03	13	4050-03	76
2004-05	87	2037-04	88	3015-04	13	4050-04	76
2005-01	88	2037-05	88	3015-05	13	4300	91
2005-02	88	2038-01	87	3015-11	13	4300/05AB	91
2005-03	88	2038-02	87	3015-12	13	4300/05WK	91
2005-04	88	2038-03	87	3015-13	13	4300K	91
2005-05	88	2038-04	87	3015-14	13	4301	91
2006-01	84	2038-05	87	3015-15	13	4302	91
2006-02	84	2039-00	89	3302D39	41	4305	91
2006-03	84	2039-01	89	3302D59	41	4305K	91
2006-04	84	2039-02	89	3302D79	41	4306	91
2006-05	84	2039-03	89	340M/LL	38	4307	91
2006-06	84	2039-04	89	340M/LR	38	4433-01	60
2006-07	84	2039-05	89	340S/ML	38	4433-02	60
2006-08	84	2039-98	89	340S/MR	38	4433-03	60
2006-09	84	2039-99	89	345M/LL	38	4433-04	60
2007-01	84	2041-00	83	345M/LR	38	4433-05	60
2007-02	84	2041-01	83	345S/ML	38	4433-10	60
2007-03	84	2041-02	83	345S/MR	38	4433-11	60
2007-04	84	2041-03	83	350LL	37	4433-22	60
2007-05	84	2041-07	83	350LR	37	4433-33	60
2007-06	84	2041-08	83	350ML	37	4442-00	15
2007-07	84	2041-09	83	350MR	37	4442-10	15
2007-08	84	2044-01	83	350SL	37	4461-02	64
2007-09	84	2044-02	83	350SR	37	4461-03	64
2009-00	89	2044-03	83	350XLL	37	4461-04	64
2020-01	84	2044-04	83	350XLR	37	4461-05	64
2020-02	84	2044-07	83	350XSL	37	4462-02	64
2020-03	84	2044-08	83	350XSR	37	4462-03	64
2020-04	84	2044-09	83	350XXLL	37	4462-04	64
2020-11	84	2045-01	85	350XXLR	37	4462-05	64
2020-12	84	2045-02	85	351LL	34	4462-06	64
2020-13	84	2045-03	85	351LR	34	4463-02	64
2020-14	84	2045-04	85	351ML	34	4463-03	64
2021-01	84	2045-05	85	351MR	34	4463-04	64
2021-02	84	3001-01	13	351SL	34	4463-05	64
2021-03	84	3001-02	13	351SR	34	4463-06	64
2021-04	84	3001-03	13	351XLL	34	4464-02	64
2021-11	84	3001-04	13	351XLR	34	4464-03	64
2021-12	84	3001-05	13	351XSL	34	4464-04	64
2021-13	84	3001-06	13	351XSR	34	4464-05	64
2021-14	84	3002-01	13	351XXLL	34	5005-01	32
2027-01	88	3002-02	13	351XXLR	34	5005-02	32
2027-02	88	3002-03	13	4004-01	77	5005-03	32
2027-03	88	3002-04	13	4004-02	77	5005-04	32
2027-04	88	3002-05	13	4004-03	77	5007-00	31
2027-05	88	3002-06	13	4004-04	77	5015-01	35
2032-00	88	3004-01	13	4005-01	78	5015-02	35
2032-01	88	3004-02	13	4005-02	78	5015-03	35
2032-02	88	3004-03	13	4005-03	78	5015-04	35
2032-03	88	3004-04	13	4005-04	78	5015-05	35
2033-01	88	3004-05	13	4013-01	76	5015-06	35
2033-02	88	3004-06	13	4013-02	76	5015-07	35

# INDEX BY PRODUCT NUMBER

5015-08	35	5066-83	32	7003-04	63	7071-11	61
5015-09	35	5066-90	32	7003-05	63	7071-12	61
5015-10	35	5066-91	32	7004-01	63	7071-13	61
5016-01	35	5066-92	32	7004-02	63	7071-15	61
5016-02	35	5066-93	32	7004-03	63	7071-17	61
5016-03	35	5073-01	33	7004-04	63	7071-55	61
5016-04	35	5073-02	33	7004-05	63	7071-57	61
5016-05	35	5073-03	33	7006-01	63	7071-59	61
5016-06	35	5073-04	33	7006-02	63	7072-02	60
5016-07	35	5073-05	33	7006-03	63	7072-03	60
5016-08	35	5073-06	33	7006-04	63	7072-04	60
5016-09	35	5073-07	33	7006-05	63	7072-05	60
5016-10	35	5073-08	33	7007-01	63	7081-00	61
5016-15	35	5073-09	33	7007-02	63	7082-00	61
5016-20	35	5073-10	33	7007-03	63	7082-22	61
5017-00	34	5074-01	33	7007-04	63	7083-00	61
5017-01	34	5074-02	33	7007-05	63	7083-33	61
5017-02	34	5074-03	33	7008-01	63	7084-00	61
5017-03	34	5074-04	33	7008-02	63	7084-44	61
5017-04	34	5074-05	33	7008-03	63	7085-00	61
5017-05	34	5074-06	33	7008-04	63	7085-55	61
5017-06	34	5074-07	33	7008-05	63	7086-00	61
5017-07	34	5074-08	33	7009-01	63	7086-66	61
5017-08	34	5074-09	33	7009-02	63	7087-00	61
5017-09	34	5074-10	33	7009-03	63	7091-03	59
5017-10	34	6001-01	29	7009-04	63	7091-04	59
5017-50	34	6001-02	29	7009-05	63	7091-05	59
5019-00	36	6001-03	29	7041-01	60	7091-06	59
5019-01	36	6001-04	29	7041-02	60	7091-07	59
5019-02	36	6003-01	28	7041-03	60	7091-08	59
5019-03	36	6003-02	28	7041-04	60	7091-50	58
5019-04	36	6003-03	28	7041-05	60	7091-51	58
5022-00	34	6003-04	28	7041-06	60	7091-52	58
5022-01	34	6004-01	27	7041-11	60	7091-55	58
5022-02	34	6004-02	27	7041-22	60	7091-56	58
5022-03	34	6004-03	27	7041-33	60	7091-57	58
5022-04	34	6004-04	27	7041-44	60	7091-58	58
5022-05	34	6005-00	26	7041-55	60	7101-00	10
5022-06	34	6101-01	9	7041-66	60	7101-20	10
5022-07	34	6101-02	9	7042-01	68	7102-00	10
5022-08	34	6101-03	9	7042-02	68	7102-22	10
5022-09	34	6102-00	9	7042-03	68	7103-00	10
5022-10	34	6103-00	9	7042-04	68	7104-02	11
5022-50	34	6104-00	10	7042-05	68	7104-22	11
5023-01	34	6105-00	10	7043-01	68	7104-23	11
5023-02	34	6114-00	9	7043-02	68	7115-00	11
5023-03	34	6117-00	12	7043-03	68	7115-20	11
5023-04	34	6157D02	24	7043-04	68	7116-01	11
5023-05	34	6157D03	24	7043-05	68	7117-01	12
5023-06	34	7001-01	63	7066-14	62	7117-02	12
5023-07	34	7001-02	63	7066-16	62	8001-01	22
5023-08	34	7001-03	63	7066-18	62	8001-02	22
5023-09	34	7001-04	63	7066-20	62	8001-03	22
5023-10	34	7001-05	63	7066-22	62	8001-04	22
5026-00	38	7002-01	63	7066-24	62	8001-05	22
5040-01	33	7002-02	63	7066-26	62	8004-01	21
5040-02	33	7002-03	63	7066-28	62	8004-02	21
5040-03	33	7002-04	63	7071-01	61	8004-03	21
5040-04	33	7002-05	63	7071-02	61	8004-04	21
5066-80	32	7003-01	63	7071-03	61	8004-05	21
5066-81	32	7003-02	63	7071-05	61	8004-06	21
5066-82	32	7003-03	63	7071-07	61	8004-11	21

# INDEX BY PRODUCT NUMBER

8004-12	21	8745-02	30	9013-03	15	9110-01	39
8004-13	21	8745-03	30	9013-04	15	9110-02	39
8004-14	21	8745-04	30	9014-00	15	9110-03	39
8004-15	21	8745-05	30	9023-00	18	9111-01	39
8004-16	21	8745-06	30	9023-01	18	9111-02	39
8006-01	23	8745-07	30	9023-02	18	9111-03	39
8006-02	23	8745-08	30	9023-03	18	9111-04	39
8006-03	23	8745-09	30	9023-04	18	9111-05	39
8006-04	23	8745-10	30	9100-00	40	9111-06	39
8006-05	23	8761-01	76	9101-01	44	9111-07	39
8011-01	25	8761-02	76	9101-02	44	9111-08	39
8011-02	25	8761-03	76	9101-03	44	9112-01	40
8011-03	25	8761-04	76	9101-04	44	9112-02	40
8013-01	24	8788-03	25	9101-05	44	9112-03	40
8013-02	24	8788-50	25	9101-06	44	9112-04	40
8017-00	22	8824-00	30	9102-01	44	9112-05	40
8017-01	22	8824-01	30	9102-02	44	9112-06	40
8017-02	22	8824-02	30	9102-03	44	9113-01	40
8017-03	22	8824-03	30	9102-04	44	9113-02	40
8017-04	22	8824-04	30	9102-05	44	9113-03	40
8017-51	22	8824-05	30	9102-06	44	9113-04	40
8017-52	22	8824-06	30	9103-01	44	9113-05	40
8017-53	22	8824-07	30	9103-02	44	9114-01	40
8017-54	22	8824-08	30	9103-03	44	9114-02	40
8020-00	21	8824-09	30	9104-01	44	9114-03	40
8020-01	21	8824-10	30	9104-02	44	9114-04	40
8023-01	23	8824-50	30	9104-03	44	9115-01	41
8023-02	23	8829-01	106	9105-01	44	9115-02	41
8023-03	23	8875-00	24	9105-02	44	9115-03	41
8023-04	23	8875-100	24	9105-03	44	9115-04	41
8023-05	23	9001-01	17	9105-04	44	9115-05	41
8103-01	72	9001-02	17	9105-05	44	9115-06	41
8103-02	72	9001-03	17	9105-06	44	9116-01	47
8103-03	72	9001-04	17	9105-07	44	9116-02	47
8103-04	72	9001-05	17	9105-08	44	9116-03	47
8103-05	72	9002-01	17	9105-20	45	9117-01	47
8103-06	72	9002-02	17	9105-21	45	9117-02	47
8105-01	72	9002-03	17	9105-22	45	9117-03	47
8105-02	72	9002-04	17	9105-24	45	9118-01	47
8105-03	72	9002-05	17	9105-25	45	9118-02	47
8105-04	72	9003-01	17	9105-26	45	9118-03	47
8105-05	72	9003-02	17	9106-01	44	9119-01	43
8105-06	72	9003-03	17	9106-02	44	9119-02	43
8106-01	71	9003-04	17	9106-03	44	9119-03	43
8106-02	71	9003-05	17	9106-04	44	9119-04	43
8106-03	71	9006-00	15	9106-05	44	9119-05	43
8106-04	71	9006-10	15	9106-06	44	9119-06	43
8106-05	71	9006-50	15	9107-01	45	9119-07	43
8106-06	71	9007-00	15	9107-02	45	9119-08	43
8108-02	73	9008-00	16	9107-03	45	9119-09	43
8108-03	73	9008-01	16	9107-04	45	9120-00	43
8108-04	73	9008-02	16	9107-05	45	9121-01	43
8108-05	73	9008-03	16	9107-06	45	9121-02	43
8108-06	73	9008-04	16	9108-01	45	9121-03	43
8120-00	72	9008-05	16	9108-02	45	9121-04	43
8121-00	71	9008-06	16	9108-03	45	9121-05	43
8130-00	71	9009-00	18	9108-04	45	9121-06	43
8140-00	71	9009-50	18	9108-05	45	9121-07	43
8726-00	102	9011-00	18	9108-06	45	9121-08	43
8726-40	102	9012-00	18	9109-01	39	9122-00	43
8726-50	102	9013-00	15	9109-02	39	9125-01	41
8745-01	30	9013-02	15	9109-03	39	9125-02	41

# INDEX BY PRODUCT NUMBER

9125-03	41	A125108	37	AB2801-20	75	BF7007-03	63
9125-04	41	A125109	37	AB2801-22	75	BF7007-04	63
9126-00	43	A125204	37	AB2801-24	75	BF7007-05	63
9127-01	47	A125205	37	AB2801-26	75	BF7007-06	63
9127-02	47	A125206	37	AB2801-28	75	BF7008-01	63
9127-03	47	A125207	37	AF1000-01	78	BF7008-02	63
9128-01	47	A125208	37	AF1000-02	78	BF7008-03	63
9128-02	47	A125209	37	AF1000-03	78	BF7008-04	63
9128-03	47	A141000	59	AF1000-04	78	BF7008-05	63
9130-00	42	A142000	59	AF1001-01	78	BF7008-06	63
9131-00	42	A142005	64	AF1001-02	78	BF7009-01	63
9132-00	42	A142006	64	AF1001-03	78	BF7009-02	63
9133-01	46	A142007	64	AF1001-04	78	BF7009-03	63
9133-03	46	A142008	64	AF1001-05	78	BF7009-04	63
9133-04	46	A142009	64	AF1001-06	78	BF7009-05	63
9133-06	46	A143000	59	AF1001-07	78	BF7009-06	63
9134-01	45	A143005	64	AF1001-08	78	BF7072-01	60
9134-02	45	A143006	64	BF1120917	58	BF7072-02	60
9135-01	45	A143007	64	BF1160917	58	BF7072-03	60
9135-02	45	A143008	64	BF1190917	58	BF7072-04	60
9136-01	45	A143009	64	BF1220917	58	BF7072-05	60
9136-02	45	A144000	59	BF1240917	58	BF7072-06	60
9148-01	46	A144005	64	BF4433-12	60	BS1000-10	53
9148-02	46	A144006	64	BF4433-16	60	BS1000-12	53
9148-03	46	A144007	64	BF4433-19	60	BS1000-14	53
9148-04	46	A144008	64	BF4433-22	60	BS1000-16	53
9148-05	46	A144009	64	BF4433-24	60	BS1000-18	53
9148-06	46	A145000	59	BF4433-26	60	BS1000-20	53
9150-01	42	A145005	64	BF5002-01	31	BS1000-21	53
9150-02	42	A145006	64	BF5002-02	31	BS1000-22	53
9151-01	46	A145007	64	BF5002-03	31	BS1000-24	53
9151-02	46	A145008	64	BF5002-04	31	BS1000-26	53
9151-03	46	A145009	64	BF5002-05	31	BS1000-28	53
9151-04	46	A146000	59	BF5002-06	31	BS1000-30	53
9151-05	46	A146005	64	BF5002-07	31	BS1000-40	53
9151-06	46	A146006	64	BF5002-08	31	BS1000-50	53
9221-00	106	A146007	64	BF5002-09	31	BS1000-52	53
9222-00	106	A146008	64	BF5002-10	31	BS1000-54	53
9307-01	62	A146009	64	BF5011-01	31	BS1000-56	53
9307-02	62	AB2000-10	75	BF5011-02	31	BS1000-58	53
9307-03	62	AB2000-12	75	BF5011-03	31	BS1000-59	53
9307-04	62	AB2000-14	75	BF5011-04	31	BS1000-99	53
9307-05	62	AB2000-16	75	BF5011-05	31	BS5000-12	53
9307-06	62	AB2340-00	74	BF5011-06	31	BS5000-14	53
9380-00	102	AB2340-06	74	BF5011-07	31	BS5000-16	53
9380-20	102	AB2341-00	74	BF5011-08	31	BS5000-18	53
9381-00	101	AB2341-06	74	BF5011-09	31	BS5000-20	53
9381-10	101	AB2350-00	74	BF5011-10	31	CA298	109
9381-50	101	AB2351-00	74	BF5030-01	31	CR1000	109
9383-00	101	AB2360-00	74	BF5030-02	31	CR2000	109
9383-10	101	AB2361-00	74	BF5066-80	32	CR3000	109
9383-50	101	AB2362-00	74	BF5066-81	32	CSG53510-A	80
9384-00	101	AB2700-10	82	BF5066-82	32	CSG53510-B	80
9384-50	101	AB2700-12	82	BF5066-83	32	CSG53510-C	80
9385-00	101	AB2700-14	82	BF7006-01	63	CSG53520-A	80
9386-00	101	AB2700-16	82	BF7006-02	63	CSG53520-B	80
9387-00	101	AB2700-18	82	BF7006-03	63	CSG53520-C	80
9387-50	101	AB2801-10	75	BF7006-04	63	CSG53550-A	80
A125104	37	AB2801-12	75	BF7006-05	63	CSG53550-B	80
A125105	37	AB2801-14	75	BF7006-06	63	CSG53550-C	80
A125106	37	AB2801-16	75	BF7007-01	63	DEB2	108
A125107	37	AB2801-18	75	BF7007-02	63	DEB3	108



# INDEX BY PRODUCT NUMBER

DEB4	108	KB5000-02	69	M60-033	98	NE7717-75	66
DEB6	108	KB6000-00	69	M60-034	97	NE7717-76	66
DEBV2	108	KB6000-01	69	M60-045	97	NE7717-77	66
DEBV3	108	M0501D01	92	M60-076-I	95	NE7718-71	67
DEBV4	108	M0501D05	92	M60-076-XS	95	NE7718-72	67
DEBV6	108	M10-070	107	M60-076-S	95	NE7718-73	67
DS1012-00	83	M10-070C	107	M60-076-M	95	NE7718-74	67
DS1013-10	83	M1128	103	M60-076-L	95	NE7718-75	67
DS1013-12	83	M1133	103	M60-078-S	95	NE7718-76	67
DS1013-14	83	M1133-38	103	M60-078-M	95	NE7718-77	67
DS1014-10	83	M1150	103	M60-078-L	95	NE7727-71	28
DS1014-12	83	M1151	103	M60-100	92	NE7727-72	28
DS1014-14	83	M1152	103	M60-200	92	NE7727-73	28
DS1014-16	83	M1153	103	M60-224	96	NE7727-74	28
DS1015-10	85	M13-055	105	M60-300	92	NE7727-75	28
DS1015-12	85	M13-100	104	M7035-S	104	NE7727-76	28
DS1015-14	85	M13-500	104	M7035-L	104	NE7728-71	28
DS1016-10	85	M190FF	93	M8125-18	104	NE7728-72	28
DS1016-12	85	M190FF-24	93	M8125-A	104	NE7728-73	28
DS1016-14	85	M190FFE	93	M8125-C	104	NE7728-74	28
DS1016-16	85	M2035	104	M8125-I	104	NE7728-75	28
EB2	107	M2035B	104	M8125-P	104	NE7728-76	28
EB2040-00	29	M3000D08	105	NE7701-01	65	NE7730-70	27
EB3	107	M3000KU	94	NE7701-02	65	NE7731-71	27
EB4	107	M3000U	94	NE7701-03	65	NE7731-72	27
EB5000-01	29	M3001U	94	NE7701-04	65	NE7731-73	27
EB5000-02	29	M3001KU	94	NE7701-05	65	NE7731-74	27
EB5000-03	29	M3002D08	105	NE7701-06	65	NE7731-75	27
EB6	107	M3006	93	NE7701-07	65	NE7731-76	27
EBV2	107	M3014D	93	NE7701-71	65	NE7732-72	77
EBV3	107	M3014D-6	93	NE7701-72	65	NE7732-73	77
EBV4	107	M3015	93	NE7701-73	65	NE7732-74	77
EBV6	107	M3015B	93	NE7701-74	65	NE7732-75	77
FD1	42	M3022	93	NE7701-75	65	NE7732-76	77
FD2	42	M3022-L	93	NE7701-76	65	NE7738-70	37
FD3	42	M3040	92	NE7701-77	65	NE7739-71	35
FF4061-00	68	M4209D06	105	NE7702-01	65	NE7739-72	35
HC1000-01	90	M4210D06	105	NE7702-02	65	NE7739-73	35
HC1000-02	90	M5001-L	103	NE7702-03	65	NE7739-74	35
HC2000-10	90	M5001-M	103	NE7702-04	65	NE7739-76	35
HC2000-12	90	M5001-S	103	NE7702-05	65	NE7739-77	35
HC2000-14	90	M5001-XL	103	NE7702-06	65	NE7739-78	35
HEB4	108	M5001-XS	103	NE7702-07	65	NE7739-79	35
HEB6	108	M5001-XXL	103	NE7702-08	65	NE7740-70	35
HS1003-08	85	M5001-XXXL	103	NE7702-71	65	NE7741-02	52
HS1003-10	85	M60-015	98	NE7702-72	65	NE7741-03	52
HS1003-12	85	M60-015A	97	NE7702-73	65	NE7741-04	52
HS1003-14	85	M60-015B	97	NE7702-74	65	NE7741-05	52
HS1003-16	85	M60-015CD	97	NE7702-75	65	NE7741-06	52
HS1005-08	86	M60-018	97	NE7702-76	65	NE7741-07	52
HS1005-10	86	M60-024-S	95	NE7702-77	65	NE7741-12	52
HS1005-12	86	M60-024-M	95	NE7707-72	67	NE7741-13	52
HS1005-14	86	M60-024-L	95	NE7707-73	67	NE7741-14	52
HS1005-16	86	M60-025-S	96	NE7707-74	67	NE7741-15	52
KB2040-00	69	M60-025-M	96	NE7707-75	67	NE7741-16	52
KB2040-10	69	M60-025-L	96	NE7707-76	67	NE7741-17	52
KB2040-40	69	M60-031-S	96	NE7710-70	66	NS2000-10	90
KB3001-00	70	M60-031-M	96	NE7710-71	66	NS2000-12	90
KB3001-10	70	M60-031-L	96	NE7717-71	66	NS2000-14	90
KB3001-40	70	M60-032	98	NE7717-72	66	SC1000-00	74
KB5000-00	69	M60-032B	98	NE7717-73	66	SC2000-10	90
KB5000-01	69	M60-032NC	98	NE7717-74	66	SC2000-12	90



# INDEX BY PRODUCT NUMBER

SC2000-14	90	T530NS	99	UDLE	57	UPT1433-25	57
SL2001-01	81	T530ST	99	UEP2031	57	UPT1433-35	57
SL2001-03	81	T531NS	99	UEP2051	57	UPT1533-15	57
SL2001-05	81	T531ST	99	UPL1330-0	56	UPT1533-25	57
SL2001-07	81	T532NS	99	UPL1330-15	56	UPT1533-35	57
SL2001-08	81	T532ST	99	UPL1330-25	56	UPT1633-15	57
SL2001-09	81	T533NS	99	UPL1430-0	56	UPT1633-25	57
SL2001-10	81	T533ST	99	UPL1430-15	56	UPT1633-35	57
SL2001-12	81	T580NS	99	UPL1430-25	56	UPT1733-15	57
SL2001-14	81	T580ST	99	UPL1530-0	56	UPT1733-25	57
SL2001-16	81	T600	99	UPL1530-15	56	UPT1733-35	57
SL5002-01	81	T6100	100	UPL1530-25	56	UPT2333-15	57
SL5002-03	81	T6110	100	UPL1630-0	56	UPT2333-25	57
SL5002-05	81	T6120	100	UPL1630-15	56	UPT2333-35	57
SL5002-07	81	T6130	100	UPL1630-25	56	UPT2433-15	57
SL5002-09	81	T6140	100	UPL1730-0	56	UPT2433-25	57
SL8001-01	81	T6160	100	UPL1730-15	56	UPT2433-35	57
SL8001-03	81	TX4442-00	16	UPL1730-25	56	UPT2533-15	57
SL8001-05	81	TX4442-02	16	UPL2330-0	56	UPT2533-25	57
SL8001-07	81	TX4442-03	16	UPL2330-15	56	UPT2533-35	57
SL8001-09	81	TX4442-04	16	UPL2330-25	56	UPT2633-15	57
SL8001-33	81	TX8020-00	23	UPL2430-0	56	UPT2633-25	57
SL8001-35	81	TX8023-21	23	UPL2430-15	56	UPT2633-35	57
SL8001-37	81	TX8023-22	23	UPL2430-25	56	UPT2733-15	57
SL8001-43	81	TX8023-23	23	UPL2530-0	56	UPT2733-25	57
SL8001-45	81	TX8023-24	23	UPL2530-15	56	UPT2733-35	57
SL8001-47	81	TX9023-00	18	UPL2530-25	56	USPK	57
T5010NS	100	TX9023-01	18	UPL2630-0	56	WLK02	109
T5010ST	100	TX9023-02	18	UPL2630-15	56	WLK02-W	109
T5020NS	100	TX9023-03	18	UPL2630-25	56		
T5020ST	100	TX9023-04	18	UPL2730-0	56		
T5030NS	100	TX9902-01	23	UPL2730-15	56		
T5030NS-10	100	TX9902-02	23	UPL2730-25	56		
T5030ST	100	TX9902-02/100	23	UPLP2330-15	56		
T5030ST-10	100	TX9902-03	23	UPLP2330-25	56		
T5031NS	100	TX9902-03/100	23	UPLP2430-15	56		
T5031NS-10	100	TX9902-04	23	UPLP2430-25	56		
T5031ST	100	TX9902-04/100	23	UPLP2530-15	56		
T5031ST-10	100	TX9902-05	23	UPLP2530-25	56		
T5032NS	100	TX9902-05/100	23	UPLP2630-15	56		
T5032NS-10	100	TX9902-06	23	UPLP2630-25	56		
T5032ST	100	TX9902-06/100	23	UPLP2730-15	56		
T5032ST-10	100	TX9903-01	31	UPLP2730-25	56		
T5033NS	100	TX9903-02	31	UPLT2333-15	57		
T5033NS-10	100	TX9903-03	31	UPLT2433-15	57		
T5033ST	100	TX9903-04	31	UPLT2533-15	57		
T5033ST-10	100	TX9903-05	31	UPLT2633-15	57		
T5040NS	100	TX9903-06	31	UPLT2733-15	57		
T5040ST	100	TX9903-07	31	UPLT2333-25	57		
T505	99	TX9903-08	31	UPLT2433-25	57		
T5050NS	100	TX9903-09	31	UPLT2533-25	57		
T5050ST	100	TX9903-10	31	UPLT2633-25	57		
T5060NS	100	TX9904-01	31	UPLT2733-25	57		
T5060ST	100	TX9904-02	31	UPMP1330-15	56		
T5070NS	100	TX9904-03	31	UPMP1430-15	56		
T5070ST	100	TX9904-04	31	UPMP1530-15	56		
T5080NS	100	TX9904-05	31	UPMP1630-15	56		
T5080ST	100	TX9904-06	31	UPMP1730-15	56		
T5090NS	100	TX9904-07	31	UPT1333-15	57		
T5090ST	100	TX9904-08	31	UPT1333-25	57		
T520NS	99	TX9904-09	31	UPT1333-35	57		
T520ST	99	TX9904-10	31	UPT1433-15	57		

# INDEX BY PRODUCT NAME

## Ab/Rib

Premium 6" Rib Belt .....	50	Sized Abdominal Binder .....	48
Premium 8" Rib Belt .....	50	Universal Abdominal Binder .....	49
Premium Sized Abdominal Binder .....	48	Universal Rib Belt .....	50
Premium Universal Abdominal Binder .....	48		

## Aluminum Splints

Baseball Finger Splints .....	39	Hand Grip Colles' Splint .....	45
Burnham Finger Splint .....	42	Long Metacarpal Splint .....	46
Burnham Thumb Splint .....	42	Mason Allen Hand Splint .....	46
Cock-Up Hand Splint .....	44	Padded Finger Splint Assortment Kit .....	40
Colles' Splint .....	44	Plastalume Splints .....	43
Elbow Splint .....	47	Plastic Colles' Splint .....	45
Finger Dressing/ Finger Cage .....	42	Posterior Leg Splint .....	47
Finger Splints .....	41	Short Metacarpal Splint .....	45
Foam Padded Finger Splints .....	41	STAX Finger Splint Kit .....	43
Foam Padded Splinting Material .....	41	Toad Finger Splints .....	39
Fold Over Finger Splints .....	39	Velcro® Finger Locks .....	43
Four Prong Finger Splints .....	40	Wrist & Forearm Splint .....	44
Full Hand Splint .....	46		
Gutter Splints .....	40		

## Ambulatory Aids

Crutches .....	109
Folding Walkers .....	109
Walking Canes .....	109

## Ankle

Air Ankle Stirrup .....	74	DeRoyal® Sports Ankle Brace .....	75
Air/Gel Ankle Stirrup .....	74	DeRoyal® Sports Ankle Brace 2 .....	75
Ankle Foot Orthosis .....	78	Elastic Ankle Sleeve .....	77
Canvas Lace-Up Ankle Brace .....	76	Figure 8 Wrap Ankle Support .....	77
Closed Heel Elastic Ankle Sleeve .....	78	Leatherette Ankle Brace .....	76
Confor® Ankle Stirrup .....	74	Neoprene Ankle Support .....	77
Deluxe Ankle Foot Orthosis .....	78	R.I.C.E. Ankle Splint .....	76
DeRoyal® Functional Ankle Brace .....	75		

## Arm

Arm Elevator/Sling w/ Ice Pockets .....	25	Narrow Pouch Arm Sling .....	22
Arm Sling w/ Foam Strap .....	23	Perforated Arm Elevator/Sling .....	25
Collar and Cuff .....	24	Premium Arm Sling .....	21
Convuluted Foam Arm Elevator/Sling .....	25	Specialty Arm Sling .....	21
CVA Sling .....	24	Strap Arm Sling .....	23
Deep Pocket Arm Sling .....	23	Triangle Bandage Arm Sling .....	24
Envelope Arm Sling .....	22	Universal Arm Sling .....	21
Light Blue Arm Sling .....	22		

# INDEX BY PRODUCT NAME

## Bandages & Wraps

Elastic Bandage ..... 108

## Bucks Traction

Bucks Traction w/ Triangle Buckle ..... 72  
Foam/Canvas Bucks Traction Boot ..... 72  
Full Foot Bucks Traction Boot..... 71  
Sized Bucks Traction Boot..... 72  
Specialty Bucks Traction Boot ..... 73  
Total Universal Bucks Traction..... 71  
Universal Bucks Traction ..... 71

## Clavicle

Acromioclavicular Splint ..... 14  
3-Way Clavicle Strap ..... 13  
Heavy Duty Clavicle Strap ..... 13  
Premium Clavicle Strap ..... 13

## Elbow

Elastic Elbow Support ..... 28  
Neoprene Elbow Sleeve ..... 28  
Neoprene Tennis Elbow Support w/ Pad ..... 27  
Pneu-Gel™ Tennis Elbow Support ..... 26  
Premium Tennis Elbow Support ..... 26  
Slimline Elbow Brace ..... 29  
Transition Elbow Brace ..... 29  
Universal Neoprene Tennis Elbow Support ..... 27  
Universal Tennis Elbow Support..... 26  
Sized Tennis Elbow Support..... 27  
Wrist & Elbow Splint ..... 29

## Foot & Heel

Aegis Post-Op Shoe ..... 83  
Black Canvas Cast Shoe ..... 87  
Canvas Cast Shoe ..... 87  
Canvas Post-Op Shoe..... 83  
Canvas Wooden Post-Op Shoe ..... 84  
Comfo-Eze Rocker Shoe ..... 87  
Darco® All Purpose Boot..... 85  
Darco® Med-Surg™ Shoe ..... 85  
Darco® SlimLine Cast Boot ..... 89  
Darco® Orthowedge™ Shoe ..... 86  
Deluxe Wooden Post-Op Shoe..... 84  
DeRoyal® Night Splint ..... 90  
E-Z Walk Cast Boot ..... 88  
Grey Spot Heel Cups..... 90  
Healing Shoe ..... 85  
Original Darby Shoe ..... 83  
Rainbow Cast Sandal ..... 88  
Silicone Heel Cups ..... 90  
Vinyl Cast Boot..... 88

## Foam Positioners

Arm Cradle ..... 97  
Body Aligner..... 98  
Concave Hip Abduction Pillow - High Density ..... 95  
Concave Hip Abduction Pillow - Wide ..... 95  
Concave Hip Abduction Pillow - Wide w/ Liner ..... 96  
Flat Hip Abduction Pillow - Narrow ..... 96  
Foot Cradle..... 97  
Foot Drop Support ..... 97  
Hip Abduction Pillow - Adjustable ..... 96  
Leg Positioner..... 98  
Leg Positioner w/ Foot Drop Strap ..... 98

# INDEX BY PRODUCT NAME

## Head/Neck

Cerv-o-Collar .....	3	Philadelphia® Cervical Collar.....	5
Comfo-Eze Cervical Collar.....	2	Philly EMT Extrication Collar .....	6
Contoured Cervical Collar.....	1	Philly One-Piece Extrication Collar .....	6
Definitive Collar .....	5	Repro Type Universal Head Halter.....	9
Deluxe Head Halter.....	9	Serpentine Cervical Collar.....	3
Discard Head Halter.....	10	Soft Foam Cervical Collar.....	4
Exo-Static/Exo-Bed Over-Door Traction Unit.....	10	STAT Neck Extrication Collar .....	7
Medium Density Cervical Collar .....	1	Spreader Bar.....	12
Medium Density Cervical Collar .....	2	Standard Over-Door Traction Unit .....	10
Medium Density Cervical Collar .....	3	Traction Cord.....	11
Medium/Firm Density Cervical Collar .....	1	Upper-7 Universal Head Halter .....	9
Over-Door Traction Kit w/ Repro Head Halter .....	11	Weight Bag .....	12
Over-Door Traction w/ Deluxe Head Halter .....	11	XTW Extended Wear Collar.....	8

## Hot & Cold Therapy

Chill-n-Wrap Ice Bag.....	102	Foam Hot/Cold Cold Therapy Wraps.....	101
DeRoyal® Cold Therapy Unit.....	99	Hot/Cold Pack.....	102
DeRoyal® Cold Therapy Unit Combos.....	99	Hot/Cold Therapy Blankets .....	100
DeRoyal® Hot/Cold Therapy Unit .....	99	Hot/Cold Therapy Wraps .....	100

## Knee

Arthroscopic Foam Knee Wrap.....	68	Knee Immobilizer w/ Elastic Straps .....	61
Closed Patella Knee Support.....	65	Open Patella Knee Support.....	65
Cutaway Knee Immobilizer.....	60	Sized Black Foam Knee Immobilizer.....	63
Cutaway Knee Immobilizer.....	60	Sized Canvas Knee Immobilizer .....	62
Cutaway Knee Immobilizer w/ Patella Strap .....	62	Sized Foam Knee Immobilizer.....	64
Cutaway Knee Immobilizer w/ Quick Fit Tabs.....	61	Sized Superlite Knee Immobilizer .....	63
Deluxe Knee Support w/ Trimmable Buttress .....	67	Sized Tietex Knee Immobilizer .....	64
Deluxe Tri-Panel Knee Immobilizer.....	58	Slimline Knee Brace .....	69
Elastic Knee Support .....	68	Slimline II Knee Brace .....	69
Fan Knee Stabilizer .....	67	Superlite Tri-Panel Knee Immobilizer .....	58
Foam Tri-Panel Knee Immobilizer.....	59	Tietex Tri-Panel Knee Immobilizer .....	59
Key Knee Brace .....	70	Transition Knee Brace .....	69
Knee Support w/ Trimmable Buttress.....	66	Universal Knee Support.....	66
Knee Immobilizer w/ Elastic Straps.....	60	Visco Elastic Knee Support .....	68

## Lumbo Sacral/Pelvis

Criss-Cross Lumbo Sacral Support .....	51	Paneled Pelvic Traction Belt .....	54
Elastic Lumbo Sacral Support.....	51	S-I Support Belt .....	54
Elastic/Neoprene Back Support.....	52	Specialty Pelvic Traction Belt.....	53
Heat & Mold Lumbo Sacral Support.....	52	Workers Back Support.....	54

## Pressure Reduction

Ankle Contracture Boot.....	91	Gelbo™ Heel & Elbow Protector.....	93
Bunny Boot .....	92	Heel/Ankle Protector - Full Foot Design .....	93
DeRoyal® Heel Suspension Boot.....	92	Heel Pillow .....	92
Foam Heel & Elbow Protector.....	93	Padded Heel & Elbow Protector.....	94

# INDEX BY PRODUCT NAME

## Specialty Nursing

Bed Rail Covers.....	105	Flexible I.V. Arm Board .....	104
Cardiac Cushion .....	104	Invalid Ring .....	108
Catheter Strap.....	103	Orthopedic Felt .....	106
Cervical Tension Pillow.....	106	Polyester Telemetry Pouch .....	105
Cerv-O-Rest Pillow.....	106	Surgical Bra .....	103
Disposable Telemetry Pouch .....	105	Trach Tube Holder.....	103
Elbow Immobilizer .....	104		

## Shoulder

Canvas Shoulder Immobilizer w/ Waist Strap .....	17	Shoulder Immobilizer w/ Foam Straps .....	18
Deluxe Swing & Swathe .....	15	Shoulder P.A.D.....	19
Elastic Shoulder Immobilizer .....	16	Shoulder P.A.D. II.....	19
Foam Sling & Swathe .....	15	Universal Foam Sling & Swathe .....	15
Foam Sling & Swathe .....	16	Universal Foam Shoulder Immobilizer .....	18
Muslin Shoulder Immobilizer w/ Waist Strap.....	17	Universal Sling Style Shoulder Immobilizer.....	18
Shoulder Abduction Pillow.....	20	Velpeau Shoulder Immobilizer.....	16
Shoulder Immobilizer w/ Canvas Swathe.....	17		

## Spine

Prolign® Lumbar Orthosis.....	55	Ultralign+® Lumbar Sacral Orthosis .....	56
Solace™ Low Profile Back Brace.....	55	Ultralign+® Thoracic Lumbar Sacral Orthosis .....	57

## Walkers

Inline Pneumatic Walker.....	79	Splintsrite™ .....	82
Inline Walker .....	80	Three-D Low Profile Walker .....	81
One-Piece Samson Walker .....	81	Three-D Lowest Profile Walker.....	81
Pacesetter™ Air Walker.....	79	Tracker EX Walker.....	80
Pacesetter™ II Walker .....	80		

## Wrist, Wrist/Forearm, Thumb

Black Foam Wrist & Thumb .....	37	Tietex Wrist & Wrist/Forearm Splint.....	31
Black Foam Wrist & Wrist/Forearm Splint .....	31	Thumb & Wrist/Thumb Spica.....	38
Black Wrist Splint .....	33	Universal Ambidextrous Wristlet .....	35
Canvas Wrist & Wrist/Forearm Splint.....	30	Universal Ambidextrous Wristlet w/ Abducted Thumb .....	37
Elastic Wrist Splint.....	35	Universal Black Foam Wrist & Wrist/Forearm Splint .....	32
Elastic Wrist Support .....	36	Universal Cutaway Wrist & Wrist/Forearm Splint .....	33
Lace Up Canvas Wrist Splint .....	34	Universal Foam Wrist & Wrist/Forearm Splint .....	32
Lace Up Suede Leatherette Wrist Splint.....	34	Universal Intera® Wrist & Wrist/Forearm Splint.....	32
Leatherette Wrist & Wrist/Forearm Splint .....	30	Universal Leatherette Wrist/Forearm Splint.....	31
Neoprene Wrist Support.....	35	Universal Leatherette Thumb Splint.....	38
Premium Wrist Splint .....	34	Universal Neoprene Wristlet .....	36
Premium Wrist & Thumb Splint.....	37	Universal Two-Panel Wrist & Wrist/Forearm Splint .....	33
Suede Leatherette Wrist&Wrist/Forearm Splint.....	30		



200 DeBusk Lane  
Powell, TN 37849  
888.938.7828 (inside U.S.)  
865.938.7828 (outside U.S.)  
[www.deroyal.com](http://www.deroyal.com)  
Reprint# 0-2098 rev. 6/08  
7.5m