


When it has to be right,
it has to be Medical Action.

Collection & Containment Solutions



Medical Action
INDUSTRIES INC.

Right. From the Beginning.

Table of Contents

Today, more than ever, medical professionals can't afford to settle for just average performance in their health care products. It's critical for you to know you can count on every product performing as expected. That's why Medical Action pays meticulous attention to every detail of its full line of Collection & Containment Solutions. Working under the most stringent quality controls ensures that every product from Medical Action is crafted to the same exacting standards of excellence, from design and manufacture to packaging and delivery. Because we know that superior performance begins with uncompromising quality, and demands that every item from Medical Action is always . . .

Right. From the Beginning.



Specialty Products	4
Patient Belongings Bags	5, 6
Lab Specimen Transport Bags	7-9
Zip Closure Bags	10, 11
Sterility Maintenance Covers	12
Medical Equipment and Cart Dust Covers	13
Biohazardous Waste Containment	14
Biohazardous Waste Bags	15-19
Autoclavable Biohazardous Waste Bags	20, 21
Laundry and Linen Bags	22-25
Water Soluble Laundry Bags	26
Hardware and Accessories	27
Institutional Trash Can Liners	28-33
Stackable/Non-Stackable Sharps Containers	34
Stackable Sharps Containers	35, 36
Non-Stackable Sharps Containers	36-38

Safeguard® and Tear Zone® are registered trademarks of BAGCO, Inc.



Bag dimensions, seam integrity and nomenclature imprinting must adhere to stringent guidelines.

Collection & Containment Solutions

Specialty Products

When it comes to Medical Waste Containment, Medical Action is the only supplier dedicated exclusively to meeting the needs of the health care provider. Offering a comprehensive selection of Specialty Bags, Biohazardous Waste Bags and Institutional Can Liners, as well as the ability to satisfy most any custom size, thickness, color or labeling preference, Medical Action is your sole source Medical Waste Containment Specialist. Our experienced Field Sales Consultants understand the unique demands of Health Care Providers and can respond with the right product to meet those needs. By performing an on-site utilization analysis they can determine where opportunities exist for cost savings and improved performance. Choosing the appropriate product specification for the intended application will save time and money, and help to optimize resources and operational efficiency.



Collection & Containment Solutions

Patient Belongings Bags

Patient Belongings Bags are a convenient and cost-effective means to collect, transport and store patient belongings from the time of admission, transfer to room, between departments and ultimate discharge. Bags come in a variety of sizes and strengths to handle most personal belongings.

- Available in a variety of standard sizes, colors and custom options designed to meet your facility's individualized needs.
- Colored film for personal privacy/security or clear for loss prevention.
- Durable polyethylene construction with built-in gusset expands capacity, is tear resistant and provides a stable platform for adding/removing clothing.
- Choice of closures include: Drawstring, Draw Tape, Rigid Handle with snap lock feature, Zip Closure and Adhesive Tabs.
- Standard labeling facilitates an easy to recognize write-on space for client name and room number. Also available without print.
- Ideal solution to help organize and manage patient belongings or set-up materials during stay.
- Set-Up bags identify specific use applications and client specific ease of control.
- **Customized labeling** provides immediate identification and **public relations/brand awareness** for your facility (name, logo, address, etc.) or specific department (ER, Same Day Surgery, Maternity, etc.). *Minimum order quantity required.*
- Typical areas of use: Admissions, ER, Outpatient Clinics, Labor and Delivery, Discharge and wherever client clothing/personal effects are collected.

Patient Belongings Bags — LLDPE Film, Flat Pack, Drawstring

Code#	Bag/Print Color	Size	Gauge/mil	Capacity/gal	Qty/cs	Qty/pkg
Print/Label: Patient Belongings, Patient Name, Room Number.						
3534	White/Blue	18 x 19 x 3	1.75	5	250	25
3539	White/Blue	20 x 19 x 4	1.75	5	250	25
3537	Clear/Blue	20 x 20 x 3	1.75	5	250	25
50-20	White/Blue	20 x 20 x 4	1.25	5	250	25
50-50	Clear/Blue	20 x 20 x 4	1.25	5	250	25
50-53*	Green/White	20 x 20 x 4	1.50	5	250	25
70-20	White/Blue	20 x 20 x 4	2.00	5	250	25
70-50	Clear/Blue	20 x 20 x 4	2.00	5	250	25
50-80	Fuschia/Black	20 x 20 x 4	2.00	5	250	25
70-52	Clear/Blue	20 x 20 x 4	2.75	5	250	25
3549	White/Blue	28 x 32 x 4	2.00	6-8	150	15
No Print/Label.						
3540	Clear/None	20 x 19	1.75	3-4	250	25
50-72	Clear/None	24 x 24	1.25	6-8	250	25

**Minimum order quantity required; please contact your distributor for details.
Measurements in inches / width x length x side gusset.*



Patient Belongings Bags



Patient Belongings Bags — LLDPE Film, Flat Pack, Rigid Handle

Code#	Bag/Print Color	Size	Gauge/mil	Capacity/gal	Qty/cs	Qty/pkg
Print/Label: Patient Belongings, Patient Name, Room Number.						
3536	White/Blue	20 x 19 x 3	1.75	5	200	20
D3535.1	Clear/Blue	20 x 20 x 3	2.00	5	100	10
50-21	White/Blue	20 x 20 x 4	1.50	5	250	25
50-51	Clear/Blue	20 x 20 x 4	1.50	5	250	25
70-21	White/Blue	20 x 20 x 4	2.50	5	250	25
70-51	Clear/Blue	20 x 20 x 4	2.50	5	250	25
D3535	Clear/Blue	21 x 20 x 3	2.00	5	200	20
70-31	Blue/Black x 2 sides	22 x 19 x 3	2.75	5	250	25

Personal Effects with Zip Closure

Code#	Bag/Print Color	Size	Gauge/mil	Capacity/qt	Qty/cs	Qty/pkg
Print/Label: Patient Name, Room Number, Eyewear, Hearing Aid, Dentures.						
3660*	Fuschia & Clear/Black	7½ x 13	2.00	1.5	1000	100

Patient Set-Up Bags with Draw Tape

Code#	Bag/Print Color	Size	Gauge/mil	Capacity/gal	Qty/cs	Qty/pkg
Print/Label: Patient Name, Room Number, Time/Date Issued, Equipment in Use, Hospital Property Do Not Remove, Return to Resp. Therapy Department (English/Spanish).						
50-30	Clear/Black	12 x 16	1.35	1-2	500	25
Print/Label: Patient Set-Up Bag, Patient Name, Room Number.						
D3602	Clear/Blue x 2 sides	12 x 16	2.00	1-2	500	–

Patient Belongings Bags — HDPE Film, Flat Pack, Bedside Bag with Blue Floral Print, Adhesive Tabs

Code#	Bag/Print Color	Size	Gauge/mic	Capacity/qt	Qty/cs	Qty/pkg
Print/Label: Personal Belongings, Patient Name, Room Number.						
50-45	White/Blue	7 x 11½ x 3½	15	1.0	2000	100

*Minimum order quantity required; please contact your distributor for details.
Measurements in inches / width x length x side gusset.

Lab Specimen Transport Bags

Laboratory Specimen Transport Bags facilitate collection, storage, identification, and transport of laboratory specimens and cultures within the hospital or surgery center, as well as off-site facilities and laboratories. Bags contain a dual pouch system which separates and protects samples from the accompanying paperwork.

- Available in a variety of standard sizes in clear or colored film with easy to use Zip Closure, Reclosable Adhesive, or Tamper Evident Adhesive closure systems.
- Leak resistant construction with Biohazard symbol for immediate identification meets federal regulations and OSHA requirements.
- Liquid tight specimen transport bags available for Pneumatic Tube Systems feature reinforced seams, Tamper Evident Adhesive Seal, scored linear tear line for easy opening and are approved shipping containers by NCCLS, CFR 49, DOT, NIDA.
- Customized labeling can identify specific use applications and departmental ease of control. *Minimum order quantity required.*
- Typical areas of use: Phlebotomy Labs, Chemistry & Hematology, Pathology, Toxicology, Microbiology, Emergency Room, Treatment Areas, Drug Screening, or wherever Routine Clinical Laboratory Specimen Collections are drawn.

Lab Specimen Transport Bags — LLDPE Film, Flat Pack, Zip Closure with Pouch, Flap

Code#	Bag/Print Color	Size	Gauge/mil	Capacity/qt	Qty/cs	Qty/pkg
Print/Label: Biohazard, Biohazard Symbol, Frozen, Refrigerate, Room Temp.						
3670*	Clear/Black & Orange	4 x 6	2.00	0.5	1000	100
4075	Clear/Black & Orange	6 x 9	2.00	1.0	1000	100
49-98	Clear/Black & Red	6 x 9	2.00	1.0	1000	100
D4078	Clear/Black & Orange	6 x 10	2.00	1.0	1000	100
C4917	Clear/Black & Orange	7 x 9½	2.00	1.5	1000	100
4915	Clear/Black & Red	8 x 10	2.00	1.5	1000	100
4915.ORG	Clear/Black & Orange	8 x 10	2.00	1.5	1000	100
4925	Clear/Black & Orange	10 x 10	2.00	1.8	1000	100
4916	Clear/Black & Orange	12 x 15	2.00	2.0	500	100

Print/Label: Biohazard, Biohazard Symbol.						
49-91	Clear/Black & Red	6 x 6	2.00	0.8	1000	100
49-97	Clear/Black & Red	6 x 9	2.00	1.0	1000	100
D4099	Clear/Black	6 x 9	2.00	1.0	1000	100
59-19	Clear/Black & Red	6 x 9	2.00	1.0	1000	100
49-95	Clear/Black & Red	8 x 8	2.00	1.5	1000	100
49-99	Clear/Black & Red	9 x 12	2.00	1.7	1000	100
59-94	Clear/Black & Red	12 x 13 x 4	2.00	2.0	1000	100

(continued on pg 8)



*Minimum order quantity required; please contact your distributor for details.
Measurements in inches / width x length x side gusset.

Lab Specimen Transport Bags

Lab Specimen Transport Bags — LLDPE Film, Flat Pack, Zip Closure with Pouch, Flap *[continued]*

Code#	Bag/Print Color	Size	Gauge/mil	Capacity/qt	Qty/cs	Qty/pkg
Print/Label: Biohazard (English/Spanish), Biohazard Symbol.						
4923	Red/Black	9 x 12 ¹ / ₂	2.00	1.7	1000	100
4930	Red Opaque/Black	12 x 12	2.00	2.0	1000	100
4921	Red/Black	12 x 15	2.00	2.0	500	100
Print/Label: STAT with Biohazard Symbol.						
3785*	Yellow/Black	6 x 9	2.00	1.0	1000	100
4001	Red/Black	6 x 9	2.00	1.0	1000	100
4002	Blue Opaque/Black	6 x 9	2.00	1.0	1000	100
Print/Label: Liver Transplant, Specimen Transport, with Biohazard Symbol.						
3723*	Red Translucent/Black	6 x 9	2.00	1.0	1000	100

Zip Closure with Pouch

Code#	Bag/Print Color	Size	Gauge/mil	Capacity/qt	Qty/cs	Qty/pkg
Print/Label: Biohazard, Biohazard Symbol, Frozen, Refrigerate, Room Temp.						
4914	Clear/Black & Red	6 x 9	2.00	1.0	1000	100
4914.ORG	Clear/Black & Orange	6 x 9	2.00	1.0	1000	100
Safeguard® Specimen Bag, Tear Zone®. Print/Label: Biohazard, Biohazard Symbol, Frozen, Refrigerate, Room Temp.						
4810	Clear/Black & Orange	8 x 10	2.00	1.5	1000	100

Slide Zip Closure with Pouch, Flap

Code#	Bag/Print Color	Size	Gauge/mil	Capacity/qt	Qty/cs	Qty/pkg
Print/Label: Biohazard, Biohazard Symbol, Frozen, Refrigerate, Room Temp.						
4510	Clear/Black & Orange	7 x 9	1.80	1.0	500	50

Reclosable Adhesive

Code#	Bag/Print Color	Size	Gauge/mil	Capacity/qt	Qty/cs	Qty/pkg
Print/Label: Biohazard, Biohazard Symbol.						
39-97	Clear/Black & Red	6 x 10	1.60	1.0	1000	100
39-99	Clear/Black & Red	8 x 11 ¹ / ₂	1.80	1.5	1000	100
Print/Label: Specimen Transport Bag, with instructions.						
47-STB1	Clear/White	8 x 12	1.50	1.5	200	–

Safeguard® and Tear Zone® are registered trademarks of BAGCO, Inc.

**Minimum order quantity required; please contact your distributor for details.
Measurements in inches / width x length x side gusset.*

Lab Specimen Transport Bags

Lab Specimen Transport Bags — LLDPE Film, Flat Pack, Permanent Adhesive Closure, Tamper Evident

Code#	Bag/Print Color	Size	Gauge/mil	Capacity/qt	Qty/cs	Qty/pkg
Print/Label: Biohazard, Biohazard Symbol.						
4065	Clear/Black & Red	6 x 10	2.00	1.0	1000	100
Perforated. Print/Label: Biohazard, Biohazard Symbol.						
59-89	Clear/Black & Red	6 x 9	2.00	1.0	1000	100
4060	Clear/Black & Red	6 x 10	2.00	1.0	1000	100
59-93	Clear/Black & Red	10 x 10	2.00	1.8	1000	100
Welded Seams, Flap, Pneumatic Tube Compatible. Print/Label: Biohazard, Biohazard Symbol.						
1130	Clear/Black & Red	6 x 9	2.00	1.0	1000	100

Economy Specimen Transport Bag — LLDPE Film, Flat Pack, Zip Closure with Pouch

Code#	Bag/Print Color	Size	Gauge/mil	Capacity/qt	Qty/cs	Qty/pkg
Print/Label: Biohazard, Biohazard Symbol.						
58-97	Clear/Black & Red	6 x 9	1.70	1.0	1000	100

Chemotherapy Specimen Transport Bags — LLDPE Film, Flat Pack, Zip Closure with Pouch

Code#	Bag/Print Color	Size	Gauge/mil	Capacity/qt	Qty/cs	Qty/pkg
Print/Label: CYTA Symbol, Chemotherapy Waste.						
9218	Yellow/Black	12 x 15	4.00	2.0	250	–
47-CBS1	Yellow/Black	12 x 15	4.00	2.0	100	–

Lab Utility Bags — LLDPE Film, Flat Pack, Zip Closure

Code#	Bag/Print Color	Size	Gauge/mil	Capacity/qt	Qty/cs	Qty/pkg
Pouch/Flap. No Print/Label.						
49-90	Clear/None	6 x 6	2.00	0.8	1000	100
49-96	Clear/None	6 x 9	2.00	1.0	1000	100
49-93	Clear/None	8 x 8	2.00	1.5	1000	100
Pouch/Flap, Gusset Bottom. No Print/Label.						
3712*	Clear/None	12 x 13 x 4	2.00	2.0	1000	100

**Minimum order quantity required; please contact your distributor for details.
Measurements in inches / width x length x side gusset.*

Zip Closure Bags provide a convenient low cost solution for collection, storage and transport of personal items or pre-packaging of individual articles for client use.

- Convenient Zip Closure system facilitates repeated easy opening and secure closing.
- Durable clear polyethylene film and reinforced seams provide excellent holding capacity.
- Available in a wide variety of sizes with and without white writing block to meet your requirements for multi-use versatility throughout the facility.
- Typical areas of use: Pharmacy, Central Supply, Dietary and throughout facility.



Zip Closure Bags — Clear, LLDPE Film, Flat Pack

Code#	Bag/Print Color	Size	Gauge/mil	Capacity	Qty/cs	Qty/pkg
White Writing Block. No Print/Label.						
Z2.0203W	Clear/None	2 x 3	2.00	–	1000	100
Z2.0304W	Clear/None	3 x 4	2.00	–	1000	100
Z2.0305W	Clear/None	3 x 5	2.00	–	1000	100
Z2.0406W	Clear/None	4 x 6	2.00	–	1000	100
Z2.0508W	Clear/None	5 x 8	2.00	–	1000	100
Z2.0609W	Clear/None	6 x 9	2.00	–	1000	100
Z2.0912W	Clear/None	9 x 12	2.00	–	1000	100
No Print/Label.						
Z2.0202	Clear/None	2 x 2	2.00	–	1000	100
Z2.0203	Clear/None	2 x 3	2.00	–	1000	100
Z2.0303	Clear/None	3 x 3	2.00	–	1000	100
Z2.0304	Clear/None	3 x 4	2.00	–	1000	100
Z2.0305	Clear/None	3 x 5	2.00	–	1000	100
Z2.0306	Clear/None	3 x 6	2.00	–	1000	100
Z2.0404	Clear/None	4 x 4	2.00	–	1000	100
Z2.0406	Clear/None	4 x 6	2.00	–	1000	100
Z2.0408	Clear/None	4 x 8	2.00	–	1000	100
Z2.0507	Clear/None	5 x 7	2.00	–	1000	100
Z2.0508	Clear/None	5 x 8	2.00	–	1000	100
Z2.0510	Clear/None	5 x 10	2.00	–	1000	100
Z2.0604	Clear/None	6 x 4	2.00	–	1000	100

(continued on pg 11)

Measurements in inches / width x length x side gusset.

Zip Closure Bags — Clear, LLDPE Film, Flat Pack (continued)

Code#	Bag/Print Color	Size	Gauge/mil	Capacity	Qty/cs	Qty/pkg
No Print/Label. (continued)						
Z2.0606	Clear/None	6 x 6	2.00	–	1000	100
Z2.0608	Clear/None	6 x 8	2.00	–	1000	100
Z2.0609	Clear/None	6 x 9	2.00	–	1000	100
Z2.0610	Clear/None	6 x 10	2.00	–	1000	100
Z2.0808	Clear/None	8 x 8	2.00	–	1000	100
Z2.0810	Clear/None	8 x 10	2.00	–	1000	100
Z2.0912	Clear/None	9 x 12	2.00	–	1000	100
Z2.1012	Clear/None	10 x 12	2.00	–	1000	100
Z2.1013	Clear/None	10 x 13	2.00	–	1000	100
Z2.1212	Clear/None	12 x 12	2.00	–	1000	100
Z2.1215	Clear/None	12 x 15	2.00	–	1000	100
Z2.1315	Clear/None	13 x 15	2.00	–	1000	100
Z2.1318	Clear/None	13 x 18	2.00	–	1000	100
Z4.0304	Clear/None	3 x 4	4.00	–	1000	100
Z4.0404	Clear/None	4 x 4	4.00	–	1000	100
Z4.0406	Clear/None	4 x 6	4.00	–	1000	100
Z4.0508	Clear/None	5 x 8	4.00	–	1000	10
Z4.0606	Clear/None	6 x 6	4.00	–	1000	100
Z4.0608	Clear/None	6 x 8	4.00	–	1000	100
Z4.0609	Clear/None	6 x 9	4.00	–	1000	100
Z4.0808	Clear/None	8 x 8	4.00	–	1000	100
Z4.0810	Clear/None	8 x 10	4.00	–	1000	100
Z4.0912	Clear/None	9 x 12	4.00	–	500	50
Z4.1012	Clear/None	10 x 12	4.00	–	500	50
Z4.1212	Clear/None	12 x 12	4.00	–	500	50
Z4.1215	Clear/None	12 x 15	4.00	–	1000	100
Z4.1318	Clear/None	13 x 18	4.00	–	500	50

Measurements in inches / width x length x side gusset.

Sterility Maintenance Covers

Sterility Maintenance Covers protect sterile packs from dust, dirt, odors, moisture and foreign materials, and are used to package, store, transport and maintain a dust free environment. Instrument covers facilitate transport and easy identification of contaminated devices.

- Clear polyethylene material provides a durable see-through barrier allowing easy identification of contents.
- Choice of Heat-Seal Closure or Self-Seal Adhesive strip secures contents and provides a tamper-evident barrier for sterile integrity.
- Easy opening linear tear strip technology allows for quick access and eliminates the need and risk associated with using scissors or blades to open.
- Available in a wide variety of sizes with and without Biohazard labeling.
- Typical areas of use: Central Processing and anywhere medical instruments and sterile products are packaged/reprocessed.

Sterility Maintenance Covers — LLDPE Film, Flat Pack, Self-Seal Adhesive Strip

Code#	Bag/Print Color	Size	Gauge/mil	Capacity	Qty/cs	Qty/pkg
Print/Label: Tear and Seal Instructions.						
810	Clear/Blue	6 x 10	2.25	–	1000	250
815	Clear/Blue	8 x 12	2.25	–	1000	250
820	Clear/Blue	10 x 15	2.25	–	1000	250
825	Clear/Blue	12 x 20	2.25	–	500	250
830	Clear/Blue	16 x 22	2.25	–	500	250
835	Clear/Blue	16 x 30	2.25	–	250	–
840	Clear/Blue	24 x 30	2.25	–	250	–
845	Clear/Blue	25 x 37	2.25	–	100	–

Heat-Seal

Code#	Bag/Print Color	Size	Gauge/mil	Capacity	Qty/cs	Qty/pkg
Print/Label: Tear and Seal Instructions.						
850	Clear/Blue	6 x 10	2.25	–	1000	250
855SMC	Clear/Blue	8 x 12	2.25	–	1000	250
860	Clear/Blue	10 x 15	2.25	–	1000	250
865	Clear/Blue	12 x 20	2.25	–	500	250
870	Clear/Blue	16 x 22	2.25	–	500	250
875	Clear/Blue	16 x 30	2.25	–	250	–
880SMC	Clear/Blue	24 x 30	2.25	–	250	–

Instrument Transport Covers — LLDPE Film, Flat Pack

Code#	Bag/Print Color	Size	Gauge/mil	Capacity	Qty/cs	Qty/pkg
Print/Label: Red Biohazard, Contaminated Materials, Instrument Transport Bag.						
52-00	Clear/Red	8 x 24 x 4	1.30	–	250	25
52-02	Clear/Red	11 x 31 x 5	1.30	–	250	25
52-04	Clear/Red	15 x 31 x 5	1.30	–	250	25

Measurements in inches / width x length x side gusset.

Medical Equipment and Cart Dust Covers

Equipment Dust Covers are designed for non-sterile protection of medical equipment, linen, and supplies during transport or while not in use. JCAHO recommends the use of protective covers for shielding clean linen from dirt and contaminants.

- Effective dust and dirt barrier helps keep contents cleaner longer.
- Flat or Gusseted covers are puncture resistant, easy to position, and come in a wide range of sizes.
- Available in flat pack or rolls for easy dispensing.
- Typical areas of use: Central Processing and anywhere medical instruments/equipment are used.

Medical Equipment and Cart Dust Covers — LLDPE Film, Flat Pack

Code#	Bag/Print Color	Size	Gauge/mil	Capacity	Qty/cs	Qty/pkg
Print/Label: Equipment Cover.						
43-01	Clear/Black	16 x 38 x 14	1.35	–	150	–
43-03	Clear/Black	16 x 56 x 14	1.35	–	100	–
43-05	Clear/Black	28 x 60 x 22	1.35	–	50	–
Print/Label: Orange circle with Biohazard Symbol.						
43-10	Clear/Black	18 x 38 x 16	1.50	–	150	–
No Print/Label.						
6089	Clear/None	7 x 8	4.00	–	1000	–
6058	Clear/None	16 x 22	2.00	–	1000	–
6059	Clear/None	24 x 30	2.00	–	500	–
48-46	Clear/None	37 x 66 x 25	1.00	–	100	–
6456	Green/None	38 x 50 x 27	1.30	–	100	–
48-65	Clear/None	42 x 68 x 24	1.25	–	100	–

Medical Equipment and Cart Dust Covers — LLDPE Film, Cored Roll

Code#	Bag/Print Color	Size	Gauge/mil	Capacity	Qty/cs	Qty/roll
No Print/Label.						
6220	Clear/None	26 x 60 x 22	1.50	–	50	50
6490	Clear/None	36 x 43 x 24	1.00	–	300	300
48-07	Clear/None	38 x 48 x 26	0.90	–	50	50
48-10	Clear/None	38 x 48 x 26	1.85	–	50	50
48-19	Clear/None	45 x 48 x 27	0.90	–	50	50
48-20	Clear/None	45 x 48 x 27	1.85	–	50	50
48-29	Clear/None	55 x 62 x 27	0.90	–	50	50
48-30	Clear/None	55 x 62 x 27	1.85	–	50	50
6625	Yellow/None	55 x 62 x 27	1.00	–	50	50
D6624	Clear/None	55 x 62 x 27	2.00	–	50	50
48-39	Clear/None	64 x 67 x 27	0.90	–	50	50
48-40	Clear/None	64 x 67 x 27	1.85	–	50	50
6092	Clear/None	64 x 67 x 27	2.00	–	50	50

Measurements in inches / width x length x side gusset.





Collection & Containment Solutions

Biohazardous Waste Bags

Biohazardous Waste Collection Bags are used to collect, store, and transport biohazardous/infectious waste within the health care facility, outpatient surgery center, alternate care clinic, physician’s office, and lab.

- Offered in both High Density/HMW and Linear Low Density resins for your specific application. Standard sizes include Flat Seal, Star Seal, and Gusseted Bottoms.
- Customized printing provides easy identification of facility or specific department. *Minimum order quantity required.*
- Biohazardous or Infectious labeling in both English and Spanish for easy identification of hazardous waste to reduce the risk of cross-contamination during collection and transport.
- International Biohazard Symbol in easy to see contrasting color helps reduce risk of cross-contamination from infectious waste handling.
- Meets 165 gram Dart Impact Resistance Test.
- Compact coreless interleaved rolls or flat pack dispensing offers complete versatility for ease of use.
- Specifications available to comply with Florida, Pennsylvania and New York state requirements.
- Accessory Hampers and wall mounted dispensers for convenient placement throughout the facility.

Collection & Containment Solutions

Biohazardous Waste Containment

Biohazardous Waste Liners are manufactured to meet strict performance standards in accordance with applicable Federal and specific State guidelines. Quality Assurance is monitored to ensure load capacity and seal integrity, as well as linear strength. As individually specified by product category, Biohazardous Waste Liners comply with A.S.T.M. standards for 165 gram dart impact resistance test (D1709-04), 480 grams Elmendorf linear tear strength (D1922-06A) and US DOT CFR 49 173-197.

- Premium resins help ensure excellent performance with superior strength film and reinforced seams that are impervious to liquids.
- Latex free and environmentally safe for incineration disposal.
- Wide variety of standard sizes, film gauges and colors provides easy availability to meet your waste handling requirements.
- Color coded to help facilitate waste stream segregation and minimize disposal costs associated with biohazardous waste.



Biohazardous Waste Bags — Meet A.S.T.M. Dart Test Requirements, LLDPE Film, Flat Pack

Code#	Bag/Print Color	Size	Gauge/mil	Capacity/gal	Qty/cs	Qty/pkg
Print/Label: Biohazardous/Infectious Waste, Biohazard Symbol.						
2302	Red/Black	14 x 19	2.00	2-3	500	–
2634	Red/Black	19 x 23	1.25	5-7	400	–
115H	Red/Black	23 x 23	1.10	7-10	500	–
116	Red/Black	23 x 23	1.50	7-10	500	–
116BX	Red/Black	23 x 23	1.50	7-10	400	100
108M	Red/Black	25 x 34	1.20	12-16	250	–
108MBX	Red/Black	25 x 34	1.20	12-16	200	50
108MP	Red/Black	25 x 34	1.20	12-16	250	25
2305	Red/Black	25 x 35	2.25	12-16	200	–
D2210	Red/Black	30 x 36	1.00	20-30	250	–
2202	Red/Black	30 x 40	3.00	30	100	–
2042	Red/Black	30½ x 41	1.10	30	250	–
117M	Red/Black	30½ x 41	1.20	30	250	–
117MBX	Red/Black	30½ x 41	1.20	30	200	50
117MP	Red/Black	30½ x 41	1.20	30	250	25
117	Red/Black	30½ x 41	1.50	30	250	–
117A	Red/Black	30½ x 41	2.00	30	250	–

(continued on pg 16)

Measurements in inches / width x length x side gusset.



Biohazardous Waste Bags

Biohazardous Waste Bags — Meet A.S.T.M. Dart Test Requirements, LLDPE Film, Flat Pack *(continued)*

Code#	Bag/Print Color	Size	Gauge/mil	Capacity/gal	Qty/cs	Qty/pkg
Print/Label: Biohazardous/Infectious Waste, Biohazard Symbol. <i>(continued)</i>						
119	Red/Black	30½ x 41	3.00	30	125	–
2020	Red/Black	31 x 38	1.85	30	100	–
104M	Red/Black	33 x 40	1.20	33	250	–
2045	Red/Black	38 x 45	1.00	44	100	–
172M	Red/Black	38 x 45	1.20	44	250	–
109	Red/Black	38 x 45	1.50	44	200	–
110	Red/Black	38 x 45	1.75	44	125	–
2204	Red/Black	38 x 45	2.00	44	100	–
169	Red/Black	38 x 45	3.00	44	100	–
D2209	Red/Black	40 x 46	1.20	40-45	250	–
173	Red/Black	43 x 48	2.00	55	100	–
2122	Red/Black	43 x 55	1.50	55	200	–
2125	Red/Black	43 x 55	3.00	55	100	–
Star Seal. Print/Label: Biohazardous/Infectious Waste, Biohazard Symbol.						
2214	Red/Black	11 x 14	1.20	1-3	1000	–
57-00	Red/Black	24 x 24	1.05	8-10	500	–
57-05	Red/Black	30½ x 41	1.05	30	250	–
47-73	Red/Black	40 x 46	1.60	40-45	100	–
Star Seal, Twist Ties. Print/Label: Biohazardous/Infectious Waste, Biohazard Symbol.						
47-50	Red/Black	24 x 24	1.20	8-10	500	–
47-71	Red/Black	31 x 41	1.60	30-32	100	–

Biohazardous Waste Bags — Meet A.S.T.M. Dart Test Requirements, LLDPE Film, Coreless Roll

Code#	Bag/Print Color	Size	Gauge/mil	Capacity/gal	Qty/cs	Qty/roll
Print/Label: Biohazardous/Infectious Waste, Biohazard Symbol.						
RD610	Red/Black	11 x 14½	1.25	1-3	200	20
RD630	Red/Black	14½ x 19	1.25	2-3	200	20
RD650	Red/Black	16 x 24	1.25	5-6	200	20
R116	Red/Black	23 x 23	1.50	7-10	300	30
R108M	Red/Black	25 x 34	1.20	12-16	300	30
R117M	Red/Black	30½ x 41	1.20	30	250	25
R117	Red/Black	30½ x 41	1.50	30	200	20
R172M	Red/Black	38 x 45	1.20	44	200	20
R109	Red/Black	38 x 45	1.50	44	150	15
R110	Red/Black	38 x 45	1.75	44	150	15

(continued on pg 17)

Measurements in inches / width x length x side gusset.

Biohazardous Waste Bags

Biohazardous Waste Bags — Meet A.S.T.M. Dart Test Requirements, LLDPE Film, Coreless Roll *(continued)*

Code#	Bag/Print Color	Size	Gauge/mil	Capacity/gal	Qty/cs	Qty/roll
Star Seal. Print/Label: Biohazardous/Infectious Waste, Biohazard Symbol.						
37-97	Red/Black	24 x 24	1.20	12-16	150	150
55445	Red/Black	24 x 24	1.50	8-10	250	50
55436	Red/Black	24 x 32	1.50	12-16	250	50

Biohazardous Waste Bags — Meet A.S.T.M. Dart/Elmendorf & DOT CFR 49 173-197, LLDPE Film, Flat Pack

Code#	Bag/Print Color	Size	Gauge/mil	Capacity/gal	Qty/cs	Qty/pkg
Print/Label: Biohazardous/Infectious Waste, Biohazard Symbol.						
50-40	Red/Black	11 x 14	1.50	1-3	500	100
50-42	Red/Black	11 x 14	1.50	1-3	500	50
F116	Red/Black	23 x 23	1.20	7-10	500	–
F116BX	Red/Black	23 x 23	1.20	7-10	400	100
F124	Red/Black	24 x 30	1.20	8-12	500	–
2635	Red/Black	25 x 35	1.35	12-16	200	–
F145	Red/Black	30 x 36	1.20	20-30	250	–
F117	Red/Black	30½ x 41	1.20	20-30	250	–
F118	Red/Black	33 x 39	1.20	33	250	–
F130	Red/Black	37 x 50	1.30	44	150	–
F110	Red/Black	37 x 50	1.50	44	150	–
F109	Red/Black	38 x 45	1.20	44	200	–
F169	Red/Black	38 x 45	3.00	44	100	–
F149	Red/Black	40 x 45	1.20	40-45	250	–
F135	Red/Black	40 x 46	1.20	40-45	100	–
45-75	Red/Black	55 x 60	2.00	60	50	–

Side Gusset, Twist Ties. Print/Label: Biohazardous/Infectious Waste, Biohazard Symbol.

47-41	Red/Black	23 x 41 x 8	2.40	30-40	100	–
47-43	Red/Black	23 x 46 x 17	2.40	40-55	100	–
47-45	Red/Black	23 x 55 x 17	2.40	45-55	100	–

Star Seal. Print/Label: Biohazardous/Infectious Waste, Biohazard Symbol.

45-45	Red/Black	40 x 55	1.25	45-55	150	–
-------	-----------	---------	------	-------	-----	---

Star Seal, Twist Ties. Print/Label: Biohazardous/Infectious Waste, Biohazard Symbol.

45-54	Red/Black	24 x 32	1.25	12-16	250	–
45-30	Red/Black	28 x 31	1.25	20	250	–

Measurements in inches / width x length x side gusset.



Right. From the Beginning.

www.Medical-Action.com ■ 1.800.645.7042

Biohazardous Waste Bags

Chemotherapy Waste Collection Bags are used for transport of anti-neoplastic medications and chemotherapy waste. Bags feature reinforced seals for leak resistance and meet labeling requirements for chemotherapy waste collection.

- Color Coded and Bold Labeling with International Biohazard Symbol for easy identification.
- Heavy Duty Construction for Puncture Resistance.
- Also available in a 12" x 15" Specimen Transport Size with Zip Closure.



Chemotherapy Waste/CYTA Symbol, LLDPE Film, Flat Pack

Code#	Bag/Print Color	Size	Gauge/mil	Capacity/gal	Qty/cs	Qty/pkg
Print/Label: Chemotherapy Waste, CYTA Symbol.						
47-CBL4	Yellow/Black	31 x 41	4.00	20-30	100	–
Print/Label: Biohazard, Biohazard Symbol, CYTA, Meets A.S.T.M. Dart Test Requirements.						
4133	White/Black	24 x 33	2.00	12-16	250	–
4138	White/Black	31 x 38	2.00	20-30	200	–
4210	Yellow/Black	30 x 43	3.00	20-30	100	–
Star Seal. Print/Label: Biohazard Chemotherapy Waste, Chemotherapy Symbol, Meets A.S.T.M. Dart Test Req.						
9230	Yellow/Black	25 x 34	2.00	12-16	250	–



Anti-Static Kick Bucket Liners comply with Operating Room requirements of NFPA Standard 56A. Bags contain a carbon additive to help prevent the buildup of electrostatic charges.

Anti-Static Kick Bucket Liner — NFPA 56A, Twist Ties, LLDPE Film, Flat Pack

Code#	Bag/Print Color	Size	Gauge/mil	Capacity/gal	Qty/cs	Qty/pkg
Print/Label: Anti-Static, Biohazard, Infectious Waste, Biohazard Symbol.						
50-02	Red/Black	12 x 24 x 8	1.30	7	500	–
50-04	Red/Black	18 x 16 x 8	1.30	7	500	–

Biomedical Waste Containers provide an effective and inexpensive means for handling large amounts of biomedical waste while controlling disposal costs.

- 3-Part System: Outer Corrugated Container, Inner Corrugated Container, Plastic Liner separation.
- Locking Lid, Puncture and Leak Resistant, Incinerator-Safe. Store flat to minimize space.



Biomedical Waste Containers — Corrugated Box with Liner, Flat Pack

Code#	Bag/Print Color	Size	Gauge	Capacity	Qty/cs	Qty/pkg
Print/Label: Biohazardous/Infectious Waste, Biohazard Symbol.						
10-2002	Red/Black & Red	6 x 10¾ x 12	–	2.5 gal/41.6 L	30	–
10-2005	Red/Black & Red	12 x 12 x 12	–	5 gal/9.45 L	25	–
10-2011	Red/Black & Red	12 x 24 x 12	–	11 gal/18.9 L	14	–

Measurements in inches / width x length x side gusset.

Biohazardous Waste Bags

Biohazardous Waste Bags — Meet A.S.T.M. Dart Test Requirements, HDPE Film, Flat Pack

Code#	Bag/Print Color	Size	Gauge/mic	Capacity/gal	Qty/cs	Qty/pkg
Print/Label: Biohazardous/Infectious Waste, Biohazard Symbol.						
HRD24246	Red/Black	24 x 24	6.00	8-10	1000	–
HRD24248	Red/Black	24 x 24	8.00	8-10	2000	–
HRD242411	Red/Black	24 x 24	11.00	8-10	1000	–
HRD334014	Red/Black	33 x 40	14.00	33	500	–
44-01	Red/Black	33 x 40	16.00	33	250	–
HRD384616	Red/Black	38 x 46	16.00	44	250	–
HRD404816	Red/Black	40 x 48	16.00	40-45	250	–
Star Seal. Print/Label: Biohazardous/Infectious Waste, Biohazard Symbol.						
44-05	Red/Black	31 x 41	16.00	30-32	250	–
Twist Ties. Print/Label: Biohazardous/Infectious Waste, Biohazard Symbol.						
44-13	Red/Black	40 x 46	16.00	40-45	200	–

Biohazardous Waste Bags — Meet A.S.T.M. Dart Test Requirements, HDPE Film, Coreless Roll

Code#	Bag/Print Color	Size	Gauge/mic	Capacity/gal	Qty/cs	Qty/pkg
Print/Label: Biohazardous/Infectious Waste, Biohazard Symbol.						
54798	Red/Black	17 x 18	11.00	4-6	1000	100
RS242411RH	Red/Black	24 x 24	11.00	8-10	1000	50
RS243311RH	Red/Black	24 x 33	11.00	12-16	1000	50
RS304314RH	Red/Black	30½ x 43	14.00	30-32	250	25
55446	Red/Black	31 x 41	17.00	30-32	250	50
RS384616RH	Red/Black	38 x 46	16.00	44	250	25
RS404816RH	Red/Black	40 x 48	10.00	40-45	250	25
Star Seal, Twist Ties. Print/Label: Biohazardous/Infectious Waste, Biohazard Symbol.						
44-03	Red/Back	17 x 18	8.00	4-6	1000	100



Wall Mounted Roll Dispensers feature:

- Durable wire frame construction with built-in lid assembly requires limited space and secures easily to the wall.

Biohazardous Waste Containment Accessories

Code#	Description	Capacity/gal	Qty/cs
A510	(1) Gallon Wall Mounted Hamper System for RD610 Bags, 9 x 12 x 6	1	6
A530	(3) Gallon Wall Mounted Hamper System for RD630 Bags, 11 x 15 x 7½	3	4
3000	Biohazard Spill Kit, powder-free gloves, PDI-Sanicleth w/germicidal, Red-Z ¾ oz, plastic scoop, 9 x 12 x 2 red zip closure bag, 24 x 33 brown bag		4

Measurements in inches / width x length x side gusset.

Autoclavable Biohazardous Waste Bags

Autoclavable Biohazardous Waste Bags are built to withstand the rigors of Steam Sterilization while meeting the requirements for the collection, storage, transport, and processing of biohazardous waste within the acute care facility. Autoclavable bags made from high strength Polypropylene (PP) Film have an operating temperature of 274°F (134°C) and can tolerate temperatures of up to 285°F (140°C). High Density Polyethylene (HDPE) Film Autoclavable bags have an operating temperature of 250°F (120°C) and can tolerate temperatures of up to 260°F (127°C). Bags contain an integral steam processing indicator that turns brown when autoclaved.

- Steam indicator identifies processed versus unprocessed bags to help reduce risk of cross-contamination from infectious waste handling.
- Constructed of High Strength Polypropylene for maximum performance. Eliminates Double Bagging Procedure to reduce cost and labor.
- Welded Seams increase durability to resist punctures, tears and leaks.
- Biohazardous or Infectious labeling in both English and Spanish or French.
- International Biohazard Symbol in easy to see contrasting color.
- Meets 165 gram Dart Impact Resistance Test.
- Available in several color-coded styles for easy identification by department.
- Twist Ties included for easy securement.



Autoclavable Biohazardous Waste Bags — Polypropylene Film with Indicator, Flat Pack, Twist Ties

Code#	Bag/Print Color	Size	Gauge/mil	Capacity/gal	Qty/cs	Qty/pkg
Max. Temp 285°F. Print/Label: Biohazard Labeling, Biohazard Symbol, English/Spanish.						
880	Orange/Black	8 x 12	1.80	1-2	400	–
5001.4	Red/Black	8 x 12	2.00	1-2	400	–
852	Clear/Black	12 x 24	1.80	2-3	200	–
5010.1	Red/Black	12 x 24	2.00	2-3	100	–
881	Orange/Black	14 x 18	1.80	3-4	400	–
5012	Red/Black	14 x 19	2.00	3-4	500	–
5119	Buff/Black	19 x 23	1.75	7	200	–
882	Orange/Black	19 x 24	1.80	7-10	200	–
842	Red/Black	19 x 24	1.80	7-10	200	–
5013	Red/Black	19 x 24	2.00	7-10	200	–
844	Red/Black	25 x 30	1.80	12-14	200	–
854	Clear/Black	25 x 30	1.80	12-14	200	–
884	Orange/Black	25 x 35	1.80	13-16	200	–
846	Red/Black	25 x 35	1.80	13-16	200	–
5237	Red/Black	25 x 35	2.00	13-16	200	–
855	Clear/Black	31 x 38	1.80	20-30	100	–
856	Red/Black	31 x 38	1.80	20-30	100	–
858	Buff/Black	31 x 38	1.80	20-30	100	–

[continued on pg 21]

Measurements in inches / width x length x side gusset.

Autoclavable Biohazardous Waste Bags

Autoclavable Biohazardous Waste Bags — Polypropylene Film with Indicator, Flat Pack, Twist Ties

[continued]

Code#	Bag/Print Color	Size	Gauge/mil	Capacity/gal	Qty/cs	Qty/pkg
Max. Temp 285°F. Print/Label: Biohazard Labeling, Biohazard Symbol, English/Spanish.						
8-300	Clear/Black	36 x 35	1.80	33	100	–
848	Red/Black	38 x 47	1.80	44	100	–
864	Buff/Black	38 x 47	1.80	44	100	–
872	Clear/Black	38 x 47	1.80	44	100	–
886	Orange/Black	38 x 47	1.80	44	100	–
5028	Red/Black	40 x 70	2.00	55-60	100	–
2120	Clear/Black	45 x 58	1.80	55	100	–
888	Clear/Black	106 x 86	1.80	60+	25	–

Print/Label: Biohazard Labeling, Biohazard Symbol, English/French.

842Y	Yellow/Black	19 x 24	2.20	7-10	200	–
5142.F	Clear/Black	24 x 36	1.75	12-16	200	–
844Y	Yellow/Black	25 x 30	2.20	12-14	200	–
846Y	Yellow/Black	25 x 35	2.20	13-16	200	–
5064.F	Yellow/Black	30 x 38	2.00	20-30	125	–
856Y	Yellow/Black	31 x 38	2.20	20-30	100	–
848Y	Yellow/Black	38 x 48	2.20	44	100	–

No Print/Label.

5023	Clear/None	38 x 47	1.80	44	100	–
5034	Clear/Black	33 x 40	2.00	33	200	–
5207	Clear/None	70 x 55	1.80	60	100	–

Autoclavable Biohazardous Waste Bags — Polypropylene Film with Indicator, Coreless Roll, Twist Ties

Code#	Bag/Print Color	Size	Gauge/mil	Capacity/gal	Qty/cs	Qty/roll
Print/Label: Biohazard Labeling, Biohazard Symbol.						
D5124	Buff/Black	25 x 35	1.75	13-16	200	20

Autoclavable Biohazardous Waste Bags — HDPE Film with Indicator, Flat Pack, Twist Ties

Code#	Bag/Print Color	Size	Gauge/mic	Capacity/gal	Qty/cs	Qty/pkg
Max Temp 260° F. Print/Label: Biohazard Labeling, Biohazard Symbol.						
8-900	Red/Black	8½ x 11	25.40	1-2	1000	–
8-904	Red/Black	14 x 19	25.40	3-4	200	–
8-908	Red/Black	19 x 24	25.40	7-10	200	–
8-912	Red/Black	25 x 35	25.40	13-16	200	–

Measurements in inches / width x length x side gusset.



A full line of printed and unprinted, color-coded **Laundry and Linen Collection Bags** provides a safe and efficient method of handling soiled or infectious linen and laundry from point of collection, through storage and transport, and into final processing.

- Wide range of standard sizes and film thickness, with and without print, provides maximum flexibility to meet your facility's needs. Custom specifications and printing available. *Minimum order quantity required.*
- Offered in both High Density and Linear Low Density resins for your specific application, bags are color-coded for easy identification by department.
- Reinforced seams add durability and superior leak resistance.
- Heavy Duty Linear Low Density Bags provide exceptional odor containment and are designed with superior tear and impact resistance for use in Multi-Story Laundry Chutes.
- Typical areas of use: Housekeeping, Laundry, and most any treatment area or room where collection and transport of linens/clothing is required.



Laundry and Linen Bags — LLDPE Film, Flat Pack

Code#	Bag/Print Color	Size	Gauge/mil	Capacity/gal	Qty/cs	Qty/pkg
Print/Label: Soiled Linen.						
47-17	White/Black	23 x 41 x 8	1.20	30-40	250	–
X2475*	White/Black	30½ x 32	1.40	20-25	250	–
211M	White/Blue	25 x 34	1.20	12-16	250	–
264	White/Blue	30½ x 41	1.40	20-30	250	–
X2276*	Mint Green/Black	29 x 45	1.20	26	240	–
X2326*	Mint Green/Black	38 x 55	1.20	60	100	–
285	Blue/Blue	30½ x 41	1.20	20-30	250	–
3095	Green/Black	40 x 46	1.00	40-45	200	–
51-40	Blue/Black	40 x 46	1.20	40-45	100	–
Print/Label: Soiled Linen (English/Spanish).						
3072	White/Black	30½ x 43	1.30	20-30	250	–
Star Seal. Print/Label: Soiled Linen.						
3102	Blue/Black	30½ x 41	1.00	20-30	250	–
Print/Label: Pediatric Soiled Linen.						
129M	Pink/Blue	30½ x 41	1.10	20-30	250	–
Print/Label: Surgery Soiled Linen.						
203D	Green/Black	30½ x 41	1.40	20-30	250	–

(continued on pg 23)

**Minimum order quantity required; please contact your distributor for details.
Measurements in inches / width x length x side gusset.*

Laundry and Linen Bags — LLDPE Film, Flat Pack (continued)

Code#	Bag/Print Color	Size	Gauge/mil	Capacity/gal	Qty/cs	Qty/pkg
Print/Label: Infectious Linen.						
208M	Yellow/Red	25 x 34	1.10	12-16	250	–
112M	Yellow/Red	30½ x 41	1.10	20-30	250	–
112	Yellow/Red	30½ x 41	1.30	20-30	250	–
170M	Yellow/Red	38 x 45	1.10	44	250	–
Print/Label: Infectious Linen, Biohazard, Biohazard Symbol.						
47-03	Red/Black	23 x 41 x 8	1.20	30-40	250	–
47-07	Yellow/Black	23 x 41 x 8	1.20	30-40	250	–
47-57	Yellow/Black	31 x 41	1.00	20-30	250	–
51-45	Yellow/Black	40 x 46	1.20	40-45	100	–
Print/Label: Reject/Torn Linen.						
215XH	Clear/Black	30½ x 41	1.10	20-30	250	–
Print/Label: Recyclable Material.						
RC101	Clear/Blue	30½ x 41	1.30	20-30	250	–
Print/Label: Contaminated Laundry.						
252XH	White/Red	30½ x 41	1.40	20-30	250	–
Print/Label: Contaminated Laundry, Biohazard Symbol.						
X2419*	Yellow/Black	38 x 50	2.00	44	100	–
No Print/Label.						
3100	White/None	20 x 22	1.50	7-10	250	–
9233	Yellow/None	25 x 34	4.00	12-16	100	–
301M	White/None	27 x 43	1.10	20-30	200	–
T662*	Mint Green/None	30 x 37	0.70	20-30	250	–
205M	Buff/None	30½ x 41	0.80	20-30	250	–
260M	Blue/None	30½ x 41	1.20	20-30	250	–
264M	White/None	30½ x 41	1.20	20-30	250	–
204XH	Clear/None	30½ x 41	1.30	20-30	250	–
201	Blue/None	30½ x 41	1.40	20-30	250	–
203	Green/None	30½ x 41	1.40	20-30	250	–
202A	White/None	30½ x 41	1.80	20-30	250	–
X2350*	Blue/None	30½ x 41	2.00	20-30	250	–
305	White/None	38 x 45	1.30	44	250	–
306	Blue/None	38 x 45	1.40	44	250	–

**Minimum order quantity required; please contact your distributor for details.
Measurements in inches / width x length x side gusset.*

Laundry and Linen Bags

Laundry and Linen Bags — LLDPE Film, Coreless Roll

Code#	Bag/Print Color	Size	Gauge/mil	Capacity/gal	Qty/cs	Qty/roll
Print/Label: Soiled Linen.						
R285HA	Light Blue/Blue	30½ x 41	1.10	20-30	250	25
R285	Light Blue/Black	30½ x 41	1.20	20-30	250	25
X2304*	Mint Green/Black	30½ x 41	1.20	20-30	250	25
R264	White/Blue	30½ x 41	1.50	20-30	250	25
Star Seal. Print/Label: Soiled Linen.						
X2029*	Dark Blue/Black	29 x 45	1.20	26	240	30
Print/Label: Infectious Linen, Biohazard, Biohazard Symbol.						
R112M	Yellow/Red	30½ x 41	1.10	20-30	250	25
R170M	Yellow/Red	38 x 45	1.10	44	150	15
No Print/Label. LDPE Film, Roll.						
R260M	Blue/None	30½ x 41	1.30	20-30	250	25

Laundry and Linen Bags — Draw Tape, LLDPE Film, Flat Pack

Code#	Bag/Print Color	Size	Gauge/mil	Capacity/gal	Qty/cs	Qty/pkg
No Print/Label.						
50-08	Blue/Black	31 x 41	1.50	20-30	250	–

Laundry and Linen Bags — Draw Tape, LLDPE Film, Coreless Roll

Code#	Bag/Print Color	Size	Gauge/mil	Capacity/gal	Qty/cs	Qty/roll
Print/Label: Soiled Linen.						
288	Blue/Black	30½ x 41	1.20	20-30	100	25

Laundry and Linen Bags — HDPE Film, Flat Pack

Code#	Bag/Print Color	Size	Gauge/mic	Capacity/gal	Qty/cs	Qty/pkg
Print/Label: Soiled Linen.						
H355*	Yellow/Black	30 x 54	14.00	30-35	200	–
X2429*	Blue/Black	30½ x 43	17.00	20-30	200	–
H563*	Blue/Black	31 x 45	16.00	20-30	250	–
Print/Label: Soiled Linen (English & Spanish); Do Not Fill Above This Line.						
3693*	Blue/Black	30 x 46	19.00	25-32	250	–
No Print/Label.						
H489*	Blue/None	38 x 46	17.00	44	250	–

*Minimum order quantity required; please contact your distributor for details.
Measurements in inches / width x length x side gusset.

Laundry and Linen Bags

Laundry and Linen Bags — HDPE Film, Coreless Roll

Code#	Bag/Print Color	Size	Gauge/mic	Capacity/gal	Qty/cs	Qty/roll
Print/Label: Soiled Linen.						
RS304314PB	Blue/Black	30 x 43	14.00	20-30	250	25
RS304316PB	Blue/Black	30 x 43	16.00	20-30	250	25
3056	Blue/Black	40 x 46	14.00	40-45	250	25
Print/Label: Infectious Linen.						
RS304314Y	Yellow/Black	30½ x 43	14.00	20-30	250	25
RS384616Y	Yellow/Black	38 x 46	16.00	44	250	25
No Print/Label.						
RS304314B	Blue/None	30½ x 43	14.00	20-30	250	25
RS404816B	Blue/None	40 x 48	16.00	40-45	250	25
Dispenser Box, Twist Ties, Star Seal. No Print/Label.						
44-27	Yellow/None	38 x 45	16.00	44	100	100



Anti-Static Kick Bucket Liners comply with Operating Room requirements of NFPA Standard 56A; bags contain a carbon additive to help prevent the buildup of electrostatic charges. Available in standard sizes with/without Soiled Linen labeling to fit most popular kick buckets.

Anti-Static Kick Bucket Liners — Meet NFPA 56A, LLDPE Film, Flat Pack

Code#	Bag/Print Color	Size	Gauge/mil	Capacity/gal	Qty/cs	Qty/pkg
Fits 18" Diameter. Print/Label: Anti-Static, Soiled Linen.						
3110AS	Blue/Black	16 x 44 x 12	1.40	30	100	–
Fits 18" Diameter. Print/Label: Anti-Static.						
51-00	Blue/Black	16 x 45 x 12	1.60	30	100	–
Fits 25" Diameter. Print/Label: Soiled Linen.						
3106AS	Blue/Black	18 x 16 x 8	1.25	5-7	500	–
Fits 25" Diameter. Print/Label: Anti-Static, Soiled Linen.						
3111AS	Blue/Black	22 x 45 x 16	1.40	35-40	100	–
Fits 25" Diameter. Print/Label: Anti-Static.						
51-20	Blue/Black	22 x 45 x 16	1.60	35-40	100	–
Print/Label: Anti-Static, Soiled Linen.						
3104AS	Blue/Black	12 x 24 x 8	1.25	4-7	500	–
Print/Label: Anti-Static.						
50-00	Blue/Black	12 x 24 x 8	1.30	4-7	500	–
Fits Oval Buckets. Print/Label: Anti-Static.						
50-01	Blue/Black	18 x 16 x 8	1.30	5-7	500	–



*Minimum order quantity required; please contact your distributor for details.
Measurements in inches / width x length x side gusset.

Water Soluble Laundry Bags

Water Soluble Laundry Bags are designed for the safe and effective collection and handling of potentially contaminated or infectious linen and clothing. The bag is impervious to bacteria and viruses and is gas impermeable, exhibiting excellent odor containment properties. Contaminated linen is placed into the bag, sealed and placed directly into the washing machine, eliminating the need for staff to come in direct contact with the contents.

- Allows for easy isolation of contaminated and soiled linens to help reduce risk of cross-contamination.
- 100% Biodegradable polyvinyl alcohol film dissolves completely into a non-toxic aqueous solution at 150°F during Hot Water machine washing and at 55°F for the Cold Water formulation.
- Superior wet strength holdout capacity.
- Integral tie strap for securing bag dissolves quickly to release contents during washing cycle.
- Bags supplied on easy to dispense coreless rolls and individually packaged zip closure bags.
- Sizes designed to fit most standard linen hampers.
- Typical areas of use: in most any treatment area or room where potentially infectious linens or clothing come in contact with bodily fluids.



Water Soluble Laundry Bags — Polyvinyl Alcohol Film, Coreless Roll

Code#	Bag/Print Color	Size	Gauge/mil	Capacity/gal	Qty/cs	Qty/pkg
Hot Water (150°F, 65°C), Heavy Duty. No Print/Label.						
1-344	Clear/None	19 x 22	1.00	7-10	100	25
1-341	Clear/None	26 x 33	1.00	20-25	100	25
1-347	Clear/None	28 x 39	1.00	30-32	100	25
1-348	Clear/None	36 x 39	1.00	40-45	100	25
Hot Water (150°F, 65°C), Standard Duty. No Print/Label.						
541-A	Clear/None	26 x 33	0.80	20-25	100	25
547-A	Clear/None	28 x 39	0.80	30-32	100	25
548-A	Clear/None	36 x 39	0.80	40-45	100	25
Cold Water (55°F, 13°C), Heavy Duty. No Print/Label.						
1-342	Clear/None	26 x 33	1.00	20-25	100	25
1-345	Clear/None	28 x 39	1.00	30-35	100	25
1-346	Clear/None	36 x 39	1.00	40-45	100	25

Measurements in inches / width x length x side gusset.

Hardware and Accessories

Hamper Stands facilitate convenient placement throughout the facility and ease of use. Stands are engineered for long-lasting performance in demanding health care environments with injection molded lid made from high impact ABS crack resistant resin to withstand repeated opening and closing.

- Hands free foot pedal actuation opens lid wide for easy access & closes automatically to help contain odors.
- Solid 7/8" Chrome plated steel tubing with lower support platform for heavyweight applications.
- Upper locking ring holds bags of at least 30 inches wide securely in place.
- Open front design provides convenient access and bag removal.
- Four 3" free turning casters—two with locking foot brakes—provide easy maneuverability.
- Hamper pouch provides convenient storage for additional bags.
- 19" wide x 21" deep x 36" tall. Available with or without labels.

Premium Chrome Hampers — Foot Pedal Activated Lid Operation

Code#	Description	Size	Capacity/gal	Qty/cs
15-9000	No Decal	19 x 36 x 21	30-33	1
15-9100	Soiled Linen Decal	19 x 36 x 21	30-33	1
15-9200	Infectious Waste Decal	19 x 36 x 21	30-33	1
15-9300	Isolation Linen Decal	19 x 36 x 21	30-33	1
15-9500	Surgery Soiled Linen Decal	19 x 36 x 21	30-33	1
15-9600	Trash Decal	19 x 36 x 21	30-33	1

Assorted Hardware and Accessories

Code#	Description	Capacity/gal	Qty/cs
15-24	Disposable Cardboard Hamper, yellow, no print, 17 x 25 x 15	20	10
15-9POUCH	Replacement Hamper Pouch, 14 1/2 x 14	—	1
15-9LIDE	Premier Lid, Plastic, 19 x 21		1



Measurements in inches / width x length x side gusset.

Institutional Trash Can Liners

Institutional Trash Can Liners provide a leak resistant containment solution for non-medical waste disposal and are used throughout the health care facility to collect and contain refuse. Available in both Linear Low Density and High Density/HMW Polyethylene film, and in flat pack or coreless rolls for easy dispensing and convenience to provide a cost-effective solution for every application.

- Available in a wide variety of standard sizes, film gauges, flat, star seal or gusseted bottom designed to meet your specific application. Custom specifications available — *minimum order quantity required.*
- Superior performance liners made from premium resins with reinforced seams are puncture resistant and impervious to liquids.
- Clear liners allow quick and easy visualization of contents before disposal while colored options hide unsightly refuse.
- Offered in both High Density/HMW and Linear Low Density resins for your specific application and budget.
- Linear Low Density liners provide maximum strength and puncture/tear resistance due to their superior elasticity, provide excellent odor containment properties and are best suited for handling a wide variety of refuse including those with straight corners/sharper edges.
- High Density/HMW liners offer excellent durability and weight capacity in a lighter gauge film that is most cost effective for refuse without sharp or jagged edges.
- Latex free and environmentally safe for incineration disposal.
- Easy dispensing in convenient flat pack or coreless rolls.
- Typical areas of use: throughout health care facility for refuse collection. Not for use with Biohazardous Waste.



Institutional Trash Can Liners — LLDPE Film, Flat Pack

Code#	Bag/Print Color	Size	Gauge/mil	Capacity/gal	Qty/cs	Qty/pkg
No Print/Label.						
3900	Clear/None	12 x 18	1.00	2-3	1000	–
550	Clear/None	18 x 24	1.50	5-7	1000	–
4621	Brown/None	23 x 29	0.50	8-10	500	–
4625	Clear/None	23 x 29	0.50	8-10	500	–
4608	Brown/None	24 x 23	0.35	8-10	1000	–
4609	Clear/None	24 x 23	0.35	8-10	1000	–
4610	Brown/None	24 x 23	0.45	8-10	1000	–
4613	Clear/None	24 x 23	0.45	8-10	1000	–
1292	Clear/None	24 x 24	0.80	8-10	500	–
X2418*	Clear/None	24 x 29	0.70	10-12	750	–
4618	Brown/None	24 x 32	0.40	12-16	1000	–

(continued on pg 29)

**Minimum order quantity required; please contact your distributor for details.
Measurements in inches / width x length x side gusset.*

Institutional Trash Can Liners

Institutional Trash Can Liners — LLDPE Film, Flat Pack (continued)

Code#	Bag/Print Color	Size	Gauge/mil	Capacity	Qty/cs	Qty/pkg
No Print/Label. (continued)						
4620	Clear/None	24 x 32	0.40	12-16	1000	–
T612*	Clear/None	24 x 32	1.00	12-16	500	–
4627	Clear/None	30 x 36	0.40	20-30	500	–
4628	Brown/None	30 x 36	0.40	20-30	500	–
3914	Brown/None	30 x 36	0.45	20-30	500	–
4637	Clear/None	30 x 36	0.65	20-30	250	–
4630	Brown/None	30 x 36	0.65	20-30	250	–
X2211*	White/None	30 x 36	0.80	20-30	250	–
4664	Brown/None	33 x 39	0.50	33	500	–
4644	Brown/None	33 x 39	0.65	33	250	–
4646	Clear/None	33 x 39	0.65	33	250	–
T578*	White/None	33 x 39	0.75	33	250	–
4648	Brown/None	33 x 39	0.90	33	250	–
1286	Clear/None	33 x 39	1.20	33	250	–
428A	Clear/None	33 x 39	2.00	33	250	–
23-35C	Clear/None	33 x 40	1.25	33	100	–
X2422*	Black/None	34 x 44	3.00	33	100	–
X2413*	Gray/None	36 x 46	1.30	33	150	–
T274*	Clear/None	36 x 53	2.00	50	100	–
T336*	Clear/None	36 x 54	5.00	50	50	–
23-55C	Clear/None	36 x 58	0.65	55	100	–
1287	Clear/None	36 x 58	1.20	55	200	–
4650	Brown/None	38 x 45	0.70	44	250	–
4651	Brown/None	38 x 45	0.95	44	250	–
4659	Clear/None	38 x 45	0.95	44	250	–
4669	Brown/None	38 x 45	1.40	44	100	–
642	Clear/None	38 x 45	1.50	44	100	–
T413*	Clear/None	38 x 45	1.65	44	40	–
1295	Clear/None	38 x 50	1.20	44	200	–
4661	Brown/None	38 x 55	0.80	60	200	–
3908	Brown/None	38 x 58	0.50	60	200	–
1291	Clear/None	38 x 58	0.80	60	200	–
3910	Brown/None	38 x 58	0.90	60	100	–
3911	Clear/None	38 x 58	0.90	60	100	–
23-65C	Clear/None	38 x 58	1.20	60	100	–
23-40C	Clear/None	40 x 46	0.45	40-45	250	–
26-15	Black/None	40 x 46	0.45	40-45	250	–
3905	Clear/None	40 x 46	0.50	40-45	250	–
3906	Clear/None	40 x 46	0.60	40-45	250	–

(continued on pg 30)

**Minimum order quantity required; please contact your distributor for details.
Measurements in inches / width x length x side gusset.*

Institutional Trash Can Liners



Institutional Trash Can Liners — LLDPE Film, Flat Pack *(continued)*

Code#	Bag/Print Color	Size	Gauge/mil	Capacity	Qty/cs	Qty/pkg
No Print/Label. <i>(continued)</i>						
4662	Brown/None	40 x 46	0.60	40-45	250	–
4660	Clear/None	40 x 46	0.70	40-45	200	–
23-50C	Clear/None	40 x 46	1.10	40-45	100	–
X2306*	Black/None	40 x 46	1.20	40-45	100	–
T598*	Clear/None	40 x 46	1.45	40-45	125	–
T706*	Gray/None	40 x 48	1.20	40-45	100	–
1111	Clear/None	40 x 48	1.35	40-45	100	–
3916	Black/None	40 x 48	1.50	40-45	200	–
1289	Brown/None	40 x 48	1.80	40-45	100	–
23-49W	White/None	43 x 47	0.65	55	100	–
3907	Brown/None	43 x 47	0.80	55	100	–
X2414*	Dark Gray/None	43 x 47	1.30	55	100	–
4666	Clear/None	43 x 48	1.50	56	100	–
4665	Brown/None	43 x 48	1.50	56	100	–
X2269*	Buff/None	43 x 48	1.65	56	100	–
X2303*	Gray/None	44 x 68	2.00	60	50	–
X2430*	Clear/None	45 x 52	1.30	55	100	
1284	Clear/None	46 x 50	0.80	55	100	–
T588*	Clear/None	49 x 50	1.30	55	150	–
X2041*	Clear/None	50 x 48	2.80	55	50	–
X1971*	Black/None	52 x 60	3.00	60	50	–
X2489*	Clear/None	52 x 60	3.00	60	50	–
Star Seal. No Print/Label.						
26-21	Black/None	46 x 50	1.00	55	100	–
Drawstring. No Print/Label.						
3659	Clear/None	40 x 40	2.00	40-45	100	–
Side Gusset. No Print/Label.						
23-01C	Clear/None	13 x 17 x 4	0.35	3-4	1000	–
T583*	Clear/None	30 x 41 x 28	1.00	35-45	100	–
T456*	Gray/None	38 x 50 x 27	1.30	65	100	–
Twist Ties. No Print/Label.						
X805*	White/None	40 x 55	3.00	55	100	–

**Minimum order quantity required; please contact your distributor for details.
Measurements in inches / width x length x side gusset.*

Institutional Trash Can Liners

Institutional Trash Can Liners — LLDPE Film, Coreless Roll

Code#	Bag/Print Color	Size	Gauge/mil	Capacity/gal	Qty/cs	Qty/pkg
No Print/Label.						
R4625	Clear/None	23 x 29	0.50	8-10	500	50
R4608	Brown/None	24 x 23	0.35	8-10	500	50
R4613	Clear/None	24 x 23	0.45	8-10	500	50
R4627	Clear/None	30 x 36	0.45	20-30	500	50
R4630	Brown/None	30 x 36	0.65	20-30	125	25
R4637	Clear/None	30 x 36	0.65	20-30	125	25
R4644	Brown/None	33 x 39	0.65	33	125	25
R4646	Clear/None	33 x 39	0.65	33	125	25
R4648	Brown/None	33 x 39	0.90	33	125	25
T498*	Clear/None	33 x 39	1.50	33	200	25
T580*	White/None	37 x 46	0.70	44	250	25
R4651	Brown/None	38 x 45	0.95	44	80	20
R4659	Clear/None	38 x 45	0.95	44	80	20
R4662	Brown/None	40 x 46	0.60	40-45	100	25
R4660	Clear/None	40 x 46	0.70	40-45	80	20
R4666	Clear/None	43 x 48	1.50	56	60	15
Cored Roll. No Print/Label.						
T412*	Clear/None	33 x 39	0.80	33	100	–
Star Seal. No Print/Label.						
1212	Clear/None	38 x 58	0.80	60	200	25

**Minimum order quantity required; please contact your distributor for details.
Measurements in inches / width x length x side gusset.*

Institutional Trash Can Liners



Institutional Trash Can Liners — HDPE Film, Flat Pack

Code#	Bag/Print Color	Size	Gauge/mic	Capacity/gal	Qty/cs	Qty/pkg
No Print/Label.						
HBF24246	Buff/None	24 x 24	6.00	8-10	2000	–
HCL24246	Clear/None	24 x 24	6.00	8-10	2000	–
HCL24248	Clear/None	24 x 24	8.00	8-10	2000	–
HCL24336	Clear/None	24 x 33	6.00	12-16	2000	–
HCL24338	Clear/None	24 x 33	8.00	12-16	1000	–
HCL303710	Clear/None	30 x 37	10.00	20-30	500	–
HWT304314	White/None	30 x 43	14.00	25-30	500	–
HBF334011	Buff/None	33 x 40	11.00	33	500	–
HCL334011	Clear/None	33 x 40	11.00	33	500	–
HCL334016	Clear/None	33 x 40	16.00	33	250	–
HCL366016	Clear/None	36 x 60	16.00	55	200	–
HCL386016	Clear/None	38 x 60	16.00	60	200	–
HCL404812	Clear/None	40 x 48	12.00	40-45	250	–
HCL404814	Clear/None	40 x 48	14.00	40-45	250	–
HCL404816	Clear/None	40 x 48	16.00	40-45	250	–
HBF434814	Buff/None	43 x 48	14.00	56	200	–
HCL434812	Clear/None	43 x 48	12.00	56	200	–
Star Seal. No Print/Label.						
29-15	Black/None	33 x 40	22.00	33	200	–
29-20	Clear/None	38 x 58	14.00	60	200	–
29-21	Clear/None	38 x 58	17.00	60	200	–
29-22	Black/None	38 x 58	22.00	60	150	–
29-11	Clear/None	40 x 48	11.00	40-45	250	–
29-45C	Clear/None	40 x 48	13.00	40-45	250	–
29-27	Black/None	43 x 48	22.00	56	150	–

Institutional Trash Can Liners — HDPE Film, Coreless Roll

Code#	Bag/Print Color	Size	Gauge/mic	Capacity/gal	Qty/cs	Qty/roll
Star Seal. No Print/Label.						
RS202206N	Clear/None	20 x 22	6.00	7-10	2000	50
RS242406N	Clear/None	24 x 24	6.00	8-10	1000	50
RS242408N	Clear/None	24 x 24	8.00	8-10	1000	50
RS243306N	Clear/None	24 x 33	6.00	12-16	1000	50
1011	Clear/None	24 x 33	6.00	12-16	1000	50
RS243308N	Clear/None	24 x 33	8.00	12-16	1000	50
H462R*	Dark Green/None	24 x 33	18.00	12-16	750	50
RS303708N	Clear/None	30 x 37	8.00	20-30	500	25

(continued on pg 33)

*Minimum order quantity required; please contact your distributor for details.
Measurements in inches / width x length x side gusset.

Institutional Trash Can Liners

Institutional Trash Can Liners — HDPE Film, Coreless Roll (continued)

Code#	Bag/Print Color	Size	Gauge/mic	Capacity/gal	Qty/cs	Qty/roll
Star Seal. No Print/Label. (continued)						
RS303710N	Clear/None	30 x 37	10.00	20-30	500	25
RS303713N	Clear/None	30 x 37	13.00	20-30	500	25
RS303716N	Clear/None	30 x 37	16.00	20-30	500	25
RS334009N	Clear/None	33 x 40	9.00	33	500	25
RS334011N	Clear/None	33 x 40	11.00	33	500	25
RS334013N	Clear/None	33 x 40	13.00	33	500	25
RS334016N	Clear/None	33 x 40	16.00	33	250	25
RS334017N	Clear/None	33 x 40	17.00	33	250	25
RS334017W	White/None	33 x 40	17.00	33	250	25
H461R*	Dark Green/Black	33 x 40	21.00	33	200	20
RS366017N	Clear/None	36 x 60	17.00	55	200	25
RS374612N	Clear/None	37 x 46	12.00	44	250	25
RS374614N	Clear/None	37 x 46	14.00	44	250	25
RS374616N	Clear/None	37 x 46	16.00	44	250	25
RS386012N	Clear/None	38 x 60	12.00	60	200	25
RS386014N	Clear/None	38 x 60	14.00	60	200	25
RS386017N	Clear/None	38 x 60	17.00	60	200	25
RS386018W	White/None	38 x 60	18.00	60	200	25
RS386019K	Black/None	38 x 60	19.00	60	200	25
RS404810N	Clear/None	40 x 48	10.00	40-45	250	25
RS404812N	Clear/None	40 x 48	12.00	40-45	250	25
RS404814N	Clear/None	40 x 48	14.00	40-45	250	25
RS404816N	Clear/None	40 x 48	16.00	40-45	250	25
RS404817N	Clear/None	40 x 48	17.00	40-45	250	25
RS404818W	White/None	40 x 48	18.00	40-45	200	25
RS404819K	Black/None	40 x 48	19.00	40-45	200	25
H570*	Clear/None	40 x 48	21.00	40-45	200	25
RS404822N*	Clear/None	40 x 48	22.00	40-45	150	25
RS404822K	Black/None	40 x 48	22.00	40-45	150	25
RS434812N	Clear/None	43 x 48	12.00	56	200	25
RS434814N	Clear/None	43 x 48	14.00	56	200	25
RS434816N	Clear/None	43 x 48	16.00	56	200	25
RS434817N	Clear/None	43 x 48	17.00	56	200	25
1729	Black/None	43 x 48	22.00	56	150	25

*Minimum order quantity required; please contact your distributor for details.
Measurements in inches / width x length x side gusset.

Stackable Sharps Containers



Our complete line of **Stackable and Non-Stackable Sharps Containers** is offered in a variety of sizes and configurations designed to give you the ability to choose the right product for your specific application. Regulated as Class II Medical devices and manufactured in accordance with FDA Quality System Regulations, as well as the International Quality System Standards ISO 13485:2003; all Medical Action Sharps Containers meet or exceed OSHA and JCAHO Standards, CDC/NIOSH guidelines, ASTM requirements and DOT Regulations including 49CFR173.24, 49CFR173.24a, and 49CFR173.134.

- **Puncture and Leak Resistant** – Strong polypropylene construction meets or exceeds industry standards.
- **Locking Lids** – Help keep contents secure for final disposal.
- **Autoclave/Incinerator Safe** – Containers are Autoclave Safe up to a maximum of 280°F (137°C). When properly incinerated, containers emit only carbon dioxide and water.
- Optional Wall Mounting Brackets and large container cart with locking casters.
- Typical areas of use: Operating Room, Emergency Room, Same Day Surgery, Primary Care Treatment areas, Patient Room, Phlebotomy Labs, Chemistry & Hematology, Pathology, Toxicology, Microbiology, Immunization, Drug Screening, or wherever Biohazardous Sharps exist.

Stackable Sharps Containers

- Stackable design provides a convenient stable platform during use, economical storage, and lower shipping cost.
- Wide horizontal drop lid models accommodate larger syringes and allow full utilization of fill line capacity.
- Tortuous Path Lids help limit hand exposure for added security.
- Transparent lids facilitate quick visual identification of fill level.
- Translucent containers facilitate easy side viewing of content level.
- Compact containers are conveniently sized for easy transport.
- Extra large containers feature a wide hinged lid and round funneled vertical drop to accommodate endoscopic accessories, elongated disposable instruments and probes.

Stackable Sharps Containers



Stackable Sharps Containers — Polypropylene

Code#	Container/Print Color	Size	Capacity	Qty/cs
Small. Print/Label: Biohazard, Biohazard Symbol.				
8702	Red/Black	3½ x 3½ x 7	1 qt	72
8702T	Translucent Red/Black	3½ x 3½ x 7	1 qt	72
8702TY	Yellow/Black	3½ x 3½ x 7	1 qt	72
Medium. Print/Label: Biohazard, Biohazard Symbol.				
8703	Red/Black	10 x 7 x 5	1 gal	24
8703T	Translucent Red/Black	10 x 7 x 5	1 gal	24
8703TY	Yellow/Black	10 x 7 x 5	1 gal	24
Medium. Tortuous Path Lid. Print/Label: Biohazard, Biohazard Symbol.				
8706	Red/Black	10 x 7 x 6¾	1 gal	12
8706T	Translucent Red/Black	10 x 7 x 6¾	1 gal	12
8706TY	Yellow/Black	10 x 7 x 6¾	1 gal	12
Large. Print/Label: Biohazard, Biohazard Symbol.				
8704	Red/Black	10 x 7 x 9½	2 gal	24
8704T	Translucent Red/Black	10 x 7 x 9½	2 gal	24
8704TY	Yellow/Black	10 x 7 x 9½	2 gal	24
Large. Tortuous Path Lid. Print/Label: Biohazard, Biohazard Symbol.				
8707	Red/Black	10 x 7 x 11¼	2 gal	10
8707T	Translucent Red/Black	10 x 7 x 11¼	2 gal	10
8707TY	Yellow/Black	10 x 7 x 11¼	2 gal	10
8710	Red/Black	12½ x 6 x 13½	3 gal	12
Large. Flat Lockable Sliding Lid. Print/Label: Biohazard, Biohazard Symbol.				
8710TF	Red/Black	12½ x 6 x 12	3 gal	12
Extra Large. Print/Label: Biohazard, Biohazard Symbol.				
8705	Red/Black	13¾ x 13¾ x 14	8 gal	10
8705T	Translucent Red/Black	13¾ x 13¾ x 14	8 gal	10
8705TY	Yellow/Black	13¾ x 13¾ x 14	8 gal	10
8716	Red/Black	13½ x 13½ x 25½	16 gal	6

Chemo Waste (White) models also feature:

- Chemotherapy Waste Labeling/CYTA Symbol for easy identification.
- Transparent top facilitates quick visualization of fill level.

Stackable Sharps Containers — Polypropylene, Chemo Waste Container

Code#	Container/Print Color	Size	Capacity/gal	Qty/cs
Extra Large. Print/Label: CYTA Symbol, Chemotherapy Waste.				
9751	White/Black & Yellow	13¾ x 13¾ x 14	8	10
9752	White/Black & Yellow	13½ x 13½ x 25½	16	6



Measurements in inches / width x depth x height.

Stackable Sharps Containers



Stackable Sharps Container Hardware — White

Code#	Description	Size	Qty/cs
8780	Wall Mount for 8702 series	–	1
8782	Locking Wall Mount for 8703, 8704 series	–	1
8783	Locking Wall Mount for 8706, 8707 series	–	1
8784	Locking Wall Mount for 8705 series	–	1
8785	Locking Wall Mount for 8706, 8707 series	–	1
8786F	Locking Bracket for 3 Gallon 8710TF	–	1
8790	Cart w/Locking Casters for 8705 & 8716	13¼ x 12½ x 15¼	1

Non-Stackable Sharps Containers



Standard (Red) models feature:

- Funneled vertical drop and Extended Neck design help minimize needle exposure and overfilling.
- Durable Polypropylene construction provides a stable base during use.
- Extended handles help keep fingers away from opening during transport.
- Needle removal port (#158, #178) provides convenient and compact disposal of sharps.

Non-Stackable Sharps Containers — Polypropylene, Red

Code#	Description	Size	Capacity/qt	Qty/cs
Print/Label: Biohazard, Biohazard Symbol.				
158	Short Tray Size w/Needle Remover	3½ x 3½ x 5	0.7	40
185S	Short Tray Size, Locking Cap	3½ x 3½ x 5	0.7	40
178	Tall Tray Size w/Needle Remover	3½ x 3½ x 10	1.7	20
185	Tall Tray Size, Locking Cap	3½ x 3½ x 10	1.7	20
184	Small Size, Locking Cap	10 x 6 x 9	4.8	12
182	Medium Size, Locking Cap	10 x 6 x 11	6.2	12
180	Large Size, Locking Cap	10 x 6 x 15	10.0	12
184R	Rectangular, Locking Cap	10¼ x 6 x 7¾	4.7	12
186W	Large Wide Opening, Locking Cap	12 x 9 x 11	12.0	12
187	Extra Large Wide Opening, Locking Cap	12 x 9 x 18	23.0	4

Measurements in inches / width x depth x height.

Non-Stackable Sharps Containers



Translucent models also feature:

- Wrap around window view label provides quick visualization of fill level.
- Leak Resistant Locking Screw Cap models can be easily opened/closed during use and provide extra measure of security when fully engaged in locked position.
- Standard Locking Cap closure on all wide mouth containers.

Non-Stackable Sharps Containers — Polypropylene, Translucent

Code#	Description	Size	Capacity/qt	Qty/cs
Print/Label: Biohazard, Biohazard Symbol.				
199	Short Tray Size, Locking Screw Cap	5 x 4 x 7½	1.3	40
195	Tall Tray Size, Locking Screw Cap	5 x 4 x 10	2.0	20
184S	Small Size, Locking Screw Cap	10 x 6 x 9	4.8	12
182S	Medium Size, Locking Screw Cap	10 x 6 x 11	6.2	12
180S	Large Size, Locking Screw Cap	10 x 6 x 15	10.0	12
186S	Large Wide Opening, Locking Cap	12 x 9 x 11	12.0	12
187S	Extra Large Wide Opening, Locking Cap	12 x 9 x 18	23.0	4

Non-Stackable Sharps Container Accessories

Code#	Description	Qty/cs
R175	Dovetail Bracket for use with 185, 178, 158, 185S	12
R275	Locking Mounting Bracket for use with 195, 199	12
R279	Basket for Holding 186S, 186W, 187, 187S	4
R299E	Basket for Holding 180, 180S, 180T, 182, 182S, 184, 184S	1
191	Cable, Lock, Key Assembly for use with R279, R299E	1

Chemo Waste (White) models feature:


- Chemotherapy Waste Labeling/ CYTA Symbol for easy identification.

Chemotherapy Waste Containment — Polypropylene, White

Code#	Description	Size	Capacity/qt	Qty/cs
Print/Label: CYTA Symbol, Chemotherapy Waste.				
930	Small Size, Locking Cap	3½ x 3½ x 10	1.7	20
931W	Standard Size, Large Opening, Locking Cap	12 x 9 x 11	12.0	12
932	Large Size, Large Opening, Locking Cap	12 x 9 x 18	23.0	4



Measurements in inches / width x depth x height.



Code#	Description	Qty/cs
3775*	Umbrella Bag, 7 x 27 x 4, 0.70 mil	1000
3802	Umbrella Stand for #3775 Bag	1



Medical Action
INDUSTRIES INC.
Right. From the Beginning.

Minor Procedure Kits & Trays
Operating Room Supplies
Patient Bedside Essentials
Collection & Containment Solutions
Laboratory Products
Medical Products

800 Prime Place
Hauppauge, NY 11788
1.800.645.7042
www.Medical-Action.com