



Urinalysis Controls

BIO-RAD



Urinalysis Controls from Bio-Rad

Versatile packaging and comprehensive analytes for today's laboratory

Urinalysis is widely documented as the first clinical laboratory test developed and dates back over 6,000 years. Even today, modern urinalysis is still one of the most widely performed tests in routine health screening.

Bio-Rad Urinalysis Controls are available in a variety of convenient packaging options, including dropper bottles, flip-top squeeze bottles and centrifuge tubes.

Covering an extensive menu of analytes to monitor physical, chemical and microscopic examination of urine specimens, Bio-Rad Urinalysis controls are well suited for manual or automated methods for dipstick or microscopic analysis of urine specimens.

Participating in the Unity™ Interlaboratory Program provides valuable peer group comparisons with distinctive, easy-to-read color coded reports.



Liquichek™ Urinalysis Control

Specifically designed for urine dipstick and microscopic testing. The control is available in 12 mL glass vials.

- Liquid, human urine based
- Designed for dipstick and microscopic analysis
- 30 day open-vial stability at 2–25°C
- 2.5 year shelf life at 2–8°C

Analytes

Bilirubin	Ketones	pH
Blood	Leukocytes	Pregnancy (hCG)
Casts	Microalbumin	Protein (Total)
Clarity	Microscopics	Protein-to-Creatinine Ratio
Color	(RBC, WBC, Crystals)	Specific Gravity
Creatinine	Nitrite	Urobilinogen
Glucose	Osmolality	

Refer to the package insert of currently available lots for specific analyte and stability claims



qUAntify® Control

Designed to monitor the precision of urine dipstick testing. The control is available in 12 mL ready-to-use dropper bottles.

- Liquid
- Designed for dipstick testing
- 2 year shelf life at 2–8°C
- 31 day open-vial stability at 2–25°C
- Color-coded peer group reports available

Analytes

Albumin	Glucose	pH
Albumin-to-Creatinine Ratio	Hemoglobin	Pregnancy (hCG)
Ascorbic Acid*	Ketones	Protein (Total)
Bilirubin	Leukocytes	Protein-to-Creatinine Ratio
Blood	Microalbumin	Specific Gravity
Creatinine	Nitrite	Urobilinogen

* Tested and found to be negative.

Refer to the package insert of currently available lots for specific analyte and stability claims



qUAntify® Plus Control

A combination urine dipstick and microscopic control. The control is available in either 12 mL centrifuge tubes or 120 mL flip-top squeeze bottles.

- Liquid
- Designed for dipstick and microscopic analysis
- 2 year shelf life at 2–8°C
- 31 day open-vial stability at 2–25°C
- Color-coded peer group reports available

Analytes

Albumin	Hemoglobin	pH
Albumin-to-Creatinine Ratio	Ketones	Pregnancy (hCG)
Ascorbic Acid*	Leukocytes	Protein (Total)
Bilirubin	Microalbumin	Protein-to-Creatinine Ratio
Blood	Microscopics	Specific Gravity
Creatinine	(RBC, WBC, Crystals)	Urobilinogen
Glucose	Nitrite	

* Tested and found to be negative.

Refer to the package insert of currently available lots for specific analyte and stability claims

Ordering Information

Cat #	Product Name	Quantity
Liquichek™ Urinalysis Control		
435	Bilevel (6 of each level)	12 x 12 mL
436	Level 1	12 x 12 mL
437	Level 2	12 x 12 mL
435X	Bilevel MiniPak (1 of each level)	2 x 12 mL

Cat #	Product Name	Quantity
qUAntify® Control		
975	Bilevel (3 of each level)	6 x 12 mL
975X	Bilevel MiniPak (1 of each level)	2 x 12 mL

qUAntify® Plus Control		
995	Bilevel (5 of each level)/Tube	10 x 12 mL
995X	Bilevel MiniPak (1 of each level)/Tube	2 x 12 mL
962	Bilevel (2 of each level)/Bottle	4 x 120 mL
963	Level 1/Bottle	4 x 120 mL
964	Level 2/Bottle	4 x 120 mL
962X	Bilevel MiniPak (1 of each level)/Bottle	2 x 120 mL



These controls are featured in the Unity™ Interlaboratory Program. Participation in the Unity™ Interlaboratory Program can help improve the quality of your laboratory test results and help you meet regulatory and accreditation requirements worldwide. You may also subscribe to one of our Unity™ QC Data Management Solutions designed to improve the effectiveness of your statistical process control. Ask your Bio-Rad Account Representative for more information.

Discover the power of Unity™ Solutions at www.QCNet.com.



**Bio-Rad
Laboratories**

**For further information, please contact the Bio-Rad office nearest
you or visit our website at www.bio-rad.com/qualitycontrol**

*Clinical
Diagnostics Group*

Website www.bio-rad.com/qualitycontrol **U.S.** 1-800-2BIO-RAD **Australia** 61-2-9914-2800 **Austria** 43-1-877-8901 **Belgium** 32-9-385-5511 **Brazil** 5521-3237-9400
Canada 1-514-334-4372 **China** 86-21-64260808 **Czech Republic** 420-241-430-532 **Denmark** +45-4452-1000 **Finland** 358-9-804-22-00 **France** 33-1-47-95-60-00
Germany +49-(0)89-318-840 **Greece** 30-210-7774396 **Hong Kong** 852-2789-3300 **Hungary** +36-1-459-6100 **India** 91-124-4029300 **Israel** 972-3-9636050
Italy +39-02-216091 **Japan** 81-3-6361-7070 **Korea** 82-2-3473-4460 **Mexico** 52(55)5200-0520 **The Netherlands** +31-318-540666 **New Zealand** 64-9-415-2280
Norway 47-23-38-41-30 **Poland** 48-22-3319999 **Portugal** 351-21-472-7700 **Russia** 7-495-721-14-04 **Singapore** 65-6415-3188 **South Africa** 27-11-442-85-08
Spain 34-91-590-5200 **Sweden** 46-8-555-127-00 **Switzerland** 41-61-717-95-55 **Thailand** 662-651-8311 **United Kingdom** +44-(0)20-8328-2000