

# **BSN** medical **Product Catalog**

The Solution That Fits Your Life



Global Leaders in Woundcare, Orthopedics and Orthopedic Soft Good Products



#### BSN medical is your One Source for Total Care by offering:

- A systematic approach to casting and splinting with industry leading products in all catagory segments for sprains, strains and fracture care.
  - Ortho-glass. America's #1 splinting system
  - Specialist. America's #1 line of plaster products
  - #2 market share position in cast tape and cast padding
- High quality wound care, bandaging, orthopedic soft goods and anti-embolism product lines
- Standardization to a single source supplier which can leverage national contracts or local agreements to drive lower product and operational costs
- Specialized sales representatives that provide commercial solutions, training, education and CEUs
- Access to products through your distributor of choice

Splinting		Bandaging	
Ortho-Glass® Splinting System	4	Co-Plus <sup>®</sup> LF Latex-Free Cohesive Bandage	2
SafetySplint <sup>™</sup> Unpadded Roll	5	Ortho-Glass* Wraps	2
Ortho-Glass* Unpadded Precut Splinting System	5	BSN Elastic Bandage Latex-Free Elastic Bandage	2
Ortho-Glass® Strap Kits	6	Isoband <sup>®</sup> Elastic Light Support Bandage	2
Ortho-Glass® Ankle Treatment System	6	Tensogrip <sup>*</sup> Tubular Support Bandage	3
OCL® Splint Roll - 2 Plaster Splinting System	7	Versagrip* Latex-Free Tubular Support Bandage	3
Specialist® J-Splint® Plaster Roll Immobilizer	7	Tensoshape <sup>*</sup> Elastic Tubular Bandage	3
Plaster Roll Immobilizer Dispenser	8	Leukotape P Rigid Strapping Tape	3
Gypsona <sup>*</sup> Endless Slab Dispenser Pack	8	Lightplast Pro Elastic Adhesive Bandage	3
		Optiplaste <sup>*</sup> Elastic Adhesive Bandage	3
Plaster		Tensoplast* Elastic Adhesive Bandage	3
Specialist Plaster of Paris Bandages & Splints	9		
Orthoflex® Elastic Plaster of Paris Bandages	9	Thermoplastic Splints	
Gypsona® S Plaster of Paris Bandages & Splints	10	Specialist Fracture Bracing	3
Gypsona® HP Plaster of Paris Bandages & Splints	10	Orthoplast* Splint	3
OCL* Plaster of Paris Bandages & Splints	11	Orthoplast* II Splint	3
Synthetic Cast Tape		Woundcare	
Delta-Lite* Conformable Fiberglass Cast Tape	12	Coverlet Fabric Adhesive Bandages	3
	13	Coverlet <sup>®</sup> Eye Occlusors	3
Delta-Lite* Plus Fiberglass Cast Tape	13 14	Elastomull <sup>*</sup> Elastic Gauze Bandage	3
Delta-Cast* Conformable Polyester Cast Tape	15	Cover-Roll Stretch Non-Woven Adhesive Bandages	3
FlashCast Elite Printed Polyester Cast Tape  Delta-Cast Soft Semi-Rigid Cast Tape	15	Cuticell* Ointment Dressing	3
Delta-Cast 30ft Settil-Rigid Cast Tape	15	Leukomed Absorbent Wound Dressing	3
Padding & Stockinette		Leukomed T Transparent Wound Dressing	3
	16	Leukomed' T Plus Absorbent Film Wound Dressing	3
Delta-Dry Water Resistant Cast Padding	17	· ·	
Delta-Rol* Synthetic Cast Padding	17		
Specialist Cotton Blend Cast Padding	17	Anti-Em / GP	
Specialist 100 Cotton Cast Padding		Anti-Em / GP* Seemless Elastic Stockings with Gradient Compression	4
Sof-Rol' Absorbent Cast Padding  Polta Not' Orthogoadia Symbolic Stockingtto	18		
Delta-Net* Orthopaedic Synthetic Stockinette	18 18	Unna Boot	
Specialist Orthopaedic Cotton Stockinette	19	Gelocast <sup>*</sup> Original Unna's Boot	4
Terry-Net* C Cotton Terry Cloth Stockinette	19	5	
Terry-Net* S Synthetic Terry Cloth Stockinette	19	Orthopaedic Soft Goods	
Terry-Net* Thumb Spica Precut Synthetic Terry Cloth Stockinette  Delta Terry-Net* Adhesive Fleece Edger	20	Upper Extremities	
· ·	20	Miami J <sup>*</sup> Cervical Collar	4
Delta Terry-Net™ Adhesive Terry Cloth / Foam Padding	20	Miami JR' Cervical Collar	4
Delta Terry-Net <sup>™</sup> Adhesive Felt Padding		Miami JTO* Thoracic Extension	4
	21		
Delta Terry-Net <sup>™</sup> Adhesive Hook & Stretch Loop Closures			4
		Philidelphia <sup>*</sup> Collar	4
Cast Removal & Accessories			4 4
Cast Removal & Accessories FCT Starter Kit Featuring Delta-Cast Conformable	22	Philidelphia <sup>°</sup> Collar Vertebrace <sup>°</sup> Extrication Collar Actimove <sup>°</sup> Cervical Collar	4 4 4
Cast Removal & Accessories  FCT Starter Kit Featuring Delta-Cast* Conformable  FCT Accessories Kit	22 22	Philidelphia <sup>*</sup> Collar Vertebrace <sup>*</sup> Extrication Collar	4 4 4 4
Cast Removal & Accessories  FCT Starter Kit Featuring Delta-Cast Conformable  FCT Accessories Kit  American Orthopedic Cast Cutter Cast Saws	22 22 23	Philidelphia Collar Vertebrace Extrication Collar Actimove Cervical Collar Actimove Cervical Comfort Cervical Collar Universal Cervical Collar	4 4 4 4 4
Cast Removal & Accessories  FCT Starter Kit Featuring Delta-Cast* Conformable  FCT Accessories Kit  American Orthopedic** Cast Cutter Cast Saws  American Orthopedic** Cast Dust Vacuum	22 22 23 23	Philidelphia Collar Vertebrace Extrication Collar Actimove Cervical Collar Actimove Cervical Comfort Cervical Collar	4 4 4 4 4 4
Cast Removal & Accessories  FCT Starter Kit Featuring Delta-Cast Conformable  FCT Accessories Kit  American Orthopedic Cast Cutter Cast Saws  American Orthopedic Cast Dust Vacuum  American Orthopedic Portable Vacuum	22 22 23 23 24	Philidelphia Collar Vertebrace Extrication Collar Actimove Cervical Collar Actimove Cervical Comfort Cervical Collar Universal Cervical Collar Sized Cervical Collar Universal Cervical Collar	4 4 4 4 4 4
Cast Removal & Accessories  FCT Starter Kit Featuring Delta-Cast' Conformable  FCT Accessories Kit  American Orthopedic" Cast Cutter Cast Saws  American Orthopedic" Cast Dust Vacuum  American Orthopedic" Portable Vacuum  TurboCare" Cast Removal System	22 22 23 23 24 24	Philidelphia Collar Vertebrace Extrication Collar Actimove Cervical Collar Actimove Cervical Comfort Cervical Collar Universal Cervical Collar Sized Cervical Collar Universal Cervical Collar Actimove Clavical Collar	4 4 4 4 4 4 4
Cast Removal & Accessories  FCT Starter Kit Featuring Delta-Cast Conformable  FCT Accessories Kit  American Orthopedic Cast Cutter Cast Saws  American Orthopedic Cast Dust Vacuum  American Orthopedic Portable Vacuum  TurboCare Cast Removal System  Cast Saw Blades	22 22 23 23 24 24 25	Philidelphia Collar Vertebrace Extrication Collar Actimove Cervical Collar Actimove Cervical Comfort Cervical Collar Universal Cervical Collar Sized Cervical Collar Universal Cervical Collar	4 4 4 4 4 4 4 4
Cast Removal & Accessories  FCT Starter Kit Featuring Delta-Cast Conformable  FCT Accessories Kit  American Orthopedic Cast Cutter Cast Saws  American Orthopedic Cast Dust Vacuum  American Orthopedic Portable Vacuum  TurboCare Cast Removal System  Cast Saw Blades  BSN Buckets	22 22 23 23 24 24 25 25	Philidelphia' Collar Vertebrace' Extrication Collar Actimove' Cervical Collar Actimove' Cervical Comfort Cervical Collar Universal Cervical Collar Sized Cervical Collar Universal Cervical Collar Universal Cervical Collar Confort' Clavicle Support Comfort' Clavicle Brace Hook & Loop Closure Clavicle Support w/Four-Way Strap, Buckle Closures	4 4 4 4 4 4 4 4 4 4
Cast Removal & Accessories  FCT Starter Kit Featuring Delta-Cast* Conformable  FCT Accessories Kit  American Orthopedic** Cast Cutter Cast Saws  American Orthopedic** Cast Dust Vacuum  American Orthopedic** Portable Vacuum  TurboCare** Cast Removal System  Cast Saw Blades  BSN Buckets  Casting Instruments For the Orthopedic Specialist	22 22 23 23 24 24 25 25 26	Philidelphia Collar Vertebrace Extrication Collar Actimove Cervical Collar Actimove Cervical Comfort Cervical Collar Universal Cervical Collar Sized Cervical Collar Universal Cervical Collar Universal Cervical Collar Confort Clavicle Support Comfort Clavicle Brace Hook & Loop Closure Clavicle Support w/Four-Way Strap, Buckle Closures Shoulder Abduction Support	4 4 4 4 4 4 4 4 4 4 4 4
Cast Removal & Accessories  FCT Starter Kit Featuring Delta-Cast* Conformable FCT Accessories Kit American Orthopedic** Cast Cutter Cast Saws American Orthopedic** Cast Dust Vacuum American Orthopedic** Portable Vacuum TurboCare** Cast Removal System Cast Saw Blades BSN Buckets Casting Instruments For the Orthopedic Specialist Zip Stick Cast Removal Aid	22 22 23 23 24 24 25 25 26 27	Philidelphia Collar Vertebrace Extrication Collar Actimove Cervical Collar Actimove Cervical Comfort Cervical Collar Universal Cervical Collar Sized Cervical Collar Universal Cervical Collar Universal Cervical Collar Comfort Clavicle Support Comfort Clavicle Brace Hook & Loop Closure Clavicle Support w/Four-Way Strap, Buckle Closures Shoulder Abduction Support Actimove Shoulder Abduction Sling	4 4 4 4 4 4 4 4 4 4 4 4 4 4 4 4 4 4 4
Cast Removal & Accessories  FCT Starter Kit Featuring Delta-Cast* Conformable  FCT Accessories Kit  American Orthopedic** Cast Cutter Cast Saws  American Orthopedic** Cast Dust Vacuum  American Orthopedic** Portable Vacuum  TurboCare** Cast Removal System  Cast Saw Blades  BSN Buckets  Casting Instruments For the Orthopedic Specialist	22 22 23 23 24 24 25 25 26	Philidelphia Collar Vertebrace Extrication Collar Actimove Cervical Collar Actimove Cervical Comfort Cervical Collar Universal Cervical Collar Sized Cervical Collar Universal Cervical Collar Universal Cervical Collar Comfort Clavicle Support Comfort Clavicle Brace Hook & Loop Closure Clavicle Support w/Four-Way Strap, Buckle Closures Shoulder Abduction Support Actimove Shoulder Abduction Sling Soft Form Posture Control Brace	4 4 4 4 4 4 4 4 4 4 4 4 4 4 4 4 4 4 4
Cast Removal & Accessories  FCT Starter Kit Featuring Delta-Cast* Conformable FCT Accessories Kit American Orthopedic** Cast Cutter Cast Saws American Orthopedic** Cast Dust Vacuum American Orthopedic** Portable Vacuum TurboCare** Cast Removal System Cast Saw Blades BSN Buckets Casting Instruments For the Orthopedic Specialist Zip Stick Cast Removal Aid	22 22 23 23 24 24 25 25 26 27	Philidelphia Collar Vertebrace Extrication Collar Actimove Cervical Collar Actimove Cervical Comfort Cervical Collar Universal Cervical Collar Sized Cervical Collar Universal Cervical Collar Universal Cervical Collar Comfort Clavicle Support Comfort Clavicle Brace Hook & Loop Closure Clavicle Support w/Four-Way Strap, Buckle Closures Shoulder Abduction Support Actimove Shoulder Abduction Sling	4 4 4 4 4 4 4 4 4 4 4 4 4 4 4 4 4 4 4







Orthopaedic Soft Goods continued		Orthopaedic Soft Goods	
Actimove® Eco Shoulder Immobilizer	50	Premium Woven Series, Latex-Free Three & Four-Panel Abdominal Binders	76
Actimove® Sling & Swathe Shoulder Immobilizer	51	Soft Form® Hernia Belt	77
Deluxe Sling & Detachable Swathe Shoulder Immobilizer	51	Soft Form® Hernia Brief	77
Elastic Shoulder Immobilizer	52	Actimove® Abdominal Support	78
Shoulder Immobilizer Sling	52	Universal Abdominal Binder 9" height	78
Actimove® Sling on a Roll	53	Soft Form® Lumbar Sacral Support 11" with Flexible Stays	79
Actimove® Comfort Arm Sling	53	Soft Form® Lumbar Sacral Support 11" with Contoured Rigid Stays	79
Actimove® Eco Arm Sling	54	7" Lumbar Sacral Support	80
Universal Arm Sling	54	Lumbar Sacral Support with Overlapping Abdominal Belt, 10" Height	80
Arm Sling Sized	55	9" Neoprene Lumbar Sacral Support	81
Universal Cradle Arm Sling	55	ProLite® Mesh•Loc™ Lumbar Support	81
Cradle Arm Sling	56	Deluxe Lumbar Support 11" Height	82
ProLite® Universal Arm Sling	56	Sacral Iliac Support with Compression Pad 6" Height	82
Universal Cuff and Collar	57	Criss-Cross Lumbar Sacral Support 10" Height	83
GelBand® Arm Band	57	Dorsal Lumbar 19" Support	83
Epi-Sport® Epicondylitis Clasp	58	Soft Form® Maternity Support Belt	84
Canvas Tennis Elbow Strap	58	Thermal Lumbar Support	84
Neoprene Elbow Strap with Loop	59	V-Loc™ LSO Spinal Orthosis	85
ProLite® Elbow Support Knitted Pullover	59		
ProLite® Commpressive Elbow Support with Viscoelastic Insert	60	Lower Extremities	
Elastic Pullover Elbow Support	60	Actimove® Genu Clip Tri-Panel Knee immobilizer with Clip Closure System	86
Neoprene Elbow Sleeve with Loop Lock	61	Actimove® Tri-Panel Knee Immobilizer	86
ProLite® Low Profile 6" Wrist Splint	61	Bel Knee Immobilizer	87
Soft Fit® Suede Finish Cock Up Wrist Brace	62	Actimove® Eco Knee Immobilizer	87
ProLite® 8" Wrist Splint	62	Universal Wrap-Around Knee Support	88
ProLite® Compression Formed Composite Wrist Splint with Adducted Thumb	63	Actimove Post Trauma Knee Dressing	88
Deluxe Universal Wrist Brace	63	Breathable Universal Tri-Panel Knee Immobilizer	88
Cinch-Lock Wrist & Forearm Brace	64	Breathable Universal Cutaway Knee Immobilizer	89
Uni-Fit® Universal Wrist Splint	65	Tri-Panel Knee Immobilizer	89
CarpalMate® Wrist Support	65	GelBand® Patella Strap	90
Elastic Pullover Wristlet	66	Safe-T-Sport® Sports Neoprene Patella Support	90
Elastic Wrist Wrap	66	ProLite® Compressive Knee Support with Viscoelastic Insert	91
Neoprene Wrist Wrap	66	Safe-T-Sport® Standard Neoprene Knee Sleeve	91
HealWell* Grip Splint Wrist / Hand / Finger Orthosis	67	Safe-T-Sport® Standard Neoprene Knee Support	92
SoftForm® Elegant Wrist Support	67	EZ-ON® Neoprene Knee Wrap	92
Perforated Suede Finish Wrist & Forearm Splint	68	Neoprene Patella Stabilizer with Removable Horseshoe	93
Thumb Spica Soft Fit Wrist Brace	68	Lateral Knee Stabilizer with "J" Shaped Buttress	94
Actimove® Gauntlet Thumb and Wrist Brace	69	SoftForm® Wrap-Around Stabilizing Knee Support	95
Actimove® Eco Wrist Brace	69	OA / Arthritis Knee Brace	96
Canvas Cock-Up Wrist Splint - Lace Up	70	Safe-T-Sport® Neoprene Hinged Knee Stabilizing Brace	97
Elastic Wrist Brace	70	Safe-T-Sport® Neoprene Hinged Knee Stabilizing Brace with Composite Hinges	97
Tether® Thumb Stabilizer	71	Safe-T-Sport® Wrap-Around Hinged Knee Brace	98
ProLite® Neoprene Pull-On Thumb Support	71	FlexLite® Hinged Knee Support	98
. C3™ Deluxe Thumb Splint	72	PowerCentric® Composite, Polycentric Hinged Knee Brace	99
HealWell* Grip Splint WHFO Wrist / Hand / Finger Orthosis	72	Adjustable ROM Knee Brace	100
STAX* Finger Splints	73	Specter® ROM Standard	101
Plastalume Finger Splints	73	Specter* ROM Cool/Lightweight	101
Baseball Finger Splint	74	Specter® ROM Short	101
Curved/Gutter Finger Splint	74	HealWell® AFO Night Wrap	102
Four Prong / Protector Finger Splint	74	HealWell* Cub™ Plantar Fasciitis Night Splint Resting Comfort Slipper™	103
Toad Finger Splint	75	HealWell* Plantar Fasciitis Night Splint	104
Fold-Over / Cot Finger Splints	75	HealWell* Multi AFO/Contracture Splint	105
Finger / Thumb Splint	75	HealWell® C3™ Multi AFO/Contracture Splint	106
• · · · · · · · · · · · · · · · · · · ·	-	HealWell® Soft Ease™ Multi-AFO / Heel Suspender	107
Trunk		Ankle Foot Orthosis/Foot Drop Splint	107
Premium Woven Series, Latex-Free Two-Panel Surgical Rib Belts	76	Pullover Elastic Anklet	108
and the second s	-		

#### Orthopaedic Soft Goods continued Orthopaedic Soft Goods 108 Safe-T-Belt® Working Lumbar Belt 136 Nylon Anklet Safe-T-Sport® Neoprene Ankle Support 108 Safe-T-Belt® Plus 136 EZ-ON® Wrap-Around Ankle Support 109 DynaBack® 136 Ankle Guard with Optional Plastic Inserts 109 Safe-T-Wrist® HD Heavy Duty 137 Canvas Lace-Up Ankle Brace Safe-T-Wrist® SD Standard Duty 137 110 Safe-T-Glove® Vibration Dampening Glove ProLite® Ankle Support Knitted Pullover 110 138 F8 Ankle Brace 111 **Preformed Splints** InnerLOK 8<sup>™</sup> Ankle Brace 111 Forearm Splint Allsport Ankle Brace 112 139 Ankle Stirrup Brace with GelFoam™ 112 Metacarpal Forearm Splint 139 ProLite® Ankle Stirrup Brace with Air Liners 113 Posterior Leg / Knee Splint 139 Actimove® Talocast Functional Ankle Brace 113 Wrist & Forearm Splint 140 Sheartex® Ankle Stirrup Brace 114 Metacarpal Hand Splint 140 Active Ankle® Trainer 114 Colles Splint 140 FlexLite® Sport Articulating Hinged Ankle Brace 115 Accessories StepLite® Easy Strider™ Walker Braces 116 StepLite® Easy Air™ Pneumatic Ankle Walker Brace 116 Poncho Style Vest 141 Patient Walker Belt Actimove® Lower Leg Walker 117 141 Catheter Strap Adjustable Walker 141 117 Specialist® Post-Op Shoe Heel and Elbow Protector 118 142 Deluxe Flexible Sole Healing Shoe 118 Actimove® Heel Protector 142 C3™ Post Op Shoe 119 Flotation Heel and Elbow Protector 143 Super Post Op Shoe Thermal Wraps 119 143 Deluxe Rigid Post-Op with Collar Shoe 120 Soft Point® Viscolas® Heel Spur Cushions 144 Vinyl Molded Sole Rocker Cast Shoe 120 Soft Point® Viscolas® Heel Pads 144 Canvas Molded Sole Rocker Soft Point® Viscolas® Full Insoles 145 121 Open-Toe Cast Shoe 121 Soft Point® Viscolas® Sports Full Insoles 145 Canvas Rocker Bottom Cast Shoe 122 Soft Point® Full Insoles 146 Walking Heels 122 Soft Point® 3/4 Insoles 146 Soft Point® Heel Pads 147 Pediatric / Youth Soft Point® Heel Cushion 147 **Pediatrics Introduction** 123 Soft Point® 3/4 Semi-Rigid Orthotics 148 Neoprene Knee Sleeve 124 Soft Point® 3/4 Flexible Cork Orthotics 149 Neoprene Patella Stabilizer Notes 124 150 Knee Immobilizer 125 Terms of Use 153 Wrist Splint 125 ProLite® Wrist Splint 126 ProLite® Wrist Splint with Abducted Thumb 126 Finger / Thumb Splint 127 Neoprene Elbow Sleeve 127 Neoprene Ankle Support 128 ProLite® Knit Ankle Support 128 Universal Arm Sling 129 Soft Form® Posture Control Brace 129 Cervical Collar 130 Clavicle Support 130 Elastic Shoulder Immobilizer 131 Canyas Rocker Bottom Cast Shoe 131 Pediatric Print Cast Shoe 132 Pediatric Ordering Information 133 Safety Products Safe-T-Lift® LX 134 Safe-T-Lift® LX2 134 Safe-T-Lift® DX 135 Custom Fit® 135







### Ortho-Glass® Splinting System







Ortho-Glass® is the most cost-effective splinting system of choice for every phase of immobilization. Features Interlocking Performance Technology (ILP). ILP eliminates the risk of delamination, ensuring a strong, durable splint.



ORTHO-GLASS <sup>®</sup> SPLINTING SYSTEM - ROLL FORM		
Dimensions	Catalog Number	Quantity
1" x 15'	OG-1L2	2 rolls / case
2" x 15'	OG-2L2	2 rolls / case
3" x 15'	OG-3L2	2 rolls / case
4" x 15'	OG-4L2	2 rolls / case
5" x 15'	OG-5L2	2 rolls / case
6" x 15'	OG-6L2	2 rolls / case
8" x 7.5'	66000687	1 roll

ORTHO-GLASS® PRE-CUT SPLINTS		
Dimensions	Catalog Number	Quantity
2" x 12"	OG-2PC	10 splints / box
3" X 12"	OG-3PC	10 splints / box
3" x 35"	OG-335PC	5 splints / box
4" x 15"	OG-4PC	10 splints / box
4" x 30"	OG-430PC	5 splints / box
5" x 30"	OG-5PC	5 splints / box
5" x 45"	OG-545PC	5 splints / box
6" x 10"	OG-610PC	10 splints / box
6" x 30"	OG-6PC	5 splints / box

ACCESSORIES		
Description	Catalog Number	Quantity
Scissors	OG-SC	1 each
Water Bottle Applicator	67400043	1 each
Splint & Cast Cart	SN4094	1 each
Ortho-Glass® 1 Unit Dispenser	OG-1CD	1 each
Ortho-Glass® 2 Unit Dispenser	OG-2CDA	1 each
Ortho-Glass® 3 Unit Dispenser	OG-3CDA	1 each
Ortho-Glass® 4 Unit Dispenser	OG-4CDA	1 each
Ortho-Glass® 5 Unit Dispenser	OG-5CDA	1 each
Ortho-Glass® 6 Unit Dispenser	OG-6CDA	1 each

FEA	TURES AND BENEFITS
•	Wrinkle-free application
•	Better conformability
•	Eliminates risk of delamination
•	No wrong side application
•	Easy stretch padding covers all edges
•	Roll form eliminates waste
	Full educational support system – workshops, in-services, manuals, wall-charts, videos.



#### EXISTING 7-LAYER TECHNOLOGY



### SafetySplint<sup>™</sup>Unpadded Roll



SafetySplint™ is a unique, strong and durable unpadded splint which allows you to control the amount of cushion between the patient and the splint. The unique liner does not require gloves and allows for easy removal for dressing changes.

SAFETYSPLINT™ UNPADDED ROLL		
Dimensions	Catalog Number	Quantity
3" x 15'	67007969	1 roll / case
4" x 15'	67007970	1 roll / case
5" x 15'	67007971	1 roll / case
6" x 15'	67007968	1 roll / case



FEATURES AND BENEFITS	
Thin profile - excellent feel-through for fracture alignment	
Unique liner allows for no gloves and easy removal for dressing changes	
Roll form - No waste	
Strong and durable	
Radiolucent	
Totally clean application	
Excellent for post-op immobilization	

ACCESSORIES	
Description	Catalog Number
Compartment Dispenser w/Scissors	SUPR3CD

### Ortho-Glass® Unpadded Precut Splinting System



The Ortho-Glass® Unpadded Precut System is a fiberglass splinting system that provides strength and durability in a padded splint. Available in various widths and lengths.

ORTHO-GLASS* UNPADDED PRECUT SPLINTING SYSTEM		
Size	Catalog Number	Quantity
3" x 12"	OG-312NP	Box / 10
3" x 38"	OG-338NP	Box / 10
4" x 15"	OG-415NP	Box / 10
4" x 24"	OG-424NP	Box / 10
4" x 30"	OG-430NP	Box / 10
4" x 38"	OG-438NP	Box / 10

FEATURES AND BENEFITS	
Excellent for post-op application	
Unit dose application	
Strong and durable	
6 sizes available	







#### Ortho-Glass® Strap Kits



The Ortho-Glass® Strap Kit is for use in the rehabilitation stage. This patient-friendly product enables long-term use without the need for a full bandage, enabling easy splint removal. Strap Kits are provided with three latex-free hook and loop, stretchable straps.



ORTHO-GLASS* STRAP KITS		
Description	Catalog Number	Quantity
Upper Extremity Kit	71415-00018-01	10 kits / box
Lower Extremity Kit	71415-00019-01	10 kits / box
Upper Extremity Kit:		Lower Extremity Kit:

3 Red Straps, 1 - 5/8" x 17" and 2 - 2" x 191/2" 3 Blue Straps, 2" x 221/2"

### Ortho-Glass® Ankle Treatment System



The Ortho-Glass® Ankle Treatment System is a complete, custom-moldable single-patient care kit that is a cost-effective treatment from acute to rehab. It is strong, lightweight yet durable and fits easily in a shoe or boot. Available as single or 10 pack.

#### Ortho-Glass® AS Treatment Kit



Kits include ice pack, elastic bandage, U-shape ankle foam cushions, rehab straps and Instruction Leaflet.

ORTHO-GLASS® AS LATEX-FREE 10 PACK		
Size	Catalog Number	Quantity
Small - Medium (< 5'7")	58000013	10pk / box
Large (> 5'7")	58000014	10pk / box

# Size Catalog Number Quantity Pediatric 58000002 1 kit Small - Medium (< 5'7")</td> 58000003 1 kit Large (> 5'7") 58000004 1 kit

# Fast and convenient 100% Latex-free Sizes: Pediatric, small-medium and large

Cost effective treatment from acute to rehab







REHAB



### OCL® Splint Roll - 2 Plaster Splinting System



OCL® is the first name in Plaster Splint Roll. OCL Splint Roll is constructed with 15 or 11 layers of plaster encased in open cell foam and flannel. The open cell foam provides padding and prevents migration while the flannel reduces moisture retention.

- DuoSet<sup>™</sup> plaster allows a quick set time with minimal heat
- · Open cell foam and flannel
- · Reinforced layer
- Available in 2", 3", 4", 5" and 6"
- 11 and 15 layers

ORDERING INFO	OLL	
Dimensions 15 Layer	Catalog Number	Quantity
2" x 20'	SR215	1Roll
3" x 20'	SR315	1Roll
4" x 20'	SR415	1Roll
5" x 20'	SR515	1Roll
6" x 20'	SR615	1Roll

Dimensions 11 Layer	Catalog Number	Quantity
2" x 20'	SR211	1Roll
3" x 20'	SR311	1Roll
4" x 20'	SR411	1Roll
5" x 20'	SR511	1Roll
6" x 20'	SR611	1Roll





### Specialist® J-Splint® Plaster Roll Immobilizer



With a foam and stockinette configuration, the Specialist® J-Splint® Plaster Roll Immobilizer allows for quick, low-cost application. It is packaged in a roll format, allowing custom-cut lengths, thereby eliminating waste.

SPECIALIST* J-S	PLINT <sup>®</sup> PLASTER ROLL II	MMOBILIZER
Dimensions	Catalog Number	Quantity
2" x 20'	7332	1 roll / box
3" x 20'	7333	1 roll / box
4" x 20'	7334	1 roll / box
5" x 20'	7335	1 roll / box
6" x 20'	7336	1 roll / box

#### FEATURES AND BENEFITS

- Utilizes Specialist® Plaster of Paris
- Economical





#### Plaster Roll Immobilizer Dispenser



The Plaster Roll Immobilizer Dispenser is a handy option for the storage of the plaster splint roll products.



Available in 3- and 5-slot designs, this dispenser has a hinged door for easy access and replenishing inventory while a cutting guide allows for quick and easy applications.

PLASTER ROLL IMMOBILIZER DISPENSER	
Description	Catalog Number
5 Station Dispenser (2", 3", 4", 5" & 6" Slots)	GPS-5CD
3 Station Dispenser (3", 4" & 6" Slots)	GPS-3CD

BSN MEDICAL EL	SN MEDICAL ELASTIC BANDAGES		
Catalog Number	Color	Size	Quantity
1037023	Beige	2" x 5 yds.	10 / box*
1037033	Beige	3" x 5 yds.	10 / box*
1037043	Beige	4" x 5 yds.	10 / box*
1037053	Beige	6" x 5 yds.	10 / box*
EB-4FLDB	Beige	4" x 10 yds.	5 rolls / pack
EB-6FLDB	Beige	6" x 10 yds.	5 rolls / pack

<sup>\*</sup> Sold by case of 5 boxes

CUTTING TEMPLATE	ESET		
Catalog Number	Description	Unit	

Plaster Roll Immobilizer Dispenser Accessories: Elastic Bandages



Plaster Roll Immobilizer Dispenser Accessories: Cutting Template



### Gypsona® Endless Slab Dispenser Pack



A unique 5-ply dispenser containing creamy, extra-fast plaster which enables improved detailed impressions and a faster set-up time. The technique used allows for faster application and a shorter casting time. This product is ideal for body jackets.



GYPSONA* ENDLESS SLAB DISPENSER PACK		
Dimensions	Catalog Number	Quantity
4" x 65'	87004913	1 / box
6" x 65'	87004915	1 / box
8" x 65'	87004916	1 / box

#### FEATURES AND BENEFITS

- Convenient 5-ply system
- Extra-fast set time



### Specialist® Plaster of Paris Bandages & Splints







Specialist® plaster of Paris uses a time-tested formula and manufacturing process that produces a consistent plaster that is uniformly coated and adhered to the gauze. This unique formula offers minimal plaster loss to ensure maximum strength and a smooth cast finish.

SPECIALIST* PLA	PECIALIST * PLASTER BANDAGES - Extra fast setting / 2-4 minutes		
Dimensions	Catalog Number	Quantity	
2" x 3yds	7362	6 doz / case	
3" x 3yds	7363	6 doz / case	
4" x 5yds	7367	6 doz / case	
5" x 5yds	7369	4 doz / case	
6" x 5yds	7368	4 doz / case	
8" x 5yds	7370	2 doz / case	

SPECIALIST PLA	PECIALIST PLASTER BANDAGES - Fast setting / 5-8 minutes		
Dimensions	Catalog Number	Quantity	
2" x 3yds	7372	6 doz / case	
3" x 3yds	7373	6 doz / case	
4" x 5yds	7374	6 doz / case	
5" x 5yds	7375	4 doz / case	
6" x 5yds	7376	4 doz / case	
8" x 5yds	7378	4 doz / case	

SPECIALIST® PLASTER SPLINTS - Extra fast setting / 2-4 minutes		
Dimensions	Catalog Number	Quantity
3" x 15"	7390	12 boxes / case
4" x 15"	7391	12 boxes / case
5" x 30"	7392	6 boxes / case

SPECIALIST	g / 5-8 minutes	
Dimensions	Catalog Number	Quantity
3" x 15"	7393	12 boxes / case
4" x 15"	7394	12 boxes / case
5" x 30"	7395	6 boxes / case
5" x 45"	7396	4 boxes / case







FEATURES AND BENEFITS
Consistency
High dry/wet plaster retention
Low exotherm
Moisture-resistant packaging
Excellent strength
Early green strength
Smooth finish

### Orthoflex® Elastic Plaster of Paris Bandages

Orthoflex® elastic plaster bandages are ideal for difficult casting procedures such as post-amputation, as they conform well to body contours and bony prominences. Orthoflex® elastic bandages are manufactured with Specialist® plaster and an elastic gauze.

ORTHOFLEX <sup>®</sup> ELA	STIC PLASTER BANDAGI	ES - Fast setting / 5-8 minutes
Dimensions	Catalog Number	Quantity
4" x 4yds	7324*	6 doz / case
5" x 4vds	7325*	4 doz / case







### Gypsona® S Plaster of Paris Bandages and Splints



Gypsona® S plaster products feature a creamy plaster that results in a smoother, stronger cast. The central plastic core prevents telescoping and improves the handling of the bandage. Available in bandages, splints and elastic plaster bandages.



FEATURES AND BENEFITS
Fast immobilization of fractures
High mass weight resulting in low plaster loss
Strong and durable casts
Central plastic core

or other or endre in or entre				
Dimensions	Catalog Number	Quantity		
3" x 15"	30-7390	50 / box		
4" x 15"	30-7391	50 / box		
FII 00II	20 7222	FO / L		

GYPSONA SPLASTER SPLINTS - Fast setting / 5-8 minutes				
Dimensions	Catalog Number	Quantity		
3" x 15"	30-7393	50 / box		
4" x 15"	30-7394	50 / box		
5" x 30"	30-7395	50 / box		
5" x 45"	30-7396	50 / box		

GYPSONA * S PLASTER BANDAGES - Extra fast setting / 2-4 minutes			
Dimensions	Catalog Number	Quantity	
2" x 3yds	30-7362	12 / box	
3" x 3yds	30-7363	12 / box	
4" x 5yds	30-7367	12 / box	
5" x 5yds	30-7369	12 / box	
6" x 5yds	30-7368	12 / box	
8" x 5yds	30-7370	12 / box	

GYPSONA SPLASTER BANDAGES - Fast setting / 5-8 minutes				
Dimensions	Catalog Number	Quantity		
2" x 3 yds	30-7372	12 / box		
3" x 3 yds	30-7373	12 / box		
4" x 5 yds	30-7374	12 / box		
5" x 5 yds	30-7375	12 / box		
6" x 5 yds	30-7376	12 / box		
8" x 5 yds	30-7378	12 / box		

GYPSONA* S PLAS	TER ELASTIC BANDAGES	- Fast setting / 5-8 minutes
Dimensions	Catalog Number	Quantity
4" x 4 yds	30-7324*	12 / box
5" x 4 yds	30-7325*	12 / box

<sup>\*</sup> Caution: This product contains latex rubber which may cause an allergic reaction

# Gypsona® HP Plaster of Paris Bandages and Splints



The traditional rich and creamy texture of Gypsona® HP allows intimate molding and creates a smooth, appealing finished cast. Unique interlocking leno weave greatly enhances stability, allowing more plaster to stay on the gauze. The result is a premium plaster of Paris.



FEATURE	S AND BENEFITS
• Leno g	gauze for extra stability
Extra-f	ast setting
Cream	y texture

GYPSONA * HP PLASTER BANDAGES - Extra-fast setting / 2-4 minutes				
Dimensions	Catalog Number	Quantity		
2" x 3 yds	30-3041	12 / box		
3" x 3 yds	30-3042	12 / box		
4" x 5 yds	30-3045	12 / box		
6" x 5 yds	30-3046	12 / box		
8" x 5 yds	30-3040	12 / box		

GYPSONA HP PLASTER SPLINTS - Extra-fast setting / 2-4 minutes				
Dimensions	Catalog Number	Quantity		
3" x 15"	30-3047	50 / box		
4" x 15"	30-3048	50 / box		
5" x 30"	30-3049	50 / box		

### OCL® Plaster of Paris Bandages & Splints



OCL plaster products feature a highly creamy texture for ease of molding. The quick wet out and consistent set time are great for making custom splints in the Emergancy Department or negative molds in the O & P office.

#### FEATURES AND BENEFITS

- Creamy texture
- Fast and Extra Fast set time
- Quick Wet out
- Smooth finish

OCL Plaster Bandage - Extra Fast - Set Time 2-4 mins. (Green Label)				
	Dimensions	Catalog Number	Quantity	
	2" x 3 yds Plaster Bandage	41923X	6 doz/cs	
	3" x 3 yds Plaster Bandage	41933X	6 doz/cs	
	3" x 5 yds Plaster Bandage	41935X	6 doz/cs	
	4" x 3 yds Plaster Bandage	41943X	6 doz/cs	
	4" x 5 yds Plaster Bandage	41945X	6 doz/cs	
	5" x 5 yds Plaster Bandage	41955X	4 doz/cs	
	6" x 3 yds Plaster Bandage	41963X	4 doz/cs	
	6" x 5 yds Plaster Bandage	41965X	4 doz/cs	
	8" x 5 yds Plaster Bandage	41985X	3 doz/cs	
	12" x 5 yds Plaster Bandage	418125X	1 doz/cs	
	18" x 5 yds Plaster Bandage	418185X	1 doz/cs	
	32" x 5 yds Plaster Bandage	418325X	1 doz/cs	





ı	Dimensions Catalog Number Quantity			
	Dimensions	Catalog Number	Quantity	
	2" x 3 yds Plaster Bandage	41823F	6 doz/cs	
	3" x 3 yds Plaster Bandage	41833F	6 doz/cs	
	3" x 5 yds Plaster Bandage	41835F	6 doz/cs	
	4" x 5 yds Plaster Bandage	41845F	6 doz/cs	
	5" x 5 yds Plaster Bandage	41855F	4 doz/cs	
	6" x 5 yds Plaster Bandage	41865F	4 doz/cs	
	8" x 5 yds Plaster Bandage	41885F	3 doz/cs	
	12" x 5 yds Plaster Bandage	418125F	1 doz/cs	
	32" x 5 vds Plaster Bandage	418325F	1 doz/cs	

OCL Elastic Plaster Bandage - Set Time 4-6 mins. (Red Label)		
Dimensions	Catalog Number	Quantity
4" x 4 yds Elastic Plaster Bandage	41344E*	6 doz/cs
5" x 4 yds Elastic Plaster Bandage	41354E*	4 doz/cs
6" x 4 yds Elastic Plaster Bandage	41364E*	4 doz/cs

 $<sup>\</sup>ensuremath{^{\star}}$  Caution: This product contains latex rubber which may cause an allergic reaction

OCL Plaster Splints - Extra Fast - Set Time 2-4 mins. (Green Label)		
Dimensions	Catalog Number	Quantity
3" x 15" Plaster Splints, 50/box	43315X	12 boxes/cs
4" x 15" Plaster Splints, 50/box	43415X	12 boxes/cs
5" x 30" Plaster Splints, 50/box	43530X	6 boxes/cs
5" x 45" Plaster Splints, 50/box	43545X	6 boxes/cs

OCL Plaster Splints - Fast - Set Time 5-8 mins. (Blue Label)		
Dimensions	Catalog Number	Quantity
3" x 15" Plaster Splints, 50/box	42315F	12 boxes/cs
4" x 15" Plaster Splints, 50/box	42415F	12 boxes/cs
5" x 30" Plaster Splints, 50/box	42530F	6 boxes/cs
5" x 45" Plaster Splints, 50/box	42545F	6 boxes/cs

# Delta-Lite® Conformable The Durable Solution FREE







Fiberglass Cast Tape

The patented fiberglass substrate of Delta-Lite® Conformable provides unsurpassed conformability, greatly reducing the need to tuck and fold during application. The finished cast is lightweight, strong, resists breakdown and is comfortable to wear.







FEATURES AND BENEFITS
Primary and secondary casting
Rigid support
Superior conformability
Exceptional strength
Excellent lamination
Three-dimensional stretch
Weight bearing in approximately 20 minutes
• 3 - 5 minute set time

DELTA-LITE*	CONFORMABLE	FIBERGLASS CAS	T TAPE
Dimensions	Color	Catalog Number	Quantity
2" x 4 yds	O White	6822A	10 rolls / box
3" x 4 yds	O White	6823A	10 rolls / box
4" x 4 yds	O White	6824A	10 rolls / box
5" x 4 yds	O White	6825A	10 rolls / box
2" x 4 yds	<ul><li>Dk Blue</li></ul>	5942	10 rolls / box
3" x 4 yds	<ul><li>Dk Blue</li></ul>	5943	10 rolls / box
4" x 4 yds	<ul><li>Dk Blue</li></ul>	5944	10 rolls / box
2" x 4 yds	Lt Blue	5902	10 rolls / box
3" x 4 yds	Lt Blue	5903	10 rolls / box
4" x 4 yds	Lt Blue	5904	10 rolls / box
2" x 4 yds	<ul><li>Purple</li></ul>	5972	10 rolls / box
3" x 4 yds	<ul><li>Purple</li></ul>	5973	10 rolls / box
4" x 4 yds	<ul><li>Purple</li></ul>	5974	10 rolls / box
2" x 4 yds	Red	5932	10 rolls / box
3" x 4 yds	Red	5933	10 rolls / box
4" x 4 yds	Red	5934	10 rolls / box
2" x 4 yds	Green	5952	10 rolls / box
3" x 4 yds	Green	5953	10 rolls / box
4" x 4 yds	<ul><li>Green</li></ul>	5954	10 rolls / box
2" x 4 yds	<ul><li>Black</li></ul>	6062	10 rolls / box
3" x 4 yds	<ul><li>Black</li></ul>	6063	10 rolls / box
4" x 4 yds	<ul><li>Black</li></ul>	6064	10 rolls / box
2" x 4 yds	Neon Pink	6052	10 rolls / box
3" x 4 yds	Neon Pink	6053	10 rolls / box
4" x 4 yds	Neon Pink	6054	10 rolls / box
2" x 4 yds	<ul><li>Neon Yellow</li></ul>	6032	10 rolls / box
3" x 4 yds	<ul><li>Neon Yellow</li></ul>	6033	10 rolls / box
4" x 4 yds	<ul><li>Neon Yellow</li></ul>	6034	10 rolls / box
2" x 4 yds	Neon Orange	6022	10 rolls / box
3" x 4 yds	Neon Orange	6023	10 rolls / box
4" x 4 yds	Neon Orange	6024	10 rolls / box
2" x 4 yds	Mixed Pack*	6055	10 rolls / box
3" x 4 yds	Mixed Pack*	6056	10 rolls / box
4" x 4 yds	Mixed Pack*	6057	10 rolls / box

Colors and availability subject to change. Contains 4 pairs of large casting gloves. \*Mixed Pack: Contains 2 dark blue, 2 light blue, 2 red, 2 purple and 2 black







### Delta-Lite® Plus The Trusted Rigid Solution








#### Fiberglass Cast Tape

Delta-Lite® Plus is our new fiberglass cast tape that combines fiberglass substrate and an advanced water activated polyurethane resin system which provides low tack properties. The low tack properties allow easy application and molding using standard examination gloves.

DELTA-LITE*	PLUS FIBERGLA	SS CAST TAPE	
Dimensions	Color	Catalog No.	Quantity
1" x 2 yds	O White	7345800	10 rolls / box
2" x 4 yds	O White	7345801	10 rolls / box
3" x 4 yds	O White	7345802	10 rolls / box
4" x 4 yds	O White	7345803	10 rolls / box
5" x 4 yds	O White	7345804	10 rolls / box
2" x 4 yds	<ul><li>Black</li></ul>	7345845	10 rolls / box
3" x 4 yds	<ul><li>Black</li></ul>	7345846	10 rolls / box
4" x 4 yds	<ul><li>Black</li></ul>	7345847	10 rolls / box
2" x 4 yds	Lt Blue	7345835	10 rolls / box
3" x 4 yds	<ul><li>Lt Blue</li></ul>	7345836	10 rolls / box
4" x 4 yds	Lt Blue	7345837	10 rolls / box
5" x 4 yds	Lt Blue	7345838	10 rolls / box
2" x 4 yds	Purple	7345860	10 rolls / box
3" x 4 yds	Purple	7345861	10 rolls / box
4" x 4 yds	Purple	7345862	10 rolls / box
2" x 4 yds	Pink	7345855	10 rolls / box
3" x 4 yds	Pink	7345856	10 rolls / box
4" x 4 yds	Pink	7345857	10 rolls / box
2" x 4 yds	Deep Red	7345830	10 rolls / box
3" x 4 yds	Deep Red	7345831	10 rolls / box
4" x 4 yds	Deep Red	7345832	10 rolls / box
2" x 4 yds	Deep Blue	7345820	10 rolls / box
3" x 4 yds	Deep Blue	7345821	10 rolls / box
4" x 4 yds	Deep Blue	7345822	10 rolls / box
5" x 4 yds	Deep Blue	7345823	10 rolls / box
2" x 4 yds	Deep Green	7345825	10 rolls / box
3" x 4 yds	Deep Green	7345826	10 rolls / box
4" x 4 yds	Deep Green	7345827	10 rolls / box
5" x 4 yds	Deep Green	7345828	10 rolls / box
2" x 4 yds	Yellow	7345865	10 rolls / box
3" x 4 yds	O Yellow	7345866	10 rolls / box
4" x 4 yds	O Yellow	7345867	10 rolls / box
2" x 4 yds	*Mixed Pack	73458-50	10 rolls / box
3" x 4 yds	*Mixed Pack	73458-51	10 rolls / box
4" x 4 yds	*Mixed Pack	73458-52	10 rolls / box



\*Mixed Pack: Contains 2 dark blue, 2 light blue, 2 red, 2 purple and 2 pink





#### FEATURES AND BENEFITS

- Low tack resin
- Primary and secondary casting
- Rigid support
- Good conformability
- High lamination
- Excellent end lay down
- Smooth surface
- 3 5 minutes set time
- Weight bearing in 20 minutes

#### Delta-Cast® Conformable







The Multi-Function Solution Polyester Cast Tape

Delta-Cast® Conformable is an easy-to-use, fiberglass-free cast tape that produces versatile, lightweight and strong casts. Provides rigid and semi-rigid casting options for primary and secondary casting applications where patient fit, comfort and compliance are critical.



\* See page 21 for FCT Starter Kit.

FEATURES AND BENEFITS		
Primary and secondary casting		
Rigid and semi-rigid support		
Excellent conformability		
Weight bearing in 20 - 30 minutes		
3 - 5 minute set time		
Fiberglass-free		
Clear radiograph viewing		
Strong and lightweight		
Soft edges		
Excellent lamination		
Smooth finish		

- Compilarioc a	re ertiloai.		
DELTA - CAST®	CONFORMABLE F	IBERGLASS-FR	EE CAST TAPE
Dimensions	Color	Catalog No.	Quantity
1" x 2yds	O White	6601	10 rolls / box
2" x 4yds	O White	6602	10 rolls / box
3" x 4yds	O White	6603	10 rolls / box
4" x 4yds	O White	6604	10 rolls / box
5" x 4yds	O White	6605	10 rolls / box
2" x 4yds	<ul><li>Dark Blue</li></ul>	6622	10 rolls / box
3" x 4yds	<ul><li>Dark Blue</li></ul>	6623	10 rolls / box
4" x 4yds	<ul><li>Dark Blue</li></ul>	6624	10 rolls / box
2" x 4yds	<ul><li>Purple</li></ul>	7002	10 rolls / box
3" x 4yds	<ul><li>Purple</li></ul>	7003	10 rolls / box
4" x 4yds	<ul><li>Purple</li></ul>	7004	10 rolls / box
2" x 4yds	<ul><li>Red</li></ul>	6662	10 rolls / box
3" x 4yds	<ul><li>Red</li></ul>	6663	10 rolls / box
4" x 4yds	<ul><li>Red</li></ul>	6664	10 rolls / box
2" x 4yds	<ul><li>Green</li></ul>	7112	10 rolls / box
3" x 4yds	<ul><li>Green</li></ul>	7113	10 rolls / box
4" x 4yds	<ul><li>Green</li></ul>	7114	10 rolls / box
2" x 4yds	<ul><li>Neon Pink</li></ul>	6612	10 rolls / box
3" x 4yds	<ul><li>Neon Pink</li></ul>	6613	10 rolls / box
4" x 4yds	<ul><li>Neon Pink</li></ul>	6614	10 rolls / box
2" x 4yds	<ul><li>Neon Yellow</li></ul>	6632	10 rolls / box
3" x 4yds	<ul><li>Neon Yellow</li></ul>	6633	10 rolls / box
4" x 4yds	<ul><li>Neon Yellow</li></ul>	6634	10 rolls / box
2" x 4yds	<ul><li>Neon Orange</li></ul>	6642	10 rolls / box
3" x 4yds	<ul><li>Neon Orange</li></ul>	6643	10 rolls / box
4" x 4yds	<ul><li>Neon Orange</li></ul>	6644	10 rolls / box
2" x 4yds	<ul><li>Black</li></ul>	7212	10 rolls / box
3" x 4yds	<ul><li>Black</li></ul>	7213	10 rolls / box
4" x 4yds	<ul><li>Black</li></ul>	7214	10 rolls / box
2" x 4yds	*Mixed Pack	6645	10 rolls / box
3" x 4yds	*Mixed Pack	6646	10 rolls / box
4" x 4yds	*Mixed Pack	6647	10 rolls / box

Colors and availability subject to change. Contains 4 pairs of large casting gloves. \*Mixed Pack: Contains 2 dark blue, 2 black, 2 red, 2 purple and 2 pink.



#### FlashCast® Elite The Colorful Solution







#### Printed Polyester Cast Tape

Add personality and color to casts with the world's only printed casting material. In addition to helping relieve patient anxiety and making casts more enjoyable to wear, the finished cast is strong, lightweight and durable. Patterns are colorfast and will not interfere with radiographic viewing.

FLASHCAST® ELITE PRINTED FIBERGLASS-FREE CAST TA			E CAST TAPE
Dimensions	Pattern	Catalog Number	Quantity
3" x 4 yds	Pastel Prints	4073	10 rolls / box
2" x 4 yds	Teddy Bears	4082	10 rolls / box
2" x 4 yds	Dinosaurs	4502	10 rolls / box
2" x 4 yds	Soccer Balls	4112	10 rolls / box
3" x 4 yds	Soccer Balls	4113	10 rolls / box
2" x 4 yds	Stars & Stripes	4142	10 rolls / box
3" x 4 yds	Stars & Stripes	4143	10 rolls / box
2" x 4 yds	Camouflage	4002	10 rolls / box
3" x 4 yds	Camouflage	4003	10 rolls / box

Patterns subject to change. Contains 4 pairs of large casting gloves.

FEATURES AND BENEFITS
Long lasting, vibrant prints
Radiolucent
Smooth cast finish
Superior lamination
Weight bearing in 20 - 30 minutes
Consistent 3 - 5 minute set-time
Available in 6 popular patterns



#### Delta-Cast® Soft The Flexible Solution



#### Semi-Rigid Cast Tape

Delta-Cast® Soft is a semi-rigid fiberglass-free cast tape, ideal for customized support. Delta-Cast® Soft provides variable support for fractures and soft tissue injuries while allowing effective stabilization for the injured area. The use of a semi-rigid cast material allows for greater range of movement, earlier mobility and reduces some of the problems of rigid immobilization, such as muscle atrophy.

DELTA-CAST	DELTA-CAST" SOFT FIBERGLASS-FREE CAST TAPE		
Dimensions	Color	Catalog Number	Quantity
1" x 2 yds	White	7345900	10 rolls/box
2" x 4 yds	O White	7345901	10 rolls/box
3" x 4 yds	O White	7345902	10 rolls/box
4" x 4 yds	O White	7345903	10 rolls/box
5" x 4 yds	O White	7345904	10 rolls/box
2" x 4 yds	<ul><li>Blue</li></ul>	7345910	10 rolls/box
3" x 4 yds	<ul><li>Blue</li></ul>	7345911	10 rolls/box
4" x 4 yds	<ul><li>Blue</li></ul>	7345912	10 rolls/box
2" x 4 yds	<ul><li>Red</li></ul>	7345915	10 rolls/box
3" x 4 yds	<ul><li>Red</li></ul>	7345916	10 rolls/box
4" x 4 yds	Red	7345917	10 rolls/box



Good conformability

· Rapid removal with scissors or by unwrapping

• Excellent for circular and removable casts





### Delta-Dry® Water Resistant Cast Padding





Delta-Dry® features an unique patented technology that out performs other water resistant cast padding. The open knit structure, along with the water resistant qualities of the fibers, allow water to quickly drain from the cast.


#### FEATURES AND BENEFITS

- · Allows patients to shower, swim and bathe like normal
- No special tools are needed for application or removal
- Conformable, enabling patient to having increased comfort and protection
- Easy to apply
- Better long term water resistance
- · Offers more under cast protection
- Eliminates the need for stockinette



#### DELTA-DRY® SYNTHETIC CAST PADDING

Dimensions	Catalog No.	Quantity
2" x 2.6 yds	73443-00000-00	12 rolls / box
3" x 2.6 yds	73443-00001-00	12 rolls / box
4" x 2.6 yds	73443-00002-00	12 rolls / box







#### Moisture Vapor Transmission Rate

- Independent laboratory testing compared the Moisture Vapor Transmission Rate (MVTR)
- Three layers of Delta-Dry<sup>®</sup> outperformed one layer of Brand X by 160% MVTR Delta-Dry<sup>®</sup> vs. Brand X

Laboratory testing results are intended to illustrate a product's performance under controlled conditions. Actual use results may vary.



### Delta-Rol® Synthetic Cast Padding



Delta-Rol® Synthetic Cast Padding is a non-absorbent synthetic material that does not hold moisture against the skin. Its conformable stretch allows narrow widths around small anatomies without cutting and tearing.



DELTA-ROL® SYNTHETIC CAST PADDING		
Dimensions	Catalog No.	Quantity
2" x 4yds	6882	24 rolls - bag / 3 bags - case
3" x 4yds	6883	12 rolls - bag / 3 bags - case
4" x 4yds	6884	12 rolls - bag / 3 bags - case
6" x 4vds	6886	6 rolls - hag / 3 hags - case

#### FEATURES AND BENEFITS Conformable stretch eliminates need for scissors Comfortable Non-absorbent synthetic fiber

#### Specialist® Cotton Blend Cast Padding

Specialist® Cotton Blend Cast Padding has a unique micropleated design for excellent conformability, protection and ease of application without tucking or tearing.

Unique micro-pleats



SPECIALIST® COTTON BLEND CAST PADDING		
Dimensions	Catalog No.	Quantity
2" x 4yds	9062	24 rolls - bag / 6 bags - case
3" x 4yds	9043	12 rolls - bag / 6 bags - case
4" x 4yds	9044	12 rolls - bag / 6 bags - case
6" x 4yds	9046	6 rolls - bag / 6 bags - case

FEATURES AND BENEFITS
Micropleated construction provides excellent conformability
Highly absorbent
Cotton rayon blend
Easy application

### Specialist® 100 Cotton Cast Padding



Specialist® 100 has the product features cotton users desire. The soft, 100% cotton material feathers and self-bonds to create a smooth undercast surface. Specialist® 100 is hand-tearable and applies easily. Excellent cohesion allows for custom padding around bony prominences.



SPECIALIST® COTTON CAST PADDING			
Dimensions	Catalog No.	Quantity	
2" x 4yds	9082	24 rolls - bag / 6 bags - case	
3" x 4yds	9083	12 rolls - bag / 6 bags - case	
4" x 4yds	9084	12 rolls - bag / 6 bags - case	
6" x 4yds	9086	6 rolls - bag / 6 bags - case	

	FE	ATURES AND BENEFITS	
	•	100% cotton	
	•	Excellent cohesion	
Ī	•	Hand tearable	



### Sof-Rol® Absorbent Cast Padding

Sof-Rol® is a highly absorbent material suitable for use around wounds. This soft natural material creates a comfortable undercast surface for either synthetic or natural casting materials.





FEATURES AND BENEFITS
Highly absorbent
Easy application

SOF-ROL® ABSORBENT CAST PADDING			
Dimensions	Catalog No.	Quantity	
2" x 4 yds	9052	24 rolls - bag / 6 bags - case	
3" x 4 yds	9033	12 rolls - bag / 6 bags - case	
4" x 4 yds	9034	12 rolls - bag / 6 bags - case	
6" x 4 yds	9036	6 rolls - bag / 6 bags - case	



### Delta-Net® Orthopedic Synthetic Stockinette

Delta-Net® Orthopedic Synthetic Stockinette stretches to more than twice its width to conform comfortably to complex anatomies. It is suitable for use with synthetic or plaster casting materials and is contained in an easy-to-use dispensing box for easy storage and handling.



F E A	ATURES AND BENEFITS
•	Synthetic fibers
•	Soft conformability

DELTA-NET* ORTHOPEDIC SYNTHETIC STOCKINETTE		
Dimensions	Catalog Number	Quantity
1" x 25 yds	6861	1 roll / case
2" x 25 yds	6862	2 rolls / case
3" x 25 yds	6863	2 rolls / case
4" x 25 yds	6864	2 rolls / case
6" x 25 yds	6866	1 roll / case
8" x 25 yds	6868	1 roll / case
10" x 25 yds	6810	1 roll / case
12" x 25 yds	68612	1 roll / case



### Specialist® Orthopaedic Cotton Stockinette

Specialist® Orthopedic Cotton Stockinette consists of a high-grade absorbent cotton that is an ideal undercast layer. It stretches to three times its width to conform to complex anatomies. It is contained in an easy-to-use dispensing box for easy storage and handling.



FEATURES AND BENEFITS
100% high-grade cotton
Convenient dispensing box

SPECIALIST <sup>®</sup> ORTHOPEDIC COTTON STOCKINETTE			
Dimensions	Catalog Number	Quantity	
1" x 25 yds	9071	1 roll / case	
2" x 25 yds	9072	9 rolls / case	
3" x 25 yds	9073	8 rolls / case	
4" x 25 yds	9074	8 rolls / case	
6" x 25 yds	9076	6 rolls / case	
10" x 25 yds	9010	1 roll / case	
12" x 25 yds	9012	1 roll / case	

### 



Terry-Net" C Cotton Terry Cloth Stockinette is constructed of high-grade absorbent cotton. The terry-net design provides additional padding and stretches to five times its normal width to conform to complex anatomies.

DELTA TERRY-NET	M C COTTON TERRY CLO	TH STOCKINETTE
Dimensions	Catalog Number	Quantity
2" x 15 yds	52012	2 rolls / case
3" x 15 yds	52013	2 rolls / case
4" x 15 yds	52014	2 rolls / case
5" x 15 yds	52015	1 roll / case
6" x 15 yds	52016	1 roll / case



FE	ATURES AND BENEFITS
•	Cotton Terry-Net cushion
•	Stretches five times normal width
•	Convenient dispensing box

### Terry-Net™S Synthetic Terry Cloth Stockinette



Terry-Net<sup>\*\*</sup>S Synthetic Terry Cloth Stockinette is constructed of non-absorbent synthetic

fibers that do not hold moisture against the skin. The Terry-Net design provides additional

padding and stretches to five times its normal width to conform to complex anatomies.

DELTA TERRY-NET™	S SYNTHETIC TERRY CLO	OTH STOCKINETTE
Dimensions	Catalog Number	Quantity
2" x 15 yds	53012	2 rolls / case
3" x 15 yds	53013	2 rolls / case
4" x 15 yds	53014	2 rolls / case
5" x 15 yds	53015	1 roll / case
6" x 15 vds	53016	1 roll / case



FEATURES AND BENEFITS		
Soft synthetic yarns		
Terry cloth design		
Stretches five times normal width		
Dispensing box		

### Terry-Net™ Thumb Spica Precut Synthetic Terry Cloth Stockinette

Terry-Net™ Thumb Spica is comprised of soft, non-absorbent fibers that do not hold moisture against the skin. A built-in thumb liner saves time, enhances fit and is available in a variety of sizes to accommodate most patients.





FΕ	ATURES AND BENEFITS
•	Built-in thumb liner
•	Soft synthetic fibers
•	Terry cloth design

Dimensions	Catalog Number	Overstitus
Difficusions	Catalog Number	Quantity
2" x 16" - Short arm	53120	10 / case
2" x 19" - Long arm	53121	10 / case
3" x 16" - Short arm	53101	10 / case
3" x 19" - Long arm	53102	10 / case

TERRY-NET™ THUMB SPICA PRECUT SYNTHETIC TERRY CLOTH STOCKINET



### Delta Terry-Net<sup>™</sup> Adhesive Fleece Edger

Delta Terry-Net<sup>™</sup> Adhesive Fleece Edger features a high quality adhesive for secure placement. Its unique stretch allows the user to pad easily around corners and contours, ensuring smooth cast edges.



DELTA TERRY-NET™ ADHESIVE FLEECE EDGER			
Dimensions	Catalog Number	Color	Quantity
1 <sup>1</sup> / <sub>4</sub> " x 15 yds	56012	White	1 roll / case
1 <sup>1</sup> / <sub>4</sub> " x 15 yds	58012	Blue	1 roll / case
FEATURES AN	D BENEFITS		
Highly conformable			
Unique stretch			
Secure adhesive			
Works great with thermoplastics			



# Delta Terry-Net<sup>™</sup> Adhesive Terry Cloth/Foam Padding

Delta Terry-Net™ Adhesive Cloth/Foam Padding provides added protection and cushioning by combining a layer of soft terry cloth bonded to ½" foam. The adhesive backing makes attaching the padding simple and easy.



DELIA TERRY-NE	T™ ADHESIVE TERRY CLO	TH/FOAM PADDIN
Dimensions	Catalog Number	Quantity
23" x 39"	55012	1 roll / case
FEATURES AND E	DENEGLES	
FEATURES AND E	BENEFIIS	
<ul> <li>Open cell foam</li> </ul>		
<ul><li>Open cell foam</li><li>Adhesive back</li></ul>		



### Delta Terry-Net<sup>™</sup> Adhesive Felt Padding



Delta Terry-Net™ Adhesive Felt Padding is an adhesive orthopedic felt, used for added patient comfort. Its adhesive backing makes attaching the padding simple and easy. It may also be used as a cast edge padding or to provide protection for bony prominences.

DELTA TERRY-NET™ ADHESIVE FELT PADDING		
Dimensions	Catalog Number	Quantity
23" x 39"	54012	1 roll / case
FEATURES AND BENEFITS		
Soft orthopedic	felt	
Adhesive back		



### Delta Terry-Net Machesive Hook & Stretch Loop Closures



This easy-to-use product simplifies the FCT technique. It provides closure without creating the discomfort caused by traditional loops' rough edges. Its high quality adhesive ensures secure placement and the simple dispensing box provides convenient storage.

DELTA TERRY-NET	「™ ADHESIVE HOOK	
Dimensions	Catalog Number	Quantity
1" x 10 yds	52100	1 roll / box
2" x 10 yds	52200	1 roll / box

DELTA TERRY-NET™ STRETCH LOOP			
Dimensions	Catalog Number	Quantity	
1" x 15 yds	52101	1 roll / box	
2" x 15 yds	52201	1 roll / box	



FEATURES AND BENEFITS	
Soft stretch loop	
High quality adhesive	
Dispensing box	





### FCT Starter Kit Featuring Delta-Cast® Conformable

The FCT Starter Kit enables innovative, functional treatment options by providing ample storage for all your FCT accessories and materials. Customized casts and splints optimize stability using revolutionary techniques to enhance patient comfort, maximize muscle function and reduce rehabilitation time.



FCT STARTER	KIT			
Description	Color	Cat. No.	Size	Qty
FCT Starter Kit	N/A	8000	N/A	1 kit

CONTENTS OF FCT STARTER KIT				
Description	Color	Cat. No.	Size	Qty
Delta-Cast® Conformable	White	6602	2" x 4 yds	1 box / 10 rolls
Delta-Cast® Conformable	White	6603	3" x 4 yds	1 box / 10 rolls
FCT Dispenser	White	51200	19" x 27" x17"	1 dispenser
*FCT Accessories Kit	N/A	72983-00015-00	N/A	1 kit
+0 1 1 6 1				

#### \*See below for contents

#### FEATURES AND BENEFITS

- Easy FCT implementation for your practice
- Functional and convenient



#### **FCT Accessories Kit**

The FCT Accessories Kit contains all the FCT accessories you require in one convenient package. Order it alone or as part of the FCT Starter Kit. Whichever option you choose, the FCT Accessories Kit is your

convenient doorway to FCT.



#### FEATURES AND BENEFITS

- Essential tools
- Convenient ordering

FCT ACCESSORIES	S KIT			
Description	Color	Cat. No.	Size	Qty
FCT Accessories Kit	N/A	72983-00015-00	N/A	1 kit

CONTENTS OF THE FCT ACCESSORIES KIT				
Description	Color	Cat. No.	Size	Qty
BSN Cast Removal Stick	White/Blue	72046-00029-00	19"	3 sticks
Co-Plus® LF Cohesive Bandage	Varied Colors	72100-00020-00	3" x 5 yds	2 rolls
Terry-Net™ 'S'	White	53012	2" x 15 yds	1 case / 2 rolls
Terry-Net™ 'S'	White	53013	3" x 17 yds	1 case / 2 rolls
Terry-Net <sup>™</sup> Thumb Spica	White	53102	3" Long arm	1 box / 10
Delta Terry-Net™ Edger	White	56012	1.25" x 15 yds	1 roll
Delta Terry-Net™ Foam/Terry Liner	White	55012	23" x 39"	1 roll
Delta Terry-Net™ Felt Liner	White	54012	23" x 39"	1 roll
Delta Terry-Net™ Hook	White	52100	1" x 10 yds	1 roll
Delta Terry-Net™ Loop	White	52101	1" x 15 yds	1 roll





#### American Orthopaedic<sup>™</sup> Cast Cutter Cast Saws

The American Orthopaedic™ Cast Cutter provides precisely what is needed when needed:

power, speed, quiet operation, durability, ease of operation and a balance of weight and power.

Includes one 2.5" diameter multi-purpose blade and wrench.

#### AMERICAN ORTHOPAEDIC™ CAST CUTTER SPECIFICATIONS

Cat. No. 0295-200		
Drive	Hex	
Voltage	115 VAC, 1.5 amps hex drive	
Supply Frequency	50-60 Hz	
Cord Length	10' with hospital grade plug	
Certification	UL Certification	

#### Cat. No. 0295-202

Drive	Pin
Voltage	115 VAC, 1.5 amps pin drive
Supply Frequency	50-60 Hz
Cord Length	10' with hospital grade plug
Certification	UL Certification



#### FEATURES AND BENEFITS

- Balance of weight and power helps reduce user fatigue.
- Internally lubricated for longer service.
- · Designed for rugged use.
- Switch is designed and placed so that it will not break if dropped.
- Hex drive provides six (6) blade rotations extending blade life by 50%
  - 10 foot power cord.
  - Hospital grade plug.

# American Orthopaedic<sup>™</sup> Cast Dust Vacuum

The American Orthopaedic™ Cast Dust Vacuum provides an elegant solution to the problem of airborne dust created during cast removal. It's based on in-depth input from orthopedic technologists tired of the problems associated with modified commercial vacuums attached to cast cutters.



#### Cat. No. 0295-400\*

Voltage	115 VAC, 8.5 amps
Supply Frequency	50-60 Hz
Certifications	ETL Certification
Hose Length	10' x 1.25" I.D.
Length	23"
Width	15"
Height	23.5"
Weight	25 lbs.

The American Orthopedic Cast Dust Vacuum is supplied with a pre-installed disposable filter cartridge, 10' x 1.25" I.D. hose, 6 hose/cord straps, vacuum nozzle adapter, owner's manual and registation card.

\* NOTE: This vacuum may also be adapted to the following cast cutters: Stryker #840, Stryker #940 and Zimmer #8909-01. Accessories required.



#### AMERICAN ORTHOPAEDIC™ PORTABLE VACUUM ACCESSORIES Description Cat. No. Unit

Disposable Filter Cartridge for TCCVS1, TCCVS2,	TCFC1	ea
or 0295-400 vaacuums		
Portable filter cartridge for 0295-300 vacuum	0295-175	ea
Straight vacuum nozzle assembly connects 0295-200	0295-283	ea
cutters to 0295-300 or 0295-400		
Vacuum nozzle assembly connects 0295-200 cutters to Stryker* vacuum with cord inside the hose	4183-122	ea





### American Orthopaedic<sup>™</sup> Portable Vacuum

The American Orthopaedic Dust Vacuum is supplied with a pre-installed disposable filter cartridge, six hose straps, storage caps for filter, nozzle and bracket for cast cutter owner's manual and registration card. Hose length 6" x 1.25".

#### AMERICAN ORTHOPAEDIC™ PORTABLE VACUUM ACCESSORIES

Description	Cat. No.	Unit
Disposable Filter Cartridge for TCCVS1, TCCVS2,	TCFC1	ea
or 0295-400 vaacuums		
Portable filter cartridge for 0295-300 vacuum	0295-175	ea
Straight vacuum nozzle assembly connects 0295-200	0295-283	ea
cutters to 0295-300 or 0295-400		
Vacuum nozzle assembly connects 0295-200 cutters to	4183-122	ea
Stryker∗ vacuum with cord inside the hose		

<sup>\*</sup> NOTE: This vacuum may also be adapted to the following cast cutters: Stryker #840, Stryker #940 and Zimmer #8909-01. Accessories required.

#### AMERICAN ORTHOPAEDIC™ PORTABLE VACUUM SPECIFICATIONS

Cat. No. 0295-300*			
Voltage	115 VAC, 6.75 amps		
Supply Frequency	50-60 Hz		
Cord Length	10'		
Airflow	66 CFM		
Waterlift	63"		
Width	7.5"		
Height	9.75"		
Length	19.5"		

#### TurboCare<sup>™</sup> Cast Removal System

The TurboCare® Cast Removal System delivers the exact level of power you need when and where you need it.



#### CAST SAW AND POWER MODULE FEATURES AND BENEFITS

- Innovative motor and saw design delivers high power in a small, light-weight package.
- Automated multi-speed control with "turbo mode" for increased power.
- The ergonomic design is comfortable for the operator.

#### CAST DUST VACUUM FEATURES AND BENEFITS

- Unique lift-out disposable filter cartridge makes waste disposal quicker, easier, safer and cleaner than messy bags.
- Remote on/off switch allows the saw and vacuum to be turned on simultaneously.

TURBOCARE® CAST DUST VACUUM ACCESSORIES		
Description	Cat. No.	Qty
Vacuum nozzle adapter connects the 0295-300 or	4183-123	ea
0295-400 vacuums to Stryker cutter		
Right angle vacuum nozzle assembly connects	4183-124	ea
0295-200 cutters to 0295-300, 0295-400, 0295-250 vacs		
Vacuum nozzle adapter connects Stryker 940 cutter to	0295-405	ea
0295-300 or 0295-400 vacuums		

#### TURBOCARE® CAST SAW with POWER MODULE SPECIFICATIONS

#### Cat. No. TCCSS1

Description	TurobCare® Cast Saw with
	1 orrei module
Power Module	115 VAC, 50-60 Hz, 3.0 amps
	w/o vacuum, 10A when vacuum
	attached
Cord/Hose Length	10' cord and hospital grade plug
	Power Module Cast Saw
Length, Width, Height	7" x 6.5" x 8", 8" x 2.5" x 7"
Weight	14 lbs., 2.5 lbs.
Includes	The TurboCare® Cast Saw is also
	supplied with one (1) 2.5" diameter
	multipurpose blade, wrench, owner's
	manual and registration card.

#### TURBOCARE® CAST VACUUM with POWER MODULE SPECIFICATIONS

#### Cat. No. TCCVS1

Description	TurobCare® Vacuum
Power Module	115 VAC, 50-60 Hz, 8.5 amps
Cord/Hose Length	8.5' x 1" I.D. hose
Length, Width, Height	25" x 15" x 25"
Weight	14 lbs., 25 lbs.
Includes	The TurboCare® Cast Dust Vacuum is supplied with a pre-installed Filter
	Cartridge, 6 hose straps vacuum
	nozzle adapter, owner's manual
	and registration card.

1

6

5

#### Cast Saw Blades

**1. Plaster Blades** are designed with a "Universal Hub" which allows them to fit most popular cast cutters. The blades are manufactured from stainless steel and hardened to increase blade life. Available in 2" and 2.5" diameters. See product numbers 4183-110 and 4183-112.

- **2. Sectioned Blades** are ideal for removing small segments or "cutting windows". The blade is manufactured from stainless steel and hardened to increase blade life. Available only with pin drive hub. See product number 4183-145. BSN also offers a sectioned blade with an attached arbor for use with 0295-270 and 0295-271 autopsy saws. See product number 0295-275.
- 3. Teflon® Coated Blades Designed for use with synthetic casting materials, the Teflon® blades are manufactured from stainless steel, undergo specialized hardening processes for prolonged life, and are coated with Teflon to reduce friction. See product numbers 4183-114 and 31-0265.
- **4. Multi-Purpose Coated Blades** All multi-purpose coated blades feature the "Universal Hub" to allow them to fit most popular cast cutters. They are capable of cutting plaster, synthetic casting materials and other demanding applications. These blades utilize Dicronite®, a unique friction reducing compound, to improve the cutting process and extend the life of the blade. Available in 2" and 2.5" diameters. See product numbers 0295-225 and 0295-226.
- **5. Titanium Nitride "Gold" Blades** BSN's longest lasting cast cutter blade. The titanium nitride coating resists wear making this blade ideal for continuous cutting of casts, plastics and composites. Available in 2.5" diameter only. See product number 31-0168.
- **6. Stainless Steel "Plus" Blade** Special sharpening process provides a sharp, long lasting cutting edge with hardened steel durability. The universal hub will fit hex and pin drive cast cutters. Available in 2" and 2.5" diameters. See product number 0295-219 and 0295-220.

Teflon® is a registered trademark of E.I. du Pont de Nemours and Company Dicronite® is a registered trademark of Lubrication Sciences, Inc.



#### **BSN** Buckets

The sturdy BSN pail provides water supply for cast activation.

CAST APPLICATION ACCESSORIES				
Description	Catalog Number	Qty		
BSN Plastic Pail	72046-00026-00	1 each		
BSN Wax Bucket	4183-155	12/case		





# Casting Instruments For the Orthopaedic Specialist

Casting instruments from BSN medical are safe for the patient and effective for cutting casting material. They are

handcrafted with high quality stainless steel to provide precision performance.



A. Improved Cast Speader

B. Replacement Spring



C. Henning Stainless Plaster Spreader



D. Small Utility Scissors



E. Large Utility Scissors



F. Bohler Serrated Edge Scissors



G. Instrument Case

Instrument Case includes the following: Stille Shear, Cast Spreader, Lister Bandage Scissor (7"), Lister Bandage Scissor (8") and Wire Cutters.





H. FCT 6" Scissors



I. FCT 8" Scissors



J. Ortho-Glass® Scissors



K. Esmarch Scissors



L. Adjustable Cast Stand

# Zip Stick Cast Removal Aid

This convenient cast aid provides additional patient protection during cast removal. The stick slides under the cast and provides a smooth path for scissors or cast saw.

Standard and metric rule adds functionality.

ZIP STICK CAST REMOVAL AID				
Color	Catalog Number	Quantity		
White	72046-00029-00	1 / each		

#### Clean Cut<sup>™</sup> Scissors

These high quality scissors easily cut through rigid casting material and soft stockinette without slipping or bunching. They are manufactured with a high quality stainless steel that is easy to clean and resists corrosion and rust.

CLEAN CUT™ SCISSORS				
Description	Catalog Number	Quantity		
Clean Cut™ Scissor, Large 8″	28230	1		
Clean Cut™ Scissor, Small 6"	28235	1		
FEATURES AND BENEFITS - C	LEAN CUT™ SCISSO	RS 6"		
High quality stainless steel				
Angular small blades				
FEATURES AND BENEFITS - C	LEAN CUT™ SCISSO	RS 8"		
Dual tungsten carbide cutting ed	dges			
Serrated cutting edge				
Double finger ring lower handle				
Protective tip				







#### Co-Plus® LF Cohesive Bandage

Co-Plus® LF is a latex-free, cohesive bandage suitable for light support or compression. Co-Plus® LF is flexible, lightweight and hand-tearable. Because of its unique construction, Co-Plus® LF adheres to itself and not to skin, aiding in patient comfort. No clips or fasteners required. Available in assorted colors and sizes. (Composition: nylon and elstomeric yarns)



CO-PLUS® LF LATEX FREE COHESIVE BANDAGE					
Dimensions	Color	Catalog Number	Quantity		
1" x 5 yds	Beige	72100-00014-00	60 / pack		
2" x 5 yds	Beige	72100-00016-00	36 / pack		
3" x 5 yds	Beige	72100-00017-00	24 / pack		
4" x 5 yds	Beige	72100-00018-00	18 / pack		
6" x 5 yds	Beige	72100-00019-00	12 / pack		
2" x 5 yds	Color Pack*	72100-00033-00	36 / pack		
3" x 5 yds	Color Pack*	72100-00020-00	24 / pack		
4" x 5 yds	Color Pack*	72100-00021-00	18 / pack		

\*Color Pack: Neon pink, blue, light purple, neon green and red.

	EATURES AND BENEFITS
•	Latex-free formula
•	Adheres only to itself
•	Non-slip support
•	Hand tearable



### Ortho-Glass® Wraps Elastic Bandage

Ortho-Glass® Wraps are excellent for mild compression support. Its controlled stretch retains strength and elasticity and is secured with hook and loop closures. These latex-free bandages are ideal for patients who are sensitive to natural rubber latex. Available in assorted colors and sizes.



ORTHO-GLASS® WRAPS	5	
Dimensions	Catalog Number	Quantity
2" x 5 yds. Asst. Pack	CEB-2	18 rolls / pack
3" x 5 yds. Asst. Pack	CEB-3	18 rolls / pack
4" x 5 yds. Asst. Pack	CEB-4	18 rolls / pack

\*Colors available: Navy blue, gold, green and purple




### **BSN Elastic Bandage®**



The BSN Elastic Bandage is latex-free, adjustable and constructed of breathable elastic. Ideal for minor sprains and strains. Each bandage is individually wrapped and may be secured by hook and loop closure or with two metal fasteners. Available in a convenient dispenser box.

BSN MEDICAL ELASTIC BANDAGES						
Dimensions Color Catalog Number Quantity						
2" x 5 yds	Beige	1037023	10 / box*			
3" x 5 yds	<ul><li>Beige</li></ul>	1037033	10 / box*			
4" x 5 yds	<ul><li>Beige</li></ul>	1037043	10 / box*			
6" x 5 yds	<ul><li>Beige</li></ul>	1037053	10 / box*			
4" x 10 yds	<ul><li>Beige</li></ul>	EB-4FLDB	5 rolls / pack			
6" x 10 yds	<ul><li>Beige</li></ul>	EB-6FLDB	5 rolls / pack			





An elastic, multi-purpose bandage that is ideal for light support, compression, fixation of splints and dressing retention. This bandage is washable and reusable and highly absorbent to secretions and perspiration.

IS	ISOBAND* ELASTIC SUPPORT BANDAGE				
D	Dimensions Color Catalog Number Quantity				
6'	" x 5.5 yds	O White	01958	50 rolls / case	
8'	" x 5.5 yds	O White	01959	35 rolls / case	







### Tensogrip® Tubular Support Bandage

Tensogrip® is an elasticated, tubular support bandage, ideal for the effective management of strains, sprains and soft tissue injuries. Provides consistent pressure and support over varying body contours. Constructed from a cotton/rayon blend with

rubber latex yarn. Washable and reusable. Available in white and natural colors.



FEATURES AND BENEFITS
Easy to apply
Exerts even pressure and compression
Washable and reusable
Wide range of sizes

Caution: This product contains latex rubber which may cause an allergic reaction

TENSOGRIP*	TUBULAR	SUPPORT	BANDAGE	
Dimensions	Size	Color	Cat. Number	Quantity
1.75" x 11 yds	Α	O White	7579	1 roll
2.5" x 11 yds	В	O White	7580	1 roll
2.75" x 11 yds	С	O White	7581	1 roll
3" x 11 yds	D	O White	7582	1 roll
3.5" x 11 yds	Е	O White	7583	1 roll
4" x 11 yds	F	O White	7584	1 roll
4.5" x 11 yds	G	O White	7585	1 roll
7" x 11 yds	J	O White	7586	1 roll
8.5" x 11 yds	K	O White	7587	1 roll
12.75" x 11 yds	L	O White	7588	1 roll
14.5" x 11yds	M	O White	7570	1 roll
2.5" x 11yds	В	<ul><li>Beige</li></ul>	7580FL	1 roll
2.75" x 11yds	С	<ul><li>Beige</li></ul>	7581FL	1 roll
3" x 11 yds	D	<ul><li>Beige</li></ul>	7582FL	1 roll
3.5" x 11 yds	Е	<ul><li>Beige</li></ul>	7583FL	1 roll
4" x 11 yds	F	<ul><li>Beige</li></ul>	7584FL	1 roll
4.5" x 11 yds	G	<ul><li>Beige</li></ul>	7585FL	1 roll

### Versagrip® Tubular Support Bandage

Versagrip® is a latex-free, tubular support bandage used for the effective management of minor sprains and strains. Easy to apply, Versagrip®s unique patented elastomer provides controlled support over a wide range of limb diameters. Three sizes cover all your tubular bandage needs.



VERSAGRIP® LATEX-FREE TUBULAR SUPPORT BANDAGE					
Limb circ.	Size	= to sizes	Color	Cat No	Qty
5" to 12.5"	1	A, B, C	<ul><li>Beige</li></ul>	67000354	1 roll
11.75" to 19.5"	2	C, D, E	<ul><li>Beige</li></ul>	67000355	1 roll
19" to 27.5"	3	F, G	<ul><li>Beige</li></ul>	67000356	1 roll

#### FEATURES AND BENEFITS

- Unique patented elastomer 3 sizes of Versagrip® cover the same range of limb diameters as 6 sizes of traditional tubular support
- Latex-free formula minimizes the risk of skin reaction
- Provides controlled pressure
- Can be used to help reduce burn scars and edema, after cast removal and dressing retention



### Tensoshape® Elastic Tubular Bandage

Tensoshape® is an elastic, tubular bandage, suitable for light or medium support. Its anatomical design provides graduated compression, promoting improved blood flow and lymphatic return. Both full leg and below-the-knee bandages available in a variety of sizes.



TENSOSHAPE* TUBU	IDAGE - BELOW THE KNEE		
Limb measurement	Size	Cat. Number	Quantity
12.5" to 14"	Small (B/C)	7589	10 / case
13.7" to 15.2"	Medium (C/D)	7590	10 / case
14.8" to 16.4"	Large (D/E)	7591	10 / case

TENSOSHAPE" TUBU	TENSOSHAPE* TUBULAR SUPPORT BANDAGE - FULL LEG		
Limb measurement	Size	Cat. Number	Quantity
12.5" to 14"	Small (B/C)	7592	10 / case
13.7" to 15.2"	Medium (C/D)	7593	10 / case
14.8" to 16.4"	Large (D/E)	7594	10 / case
16" to 17.6"	X Large (D/E)	7595	10 / case



Caution: This product contains latex rubber which may cause an allergic reaction



### Leukotape® P Rigid Strapping Tape

Leukotape® P is a high strength, rigid strapping tape, with a very strong, zinc oxide adhesive. It is porous and hand-tearable, yet offers extra strong support for sprains and strains. Ideal for patellofemoral taping. Use in conjunction with Cover Roll Stretch.

LEUKOTAPE <sup>®</sup> P RIGID STRAPPING TAPE				
Dimensions	Cat. Number	Box Contents	Case contents	
1.5" x 15 yds	76168	1 roll	30 rolls / case	
FEATURES A	ND BENEFITS			
• Porous				
Hand-tearable				
High strength adhesive				
Ideal for patellofemoral taping				
Caution: This product contains latex rubber which may cause an alleroic reaction				







### Lightplast® Pro Elastic Adhesive Bandage

Lightplast® Pro is a tearable, elastic adhesive bandage, suitable for light support or compression. The unique fabric backing helps to conform to body movement. Ideal for all-purpose taping.





LIGHTPLAST	PRO ELASTI	C ADHESIVE BANDA	G E
Dimensions	Color	Catalog Number	Quantity
1" x 5 yds	O White	76952	48 rolls / case
1.5" x 5 yds	White	76953	32 rolls / case
2" x 5 yds	White	76954	24 rolls / case
3" x 5 yds	White	76955	16 rolls / case
2" x 7.5 yds	White	76956	24 rolls / case
3" x 7.5 yds	White	76957	16 rolls / case
2" x 5 yds	<ul><li>Black</li></ul>	72880	24 rolls / case
3" x 5 yds	<ul><li>Black</li></ul>	72881	16 rolls / case
3" x 7.5 yds	<ul><li>Black</li></ul>	72883	16 rolls / case

- Hand-tearable
- Excellent adhesion
- Air permeable
- · Light and close fitting

### Optiplaste® Elastic Adhesive Bandage

An elastic adhesive bandage, suitable for medium support or compression. Ideal for securing leg/limb support and wound dressings with no slippage. Patented weft-inserted construction provides strength and stability with improved softness and comfort - unique fabric backing helps to conform to body movement.



OPTIPLASTE*	ELASTIC AD	HESIVE BANDAGE	
Dimensions	Color	Catalog Number	Quantity
1" x 5 yds	Beige	72056-00060-00	12 rolls / box
2" x 5 yds	<ul><li>Beige</li></ul>	72056-00061-00	6 rolls / box
3" x 5 yds	<ul><li>Beige</li></ul>	72056-00062-00	4 rolls / box
4" x 5 yds	<ul><li>Beige</li></ul>	72056-00063-00	6 rolls / box

#### FEATURES AND BENEFITS

- Patented weft-inserted construction
- High quality extensible cloth
- Porous adhesive
- Fluffy non-fray edges

Caution: This product contains latex rubber which may cause an allergic reaction



# Tensoplast® Elastic Adhesive Bandage





Tensoplast®, formerly Elastoplast®, is an elastic adhesive bandage which is suitable for medium support or compression. Superior stretch and regain provides firm, controlled compression and support. The porous adhesive helps prevent skin maceration. Fluffy non-fray edges prevent the bandage from rucking or marking.

TENSOPLAST® ELASTIC ADHESIVE BANDAGE - INDIVIDUALLY BOXED			
Dimensions	Color	Catalog Number	Quantity
1" x 5 yds	<ul><li>Beige</li></ul>	02598002	36 rolls / case
2" x 5 yds	<ul><li>Beige</li></ul>	02599002	36 rolls / case
3" x 5 yds	<ul><li>Beige</li></ul>	02600002	36 rolls / case
4" x 5 yds	<ul><li>Beige</li></ul>	02601002	36 rolls / case
6" x 5 yds	<ul><li>Beige</li></ul>	02602002	12 rolls / case
1" x 5 yds	O White	02593002	36 rolls / case
2" x 5 yds	O White	02594002	36 rolls / case
3" x 5 yds	O White	02595002	36 rolls / case
4" x 5 yds	O White	02596002	36 rolls / case
6" x 5 yds	O White	02597002	12 rolls / case

TENSOPLAST	* ELASTIC A	DHESIVE BANDAGE -	BOXED PACKS
Dimensions	Color	Catalog No.	Quantity
2" x 5 yds	<ul><li>Beige</li></ul>	04415001	72 rolls / case
3" x 5 vds	<ul><li>Beige</li></ul>	04416001	48 rolls /case

TENSOPLAST	* ELASTIC ADF	IESIVE BANDAGE -	BULK PACKS
Dimensions	Color	Catalog No.	Quantity
1" x 5 yds	<ul><li>Beige</li></ul>	04410001	32 rolls / case
2" x 5 yds	<ul><li>Beige</li></ul>	04412001	24 rolls / case
3" x 5 yds	<ul><li>Beige</li></ul>	04413001	16 rolls / case
4" x 5 yds	<ul><li>Beige</li></ul>	04414001	16 rolls / case





FEATURES AND BENEFITS		
High adhesive properties		
Provides strong compression and long term support		
Water repellent		
Air permeable		
Caution: This product contains latex rubber which may cause an allergic reaction		





# Specialist® Fracture Bracing

These bracing products are fabricated from durable, lightweight thermoplastics. Simple trimming of these products assures a custom fit while designed to provide rigid, cast-like support with a comfortable anatomical shape.













A. SPECIALIS	ST* WRIST-HAND ORTHOSIS	
Size	Catalog Number	Catalog Number
Small	61531 - Right	61532 - Left
Medium	61541 - Right	61542 - Left
Large	61551 - Right	61552 - Left

B. SPECIALIS	B. SPECIALIST WRIST-HAND-THUMB ORTHOSIS			
Size	Catalog Number	Catalog Number		
Small	61831 - Right	61832 - Left		
Medium	61841 - Right	61842 - Left		
Large	61851 - Right	61852 - Left		

C. SPECIALIS	ST° THUMB ORTHOSIS	
Size	Catalog Number	Catalog Number
Small	61931 - Right	61932 - Left
Medium	61941 - Right	61942 - Left
Large	61951 - Right	61952 - Left

D. SPECIALIST* ULNAR FRACTURE ORTHOSIS					
Distal Circ.	Proximal Circ.	Size	Catalog Number		
6" - 7"	8 - 11"	Small	62033		
6.75" - 7.75"	9.5" - 12.5"	Medium	62043		
7.5" - 8.5"	11" - 14"	Large	62053		

E. SPECIALIST <sup>®</sup> HUMERAL FRACTURE ORTHOSIS					
Length	Bicep Circ.	Size	Catalog Number		
4.75"	7.5" - 11"	Small	62533		
5.25"	10" - 13"	Medium	62543		
5.75"	12" - 15"	Large	62553		
6.25"	14" - 17"	X Large	62563		

F. SPECIALIST* HUMERAL FRACTURE ORTHOSIS W/DELTOID EX.					
Length	Bicep Circ.	Size	Catalog Number		
7.5"	7.5" - 11"	Small	62333		
81/10"	10" - 13"	Medium	62343		
84/5"	12" - 15"	Large	62353		
9 2/5"	14" - 17"	X Large	62363		



### Orthoplast® Splint



Orthoplast® is ideal for a variety of rehabilitation and corrective prosthesis, consisting of a material that can be easily cut and molded. It is self-adhering when heated, with edges that trim and finish nicely. Ideal for large orthosis and has a working time of 5 - 8 minutes.

ORTHOPLAST* PLAIN SHEETS					
Dimensions	Color	Catalog Number	Quantity		
18" x 24"	White	9021	4 / case		
ORTHOPLAST*	PERFORATED S	HEETS			
Dimensions	Color	Catalog Number	Quantity		
24" x 36"	White	9022	2 / case		
FEATURES A	ND BENEFITS				
Excellent strength and mobility					
Easily cut and molded to body contours					
Self-adhering when heated					





## Orthoplast® II splint

FORMULA

Orthoplast® II is specially formulated to enhance upper extremity applications. It consists of a material that can be easily cut and molded and is self-adhering when heated, with edges that trim and finish nicely. Working time of 3 - 4 minutes.

ORTHOPLAST	* II PLAIN SHEET	S	
Dimensions	Color	Catalog Number	Quantity
18" x 24"	O White	9121	4 / case
ORTHOPLAST	* II PERFORATED	SHEETS	
Dimensions	Color	Catalog Number	Quantity
Dimensions 24" x 36"	O White	Catalog Number 9122	Quantity 2 / case
2	-		-
24" x 36"	-		-
24" x 36"  FEATURES A	O White	9122	-
24" x 36"  FEATURES A  Excellent of	O White	9122 remities	-







#### Coverlet® Fabric Adhesive Bandages

Coverlet® Fabric Adhesive Bandages are made to absorb liquid fast - faster than any of the competitors tested. Featuring an extra-large wound pad capable of absorbing 10 times its own weight in water. The 360° adhesive surrounding the island pad helps seal the afflicted area, keeping the wounds clean and providing an environment conducive to healing.



COVERL	COVERLET® LATEX-FREE DIGITS				
Description		Catalog No	Qty	Case Qty	
Sm. fing	er tip	퉈	01306	100 / box	12 boxes
Lg. finge	er tip		01307	50 / box	12 boxes

COVERLET LATEX-FREE FOUR WING					
Dimensions	Catalog No	Qty	Case Qty		
3" X 3"	00385	50 / box	12 boxes		

# Description Catalog No Coverlet\* Filled\* - 18 Drawers 70550168 Coverlet\* Empty\* - 18 Drawers 70550167

FEATURES AND BENEFITS
Available in 13 different shapes and sizes
Perfect for use on knees, elbows, chins, noses and knuckles
Elastic fabric moves with the curves and contours of skin and body
100% latex-free

COVERLET® LATEX-FREE STRIPS				
Dimensions		Catalog No	Qty	Case Qty
0.75" x 3"	10103	00230	100 / box	12 boxes
1" x 3"		00231	100 / box	12 boxes
1" x 3"	脚一脚	76448000	300 / box	6 boxes

COVERLET <sup>®</sup> LATEX-FREE SPOTS				
Dimensions		Catalog No	Qty	Case Qty
0.875" round		00301	100 / box	12 boxes
1.25" oval		00303	100 / box	12 boxes

COVERLET® LATEX-FREE PATCHES				
Dimensions		Catalog No	Qty	Case Qty
1.5" x 3"		00330	100 / box	12 boxes
2" x 3"	100	00340000	50 / box	12 boxes

COVERLET® L	COVERLET® LATEX-FREE KNUCKLE					
Dimensions		Catalog No	Qty	Case Qty		
1.5" x 3"	516	01390	100 / box	12 boxes		

COVERLET® TOE SHIELD (CONTAINS LATEX)				
Dimensions Catalog No Qty Case Qty				
2.5" x 1.375"	00801	100 / box	12 boxes	

#### Coverlet® Eye Occlusors



The Coverlet<sup>®</sup> is a comfortable solution for the treatment of lazy eye. This bandage is composed of an inlay wound pad that forms a complete cover over the eye.



COVERLET <sup>®</sup> EYE OCCLUSORS				
Dimensions	Color	Catalog No	Qty	Case Qty
1.8" x 2.5" Jr	Beige	46429	20 / box	12 boxes
2" x 3" Reg	Beige	46430	20 / box	12 boxes

#### FEATURES AND BENEFITS

- Non-irritating adhesive for easy removal
- Complete coverage of the eye
- Elastic & comfortable

Caution: This product contains latex rubber which may cause an allergic reaction



### Elastomull® Elastic Gauze Bandage

Elastomull® offers a precision fit, even in areas where there is constant movement such as the elbow or knee. Superior performance and stretchability in controlling bleeding or edema. It consists of rayon and polyester that is virtually lint-free, reducing risk of particle contamination.

ELASTOMUL® ELASTIC GAUZE BANDAGE - STERILE							
Dimensions Catalog No Quantity Case Quantity							
1" x 4.1 yds	02075001	12 rolls	8 bags				
2" x 4.1 yds	02076001	12 rolls	8 bags				
3" x 4.1 yds	02070001	12 rolls	8 bags				
4" x 4.1 yds	02071001	12 rolls	8 bags				
6" x 4.1 yds	02072001	6 rolls	8 bags				

ELASTOMUL® ELASTIC GAUZE BANDAGE - NON-STERILE			
Dimensions	Catalog No	Quantity	Case Quantity
1" x 4.1 yds	02088000	24 rolls	8 bags
2" x 4.1 yds	02089000	12 rolls	8 bags
3" x 4.1 yds	02101000	12 rolls	8 bags
4" x 4.1 yds	02102000	12 rolls	8 bags
6" x 4.1 yds	02103000	8 rolls	8 bags



#### FEATURES AND BENEFITS

- · High density fabric for a soft, rich texture to ensure patient comfort
- Superior performance and stretchability
- Precision fit for areas with constant movement (elbows and knees)
- Helps prevent slippage once bandage is in place
- Available in sterile and non-sterile
- Non-sterile bandage is suitable for securing primary dressing and sterile bandage may be placed directly on wound

#### Cover-Roll® stretch Non-woven Adhesive Bandage

Cover-Roll® stretch is a wide-area fixation dressing, widely used in place of surgical tapes to secure and retain primary dressings on surgical or wound sites. This latex-free bandage provides light compression without causing constriction for joints and highly contoured parts of the body. Convenient and easily cut to size, eliminating unnecessary waste.

#### FEATURES AND BENEFITS

- · Soft, stretchable non-woven polyester adapts well to the body
- Provides light compression without causing constriction
- Polyacrylate adhesive holds dressing safely, even in the presence of moisture
- 100% latex-free
- · Easily cut to eliminate waste

#### INDICATIONS

- · Retains primary dressing on small to large surgical or wound sites
- Ideal use for areas of body requiring flexibility such as elbows and knees
- Ideal base coverage for areas applying the patella/femoral taping technique

COVER-ROLL* NON-WOVEN ADHESIVE BANDAGE			
Dimensions	Catalog No	Quantity	Case Quantity
2" x 10 yds	2034	1 roll	12 boxes
4" x 10 yds	2035	1 roll	12 boxes
6" x 10 yds	2041	1 roll	12 boxes
8" x 10 yds	2042	1 roll	12 boxes
12" x 10 yds	2043	1 roll	12 boxes





Dimensions	Catalog No	Quantity	Case Quantity
2" x 2 yds	45547	1 roll	12 boxes
4" x 2 yds	45548	1 roll	12 boxes
6" x 2 yds	45549	1 roll	12 boxes
8" x 2 yds	45550	1 roll	12 boxes
12" x 2 yds	45551	1 roll	12 boxes
2" x 10 yds	45552	1 roll	12 boxes
4" x 10 yds	45553	1 roll	12 boxes
6" x 10 yds	45554	1 roll	12 boxes
8" x 10 yds	45555	1 roll	12 boxes
12" x 10 yds	45556	1 roll	12 boxes

COVER-ROLL\* STRETCH NON-WOVEN ADHESIVE BANDAGE

#### Cuticell® Ointment Dressing

Cuticell

2 ms 23 ms / 20 m, 23 ms / 20 ms /

Cuticell® is a sterile ointment dressing designed to protect and minimize pain at dressing change. The primary wound contact layer is made of smooth acetate fabric, impregnated with a non-medicated ointment. The result is an ointment dressing that ensures the wound is not disturbed and maintains a moist environment that supports wound healing.



CUTICELL* OINTMENT DRESSING			
Dimensions	Catalog Number	Quantity	
3" x 3"	7253900	10 / box	
3" x 3"	7253901	50 / box	
3" x 8"	7253902	10 / box	
3" x 8"	7253903	50 / box	

#### FEATURES AND BENEFITS

- The defined, small mesh size prevents tissue growth and does not disturb wound healing
- The smooth acetate fabric with knit structures minimizes pain at dressing change
- The dressing is permeable to air and allows exudate to drain freely into an absorbent secondary dressing
- Ointment impregnation helps maintain a moist healing environment

#### INDICATIONS

- First and second degree burns
- Cuts and abrasions
- Lacerations
- Radiation injuries
- Donor and recipient skin graft sites



### Leukomed® Absorbent Wound Dressing

Leukomed® Adhesive Island Dressing is composed of an acrylic adhesive and non-woven substrate with a highly absorbent wound pad for secreting wounds. The dressing provides a latex-free adhesive which enables easy removal,

while the non-adherent pad helps prevent adhesion to the wound.



- Post surgical incisions
- Moderate drainage wounds
- Pressure ulcers

INDICATION

- Partial and full thickness wounds
- As a secondary dressing over gels, ointments



LEUKOMED* ADHESIVE DRESSING			
Dimensions	Catalog No	Quantity	Case Quantity
2.8" x 2"	7238014	50 / box	10 boxes
3.1" x 4"	7238001	50 / box	8 boxes
3.1" x 6"	7238002	50 / box	8 boxes
4" x 7.8"	7238003	50 / box	8 boxes
4" x 9.8"	7238004	50 / box	8 boxes
4" x 11.8"	7238005	50 / box	8 boxes
4" x 13.8"	7238006	50 / box	8 boxes

#### FEATURES AND BENEFITS

- Good absorption capacity
- Non-adherent layer helps prevents wound pad from sticking to wound
- Latex-free acrylic adhesive provides safe and secure fit
- Skin friendly for easy removal
- Rounded edges
- Highly conformable with two-way stretch
- · Permeable to air
- · Radio transparent

#### Leukomed® T Transparent Wound Dressing

Leukomed® T is a sterile, breathable, waterproof, transparent film dressing designed to cover catheter or wound sites while providing a barrier to external contaminants such as bacteria, blood and body fluids.

LEUKOMED° T	TRANSPARENT FI	LM DRESSING	
Dimensions	Catalog No	Quantity	Case Quantity
2.8" x 2"	7238106	50 / box	10 boxes
3.1" x 4"	7238101	50 / box	10 boxes
4.3" x 5.5"	7238102	50 / box	10 boxes

#### FEATURES AND BENEFITS

- Permeable to water vapor and oxygen
- Sterilized film provides easy wound inspection without removing the dressing
- · Waterproof and bacteria proof
- Rounded edges prevents dressing from curling off
- Patented easy-to-use delivery system



• For use over catheter sites and devices

#### Leukomed®T plus

#### Absorbent Film Wound Dressing

Leukomed® T plus is a sterile, breathable, waterproof, transparent film dressing with a highly absorbent wound pad for exuding wounds.

LEUKOMED T PLUS ABSORBENT FILM WOUND DRESSING			
Dimensions	Catalog No	Quantity	Case Quantity
2.8" x 2"	7238212	50 / box	10 boxes
3.1" x 4"	7238201	50 / box	8 boxes
3.1" x 6"	7238202	50 / box	8 boxes
4" x 9.8"	7238203	50 / box	6 boxes
4" x 11.8"	7238204	50 / box	5 boxes
4" x 13.8"	7238205	50 / box	4 boxes

#### FEATURES AND BENEFITS

- Convenient, all-in-one sterile dressing
- · Water and bacteria proof allowing patients to shower
- Good absorption capacity
- · Coated pad prevents adhesion to wound
- Provides moist wound environment to enhance healing
- Easy to apply
- Conformable





- Post-surgical dressings to protect closed wounds
- Moderate drainage wounds
- Acute wounds i.e. cuts, burns and abrasions
- Chronic wounds



### $\textbf{Anti-Em/GP}^{\text{TM}} \text{ Seamless Elastic Stockings with Gradient Compression}$



In cooperation with leading vascular specialists, JOBST engineers developed the JOBST® Gradient Compression Seamless Anti-Embolism Elastic Stocking specifically designed to combat thromboembolism.

The JOBST Anti-Em/GP Stocking applies a gradient compression, greatest at the ankle, and gradually decreasing up the length of the leg without tourniqueting or impeding the venous return. The use of JOBST Anti-Em/GP Seamless Stockings is an accepted prophylaxis for minimizing the risk of thromboembolism in hospitalized patients requiring bed rest or prolonged periods of immobilization.

- JOBST Anti-Em/GP stockings are latex free
- A convenient inspection window is located on top of the foot so the patients skin will never touch the floor
- Toe color coordinates with sizing chart
- Seamless elastic stockings with gradient compression



ANTI EM / GP	COMPRESSION S	TOCKINGS SPE	CIFICATIONS
Description	Catalog No	Size	Quantity
Knee Length	111402	Sm Reg	12 pr / box
Knee Length	111403	Sm Long	12 pr / box
Knee Length	111406	Med Reg	12 pr / box
Knee Length	111407	Med Long	12 pr / box
Knee Length	111410	Lg Reg	12 pr / box
Knee Length	111411	Lg Long	12 pr / box
Knee Length	111414	X-Lg Reg	12 pr / box
Knee Length	111415	X-Lg Long	12 pr / box
Thigh Length	111450	Sm Short	6 pr / box
Thigh Length	111451	Sm Reg	6 pr / box
Thigh Length	111452	Sm Long	6 pr / box
Thigh Length	111454	Med Short	6 pr / box
Thigh Length	111455	Med Reg	6 pr / box
Thigh Length	111456	Med Long	6 pr / box
Thigh Length	111458	Lg Short	6 pr / box
Thigh Length	111459	Lg Reg	6 pr / box
Thigh Length	111460	Lg Long	6 pr / box
Thigh Length	111462	X-Lg Reg	6 pr / box
Thigh Length	111463	X-Lg Long	6 pr / box
Thigh Length	111464	XX-Lg Reg	6 pr / box
Thigh Length	111465	XX-Lg Long	6 pr / box
Waist Length	111625	Sm Reg	6 pr / box
Waist Length	111626	Sm Long	6 pr / box
Waist Length	111627	Med Reg	6 pr / box
Waist Length	111628	Med Long	6 pr / box
Waist Length	111629	Lg Reg	6 pr / box
Waist Length	111630	Lg Long	6 pr / box



### Gelocast® Original Unna's Boot

Gelocast® the original Unna's boot is impregnated with zinc oxide and calamine to soothe irritated skin. The non-raveling gauze is uniformly coated so that it is easy to handle and mold to the limb providing a semi-rigid cast that is ideal for immobilization of strains and sprains.

GELOCAST* ORIGINAL UNNA'S BOOT			
Dimensions	Catalog No	Quantity	Case Quantity
3" x 10 yds	01052	1 per box	12 boxes
4" x 10 yds	01053	1 per box	12 boxes





INI	DICATION
•	Venous stasis ulcers
•	Phlebitis
•	Strains and sprains
•	Dislocations
•	As an after-cast

FEATURES AND BENEFITS
Dries quickly providing a semi-rigid cast
Soothing calamine reduces skin irritation
100% cotton base does not unravel reducing wastage
Provides high working pressure and lower resting pressure
Uniformly impregnated, Gelocast is easy to apply and remove
<ul> <li>Each bandage is sealed in a ethylene vinyl acetate (EVA) bag to prolong shelf life</li> </ul>
Available in two sizes: 3" and 4"

NOTE: Please refer to the product label and/or package insert for full instructions on the safe use of these products.



#### Miami J® Cervical Collar

The Miami J<sup>®</sup> Cervical Collar is especially designed to meet the needs of ambulatory patients over and extended time. Technologically advanced Sorbatex padding adds long term comfort by quickly absorbing moisture away from the skin.



Description	Size	Cat. No.	Replacement Pads
Collar w/pads	X Small	MJ 250	MJP 250
Collar w/pads	Small	MJ 300	MJP 100
Collar w/pads	Medium	MJ 400	MJP 100
Collar w/pads	Large	MJ 500	MJP 100
Collar w/pads	X Short	MJ 200S	MJP 100
Collar w/pads	Stout	MJ 200L	MJP 200L

MIAMI J\* CERVICAL COLLAR WITH PADS

FEATURES AND BENEFITS
<ul> <li>Polyurethane shell provides rigid support and immobilization</li> </ul>
Two piece design allows symmetrical application
Angulation of chin and occipital portion of collar is adjustable
Sorbatex padding provides comfort through moisture control

MIAMI J <sup>®</sup> CERVICAL	COLLAR	WITH EXTRA	PADS
Description	Size	Cat. No.	Replacement Pads
Collar w/extra pads	X Small	MJR 250	MJP 250
Collar w/extra pads	Small	MJR 300	MJP 100
Collar w/extra pads	Medium	MJR 400	MJP 100
Collar w/extra pads	Large	MJR 500	MJP 100
Collar w/extra pads	X Short	MJR 200S	MJP 100
Collar w/extra pads	Stout	MJ 200L	MJP 200L

#### Miami JR.® Cervical Collar

The Miami JR° Cervical Collar is shaped to fit the developing anatomy of infants, toddlers and adolescents. This fit provides superior immobilization and comfort.



MIAMI JR* CERVICAL COLLAR WITH PADS			
Description	Size	Cat. No.	
Collar w/pads	0 - 6 months	MJ-P0	
Collar w/pads	6 months - 2 years	MJ-P1	
Collar w/pads	2 - 6 years	MJ-P2	
Collar w/pads	6 - 12 years	MJ-P3	

# DescriptionSizeCat. No.Collar w/extra pads0 - 6 monthsMJR-P0Collar w/extra pads6 months - 2 yearsMJR-P1Collar w/extra pads2 - 6 yearsMJR-P2Collar w/extra pads6 - 12 yearsMJR-P3



#### Miami JTO® Thoracic Extension

The Miami JTO® Thoracic Extension extends the control of any adult Miami J® collar. The modular chin strut enhances flexion control by bridging the chin and sternum while the floating back improves fit and allows patients to sleep comfortably without bars.

MIAMI JTO <sup>®</sup> THORACIC EXTENSION	
Description	Catalog Number
Thoracic Extension w/replacement pads	MTR-900



#### Philadelphia<sup>®</sup> Collar

Lightweight molded foam with rigid plastic occipital and mandible posts for greater stabilization. Two piece construction for easy application with hook and loop closure. All models are fully radiolucent. Tracheostomy opening in front.

PHILADELPHIA*	COLLAR			
2.25" Tall	3.25" Tall	4.25" Tall	Size	Inches
10-622435	10-632432	10-642439	Small	10 - 12"
10-622532	10-632539	10-642536	Medium	13 - 15"
10-622639	10-632636	10-642633	Large	16 - 19"



#### Vertebrace® Extrication Collar

The Vertebrace Extrication Collar is a one-piece construction with an integral flip-up chin support and a large trach opening. It contains a reinforced back which provides firm, cervical support and is X-ray translucent. Stack flat or roll tight for easy storage.

VERTEBRACE® EXTRICATION COLLAR			
Dimensions	Size	Catalog Number	
2.5" - 17"	Pedi-Short	79300-00000-01	
3.5" - 17.75"	Pediatric	79305-00000-01	
2.5" - 22.5"	X-Short*	79301-00000-01	
3.5" - 22.5"	Short	79302-00000-01	
4.5" - 2.5"	Regular	79303-00000-01	
5.5" - 22.5"	Tall	79304-00000-01	

\* For patients with little or no neck

#### FEATURES AND BENEFITS

- Quick and easy to apply
- Enables transportation of the patient
- Does not obstruct other emergency treatment
- Easy to store and carry
- Apply without moving the patient



#### Actimove® Cervical Collar



The Actimove® Cervical Collar is used for acute neck injuries or chronic cervical syndromes, maintaining neck stability and range-of-motion restriction. It is available in medium and firm densities and is comprised of a pre-shaped foam for a good

ACTIMOVE° CE	RVICAL CO	LLAR - MEDIUM DE	NSITY
Neck Circ.	Neck	Size	Catalog Number
10.4" - 13.2"	2.5"	X Small	72819-00025-00
12" - 14.4"	3"	Small	72819-00026-00
13.2" - 15.6"	3.5"	Medium	72819-00027-00
14.4" - 16.8"	4"	Large	72819-00028-00
15.2" - 19.2"	4.5"	X Large	72819-00029-00
15.2" - 19.2"	2.75"	X Large (short)	72819-00030-00

FEA	ATURES AND BENEFITS
•	Anatomically contoured foam
•	Long foam and long Velcro® overlap
•	Medium density and firm density
•	Extra long and low collar size
•	Additional stockinette to cover the collar

ACTIMOVE* CE	RVICAL CO	LLAR - MEDIUM DE	NSITY
THE	1		
Neck Circ.	Neck	Size	Catalog Number
10.4" - 13.2"	2.5"	X Small	72819-00031-00
12" - 14.4"	3"	Small	72819-00032-00
13.2" - 15.6"	3.5"	Medium	72819-00033-00
14.4" - 16.8"	4"	Large	72819-00034-00
15.2" - 19.2"	2.75"	X Large (short)	72819-00035-00

#### Actimove® Cervical Comfort Cervical Collar



The Actimove® Comfort Cervical Collar is a convenient, anatomically shaped collar that maintains neck stability and restricts range of motion. Medium density with a blue stockinette cover for improved aesthetics

ACTIMOVE <sup>®</sup> CERVICAL - COMFORT CERVICAL COLLAR			
Neck Circ.	Size	Catalog Number	
10" - 13"	X Small	72859-00036-00	
11.75" - 14"	Small	72859-00037-00	
13" - 15.5"	Medium	72859-00038-00	
14" - 16.5"	Large	72859-00039-00	
15.3" - 18"	X Large	72859-00040-00	
15" - 19"	X Large (short)	72859-00041-00	

# Anatomically contoured foam Medium density foam Extra long and low collar size Additional stockinette to cover the collar Navy blue collar



#### **Universal Cervical Collar**

Unique collar design with an extended hook and loop closure adjusts to ensure a proper fit for neck sizes 13" to 19". Polyurethane foam is covered with a soft beige stockinette. Available in three widths and two densities. Traditional low contour design maintains the head in a neutral position.

UNIVERSAL CERVICAL COLLAR				
Low Contour Reg. Density	Low Contour Firm Density	Width	Inches	
10-121033	10-122032	Narrow	2.5" x 21"	
10-131030	10-132039	Regular	3.25" x 21"	
10-141037	-	Wide	4.25" x 21"	



#### Sized Cervical Collar

Polyurethane foam collars are die cut to specific sizes for improved fitting. Foam is covered with a soft beige stockinette. Available in two widths and two densities. Traditional low contour design maintains the head in a neutral position.

SIZED CERVICAL COLLAR - NARROW WIDTH				
	Reg. Density	Size	Inches	
NARROW	10-121431	Small	2.5" x 15"	
WIDTH 2.5"	10-121538	Medium	2.5" x 17"	
LOW CONTOUR	10-121635	Large	2.5" x 19"	
	10-121732	X-Large	2.5" x 21"	

SIZED CERVICAL COLLAR - REGULAR WIDTH			
	Reg. Density	Size	Inches
REGULAR	10-131438	Small	3.25" x 15"
WIDTH 3.25" LOW CONTOUR	10-131535 10-131632	Medium Large	3.25" x 17" 3.25" x 19"
	10-131739	X-Large	3.25" x 21"

#### Universal Cervical Collar with Microban®

The soft stockinette is treated with Microban<sup>®</sup>. Unique universal design with extended hook and loop closure. Available in three widths in regular density. Traditional low contour design maintains the head in a neutral position.

UNIVERSAL CERVICAL COLLAR		
Catalog Number	Width Size	Size
11-121003	Narrow	2.5" x 21"
11-131000	Regular	3.25" x 21"
11-141004	Wide	4.25" x 21"



#### Actimove® Clavicle Support

The Actimove® Clavicle Support is used for the management of clavicular fractures, using a figure-of-eight harness that maintains shoulder position. It is comprised of rigid straps which maintain support and soft foam padding for patient comfort.



ACTIMOVE® CLAVICLE SUPPORT		
Dimensions	Size	Catalog Number
< 25.5"	Pediatric	72819-00018-00
25.5" - 35.5"	Small	72819-00019-00
35.5" - 43.5"	Medium	72819-00020-00
43" - 51"	Large	72819-00093-00

	EATURES AND BENEFITS
•	Four sizes, pediatric size inclusive
•	Removable and washable pads for armpit area
•	Non-stretchable straps
•	Straps covered with soft foam
•	Velcro® hook and loop closure system
•	Re-positional D-rings with Velcro® at back

Velcro® is a registered trademark of Velcro U.S.A., Inc.

## Comfort® Clavical Brace Hook & Loop Closure

The Comfort® Clavical Brace provides easy application with hook and loop closures and comfort with wide, stockinette covered padding.

# Hook and loop closure for easy application Stockinette covered felt Comfortable, wide padding for proper immobilization Felt back pad helps prevent irritation

COMFORT" CLAVICAL BRACE		
Description	Size	Catalog Number
Infant	XX Small	501-XXS
Child	X Small	501-XS
Child/Youth	Small	501-S
Youth/Small Adult	Medium	501-M
Adult	Large	501-L
Large Adult	X Large	501-XL



#### Clavical Support with Four-Way Strap, Buckle Closures

Designed to help hold shoulders back and limit clavicle movement during treatment of clavicular fractures. Large felt back panel for greater support. Four separate buckle adjustments provide an excellent fit and maintain even support. Foam padded shoulder straps for greater comfort.

CLAVICLE SUPPO	RT with FOUR-W	AY STRAP, BUCKLE	CLOSURES
Catalog Number	Size	Fits Chest	
16-101106	XX-Small	15" - 20"	Infant
16-101307	X-Small	20" - 24"	Pediatric
16-101404	Small	24" - 30"	Youth
16-101501	Medium	30" - 36"	Small Adult
16-101608	Large	36" - 42"	Adult
16-101705	X-Large	42" - 48"	Large Adult



Measure around chest at the base of the sternum

## **Shoulder Abduction Support**

Constructed of sculpted foam to fit body contour. Designed to provide 45 to 70 degrees of abduction based on positioning. Cushion can be placed in the anterior, posterior or central position. Foam and elastic support is universally designed to fit a range of sizes. Cuffs are adjustable with hook and loop closures. Wide shoulder strap helps prevent migration.

SHOULDER ABDUCTION SUPPORT		
Catalog Number Size		
16-801009	Universal	





#### Actimove® Shoulder Abduction Sling

The Actimove® Shoulder Abduction Sling provides immobilization for rotator cuff repairs, capsular shifts, Bankhart repairs, Glenohumeral dislocation/subluxation and soft tissue strains/repairs.



ACTIMOVE* SHOULDER ABDUCTION SLING			
Size	Width	Length	Catalog Number
Small / Medium	up to 61"	< 15"	7285925
Large / X Large	up to 61"	> 15"	7285926
FEATURES AND	BENEFITS		
Allows arm to be positioned to individual requirement			
Positions the shoulder and arm in a neutral position			
Prevents post-operative internal rotation contractures			
Breathable fabric			
Quick release buckle on waist strap for fast and easy application			
Easy open sling pouch for access to arm			
Soft thumb loop prevents wrist drop			
Fits left or right arm			

#### Soft Form®



#### **Posture Control Brace**

The Soft Form® Posture Control Brace is designed to correct poor posture by gently holding the shoulders back in proper position. Elastic side panels provide support compression to stabilize the abdominal and lumbar regions for improved posture and alignment. The soft, laminate body material is comfortable against the skin and ventilated elastic side panels are breathable to keep the skin cool. May be worn underneath clothing, while sitting, or during activity to correct and control poor positioning or posture. Ideal for limiting forward and downward movement of the shoulders (slouching). Helps improve poor posture and relieves pain caused by arthritis. Minimizes a "stooped" posture often associated with osteoporosis. Encourages proper posture while seated for prolonged periods. Color: Beige.

Encourages proper posture
 Gently keeps shoulders back for posture control
 Soft, unique laminate body material
 Support compression to stabilize the abdominal and lumbar regions
 Two stays for additional support
 Cool, comfortable ventilated elastic side panels
 Criss-cross, plush elastic shoulder straps for increased adjustability
 May be worn under clothing, while sitting or during activity

SOFT FORM* POS	TURE CONTROL B	RACE
Catalog Number	Size	Fits
16-900400	Small	26" - 32"
16-900500	Medium	30" - 36"
16-900600	Large	34" - 40"
16-900700	X-Large	38" - 44"
16-900800	XX-Large	42" - 48"

Measure around the fullest part of the abdomen.

#### Posture Control Shoulder Brace

Brace is designed to gently keep the shoulders back for proper posture control. High quality canvas back with a soft cotton lining. Four plastic stays on the back panel for extra support. Foam padded shoulder straps for greater comfort.

POSTURE CONTRO	L SHOULDER BR	ACE	
Catalog Number	Size	Fits	
16-420406	Small	24" - 30"	Youth
16-420503	Medium	30" - 36"	Small Adult
16-420600	Large	36" - 42"	Adult
16-420707	X-Large	42" - 48"	Large Adult



Measure around chest at the base of the sternum

# Safe-T-Sport® Neoprene Shoulder Support

Stabilizes and protects the shoulder to prevent reinjury, while still allowing for flexibility. The sports neoprene material gives comfortable compression and four-way stretch for a proper fit and comfort. Neoprene retains the body's natural heat for therapeutic and soothing warmth to the shoulder joints and surrounding muscles. Adjustable chest and arm closures allow for a custom fit. Built-in pocket holds a resuable hot/cold gel pack (sold separately) in place over the tender joint. Easy to wear under clothing. For correct sizing, measure around the chest. Color: Black.



SAFE-T-SPORT* N	NEOPRENE SHOUL	DER SUPPORT
Catalog Number	Size	Fits
16-500400	Small	up to 40"
16-500500	Medium	41" - 50"
16-500600	Large	51" - 60"

Measure chest circumference

FEATURES AND BENEFITS
Protects and supports the shoulder while allowing for flexibility
Sports neoprene with a soft lining surrounds the shoulder for therapeutic warmth and comfortable compression
Built-in pocket holds a reusable hot/cold gel pack in place over the tender joint (gel packs sold separately)
Use an optional thermal gel pack to provide hot or cold therapy
Adjustable chest and arm closures

Inconspicuous wear - easy to wear under clothing



#### Actimove® Comfort Shoulder Immobilizer

The Actimove® Comfort Shoulder Immobilizer provides comfortable and safe shoulder support with optimum wrist mobility. It is used following surgery, fracture and for surrounding soft tissue injury. Adjustable straps plus an integrated pad provide additional patient comfort.



ACTIMOVE® COMFORT SHOULDER IMMOBILIZER		
Dimensions	Size	Catalog Number
13" <	X Small	72819-00094-00
13.5" - 15"	Small	72819-00039-00
14.5" - 17"	Medium	72819-00040-00
16.5" - 19"	Large	72819-00041-00
> 19"	X Large	72819-00042-00

FEATURES AND BENEFITS
Waist strap attachment at wrist level
Envelope opening system
Neck pad
Y-tab closure system
Thumb loop
Large envelope

#### Actimove® Eco Shoulder Immobilizer

The Actimove® Eco Shoulder Immobilizer limits abduction, extension and internal and external rotation of the humerus following surgery or fractures. A soft foam neck strap provides patient comfort and a hook and loop closure ensures easy application and removal.

**Dimensions** 

Thumb loop
Waist strap

13.5" - 15"



14.5" - 17"	Medium	72819-00066-00	
16.5" - 19"	Large	72819-00067-00	
> 19"	X Large	72819-00068-00	
FEATURES AND BENEFITS			
Large envelope			
Soft foam neck strap			
Non-stretch foan	n neck strap		
Hook fastener or	Hook fastener on strap		

**Catalog Number** 

72819-00065-00

Size

Small



#### Actimove® Sling & Swathe Shoulder Immobilizer

The Actimove® Sling and Swathe provides comfortable arm immobilization. Consisting of a blue hook

receivable material that allows for custom fittings. It is used following surgery, fracture and

for surrounding soft tissue injury.

ACTIMOV	'E <sup>*</sup> SLING & SWATHE		
Size	Length	Width	Catalog Number
Small	10" - 13"	up to 43"	72859-00001-00
Medium	13" - 16"	up to 55"	72859-00002-00
Large	16" +	up to 67"	72859-00003-00
FEATUR	ES AND BENEFITS		
• Long	chest strap/wrap		
Additional elbow strap			
Thin foam for envelope			
Soft foam and terry towel lining			
• 3 sizes			
Firm support			
Extremely soft and comfortable			
Wide	foam straps		



# Deluxe Sling & Detachable Swathe

Extra long sling gives greater range of fit. Made with a long hand catch to prevent wrist drop. Comfort foam sleeve on strap. Easy slide buckle adjustment. Fits right or left. Denim material. Color: Navy/White.

DELUXE SLING	AND DETACHABLE :	SWATHE
Catalog Number	Size	Length
16-602259	X-Small	11.75"
16-602453	Small	14"
16-602550	Medium	15.5"
16-602657	Large	17.75"
16-602754	X-Large	19.25"

Measure forearm from elbow to base of little finger





#### Elastic Shoulder Immobilizer



Popular lightweight, low profile support provides shoulder immobilization and reduces abduction and arm rotation. Made of firm elastic with a humeral and wrist cuff and loop lock closures. Color: White.

ELASTIC SHOULDER IMMOBILIZER			
Catalog Number	Size	Fits Chest	
16-701104	XX-Small	16" - 22"	Infant
16-701301	X-Small	20" - 26"	Pediatric
16-701408	Small	24" - 30"	Youth
16-701505	Medium	28" - 34"	Small Adult
16-701602	Large	32" - 38"	Adult
16-701709	X-Large	36" - 42"	Large Adult

### Shoulder Immobilizer Sling



Designed to immobilize the arm and shoulder. Envelope arm sling incorporates a body strap to prevent shoulder rotation. Slide buckle closures. Navy poly/cotton material. Color: Navy/White.

SHOULDER IMMOBILIZER SLING		
Catalog Number	Size	Length
28-802401	Small	11"
28-802508	Medium	13"
28-802605	Large	14"
28-802702	X-Large	15"

Measure forearm from elbow to base of little finger.



#### Actimove® Sling on a Roll

The Actimove® Sling on a Roll is a convenient immobilization system for the upper extremity. The combination of a strap with the Y-tab hook closures offers users and patients an easy-to-adjust arm support. Custom fit for pediatric or adult with a single product.

ACTIMOVE® SLING ON A ROLL		
Contents	Catalog Number	Quantity
Strap: 2.25" x 13 yds	72859-00018-00	2 rolls / box
Hook fasteners: 14 Y-tabs per set	Included	2 sets / box
FEATURES AND BENEFITS		
13 yards of strapping		

Y-tab hook fasteners allow easy positioning adjustments of the arm
 Strap width remains constant during use which guarantees less

• Stretch resistant foam keeps arm in desired position

compression in the neck





#### Actimove® Comfort Arm Sling

The Actimove® Comfort Arm sling provides effective vertical, lightweight support to the arm. Its additional neck pad provides extra patient comfort and a hook and loop closure system allows easy application and removal.

ACTIMOVE COMFORT	ARM SLING	
Dimensions	Size	Catalog Number
13.5" - 15"	Small	72819-00061-00
14.5" - 17"	Medium	72819-00062-00
16.5" - 19"	Large	72819-00063-00
> 19"	X Large	72819-00064-00

# Neck pad cushion Y-tab hook closure system Soft thumb loop Large envelope





#### Actimove® Eco Arm Sling

The Actimove® Eco Arm sling provides effective vertical, lightweight support to the arm.

Its soft neck strap provides patient comfort and a hook and loop closure system allows easy application and removal.



ACTIMOVE* ECO ARM SLING		
Dimensions —	Size	Catalog Number
13.5" - 15"	Small	72819-00057-00
14.5" - 17"	Medium	72819-00058-00
16.5" - 19"	Large	72819-00059-00
> 19"	X Large	72819-00060-00
FEATURES AND BENE	EFITS	
Large envelope		
Soft foam neck strap		

## Universal Arm Sling

The Universal Arm Sling can be folded to the proper length to accommodate most patients.

No trimming is necessary and one code covers most patients' needs.

Non-stretch foam neck strap

Hook fastener on strap

Thumb loop



UNIVERSAL ARM SLING	
Size	Catalog Number
Universal	1802-U
FEATURES AND BENEFITS	
Envelope style	
Slide buckle on strap for easy adjustn	nents
Hook and loop closure	



### Arm Sling Sized

Extra deep pouch for bulky castsNeck pad for patient comfort

• Thumb loop allows patient to relax hand

• Adjustable shoulder strap with hook and loop closure

This arm sling offers a variety of sizes, providing a great fitting product. The pouch is extra-deep to accommodate bulky casts and has a thumb loop to keep arm positioned. For comfort, it comes with an adjustable strap and neck pad.

ARM SLING SIZED		
Catalog Number		
1804-S		
1804-M		
1804-L		
1804-XL		





### **Universal Cradle Arm Sling**

Traditional cradle arm sling in an adjustable universal design to fit most patients. Easy, fast and simple application with double D-rings and loop lock closure. Navy poly / cotton material.

UNIVERSAL CRAD	LE ARM SLING	
Catalog Number	Size	Fits
28-502006	Universal	11" - 16"
Measure forearm from elbow to base of little finger		





#### Cradle Arm Sling



Traditional envelope arm sling comfortably distributes the weight of the arm following surgical or non-surgical treatment of the arm, hand or shoulder. Made in a strong, seamless quality design. Easy, fast and simple application with slide buckle adjustment. Comfort padding on shoulder strap. Poly/cotton material.

CRADLE ARM SLING		
Poly / Cot Catalog Number	Size	Sling Length
28-302406	Small	11"
28-302503	Medium	13"
28-302600	Large	14"
28-302707	X Large	15"

Measure forearm from elbow to base of little finger

#### ProLite® Universal Arm Sling



Soft, durable sling material comfortably supports the weight of the arm following injury or surgical treatment of the arm, hand or shoulder and can be worn with or without a cast. Features a padded, adjustable shoulder strap and thumb loop that limits wrist drop and prevents migration. Universal design fits a range of sizes and right or left arm. Black color hides dirt. Color: Black.

PROLITE* UNIVERSAL ARM SLING	
Catalog Number	Size
28-402000	Universal



#### **Universal Cuff & Collar**

Simple cuff and collar design comfortably holds weight of arm without a lot of bulk. Adjustable elastic wrist band. Heavyweight web strap with slide buckle adjustment. Comfort padding on strap. Fits a range of sizes. Color: White.

UNIVERSAL CRADLE	ARM SLING	
Catalog Number	Size	
28-911000	Universal	



#### GelBand® Arm Band

With 30 Durometer Thermal Viscoelastic Conforming Compression, this arm band provides the ideal in therapeutic support - conforming compression without restricting circulation. Use this support to treat common injuries like tennis or golfer's elbow. Provides relief from arthritis, Carpal Tunnel Syndrome, strains and sprains. Unique GelCell® can be cooled or heated for therapy on acute or chronic injuries. Colors: Black or Beige.

GELBAND* ARM BAND			
Beige Catalog Number	Black Catalog Number	Size	Fits Forearm
19-500430	19-500440	X-Small	7" - 11"
19-500630	19-500640	Universal	11" - 16"

Measure around forearm





### Epi-Sport®

# **Epicondylitis Clasp**

This specially designed clasp provides targeted compression to the long tendons of the forearm extensor and flexor muscles. Dual action compression works simultaneously to treat Lateral and/or Medial Epicondylitis (forehand and backhand tennis elbow) without restricting circulation. This low profile product is extremely lightweight and comfortable to wear without any pinch points. The shock absorbing liner also helps wick moisture away from skin. Color: Black.

EPI-SPORT® EPICONDYLITIS CLASP			
Catalog Number	Size	Fits	
19-600340	X-Small	8" - 9"	
19-600440	Small	9" - 10"	
19-600540	Medium	10" - 11"	
19-600640	Large	11" - 12"	
19-600740	X-Large	12" - 13.5"	

Measure around forearm 1.5" - 2" below elbow

## Canvas Tennis Elbow Strap



CANVAS TENNIS	ELBOW STRAP	
Catalog Number	Size	Fits
19-201002	Universal	11" - 16"
19-201400	X-Small	7" - 11"

Measure around upper forearm.

FEATURES AND BENEFITS		
Canvas outer with cool cotton lining		
Tapered design		
Loop lock closure/plastic guide		



# Neoprene Elbow Strap with Loop Lock

NEOPRENE ELBOV	W STRAP WITH LO	OOP LOCK	
Catalog Number	Size	Fits	
19-602609	Universal	11" - 16"	
19-602405	Small	7" - 11"	

Measure 3" below center of elbow
FEATURES AND BENEFITS
Neoprene for therapeutic warmth
Double compression pad
Tapered design
Loop lock closure/plastic guide



#### ProLite® Elbow Support Knitted Pullover

Provides support and compression to strained tendons/muscles or injured elbows. Made of a stretch knitted elastic that is a lightweight and breathable alternative to neoprene. Gives compressive support while maintaining flexibility. Color: Beige.

PROLITE* ELBOW SUPPORT KNITTED PULLOVER		
Catalog Number	Size	Fits
19-400430	Small	8" - 9"
19-400530	Medium	9" - 10"
19-400630	Large	10"- 11.5"
19-400730	X-Large	11.5" - 13"

#### FEATURES AND BENEFITS • Comfortable, firm support for weak or injured elbows

- Unique knit design for precise compression
- Measure around the arm just above elbow.











# ProLite® Compressive Elbow Support with Viscoelastic Insert

For swollen or tender elbows resulting from sprains, strains, overuse or sports injuries like tennis and golfer's elbow. Contains one viscoelastic insert positioned on the tendon below the elbow joint. Made of lightweight, breathable knit material. Will not retain excessive body heat. Color: White/blue.

PROLITE* COMPRESSI	VE ELBOW SUPPORT	with VISCOELASTIC INSERT
Catalog Number	Size	Fits
19-450300	X-Small	9" - 10"
19-450400	Small	10" - 11"
19-450500	Medium	11" - 12"
19-450600	Large	12" - 13"
19-450700	X-Large	13" - 14"

Measure 3" above elbow joint



#### Elbow Support Elastic Pullover



ELBOW SUPPORT ELA	STIC PULLOVER	
Catalog Number	Size	Fits
19-501407	Small	8.5" - 10"
19-501504	Medium	10" - 11.5"
19-501601	Large	11.5" - 13"
19-501708	X-Large	13.5" - 15"

Measure around upper forearm.


FEATURES AND BENEFITS
Durable lightweight knitted elastic
Support with flexibility
Seamless construction won't irritate skin



#### Neoprene Elbow Sleeve with Loop Lock

NEOPRENE ELBOW SLEEVE with LOOP LOCK			
Catalog Number	Size	Fits	
19-601406	Small	7" - 9"	
19-601503	Medium	9" - 10.5"	
19-601600	Large	10.5" - 13.5"	
19-601707	X-Large	13.5" - 15"	

Measure 3" below center of elbow





### ProLite® 6" Low Profile Wrist Splint



Provides stabilization and compression to weak or injured wrists. Ideal for sprains, strains, Carpal Tunnel Syndrome and after cast removal.

Removable palmar stay and two medial/lateral stays provide maximum support and stabilization, holding the wrist in a neutral position.

Material provides comfortable, continuous compression. Low contour at the palmar crease allows for full finger dexterity and improved grip strength.

Breathable spandex insert ensures a proper fit and easy application. Low profile and shorter 6" length for less bulk. Color: Black (helps hide soil).

PROLITE* 6" LOW PROFILE WRIST SPLINT			
Right Catalog Number	Left Catalog Number	Size	Fits Wrists
22-470410	22-470420	Small	5.5" - 6.5"
22-470510	22-470520	Medium	6.5" - 7.5"
22-470610	22-470620	Large	7.5"- 8.5"
22-470710	22-470720	X-Large	8.5" & up

Measure around the center of the wrist. Specify right or left.









#### Soft Fit Suede Finish Cock-Up Wrist Brace



Made of comfortable, perforated suede finish material with a soft, padded liner that absorbs moisture. Metal palmar stay holds the wrist in the neutral, "cock-up" position and prevents flexion. Two metal dorsal stays on the back of the brace help prevent extension. Recommended to support weak or injured wrists, treat wrist strains (slight to moderate wrist sprains), Carpal Tunnel Syndrome or after cast removal.

Easy, one-hand lace pull closure system. Colors: Beige or Black.

SOFT FIT SUEDE FINISH COCK-UP WRIST BRACE					
Beige Right Catalog Number	Beige Left Catalog Number	Black Right Catalog Number	Black Left Catalog Number	Size	Fits Wrists
22-150310	22-150320	22-150340	22-151340	X-Small	Under 4.5"
22-150410	22-150420	22-150440	22-151440	Small	4.5" - 5.5"
22-150510	22-150520	22-150540	22-151540	Medium	5.5"- 6.5"
22-150610	22-150620	22-150640	22-151640	Large	6.5" - 7.5"
22-150710	22-150720	22-150740	22-151740	X-Large	7.5" & up

Measure around the center of the wrist. Specify right or left.

#### ProLite® 8" Wrist Splint



Best Fitting Splint Available! Provides immobilization for weak or injured wrists.

Ideal for sprains, strains or after cast removal. Made of a cool and breathable open cell foam material with interlock poly cover. Provides comfortable, continuous compression. Low contour at palmar crease allows for full finger dexterity and improved grip strength. Removable palmar stay and two medial/lateral stays provide maximum support and immobilization. Loop lock closures and elastic insert assure proper fit and easy one hand application. Black color helps hide soil. Specify right or left.

PROLITE® 8" WRIST SPLINT			
8" Splint Right Catalog Number	8" Splint Left Catalog Number	Size	Fits Wrists
22-450311	22-450321	X-Small	4.5" - 5.5"
22-450412	22-450422	Small	5.5"- 6.5"
22-450513	22-450523	Medium	6.5" - 7.5"
22-450614	22-450624	Large	7.5" - 8.5"
22-450715	22-450725	X-Large	8.5" & up

Measure around the center of the wrist. Specify right or left.

# ProLite® Compression-Formed Composite Wrist Splint with Abducted Thumb

Recommended when immobilization of the wrist and/or thumb is needed.

Semi-rigid Thumb Spica portion is contoured to provide intimate support, immobilization and thumb abduction. The unique design of the foam spica protects the thumb from impact or shock without any uncomfortable binding or stitch seams. The splint has a metal palmar and a dorsal stay to prevent flexion and extension. Recommended for post-surgical treatment of wrist and thumb injuries. Especially for treatment of Basal Joint Arthritis, DeQuervain's Tenosynovitis, tendinitis, accumulative trauma disorders of the thumb and wrist, severe wrist sprains, strains / sprains of the thumb joint (MCP), advanced Carpal Tunnel Syndrome, or after cast removal. Color: Black. Specify right or left.



#### PROLITE® COMPRESSION-FORMED COMPOSITE WRIST SPLINT w/ABDUCTED THUMB

Abducted Thumb Right Catalog Number	Abducted Thumb Left Catalog Number	Size	Fits Wrists
22-460310	22-460320	X-Small	4.5" - 5.5"
22-460410	22-460420	Small	5.5"- 6.5"
22-460510	22-460520	Medium	6.5" - 7.5"
22-460610	22-460620	Large	7.5" - 8.5"
22-460710	22-460720	X-Large	8.5" & up

Measure around the center of the wrist. Specify right or left.

#### Deluxe Universal Wrist Brace

The unique, two part ambidextrous design of the Deluxe Universal Wrist Brace accommodates most patients in comfortable, breathable padded material.

DELUXE UNIVERSAL WRIST BRACE		
Length	Size	Catalog Number
8"	Universal	2008-U
11"	Universal	2007-U
8"	Universal - Dispenser of 12	D2008-U
11"	Universal - Dispenser of 12	D2007-U

#### FEATURES AND BENEFITS

- Removable volar aluminum stays allow for custom contouring
- Dorsal aluminum stays allow for firm support
- D-ring wrist strap offers additional support to limit wrist flexion and extension





#### Cinch-Lock Wrist & Forearm Brace

The Cinch-Lock Wrist and Forearm Brace provides added immobilization and comfort through its wrist cinch and generously padded thumb area. This brace is constructed from durable, breathable material for a long-lasting patient comfort. This brace is available in two conventional lengths, 8 inches and 11 inches.



8" CINCH-LOCK WRIST AND FOREARM BRACE - LEFT			
Wrist Circumference	Size	Catalog Number	
4.5" - 5.5"	X Small	2015-XS	
5.5" - 6.5"	Small	2015-S	
6.5" - 7.25"	Medium	2015-M	
7.25" - 8"	Large	2015-L	
8" - 8.75"	X Large	2015-XL	

# Malleable palmer and dorsal stay for maximum support Additional foam padding at thumb fold for greater patient comfort Wrist cinch strap provides additional immobilization

8" CINCH-LOCK WRIST AND FOREARM BRACE - RIGHT			
Wrist Circumference	Size	Catalog Number	
4.5" - 5.5"	X Small	2016-XS	
5.5" - 6.5"	Small	2016-S	
6.5" - 7.25"	Medium	2016-M	
7.25" - 8"	Large	2016-L	
8" - 8.75"	X Large	2016-XL	

11" CINCH-LOCK WRIST AND FOREARM BRACE - RIGHT			
Wrist Circumference	Size	Catalog Number	
4.5" - 5.5"	X Small	2002-XS	
5.5" - 6.5"	Small	2002-S	
6.5" - 7.25"	Medium	2002-M	
7.25" - 8"	Large	2002-L	
8" - 8.75"	X Large	2002-XL	

11" CINCH-LOCK WRIST AND FOREARM BRACE - LEFT			
Wrist Circumference	Size	Catalog Number	
4.5" - 5.5"	X Small	2001-XS	
5.5" - 6.5"	Small	2001-S	
6.5" - 7.25"	Medium	2001-M	
7.25" - 8"	Large	2001-L	
8" - 8.75"	X Large	2001-XL	



U.S. Patent:

U.S. PATENTS

D468,434; 6,517.501

#### Uni-Fit® Universal Wrist Splint

Made of soft foam with a wide palmar and movable dorsal stays. Loop lock closures ensure proper fit and easy application. Extremely low profile, ultra thin, breathable foam with smooth lycra lining. Unique design fits a wide range of sizes and left or right arm. Dorsal stays can be removed to give less immobilization as healing progresses. Eight inch length. Color: Black. Universal Size - Fits Right or Left.

# Catalog Number Fits Wrist 22-602065 5" - 11"

Measure around the center of the wrist

## CarpalMate® Wrist Support

CarpalMate® is a lightweight wrist support designed to maintain the wrist in a neutral position while allowing full finger dexterity. Prevents the wrist from dropping and prevents excess thumb and little finger opposition. Controlling these motions will prevent excess flexion of the wrist and finger flexor tendons thus avoiding overuse. The support is held in place by a wide elastic strap that provides comfortable compressive support at the wrist. Ideal for Carpal Tunnel Syndrome. Colors: Black or Beige.

CARPALMATE* WRIST SUPPORT					
Black	Beige				
Catalog Number	Catalog Number	Size			
22-140100	22-140000	Universal			

Fits right or left

CarpalMate\* is ideal for activities that require repetitive, stressful motions such as:

Typing

Data Entry

• Small Parts Assembly







#### **Elastic Pullover Wristlet**



Durable lightweight knitted elastic provides support with flexibility. Seamless construction will not irritate skin. Easy slip on style. Color: Beige.

ELASTIC PULLOV	ER WRISTLET	
Catalog Number	Size	Fits
22-103403	Small	5" - 6"
22-103500	Medium	6" - 7"
22-103607	Large	7"- 8.5"
22-103704	X-Large	8.5" & up

Measure around center of wrist.



#### **Elastic Wrist Wrap**



Lightweight elastic wrap provides support with flexibility. Wrap design fits a wide range of sizes. Available with or without thumb loop. Color: White.

ELASTIC WRIST	WRAP	
Elastic Wrist Catalog Number	Elastic w/Thumb Loop Catalog Number	Fits
22-105003	22-106002	4" - 11"

Measure around center of wrist.



### Neoprene Wrist Wrap

Neoprene wrist wrap with adjustable closure fits a wide range of sizes and fits right or left wrist. Provides lightweight support and soothing warmth to aching, fatigued or weak wrists. Made of lightweight neoprene with a soft nylon loop lining. Universal size. Color: Navy.

NEOPRENE WRIST	WRAP
Catalog Number	Fits
22-109009	4" - 11"
Measure around the center	er of the wrist.



# Soft Form® Elegant Wrist Support

Provides light support to weak or injured wrists. It is ideal for treatment of Carpal Tunnel Syndrome, arthritis, mild wrist strains or sprains and weak or aching wrists. Has one metal palmar stay to help maintain the wrist in a neutral position. The base is made of a unique, soft fleece-lined Spandex material. The sleek, low-profile design also has a low palmar crease for full finger dexterity. Color: Beige. Specify right or left.

SOFT FORM* ELEGANT WRIST SUPPORT				
Right Catalog Number	Left Catalog Number	Size	Fits Wrists	
22-560310	22-560320	X-Small	4.5" - 5.5"	
22-560410	22-560420	Small	5.5" - 6.25"	
22-560510	22-560520	Medium	6.25"- 7"	
22-560610	22-560620	Large	7" - 7.75"	
22-560710	22-560720	X-Large	7.75" & up	

Measure around the center of the wrist. Specify right or left.





# Perforated Suede Finish Wrist & Forearm Splint

Recommended when greater immobilization is needed. Stabilizes the wrist following minor fractures, sprains or after cast removal. Made of a durable perforated suede finish material with a cool cotton lining. One wide contoured palmar stay and two dorsal stays. Loop lock closures. Available in eight or ten inch lengths. Specify right or left. Color: Beige.

PERFORATED SUEDE FINISH WRIST & FOREARM SPLINT					
8" Right Catalog Number	8" Left Catalog Number	10" Right Catalog Number	10" Left Catalog Number	Size	Fits
22-442418	22-442427	22-441419	22-441428	Small	4.5" - 5.5"
22-442515	22-442524	22-441516	22-441525	Medium	5.5"- 6.5"
22-442612	22-442621	22-441613	22-441622	Large	6.5" - 7.5"
22-442719	22-442728	22-441710	22-441729	X-Large	7.5" & up

Measure around the center of the wrist. Specify right or left.



### Thumb Spica Soft Fit



Anatomically designed to ensure a comfortable fit. Perspiration is absorbed by a self wicking lining and increased ventilation is provided by perforated blue foam. Ideal for arthritis, DeQuervain's Syndrome, tendinitis and post thumb injuries. Closure straps are stretchable for ease of application and removal. Universal size fits either left or right.

THUMB	SPICA	SOFT	FIT	
Catalog	Numbe	r	Size	
25-12000	03		Universal	



#### Actimove® Gauntlet Thumb and Wrist Brace

The Actimove® Gauntlet is a functional thumb and wrist support for the treatment of light to moderate injuries or chronic instability problems at the thumb.

ACTIMOVE * GAUNTLET THUMB AND WRIST BRACE			
Size	Wrist Circ.	Length	Catalog Number
Small	< 6"	7"	72859-00010-00 (left)
Medium	6" - 7.5"	7"	72859-00011-00 (left)
Large	7.5" - 8.5"	8"	72859-00012-00 (left)
X Large	> 8.5"	8"	72859-00013-00 (left)
Small	< 6"	7"	72859-00014-00 (right)
Medium	6" - 7.5"	7"	72859-00015-00 (right)
Large	7.5" - 8.5"	8"	72859-00016-00 (right)
X Large	> 8.5"	8"	72859-00017-00 (right)

FEATURES AND BENEFITS
Soft Neoprene thumb cover
Pre-shaped and removable metal stays for wrist and thumb
Distal brace edge situated under the palmar crease
Brace is contoured at the ulnar aspect of the hand
Elastic material
High coverage of dorsal aspect of the hand
Reliable thumb support
Convenient application



#### Actimove® Eco Wrist Brace

The Actimove® Eco Wrist Brace is used for stabilization of the wrist without restriction of the hand function. It is designed to allow full function of the metacarpal joints and has two removable and adjustable stays.

ACTIMOVE * ECO WRIST BRACE				
Size	Wrist Circ.	Length	Catalog Number	
Small	< 6.5"	7"	72819-00043-00 (left)	
Medium	6.5" - 7.5"	7"	72819-00044-00 (left)	
Large	7.5" - 8.5"	8"	72819-00045-00 (left)	
X Large	> 8.5"	8"	72819-00046-00 (left)	
Small	< 6.5"	7"	72819-00047-00 (right)	
Medium	6.5" - 7.5"	7"	72819-00048-00 (right)	
Large	7.5" - 8.5"	8"	72819-00049-00 (right)	
X Large	> 8.5"	8"	72819-00050-00 (right)	

# Distal brace edge situated under palmar crease Brace is contoured at the ulnar aspect of the hand Elastic material Pre-shaped and removable stay High coverage of dorsal aspect of the hand







## Canvas Cock-Up Wrist Splint - Lace Up



CANVAS COCK-	UP WRIST SPLINT	- LACE UP	
6" Right Catalog Number	6" Left Catalog Number	Size	Fits
22-111314	22-111323	X-Small	4.5" - 5.5"
22-111411	22-111420	Small	5.5" - 6.5"
22-111518	22-111527	Medium	6.5"- 7.5"
22-111615	22-111624	Large	7.5" - 8.5"

Measure around the center of the wrist. Specify right or left.

#### FEATURES AND BENEFITS

- Durable canvas material with a soft cotton lining
- Palmar stay holds wrist in cock-up position and is malleable to increase or decrease the degree of cock-up
- Lace-up front has two pull-type hook and loop straps for easy one hand closure
- · Color: Light Beige

#### **Elastic Wrist Brace**



ELASTIC WRIST	T BRACE		
6" Right Catalog Number	6" Left Catalog Number	Size	Fits
22-201316	22-201325	X-Small	Under 4.5"
22-201413	22-201422	Small	4.5" - 5.5"
22-201510	22-201529	Medium	5.5" - 6.5"
22-201617	22-201626	Large	6.5"- 7.5"
22-201714	22-201723	X-Large	7.5" & up

Measure around the center of the wrist. Specify right or left.

#### FEATURES AND BENEFITS

- Lightweight elastic support with stretch compression
- Wide, malleable and removable palmar stay
- Stay positions wrist in cock-up or cock-down position
- Low profile design maximum grip strength and finger dexterity
- One-hand hook and loop closures
- Color: Beige



### Tether® Thumb Stabilizer

Helps support and limit the motion of the MCP joint of the thumb. Ideal for use on soft tissue injuries, ligament strains, Gamekeeper's thumb, osteoarthritis, and degenerative joint disease. Allows for full finger dexterity. Thermal lining provides warmth to the thumb area and joint. Wrist wrap is made of a soft, cool and breathable material. Made of conforming genuine leather. Color: Black.

TETHER® THUME	3 STABILIZER		
Right Catalog Number	Left Catalog Number	Size	Fits
25-160410	25-160420	Small	2.25" - 2.75"
25-160510	25-160520	Medium	2.75" - 3.25"
25-160610	25-160620	Large	3.25" - 3.75"





### ProLite® Neoprene Pull-On Thumb Support

Provides warmth and uniform compression to aching or injured

thumbs and wrists. Compression helps reduce pain, inflammation, and swelling while providing support to soft tissues. Indicated for arthritis, bursitis, tendinitis, repetitive motion or soft tissue injuries, and general muscle soreness. Colors: Black or Beige.

PROLITE® NEOPRENE PULL-ON THUMB SUPPORT					
Black Catalog Number	Beige Catalog Number	Size	Fits		
25-130400	25-130450	Small	5.5" - 6"		
25-130500	25-130550	Medium	6" - 6.5"		
25-130600	25-130650	Large	6.5" - 7"		
25-130700	25-130750	X Large	7" - 8"		

Measure circumference of wrist. Specify right or left.

- Full thumb coverage Neoprene applies soothing warmth to all the joints of the thumb.
- Contoured shape Allows for full use of hand.
- Soft nylon terry lining Works with soft stitching to give a seamless feel.



### C3<sup>™</sup> Deluxe Thumb Splint



Supports and limits the motion of the thumb joints. Also comfortably supports the wrist and prevents rotation. Ideal for DeQuervain's Syndrome, tendinitis, arthritis and sports injuries. Made of unique C3<sup>TM</sup> Tri-Laminate material making it cool and comfortable to wear. Allows for changes in the amount of swelling. The anatomically conforming stay and closures help restrict movement and provide maximum support. Universal size. Fits either left or right. Color: Black.

C3™ DELUXE THU	MB SPLINT	
Catalog Number	Size	Fits
25-170002	5" - 11"	Universal

Measure around the center of the wrist.

### HealWell<sup>™</sup> Grip Splint whfo

Wrist / Hand / Finger Orthosis



Maintains a neutral position of the hand, reducing wrist flexion and increasing extension of the fingers. Holds the hand open, preventing finger flexion and damage to the palmar skin.

#### Indications:

- Mild/moderate non-fixed contractures of the hand, wrist and fingers
- Therapeutic positioning of paralyzed or injured hand or wrist
- Particularly useful after injury, stroke or disease to prevent contracture and protect the hand and wrist

HEALWELL™ GRIP	SPLINT WHFO	
Catalog Number	Size	Fits
25-800000	Universal	left or right

The HealWell\* Grip Splint (WHFO) is supplied with two palmar rolls in 1.5" and 1.75" diameters to accommodate various levels of finger contracture or hand sizes and/or offer progressive extension of the fingers.



### **STAX**® Finger Splints

STAX® Finger Splints are designed to support the distal joint of the finger extension while permitting unrestricted movement of the proximal inter-phalangeal joint.

STAX®	FINGER SPI	LINTS		
Size	Color	Catalog Number	Unit	Kit Quantity*
1	<ul><li>Beige</li></ul>	PS6-1	1	4 splints / kit
2	<ul><li>Beige</li></ul>	PS6-2	1	4 splints / kit
3	<ul><li>Beige</li></ul>	PS6-3	1	4 splints / kit
4	<ul><li>Beige</li></ul>	PS6-4	1	5 splints / kit
5	<ul><li>Beige</li></ul>	PS6-5	1	4 splints / kit
5.5	<ul><li>Beige</li></ul>	PS6-5.5	1	4 splints / kit
6	<ul><li>Beige</li></ul>	PS6-6	1	4 splints / kit
7	<ul><li>Beige</li></ul>	PS6-7	1	1 splint / kit



STAX®	FINGER SP	LINTS		
Size	Color	Catalog Number	Unit	Kit Quantity*
1	O Clear	PS6C-1	1	4 splints / kit
2	O Clear	PS6C-2	1	4 splints / kit
3	O Clear	PS6C-3	1	4 splints / kit
4	O Clear	PS6C-4	1	5 splints / kit
5	O Clear	PS6C-5	1	4 splints / kit
5.5	O Clear	PS6C-5.5	1	4 splints / kit
6	O Clear	PS6C-6	1	4 splints / kit
7	O Clear	PS6C-7	1	1 splint / kit

#### FEATURES AND BENEFITS

- Automatically elevates the distal joint to its proper extended position
- The kit is ideal for any location where a multitude of splints in a variety of sizes needs to be kept

STAX* FINGER	SPLINTS	KIT - AS	SHOWN	ABOVE
Color			(	Catalog Number
Clear			F	PS5
Beige			F	PS5C

<sup>\*</sup> The PS5 and the PS5C splint kits contain the assortments listed under the kit quantity

### Plastalume® Finger Splints

With their large base for stabilization, Plastalume Finger Splints provide treatment for trauma, dislocations and fractions on the extensive surface of fingers and thumbs.

6 each of A-1 thru AB-8

PLASTALUM	E FINGER SPLINTS	
Length	Catalog Number	Quantity
1.25"	A-1	12 splints / pack
2.25"	A-2	12 splints / pack
3.25"	A-3	12 splints / pack
4.25"	A-4	12 splints / pack
3.25"	AB-4	12 splints / pack
4.25"	AB-5	12 splints / pack
5.25"	AB-6	12 splints / pack
6.25"	AB-7	12 splints / pack
7.25"	AB-8	12 splints / pack
8.25"	AB-9	12 splints / pack

#### Catalog Number Quantity 27 3 each of A-1 thru AB-8



- Malleable aluminum with soft live foam that can be shaped easily
- Ball-shaped end prevents splint from rolling out of position and helps comform to finger
- Vinyl covering covers aluminum and cleans easily
- Latex-free and radiolucent, eliminating need for resplinting
- AB series designed with enlarged MP joint pads for additional patient comfort

### **Finger Splints**

A variety of standard splint styles are available. Each splint is constructed with malleable aluminum and foam padding for patient comfort unless otherwise indicated. Most finger splints are available in two packaging configurations - either individually bagged in 6-pack box or bulk packed in bag of 12.



Baseball Finger Splint
Stabilizes fractures of the distal phalanx.

BASEBALL FING	ER SPLINT		
Individually Bagged 6-Pack Catalog Number	Size	Length	Bulk Pack Bag of 12 Catalog Number
25-2014PK	Small	3.5"	25-201A405
25-2014FK	Jillali	3.5	25-201A405
25-2014PK 25-2015PK	Medium	4"	25-201A507



**Curved / Gutter Finger Splint** 

For treatment of transverse fractures of the phalanges.

CURVED/GUTTER	FINGER SPLINT		
Individually Bagged 6-Pack Catalog Number	Size	Length	Bulk Pack Bag of 12 Catalog Number
25-4014PK	Small	2.5"	25-401A404
25-4015PK	Medium	4"	25-401A507
25-4016PK	Large	6"	25-401A609



#### Four Prong / Protector Finger Splint

For total finger immobilization and protection with easy application.

FOUR PRONG/PRO	OTECTOR FINGER	SPLINT	
Individually Bagged 6-Pack Catalog Number	Size	Length	Bulk Pack Bag of 12 Catalog Number
25-7014PK	Small	1.5"	25-701A403
25-7015PK	Medium	2"	25-701A505
25-7016PK	Large	3"	25-701A607

### Finger Splints

### **Toad Finger Splint**

Maintains phalangeal joints in position without use of tape.

TOAD FINGER SE	PLINT		
Individually Bagged 6-Pack Catalog Number	Size	Length	Bulk Pack Bag of 12 Catalog Number
25-3014PK	Small	2.125"	25-301A403
25-3014PK 25-3015PK	Small Medium	2.125" 2.75"	25-301A403 25-301A505



#### Fold-Over / Cot Finger Splints

For total finger immobilization and protection with easy application. Good finger tip protection.

FOLD	FOLD OVER / COT FINGER SPLINT				
Ba 6-I	ridually gged Pack g Number	Size	Length	Bulk Pack Bag of 12 Catalog Number	
25-5	014PK	Small	2"	25-501A403	
25-5	015PK	Medium	3"	25-501A507	



#### Finger / Thumb Splint

Lewin style splint for finger/thumb permits any degree of flexion/hyperextension. Made of a malleable metal stay that can be contoured to the finger/thumb. Soft foam liner for increased comfort. Hook and loop closures for stability. Moveable tabs can be moved anywhere on the splint for greater adjustability. Available in 6-Packs Only.

FINGER/THUMB	SPLINT	
Individually Bagged 6-Pack Catalog Number	Size	Length
25-1010PK	Universal	10"
25-1022PK	Pediatric	5"
25-1024PK	Small	7"
25-1025PK	Medium	8"
25-1026PK	Large	10"





### Premium Woven Series, Latex-Free





### Two-Panel Surgical Rib Belts

#### For Women and Men

Two-panel surgical rib belts help stabilize rib and sternum fractures by limiting expansion through compression. They provide comfort and encourage more controlled breathing to help reduce pain. Rib belts also provide support and compression to the muscles and soft tissues of the rib cage weakened by strain, trauma, overuse, inactivity, or surgery. Female rib belt is designed with a contoured fit around the bust line. Universally sized. Color: White.

PREMIUM WOVEN SERIES, LATEX-FREE TWO-PANEL SURGICAL RIB BELTS				
Male Catalog Number	Female Catalog Number	Size	Fits Chest	
34-4100SL	34-4200SL	Universal	24" - 48"	

Measure around chest at base of the sternum

### Premium Woven Series, Latex-Free Three



For Women and Men

THREE-PANEL SURGICAL BINDERS (9" HEIGHT) provide support and compression to abdominal muscles weakened by strain without restricting breathing or muscle redevelopment.

**FOUR-PANEL SURGICAL BINDERS (12" HEIGHT)** provide support and compression to abdominal muscles weakened by strain without restricting breathing or muscle redevelopment. Tall height for taller patients or long line surgeries.



Measure around fullest part of the abdomen.





### Soft Form® Hernia Belt

Provides gentle relief from reducible inguinal hernias. Manufactured using lightweight materials with no uncomfortable or hard to use buckles or snaps. Provides constant, comfortable and adjustable pressure to the hernia. Two uniquely shaped foam compression pads provide gradual pressure and support to the weakened muscles with focused compression on the hernia. Either pad can be removed to treat a single hernia. Color: Beige.

SOFT FORM® HERNIA BELT				
Catalog Number	Size	Waist Measurement		
67-350400	Small	30" - 35"		
67-350500	Medium	35" - 41"		
67-350600	Large	41" - 46"		

FEATURES AND BENEFITS
Easy hook and loop closure
Adjustable waist band
Adjustable perineal straps
Flexible foam cushions
Neutral beige color can be worn inconspicuously under most clothing





6,422,242B1

### Soft Form® Hernia Brief

Offers gentle relief from reducible inguinal hernias by providing adjustable and comfortable compression. Low in profile and white in color, it can be worn inconspicuously under most clothing. May be worn in lieu of a regular undergarment and, with hernia support built into the brief, eliminates the need to wear both a regular undergarment and a bulky truss. The removable foam cushions are designed to provide gradual compression around the affected area with focused compression on the hernia. Either pad can be removed to treat a single hernia.





U.S. Patents: D482.512 6,622,719B1

SOFT FORM® HERNIA BRIEF				
Catalog Number	Size	Waist Measurement		
67-500400	Small	28" - 30"		
67-500500	Medium	32" - 34"		
67-500600	Large	36" - 38"		
67-500700	X-Large	40" - 42"		
67-500800	XX-Large	44" - 46"		

FEATURES AND BENEFITS	
Built-in hernia support	
Removable foam cushions	
Easily adjustable for the ideal amount of compression	
Treats right, left or double hernias	
Advanced plush, micro-soft materials for extra comfort	
Anatomically designed - creates compression only where needed	
Easy and comfortable to wear	



### Actimove® Abdominal Support



The Actimove® Abdominal Support is used following surgery as it gives support to the soft tissue around the abdomen. This lightweight product is washable and comfortable to wear, with an elastic strap that allows air circulation between the skin and product.

ACTIMOVE® ABDOMINAL SUPPORT				
Waist Circ.	Height	Size	Catalog Number	
22" - 47"	9"	Small Short	72819-00000-00	
40" - 63"	9"	Large Short	72819-00001-00	
22" - 47"	12"	Small Long	72819-00002-00	
40" - 63"	12"	Large long	72819-00003-00	

#### FEATURES AND RENEELTS

- V-shaped elastic strap
- Four sizes: Two different heights and two strap lengths available
- Lightweight elastic strap with ventilation sections
- Soft, flexible elastic strap edges
- Velcro® closure system

### **Universal Abdominal**

Binder 9" height



## UNIVERSAL ABDOMINAL BINDER 9" HEIGHT Universal Catalog Number Size Fits 34-603005 Long 36" - 60" 34-604004 Short 26" - 48"

Measure around fullest part of the abdomen.

## FEATURES AND BENEFITS 9" height Elastic/foam construction Plastic stays for increased stability Fits range of sizes



### Soft Form® Lumbar Sacral Support

#### 11" with Flexible Stays

Indicated for muscle strain, lower back pain and discomfort, lumbar instability, lumbar spine sprain/strain and sacral contusion.

Stabilizes and supports the lumbar sacral region to help provide relief from lower back pain and discomfort. Unique multi-layered soft laminate provides excellent support and comfort with a soft cotton/nylon lining. Sewn in a five-piece design for an ultimate tapered fit without binding or rolling. Elastic compression straps in a criss-cross design give compression to the entire lumbar region. Lightweight stays in the back provide stabilizing support with flexibility. Color: Beige.

SOFT FORM® LUMBAR	SACRAL SUPPOR	RT 11" W/FLEXIBLE STAYS
Catalog Number	Size	Fits
31-550400	Small	28" - 32"
31-550500	Medium	33" - 36"
31-550600	Large	37" - 40"
31-550700	X-Large	41" - 44"
31-550800	XX-Large	45" - 48"
31-550900	XXX-Large	49" - 52"
31-550950	XXXX-Large	53" - 56"





### Soft Form® Lumbar Sacral Support

#### 11" with Contoured Rigid Stays

#### For Men and Women

Indicated for muscle strain, lower back pain and discomfort, lumbar instability, lumbar spine sprain/strain, sacral contusion and lumbar disc dysfunction.

Similar to the above product with two removable contoured metal back stays which can be easily shaped for a customized fit. Provides more rigid and firm support.

Color: Beige.

SOFT FORM® LUMBAR	SACRAL SUPPORT 11"	w/CONTOURED RIGID STAYS
Catalog Number	Size	Fits
31-560400	Small	28" - 32"
31-560500	Medium	33" - 36"
31-560600	Large	37" - 40"
31-560700	X-Large	41" - 44"
31-560800	XX-Large	44" - 48"
31-560900	XXX-Large	49" - 52"
31-560950	XXXX-Large	53" - 56"

 $\label{eq:measure around the fullest part of the abdomen.}$ 





## 7" Lumbar Sacral Support®





### Cool - Lightweight

Ideal for Lightweight Lumbar Support

Lower height is great for shorter and smaller patients. Made of high quality seven inch high, latex free elastic. Four flexible plastic stays in the back provide support to the lumbar region. Four-inch wide tapered breathable elastic bands overlap and give a contouring shape to most body types. Soft plush finish. Easy hook and loop closure.

Color: White

7" LUMBAR SACRAL SUPPORT				
Catalog Number	Size	Fits		
31-701400	Small	28" - 32"		
31-701500	Medium	33" - 36"		
31-701600	Large	37" - 40"		
31-701700	X-Large	41" - 44"		


Measure around hips just below waistline

### **Lumbar Sacral Support**









Measure around hips just below waistline



### 9" Neoprene Lumbar Sacral Support

This support provides compression to the lumbar sacral and abdominal regions, stabilizing and supporting these areas to provide relief from lower back pain and discomfort.

#### Indications:

- Occasional or Chronic Back Pain
- Muscle Strain
- Lumbar Sprain/Strain
- Sacral Contusion





9" NEOPRENE	LUMBAR SACRAL	SUPPORT
Catalog Number	Size	Fits
31-722400	Small	28" - 32"
31-722507	Medium	33" - 36"
31-722604	Large	37" - 40"
31-722701	X-Large	41" - 44"
31-722808	XX-Large	45" - 48"

Measure around hips just below the waistline.

### ProLite<sup>®</sup> Mesh•Loc<sup>™</sup> Lumbar Support

This support provides abdominal compression and lift to the lumbar region, stabilizing and supporting these areas to relieve lower back pain and discomfort. Ideal for occasional or chronic back pain, lumbar sprains/strains, lumbar instability, muscle weakness and postural support.



#### Indications:

- · Occasional or Chronic Back Pain
- Lumbar Sprain/Strains
- Lumbar Instability
- Muscle Weakness and Postural Support

FE	EATURES AND BENEFITS
•	Patented strapping method provides hydrostatic lift to instantly relieve the load and stress on the spine
•	Tapered panels contour the body for excellent compression and support
•	Rigid stays provide strength and stability to the lumbar region
•	Easy application with simple hook and loop closures
•	Black color for extended wear

PROLITE® MESH•LOC™ LUMBAR SUPPORT					
Catalog Number	Size	Fits			
31-330400	Small	24" - 28"			
31-330500	Medium	29" - 33"			
31-330600	Large	34" - 39"			
31-330700	X-Large	40" - 44"			
31-330800	XX-Large	45" - 49"			
31-330900	XXX-Large	50" - 55"			

Measure around hips just below the waistline.



### Deluxe Lumbar



Support 11" Height

Alleviates pain caused from muscle strain by relieving stress on the lower spine. This popular support is exceptionally strong and durable. Back panel is made of perforated vinyl with a soft cotton lining. Strength is achieved with four contoured rigid stays that can be shaped for a customized fit. Double front closures are made of high quality elastic and soft foam. Color: White.


DELUXE LUMBAR SUPPORT 11" HEIGHT						
Catalog Number	Size	Fits				
31-101403	Small	25" - 34"				
31-101607	Medium - Large	32" - 42"				
31-101704	X-Large	38" - 48"				

Measure around hips just below waistline



### Sacral Iliac Support

with Compression Pad 6" Height





SACRAL ILIAC SUPPORT W/COMPRESSION PAD 6" HEIGHT				
Catalog Number	Size	Fits Hips		
31-801008	Universal	33" - 48"		
31-801406	Small	28" - 32"		
31-801503	Medium	33" - 36"		
31-801600	Large	37" - 40"		
31-801707	X-Large	41" - 44"		
31-801804	XX-Large	45" - 48"		

Measure around hips just below waistline.



## Criss-Cross Lumbar Sacral Support 10" Height

Criss-cross wrap around ten inch high lumbar sacral support. All elastic construction for comfortable compression. Six-inch wide overlapping front elastic bands. Three-inch wide elastic criss-cross bands in back. Flexible plastic stays prevent the support from rolling. Easy hook and loop closure. Universal design. Color: White.

CRISS-CROSS	LUMBAR SACRAL	SUPPORT 10" HEIGHT
Catalog Number	Size	Fits
31-410000	Universal	33" - 48"

Measure around hips just below waistline.





### Dorsal Lumbar 19" Support

Provides compression to the dorsal lumbar and abdominal areas. Back panel is 19" high with four rigid contoured stays for extra support. Six-inch wide overlapping front elastic bands. Easy hook and loop closure, with adjustable non-elastic web straps that close with locking buckles. Adjustable shoulder straps. Color: White.

DORSAL LUMBAR 19" SUPPORT					
Catalog Number	Size	Fits			
31-950401	Small	28" - 32"			
31-950508	Medium	33" - 36"			
31-950605	Large	37" - 40"			
31-950702	X-Large	41" - 44"			
31-950809	XX-Large	45" - 48"			

Measure around hips just below waistline.







### Soft Form®

### **Maternity Support Belt**

Comfortably supports the lumbar and abdominal regions to help provide relief of lower back pain and discomfort during pregnancy. Made of a unique multi-layered soft laminate that provides excellent support and comfort with a soft cotton/nylon lining. Abdominal lift attachment provides lift without undesirable pressure and can expand to accommodate growth during various stages of pregnancy. Easy to wear under most clothing. Color: Beige.



U.S. Patents D437.416



SOFT FORM* MATERNITY SUPPORT BELT					
Catalog Number	Size	Dress Size Before Pregnancy			
31-570400	Petite	2 - 6			
31-570500	Universal (Medium)	6 - 14			
31-570600	Large	14 - 20			

### Thermal Lumbar Support



U.S. Patent D496.108

Reusable hot/cold gel pack for treatment of minor back aches and pains resulting from sprains, muscle soreness and strains, cramping, or contusion (bruising). Stabilizes and supports the lumbar sacral region to help provide relief from lower back pain and discomfort. Color: Black.

THERMAL LUMBAR SUPPORT				
Catalog Number	Size	Fits		
31-7504	Small	28" - 32"		
31-7505	Medium	33" - 36"		
31-7506	Large	37" - 40"		
31-7507	X-Large	41" - 44"		

Measure around the fullest part of the abdomen.







- Reusable thermal gel provides penetrating hot or cold therapy.
- Mesh pocket holds reusable hot/cold thermal gel insert over lower lumbar region.
- Support gives an ultimate tapered fit without binding or rolling.
- Plush, ventilated elastic panels lightweight elastic with a unique ventilated filament stitch provides cool and comfortable compression. Panels overlap and give a contouring shape to most body types.
- Easy hook and loop closure.

### V-Loc<sup>™</sup> LSO Spinal Orthosis

When the strapping system is closed over the V-Shaped rigid anterior panel, it provides intra-abdominal muscle compression and provides a lift to hold the weight of the abdominal area. By lifting and holding the weight, pressure is relieved from lumbar discs reducing the associated pain. The posterior and anterior panels work together, giving opposing support to stabilize the lumbar region. The result is a supported and pressure free lumbar spine held in a healthy, correct posture.

U.S. Patents: 5,207,636; 5,267,948

#### Indications:

- Chronic Lower Back Pain
- Lumbar Strains or Sprains
- Muscle Weakness
- Lumbar Instability
- Post Operative Rehabilitation
- Postural Support

V·LOC™ LSO SPI	INAL ORTHOSIS	
Catalog Number	Size	Fits
32-1504	Small	24" - 28"
32-1505	Medium	29" - 33"
32-1506	Large	34" - 39"
32-1507	X-Large	40" - 44"
32-1508	XX-Large	45" - 49"
32-1509	XXX-Large	50" - 55"

Measure around waist at the navel.

- Unloads pressure from lumbar discs to alleviate lower back pain and discomfort
- V-Shaped anterior panel creates true hydrostatic lift of lower abdominal muscles
- Rigid posterior and anterior panels stabilize the lower spine and control rotation and pelvic tilt
- 10" high back panel covers L1-L5 vertebrae
- Elastic insets provide a conforming fit around the hips
- Soft, padded lining for maximum comfort and compliance
- Ultra-thin low profile for comfortable, inconspicuous wear
- Fully adjustable, patented closure system does not rely on elastic for extended product life
- Provides pain relief during activity and allows for immediate participation in rehabilitation exercises
- For men and women



### Actimove® Genu Clip Knee Immobilizer with Clip Closure System

The Actimove® Genu Clip Knee Immobilizer with clip closure system is convenient and easy to apply. It stabilizes the knee joint in a comfortable way after acute injury or surgical interventions.



ACTIMOVE <sup>®</sup> GENU CLIP - SMALL					
Length	Size	Height 🔐	Thigh Circ.	Catalog Number	
12"	P1	< 3'10"	< 12.5"	72859-00045-00	
16"	P2	< 4'7"	< 16.5"	72859-00046-00	

ACTIMOVE ®	ACTIMOVE* GENU CLIP - LARGE				
Length	Size	Height	Thigh Circ.	Catalog Number	
20"	A1-S	≤ 5'3"	Max 20.5"	72859-00047-00	
20"	A2-L	≤ 5'3"	Max 20.5"	72859-00048-00	
24"	B1-S	≤ 5'10"	Max 20.5"	72859-00049-00	
24"	B2-L	≤ 5'10"	Max 25.5"	72859-00050-00	
28"	C1-S	≥ 5'10"	Max 20.5"	72859-00051-00	
28"	C2-L	≥ 5'10"	Max 25.5"	72859-00052-00	



#### FEATURES AND BENEFITS

- Lightweight and slim design
- Preshaped metal stay
- Positionable lateral stays
- Terry net liner
- Convenient straps with a mechanical clip system
- Easy application and adjustment
- Quick and easy removal

### Actimove® Tri-Panel Knee Immobilizer



The Actimove® Tri-Panel Knee Immobilizer supports the knee and helps to reduce pain following trauma or surgery. It has adjustable panels and detachable hook-and-loop closures for easy application and removal.

ACTIMOVE® TRI-PANEL KNEE IMMOBILIZER					
Length	Size	Thigh Circ.	Catalog Number		
20"	A1-S	<25"	72819-00014-00		
20"	A2-L	>25"	72819-00015-00		
24"	B1-S	<25"	72819-00016-00		
24"	B2-L	>25"	72819-00017-00		

- Three panel design
- Adjustable side panels
- Terry net liner
- Wide lateral straps
- Two T-shaped posterior straps
- Hook-and-loop closure

### Bel Knee Immobilizer

The versatile and economic Bel Knee Immobilizer uses a universal closure system and removable, repositionable aluminum stays to fit most sizing needs.

BEL KNEE IMMOBILIZER				
Length	Catalog Number			
13"	BK13-U			
16"	BK16-U			
19"	BK19-U			
23"	BK23-U			

FEATURES AND BENEFITS
Lightweight and slim design
Preshaped metal stay
Positionable lateral stays
Terry net liner
Foam exterior
Convenient straps with a mechanical clip system
Easy application and adjustment
Quick and easy removal



### Actimove® Eco Knee Immobilizer

The Actimove® Eco Knee Immobilizer supports the knee and helps to reduce pain following trauma or surgery. It contains removable stays that provide a good anatomical fit, and detachable hook and loop closures with D-rings allow patient adjustments.

ACTIMOVE® ECO KNEE IMMOBILIZER - SMALL				
Length	Size	Height	Thigh Circ.	Catalog Number
12"	P1	< 3'10"	< 12.5"	72819-00072-00
16"	P2	< 4'7"	< 16.5"	72819-00073-00

ACTIMOVE <sup>®</sup> ECO KNEE IMMOBILIZER - LARGE				
Length	Size	Height	Thigh Circ.	Catalog Number
20"	Small	≤ 5'3"	Max 20.5"	72819-00008-00
20"	Large	≤ 5'3"	Max 25.5"	72819-00009-00
24"	Small	≤ 5'10"	Max 20.5"	72819-00010-00
24"	Large	≤ 5'10"	Max 25.5"	72819-00011-00
28"	Small	≤ 5'10"	Max 20.5"	72819-00012-00
28"	Large	≤ 5'10"	Max 25.5"	72819-00013-00

ACTIMOVE <sup>®</sup> ECO KNEE IMMOBILIZER - EXTENSION TABS				
Size	Actimove Eco	Sales Unit	Catalog Number	
Small	A1 - A2	10 tabs / bag	72819-00074-00	
Medium	B1 - B2	10 tabs / bag	72819-00075-00	
Large	C1 - C2	10 tabs / bag	72819-00076-00	



FE	ATURES AND BENEFITS		
•	Pre-shaped metal stays		
•	Positionable lateral stays with Velcro®		
•	Lightweight and slim design		
•	Terry net liner		
•	Four detachable straps		
•	Foam extension tabs for larger sizes		

### Universal Wrap-Around Knee Support



Unique gel pads absorb heat and cold, creating a comfortable environment for the open patella while the neoprene, universal, closure construction allows the Universal Wrap-Around Knee Support to fit most patients.

UNIVERSAL WRAP-AROUND KNEE SUPPORT				
Size	Catalog Number			
Universal	E202000			
FEATURES AND BENEFITS				
1/s" thick neoprene with brushed nylon for universal positioning of gelpads and easy application				
Gel surrounds and conforms to the	affected area, making an ideal padding			

### Actimove® Post Trauma Knee Dressing



The Actimove® Post Trauma Knee Dressing replaces the multi-layer bandaging system used for the traditional "Jones Dressing". Ideal for arthroscopic interventions or trauma of the knee, providing compression and moderate support where needed.

Dimensions	Catalog Number	Quantity	
Universal	72859-00000-00	1 each	
FEATURES AND BENEFITS			
Soft foam and to	erry net liner		
<ul><li>Soft foam and te</li><li>Hook and loop c</li></ul>			
Hook and loop c			

### **Breathable Universal**





Popular universal knee immobilizer has adjustable side panels to fit a range of sizes. Made of a blue perforated mesh material laminated to a soft, light gray medical grade foam. Lateral and medial stays increase stability. Immobilization created with contoured posterior stays. Easy loop lock closures. Universal size.

BREATHABLE UNIVERSAL TRI-PANEL FOAM KNEE IMMOBILIZER					
14" Length Catalog Number	16" Length Catalog Number	18" Length Catalog Number	20" Length Catalog Number	22" Length Catalog Number	24" Length Catalog Number
37-614008	37-616006	37-618004	37-620009	37-622007	37-624005

Universal Size fits 13" - 26". Measure around leg 4" above center of kneecap.

## Breathable Universal Cutaway Knee Immobilizer

Practical universal knee immobilizer can be cut to fit a range of sizes. Movable medial and lateral stays can be placed in desired position. Made of a blue perforated mesh material laminated to a soft, light gray medical grade foam. Immobilization created by posterior stays. Easy loop lock closures. Includes a removable elastic patella closure. Universal size - trimmable.

BREATHABLE UNIVER	SAL CUTAWAY KNEE II	MMOBILIZER
18" Length Catalog Number	20" Length Catalog Number	24" Length Catalog Number
37-718003	37-720008	37-724004

Universal Size fits up to 23". Measure around leg 4" above center of kneecap.



### Tri-Panel Knee Immobilizer

Popular universal knee immobilizer has adjustable side panels to fit a range of sizes. Constructed of a medical grade foam with a soft loop liner for patient comfort. Lateral and medial stays increase stability. Immobilization created with contoured posterior stays. Easy loop lock closures. Universal size.

TRI-PANEL KNEE IMMOBILIZER					
14" Length Catalog Number	16" Length Catalog Number	18" Length Catalog Number	20" Length Catalog Number	22" Length Catalog Number	24" Length Catalog Number
38-140002	38-141004	38-142006	38-143008	38-144000	38-145002

Universal Size fits 13" - 26". Measure around leg 4" above center of kneecap.



### GelBand® Patella Strap



U-Point™ Gel Insert focuses pressure on the tendon below the kneecap

This patella strap provides the ideal in therapeutic support - conforming compression without restricting circulation. Use this support to treat symptomatic pain from Chondromalacia Patella Syndrome, Patellar Tendinitis, and Osgood-Schlatter's disease. Strap applies uniform compression which helps guide the patella and improve tracking. Includes a uniquely designed viscoelastic insert that provides consistent and focused compression.

Angulated buckles allow strap to reside properly for minimal migration in either inferior or superior applications. Colors: Black or Beige.

GELBAND* PATELLA STRAP				
Black Catalog Number	Beige Catalog Number	Size	Fits	
37-500040	37-500030	Universal	11" - 18"	

Measure around kneecap.

### Safe-T-Sport® Sports Neoprene Patella Support



Neoprene retains the body's natural heat to provide therapeutic and soothing warmth to joints and muscles. The inferior patella pad focuses compression on the patellar tendon to relieve discomfort, increase stabilization and improve patella tracking. Indicated for Chondromalacia Patella, fat pad impingement, Osgood-Schlatter Disease, patella instability, patella tendonitis, patellar subluxation, degeneration, and arthritis. Color: Black.

SAFE-T-SPORT*	SPORTS NEOPRENE	PATELLA SUPPORT
Catalog Number	Size	Fits
37-300400	Small	13" - 14"
37-300500	Medium	14.5" - 15.5"
37-300600	Large	16" - 17"
37-300700	X-Large	17.5" - 18.5"

Measure circumference around the center of the patella (kneecap).

Note - If measurement falls between two sizes, select the next larger size.

FEATURES AND BENEFITS
Sports Neoprene provides therapeutic warmth and support around the joint to speed healing.
<ul> <li>Inferior Patella Pad applies uniform compression to reduce the load on the patellar tendon to increase stability, improve patella tracking, and help alleviate pain and pressure.</li> </ul>
Adjustable hook and loop straps for an improved fit and enhanced compression control.

Simple wrap-around design for quick, easy application.



## ProLite® Compressive Knee Support with Viscoelastic Insert

For swollen or tender knees resulting from sprains, strains, sports injuries, arthritis or post-operative rehabilitation. Contains a donut-shaped viscoelastic insert that surrounds the patella and provides a gentle massaging effect as the knee is flexed and extended. For added side-to-side stability, the support includes flexible stays. Made of lightweight, breathable knit material. Will not retain excessive body heat. Color: White/Blue.

PROLITE® COMPRESSIVE KNEE SUPPORT			
Catalog Number	Size	Fits	
37-850400	Small	14" - 16"	
37-850500	Medium	16" - 18"	
37-850600	Large	18" - 20"	
37-850700	X-Large	20" - 22"	

Measure 4" above center of the kneecap.



## Safe-T-Sport® Standard Neoprene Knee Sleeve

SAFE-T-SPORT® STANDARD NEOPRENE KNEE SLEEVE			
Open Patella Catalog Number	Closed Patella Catalog Number	Size	Fits
37-373305	37-374304	X-Small	12" - 13"
37-373402	37-374401	Small	14" - 15"
37-373509	37-374508	Medium	16" - 17"
37-373606	37-374605	Large	18" - 19"
37-373703	37-374702	X-Large	20" - 21"
37-373804	37-374806	XX-Large	22" - 23"
37-373907	37-374903	XXX-Large	24" - 25"

Measure 4" above center of the kneecap.





## Safe-T-Sport® Standard Neoprene Knee

Support



SAFE-T-SPORT*	STANDARD NEOPR	ENE KNEE SUPPORT
Catalog Number	Size	Fits
37-201309	X-Small	12" - 13"
37-201406	Small	14" - 15"
37-201503	Medium	16" - 17"
37-201600	Large	18" - 19"
37-201707	X-Large	20" - 21"
37-201800	XX-Large	22" - 23"
37-201902	XXX-Large	24" - 25"

Measure 4" above center of the kneecap

### FEATURES AND BENEFITS

- Double layer front for added compression
- Longer length
- Sports neoprene for therapeutic warmth
- Soft nylon loop lining
- Open patella
- Navy

## EZ-ON<sup>®</sup> Neoprene Knee Wrap



EZ-ON* NEOPREN	E KNEE WRAP		
Catalog Number	Size	Fits	
37-307001	Universal	12" - 21"	
37-307700	Plus Size	20" - 30"	

Measure 4" above center of the kneecap.

- Wrap design fits a range of sizes and is easy to apply to swollen or injured knees
- Sports neoprene for therapeutic warmth
- Soft, terry material against skin
- Rubber patella support donut can be removed or cut to horseshoe
- Open behind knee
- Navy

### Neoprene Patella Stabilizer

#### with Removable Horseshoe

This stabilizer comes with a removable horseshoe that can be adjusted superiorly, inferiorly, medially or laterally depending on the stability needs of the knee. The horseshoe provides increased pressure around the patella for patellar tracking problems. Medial/lateral stability is provided by the spiral stays.

#### INDICATIONS:

- Chondromalacia Patella
- Infra and Subpatellar Bursitis
- Lateral Tracking Patella
- Osgood-Schlatter Disease
- Patellar Instability
- Patellar/Pain Syndrome
- Patellar Subluxation
- Patellar Tendonitis

NEOPRENE PATE	LLA STABILIZER W	ith REMOVABLE	HORSESHOE
Navy Catalog Number	Black Catalog Number	Size	Fits
37-103308	37-103340	X-Small	12" - 13"
37-103405	37-103440	Small	14" - 15"
37-103502	37-103540	Medium	16" - 17"
37-103609	37-103640	Large	18" - 19"
37-103706	37-103740	X-Large	20" - 21"
37-103801	37-103840	XX-Large	22" - 23"
37-103903	37-103940	XXX-Large	24" - 25"

Measure 4" above center of the kneecap.

- Four way stretch for proper fit and comfort
- Sports neoprene for therapeutic warmth
- Adjustable support straps
- Horseshoe support, spiral stays
- Easy slip-on style





### Lateral Knee Stabilizer



#### with "J" Shaped Buttress

The Safe-T-Sport® Lateral Knee Stabilizer with "J" Shaped Buttress provides excellent side-to-side support of the knee and stabilization of the kneecap. Indicated for conditions such as: Chondromalacia, pain from lateral pressure syndrome, patella instabilities, lateral patellofemoral malalignment, and lateral dislocation.

#### Sports Neoprene

Neoprene has soft, long loop nylon lining for comfortable wear. Durable enough for wear during sports. Sports neoprene provides consistent compression and therapeutic warmth to joints and muscles.

#### Tubular "J" Buttress

Maintains proper alignment of the patella (kneecap). Provides lateral side support and stabilization.

#### Adjustable Support Straps

Straps adjust medially for an adjustable amount of compression and increased stability. The unique angle of the straps creates extra support compression across the lateral "J" buttress for increased stabilization.

#### Flexible Stays

Spiral stays provide medial and lateral stability while maintaining full range of motion.

LATERAL KNEE STABILIZER with "J" SHAPED BUTTRESS			
Right Catalog Number	Left Catalog Number	Size	Fits
37-250310	37-250320	X-Small	12" - 13"
37-250410	37-250420	Small	14" - 15"
37-250510	37-250520	Medium	16" - 17"
37-250610	37-250620	Large	18" - 19"
37-250710	37-250720	X-Large	20" - 21"
37-250810	37-250820	XX-Large	22" - 23"
37-250910	37-250920	XXX-Large	24" - 25"

Measure 4" above center of the kneecap. Specify right or left.



### Soft Form Wrap-Around Stabilizing Knee Support

#### **INDICATIONS:**

- Slight to moderate ligament sprains and strains
- Tendinitis of the knee
- Arthritis
- Knee joint instability

SOFI FORM WRAP-AROUND STABILIZING KNEE SUPPORT			
Catalog Number	Size	Fits	
37-303300	X-Small	12" - 13"	
37-303400	Small	14" - 15"	
37-303500	Medium	16" - 17"	
37-303600	Large	18" - 19"	
37-303700	X-Large	20" - 21"	
37-303800	XX-Large	22" - 23"	
37-303900	XXX-Large	24" - 25"	

Measure circumference of leg, 4" above center of the kneecap.

- Unique Breathable Foam Material is latex free and neoprene free. Provides compression with a soft plush fabric against the skin. Non-irritating for those with sensitive skin.
- Wrap-Around Design is ideal for swollen or injured knees quick and easy to apply.
- Flexible Spiral Stays provide medial and lateral (side-to-side) stability of the knee joint.
- Open Popliteal eliminates bunching and irritation behind the knee.
- Condyle Pads offer additional padding and stabilization to the patella.
- Adjustable Support Straps for an improved fit that stays in place.
- Latex free and Neoprene free.





### OA / Arthritis Knee Brace



This lightweight, low profile brace is ideal for daily activities. The sports neoprene provides therapeutic warmth to the joint and surrounding tissues to help soothe aches and pains.

#### Indications:

- Mild Unicompartmental Osteoarthritis (OA)
- Knee Pain or Inflammation Due to Cartilage or Meniscus Deterioration
- Mild Ligament Instabilities
- Medial or Lateral Instability
- Post-Op (Meniscal or Articular Cartilage Repair)

OA / ARTHRITIS	KNEE BRACE		
Right Catalog Number	Left Catalog Number	Size	Fits
37-150310	37-150320	X-Small	12" - 13"
37-150410	37-150420	Small	14" - 15"
37-150510	37-150520	Medium	16" - 17"
37-150610	37-150620	Large	18" - 19"
37-150710	37-150720	X-Large	20" - 21"
37-150810	37-150820	XX-Large	22" - 23"
37-150910	37-150920	XXX-Large	24" - 25"

Specify Medial Right / Lateral Left or Medial Left / Lateral Right. \*Measure circumference of leg, 4" above center of the kneecap.

- Provides Excellent Support and Targeted Pain Relief
- May Relieve the Pain Associated with Mild Osteoarthritis (OA) of the Knee
- Unilateral Hinge and Strapping Configuration Apply Corrective Force to Offload the Affected Side of the Joint
- Wrap-Around, Patient-Friendly Design Easy to Apply to Swollen Tender Knees
- Neoprene Provides Therapeutic Warmth to Soothe the Aching Joint
- Support Straps Allow for Adjustable Compression
- Available in Medial and Lateral Designs for Left and Right Knee Applications
- Interchangeable Condyle Pads for Day-to-Day Adjustment



# Safe-T-Sport® Neoprene Hinged Knee Stabilizing Brace

Neoprene knee sleeve with durable steel hinged reinforcements allow for flexion and limit medial/lateral movement with a 180° hyperextension stop. Hinges in neoprene covered pockets for added comfort. Brace includes a superior horseshoe for added patella stability. Two adjustable straps prevent brace migration. Sports neoprene with nylon terry liner has a soft, seamless feel against the skin. Open patella and easy slip on style. Colors: Navy or Black.

SAFE-T-SPORT® NEOPRENE HINGED KNEE STABILIZING BRACE			
Navy Catalog Number	Black Catalog Number	Size	Fits
37-104307	37-104340	X-Small	12" - 13"
37-104404	37-104440	Small	14" - 15"
37-104501	37-104540	Medium	16" - 17"
37-104608	37-104640	Large	18" - 19"
37-104705	37-104740	X-Large	20" - 21"
37-104802	37-104840	XX-Large	22" - 23"
37-104906	37-104940	XXX-Large	24" - 25"



### Safe-T-Sport® Neoprene Hinged Knee Stabilizing Brace with Composite Hinges

Neoprene knee sleeve with lightweight composite hinges. Lightweight, thin neoprene with nylon terry liner for a soft, seamless feel against the skin. Ideal for use on swollen or tender knees resulting from sprains, strains, injuries, arthritis or post-operative rehabilitation. Color: Beige.

SAFE-T-SPORT" NEOPRENE HINGED KNEE STABILIZING BRACE with COMPOSITE HINGES			
Catalog Number	Size	Fits	
37-107330	X-Small	12" - 13"	
37-107430	Small	14" - 15"	
37-107530	Medium	16" - 17"	
37-107630	Large	18" - 19"	
37-107730	X-Large	20" - 21"	
37-107830	XX-Large	22" - 23"	

- Lightweight neoprene for consistent compression and therapeutic warmth
- Interior horseshoe support around knee joint for patellar tracking problems
- Lightweight composite hinges allow for flexion while limiting medial/lateral movement and preventing hyperextension
- Pull-on style with two loop lock adjustable straps





## Safe-T-Sport® Wrap-Around

**Hinged Knee Brace** 

Sports neoprene provides consistent compression and four-way stretch for proper fit and comfort. Neoprene retains the body's natural heat to provide therapeutic and soothing warmth to joints and muscles. Indicated for slight to moderate ligament sprains and tears, medial/lateral meniscal tears, arthritis and knee joint instability. Durable hinges give medial and lateral stabilization to the knee and allow for full knee flexion with hyperextension stops. Color: Black.

SAFE-T-SPORT*	WRAP-AROUND	HINGED KNEE BRACE
Catalog Number	Size	Fits
37-350300	X-Small	12" - 13"
37-350400	Small	14" - 15"
37-350500	Medium	16" - 17"
37-350600	Large	18" - 19"
37-350700	X-Large	20" - 21"
37-350800	XX-Large	22" - 23"
37-350900	XXX-Large	24" - 25"

FEATURES AND BENEFITS
Sports neoprene with four-way stretch compression
<ul> <li>Front wrap closures for quick, easy application - especially easy to apply over tender or swollen knees</li> </ul>
Medial/lateral metal hinges allow for full flexion and prevent hyperextension
Open behind the knee to help prevent bunching and improve comfort
Two adjustable support straps

Measure 4" above center of the kneecap

### FlexLite® Cool Walking Knee **Hinged Knee Support**

#### with Neoprene Knee Support

This practical hinged knee support is extremely lightweight and recommended for mild ligament injury, knee instability or degenerative joint disease. Ideal for walking or light activities. Provides medial and lateral stabilization. Hinges are made of an ultra light space age composite material that is strong and durable. Hinges include stops to help prevent the knee from hyperextending. Unlike neoprene, the base foam material is breathable and cool with less bulk. Wrap style front closures and tension control straps. Color: Black.

FLEXLITE" COOL WALKING KNEE HINGED KNEE SUPPORT					
Catalog Number	Size	Thigh	Calf		
37-108440	Small	14" - 16"	11.5" - 13"		
37-108540	Medium	16" - 18"	13" - 14.5"		
37-108640	Large	18" - 21"	14.5" - 16"		
37-108740	X-Large	21" - 24"	16" - 18"		
37-108840	XX-Large	24" - 27"	18" - 20"		



### PowerCentric® Composite, Polycentric

### **Hinged Knee Brace**

Offers intermediary support for mild to moderate ligament strains, sprains (ACL/PCL), meniscal injuries, medial/lateral co-ligament instabilities and tears, and hyperextension or chronic knee instability problems. The polycentric hinges allow for full range of motion and include a built-in "stop" to prevent hyperextension of the knee. Contoured condyle shells provide additional protection against lateral impacts during activity. Color: Black.

#### POWERCENTRIC\* COMPOSITE, POLYCENTRIC HINGED KNEE BRACE

Black Catalog Number	Size	Thigh	Calf
37-109400	Small	15" - 18"	11" - 13"
37-109500	Medium	18" - 21"	13" - 16"
37-109600	Large	21" - 25"	16" - 19"
37-109700	X-Large	25" - 29"	19" - 22"
37-109800	XX-Large	29" - 33"	22" - 24"

Measure your thigh 4 inches above the center of your kneecap. Measure around your calf 4 inches below the center of your kneecap. If measurements differ slightly, use the thigh measurement

- Polycentric hinges allow for full range of motion
- Lightweight and durable composite hinge material
- Hinge stop prevents knee hyperextension
- Adjustable straps
- Excellent patella stabilization
- Open popliteal design prevents bunching and irritation at the back of the knee
- Sports neoprene for consistent support and therapeutic warmth
- Easy slip-on style





### Adjustable ROM Knee Brace



#### **INDICATIONS:**


- Mild to moderate ligament strains and sprains (ACL/PCL/MCL/LCL)
- · Meniscal strains and tears
- Medial/lateral knee instability
- Chronic knee instability
- Post-op range of motion control

ADJUSTABLE ROM KNEE BRACE				
Size	Catalog Number			
Universal Regular	12" - 20"	37-450000		
Universal Plus	18" - 26"	37-450700		

Measure around center of kneecap. Brace fits right or left knee.

- Polycentric Range of Motion Hinges allow for progressive rehabilitation. Flexion and extension settings can be adjusted in 20 degree increments to control movement of the joint, preventing further injury and speeding healing. Dual axis, polycentric hinges track more closely with the knee motion and provide medial and lateral stabilization.
- Tri-Permalon™ Material is latex free and neoprene free and is ideal
  for those with sensitive skin or who are susceptible to heat build-up.
  Durable outer material is hook compatible and allows wrap closures
  to be placed anywhere on the support and gives excellent strength
  and support. Inner breathable foam material does not retain body
  heat. Allows for cool and comfortable wear. Smooth lycra lining is
  comfortable against the skin and provides power stretch for
  compressive support.
- Universal Two-Way Adjustable Straps fit a wide range of sizes, create a custom-like fit, and gives an improved fit that stays in place.
- Open Popliteal eliminates bunching and irritation behind the knee.
- Wrap-Around Design for easy application over swollen, tender joints.







### Specter® ROM Standard

Standard Specter® ROM provides post-operative immobilization and range of motion needed after ACL, PCL, MCL and LCL surgeries.

#### SPECTER® HINGE

All Specter® ROM Braces have a spring loaded snap-lock hinge to make ROM adjustments quick and easy. For patient lock-out, use the optional hex screw. Range of motion increments are in 7.5 degrees.

Use with or without optional Patient Lock-Out

Spring loaded snap-lock for easy setting



SPECTER® ROM STANDARD KNEE BRACE	
Size	Catalog Number
Universal	37-0170

#### FEATURES AND BENEFITS

- Strong durable hinge with snap-lock design and optional patient lock-out
- 7.5 degree incremental flexion/extension settings
- Integrated strap/strut design
- Conformable aluminum struts for varus/valgus position
- Breakdown conversion to shorter 4-strap design



## Specter® ROM Cool/Lightweight

Cool and lightweight variation of the Standard Specter® ROM. Ideal for warmer climates and long term wear.

SPECTER* ROM	COOL / LIGHTWEIGHT KNEE BRACE
Size	Catalog Number
Universal	37-0174



### Specter® ROM Short

Short version of Standard Specter® ROM. Ideal for post-operative patients during rehabilitation. Shorter length is more comfortable for ambulatory patients.

SPECTER* ROM SHORT KNEE BRACE	
Size	Catalog Number
Universal	37-0176



### HealWell® AFO Night Wrap



Helps prevent nighttime contracture with countervailing tension to hold the foot neutral. This helps relieve pain and discomfort and speed healing in new or mild conditions such as Achilles tendinitis, plantar fasciitis or arthritis. Designed to wear during rest or sleep.

#### Indications:

- Arthritis
- Achilles Tendinitis
- Plantar Fasciitis
- Calf Muscle Cramps
- Equinus (Muscle Tightening)
- Heel & Arch Pain
- Achilles Paratenonitis
- Post Surgical Treatment of the Foot
- Fits right or left

HEALWELL® AFC	NIGHT WRAP		
Catalog Number	Size	Fits Woman's Shoe	Fits Men's Shoe
58-6004	Small/Medium	6 - 9	up to 8
58-6006	Large/X-Large	9.5 and up	8.5 to 11.5

For shoe size greater than 11.5 it is recommended to use the HealWell\*  $\text{Cub}^{\text{TM}}$  or Standard Night Splint.

- SOFT LAMINATED FOAM Breathable, lightweight material is comfortable for night-time wear with a smooth lining.
- SLIM, LOW PROFILE No bulk, compact wrap design is easy to sleep in.
- NON-SLIP PAD For safe, occasional walking.
- COMPRESSION FORMED BACK Holds the wrap in place and conforms around the ankle and heel for improved foot positioning.
- DUAL TENSION STRAPS Countervailing tension straps help hold the foot in position and prevent plantar flexion. Straps criss-cross for less hulk
- PLANTAR SURFACE SUPPORT Tension straps work with a low profile, rigid stay on the plantar surface, functioning as a lever to maintain foot position.



### HealWell® Cub™ Plantar Fasciitis Night Splint Resting Comfort Slipper™

The HealWell® Cub™ Resting Comfort Slipper™ is designed to comfortably position the foot in the controlled amount of dorsiflexion to provide a gentle stretch across the plantar fascia, Achilles tendon and calf muscles. This gentle stretch helps reduce the muscle and ligamentous contracture, inflammation and associated pain - clinically proven to reduce symptoms of plantar fasciitis. This action also helps relieve many other conditions of the lower extremities including foot drop, Achilles tendinitis, and painful night calf muscle cramping. Color: Black. Specify right or left.

HEAL	HEALWELL® PLANTAR FASCIITIS NIGHT SPLINT					
Siz	:e	Right Catalog Number	Left Catalog Number	Fits Woman's Shoe	Fits Men's Shoe	
Sm	all	58-5004	58-5014	7.5 and under	6.5 and under	
Medi	ium	58-5005	58-5015	8 -11.5	7 -10.5	
Larç	ge	58-5006	58-5016	12 and up	11 -13	

#### FEATURES AND BENEFITS

- · Low profile, lightweight with less bulk for superior comfort.
- Open toe "slipper" is comfortable against the foot and easy to sleep in. Cool and comfortable to wear.
- Lateral side support with opening ensures no pressure on the sensitive heel, Achilles tendon, or bony ankle malleolus.
- Adjustable hinge can be positioned from 0 to 10 degrees depending on the amount of stretch required or desired.
- Built in toe wedge for more effective stretch across the plantar fascia.
- Ankle strap is perfectly positioned to hold the foot and heel in place for a more effective stretch.
- · Easy and quick to apply.
- · Cool, breathable foam material.



US Patents D497,207: D495,058 D481,798 & 7,182,743



#### ADJUSTABLE DORSIFLEXION Adjustable hinge can be positioned from 0 to 10 degrees, depending on the amount of stretch required or desired.

Preset at most common 5 degree setting.

### HealWell® Plantar Fasciitis

Night Splint



#### Our Standard Plantar Fasciitis Night Splint

Comfortably positions the foot in a controlled amount of dorsiflexion. Adjustable bilateral dorsiflexion straps enable the application of prolonged stretch to the soft tissues of the Plantar Fascia and the Achilles tendon. Both the heel orientation and tibial straps are easy to adjust and use with push lock buckles. Padded straps for increased comfort. Made of high density plastic with ventilated calf area for patient comfort and breathability. Removable, soft liner. Removable foam wedge to provide additional toe extension.

Universal design fits both left and right.

#### Indicated to Treat:

Plantar Fasciitis • Achilles Tendinitis • Heel Spurs Muscle Tightening • Runner's Cramp



HEALWELL® PLANTA	R FASCIITIS	NIGHT SPLINT	
Catalog Number	Size	Men's Shoe Size	Women's Shoe Size
58-409400	Small	Under 6	5 - 7
58-409500	Medium	6 - 8	7.5 to 9.5
58-409600	Large	8.5 - 10.5	7.5 - 9.5
58-409700	X-Large	11 - 13	Over 12
58-409800	XX-Large	Over 13	_



# HealWell<sup>®</sup> Multi AFO / Contracture Splint

The HealWell® Multi AFO is designed to help prevent and/or treat many conditions of the lower extremities. The shell is made of lightweight plastic and has a spring-like action. It has a soft and comfortable liner that allows air to circulate, absorbs perspiration and also keeps the foot warm.

Padded adjustable toe post and hip rotation bar. Liner is removable for

laundering. Available with or without transfer pad.

#### Indications:

- Non-Fixed Ankle Contractures
- Foot Drop
- Pressure Ulcers
- Hip Abduction
- Ankle Deformity
- Hip or Leg Rotation
- Post Surgery Knee or Hip

Available with or without transfer pad.

Pad assists the patient in transfers and in therapy treatment plans.



HEALWELL* M	HEALWELL* MULTI AFO/CONTRACTURE SPLINT						
Size	Men's Shoe Size	Women's Shoe Size	Calf Circumference	Catalog Number	Catalog Number		
Small	4 - 5	3 - 6	12" - 16"	58-100409	58-120400*		
Medium	6 - 9	7 to 10	14" - 18"	58-100506	58-120507*		
Large	10	11	16" - 20"	58-100603	58-120604*		
X-Large	Over 10	Over 11	18" - 22"	58-100700	58-120701*		

<sup>\*</sup> with Transfer Pad



### HealWell® C3™ Multi AFO/Contracture Splint

The HealWell® Tri-Laminate Multi AFO is similar to the standard Multi AFO, except with a Tri-Laminate material liner. The supple material allows for changes in the amount of swelling. Available with or without transfer pad.

Color: Black (hides soil).



HEALWELL® TRI-LAMINATE MULTI AFO/CONTRACTURE SPLINT						
Catalog Number	With Transfer Pad Catalog Number	Size	Men's Shoe	Women's Shoe	Calf Circum.	
58-220401	58-225400	Small	4 - 5	3 - 6	12" - 16"	
58-220503	58-225500	Medium	6 - 9	7 - 10	14" - 18"	
58-220604	58-225600	Large	10	11	16" - 20"	
E0 22070E	E0 22E700	V Lorgo	Over 10	Over 11	10" 22"	

- Completely suspends the heel, protecting it from pressure or friction
- Liner is made of Tri-Laminate foam material that is cool and comfortable to wear
- Tri-Laminate material has soft liner to help wick moisture into the foam inner layer
- Supple Tri-Laminate material allows for changes in the amount of swelling
- Padded toe post is adjustable, protecting the toes from contact with bed linens
- Hip rotation control bar helps position the foot correctly by preventing internal or external rotation of the hip, foot and/or ankle
- Ankle strap reinforces foot position
- Available with or without transfer pad



### HealWell<sup>®</sup> Soft Ease<sup>™</sup> Multi-AFO / Heel Suspender

#### with Soft Replaceable Foam Cover

The HealWell® Soft Ease™ Multi-AFO/Heel Suspender is designed to position the foot with the heel completely suspended, eliminating pressure and preventing ulcers. The unique shape allows the heel to float in an air cavity with absolutely no contact. This product is an economical, effective solution for heel ulcer prevention for any at-risk patient.

By suspending the heel, the Soft Ease promotes air circulation around the foot, aiding in treatment and healing of existing pressure ulcers. The resilient plastic shell base correctly positions the foot for treatment of foot drop and non-fixed ankle/foot contractures. It also maintains alignment of the foot with the lower leg, minimizing inversion and eversion of the foot.

HEALWELL® SOFT EASE™ MULTI-AFO / HEEL SUSPENDER					
Catalog Number	Size	Men's Shoe	Women's Shoe		
58-240400	Small / Medium	4 - 9	5 - 10		
58-240600	Large / X-Large	10 - 13	11+		

Measure according to U.S. shoe size. Fits Right or Left

FEATURES AND BENEFITS		
<ul> <li>Super soft material wicks moisture away to keep skin dry and reduces friction and shear</li> </ul>		
<ul> <li>Integrated padded toe post protects the toes from contact with bed linens to prevent skin breakdown</li> </ul>		
Fully adjustable straps accommodate changes in edema		
Extra padding at sensitive Achilles area		

#### INDICATIONS:

- Heel Ulcer Treatment and Prevention
- Diabetic Foot Care Management
- Foot Drop
- Post Surgery
- Non-Fixed Ankle Contractures

# Ankle Foot Orthosis/Foot Drop Splint

Indicated to correct foot drop during ambulation in patients with neurological deficit. Made of durable and resilient polypropylene. Fits into standard shoe. Can be easily trimmed or molded with heat. Forms at 310°-325°F. Adjustable calf band with padding and adjustable closure. Specify right or left.

ANKLE FOOT ORTHOSIS/FOOT DROP SPLINT					
Right Catalog Number	Left Catalog Number	Size	Men's Shoe	Women's Shoe	
58-320412	58-320421	Small	3 - 6	5 - 8	
58-320519	58-320528	Medium	6 - 9	8 - 11	
58-320616	58-320625	Large	9 - 11+	-	





### **Pullover Elastic Anklet**





PULLOVER ELASTIC ANKLET				
Catalog Number	Size	Fits		
40-101401	Small	8" - 9.5"		
40-101508	Medium	9.5" - 11"		
40-101605	Large	11" - 12.5"		
40-101702	X-Large	12.5" - 14.5"		

Measure around ankle.

FEATURES AND BENEFITS	
Lightweight	
Knit elastic	
Seamless, slip-on style	
White	

### **Nylon Anklet**





NYLON ANKLET		
Catalog Number	Size	Fits
40-401408	Small	7" - 8.5"
40-401505	Medium	8.5" - 10.5"
40-401602	Large	10" - 11"
40-401709	X-Large	11.5" - 12.5"

Measure around ankle.

FEATURES AND BENEFITS				
oe				
control swelling and inflammation				

Safe-T-Sport® Neoprene
Ankle Support



SAFE-T-SPORT® NEOPRENE ANKLE SUPPORT					
Size	Fits				
Pediatric	6.5" - 7.25"				
Youth	7.25" - 8.5"				
Small	8" - 9.5"				
Medium	9.5" - 11"				
Large	11" - 12.5"				
X-Large	12.5" - 14.5"				
	Size Pediatric Youth Small Medium Large				

Measure around ankles


- Neoprene provides therapeutic warmth and compression to ankle joint
- Open heel, slip-on style
- Navy



# EZ-ON® Wrap-Around Ankle Support

Double elastic strap offers adjustable compression to both the medial and lateral sides of the ankle and arch of the foot. The wide strap simulates taping and wraps in a figure-eight design with a simple hook closure. Indicated for swollen, tender, or weak ankles, sprains, strains, sports injuries, and arthritis. Color: Black.

EZ-ON° WRAP-AROUN	ND ANKLE SUPPOR	LE SUPPORT		
Catalog Number	Size	Fits		
40-550400	Small	8" - 10.5"		
40-550500	Medium	10.5" - 12.5"		
40-550600	Large	12.5" - 14.5"		
40-550700	X-Large	14.5" - 16.5"		
40-550800	XX-Large	16.5" - 18.5"		
40-550900	XXX-Large	18.5" - 20.5"		





- SUPPORTS AND STABILIZES the ankle while allowing for full flexibility.
- MADE WITH SPORTS NEOPRENE AND TRI-PERMALON™
   The foot and ankle portion is made of neoprene which provides therapeutic warmth and compression to the joint and soft tissues.
   The top portion is made of a breathable Tri-Permalon™ foam to let the heat dissipate and air circulate.
- WIDE, DOUBLE ELASTIC STRAP wraps around the joint for increased stability. Figure-eight strap design simulates taping of the ankle and provides arch support.
- SIMPLE HOOK CLOSURE allows for quick and easy application.
- OPEN HEEL comfortable wrap support.
- MAY BE WORN UNDER OR OVER SOCKS.

### Ankle Guard with Optional Plastic Inserts

Soft, flexible vinyl ankle brace provides medial/lateral support to weak or injured ankles. Lace-up front allows for an intimate fit. Light enough to wear during athletic activity to help reduce the chance of sprains and strains. Medial/lateral support can be increased with optional plastic inserts. Plastic inserts are included. Fits left or right. Color: White.

ANKLE GUARD with OPTIONAL PLASTIC INSERTS				
Catalog Number	Size	Men's Shoe Size	Women's Shoe Size	
40-411201	XX-Small	Under 6	Under 7	
40-411308	X-Small	6 - 7	7 - 8	
40-411405	Small	8 - 9	9 - 10	
40-411502	Medium	9 - 11	10 - 12	
40-411609	Large	11 - 13	12 - 14	
40-411706	X-Large	14 & up	-	

Size according to shoe size.



With Optional Plastic Support Inserts



# Canvas Lace-Up ■ Ankle Brace



LITE

Indicated when support of ankle joint is needed. Ideal for after cast removal or for mild/moderate sprains. Made of a breathable canvas material with soft flannel lining. Has metal spring stays for immobilization and added stability. Lace-up closure allows for an intimate medial/lateral fit. Fits left or right. Color: White.

CANVAS LACE-UP A	NKLE BRACE		
Catalog Number	Size	Men's Shoe Size	Women's Shoe Size
40-311202	X-Small	6 - 7	Under 7-8
40-311406	Medium	8 - 9	9 - 10
40-311503	Medium	9 - 11	10 - 12
40-311600	Large	11 - 13	12 - 14
40-311707	X-Large	13 - 14	14 - 15

Size according to shoe size.

# ProLite® Ankle Support Knitted Pullover

Provides excellent compression and stabilization for weak or injured ankles. Made of a stretch knitted elastic that is a lightweight and breathable alternative to neoprene. Elastic provides compressive support while maintaining flexibility.

Colors: Blue or Beige.



PROLITE* ANKL	E SUPPORT KNITTED	PULLOVER	
Blue Catalog Numbe	Beige er Catalog Number	Size	Fits Ankles
40-400400	40-400430	Small	7" - 8.5"
40-400500	40-400530	Medium	8.5" - 9.5"
40-400600	40-400630	Large	9.5" - 11"
40-400700	40-400730	X-Large	11" - 12"

Measure around ankle just above the ankle bone.

- Comfortable, firm support for weak or injured ankles
- Made of a lightweight material for a breathable alternative to neoprene
- Unique knit design for precise compression
- Available in blue or beige



### F8 Ankle Brace

The figure eight straps of the F8 Ankle Brace eliminates the need for taping and the lightweight, ballistic nylon provides a comfortable fit without bulk.

F8 ANKLE BRACE (MALE)				
Size	Shoe Size	Catalog Number		
X Small	3 - 5	2231		
Small	6 - 7	Q2232		
Medium	8 - 10	2233		
Large	11 - 12	2234		
X Large	13 - 14	2235		

F8 ANKLE BRACE (FEMALE)			
Size	Shoe Size	Catalog Number	
X Small	4 - 6	2231	
Small	7 - 9	Q2232	
Medium	10 - 11	2233	
Large	12 - 13	2234	
X Large	14 +	2235	



### InnerLOK 8™ Ankle Brace

Ideal for sports and to support weak or injured ankles. Made of lightweight, low-profile nylon material for long-lasting wear. Extra ankle stabilization with internal figure eight straps that are preconfigured to ensure correct placement every time. Application is quick and easy. Comfortably fits in athletic shoes. Elastic back. Fits left or right. Color: Black.

INNERLOK 8™ ANKLE BRACE				
Catalog Number	Size	Men's Shoe Size	Women's Shoe Size	
40-511400	Small	6 - 7	7 - 9	
40-511500	Medium	8 - 10	10 - 11	
40-511600	Large	11 - 12	12 - 13	
40-511700	X-Large	13 - 14	14 & up	

Size according to shoe size.



Figure-Eight Strapping



### Allsport Ankle Brace



ALLSPORT ANKLE BRACE	
Size	Catalog Number
Universal Left	144ALGV
Universal Right	144ARGV

# Ankle Stirrup Brace with GelFoam™

Gel filled liner conforms to the malleoli of the ankle for intimate support and a custom fit. GelFoam™ liner can be removed and placed in the refrigerator for cold therapy. Comfortable, soft terry lining. Shell design fits left or right ankle. Color: Navy.



- Gel/Foam bladder with soft terry liner
- Liners can be placed in the refrigerator for cold therapy
- Adjustable straps



# ProLite® Ankle Stirrup Brace with Air Liners

This brace provides medial and lateral support to stabilize the ankle joint for faster recovery and healing of tender soft tissues. Ideal brace for mild or moderate sprains, strains, or after cast removal. Provides stabilization to chronically weak ankle joints as well.

Color: White/Blue.





### Actimove® Talocast Functional Ankle Brace

The Actimove Talocast Ankle Brace is anatomically shaped to provide an ideal fit. The low profile design allows an easy fit within an athletic shoe. Features a foam, interior pad which provides patient comfort.

ACTIMOVE® TALOCAST FUNCTIONAL ANKLE BRACE			
Size	Quantity	Catalog Number	
Left - Standard	Single	72819-00083-00	
Right - Standard	Single	72819-00084-00	
Left - Junior	Single	73088-00001-00	
Right - Junior	Single	73088-00000-00	

- Strong and lightweight
- Medial and lateral anatomically shaped shells and large malleolar indent
- Straps with D-rings allows close, precise placement
- Foam interior padding
- Additional protection padding





U.S. Patent 5,014,691

# Sheartex® Ankle Stirrup Brace

Sheartex\* padding is extra soft for greater patient comfort. Liner is completely air permeable to keep heat and perspiration to a minimum - can easily be removed for laundering. Shell design fits left or right ankle. Color: Navy.

Size
Universal



### Active Ankle® Trainer

The Active Ankle\* Trainer's U-shaped stirrup relieves pressure from the ankle joint and offers superior inversion/eversion protection. Its bi-lateral hinge allows for freedom of motion in dorsiflexion and plantarflexion. Neoprene padding and dual strapping provide comfort and a snug fit and will work on either the left or right ankle.

ACTIVE ANKLE® TRAINER				
Size	Male Shoe Size	Female Shoe Size	Catalog Number	
X Small	n/a	5 - 6	AAT-XS	
Small	5 - 8.5	6.5 - 9.5	AAT-S	
Medium	9 - 12	10 - 13	AAT-M	
Large	12.5 - 16	13.5 - 16	AAT-L	



FlexLite<sup>®</sup> Sport Articulating Hinged Ankle Brace

Provides ankle stability to weak or injured ankles. Ideal for treatment of slight to moderate ankle sprains, strains, and ankle instability or weakness. Makes a great preventative brace during activity to help reduce the chance of sprains or strains. Color: Black.

FLEXLITE SPORT A	RTICULATING H	INGED ANKLE BR	ACE
Catalog Number	Size	Men's Shoe Size	Women's Shoe Size
40-800400	Small	5 - 8.5	6.5 - 9.5
40-800500	Medium	9 - 12	10 - 13
40-800600	Large	12.5 - 16	13.5 +

Size according to shoe size.

- Bi-lateral articulating hinges glide smoothly allowing for full, unrestricted foot flexion.
- Foot plate transfers weight-bearing forces from the ankle joint and increases foot and ankle stability.
- Solid u-shape provides superior inversion/eversion protection.
- Lightweight support shells molded from a lightweight plastic material in a design tailored to the ankle for maximum stabilization of the joint.
- Padded, formed foam liners. Ultra-lightweight liners cushion and hug the ankle with less bulk for a sleek, low profile fit.
- Integrated elastic straps. Two closure straps are integrated for a secure fit with no movement.
- Fits Right or Left



FULL, UNRESTRICTED FOOT FLEXION



# StepLite® Easy Strider™ Walker Braces

#### Indications:

- Acute or severe ankle sprains
- Lower extremity soft tissue injuries
- Stable fractures of the foot or ankle
- After cast removal

- · Severe ankle or foot strains
- Stress fractures of the foot or ankle
- Post achilles tendon repairs

FEATURES AND BENEFITS
Lightweight and low profile
Rocker bottom with no skid sole
Criss-cross ankle straps for improved foot stabilization
Adjustable top straps
Breathable foam liner with soft wicking inner material to keep skin cool and dry
Extra comfort padding
Compact packaging for easy handling and storage
High and low height walkers available

STEPLITE* EASY	STRIDER™ WALKER	BRACES	- LOW & HIGH	HEIGHT
Low Height Catalog Number	High Height Catalog Number	Size	Fits Men's Shoe	Fits Women's Shoe
43-420400	43-430400	Small	6 - 7.5	5 - 8.5
43-420500	43-430500	Medium	8 - 9.5	9 - 10.5
43-420600	43-430600	Large	10 - 12.5	11 - 12.5

StepLite<sup>®</sup> Easy Air<sup>™</sup> Pneumatic Ankle Walker Brace

The built-in air bladder contains air chambers that can be inflated or deflated for a customized fit or to accommodates changes in swelling. The bladder and pump are built right into the liner for easy adjustments on the go -- no pump attachments required! The inflated air bladder gives a massaging effect as patient walks. The soft, breathable liner wicks moisture away from the skin. Also features a non-slip rocker bottom. Available in low or high height. Fits right or left. Color: Black/Gray.

STEPLITE" EASY	AIR™ PNEUMATIC A	NKLE WALKER	BRACES - LOW	& HIGH HEIGHT
Low Height Catalog Number	High Height Catalog Number	Size	Fits Men's Shoe	Fits Women's Shoe
43-440300	43-450300	X-Small	2 - 4	3.5 - 5.5
43-440400	43-450400	Small	6 - 7.5	5 - 8.5
43-440500	43-450500	Medium	8 - 9.5	9 - 10.5
43-440600	43-450600	Large	10 - 12.5	11 - 12.5
43-440700	43-450700	X-Large	12.5+	14+





### Actimove® Lower Leg Walker

The Actimove Lower Leg Walker combines comfort and support at an affordable price. Its low profile, contoured construction ensures anatomical fit. The rocker bottom is designed to increase stability of the patient while allowing him to resume natural gait.

ACTIMOVE*	ACTIMOVE® LOWER LEG WALKER				
Size	Shoe Si Male	<sup>ze</sup> Female	Catalog Number		
X Small	2 - 4	3.5 - 5.5	72819-00078-00		
Small	4.5 - 7	6 - 8.5	72819-00079-00		
Medium	7.5 - 10.5	9 - 12	72819-00080-00		
Large	10.5 - 12.5	12 - 14	72819-00081-00		
X Large	12.5 +	14 +	72819-00082-00		

FEATURES AND BENEFITS	
Five sizes	
Quick and easy application	
Vented liner	
Superior comfort	
Replacement liners also available	



### Adjustable walker

The adjustable, low-profile walker is designed to promote a normal gait while offering 45 degrees of ROM adjustment. It has a three-way compression system with telescoping medial-lateral stays and extra Achilles heel pads, providing a custom fit, enhanced stability and support. Its rigid, slip resistant sole, washable tricot foam liner and fully cushioned insole provide added comfort.


ADJUSTABLE W	VALKER			
Dimensions	Sho Male	e Size Female	Size	Catalog Number
3.75" x 8.5"	1 - 4	1 - 5	X Small	SWXS
3.75" x 9.75"	4 - 8	5.5 - 8.5	Small	SWS-1
4.25" x 10.75"	7 - 9	8.5 - 10.5	Medium	SWM-1
4.5" x 11.5"	9 - 13	10.5 - up	Large	SWL-1

- · Low profile design promotes normal gait
- Three-way compression system telescoping medial-lateral and posterior Achilles pad provide custom fit, enhanced stability, support and comfort
- Adjustable ankle joints allow up to 45 degrees total range of motion or fixed setting in dorsiflexion, plantarflexion and neutral positions
- Rigid sole features a computer-aided design for natural walking gait and a slip-resistant tread for safety
- Washable, tricot foam liner and a fully cushioned insole for patient comfort
- Additional pads included to help contour the ankle for better control
- Easy step-in/step-out application saves time



### Specialist® Post-Op Shoe

The Specialist Post-Op Shoe consists of a large wrap with hook and loop closures and an open toe to accommodate post-surgical dressings. A semi-rigid flared sole in the heel provides enhanced medial lateral stability. The padded collar and reinforced heel provide extra patient comfort.



SPECIALIST® POST-OP SHOE (FEMALE)					
Description	Color	Catalog Number	Quantity		
Female Small	Black	9310	6 / case		
Female Medium	Black	9311	6 / case		
Female Large	Black	9312	6 / case		

### Deluxe Flexible Sole Healing Shoe

The insole construction of dense foam in the Deluxe Flexible Sole Healing Shoe provides exceptional comfort while helping to absorb the shock of heel strike.

Width

3.75"

DELUXE, FLEXIBLE SOLE HEALING SHOE (MALE)



Padded heel counter provides comfort and prevents heel slippage

Adjustable and placeable hook and loop closure provides even

	12.5"	4.5"	X Large	800
	DELUXE,	FLEXIBLE SOLE HI	EALING SHOE (F	EMALE)
ENEFITS	Length	Width	Size	Cat
Sole	8.25"	3"	X Small	800

Length

11"

12"

11.5"



Size

Small

Large

Medium

**Catalog Number** 

800-MS

800-MM

800-ML

800-MXL



FEATURES AND B

Flexible, non-skid

compression and secure fit

### C3<sup>™</sup> Post Op Shoe

Designed for wear after soft tissue procedures or post-operative care. Especially useful for geriatric and diabetic patients with insensitive feet. C3<sup>TM</sup> Tri-Laminate upper makes this product cool and comfortable to wear. Adjustable straps accommodate the bulkiest dressings. Removable trim-fit tongue provides a better fit. Semi-rigid molded non-skid sole absorbs shock during ambulation. Fits right or left. Sold as each. Color: Black.



C3™ POST OP SHO	E			
Men's Catalog Number	Women's Catalog Number	Size	Men's Shoe	Women's Shoe
-	-	X-Small	-	Under 4
43-821400	43-820400	Small	6 - 8	4 - 6
43-821500	43-820500	Medium	8.5 - 10	6.5 - 8
43-821600	43-820600	Large	10.5 - 12	8.5 - 10
43-821700	-	X-Large	12.5 - 14	-

### Super Post Op Shoe

Designed for wear after soft tissue procedures. Especially useful for geriatric and diabetic patients with insensitive feet. Durable open weave nylon upper with soft inner lining. Molded heel provides a smooth and stable contour. Removable tongue can be trimmed for proper fit. Semi-rigid molded non-skid sole absorbs shock during ambulation. Fits right or left. Sold as each. Color: Navy.



SUPER POST OP SI	1 O E			
Men's Catalog Number	Women's Catalog Number	Size	Men's Shoe	Women's Shoe
-	43-890308	X-Small	-	Under 4
43-891404	43-890405	Small	6 - 8	4 - 6
43-891501	43-890502	Medium	8.5 - 10	6.5 - 8
43-891608	43-890609	Large	10.5 - 12	8.5 - 10
43-891705	-	X-Large	12.5 - 14	-



### Deluxe Rigid Post-Op with Collar Shoe

The Deluxe Rigid Post-Op Shoe is constructed with a rigid sole for total stabilization, a heel counter for support and a padded collar to reduce heel slippage.



DELUXE RIC	SID AND POST-O	P WITH COLLAR	SHOE (MALE)
Length	Width	Size	Catalog Number
11"	3.5"	Small	445-MS
11.5"	4"	Medium	445-MM
12"	4"	Large	445-ML

DELUXE RIC	SID AND POST-O	P WITH COLLAR :	SHOE (FEMALE)
Length	Width	Size	Catalog Number
9.5"	3.25"	Small	445-WS
10"	3.5"	Medium	445-WM
10.5"	3.5"	Large	445-WL

#### EATURES AND BENEFITS

- Durable, padded, breathable upper with soft, white inner lining and padded tongue
- Loop-lock / D-ring for best fit

### Vinyl Molded Sole Rocker Cast Shoe

The molded rocker sole of the Vinyl Molded Rocker Cast Shoe provides a more natural gait pattern and the rigid insole design helps prevent slippage inside the cast shoe. The closed toe and water resistant vinyl construction help protect the cast against the elements. An open toe version is also available.



VINYL MOLE	DED ROCKER CA	ST SHOE	
Length	Width	Size	Catalog Number
7.25"	3.5"	X Small	VRCT-XS
9.5"	4.25"	Small	VRCT-S
11.25"	5"	Medium	VRCT-M
12.75"	5.5"	Large	VRCT-L
13.5"	5.5"	X Large	VRCT-XL

- Keeps foot warmer in cool weather
- · Hook and loop closures for a more secure fit

VINYL MOLD	ED ROCKER CA	ST OPEN TOE CA	ST SHOE
Length	Width	Size	Catalog Number
7.25"	3.5"	X Small	VROT-XS
9.5"	4.25"	Small	VROT-S
11.25"	5"	Medium	VROT-M
12.75"	5.5"	Large	VROT-L
13.5"	5.5"	X Large	VROT-XL



### Canvas Molded Sole Rocker

The durable canvas duck design of the Canvas Molded Sole Rocker allows sizing versatility with its open toe and heel design. Hook and loop closures secure the cast shoe against slippage.

CANVAS MOLDED SOLE ROCKER					
Length	Width	Size	Catalog Number		
7.25"	3.5"	X Small	CMR-XS		
9.5"	4"	Small	CMR-S		
11"	4.75"	Medium	CMR-M		
12.5"	5.5"	Large	CMR-L		
13.5"	5.5"	X Large	CMR-XL		



The Open-Toe Cast Shoe is constructed with a heavy canvas upper for increased strength and a tri-level EVA rocker bottom that provides shock absorption and smooth, roll-through gait for greater walking stability. Affordably priced, this premium cast shoe can be worn on either the left or right foot.

OPEN-TOE CAST SHOE					
Dimensions	Color	Catalog Number	Quantity		
6.5" - Pediatric	Blue	58-35007	1 / each		
7.5" - X Small	Blue	58-35008	1 / each		
8.5" - Small Adult	Blue	58-35010	1 / each		

OPEN-TOE CAST SHOE					
Dimensions	Color	Catalog Number	Quantity		
10" - Medium Adult	Blue	58-35011	1 / each		
11.5" - Large Adult	Blue	58-35012	1 / each		
12.5" - X Large Adult	Blue	58-35013	1 / each		






### Canvas Rocker Bottom Cast Shoe

Provides protection for walking casts. Made of a durable canvas material with two adjustable closures. Has open toe and open heel with non-skid sole. Color: Navy (infant size/pediatric print). Fits right or left





Sole Dimensions

CANVAS ROCKER BOTTOM CAST SHOE						
Catalog Number	Size	Α	В			
43-101107	Infant	5.5"	3"			
43-101204	Pediatric	7.25"	3.625"			
43-101408	Small	9.375"	4.25"			
43-101505	Medium	11"	4.875"			
43-101602	Large	12.625"	5.375"			
43-101709	X-Large	13.5"	5.625"			

### Walking Heels

**Pedestal** 



The deep tread of our Walking Heels provides good patient comfort and stability by allowing regular heel-toe walking. The specialized synthetic rubber is resilient and durable while the center slot permits easy and secure application.


WALKING HEELS					
Description	Catalog Number	Quantity			
Pedestal Walking Heel	4183-133	12 / box			
Cushion Tread Walking Heel	4183-135	12 / box			
Adult Walking Heel	4183-130	12 / box			
Child Walking Heel	4183-132	12 / box			



### Orthopedic Supports for Pediatric / Youth

#### **Anatomically Designed...**

A full range of orthopedic supports and braces specifically designed for pediatric and youth sizes. These pediatric products have been anatomically designed using anthropometric data compiled from two studies to accommodate 95% of children, providing a more precise fit and better overall support.



#### ...Specifically for Kids

Anatomically designed using anthropometric data compiled from studies performed by the Consumer Product Safety Commission and the Information Technology Laboratory at the National Institute of Standards and Technology. These studies collected and recorded average body measurements of each age range to aid in the design of products specifically for children. Using data from the resulting report, Physical Characteristics of Children\*, our designed its pediatric products have been designed to accommodate 95% of children.

Examples of Sizing Data

#### **Quality Materials**

The high quality, soft materials used in our pediatric line provide optimal comfort, support, and durability that both kids and parents can count on. All infant and pediatric sized products are child friendly – latex free and free of metal buckles and harsh closures, reducing the chance of irritation.

#### **Approximate Age Ranges**

- Infant Fits 0-24 months
- Pediatric Fits 2-7 years
- Youth Fits 8-12 years
- Small Adult Fits 13+ years
- \* Age ranges for each size are approximate and should be used as a general guide. Measurements should be taken to ensure proper fit and function.

#### Informative Color Packaging & Instruction Sheet

The infant and pediatric sized products in this offering are supplied in consumer friendly packaging with complete care and use instructions on the carton or on the enclosed insert sheet. Each box is labeled with measuring information and a provider peel-off label for easy transfer, billing and reordering. Many of the x-small and small sizes in our regular line fit youth and teens as well.



tems indicated with this icon are only available in our Pediatric packaging.







insert

### Neoprene Knee Sleeve



FEATURES AND BENEFITS
Sports neoprene provides therapeutic warmth and compression
Open patella, easy slip-on style
Supports the knee
. Limite avalling and degrapes noin

#### Indications:

- Joint Effusion (swelling)
- Knee Contusion (bruising)
- Knee Strains/Sprains
- Patella Instability
- Muscle Strains
- Juvenile Arthritis

NEOPRENE KNEE SLEEVE - PEDIATRIC LINE					
Catalog Number	Size	Approx. Age	Circumference at Kneecap		
37-373225	Ped / Sm	2 - 3 yrs	6" - 8"		
37-373250	Ped / Md	4 - 5 yrs	8" - 10"		
37-373275	Ped / Lg	6 - 7 yrs	10" - 12"		

NEOPRENE KNEE SLEEVE - REGULAR LINE					
Catalog Number	Size	Approx. Age	Circumference 4" Above Knee*		
37-373305	Youth / XS	8 - 12 yrs	12" - 13"		
37-373402	Small	13+ yrs	14" - 15"		

<sup>\*</sup>Measure circumference 4" above center of knee.

### Neoprene Patella Stabilizer



#### Indications:

- Chondromalacia Patella
- Lateral Tracking Patella
- Osgood-Schlatter Disease
- Knee Sprains
- Patella Instability
- Patella Subluxation

NEOPREN	E PATELLA	STABILIZER	- PEDIATRIC	LINE
Catalog Number		Size	Approx. Age	Circumference at Kneecap
37-10320	7 1	Ped / Lg	6 - 7 yrs	10" - 12"

FEATURES AND BENEFITS					
Sports neoprene provides therapeutic warmth and compression	n				
Protects and supports the knee					
Removable horseshoe can be placed superior, inferior, lateral					

- or medial, depending on where support is needed

  Two straps provide adjustable compression and prevent migration
- Open patella, easy slip-on style

NEOPRENE F	PATELLA STABILIZER	R - REGULAR	LINE
Catalog Number	Size	Approx. Age	Circumference 4" Above Knee*
37-103308	Youth / XS	8 - 12 yrs	12" - 13"
37-103405	Small	13+ yrs	14" - 15"

<sup>\*</sup>Measure circumference 4" above center of knee.

### Knee Immobilizer

#### Indications:

- Post Surgery
- After Cast Removal
- Knee Stability After Injury

	KNEE IMMOBILIZER - REGULAR LINE						
С	atalog Number Each	Size	Approx. Length*	Fits**			
	38-140002	Youth / XS	12"	12" - 13"			
	38-141004	Youth / Lg	14"	14" - 15"			

<sup>\*</sup>Measure inseam to determine appropriate length.

#### FEATURES AND BENEFITS

- Popular tri-panel design with adjustable side panels
- Made of a soft foam treated with Antimicrobial Protection for a fresher, more hygienic product
- Lateral and medial stays increase stability
- Contoured posterior stays immobilize the knee
- Easy loop lock closures



### Wrist Splint

#### Indications:

- Sprains
- Strains
- After Cast Removal

WRIST SPLINT	- PEDIATRIC LINE				
Catalog Number Right	Catalog Number Left	Size	Approx. Age	Wrist Circumference	
22-300110	22-300120	Infant	0 -24 mos	4.5" and under	
22-300210	22-300220	Pediatric	2 - 7 yrs	4.5" - 5.5"	
22-300310	22-300320	Youth	8 - 12 yrs	5.5" - 6.5"	

Measure circumference at wrist. Specify right or left.

- Soft, breathable foam for comfort
- Antimicrobial Protection for a fresher, more hygienic product
- Palmar/dorsal stays immobilize and support the wrist, prevent extension and flexion
- Low contour at palmar crease allows for full finger dexterity
- Simple hook and loop closures



<sup>\*\*</sup>Measure circumference 4" above center of knee.

### ProLite® Wrist Splint



#### Indications:

- Sprains
- Strains
- After Cast Removal

PROLITE" WRIST SPLINT - REGULAR LINE						
Catalog Number Right	Catalog Number Left	Size	Approx. Age	Wrist Circumference		
22-450412	22-450422	Youth / Sm	8 - 12 yrs	5.5" - 6.5"		
22-450513	22-450523	Sm Adult / Med	13+ yrs	6.5" - 7.5"		

Specify right or left

#### FEATURES AND BENEFITS

- · Cool, breathable open cell foam for comfort
- Provides continuous compression
- Low contour at palmar crease allows for full finger dexterity
- Removable palmar stay and medial/lateral stays for maximum support and immobilization
- Easy hook and loop closures

### ProLite<sup>®</sup> Wrist Splint





- Post Surgery of Wrist/Thumb
- Severe Wrist Sprains
- Strains/Sprains of the MCP Joint
- After Cast Removal



PROLITE" WRIST SPLINT with ABDUCTED THUMB - REGULAR LINE						
Catalog Number Right	Catalog Number Left	Size	Approx. Age	Wrist Circumference		
22-460410	22-460420	Youth / Sm	8 - 12 yrs	5.5" - 6.5"		
22-460510	22-460520	Sm Adult / Med	13+ yrs	6.5" - 7.5"		

Specify right or left

- Semi-rigid thumb spica, contoured for immobilization and thumb abduction
- Protects thumb from impact and shock
- Palmar and dorsal stays prevent flexion and extension
- Cool, breathable open cell foam for comfort
- Easy hook and loop closures

### Finger / Thumb Splint

ı	FINGER / THUMB	SPLINT - REGUI	LAR LINE		
	Catalog Number Box of 6	Size	Approx. Age	Splint Length	
	25-1022PK	Ped	2 - 7 yrs	5"	
	25-1024PK	Youth / Sm	8 - 12 yrs	7"	
	25-1025PK	Sm Adult / Med	13+ yrs	8"	

FEATURES AND BENEFITS
Permits any degree of flexion/hyperextension
Malleable stay can be contoured to finger/thumb
Soft foam liner for increased comfort
Movable hook and loop closures for greater adjustability
Individually packaged in a convenient 6-pack box



### Neoprene Elbow Sleeve

#### Indications:

- Tendonitis
- Tennis Elbow
- Juvenile Arthritis
- Muscle Strains

NEOPRENE ELBO	OW SLEEVE - PE	DIATRIC LINE	
Catalog Number Each	Size	Approx. Age	Forearm Circumference
19-601304	Ped / Lg	5 - 7 yrs	6.5 - 7.5"

NEOPRENE ELBOW SLEEVE - REGULAR LINE						
Catalog Number Each	Size	Approx. Age	Forearm Circumference			
19-601406	Youth / Sm	6 - 7 yrs	7" - 9"			
19-601503	Sm Adult / Med	13+ yrs	9" - 10.5"			



- Sports neoprene provides therapeutic warmth and compression around entire joint
- Soft nylon loop lining
- Adjustable strap for increased support

### Neoprene Ankle Support



#### Indications:

- Mild Sprains
- Ankle Strains
- Juvenile Arthritis
- Weak Ankles
- Tendinitis

FEATURES AND BENEFITS					
Sports neoprene for compression and therapeutic warmth					
Supports the ankle					
Limits swelling and decreases pain					
Allows for a full range of motion					
Easy slip-on style, open heel					
Can be worn over or under a sock					

NEOPRENE ANKLE SUPPORT - PEDIATRIC LINE						
Catalog Number	Size	Approx. Age	Ankle Circumference			
40-701204	Ped / Lg	6 - 7 yrs	6.5" - 7.25"			
40-701307	Youth	8 - 12 yrs	7.25" - 8.5"			

NEOPRENE ANK	LE SUPPORT -	REGULAR LINE	
Catalog Number	Size	Approx. Age	Ankle Circumference
40-701405	Small	13+ yrs / Sm Adult	8 - 9.5"

### ProLite® Knit Ankle Support



#### Indications:

- Weak or Injured Ankles
- Ankle Strains
- Tendinitis

Catalog Number Beige	Catalog Number Blue	Size	Approx.	Ankle Circumference
40-400430	40-400400	Youth / Sm	8 - 12 yrs	7" - 8.5"
40-400530	40-400500	Sm Adult / Med	13+ yrs	8.5" - 9.5"

# Comfortable firm support for weak or injured ankles Lightweight knit material is a breathable alternative to neoprene Unique knit design for precise compression Available in blue or beige Latex free



### **Universal Arm Sling**

UNIVERSAL A					
Catalog Number	Size	Approx. Age	Sling Length	Forearm Length	Color
28-502102	Infant	0 - 24 mos	5"	3.5" - 5.5"	Ped Print
28-502204	Ped	2 - 7 yrs	7"	5" - 8"	Ped Print
28-502307	Youth	8 - 12 yrs	9"	7" - 10"	Navy Poly / Ctn

UNIVERSAL ARM SLING - REGULAR LINE					
Catalog Number	Size	Approx. Age	Sling Length	Forearm Length	Color
28-502006	Universal	13+ yrs	17"	11" - 16"	Navy Poly / Ctn



#### FEATURES AND BENEFITS

- Traditional cradle arm sling
- Fast, simple application
- Hook and loop closures for increased adjustability

# Soft Form<sup>®</sup> Posture Control Brace

#### Indications:

- For limiting forward and downward movement of the shoulders (slouching)
- Helps improve poor posture and relieves pain caused by juvenile arthritis
- To encourage proper posture

#### Contraindications:

- Chronic pain associated with juvenile arthritis
- Severe developmental abnormalities or permanent anatomical deformities

SOFT FORM*	POSTURE CONTRO	L BRACE - RE	GULAR LINE
Catalog Number	Size	Approx. Age	Abdomen Circumference
16-900100	Ped / Sm	2 - 4 yrs	20" - 24"
16-900200	Ped / Lg	5 - 7 yrs	23" - 27"
16-900300	Youth	8 - 12 yrs	25" - 29"



# Promotes proper posture Gently holds shoulders back Elastic side panels provide support compression to abdominal and lumbar regions for improved posture and alignment Soft, breathable materials keep skin cool Flexible plastic stays for additional support Criss-cross shoulder straps prevent pinching of the skin and increase adjustability



### **Cervical Collar**



CERVICAL COLLAR - PEDIATRIC LINE					
Catalog Number	Size	Approx. Age	Neck Circumference	Width x Length	
11-111101	Infant	0 - 24 mos	9" - 10.5"	2" x 13"	
11-121202	Ped	2 - 7 yrs	10" - 12"	2.5" x 15"	
11-131306	Youth	8 - 12 yrs	12" - 14"	2.5" x 17"	

#### FEATURES AND BENEFITS

- Low contour maintains head in neutral position
- Simple application with hook and loop closures

### Clavicle Support



CLAVICAL SUPPORT - PEDIATRIC LINE					
Catalog Number	Size	Approx. Age	Fits Chests		
16-101106	Infants	0 - 24 mos	15" - 20"		
16-101307	Ped / XS	2 - 7 vrs	20" - 24"		

CLAVICAL SUP	PORT - REGULAR	LINE	
Catalog Number	Size	Approx. Age	Fits Chests
16-101404	Youth / Sm	8 - 12 yrs	24" - 30"

- Holds shoulders back and limits clavicle movement during treatment of clavicular fractures
- Felt back panel for greater support
- Adjustable top straps for an excellent fit
- Padded straps for comfort
- Treatment of Clavicular Fractures



### Elastic Shoulder Immobilizer

#### ELASTIC SHOULDER IMMOBILIZER - PEDIATRIC LINE Catalog Number Approx. Size Fits Age Chests 16-701104 Infants 0 - 24 mos 16" - 22" 16-701301 Ped / XS 2 - 7 yrs 20" - 26"

ELASTIC SHOU	LDER IMMOBILIZ	ER - REGULAR	LINE	
Catalog Number	Size	Approx. Age	Fits Chests	
16-701408	Youth / Sm	8 - 12 yrs	24" - 30"	

FEATURES AND BENEFITS
Lightweight, low profile
Immobilizes shoulder, reduces abduction and arm rotation
Loop lock closures for easy application
Fits right or left




### Canvas Rocker Bottom Cast Shoe

CANVAS ROCE	KER BOTTOM C	AST SHOE - PEDI	ATRIC LINE	
Catalog		Sole Dim	nensions	
Number	Size	Α	В	Color
43-101107	Infant	5.5"	3"	Pediatric Print
43-101204	Ped	7.25"	3.625"	Navy

CANVAS ROC	CKER BOTTOM CA	ST SHOE - REG	ULAR LINE	
Catalog		Sole Dir	mensions	
Number	Size	Α	В	Color
43-101408	Youth / Sm	9.375"	4.25"	Navy

FEATURES AND BENEFITS
Provides protection for walking casts
Durable canvas material with adjustable closures
Open toe and heel
Non-skid sole
Fits right or left







### Pediatric Print Cast Shoe



The Pediatric Print Cast Shoe, sized and designed for small children. The varying bright prints add an element of fun while the three hook and loop closures provide security for little feet.

PEDIATRIC PRINT CAST SHOE					
Length	Width	Size	Catalog Number		
4.75"	2.5"	XXX Small	PPCR-XXXS		
5.5"	3"	XX Small	PPCR-XXS		

#### FEATURES AND BENEFITS

Non-skid rocker bottom



Part		,	Appropriate Fo	or
Series	Description	Infant	Ped.	Youth
10-622	Philadelphia <sup>®</sup> Collar [Pediatric = 2.25" Tall (8 - 10"), Small = 2.25" Tall (10 - 12")]		Х	Х
11-1X1	Cervical Collar	Χ	X	X
16-101	Clavicle Support	X	X	X
16-420	Posture Control Shoulder Brace (Pediatric = 20 - 24", Youth = 24 - 30" Chest)			Х
16-602	Sling and Swathe			X
16-701	Elastic Shoulder Immobilizer	Χ	X	X
16-900	Soft Form® Posture Control Brace	X	X	
19-400	ProLite® Knitted Pullover Elbow Support (Small = 8 - 9" Forearm Circumference)			X
19-500	GelBand® Arm Band (X-Small = 7 - 11" Forearm Circumference)			X
19-600	EpiSport® Epicondylitis Clasp (X-Small = 8 - 9" Forearm Circumference)			X
19-601	Neoprene Elbow Sleeve	X	X	
19-602	Neoprene Elbow Strap w/Loop Lock (Small = 7 - 11" Forearm Circumference)			X
22-103	Elastic Pullover Wristlet (Small = 5 - 6" Wrist Circumference)		X	
22-106	Elastic Wrist Wrap w/Thumb Loop (Universal = 4 - 11" Wrist Circumference)		X	X
22-109	Neoprene Wrist Wrap (Universal = 4 - 11" Wrist Circumference)		X	X
22-150	Soft Fit Suede Finish Wrist Brace (Medium = 5.5 - 6.5" Wrist Circumference, Specify Right or Left)			X
22-201	Elastic Wrist Brace 6" Length (Small = 4.5 - 5.5", Medium = 5.5 - 6.5" Wrist Circumference, Specify Right or Left)		X	X
22-300	Wrist Splint (Specify Right or Left)  Prol the Wrist Connect Knitted Bullouer (Small   E.E. 6.5" Wrist Circumference)	X	X	X
22-400	ProLite® Wrist Support Knitted Pullover (Small = 5.5 - 6.5" Wrist Circumference)			X
22-442	8" Perforated Suede Finish Wrist Brace (Medium = 5.5 - 6.5" Wrist Circumference, Specify Right or Left)			
22-450	ProLite® 8" Wrist Splint (Small = 5.5 - 6.5" Wrist Circumference, Specify Right or Left)  ProLite® Wrist Splint w/Abducted Thumb (Small = 5.5 - 6.5" Wrist Circumference, Specify Right or Left)			X
22-460				X
22-470	ProLite® 6" Wrist Splint (Small = 5.5 - 6.5" Wrist Circumference, Specify Right or Left)			X
22-560 25-102	Soft Form® Elegant Wrist Support (Small = 5.5 - 6.25" Wrist Circumference, Specify Right or Left)	Х	X	X
	Finger / Thumb Splint w/Moveable Tabs  Tothor® Thumb Stabilizor (Small   2.35   2.75" Hand Width)	^	X	
25-160	Tether® Thumb Stabilizer (Small = 2.25 - 2.75" Hand Width)  Tood Finger Splint (Small = 2.105" Legeth)		X	
25-301 25-401	Toad Finger Splint (Small = 2.125" Length)  Curved / Gutter Finger Splint (X-Small = 1.5", Small = 2.5" Length)		X	X
25-401	Fold-Over / Cot Finger Splint (X-Small = 1.5", Small = 2" Length)		X	X
25-701	Four Prong Finger Splint (X-Small = 1.5", Small = 2" Length)		X	X
25-773	Aluminum Colles Splint, Small (Specify Right or Left)		X	
25-778	Aluminum Wrist & Forearm Splint, Small (Specify Right or Left)		Α	X
25-779	Aluminum Wrist & Foream Splint, Small (Specify Right or Left)  Aluminum Metacarpal Hand Splint, Small (Specify Right or Left)			X
28-502	Universal Arm Sling	Х	Х	X
28-802	Shoulder Immobilizer Sling, Navy/Poly/Cotton (Pediatric = 7" Sling Length, Small = 11" Sling Length)			X
31-701	7" Lumbar Sacral Support (Small = 28 - 32" Circumference at Hips)			X
34-410	Two Panel Surgical Rib Belt, Male (Universal = 24 - 48" Chest Circumference at Base of Sternum)			X
34-420	Two Panel Surgical Rib Belt, Female (Universal = 24 - 48" Chest Circumference at Base of Sternum)			X
37-103	Neoprene Patella Stabilizer		Х	X
37-104	Neoprene Hinged Knee Stabilizing Brace [X-Small fits 12 - 13" (Circumference 4" Above Center of Kneecap)]			X
37-107	Neoprene Stabilizing Knee Brace [X-Small fits 12 - 13" (Circumference 4" Above Center of Kneecap)]			X
37-350	Wrap-Around Hinged Knee Stabilizing Brace [X-Small fits 12 - 13" (Circumference 4" Above Center of Kneecap)]			X
37-373	Standard Neoprene Knee Sleeve w/Open Patella		Х	X
37-374	Standard Neoprene Knee Sleeve w/Closed Patella [X-Small fits 12 - 13" (Circumference 4" Above Center of Kneecap)]			X
37-400	ProLite® Knitted Pullover Knee Support [Small fits 12.5 - 14" (Circumference 4" Above Center of Kneecap)]			X
37-500	GelBand® Patella Strap (Universal = 11-18" Circumference at Kneecap)			X
37-71X	Universal Cutaway Knee Immobilizer, 10", 12" & 14" Lengths, Fits up to 23" (Circumference 4" Above Center of Kneecap)		Х	
38-140	Knee Immobilizer, 12" & 14" Lengths		Х	
40-400	ProLite® Knit Ankle Support		Х	
40-401	Nylon Anklet (Small = 7 - 8.5" Ankle Circumference)		Х	
40-411	Ankle Guard Lace-Up (XX-Small = Men's Shoe Under 6, Women's Shoe Under 7)			Х
40-450	Compressive Ankle Support with Viscoelastic (X-Small = 5 - 7.5" Ankle Circumference)			X
40-701	Neoprene Ankle Support	Х	Х	
40-800	FlexLite® Sport Articulating Hinged Ankle Brace (Small = Men's Shoe 5 - 8.5, Women's Shoe 6.5 - 9.5)			Х
43-101	Canvas Rocker Bottom Cast Shoe ( = Infant & Ped)	Х	Х	X
43-820	C3™ Post-Op Shoe (Small = Women's Shoe 4 - 6)		Х	
43-821	C3™ Post-Op Shoe (Small = Men's Shoe 6 - 8)		Х	
43-890	Super Post-Op Shoe (X-Small = Women's Shoe Under 4, Small = Women's Shoe 4 - 6)			Х
43-891	Super Post-Op Shoe (Small = Men's Shoe 6 - 8)		X	
55-3X0	Cast Protectors (Short Arm, Full Arm, Short Leg, Full Leg)	X	X	





#### USE THIS BACK SIZING CHART FOR ALL BACK SUPPORTS.

To find your proper size, measure your waist just below the navel, and hips at their widest point.

SIZE	WAIST	HIPS	SIZE	WAIST	HIPS	SIZE	WAIST	HIPS
XX-Small	22 - 25"	26 - 29"	Medium	34 - 37"	38 - 41"	XX-Large	46 - 49"	50 -53"
X-Small	26 - 29"	30 - 33"	Large	38 - 41"	42 - 45"	XXX-Large	50 - 53"	54 - 57"
Small	30 - 33"	34 - 37"	X-Large	42 - 45"	46 - 49"	XXXX-Large	54 - 57"	58 - 61"

Certain Back Supports on are manufactured under one or more of the following U.S. Patents: 5,147,261; 5,399,151; 5,499,965; 5,656,020; 5,693;006.



### Safe-T-Lift® LX

Designed for high intensity material handling or when maximum longevity is required. Wide comfortable base, with four coated,

spring steel stays provide excellent support and durability. High grade spandex base allows for long product life, and tension straps protect product from overstretching. Non-slip webbing and extra stay pockets prevent rolling and keep belt firmly in place. Snag resistant

SAFE-T-LIFT* LX				
Size	Catalog Number			
X-Small	70-110349			
Small	70-110446			
Medium	70-110543			
Large	70-110640			
X-Large	70-110747			
XX-Large	70-110844			
XXX-Large	70-110941			

woven suspenders are rugged, supportive and comfortable.

Color: Black.



### Safe-T-Lift® LX2

#### High Visibility Meets OSHA Reflective Guidelines!

This premium product provides support to the lower back and helps promote proper posture and lifting techniques. Wide support base is made of high grade, breathable spandex for extended product life. Color: Orange/Black.

SAFE-T-LIFT® LX2		
Size	Catalog Number	
Small	70-1004B2	
Medium	70-1005B9	
Large	70-1006B6	
X-Large	70-1007B3	
XX-Large	70-1008B0	
XXX-Large	70-1009B7	

- High visibility mylar reflective stitching in shoulder straps/compression panels
- Four coated and plastic tipped spring steel support stays
- Rubberized gripper strips hold belt in place
- Durable woven elastic power compression panels
- Adjustable, snag resistant woven suspenders with snap-away safety buckles



### Safe-T-Lift® DX

Features four polypropylene stays for comfortable support.

Breathable spandex base with non-slip webbing to prevent rolling back and keep the belt firmly in place. Snag resistant woven suspenders are rugged, supportive and comfortable. Durable elastic side pulls provide long product life. Color: Black.

SAFE-T-LIFT* DX		
Size	Catalog Number	
X-Small	70-120346	
Small	70-120443	
Medium	70-120540	
Large	70-120647	
X-Large	70-120744	
XX-Large	70-120841	
XXX-Large	70-122946	



### **Custom Fit®**

Designed to stay in place on a wide variety of body types and in a wide variety of situations. Patented contour design prevents belt from riding up when sitting or bending. Ideal for females, unique physiques or heavier body types. Adjustable shaping strap provides a tailored fit. Spandex base and non-slip webbing keeps belt firmly in place. Snag resistant woven suspenders, break away buckles and side pulls for greater safety and durability.

CUSTOM FI	T °
Size	Black Catalog Number
X-Small	70-160344
Small	70-160441
Medium	70-160548
Large	70-160645
X-Large	70-160742
XX-Large	70-160849
XXX-Large	70-160946





### Safe-T-Belt® Working Lumbar Belt



Ideal choice when shoulder straps are not desired. Excellent choice for workers who spend time sitting or need to be unencumbered by shoulder straps. Contoured base prevents riding up when sitting or bending. Two polypropylene stays provide excellent support and durability. Compression pad provides focused support to the lumbar region. Overlapping elastic side pulls provide excellent support and freedom of movement. Color: Black.

SAFE-T-BELT*		
Size	Catalog Number	
X-Small	70-910341	
Small	70-910448	
Medium	70-910545	
Large	70-910642	
X-Large	70-910749	
XX-Large	70-910846	

### Safe-T-Belt® Plus



Contoured base prevents riding up when sitting or bending. Two polypropylene stays provide excellent support and durability. Compression pad provides focused support to the lumbar region. Overlapping elastic side pulls provide excellent support and freedom of movement. Snag resistant Safe-T-Straps® with woven suspenders and break away buckles for greater safety and durability. Color: Black.

SAFE-T-BELT® PLUS	
Size	Catalog Number
X-Small	70-920348
Small	70-920445
Medium	70-920542
Large	70-920649
X-Large	70-920746
XX-Large	70-920843
XXX-Large	70-920946

### DynaBack<sup>®</sup>



Universal sizing fits 24 to 46 inch waists. Back panel is made of durable elastic with polypropylene stays for added support. Shoulder straps are adjusted or removable for flexibility and comfort. Color: Black.

DYNABACK*	
Size	Order #
Universal	73-130003

### Safe-T-Wrist® HD Heavy Duty

Designed for tough, highly repetitive tasks, as well as for use with Carpal Tunnel Syndrome. Flexible palmar and dorsal stays provide firm support while the contoured palm design allows for full finger dexterity. Plush elastic base material for long life and high durability, vinyl palm pad. Full elastic wrist wrap provides excellent compression. Color: Black. Specify right or left.

SAFE-T-WRIST*	HD HEAVY DUTY			
Right Catalog Number	Left Catalog Number	Size	Fits Wrist	Palm Width
71-111444	71-112443	Small	5" - 6"	2.75" - 3.5"
71-111541	71-112540	Medium	6" - 7"	3.5" - 4"
71-111648	71-112647	Large	7" - 8"	4" - 5.5"
71-111745	71-112744	X-Large	8" - 9"	5.5" - 6"



#### FEATURES AND BENEFITS

- Stabilizes the wrist with full finger dexterity
- · Limits wrist flexion and extension
- Palmar and dorsal flexible plastic stays stabilize both sides of the wrist
- Full elastic wrist wrap for greater stability, minimizes wrist rotation
- Durable material and vinyl covered palm area
- Palm area padded inside with moisture wicking material for comfortable, extended wear
- Sleek, low profile design

### Safe-T-Wrist® SD Standard Duty

Provides moderate support. Ideal for tasks that require extension, flexion and rotation of the wrist. High quality neoprene absorbs perspiration while providing therapeutic warmth. Hook and loop closure is removable. Fits left or right. Color: Black.

SAFE-T-WRIST* SD STANDARD DUTY				
Right Catalog Number	Size	Fits Wrists	Palm Width	
71-210444	Small	Under 5.5"	Under 3"	
71-210541	Medium	5.5" - 6.5"	3" - 3.75"	
71-210648	Large	6.5" - 7.5"	3.75" - 4.5"	
71-210745	X-Large	Over 7.5"	Over 4.5"	

- Sports neoprene provides soothing warmth
- Removable, adjustable support strap for increased stability
- Open palm
- Sleek, low profile design
- Fits right or left







Durable, comfortable protection for the hands. GelShield® insert redistributes force away from the hands. Constructed of high quality leather with an absorbent terry cloth wipe pad. Flexible, breathable spandex. Hook and loop closure for easy fit. Washable for continued wear. Available with or without elastic wrist wrap. Packaged as a pair. Color: Black.

SAFE-T-GLOVE* VIBRATION DAMPENING GLOVE			
Catalog Number	w/Wrist Strap Catalog Number	Size	Fits Palm
71-610440	71-620447	Small	3"
71-610547	71-620544	Medium	3.5"
71-610644	71-620641	Large	4"
71-610741	71-620748	X-Large	4.5"

Measure the widest part of the palm.

138



### Forearm Splint

The Forearm Splint can be used for shaft fractures of the ulna and radius and as synovitis of the wrist, sprains and colles fractures where an abducted hand position is not desired.

FOREARM SPLINT	
Size	Catalog Number
Adult - Right	Q101501
Adult - Left	101502
Youth - Right	101503
Youth - Left	101504
Child - Right	Q101505
Child - Left	101506
Infant - Right	Q101507
Infant - Left	Q101508



### Metacarpal Forearm Splint

The Metacarpal Forearm Splint is useful in the reduction of radius, ulna, metacarpal and phalange fractures.

Size C	atalaa Numbar
	atalog Number
Adult - Right Q	2101701
Adult - Left 10	01702
Youth - Right	01703
Youth - Left 10	01704
Child - Right	01705
Child - Left 10	01706
Infant - Right 10	01707
Infant - Left 10	01708



### Posterior Leg / Knee Splint

The Posterior Leg / Knee Splint is used with fractures of tibia / fibula, ankle fractures also used to immobilize the lower leg during dressing application.

POSTERIOR LEG/KNEE SPLINT	
Size	Catalog Number
Adult	102601
Youth	102602
Child	102603





### **Aluminum Hand Splints**

Made of durable coated aluminum with a soft foam lining. Ventilation holes promote air circulation for greater comfort. Hand splints are packaged in durable, hang tag ready polybags.

#### Wrist & Forearm Splint

Provides temporary support for wrist and forearm injuries prior to surgery or casting. Specify right or left.



WRIST & FOREARI	M SPLINT	
Right Catalog Number EA	Left Catalog Number EA	Size
25-778417	25-778426	Small
25-778514	25-778523	Medium
25-778611	25-778620	Large

#### **Metacarpal Hand Splint**

Designed to provide support for metacarpal and phalangeal fractures. Specify right or left.



METACARPAL HAN	D SPLINT	
Right Catalog Number EA	Left Catalog Number EA	Size
25-779416	25-779425	Small
25-779610	25-779629	Large

#### **Colles Splint**

Provides temporary support of colles injuries prior to surgery or casting. Specify right or left.



COLLES SPLINT		
Right Catalog Number EA	Left Catalog Number EA	Size
25-773412	25-773421	Small
25-773519	25-773528	Medium
25-773616	25-773625	Large

### Poncho Style Vest

Prevent patients from falling or slipping out of bed or wheelchair. Comfortable over-the-head poncho style vest prevents shoulders from slipping out. Vest closes in front and ties are secured out of patient's reach. Cool open weave mesh. Binding is color coded by size.

PONCHO STYL	E VEST		
Catalog Number	Size	Fits Chests	Binding Color
49-104409	Small	30" - 36"	Grey
49-104506	Medium	37" - 42"	Yellow
49-104603	Large	43" - 48"	White
49-104700	X-Large	49" - 52"	Black





### Patient Walker Belt

When worn around the patient's waist, this belt assists care givers by allowing for a secure hold on the patient during walking or exercising. The belt is made of a soft material which is extremely strong and durable. Closure is fast, simple and secure with the two-sided sure-lock buckle. Length is easily adjusted to fit a wide range of sizes. Two convenient hand holds are located on back of belt for added stability.

PATIENT WALKER BELT				
Catalog Number	6-Pack Catalog Number	Size	Chests	
52-150403	_	Small	22" - 30"	
52-150005	52-1500PK	Universal	30" - 48"	
52-150607	_	Large	38" - 60"	



### Catheter Strap

Helps provide security against irritation by reducing catheter movement. Made of two inch wide elastic with hook closure. Universal size.

CATHETER STRAP		
Catalog Number	Size	
49-240002	Universal	



### **Heel and Elbow Protector**

Protector helps reduce pressure and chance of decubitus ulcers on injured or sensitive heels or elbows. A soft foam pad is encased in a seamless elastic sock which stretches and conforms for a perfect fit. Made of lightweight mesh material. Universal size fits a range of sizes. Color: White.



HEEL and	ELBOW	PROTECTOR	
Catalog N	lumber	Size	

### Actimove® Heel Protector



The Actimove'Heel protector is ideal protection for bed-bound patients. It consists of a synthetic sheep skin for comfortable heel cushioning and easy-to-apply hook and loop straps. Universal size and sold as a pair.

ACTIMOVE® HEEL	PROTECTOR	
Dimensions	Catalog Number	Quantity
Universal	72819-00070-00	1 pair



# Flotation Heel and Elbow Protector

Helps reduce pressure and chance of decubitus ulcers on injured or sensitive heels or elbows. Unique memory foam pad is encased in an elastic support for maximum protection. Universal size fits a range of sizes. Color: Light Blue.

#### FLOTATION HEEL and ELBOW PROTECTOR

Catalog Number	Size	
61-130008	Universal	



### Thermal Wraps

For either hot or cold therapy. Ideal for treatment of minor injuries such as scrapes, bruises, sprains, strains and chronic aches and pains. Reusable hot/cold compresses which consist of one hot/cold gel pack and a universal wrap cover with a soft foam material that is hook compatible with a soft, absorbent terry cloth liner. Smaller size  $(4 \times 15)$  is ideal for the neck, leg and ankle. Larger size  $(6 \times 10)$  is ideal for the back, shoulders and thighs.

#### THERMAL WRAPS

4" x 15" 6" x 10"

Neck / Leg Back / Shoulder

Catalog Number

53-120061 53-121066





Viscolas® is the only viscoelastic polymer expressly formulated for orthopedic applications. Reproducing the shock attenuating properties of the fatty pads on the bottom of the foot, Viscolas® effectively dissipates the skeletal shock and reverberations that occur with each step, and equalizes underfoot pressures. Use of Viscolas® has proven to speed post-operative recovery of the foot, leg, and lower-back; reduces symptoms of heel strike injuries such as tendinitis and shin splints; and provides effective soft tissue cushioning under the foot. In trial after trial, Viscolas® has been proven to provide significantly better shock absorbing efficiency than any other viscoelastic polymer. Insoles made of Viscolas® are the ones most prescribed by orthopedic surgeons, therapists, and other health professionals... and they are the insoles most chosen for personal use by the staffs of orthopedic facilities.

Viscolas is a registered trademark of the CABOT CORPORATION

### Soft Point<sup>®</sup> Viscolas<sup>®</sup> Heel Spur Cushions

Removable plug creates a cavity beneath the area of greatest discomfort. They also provide extra cushioning on the medial side of the plantar fascia for additional relief. Covered with super-soft, ultra-absorbent, abrasion-resistant suede fabric to keep the foot dry and comfortable.

#### Indicated for:

- Heel Pain
- Calcaneal Heel Spurs
- Plantar Fasciitis

"VISCOLAS"	HEEL SPUR CUSHI	ONS
Size	Men's Shoe Size	Women's Shoe Size
А	2 - 4	4 - 6
В	5 - 7	7 - 9
С	8 - 10	10 - 12
D	11 - 13	-
	<b>Size</b> A B	Size         Shoe Size           A         2 - 4           B         5 - 7           C         8 - 10

### Soft Point® Viscolas® Heel Pads

This low profile pad is ideal for most shoe types. It is extremely easy and comfortable to wear. Natural tacky surface to keep heel pad in place, covered with super-soft, ultra-absorbent, abrasion-resistant suede fabric to keep the foot dry and comfortable.

#### Indicated for:

- Heel Pain
- Fat Pad Atrophy
- Leg Length
   Discrepancies

VISCOLAS*	HEEL PADS	
Size	Men's Shoe Size	Women's Shoe Size
А	2 - 4	4 - 6
В	5 - 7	7 - 9
С	8 - 10	10 - 12
D	11 - 13	-
	Size A B C	Size         Shoe Size           A         2 - 4           B         5 - 7           C         8 - 10



# Soft Point® Viscolas® Full Insoles

Viscolas® Full Insoles have a thin, low profile that is ideal for a wide variety of shoe types and styles. Slim toe area for a comfortable, non-binding fit. Extra padding at the heel for increased shock absorption. Natural tacky surface to keep insole in place. Covered with super-soft, ultra-absorbent, abrasion-resistant suede fabric to keep the foot dry and comfortable.

#### Indicated for:

- Lower Back, Foot & Leg Fatigue
- Fat Pad Atrophy
- Post-Operative Conditions
- Shin Splints

SOFT POINT®	VISCOLAS*	FULL INSOLES		
Catalog Number	Size	Men's Shoe Size	Women's Shoe Size	
F251	А	2 - 4	4 - 6	
F252	В	5 - 7	7 - 9	
F253	С	8 - 10	10 - 12	
F254	D	11 - 13	-	



### Soft Point® Viscolas® Sports Full Insoles

The Soft Point® Viscolas® Sports Full Insoles are made with original Viscolas® viscoelastic polymer inserts for superior cushioning at high-impact areas. The shape of the insoles also provides arch support and heel stabilization. This insole is ideal for athletic shoes.

#### Indicated for:

- Athletic Activity
- Foot & Leg Pain
- Shin Splints
- Medial Ankle Pain

SOFT POINT	VISCOLAS*	FULL INSOLES	
Catalog Number	Size	Men's Shoe Size	Women's Shoe Size
F261	А	2 - 4	4 - 6
F262	В	5 - 7	7 - 9
F263	С	8 - 10	10 - 12
F264	D	11 - 13	-





SOFT POINT® Footcare Products are made of medical grade silicone for long lasting support and comfort. Silicone gives optimal absorption of the shock generated from everyday and high-impact heel strike. The impact and force on the heel, joints and back are thereby reduced helping to relieve pain and discomfort. Non-irritating and odor free, silicone will always remain soft and retain its shape.

Many of the SOFT POINT® Footcare Products have blue inserts which have a softer durometer profile for enhanced comfort. The blue inserts are scientifically placed to help relieve pressure and pain on the most sensitive areas of the foot.

Soft Point is a registered trademark of FLA

### Soft Point® Full Insoles



The blue inserts in the Full Insoles have a softer durometer profile in the heel and metatarsal area for enhanced comfort. The positioning of the blue inserts helps relieve pressure and pain on the most sensitive areas of the foot. The shape of the Full Insole gives slight arch support and has a tapered, lower profile at the toes for an improved fit.

SOFT POINT	FULL INSOLE	s	
Catalog Number	Size	Men's Shoe Size	Women's Shoe Size
F151	Small	4 - 5	6 - 7
F152	Medium	6 - 7	8 - 9
F153	Large	8 - 10	10 - 11
F154	X-Large	10.5 - 12	-

### Soft Point® 3/4 Insoles



The blue inserts in the ¾ Insoles have a softer durometer profile in the heel and metatarsal area for enhanced comfort. The positioning of the softer, blue inserts helps relieve pressure and pain on the most sensitive areas of the foot. The shape of the ¾ Insole gives additional arch support and heel stabilization.

SOFT POINT	* 3/4 INSOLES		
Catalog Number	Size	Men's Shoe Size	Women's Shoe Size
F140	X-Small	4 - 5	5.5 - 6.5
F141	Small	5.5 - 6.5	7 - 8
F142	Medium	7 - 8.5	8.5 - 10
F143	Large	9 - 10	10.5 - 11
F144	X-Large	10.5 - 12	-



### Soft Point® Heel Pads

This low profile pad is ideal for most shoe types. It is extremely easy and comfortable to wear with a hollowed center to relieve pressure on the painful heel area – ideal for heel spurs.

SOFT POINT	* HEEL PADS		
Catalog Number	Size	Men's Shoe Size	Women's Shoe Size
F111	Small	1 - 5	5 - 6.5
F112	Medium	5.5 - 9.5	7 - 11.5
F113	Large	10 - 15	-



### Soft Point® Heel Cushion

The blue spot in the Heel Cushion has a softer durometer profile at the heel for enhanced comfort. The shape of the cushion and the medial position of the soft, blue spot help relieve pressure and pain on the most sensitive area of the plantar fascia.

SOFT POINT	* HEEL CUSHI	O N	
Catalog Number	Size	Men's Shoe Size	Women's Shoe Size
F131	Small	3.5 - 7.5	5 - 8
F133	Large	8 - 13	8.5 & up





# Soft Point® 3/4 Semi-Rigid Orthotics



These <sup>3</sup>/<sub>4</sub> orthotics provide superior support to the arch of the foot and relieve the pain associated with many foot conditions. Deep contours cradle the heel and create full arch contact to comfortably stabilize the foot. The flexible, semi-rigid base gives spring-action flex with each step providing continuous, comfortable arch support.

#### Ideal for:

Relief from Plantar Fasciitis; Forefoot, Arch and Heel Pain; Back and Joint Ache; Weak or Fallen Arches; Aching and Fatigued Feet; and Prevention of Pronation, or the inward rolling of the ankles.

SOFT POINT	3/4 SEMI-R	IGID ORTHOTICS	
Catalog Number	Size	Men's Shoe Size	Women's Shoe Size
F541	В	-	5 - 7
F542	С	5 - 7	7 - 9
F543	D	7 - 9	9 - 11
F544	Е	9 - 11	-
F545	F	11 - 13	-

Size according to US shoe size.

- Shock Absorbing Poron® Liner Soft and flexible liner absorbs the shock
  of each step, removing the shock effects from the feet and joints. Poron's
  super-soft quality means maximum comfort.
- Breathable Fabric reduces perspiration, allowing for extended, comfortable wear and reduces shear to minimize blisters or calluses
- 3/4 Length and Moderate Width relieves pressure on the forefoot and accommodates a broad range of shoes – from walking and athletic shoes to casual and dress shoes. Can easily be moved from shoe to shoe.
- Flexible, Semi-Rigid Polypropylene Base is strong enough to provide solid foot support and stability. Unique spring action flexes with each step to continuously support the arch.
- High Side Rims allow for full contact with the arch, stabilizing the feet and preventing ankles from rolling inward.
- Deep Heel Cup cradles the heel and helps stabilize the foot.
- Biomechanically Designed to provide longitudinal support across the foot and full support to the arch.



# Soft Point® 3/4 Flexible Cork Orthotics

These <sup>3</sup>/<sub>4</sub> orthotics provide superior support to the arch of the foot and relieve the pain associated with many foot conditions. Deep contours cradle the heel and create full arch contact to comfortably stabilize the foot. The extremely lightweight and durable Cork/EVA composite material molds to the shape of the foot with the body's natural heat and weight, creating a customized feel.

#### Ideal for:

Relief from Plantar Fasciitis; Forefoot, Arch and Heel Pain; Back and Joint Ache; Weak or Fallen Arches; Aching and Fatigued Feet; and Prevention of Pronation, or the inward rolling of the ankles.

SOFT POINT	Γ° 3/4 FLEXIBL	E CORK ORTHOTIC	cs
Catalog Number	Size	Men's Shoe Size	Women's Shoe Size
F640	Small	-	5 - 7
F641	Medium	6 - 7	8 - 10
F642	Large	8 - 10	-
F643	X-Large	11 - 13	-

Size according to US shoe size.

- Extremely thin, flexible, lightweight base is super lightweight, durable and shock absorbing.
- Cork / EVA material molds to natural contours of the foot to give a comfortable, custom-like fit and absorb shock.
- Biomechanically designed to provide longitudinal support across the foot and full support to the arch.
- Breathable fabric reduces perspiration, allowing for extended, comfortable wear and reduces shear to minimize blisters or calluses.
- Deep heel cup cradles the heel and helps stabilize the foot.
- High side rims allow for full contact with the arch, stabilizing the feet and preventing ankles from rolling inward.
- 3/4 length and moderate width relieves pressure on the forefoot and accommodates a broad range of shoes – from walking and athletic shoes to casual and dress shoes. Can easily be moved from shoe to shoe.











#### BSN medical Inc.

#### TERMS AND CONDITIONS OF SALE AND DELIVERY

BSN medical Inc., is the "Seller" for all transactions contemplated below. Any person buying or offering to buy goods or services from Seller is below referred to as the "Buyer." Products" as used below means the products and any related services sold or supplied, including any replacement products or services.

#### I. Modification of Material Furnished by Seller

- 1.1 Catalogs and Written Materials. All catalogs and other material of Seller are subject to modification by Seller and are not binding unless so stated in writing by Seller. Seller may correct clerical errors and typos in any catalogs or other materials, including pricing contracts, at any time.
- 1.2 All prices listed in any catalogs, price lists or other materials of Seller are subject to modification by Seller and are not binding unless so stated in writing by Seller.
- 1.3 Product Discontinuation. Seller reserves the right to modify or discontinue Products at any time without notice.

#### II. Prices

- 2.1 Unless Seller has specified in writing that any quotation is binding for a specified period not yet expired, price quotes may be changed at any time prior to shipment of the Products. Unless otherwise indicated in writing by Seller, all pricing is freight pre-paid to Buyer's U.S. warehouse or other continental U.S. shipping point designated by Buyer. The invoice price shall be the price in effect when Seller delivers the Products pursuant to Article III hereof.
- 2.2 Other Charges: Taxes. Unless otherwise agreed by Seller in writing, sales, use, value added, excise, or similar taxes, as well as insurance (which Seller does not provide) are not included in the price.

#### III. Delivery

- 3.1 Delivery Terms; Errors. Unless otherwise agreed by the parties in writing, delivery is deemed made and risk of loss passes at the point to which Seller has paid freight; otherwise delivery is made and risk of loss passes at Seller's warehouse or other shipping point. Any shipping or pricing errors are to be reported by Buyer to Seller within thirty (30) days after delivery, or shall otherwise be deemed waived.
- 3.2 Delivery Dates. Any delivery dates indicated herein or otherwise set by Seller do not bind Seller unless specifically agreed in writing by Seller.
- 3.3 Separate and Partial Shipments. Each shipment by Seller is a separate and independent transaction and may be invoiced separately by Seller. In addition, Seller may make partial shipments and submit separate invoices for each such partial shipment.
- 3.4 Shipping Address. It shall be Buyer's sole responsibility to ensure that the correct shipping address is promptly communicated to Seller. Seller is not required to contact Buyer to obtain or confirm any shipping destination.

#### IV. Terms

- 4.1 Payment Terms. Payments shall be made in U.S. Dollars, unless Seller agrees in writing to invoice in another currency. Time is of the essence in regard to payment of any amounts.
- 4.2 Bank Charges. All bank charges in respect of payment instruments shall be paid by Buyer immediately.
- 4.3 Interest. If the amounts owed are in U.S. dollars, interest will be charged daily on past due amounts at a rate of the lesser of (i) one and one-half percent (1.5%) per month and (H) the maximum lawful rate of interest. If the amounts owed are in a currency other than U.S. dollars, interest on past due accounts will be charged at the rate of interest customarily charged by Seller on sales in the currency specified.
  4.4 Change in Terms. If Seller concludes in its sole discretion that the financial condition of Buyer at any
- 4.4 Change in Terms. If Seller concludes in its sole discretion that the financial condition of Buyer at any time jeopardizes its ability to pay, Seller may require cash payments or otherwise change payment terms or require additional security satisfactory to Seller before further performance by Seller. Buyer's failure to timely pay any invoice under any Contract shall make all other invoices of Seller immediately due and payable and, at the discretion of Seller, shall be grounds for suspension or cancellation of any further performance by Seller. The receipt by Seller of part payment shall not constitute a waiver of any of Seller's rights herein or at law, including the right to cancel.

#### V. Limitation of Warranty

- 5.1 Warranty. Seller warrants that the Products are manufactured and function in accordance with the Seller's ratings and specifications. All Products are to be inspected fully and all obvious non-conformities reported in writing to Seller within thirty (30) days (in the case of custom-made Products, fifteen (15) days) after receipt by the end user, even if samples were previously sent. Seller, upon being satisfied of the existence of such non-conformity for which it is liable, will correct the same by delivering replacement Product to Buyer. If Seller cannot correct such non-conformity by delivering replacement Product, Seller will return to Buyer the purchase price thereof.

  The foregoing shall be Buyer's sole and exclusive remedies, and the liability of Seller hereunder is expressly limited to replacement of nonconforming Products or the repayment of the purchase price, as the case may be. The warranties set forth apply both to original and replacement Products.
- 5.2 No Informal Warranties. Any advice, information, or suggestions given by Seller to Buyer in regard to the application or use of the Products are not warranties unless given in writing and designated specifically by Seller as warranties. Statements made in Seller's general advertising or printed material do not constitute warranties. Samples are not warranted as conforming to Seller's ratings and specifications, nor is any warranty made that Products conform to any sample, unless Seller specifically so warrants in writing. Seller makes no warranty that any Product delivered under this Contract is identical in nature or quality to any Product previously sold or delivered to Buyer.
- 5.3 Warranty Exclusions. THE FOREGOING WARRANTIES SHALL BE IN LIEU OF ALL OTHER WARRANTIES, AND THE WARRANTIES OF MERCHANTABILITY AND FITNESS FOR PURPOSE AND OTHERWISE, EXPRESS OR IMPLIED, ARE EXCLUDED.
- 5.4 Return or Destruction. No Products shall be returned or destroyed without Seller's prior consent

#### VI. Confidentiality; Property Rights

- 6.1 Confidentiality and Use. Buyer agrees that all specifications, data and other technical, business, or financial information furnished by Seller to Buyer are the property of Seller, are furnished solely for Buyer's performance hereunder, and may not be copied or made accessible to third parties without Seller's prior written consent. Buyer shall promptly return such specifications, data and other technical, business, or financial information and all copies thereof to Seller upon Seller's request. Such request may be made at any time prior to or after delivery of the Products. The foregoing shall survive cancellation or completion of the Contract.
- 6.2 Identification. Seller reserves the right to place its firm name, symbol, and relevant manufacturing and product identification numbers on Products delivered pursuant hereto.

#### Cancellation

7.1 No Cancellation by Buyer. The Contract may not be cancelled by Buyer without the prior written permission of Seller, and Buyer shall hold Seller harmless from any and all costs and damages flowing from any wronful cancellation.

- 7.2 Products Made to Order. If Buyer purports to cancel the Contract in violation of Paragraph 7.1 in respect of any Products specially manufactured for Buyer, Seller shall also be reimbursed by Buyer for all direct costs of labor, raw materials and parts purchased or contracted to be purchased and overhead costs incurred by Seller in the manufacture of said Products. Buyer agrees that Seller's cost calculation shall be accepted by the parties hereto. In addition, Seller may complete and sell to a third party any specially manufactured Products for which the order has been purportedly cancelled by Buyer, without any liability whatsoever to Buyer.
- 7.3 Seller's Cancellation. If Buyer shall (I) become insolvent or otherwise fail to pay its debts on time; (II) be placed in receivership; (III) be the subject of any voluntary or involuntary bankruptcy petition; or (IV) default in the performance of its obligation under this or any other Contract and fail to cure such default within thirty (30) days after written notice thereof, Seller may cancel all or any portion of Buyer's pending orders under this Contract.

#### VIII. Force Maieure

- 8.1 Definition. "Force Majeure" means any cause not within the reasonable control of the party affected, but no payment of monies can be excused by Force Majeure. Events of Force Majeure include, without limitation, acts of God, lockouts or other labor disturbances, wars, blockades, quarantine embargoes, riots, fires, explosions, failures of production facilities, shortages of fuel, transportation, utilities, or raw materials, and governmental laws and regulations.
- 8.2 Consequences. The party prevented from performing shall promptly so notify the other party hereto and shall provide notice of the termination thereof. Upon termination of the Force Majeure the performance of any suspended obligation shall recommence. Either party may terminate the Contract if any Force Majeure lasts more than ninety (90) days. Buyer's insolvency may not be asserted as the basis for Force Majeure.

#### IX. Miscellaneous

- Buyer's Compliance with Medical Devices Act. Products manufactured by or for Seller are medical devices regulated by the Food and Drug Administration under the 1976 Medical Device Amendments. Some of these Products are labeled with a statement required for Prescription Devices, namely, "Caution: Federal (U.S.) Law Restricts This Device to Sale by or on the Order of Physician." Pursuant to federal and state law, such Products must only be sold to or on the prescription or order of a duly licensed healthcare professional. In addition, state law may prohibit substitution of a medical device to another brand without receiving permission from the prescribing healthcare professional. Buyer agrees to indemnify and hold Seller harmless from, and to reimburse it for, all liabilities, claims, suits, judgments, costs and expenses, including reasonable attorney's fees and disbursements, arising out of or relating to any failure of Buyer to comply with any such laws, rules, or regulations.
- 9.2 Trademarks. The term BSN medical Inc. is a trademark of Seller's parent company and its logo, as well as other trademarks, are U.S. registered trademarks and service marks. They are, and shall remain, the exclusive property of Seller's parent. It is illegal to use these trademarks, service marks, and any other trademarks or copyrighted marketing material belonging to Seller or Seller's parent in advertising, promotion, or in association with the name of a dealership/distributor or in any other way, without receiving prior written permission from Seller or Seller's parent to do so.
- 9.3 No Waiver. The failure of either party to insist upon strict performance by the other party of any term hereof or to exercise any right hereunder shall not be deemed to be a modification of or a waiver of the future performance of such term.
- 9.4 Governing Law. This form and the Contract shall be governed by North Carolina law
- 9.5 Arbitration. Any controversy or claim arising out of or relating to this form or the Contract, or the negotiation or breach hereof or thereof, or any sales pursuant hereto or thereto, shall be settled by arbitration in accordance with the Commercial Arbitration Rules of the American Arbitration Association, and judgment upon the award rendered by the arbitrator may be entered in any court having jurisdiction thereof. The arbitration shall be held in Charlotte, North Carolina, before a single arbitrator mutually agreeable to Seller and Buyer, or if no agreement can be reached, then selected by the American Arbitration Association. The arbitrator shall render a reasoned opinion and shall award reimbursement of attorneys' and other experts' fees and disbursements and other costs of arbitration to the prevailing party, in such manner as the arbitrator shall deem appropriate.

#### X. Damages. Limits and Exclusions

- 10.1 In no event shall Seller's liability for any damages arising out of the sale of Products, regardless of the legal theory on which such damages may be based, exceed the amount that Seller has been paid for the Products under this Contract.
- 10.2 IT IS EXPRESSLY UNDERSTOOD AND AGREED THAT SELLER SHALL NOT BE LIABLE TO BUYER OR ANY OTHER PERSON, WHETHER BY WAY OF INDEMNIFICATION OR CONTRIBUTION OR OTHERWISE, FOR SPECIAL, CONSEQUENTIAL OR PUNITIVE DAMAGES, WHETHER ARISING FROM DEFECTIVE WORKMANSHIP OR MATERIALS OR DESIGN, BREACH OF WARRANTY, DELAYS IN DELIVERY OR OTHER BREACH OF CONTRACT, FORCE MAJEURE, OR FROM ANY OTHER CAUSE WHATSOEVER, INCLUDING THE NEGLIGENCE OR MISCONDUCT OF SELLER OR ITS AGENTS OR EMPLOYEES.