



Latex-Free Mobile Aneroid Blood Pressure Monitor
ANEROID, MOBILE, ADULT CUFF

- Large, raised numbers and a luminescent coating make our mobile aneroid easy to read in dim light
- Eight-foot coiled tubing and pivoting bracket extend service area
- Inflation system and bulb assembly store easily and safely behind gauge
- Includes Medline's nylon range finder cuff
- Neoprene, latex-free

Packaging	1 Case / Each
Manuf / Supplier	Medline
Manuf / Supplier #	MDS9407LF

SKU **MDS9407LF**

[Material Data Sheet Documents](#)

No MSDS content is available for this product.

Specifications

Specifications	Values
Aneroid BP Type	Mobile
Combination Kit	N
Cuff Included	Yes
HPIS Classification	650_10_40_20
Latex Free	Yes
Mobile Unit BP	Y
Product Size	Adult
Reusable	Yes
Single Patient Use	No
UNSPSC code	42181601

Packaging Details

Unit of Measure	Conversion	Net / Gross Weight (lbs)	Volume (cubic ft)	Dimensions (inch) L x W x H	GTIN
Case (CS)	1.0 Each	0.0 / 12.0	1.885	20.5 x 20.5 x 7.75	40080196771340
Each (EA)	1.0 Each	0.0 / 12.0	1.885	20.5 x 20.5 x 7.75	80196771342

HCPCS Information

HCPCS Code: NO CODE - PRODUCT DOES NOT MEET HCPCS SPECIFICATIONS

HCPCS Disclaimer: HCPCS codes and Home Health Consolidated Billing codes provided by Medline are intended as general guidelines only. Medline does not guarantee coverage or reimbursement of any products. You must address all coverage and reimbursement issues (including the correctness and accuracy of codes) with your individual payers. It is your responsibility to ensure the accuracy and appropriateness of each claim you submit, in accordance with all applicable payer requirements.