



Enteral Feeding, Access & GI Products

1-800-962-9888

www.SafeEnteralConnections.com



COVIDIEN

Table of Contents

Enteral Feeding, Access & GI Products

ENTERAL FEEDING

Enteral Feeding Pumps/Pump Sets

Kangaroo™ ePump	1
Kangaroo™ Joey	2
Kangaroo™ 924 Feeding Sets	3

Screw Spike Compatible Feeding Sets

Kangaroo™ Feeding Sets for use with Screw Spike System	3
--	---

Pump Sets

Kangaroo™ DEHP-Free	4
Kangaroo™ DEHP-Free Anti-Free Flow	4
Kangaroo™ Anti-Free Flow	4
Kangaroo™ Pump Sets	5
Kangaroo™ Gravity	6

Enteral Irrigation Kits

Kangaroo™ I.V. Pole Pouch and Kits	6
--	---

ENTERAL ACCESS

Nasogastric Feeding Tubes Neonatal and Pediatric Feeding System

Kangaroo™ Purple Polyurethane Feeding Tubes	7
Kangaroo™ PVC with Sentinel Line™	7
Kangaroo™ Extension Sets	7
Monoject™ Purple Oral Syringes	8
Kangaroo™ Polyurethane Feeding Tubes	9
Kangaroo™ PVC Feeding Tubes with Radiopaque Line	9
Curity™ PVC Feeding Tubes	7

Nasogastric Feeding Tubes Adult

Kangaroo™ with Dobbhoff™ Tip	10
Kangaroo™ with Standard Tip	11
Kangaroo™ Non-Weighted Feeding Tubes	11

CO₂ Detector

Kangaroo™ CO ₂ Detector	12
--	----

Specialty Tubes

Kangaroo™ Naso-Jejunal Feeding and Gastric Decompression Tube	11
Kangaroo™ Naso-Jejunal Feeding Tubes	13
Kangaroo™ Jejunal Feeding and Gastric Decompression Tube	13

Safety P.E.G. Systems

Kangaroo™ Safety P.E.G. Kits	13
Kangaroo™ Safety Single Pass P.E.G. Kits	14
Kangaroo™ Y-Port-G-Tube/P.E.G. Adapter	14

Gastrostomy and Skin Level Tubes

Kangaroo™ Gastrostomy Tubes	15
Kangaroo™ Skin Level Balloon Gastrostomy Kits 16–17	

GI PRODUCTS

Salem Sump™ Tubes

Kangaroo™ Multifunction Port	18
Kangaroo™ Multifunction Silicone Port	18
Salem Sump™ PVC Tube	19
Salem Sump™ PVC Tube w/Anti-Reflux Valve	19
Salem Sump™ Anti-Reflux Valve	19
Salem Sump™ Silicone Tube	20
Salem Sump™ Silicone Tube w/Anti-Reflux Valve..	20

Stomach Tubes

Argyle™ Stomach Tube (Levin type)	21
Stomach Tube (Levin type)	21
Argyle™ Edlich Gastric Lavage Trays	21

Index 22–25

Enteral Feeding Pumps/Pump Sets



382400

Kangaroo™ ePump Enteral Feeding Pump

- Attitude independent
- "Smart" pump set technology
- "ATM-style" ease of operation
- Feeding history up to 72 hours
- DEHP-Free pump sets
- Anti-Free-flow pump sets
- Flushing/hydration capabilities
- Features safe enteral connection system

Specifications

- 15 hour battery
- Weight – 2.34 lbs without pole clamp
2.69 lbs with pole clamp
- Feed Rate – 1-400 in 1 mL increments
- Flush Volume – 10-500 in 1 mL increments
- Flush Interval – 1-24 in 1 hr increments
- Multilingual – 15 languages
- Memory – 72 hours



773656, 772055

Code	Description	Ship Case
382400	Kangaroo ePump with Pole Clamp	1
382491	Kangaroo ePump Power Adapter	1
382492	Kangaroo ePump Pole Clamp	1
1041411	Kangaroo ePump Battery	1

Pump sets

772055	500 mL Pump Set	30
773656	1000 mL Pump Set	30
776255	600 mL Rigid Container Pump Set	30
776250	1200 mL Rigid Container Pump Set	30
674668	500 mL Feed Set with 500 mL Flush Bag	30
773662	1000 mL Feed Set with 1000 mL Flush Bag	30
712054	Sterile 500 mL Pump Set	30
713650	Sterile 1000 mL Pump Set	30
716154	Sterile 100 mL Burette Set	30

Recertification set

776150	100 mL Burette Recertification Set	5
--------	------------------------------------	---

Kangaroo™ Feeding Sets for Safety Screw Spike System

775659	Kangaroo ePump Safety Screw Spike	30
775100	Kangaroo ePump Safety Screw Spike with 1000 mL Flush Bag	30



775659



775100

Enteral Feeding Pumps/Pump Sets (Continued)



383400

Kangaroo™ Joey Enteral Feeding Pump

- Unique ambulatory feeding and flushing/hydrating capabilities
- Attitude Independent
- “Smart” pump set technology
- “ATM-style” ease of operation
- Feeding history up to 72 hours
- DEHP-Free pump sets
- Anti-Free-flow pump sets
- Fully portable design with back-pack
- Compact, lightweight construction
- Features the safe enteral connection system
- New ambulatory-specific feed set design

Specifications

- 18 hour battery
- Weight – 1.7 lbs., with pole clamp – 2.3 lbs
- Feed Rate – 1-400 in 1 mL increments
- Flush Volume – 10-500 in 1 mL increments
- Flush Interval – 1-24 in 1 hr intervals
- Memory – 72 hours
- Multilingual – 15 languages
- Tubing length: 9" (From pump to feeding container).

Code	Description	Ship Case
383400	Kangaroo Joey Enteral Feeding Pump w/ Pole Clamp	1
383493	Kangaroo Joey Pole Clamp	1
383494	Kangaroo Joey Flexible Pole Clamp	1
383491	Kangaroo Joey Power Cord	1
F010506	Kangaroo Joey Battery	1
770026	Kangaroo Joey Blue Mini Backpack	1
770028	Kangaroo Joey Blue Super Mini Backpack	1
770034	Kangaroo Joey Pink Mini Backpack	1
770029	Kangaroo Joey Pink Super Mini Backpack	1
770027	Kangaroo Joey Green Mini Backpack	1
770032	Kangaroo Joey Green Super Mini Backpack	1
770025	Kangaroo Joey Black Mini Backpack	1
770031	Kangaroo Joey Black Super Mini Backpack	1

Pump sets

762055	Kangaroo Joey 500 mL Pump Set	30
763656	Kangaroo Joey 1000 mL Pump Set	30
763662	Kangaroo Joey 1000 mL Pump Set with 1000 mL Flush Bag	30

Recertification set

776150	100 mL Burette Recertification Set	5
--------	------------------------------------	---

Kangaroo™ Feeding Sets for use with NEW Safety Screw Spike System

765559	Kangaroo Joey Safety Screw Spike Set	30
765100	Kangaroo Joey Safety Screw Spike with 1000 mL Flush Bag	30



763662

Enteral Feeding Pumps/Pump Sets (Continued)




773621

Kangaroo™ 924 Pump Sets

Code	Description	Ship Case
Pump sets		
774625	Kangaroo 924 Proximal Spike Set	30
772025	Kangaroo 924 500 mL Pump Set	30
773621	Kangaroo 924 1000 mL Pump Set (Individual Pack)	30
703621	Kangaroo 924 1000 mL Pump Set (Bulk Pack)	30
Kangaroo™ Feeding Sets for use with NEW Safety Screw Spike System		
775759	Kangaroo 924 Safety Screw Spike	30

Safety Screw Spike Compatible Feeding Sets



765559

Kangaroo™ Feeding Sets for use with Safety Screw Spike

- Covidien provides a safety spike compatible with the Nestlé SpikeRight system
- The new enteral-specific SpikeRight™* port is threaded and connects with screw-top enteral spike sets

Code	Description	SHIP Case
765559	Kangaroo Joey Safety Screw Spike Set	30
765100	Kangaroo Joey Safety Screw Spike with 1000 mL Flush Bag	30
775659	Kangaroo ePump Safety Screw Spike Set	30
775100	Kangaroo ePump Safety Screw Spike with 1000 mL Flush Bag	30
775759	Kangaroo 924 Safety Screw Spike Set	30

Pump Sets for use with all Kangaroo™ 224, Kangaroo™ 324, Kangaroo™ PET and Kangaroo™ Control Enteral Feeding Pumps



773630

Kangaroo™ DEHP-Free Pump Sets

- New PVC material eliminates the use of the compound/chemical DEHP

Code	Type	Description	Ship Case
773630	Non-sterile	1000 mL DEHP-Free Pump Set	30
713006	Sterile	1000 mL DEHP-Free Pump Set	36
702035	Non-sterile	500 mL DEHP-Free Pump Set	30
713100	Sterile	DEHP-Free 100 mL Burette Pump Set	36
702520	Non-sterile	1000 mL DEHP-Free Gravity Set	30



773620

Kangaroo™ DEHP-Free Pump Sets with Anti-Free-Flow (AFF)

- New PVC material eliminates the use of the compound/chemical DEHP
- Incorporates the internal Anti-Free-Flow (AFF) device for continuous protection for the patient from the possibility of a free flow incident

Code	Type	Description	Ship Case
773620	Non-sterile	1000 mL DEHP-Free with AFF Pump Set	30
702025	Non-sterile	500 mL DEHP-Free with AFF Pump Set	30



773610

Kangaroo™ Anti-Free-Flow Pump Sets

- Anti-Free-Flow (AFF) device protects the patient from the possible free flow of formula that may result from the pump set becoming dislodged or loaded improperly
- Unique internal design ensures the safety mechanism may not be tampered with by either the patient or caregiver, providing continuous protection

Code	Type	Description	Ship Case
773610	Non-sterile	1000 mL AFF Pump Set	30
702015	Non-sterile	500 mL AFF Pump Set	30

Pump Sets for use with all Kangaroo™ 224, Kangaroo™ 324, Kangaroo™ PET and Kangaroo™ Control Enteral Feeding Pumps (Continued)



8884773600

Kangaroo™ Pump Sets

- Overinfusion safeguard (when used with Kangaroo 224, Kangaroo 324, Kangaroo PET, or Kangaroo Control Pumps)
- Non-IV compatible connectors

Code	Type	Description	SHIP Case
8884712004	Sterile	500 mL Bag with Easy-Cap Closure and Ice Pouch	36
8884702005	Non-sterile	500 mL Bag with Easy-Cap Closure	30
8884702105	Non-sterile	Deluxe 500 mL Graduated Bag with Easy-Cap Closure	30
8884713600	Sterile	1000 mL Bag with Easy-Cap Closure and Ice Pouch	36
8884773600	Non-sterile	1000 mL Bag with Easy-Cap Closure and Ice Pouch	30
8884712100	Sterile	Deluxe 500 mL Graduated Bag with Easy-Cap Closure	36
8884714300	Sterile	Deluxe 1200 mL Graduated Bag with Easy-Cap Closure	36
8884704305	Non-sterile	Easy-Cap 1200 mL	30
8884714800	Sterile	Deluxe 1600 mL Graduated Bag with Easy-Cap Closure	36
8884704805	Non-sterile	1600 mL Bag with Easy-Cap Closure	30
8884704705	Non-sterile	Screw Cap Pump Set — Screw Cap Designed to Fit Glass Enteral Formula Bottles	30
8884706204	Non-sterile	1200 mL Rigid Container with Preattached Pump Set	30
8884706212	Non-sterile	1200 mL Rigid Container without Pump Set	30
8884706220	Non-sterile	600 mL Rigid Container with Preattached Pump Set	30
8884706238	Non-sterile	600 mL Rigid Container without Pump Set	30
8884716100	Sterile	100 mL Rigid Burette	36



8884706204

Kangaroo™ Pump Set Accessories

8884705008	Non-sterile	Y-site Extension Set	30
8884700207	Non-sterile	Pump Set 4 ft. Tubing Extension	30

Kangaroo™ Gravity Sets



702520

Kangaroo™ Gravity Sets

Code	Type	Description	Ship Case
702505	Non-sterile	1000 mL Large Bore Gravity Set	30
702520	Non-sterile	1000 mL DEHP-Free, Large Easy-Cap Opening	30
8884702500	Non-sterile	1000 mL Gravity Bag	30
8884712300	Sterile	1000 mL Gravity Bag	24

Kangaroo™ I.V. Pole Pouch and Kits




84064

Kangaroo™ I.V. Pole Pouch and Kits

- Convenient tear-to-open packaging, resealable flap
- Transparent, polyethylene bag hangs easily, clearly labeled
- Non-sterile

Code	Description	Ship Case
84064	I.V. Pole Pouch, includes 60 cc Flat-Top Syringe	30
84065	I.V. Pole Pouch, includes 60 cc Piston Syringe	30
84066	I.V. Pole Kit, includes 60 cc Piston Syringe w/500 cc Container	30
8884700108	60 cc Piston Syringe w/ 500 cc Container and Basin	20
8884700116	60 cc Piston Syringe w/ 500 cc Container	30
8884700124	60 cc Piston Syringe	30

Nasogastric Feeding Tubes — Neonatal and Pediatric Feeding Systems



461701

Kangaroo™ Purple Polyurethane Feeding Tubes with Safe Enteral Connections

- Enteral connector engineered to reduce risk of tubing misconnections
- Addresses 2006 JCAHO Tubing Misconnections Sentinel Event alert[†]
- Purple color to denote enteral feeding — oral tip syringe only
- Transparent tube with purple radiopaque Sentinel Line™ with Sentinel Eye™ tubing
- Sterile
- DEHP-free

Code	Markings	French Size	Length	Ship Case
461355	From 5 cm to 25 cm in 1 cm increments	3.5 (1.2 mm)	12" (31 cm)	10
461412		5 (1.7 mm)	20" (51 cm)	10
461503		5 (1.7 mm)	36" (91 cm)	10
461420		6.5 (2.2 mm)	20" (51 cm)	10
461438		6.5 (2.2 mm)	36" (91 cm)	10
461701		8 (2.7 mm)	20" (51 cm)	10
461800		8 (2.7 mm)	42" (107 cm)	10
461042		10 (3.3 mm)	42" (107 cm)	10



460802

Kangaroo™ Purple PVC Neonatal and Pediatric Feeding Tubes with Safe Enteral Connections

Code	French Size	Length	Ship Case
461206	3.5 (1.2 mm)	12" (31 cm)	10
460802	5 (1.7 mm)	16" (41 cm)	50
461008	5 (1.7 mm)	36" (91 cm)	50
460711	6.5 (2.2 mm)	16" (41 cm)	50
460604	8 (2.7 mm)	16" (41 cm)	50
460406	8 (2.7 mm)	42" (107 cm)	50
460208	10 (3.3 mm)	42" (107 cm)	50



35ES

Kangaroo™ Feeding Tube Extension Sets

- PVC tubing
- With or without bifurcation (med port)
- DEHP-free
- Latex free
- Accepts oral tip syringe

Code	Description	Ship Case
35ES	35" Extension Set	50
60ES	60" Extension Set	50
35EY	35" Bifurcated Extension Set	50
60EY	60" Bifurcated Extension Set	50

Nasogastric Feeding Tubes — Neonatal and Pediatric Feeding Systems (Continued)



Monoject™ Purple Oral Syringes

Code Number	Description	Case Qty
Bulk		
401S	1 mL Non-sterile Purple Oral Syringe	500
403S	3 mL Non-sterile Purple Oral Syringe	500
406S	6 mL Non-sterile Purple Oral Syringe	500
410S	10 mL Non-sterile Purple Oral Syringe	500
Sterile		
401SG	1 mL Sterile Purple Oral Syringe	240
403SG	3 mL Sterile Purple Oral Syringe	480
406SG	6 mL Sterile Purple Oral Syringe	400
412SG	12 mL Sterile Purple Oral Syringe	192
435SG	35 mL Sterile Purple Oral Syringe	160
460SG	60 mL Sterile Purple Oral Syringe	120

Nasogastric Feeding Tubes — Neonatal and Pediatric Feeding Systems (Continued)




8888261701

Kangaroo™ Polyurethane Feeding Tubes

- Polyurethane material remains flexible for prolonged use
- Transparent
- Colored radiopaque Sentinel Line™ with Sentinel Eye™ tubing
- Proximal end accepts male luer connector
- Attached cap to cover proximal end
- Sterile

Code	Markings	French Size	Length	Ship Case
8888261355	From 5 cm to 25 cm in 1 cm increments	3.5 (1.2 mm)	12" (31 cm)	10
8888261412		5 (1.7 mm)	20" (51 cm)	10
8888261503		5 (1.7 mm)	36" (91 cm)	10
8888261420		6.5 (2.2 mm)	20" (51 cm)	10
8888261438		6.5 (2.2 mm)	36" (91 cm)	10
8888261701		8 (2.7 mm)	20" (51 cm)	10
8888261800		8 (2.7 mm)	42" (107 cm)	10
8888261042		10 (33 mm)	42" (107 cm)	10



8888260604

Kangaroo™ Polyvinyl Chloride Feeding Tubes with Radiopaque Line

- Thermosensitive material with Slide-TEX finish
- Colored radiopaque Sentinel Line™ tubing
- Two eyes
- Proximal end accepts male Luer connector
- Funnel connector
- Attached cap to cover proximal end
- Sterile

Code	Markings	French Size	Length	Ship Case
8888261206	From 5 cm to 25 cm in 1 cm increments	3.5 (1.2 mm)	12" (31 cm)	10
8888260802		5 (1.7 mm)	16" (41 cm)	50
8888261008		5 (1.7 mm)	36" (91 cm)	50
8888260711		6.5 (2.2 mm)	16" (41 cm)	50
8888260604		8 (2.7 mm)	16" (41 cm)	50
8888260406		8 (2.7 mm)	42" (107 cm)	50
8888260208		10 (3.3 mm)	42" (107 cm)	50

Curity™ PVC Feeding Tubes

155720*	8" (20 cm)	5 (1.7 mm)	15" (38 cm)	50
155721*	8" (20 cm)	5 (1.7 mm)	36" (91 cm)	50
155722*	8" (20 cm)	8 (2.7 mm)	15" (38 cm)	50
155723*	8" (20 cm)	8 (2.7 mm)	42" (107 cm)	50

* Does not contain radiopaque Sentinel Line or cm markings.

Nasogastric Feeding Tubes — Adult



Kangaroo™ Feeding Tubes with Dobbhoff™ Tip

- Dual port enteral connector engineered to reduce risk of tubing misconnections
- Addresses 2006 JCAHO Tubing Misconnections Sentinel Event alert[†]
- Purple color to denote enteral feeding
- Weighted feeding tubes with and without flow-through stylets
- Radiopaque polyurethane material
- CM markings
- MRI conditional once stylet is removed
- DEHP-Free

Code	Includes Stylet	Rigid Outlet Port	Weight	French Size	Length	SHIP Case
8884710867	Yes	Yes	5g	8	43" (109 cm)	10
8884710842	No	Yes	7g	8	43" (109 cm)	10
8884710859	Yes	Yes	7g	8	43" (109 cm)	10
8884710826	Yes	Yes	7g	8	55" (140 cm)	10
8884711006	Yes	Yes	7g	10	43" (109 cm)	10
8884711246	No	No	7g	12	43" (109 cm)	10
8884711253	Yes	Yes	7g	12	43" (109 cm)	10

Nasogastric Feeding Tubes — Adult (Continued)



8884721211

Kangaroo™ Feeding Tubes with Standard Tip

- Dual port enteral connector engineered to reduce risk of tubing misconnections
- Addresses 2006 JCAHO Tubing Misconnections Sentinel Event alert[†]
- Purple color to denote enteral feeding — oral/cath tip syringe only
- Weighted feeding tubes with and without flow-through stylets
- Radiopaque polyurethane material
- CM markings
- MRI conditional once stylet is removed
- DEHP-Free

Code	Includes Stylet	Rigid Outlet Port	Weight	French Size	Length	Ship Case
8884720817	No	Yes	3g	8	36" (91 cm)	10
8884720841	No	Yes	3g	8	43" (109 cm)	10
8884720825	Yes	Yes	3g	8	36" (91 cm)	10
8884720858	Yes	Yes	3g	8	43" (109 cm)	10
8884720855	Yes	Yes	5g	8	55" (140 cm)	10
8884721088	Yes	Yes	5g	10	43" (109 cm)	10
8884721096	No	Yes	5g	10	43" (109 cm)	10
8884721055	Yes	Yes	5g	10	55" (140 cm)	10
8884721211	No	No	5g	12	36" (91 cm)	10
8884721237	No	No	5g	12	43" (109 cm)	10
8884721252	Yes	Yes	5g	12	43" (109 cm)	10
8884721255	Yes	Yes	5g	12	55" (140 cm)	10



8884711501

Kangaroo™ Non-Weighted Feeding Tubes


- Clear polyurethane material
- Proximal Y-site connector with catheter tip adapter
- Large exit holes and open distal end
- Radiopaque stripe
- DEHP-Free

Code	French Size	Length	Ship Case
8884711568	8	36" (91 cm)	10
8884711550	10	36" (91 cm)	10
8884711501	12	36" (91 cm)	10
8884711519	14	36" (91 cm)	10

CO₂ Detector



777702




7771610CN

Kangaroo™ CO₂ Detector For Gastric Tube Placement

- Proprietary colorimetric paper technology
- Detects presence of CO₂ within seconds
- Reduces unnecessary waiting time and minimizes the number of X-rays needed for procedure
- Ends the guesswork about proper gastric tube placement
- Convenient bellows device:
 - Clears obstructions from distal tip
 - Prevents CO₂ from entering tube prematurely
 - Facilitates airflow for quick results
- Easily attaches to either a bellows or syringe
- Requires no change in insertion technique
- Universally adapts to all feeding and NG tubes

Code	Description	Ship Case
777702	Kangaroo CO ₂ Detector, with Bellows	10
710859CN	Kangaroo 8Fr; 43" w/stylet, Safe Enteral Connections and Kangaroo CO ₂ Detector	10
711006CN	Kangaroo 10Fr; 43" w/stylet, Safe Enteral Connections and Kangaroo CO ₂ Detector	10
711253CN	Kangaroo 12Fr; 43" w/stylet, Safe Enteral Connections and Kangaroo CO ₂ Detector	10
720858CN	Kangaroo 8Fr; 43" w/stylet, Safe Enteral Connections and Kangaroo CO ₂ Detector	10
721088CN	Kangaroo 10Fr; 43" w/stylet, Safe Enteral Connections and Kangaroo CO ₂ Detector	10
721252CN	Kangaroo 12Fr; 43" w/stylet, Safe Enteral Connections and Kangaroo CO ₂ Detector	10
266130CN	16Fr Salem Sump Tube w/ ARV and Kangaroo CO ₂ Detector	50
266148CN	18Fr Salem Sump Tube w/ ARV and Kangaroo CO ₂ Detector	50
7771610CN	16Fr Salem Sump Tube w/ Multifunctional Port and CO ₂ Detector	10
7771810CN	18Fr Salem Sump Tube w/ Multifunctional Port and CO ₂ Detector	10

Specialty Tubes




Kangaroo™ Naso-Jejunal Feeding and Gastric Decompression Tube

- For "push" insertion over guidewire (kit included)
- Radiopaque polyurethane material
- Unique adjustable length feeding tube design
- Spiral feeding outlet holes
- Jejunal tube can be retracted within gastric portion
- MRI conditional once stylet is removed

Code	Overall Diameter	Jejunal Tube Diameter	Overall Length	Ship Case
8884752513	16Fr	9Fr.	67" (170 cm)	2

Specialty Tubes (Continued)



8884752505

Kangaroo™ Naso-Jejunal Feeding Tubes

- Designed for precision placement in the small bowel using either endoscopy and/or fluoroscopy (includes guidewire and transfer cannula)
- Unique hollow bolus weighted tip
- Radiopaque polyurethane material
- MRI conditional once stylet is removed

Code	Weight	French Size	Length	Ship Case
8884752505	7g	12	60" (152 cm)	2



8884750095

Kangaroo™ Jejunal Feeding and Decompression Tube

- For use with 20Fr EntriStar single pass P.E.G. and Dobbhoff™ safety P.E.G. kits
- Guidewire insertion
- Simultaneous gastric decompression and jejunal feeding through a P.E.G.
- Sterile
- MRI conditional once stylet is removed

Code	French Size of Jejunal Tube	French Size of P.E.G.	Length of Jejunal Tube	Ship Case
8884750095	9	20	35" (89 cm)	2

Safety P.E.G. Systems



761610

Kangaroo™ Safety P.E.G. Kits with Safe Enteral Connections


- Enteral connector engineered to reduce risk of tubing misconnections
- Addresses 2006 JCAHO Tubing Misconnections Sentinel Event alert[†]
- Purple color to denote enteral feeding
- Featuring the Grant 90° bend which straightens easily for P.E.G. placement and insertion of accessory feeding tubes
- Endoscopic removal
- Radiopaque polyurethane material
- Accessory jejunal feeding/gastric decompression tube is available for 20Fr P.E.G.
- Sterile
- DEHP-Free

Code		Size	Procedural Aids	Ship Case
------	--	------	-----------------	-----------

Safety Percutaneous Endoscopic Gastrostomy "Pull" Kits

761610	Kit Contains: 21ga x 1 1/2" Monoject Magellan™ Safety Needle, Dual Port Enteral Connector, 6cc Syringe, Pinch Clamp, #11 Safety Scalpel, Catheter Assembly w/90 Grant Bend, 18ga Safety Introducer Needle, External Retention Band, Hemostat, PVP Iodine Swabstick, Obturator Removal Kit, Gauze Sponges, Looping Wire/Guidewire, Lubrication Jelly Packet, Fenestrated Drape, Disposable Snare	16Fr	Yes	2
762010	Kit Contains: 21ga x 1 1/2" Monoject Magellan™ Safety Needle, Dual Port Enteral Connector, 6cc Syringe, Pinch Clamp, #11 Safety Scalpel, Catheter Assembly w/90 Grant Bend, 18ga Safety Introducer Needle, External Retention Band, Hemostat, PVP Iodine Swabstick, Obturator Removal Kit, Gauze Sponges, Looping Wire/Guidewire, Lubrication Jelly Packet, Fenestrated Drape, Disposable Snare	20Fr	Yes	2

Safety P.E.G. Systems (Continued)



761600

Kangaroo™ Safety Single Pass P.E.G. Kits with Safe Enteral Connections

- Enteral connector engineered to reduce risk of tubing misconnections
- Addresses 2006 JCAHO Tubing Misconnections Sentinel Event alert[†]
- Purple color to denote enteral feeding
- Radiopaque internal retention bolster
- Non-endoscopic removal
- Easy to read placement marks
- Hydromer™* hydrophilic grafting on the “Pull” dilator tip
- Accessory jejunal feeding/gastric decompression tube is available for 20Fr P.E.G.
- Available in both “Pull” and “Push” over the guidewire methods
- Sterile
- DEHP-Free

Code		Size	Procedural Aids	Ship Case
Safety Percutaneous Endoscopic Gastrostomy “Pull” Kits				
761600	<u>Kit Contains:</u> 21ga x 1 1/2" Monoject Magellan™ Safety Needle, Dual Port Enteral Connector, 6cc Syringe, Pinch Clamp, #11 Safety Scalpel, Catheter Assembly, 18ga Safety Introducer Needle, External Retention Band, Hemostat, PVP Iodine Swabstick, Obturator Removal Kit, Gauze Sponges, Looping Wire/Guidewire, Lubrication Jelly Packet, Fenestrated Drape, Disposable Snare	16Fr	Yes	2
762000	<u>Kit Contains:</u> 21ga x 1 1/2" Monoject Magellan™ Safety Needle, Dual Port Enteral Connector, 6cc Syringe, Pinch Clamp, #11 Safety Scalpel, Catheter Assembly, 18ga Safety Introducer Needle, External Retention Band, Hemostat, PVP Iodine Swabstick, Obturator Removal Kit, Gauze Sponges, Looping Wire/Guidewire, Lubrication Jelly Packet, Fenestrated Drape, Disposable Snare	20Fr	Yes	2



8884751622

Kangaroo™ Y-Port-G-Tube/P.E.G. Adapter with Safe Enteral Connections

- Dual port enteral connector engineered to reduce risk of tubing misconnections
- Addresses 2006 JCAHO Tubing Misconnections Sentinel Event alert[†]
- Purple color to denote enteral feeding
- For use with the Kangaroo™ P.E.G. Kits
- The 20Fr Y-port can be used as the replacement Y-port on all Kangaroo™ gastrostomy tubes

Code	Description	Ship Case
8884751622	16Fr Y-port/P.E.G. Adapter, Non-sterile	20
8884752042	20Fr Y-port/P.E.G. Adapter, Non-sterile	20

Gastrostomy and Skin Level Tubes



Kangaroo™ Gastrostomy Feeding Tubes with Y-Ports with Safe Enteral Connections

- Dual port enteral connector engineered to reduce risk of tubing misconnections
- Addresses 2006 JCAHO Tubing Misconnections Sentinel Event alert†
- Purple color to denote enteral feeding — oral/cath tip syringe only
- Medical-grade silicone tubing
- Rounded skin disc for external retention with feet
- Graduated shaft
- Gastric balloon for internal retention
- Rounded flush tip with open distal end
- DEHP-Free

Code	French Size	Balloon Size	Ship Case
8884715122	12Fr	5 cc	5
8884715148	14Fr	5 cc	5
8884720163	16Fr	20 cc	5
8884720189	18Fr	20 cc	5
8884720205	20Fr	20 cc	5
8884720221	22Fr	20 cc	5
8884720247	24Fr	20 cc	5
8884720265	26Fr	20 cc	5
8884720285	28Fr	20 cc	5

Gastrostomy and Skin Level Tubes (Continued)



Kangaroo™ Skin Level Balloon Gastrostomy Kits with Safe Enteral Connections

- Enteral connector engineered to reduce risk of tubing misconnections
- Addresses 2006 JCAHO Tubing Misconnections Sentinel Event alert[†]
- Purple color to denote enteral feeding
- 100% silicone material
- Durable balloon construction
- Raised feet design
- ClipStar™ enteral security clip
- Large range of sizes available
- DEHP-Free

Contents of NutriPort Kit:

Skin Level Gastrostomy Device, Step-Up Shim, Bolus Feeding Set, Continuous Feeding Set with ClipStar Enteral Security Clip, 60 cc Irrigation Syringe, 12cc Syringe, Gauze Pads

Code	Description	Ship Carton	Code	Description	Ship Carton
712080	Balloon SLG 12Fr x 0.8 cm	1	714080	Balloon SLG 14Fr x 0.8 cm	1
712100	Balloon SLG 12Fr x 1.0 cm	1	714100	Balloon SLG 14Fr x 1.0 cm	1
712120	Balloon SLG 12Fr x 1.2 cm	1	714120	Balloon SLG 14Fr x 1.2 cm	1
712150	Balloon SLG 12Fr x 1.5 cm	1	714150	Balloon SLG 14Fr x 1.5 cm	1
712170	Balloon SLG 12Fr x 1.7 cm	1	714170	Balloon SLG 14Fr x 1.7 cm	1
712200	Balloon SLG 12Fr x 2.0 cm	1	714200	Balloon SLG 14Fr x 2.0 cm	1
712230	Balloon SLG 12Fr x 2.3 cm	1	714230	Balloon SLG 14Fr x 2.3 cm	1
712250	Balloon SLG 12Fr x 2.5 cm	1	714250	Balloon SLG 14Fr x 2.5 cm	1
712270	Balloon SLG 12Fr x 2.7 cm	1	714270	Balloon SLG 14Fr x 2.7 cm	1
712300	Balloon SLG 12Fr x 3.0 cm	1	714300	Balloon SLG 14Fr x 3.0 cm	1
712350	Balloon SLG 12Fr x 3.5 cm	1	714350	Balloon SLG 14Fr x 3.5 cm	1
712400	Balloon SLG 12Fr x 4.0 cm	1	714400	Balloon SLG 14Fr x 4.0 cm	1
712450	Balloon SLG 12Fr x 4.5 cm	1	714450	Balloon SLG 14Fr x 4.5 cm	1
712500	Balloon SLG 12Fr x 5.0 cm	1	714500	Balloon SLG 14Fr x 5.0 cm	1


Gastrostomy and Skin Level Tubes (Continued)

Kangaroo™ Skin Level Balloon Gastrostomy Kits (continued)

Code	Description	Ship Carton	Code	Description	Ship Carton
716080	Balloon SLG 16Fr x 0.8 cm	1	720080	Balloon SLG 20Fr x 0.8 cm	1
716100	Balloon SLG 16Fr x 1.0 cm	1	720100	Balloon SLG 20Fr x 1.0 cm	1
716120	Balloon SLG 16Fr x 1.2 cm	1	720120	Balloon SLG 20Fr x 1.2 cm	1
716150	Balloon SLG 16Fr x 1.5 cm	1	720150	Balloon SLG 20Fr x 1.5 cm	1
716170	Balloon SLG 16Fr x 1.7 cm	1	720170	Balloon SLG 20Fr x 1.7 cm	1
716200	Balloon SLG 16Fr x 2.0 cm	1	720200	Balloon SLG 20Fr x 2.0 cm	1
716230	Balloon SLG 16Fr x 2.3 cm	1	720230	Balloon SLG 20Fr x 2.3 cm	1
716250	Balloon SLG 16Fr x 2.5 cm	1	720250	Balloon SLG 20Fr x 2.5 cm	1
716270	Balloon SLG 16Fr x 2.7 cm	1	720270	Balloon SLG 20Fr x 2.7 cm	1
716300	Balloon SLG 16Fr x 3.0 cm	1	720300	Balloon SLG 20Fr x 3.0 cm	1
716350	Balloon SLG 16Fr x 3.5 cm	1	720350	Balloon SLG 20Fr x 3.5 cm	1
716400	Balloon SLG 16Fr x 4.0 cm	1	720400	Balloon SLG 20Fr x 4.0 cm	1
716450	Balloon SLG 16Fr x 4.5 cm	1	720450	Balloon SLG 20Fr x 4.5 cm	1
716500	Balloon SLG 16Fr x 5.0 cm	1	720500	Balloon SLG 20Fr x 5.0 cm	1
718080	Balloon SLG 18Fr x 0.8 cm	1	724080	Balloon SLG 24Fr x 0.8 cm	1
718100	Balloon SLG 18Fr x 1.0 cm	1	724100	Balloon SLG 24Fr x 1.0 cm	1
718120	Balloon SLG 18Fr x 1.2 cm	1	724120	Balloon SLG 24Fr x 1.2 cm	1
718150	Balloon SLG 18Fr x 1.5 cm	1	724150	Balloon SLG 24Fr x 1.5 cm	1
718170	Balloon SLG 18Fr x 1.7 cm	1	724170	Balloon SLG 24Fr x 1.7 cm	1
718200	Balloon SLG 18Fr x 2.0 cm	1	724200	Balloon SLG 24Fr x 2.0 cm	1
718230	Balloon SLG 18Fr x 2.3 cm	1	724230	Balloon SLG 24Fr x 2.3 cm	1
718250	Balloon SLG 18Fr x 2.5 cm	1	724250	Balloon SLG 24Fr x 2.5 cm	1
718270	Balloon SLG 18Fr x 2.7 cm	1	724270	Balloon SLG 24Fr x 2.7 cm	1
718300	Balloon SLG 18Fr x 3.0 cm	1	724300	Balloon SLG 24Fr x 3.0 cm	1
718350	Balloon SLG 18Fr x 3.5 cm	1	724350	Balloon SLG 24Fr x 3.5 cm	1
718400	Balloon SLG 18Fr x 4.0 cm	1	724400	Balloon SLG 24Fr x 4.0 cm	1
718450	Balloon SLG 18Fr x 4.5 cm	1	724450	Balloon SLG 24Fr x 4.5 cm	1
718500	Balloon SLG 18Fr x 5.0 cm	1	724500	Balloon SLG 24Fr x 5.0 cm	1

Accessories

Code	Description	Ship Case
8884741805	Stoma Measuring Device, Sterile	5
8884741813	Bolus Feeding Set, Non-sterile	10
8884741821	Continuous Feeding Set, Non-sterile	10



8884741805

Salem Sump™ Tubes



Kangaroo™ Multifunctional Port with Safe Enteral Connections

- Addresses 2006 JCAHO Tubing Misconnections Sentinel Event alert†
- Purple color to denote enteral feeding
- All in one closed port system
- Allows for feeding, suction, irrigation and medication
- Internal ARV design
- DEHP-Free
- Reduce clinician exposure
- Improve clinician efficacy
- Radiopaque tubing

Code	Description	Ship Case
7771111	Kangaroo Multifunctional Port with Single Lumen Adaptor	10

Salem Sump

7771010	10Fr Salem Sump with Kangaroo Multifunctional Port	10
7771210	12Fr Salem Sump with Kangaroo Multifunctional Port	10
7771410	14Fr Salem Sump with Kangaroo Multifunctional Port	10
7771610	16Fr Salem Sump with Kangaroo Multifunctional Port	10
7771810	18Fr Salem Sump with Kangaroo Multifunctional Port	10

Silicone Salem Sump w/Anti-Reflux Valve Irrigation

7771011	10Fr Silicone Salem Sump w/Anti-Reflux Valve Irrigation Kangaroo Multifunctional Port	10
7771211	12Fr Silicone Salem Sump w/Anti-Reflux Valve Irrigation Kangaroo Multifunctional Port	10
7771411	14Fr Silicone Salem Sump w/Anti-Reflux Valve Irrigation Kangaroo Multifunctional Port	10
7771611	16Fr Silicone Salem Sump w/Anti-Reflux Valve Irrigation Kangaroo Multifunctional Port	10
7771811	18Fr Silicone Salem Sump w/Anti-Reflux Valve Irrigation Kangaroo Multifunctional Port	10

Salem Sump™ Tubes (Continued)



8888266106

Salem Sump™ Tube with Anti-Reflux Valve

- New cm markings allow clinician to easily determine tube depth during placement
- Gastric reflux prevention valve
- Double lumen PVC tube – one for suction drainage; one for sump vent
- Radiopaque Sentinel Line™ with Sentinel Eye™ tubing
- Integral irrigation funnel for irrigating through sump vent lumen
- Funnel end of vent lumen can be used to cap 5-in-1 adapter
- Sterile

Code	French Size	Length	Ship Case
8888266106	10 (3.3 mm)	36" (91 cm)	50
8888266114	12 (4.0 mm)	48" (122 cm)	50
8888266122	14 (4.7 mm)	48" (122 cm)	50
8888266130	16 (5.3 mm)	48" (122 cm)	50
8888266148	18 (6.0 mm)	48" (122 cm)	50




8888266197

Salem Sump™ Anti-Reflux Valve

- Gastric reflux prevention valve
- Allows for maximum air flow to assist with stomach ulcer prevention
- Sterile

Code	Description	Ship Case
8888266197	Valve only	10



8888268060

Salem Sump™ Tube

- New cm markings allow clinician to easily determine tube depth during placement
- Double lumen PVC tube – one for suction drainage; one for sump vent
- Radiopaque Sentinel Line™ with Sentinel Eye™ tubing
- Integral irrigation funnel for irrigating through sump vent lumen
- Funnel end of vent lumen can be used to cap 5-in-1 adapter
- Sterile

Code	French Size	Length	Ship Case
8888268060	6 (2.0 mm)	24" (61 cm)	10
8888268086	8 (2.7 mm)	24" (61 cm)	10
8888264911	10 (3.3 mm)	36" (91 cm)	50
8888264929	12 (4.0 mm)	48" (122 cm)	50
8888264945	14 (4.7 mm)	48" (122 cm)	50
8888264960	16 (5.3 mm)	48" (122 cm)	50
8888264986	18 (6.0 mm)	48" (122 cm)	50

Salem Sump™ Tubes (Continued)



8888266213

Salem Sump™ Silicone Tube with Anti-Reflux Valve

- Gastric reflux prevention valve
- Double lumen silicone tube – one for suction drainage; one for sump vent
- Radiopaque Sentinel Line™ with Sentinel Eye™ tubing
- Removable 5-in-1 connector
- Capping system with sump vent lumen cap, suction lumen adapter and suction lumen adapter cap
- Sump vent lumen can also be used for irrigation
- Sterile

Code	French Size	Length	Ship Case
8888266213	10 (3.3 mm)	36" (91 cm)	10
8888266221	12 (4.0 mm)	48" (122 cm)	10
8888266247	14 (4.7 mm)	48" (122 cm)	10
8888266262	16 (5.3 mm)	48" (122 cm)	10
8888266288	18 (6.0 mm)	48" (122 cm)	10



8888265108

Salem Sump™ Silicone Tubes

- Double lumen silicone tube – one for suction drainage; one for venting action
- Radiopaque Sentinel Line™ with Sentinel Eye™ tubing
- Removable 5-in-1 connector
- Capping system with sump vent lumen cap, suction lumen adapter and suction lumen adapter cap
- Sump vent lumen can also be used for irrigation
- Sterile

Code	French Size	Length	Ship Case
8888265108	10 (3.3 mm)	36" (91 cm)	10
8888265116	12 (4.0 mm)	48" (122 cm)	10
8888265124	14 (4.7 mm)	48" (122 cm)	10
8888265132	16 (5.3 mm)	48" (122 cm)	10
8888265140	18 (6.0 mm)	48" (122 cm)	10

Stomach Tubes



8888265314

Argyle™ Stomach Tubes (Levin Type)

- Clear, thermosensitive PVC tube with four eyes
- Radiopaque Sentinel Line™ with Sentinel Eye™ tubing
- Open distal end
- Integral funnel connector
- Markings 18" (46 cm), 22" (56 cm), 26" (66 cm) and 30" (76 cm) From distal end
- Non-sterile

Code	French Size	Length	Ship Case
8888265314	12 (4.0 mm)	48" (122 cm)	50
8888265322	16 (5.3 mm)	48" (122 cm)	50
8888265330	18 (6.0 mm)	48" (122 cm)	50



155712

Stomach Tubes (Levin Type)

- Tube has four staggered eyes and center opening to reduce trauma
- Depth marks begin 20" from distal tip with three additional marks at 4" intervals
- Smooth, rounded tip reduces insertion trauma
- A variety of French sizes, 48" long
- Non-sterile

Code	French Size	Ship Case
155710	12Fr	50
155711	14Fr	50
155712	16Fr	50
155713	18Fr	50



8888750018

Argyle™ Edlich Gastric Lavage Trays

- Large 140 cc (5 oz.) graduated Monoject™ syringe
- Transparent PVC gastric tube with closed end and 4 eyes
- Depth marks at 18" (46 cm), 22" (56 cm), 26" (66 cm) and 30" (76 cm)
- Clamp

Code	French Size	Length	Ship Case
8888750018	Large, 34Fr (11.3 mm)	36" (91 cm)	20
8888750513	Medium/Pediatric, 24Fr (8.0 mm)	36" (91 cm)	20

INDEX

Code	Product	Page	Code	Product	Page
1041411	Kangaroo ePump Battery	1	461420	Kangaroo Purple Polyurethane Feeding Tubes	7
155710	Stomach Tubes, PVC (Levin Type)	21	461438	Kangaroo Purple Polyurethane Feeding Tubes	7
155711	Stomach Tubes, PVC (Levin Type)	21	461503	Kangaroo Purple Polyurethane Feeding Tubes	7
155712	Stomach Tubes, PVC (Levin Type)	21	461701	Kangaroo Purple Polyurethane Feeding Tubes	7
155713	Stomach Tubes, PVC (Levin Type)	21	461800	Kangaroo Purple Polyurethane Feeding Tubes	7
155720	Curity PVC Feeding Tubes.....	9	60ES	Kangaroo Feeding Tube Bifurcated Extension Set	7
155721	Curity PVC Feeding Tubes.....	9	60EY	Kangaroo Feeding Tube Bifurcated Extension Set	7
155722	Curity PVC Feeding Tubes.....	9	674668	Kangaroo ePump 500 mL Pump Set, 500 mL Flush Bag.....	1
155723	Curity PVC Feeding Tubes.....	9	702015	Kangaroo 500 mL AFF Pump Set	4
266130CN	Salem Sump Tube w/ ARV and CO ₂ Detector.....	12	702025	Kangaroo 500 mL DEHP-Free with AFF Pump Set.....	4
266148CN	Salem Sump Tube w/ ARV and CO ₂ Detector.....	12	702035	Kangaroo 500 mL DEHP-Free Pump Set	4
35ES	Kangaroo Feeding Tube Extension Set.....	7	702505	Kangaroo 1000 ml Large Bore Gravity Set.....	6
35EY	Kangaroo Feeding Tube Extension Set.....	7	702520	Kangaroo 1000 ml DEHP Free, Large Opening.....	6
382400	Kangaroo ePump with Pole Clamp	1	702520	Kangaroo 1000 ml DEHP-Free Gravity Set	4
382491	Kangaroo ePump Power Adapter.....	1	703621	Kangaroo 924 1000 mL Pump Set (Bulk Pack)	3
382492	Kangaroo ePump Pole Clamp	1	704645	Kangaroo DEHP-Free Proximal Spike Set.....	4
383400	Kangaroo Joey Enteral Feeding Pump.....	2	710859CN	Kangaroo, Safe Enteral Connections & CO ₂ Detector	12
383491	Kangaroo Joey Power Cord.....	2	711006CN	Kangaroo, Safe Enteral Connections & CO ₂ Detector	12
383493	Kangaroo Joey Pole Clamp	2	711253CN	Kangaroo, Safe Enteral Connections & CO ₂ Detector	12
383494	Kangaroo Joey Flexible Pole Clamp.....	2	712054	Kangaroo ePump Sterile 500 mL Pump Set.....	1
401S	Monoject Purple Oral Syringes.....	8	712080	Kangaroo Skin Level Balloon Gastrostomy Kit.....	16
401SG	Monoject Purple Oral Syringes.....	8	712100	Kangaroo Skin Level Balloon Gastrostomy Kit.....	16
403S	Monoject Purple Oral Syringes.....	8	712120	Kangaroo Skin Level Balloon Gastrostomy Kit.....	16
403SG	Monoject Purple Oral Syringes.....	8	712150	Kangaroo Skin Level Balloon Gastrostomy Kit.....	16
406S	Monoject Purple Oral Syringes.....	8	712170	Kangaroo Skin Level Balloon Gastrostomy Kit.....	16
406SG	Monoject Purple Oral Syringes.....	8	712200	Kangaroo Skin Level Balloon Gastrostomy Kit.....	16
410S	Monoject Purple Oral Syringes.....	8	712230	Kangaroo Skin Level Balloon Gastrostomy Kit.....	16
412SG	Monoject Purple Oral Syringes.....	8	712250	Kangaroo Skin Level Balloon Gastrostomy Kit.....	16
435SG	Monoject Purple Oral Syringes.....	8	712270	Kangaroo Skin Level Balloon Gastrostomy Kit.....	16
460208	Kangaroo Purple PVC Feeding Tubes	7	712300	Kangaroo Skin Level Balloon Gastrostomy Kit.....	16
460406	Kangaroo Purple PVC Feeding Tubes	7	712350	Kangaroo Skin Level Balloon Gastrostomy Kit.....	16
460604	Kangaroo Purple PVC Feeding Tubes	7	712400	Kangaroo Skin Level Balloon Gastrostomy Kit.....	16
460711	Kangaroo Purple PVC Feeding Tubes	7	712450	Kangaroo Skin Level Balloon Gastrostomy Kit.....	16
460802	Kangaroo Purple PVC Feeding Tubes	7	712500	Kangaroo Skin Level Balloon Gastrostomy Kit.....	16
460SG	Monoject Purple Oral Syringes.....	8	713006	Kangaroo Sterile 1000 mL DEHP-Free Pump Set	4
461008	Kangaroo Purple PVC Feeding Tubes	7	713100	Kangaroo Sterile DEHP-Free Burette Pump Set	4
461042	Kangaroo Purple Polyurethane Feeding Tubes	7	713650	Kangaroo ePump Sterile 1000 mL Pump Set.....	1
461206	Kangaroo Purple PVC Feeding Tubes	7	714080	Kangaroo Skin Level Balloon Gastrostomy Kit.....	16
461355	Kangaroo Purple Polyurethane Feeding Tubes	7	714100	Kangaroo Skin Level Balloon Gastrostomy Kit.....	16

INDEX

Code	Product	Page	Code	Product	Page
714120	Kangaroo Skin Level Balloon Gastrostomy Kit.....	16	718450	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17
714150	Kangaroo Skin Level Balloon Gastrostomy Kit.....	16	718500	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17
714170	Kangaroo Skin Level Balloon Gastrostomy Kit.....	16	720080	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17
714200	Kangaroo Skin Level Balloon Gastrostomy Kit.....	16	720100	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17
714230	Kangaroo Skin Level Balloon Gastrostomy Kit.....	16	720120	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17
714250	Kangaroo Skin Level Balloon Gastrostomy Kit.....	16	720150	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17
714270	Kangaroo Skin Level Balloon Gastrostomy Kit.....	16	720170	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17
714300	Kangaroo Skin Level Balloon Gastrostomy Kit.....	16	720200	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17
714350	Kangaroo Skin Level Balloon Gastrostomy Kit.....	16	720230	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17
714400	Kangaroo Skin Level Balloon Gastrostomy Kit.....	16	720250	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17
714450	Kangaroo Skin Level Balloon Gastrostomy Kit.....	16	720270	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17
714500	Kangaroo Skin Level Balloon Gastrostomy Kit.....	16	720300	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17
716080	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17	720350	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17
716100	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17	720400	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17
716120	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17	720450	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17
716150	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17	720500	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17
716154	Kangaroo ePump Sterile 100 mL Burette Set	1	720858CN	Kangaroo, Safe Enteral Connections & CO ₂ Detector	12
716170	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17	721088CN	Kangaroo, Safe Enteral Connections & CO ₂ Detector	12
716200	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17	721252CN	Kangaroo, Safe Enteral Connections & CO ₂ Detector	12
716230	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17	724080	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17
716250	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17	724100	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17
716270	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17	724120	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17
716300	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17	724150	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17
716350	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17	724170	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17
716400	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17	724200	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17
716450	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17	724230	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17
716500	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17	724250	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17
718080	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17	724270	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17
718100	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17	724300	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17
718120	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17	724350	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17
718150	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17	724400	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17
718170	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17	724450	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17
718200	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17	724500	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17
718230	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17	761600	Kangaroo Safety Single Pass P.E.G. "Pull" Kit.....	14
718250	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17	761610	Kangaroo Safety P.E.G. Kit.....	13
718270	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17	762000	Kangaroo Safety Single Pass P.E.G. "Pull" Kit.....	14
718300	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17	762010	Kangaroo Safety P.E.G. Kit.....	13
718350	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17	762055	Kangaroo Joey 500 mL Pump Set	2
718400	Kangaroo Skin Level Balloon Gastrostomy Kit.....	17	763656	Kangaroo Joey 1000 mL Pump Set	2

INDEX

Code	Product	Page	Code	Product	Page
763662	Kangaroo Joey 1000 mL Pump Set, 1000 mL Flush Bag.....	2	84065	Kangaroo I.V. Pole Pouch, 60 cc Piston Syringe.....	6
765100	Kangaroo Joey Safety Screw Spike, Flush bag.....	2, 3	84066	Kangaroo I.V. Pole Kit; 60 cc Piston Syringe; Container.....	6
765559	Kangaroo Joey Safety Screw Spike Set.....	2, 3	8884700108	60 cc Piston Syringe, 500 cc Container and Basin.....	6
770025	Kangaroo Joey Backpack.....	2	8884700116	60 cc Piston Syringe, 500 cc Container	6
770026	Kangaroo Joey Backpack.....	2	8884700124	60cc Piston Syringe	6
770027	Kangaroo Joey Backpack.....	2	8884700207	Kangaroo Pump Set.....	5
770028	Kangaroo Joey Backpack.....	2	8884702005	Kangaroo 500 ml Bag with Easy-Cap Closure.....	5
770029	Kangaroo Joey Backpack.....	2	8884702105	Kangaroo Deluxe 500 ml Graduated Bag.....	5
770031	Kangaroo Joey Backpack.....	2	8884702500	Kangaroo 1000 mL Gravity Bag.....	6
770032	Kangaroo Joey Backpack.....	2	8884704305	Kangaroo Easy-Cap 1200 mL	5
770034	Kangaroo Joey Backpack.....	2	8884704705	Kangaroo Screw Cap Pump Set for Glass Bottles	5
772025	Kangaroo 924 500 mL Pump Set.....	3	8884704805	Kangaroo 1600 mL Bag with Easy-Cap Closure	5
772055	Kangaroo Joey 500 mL Pump Set	1	8884705008	Kangaroo Y-site Extension Set.....	5
773610	Kangaroo 1000 mL AFF Pump Set	4	8884706204	Kangaroo 1200 mL Rigid Container, Pump Set	5
773620	Kangaroo 1000 mL DEHP-Free, AFF Pump Set.....	4	8884706212	Kangaroo 1200 mL Rigid Container	5
773621	Kangaroo 924 1000 mL Pump Set (Indiv. Pack)	3	8884706220	Kangaroo 600 mL Rigid Container, Pump Set	5
773630	Kangaroo 1000 mL DEHP-Free Pump Set.....	4	8884706238	Kangaroo 600 mL Rigid Container.....	5
773656	Kangaroo ePump 1000 mL Pump Set	1	8884710826	Kangaroo Nasogastric Feeding Tube, Dobbhoff Tip	10
773662	Kangaroo ePump 1000 mL Feed set, 1000 mL Flush Set.....	1	8884710842	Kangaroo Nasogastric Feeding Tube, Dobbhoff Tip	10
774625	Kangaroo 924 Proximal Spike Set.....	3	8884710859	Kangaroo Nasogastric Feeding Tube, Dobbhoff Tip	10
775100	Kangaroo ePump Safety Screw Spike, Flush Bag.....	1	8884710867	Kangaroo Nasogastric Feeding Tube, Dobbhoff Tip	10
775655	Kangaroo ePump Safety Screw Spike Set	3	8884711006	Kangaroo Nasogastric Feeding Tube, Dobbhoff Tip	10
775659	Kangaroo ePump Safety Screw Spike.....	1	8884711246	Kangaroo Nasogastric Feeding Tube, Dobbhoff Tip	10
775662	Kangaroo ePump Safety Screw Spike Set, Flush Bag.....	3	8884711253	Kangaroo Nasogastric Feeding Tube, Dobbhoff Tip	10
775759	Kangaroo 924 Pump Safety Screw Spike Set.....	3	8884711501	Kangaroo Non-weighted Feeding Tubes.....	11
776150	100 mL Burette Recertification Set	1, 2	8884711519	Kangaroo Non-weighted Feeding Tubes.....	11
776250	Kangaroo ePump 1200 mL Rigid Container Pump Set.....	1	8884711550	Kangaroo Non-weighted Feeding Tubes.....	11
776255	Kangaroo ePump 600 mL Rigid Container Pump Set	1	8884711568	Kangaroo Non-weighted Feeding Tubes.....	11
7771210	12 Fr. Salem Sump Tube with Multifunctional Port.....	18	8884712004	Kangaroo Sterile; Bag, Easy-Cap, Ice Pouch.....	5
7771211	Multifunctional Silicone Port, Safe Enteral Connections	18	8884712100	Kangaroo Sterile Deluxe Bag, Easy-Cap.....	5
7771410	14 Fr. Salem Sump Tube with Multifunctional Port.....	18	8884712300	Kangaroo Sterile Gravity Bag.....	6
7771610	16 Fr. Salem Sump Tube with Multifunctional Port.....	18	8884713600	Kangaroo Sterile Bag, Easy-Cap, Ice Pouch.....	5
7771610CN	Salem Sump Tube, Multifunctional Port and CO ₂ Detector .	12	8884714300	Kangaroo Sterile Deluxe Bag, Easy-Cap	5
7771611	Multifunctional Silicone Port, Safe Enteral Connections	18	8884714800	Kangaroo Sterile Deluxe Bag, Easy-Cap	5
7771810	18 Fr. Salem Sump Tube with Multifunctional Port.....	18	8884715122	Kangaroo Gastrostomy Feeding Tubes, Y-Ports	15
7771810CN	Salem Sump Tube, Multifunctional Port and CO ₂ Detector .	12	8884715148	Kangaroo Gastrostomy Feeding Tubes, Y-Ports	15
7771811	Multifunctional Silicone Port, Safe Enteral Connections	18	8884716100	Kangaroo Sterile 100 mL Rigid Burette	5
777702	Kangaroo CO ₂ Detector	12	8884720163	Kangaroo Gastrostomy Feeding Tubes, Y-Ports	15
84064	Kangaroo I.V. Pole Pouch, 60 cc Flat-Top Syringe.....	6	8884720189	Kangaroo Gastrostomy Feeding Tubes, Y-Ports	15

INDEX

Code	Product	Page	Code	Product	Page
8884720205	Kangaroo Gastrostomy Feeding Tubes, Y-Ports	15	8888261701	Kangaroo Polyurethane Feeding Tube	9
8884720221	Kangaroo Gastrostomy Feeding Tubes, Y-Ports	15	8888261800	Kangaroo Polyurethane Feeding Tube	9
8884720247	Kangaroo Gastrostomy Feeding Tubes, Y-Ports	15	8888264911	Salem Sump Tube.....	19
8884720265	Kangaroo Gastrostomy Feeding Tubes, Y-Ports	15	8888264929	Salem Sump Tube.....	19
8884720285	Kangaroo Gastrostomy Feeding Tubes, Y-Ports	15	8888264945	Salem Sump Tube.....	19
8884720817	Kangaroo Nasogastric Feeding Tubes.....	11	8888264960	Salem Sump Tube.....	19
8884720825	Kangaroo Nasogastric Feeding Tubes.....	11	8888264986	Salem Sump Tube.....	19
8884720841	Kangaroo Nasogastric Feeding Tubes.....	11	8888265108	Salem Sump Silicone Tube	20
8884720855	Kangaroo Nasogastric Feeding Tubes.....	11	8888265116	Salem Sump Silicone Tube	20
8884720858	Kangaroo Nasogastric Feeding Tubes.....	11	8888265124	Salem Sump Silicone Tube	20
8884721055	Kangaroo Nasogastric Feeding Tubes.....	11	8888265132	Salem Sump Silicone Tube	20
8884721088	Kangaroo Nasogastric Feeding Tubes.....	11	8888265140	Salem Sump Silicone Tube	20
8884721096	Kangaroo Nasogastric Feeding Tubes.....	11	8888265314	Argyle Stomach Tube.....	21
8884721211	Kangaroo Nasogastric Feeding Tubes.....	11	8888265322	Argyle Stomach Tube.....	21
8884721237	Kangaroo Nasogastric Feeding Tubes.....	11	8888265330	Argyle Stomach Tube.....	21
8884721252	Kangaroo Nasogastric Feeding Tubes.....	11	8888266106	Salem Sump Tube w/ARV	19
8884721255	Kangaroo Nasogastric Feeding Tubes.....	11	8888266114	Salem Sump Tube w/ARV	19
8884741805	Stoma Measuring Device, Sterile.....	17	8888266122	Salem Sump Tube w/ARV	19
8884741813	Bolus Feeding Set.....	17	8888266130	Salem Sump Tube w/ARV	19
8884741821	Continuous Feeding Set.....	17	8888266148	Salem Sump Tube w/ARV	19
8884750095	Kangaroo Jejunal Feeding and Decompression Tube.....	13	8888266197	Salem Sump Anti-Reflux Valve.....	19
8884751622	Kangaroo 16 Fr. Y-port/P.E.G. Adapter	14	8888266213	Salem Sump Silicone Tube w/ARV.....	20
8884752042	Kangaroo 20 Fr. Y-port/P.E.G. Adapter	14	8888266221	Salem Sump Silicone Tube w/ARV.....	20
8884752505	Kangaroo Naso-jejunal Feeding Tube	13	8888266247	Salem Sump Silicone Tube w/ARV.....	20
8884752513	Kangaroo Naso-Jejunal Tube	12	8888266262	Salem Sump Silicone Tube w/ARV.....	20
8884773600	Kangaroo 1000 mL Bag, Easy-Cap, Ice Pouch	5	8888266288	Salem Sump Silicone Tube w/ARV.....	20
8888260208	Kangaroo PVC Feeding Tube	9	8888268060	Salem Sump Tube.....	19
8888260406	Kangaroo PVC Feeding Tube	9	8888268086	Salem Sump Tube.....	19
8888260604	Kangaroo PVC Feeding Tube	9	8888750018	Argyle Edlich Gastric Lavage Trays.....	21
8888260711	Kangaroo PVC Feeding Tube	9	8888750513	Argyle Edlich Gastric Lavage Trays.....	21
8888260802	Kangaroo PVC Feeding Tube	9	F010506	Kangaroo Joey Battery	2
8888261008	Kangaroo PVC Feeding Tube	9			
8888261042	Kangaroo Polyurethane Feeding Tube	9			
8888261206	Kangaroo PVC Feeding Tube	9			
8888261355	Kangaroo Polyurethane Feeding Tube	9			
8888261412	Kangaroo Polyurethane Feeding Tube	9			
8888261420	Kangaroo Polyurethane Feeding Tube	9			
8888261438	Kangaroo Polyurethane Feeding Tube	9			
8888261503	Kangaroo Polyurethane Feeding Tube	9			

The products shown in this catalog may not be licensed or registered in all parts of the world.
Please contact your local sales representative for availability

† Sentinel Event Alert, Joint Commission on Accreditation of Healthcare Organizations, Issue 36, April 3, 2006.

Cover photo credit: Thinkstock.

COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and internationally registered trademarks of Covidien AG.

™ * Hydromer is a trademark of Hydromer Inc. and is used under license.

™ * SpikeRight is a trademark of the Société des Produits Nestlé S.A.

Other brands are trademarks of a Covidien company.

©2013 Covidien.



H6200ENAF – xM – 0613