




# Exos Rapid Splint™ System

*It's About Time...*



## *...Applying Advanced Technology to the Oldest Treatment in Orthopedics.*

The Exos Rapid Splint System advances the art of immobilization with an easier, cleaner, and more efficient way to splint. Its reinforced composite polymer creates a lightweight, stable structure that quickly becomes moldable using dry heat with no water and no mess. Rapid Splint goes on smooth and dry, without sharp fiberglass edges, enhancing patient safety and comfort. It reaches stability quickly, is reformable and offers superior radiolucency. The Exos Rapid Splint System — Applying advanced technology to the oldest treatment in Orthopedics...It's about time!



*Easier and Simpler*



*Time and Efficiency*



*Patient Safety  
and Comfort*



## Dry Heat Forming

**No water needed** and no plaster or resin mess

**Easy preparation** and fewer steps save time

- ❶ Heat
- ❷ Apply
- ❸ Mold

**Adjustable curing time** to fit individual patient requirements and optimize patient throughput

**Reformable** to reduce waste and ensure best fit

**Full rigidity** in 10-15 minutes

## Exothermix Splint

**Patented Exothermix Heating**  
for all-in-one splinting to be  
done anywhere



## Rapid Splint Warmer

**Silicone Heating** for fast and easy preparation

## Patient Safety and Comfort

**Dry application** eliminates introducing moisture to the skin

**Smooth edges** without sharp fiberglass enhance patient comfort and safety

**Trimnable warm or cold** for best fit, comfort and function

**Integrated padding** for comfort and protection

**Wrinkle-free molding** for a smooth interface, better fit and less skin pressure



Posterior Short Leg



Volar Wrist



## Revolutionary Technology

Never before have such advanced technologies been used to deliver unmatched simplicity and elegance in splint immobilization. Rapid Splint's dry heat, thermoformable polymer approach is easier, cleaner, and faster for the clinician, and provides an overall better patient experience.

### **Technologically advanced reinforced composite polymer**

Thermoformable splints are high-strength, lightweight and easily formed



### **Integrated open-cell foam padding**

Enhanced comfort and protection; eliminates or reduces the need for cast padding



WITHOUT SPLINT



RAPID SPLINT



PLASTER

**Superior radiolucency** for clearer visualization of injury or fracture



Part #	Description	Part #	Description
<b>Exos Rapid Splint</b>		<b>Exos Rapid Splint with Padding</b> <i>(continued)</i>	
100-001-210	Exos Rapid Splint, 2" x 10" 10 Pack	100-001-P513	Exos Rapid Splint with Padding, 5" x 13" 5 Pack
100-001-312	Exos Rapid Splint, 3" x 12" 10 Pack	100-002-P326	Exos Rapid Splint Extra with Padding, 3" x 26" 5 Pack
100-001-336	Exos Rapid Splint, 3" x 36" 5 Pack	100-002-P426	Exos Rapid Splint Extra with Padding, 4" x 26" 5 Pack
100-001-413	Exos Rapid Splint, 4" x 13" 10 Pack	100-002-P526	Exos Rapid Splint Extra with Padding, 5" x 26" 5 Pack
100-001-436	Exos Rapid Splint, 4" x 36" 5 Pack	<b>Exos Rapid Splint with Exothermix Heater Pack</b>	
100-001-513	Exos Rapid Splint, 5" x 13" 10 Pack	100-001-E312	Exos Rapid Splint with Exothermix Heater Pack, 3" x 12" 10 Pack
100-002-326	Exos Rapid Splint Extra, 3" x 26" 5 Pack	100-001-E413	Exos Rapid Splint with Exothermix Heater Pack, 4" x 13" 10 Pack
100-002-426	Exos Rapid Splint Extra, 4" x 26" 5 Pack	100-001-E513	Exos Rapid Splint with Exothermix Heater Pack, 5" x 13" 10 Pack
100-002-526	Exos Rapid Splint Extra, 5" x 26" 5 Pack	<b>Accessories</b>	
<b>Exos Rapid Splint with Padding</b>		800-HEATERBAG	Exos Rapid Splint Warmer
100-001-P210	Exos Rapid Splint with Padding, 2" x 10" 5 Pack	79-98233	Elastic Bandages - Self Closure, 5 yards x 2" 10 Pack
100-001-P312	Exos Rapid Splint with Padding, 3" x 12" 5 Pack	79-98235	Elastic Bandages - Self Closure, 5 yards x 3" 10 Pack
100-001-P336	Exos Rapid Splint with Padding, 3" x 36" 5 Pack	79-98237	Elastic Bandages - Self Closure, 5 yards x 4" 10 Pack
100-001-P413	Exos Rapid Splint with Padding, 4" x 13" 5 Pack	79-98238	Elastic Bandages - Self Closure, 5 yards x 6" 10 Pack
100-001-P436	Exos Rapid Splint with Padding, 4" x 36" 5 Pack	800-40-STRAP	Accessory Strap 4" (10.2 cm), 10 Pack



DJO, LLC | A DJO Global Company

T 800.336.6569 D 760.727.1280 F 800.936.6569  
1430 Decision Street | Vista, CA 92081-8553 | U.S.A.

[DJOglobal.com](http://DJOglobal.com)

*Together in Motion™*

