



18-2



18-3



18-4



18-8



18-8SR

- 18-2 DERF Needle Holder 4 3/4" (12.1 cm)
- 18-3 WEBSTER Needle Holder 4 3/4" (12.1 cm), extra delicate, smooth jaws
- 18-4 WEBSTER Needle Holder 5" (12.7 cm), delicate, smooth jaws, standard pattern
- 18-8 HALSEY Needle Holder 5" (12.7 cm), smooth jaws
- 18-8SR HALSEY Needle Holder 5" (12.7 cm), serrated jaws



18-10	<b>BAUMGARTNER</b> Needle Holder 5 1/2" (14 cm)
18-14	<b>COLLIER</b> Needle Holder 5" (12.7 cm), fenestrated jaws
18-16	<b>BROWN</b> Needle Holder 5 1/4" (13.3 cm), convex jaws, fine serrations
18-18	<b>BROWN</b> Needle Holder 6 3/4" (13.3 cm), convex jaws, fine serrations
<b>OLSEN-HEGAR NEEDLE HOLDERS</b> — Needle holder and Suture Scissors Combined	
18-20	4 3/4" (12.1 cm) extra delicate pattern, smooth jaws
18-22	4 3/4" (12.1 cm) extra delicate pattern, serrated jaws
18-24	5 1/2" (14 cm), serrated jaws
18-26	6 1/2" (16.5 cm), serrated jaws
18-28	7 1/2" (18.4 cm), serrated jaws



18-48	BOYNTON Needle Holder 5" (12.7 cm), delicate
18-50	MATHIEU Needle Holder 5 1/2" (14 cm), with spring handles
18-60	MAYO-HEGAR Needle Holder 5" (12.7 cm)
18-62	MAYO-HEGAR Needle Holder 6" (15.2 cm)
18-64	MAYO-HEGAR Needle Holder 7" (17.8 cm)
18-66	MAYO-HEGAR Needle Holder 8" (20.3 cm)
18-70	MAYO-HEGAR (MASSON) Needle Holder 10 1/2" (26.7 cm)
18-80	CRILE-WOOD Needle Holder 6" (15.2 cm)
18-82	CRILE-WOOD Needle Holder 7" (17.8 cm)
18-84	CRILE-WOOD Needle Holder 8" (20.3 cm)
18-86	CRILE-WOOD Needle Holder 9" (22.9 cm)



18-120 SAROT Needle Holder 7 1/8" (18.1 cm) light model

18-122 SAROT Needle Holder 10 1/2" (26.7 cm)

18-126 HEANEY Needle Holder 8 1/2" (21.6 cm), curved jaws

18-128 STRATTE Needle Holder 9" (22.9 cm), double curve jaws and bent shanks



18-200  
18-202



18-204  
18-206



18-308

18-200	CASTROVIEJO Needle Holder 5 1/2" (14 cm), straight, without lock, smooth jaws
18-202	CASTROVIEJO Needle Holder 5 1/2" (14 cm), curved, without lock, smooth jaws
18-204	CASTROVIEJO Needle Holder 5 1/2" (14 cm), straight, with lock, smooth jaws
18-206	CASTROVIEJO Needle Holder 5 1/2" (14 cm), curved, with lock, smooth jaws
18-308	BARRAQUER Needle Holder 5 1/8" (13 cm), curved, delicate tapered jaws, 8 mm wide round handle

## TUNGSTEN CARBIDE NEEDLE HOLDERS

Tungsten Carbide Needle Holders grasp surgical needles far more securely than conventional stainless steel needle holders. The Tungsten Carbide inserts, the hardest steel made for surgical instruments, are bonded to the finest quality stainless steel needle holders. These highly durable inserts securely grip the needle, preventing it from rotating or slipping.

While this needle holder has a higher initial cost than standard stainless steel needle holders, it will provide the user with superior performance and longer life.

TUNGSTEN  
CARBIDE



18-3C  
18-3SRC

18-4C

**WEBSTER NEEDLE HOLDER**—Extra delicate jaws ideal for facial, pediatric and ophthalmic procedures. Recommended for 7/0 and finer suture.

18-3C 4 3/4" (12.1 cm), smooth jaws  
18-3SRC 4 3/4" (12.1 cm), serrated jaws

**WEBSTER Needle Holder**—Delicate smooth jaws ideal for facial, pediatric and ophthalmic procedures. Recommended for 7/0 and finer suture.

18-4C 5 1/4" (13.3 cm), smooth jaws

**TUNGSTEN  
CARBIDE**



**PAR Needle Holder**—Extra Delicate, serrated jaws for improved needle control. Ideal for facial, pediatric and hand procedures. Recommended for 7/0 and finer suture.

18-5C 4 1/2" (11.4 cm)  
18-6C 5" (12.7 cm)

**HALSEY Needle Holder**—General purpose plastic surgery needle holder. Suggested for 6/0 and finer suture.

18-8C 5" (12.7 cm), smooth jaws  
18-8SRC 5" (12.7 cm), serrated jaws

**BAUMGARTNER Needle Holder**—Heavier pattern ideal for minor procedures and routine skin closures with serrated jaws. Recommended for 3/0 to 6/0 suture.

18-10C 5 1/2" (14 cm)

**OLSEN-HEGAR Needle Holder**—Combination needle holder and scissors with delicate jaws. Recommended for 6/0 and finer suture.

18-20C 4 3/4" (12.1 cm), smooth jaws  
18-22C 4 3/4" (12.1 cm), serrated jaws

TUNGSTEN  
CARBIDE



**OLSEN-HEGAR Needle Holder**—Heavier pattern combination needle holder and scissors with serrated jaws. Recommended for 4/0 to 2 suture.  
 18-24C 5 1/2" (14 cm)    18-26C 6 1/2" (16.5 cm)    18-28C 7 1/2" (18.4 cm)

**MAYO-HEGAR Needle Holder**—Heavier general purpose needle holder with serrated jaws. Recommended for 4/0 to 4 suture.  
 18-60C 5" (12.7 cm)    18-64C 7" (17.8 cm)    18-70C 10 1/2" (26.7 cm)  
 18-62C 6" (15.2 cm)    18-66C 8" (20.3 cm)

**CRILE-WOOD Needle Holder**—General purpose needle holder with serrated jaws. Recommended for 3/0 to 6/0 suture.  
 18-80C 6" (15.2 cm)    18-82C 7" (17.8 cm)    18-84C 8" (20.3 cm)    18-86C 9" (22.9 cm)

**COOLEY Needle Holder**—Extra delicate with finely tapered serrated jaws ideal for vascular procedures. Recommended for 2/0 to 6/0 suture.  
 18-90C 6" (15.2 cm)    18-92C 6 3/4" (17.1 cm)

TUNGSTEN  
CARBIDE



**SAROT Needle Holder**—Bend in shank enables jaws to open with minimal movement of handles and improves tip visualization. Ideal for small incisions with limited exposure with serrated jaws  
 18-120C 7 1/4" (18.4 cm)      18-122C 10 1/2" (26.7 cm)

**HEANEY Needle Holder**—Curved heavy design ideal for placement of sutures in difficult to reach places with serrated jaws. Recommended for 3/0 to 4  
 18-126C 8 1/2" (21.6 cm)

**STRATTE Needle Holder**—Double curved jaws and shanks allows placement of needles at difficult angles with serrated jaws  
 18-128C 9" (22.9 cm)

**RYDER-Type Intracardiac Needle Holder**—Fine parallel jaws prevent flattening of needle ideal for cardiovascular procedures with serrated jaws. Recommended for 4/0 to 6/0  
 18-150C 5" (12.7 cm)      18-152C 6" (15.2 cm)      18-154C 7" (17.8 cm)      18-156C 9" (22.9 cm)

TUNGSTEN  
CARBIDE



18-204C  
18-204SRC



18-206C  
18-206SRC



18-220C  
18-222C  
18-225C



18-240C  
18-246C  
18-250C

**CASTROVIEJO Needle Holder** — Flat handle provides non-slip control ideal for facial and hand surgery. Recommended for 6/0 and finer suture.

18-204C 5 1/2" (14 cm), straight smooth jaws with lock  
18-204SRC 5 1/2" (14 cm), straight serrated jaws with lock

**CASTROVIEJO Needle Holder**— Flat handle provides non-slip control ideal for facial and hand surgery. Recommended for 6/0 and finer suture

18-206C 5 1/2" (14 cm), curved smooth jaws with lock  
18-206SRC 5 1/2" (14 cm), curved serrated jaws with lock

**JACOBSON Microvascular Needle Holder** — Fine serrated jaws with light action ideal for vascular surgery. Recommended for 7/0 and finer suture

18-220C 7" (17.8 cm), straight, serrated jaws with lock  
18-222C 7" (17.8 cm), curved, serrated jaws with lock  
18-225C 8 1/2" (21.6 cm), serrated jaws with lock

**JACOBSON Microvascular Needle Holder** — Round handles reduce hand movement and enhance finger control. Ideal for vascular surgery. Recommended for 8/0 and finer suture

18-240C 6" (15.2 cm), straight jaws with lock  
18-246C 7" (17.8 cm), straight jaws with lock  
18-250C 8" (20.3 cm), straight jaws with lock



18-504C

18-562C  
18-564C  
18-566C

18-580C  
18-582C  
18-584C  
18-586C

## AUTO-RELEASE NEEDLE HOLDERS

Opens with just a squeeze of the palm. Suitable for left handed use.

18-504C	WEBSTER Auto-Release Needle Holder 5 1/4" (13.3 cm), smooth jaws, carbide
18-562C	MAYO-HEGAR Auto-Release Needle Holder 6" (15.2 cm) serrated jaws, carbide
18-564C	MAYO-HEGAR Auto-Release Needle Holder 7" (17.8 cm), serrated jaws, carbide
18-566C	MAYO-HEGAR Auto-Release Needle Holder 8" (20.3 cm), serrated jaws, carbide
18-580C	CRILE-WOOD Auto-Release Needle Holder 6" (15.2 cm), serrated jaws, carbide
18-582C	CRILE-WOOD Auto-Release Needle Holder 7" (17.8 cm), serrated jaws, carbide
18-584C	CRILE-WOOD Auto-Release Needle Holder 8" (20.3 cm), serrated jaws, carbide
18-586C	CRILE-WOOD Auto-Release Needle Holder 9" (22.9 cm), serrated jaws, carbide

TUNGSTEN  
CARBIDE



## LEFT HANDED NEEDLE HOLDERS

18-1004C WEBSTER Needle Holder 5 1/4" (13.3 cm) smooth jaws

18-1010C BAUMGARTNER Needle Holder 5 1/2" (14 cm) serrated jaws

18-1062C MAYO-HEGAR Needle Holder 6" (15.2 cm) serrated jaws

18-1064C MAYO-HEGAR Needle Holder 7" (17.8 cm) serrated jaws

18-1080C CRILE-WOOD Needle Holder 6" (15.2 cm) serrated jaws

18-1082C CRILE-WOOD Needle Holder 7" (17.8 cm) serrated jaws