

FOOD MICROBIOLOGY

PRODUCT NAME AND DESCRIPTION	TYPE/FILL	PACK SIZE	PRODUCT CODE
ANTISERA (Please refer to the Diagnostic Tests Section)			
PREPARED MEDIA			
Acetate Differential Agar, Slant	15 x 103mm	20/Pk.	R060022
Brain Heart Infusion (BHI) Broth	500mL	10/Pk.	R112022
Brilliant Green Bile (2%) Broth For detection of coliform microorganisms in water, wastewater, foods, and dairy products. Fermentation is detected by gas in Durham tube. Industrial Reference: AOAC, BAM, EPA, IDF, SMWW.			
Brilliant Green Bile Broth 2%	E - DT	20/Pk.	R07022
Brilliant Green Bile Broth 2%	E - DT	100/Pk.	R07020
Decarboxylase Broth Lysine, 5mL	15 x 103mm, tube	20/Pk.	R060760
EMB Agar Levine For isolation and differentiation of gram negative enteric bacilli from clinical and industrial samples. Industrial Reference: BAM, COMPF, IDF, USP.			
EMB Agar, Levine	Monoplate	10/Pk.	R01400
EMB Agar, Levine	Monoplate	100/Pk.	R01402
EMB Agar, Levine	Spacesaver Plate	15/Pk.	R01401
E.C. Medium For detection of coliform bacteria at 37°C and <i>Escherichia coli</i> at 44.5°C. Industrial Reference: BAM, EPA, SMWW. Not intended for IVD use.	E - 10mL DT	20/Pk.	R07102
E.C. Medium w/MUG For detection of <i>Escherichia coli</i> in water and food samples by a fluorogenic assay procedure. Industrial Reference: BAM. Not intended for IVD use.			
EC Medium w/MUG	E - 10mL	20/Pk.	R07655
EC Medium w/MUG	E - 10mL DT	20/Pk.	R07657
MacConkey Agar For selective and differential isolation of gram-negative bacilli on the basis of lactose fermentation. Industrial Reference: BAM, COMPF, SMWW, USP.			
MacConkey Agar	Monoplate	10/Pk.	R01550
MacConkey Agar	Monoplate	100/Pk.	R01552
MacConkey Agar	Spacesaver Plate	15/Pk.	R01551
MacConkey Agar	Spacesaver Plate	150/Pk.	R01562
MacConkey Agar	D - slant	20/Pk.	R061322
MacConkey Agar w/MUG For isolation and presumptive identification of <i>Escherichia coli</i> by a fluorogenic assay procedure.	Monoplate	10/Pk.	R01554
MacConkey Agar w/Sorbitol (SMAC) For selective and differential isolation of sorbitol-negative <i>Escherichia coli</i> associated with verotoxin production, especially O157:H7. Industrial Reference: BAM, COMPF.			
MacConkey Agar w/Sorbitol	Monoplate	10/Pk.	R01556
MacConkey Agar w/Sorbitol	Monoplate	100/Pk.	R01561
MacConkey Agar w/Sorbitol	G - PT	20/Pk.	R09548
MacConkey Agar w/Sorbitol, Cefixime, Tellurite (CT-SMAC) For selective and differential isolation of sorbitol-negative <i>Escherichia coli</i> associated with verotoxin production, especially O157:H7.	Monoplate	10/Pk.	R110241
MacConkey Agar w/Sorbitol, MUG For isolation of Gram-negative bacilli and for differentiation of those microorganisms based on lactose fermentation.	Spacesaver Plate	15/Pk.	R01563



Sorbitol MacConkey (SMAC) with BCIG Agar Plate