



Bonded Syringe

Needle-Free Syringe with Bonded Texium Closed Male Luer. Not made with DEHP. PV: 30 mL

Summary

SKU	MY8030
Product name	Bonded Syringe
Product Type	

Weight (lbs)	3.7500	Country of Manufacture		Approximate Priming Volume (mL)	30
Case Dimensions	10.687in x 10.187in x 19.875in 3.75lb	Not Made with Natural Rubber	Yes	Luer Lock	Texium Closed Male Luer
Manifold	None	Number of Back Check Valves	0	Number of Connectors	1
Number of Fuses	0	Removable or Bonded Connectors	Bonded	Sterility	Fluid Path Sterile
Stopcock	None	Texium	Yes	Color of Tubing	None
T-Connector Type	None	Number of y-sites	0	Y-Site Port Type	None
Type of Connector	Texium	Filter	None		