

## COMPRESSION RECOVERY SOCKS

### MATERIALS

- Advanced Lycra/nylon knit structure
- Bamboo charcoal fibers at the toe and heel
- COOLMAX® interwoven fabric
- Spandex 25%, COOLMAX® 35%, Nylon 30%
- White Socks: 10% Cotton from Bamboo Charcoal, Black Socks: 10% Rayon from Bamboo Charcoal

### FEATURES & BENEFITS

- Gradient 20-30mm Hg compression to enhance circulation for faster recovery
- Ventilated design with COOLMAX® moisture-wicking fibers for cool, comfortable wear
- Comfortable top and bottom bands hold sleeve in place without constricting
- Figure-8 support to increase stabilization of the arch and ankle
- Bamboo charcoal fibers to help regulate temperature, bacteria and perspiration
- Embedded Achilles and arch support to alleviate aches and pains

**IDEAL FOR:** Athletes and runners to speed recovery and provide support during or after running or athletic activity.



Part No	Size	Shoe Size	Fits Calf	Fits Ankle	UPC Code
6020-01MS	Men's Small, White/Gray	7 & under	11 – 15"	7.5 – 8.5"	749852111198
6020-01MM	Men's Medium, White/Gray	7.5 – 10	13 – 17"	8.5 – 9.5"	749852111204
6020-01ML	Men's Large, White/Gray	10.5 – 12	14 – 18"	9 – 11.25"	749852111211
6020-01MX	Men's X-Large, White/Gray	12.5 & over	15 – 20"	11.5 – 12.5"	749852111228
6020-11MS	Men's Small, Black/Gray	7 & under	11 – 15"	7.5 – 8.5"	749852111235
6020-11MM	Men's Medium, Black/Gray	7.5 – 10	13 – 17"	8.5 – 9.5"	749852111242
6020-11ML	Men's Large, Black/Gray	10.5 – 12	14 – 18"	9 – 11.25"	749852111259
6020-11MX	Men's X-Large, Black/Gray	12.5 & over	15 – 20"	11.5 – 12.5"	749852111266
6020-15MS	Men's Small, Black/Green	7 & under	11 – 15"	7.5 – 8.5"	749852111273
6020-15MM	Men's Medium, Black/Green	7.5 – 10	13 – 17"	8.5 – 9.5"	749852111280
6020-15ML	Men's Large, Black/Green	10.5 – 12	14 – 18"	9 – 11.25"	749852111297
6020-15MX	Men's X-Large, Black/Green	12.5 & over	15 – 20"	11.5 – 12.5"	749852111303
6025-01WS	Women's Small, White/Gray	4 – 5	10.5 – 14.5"	6.5 – 8"	749852111310
6025-01WM	Women's Medium, White/Gray	5.5 – 7.5	11.5 – 15.5"	8 – 9.5"	749852111327
6025-01WL	Women's Large, White/Gray	8 – 10.5	12.5 – 17"	9.5 – 11"	749852111334
6025-01WX	Women's X-Large, White/Gray	10.5 & over	13.5 – 17.5"	11 – 12.5"	749852111341
6025-15WS	Women's Small, Black/Green	4 – 5	10.5 – 14.5"	6.5 – 8"	749852111358
6025-15WM	Women's Medium, Black/Green	5.5 – 7.5	11.5 – 15.5"	8 – 9.5"	749852111365
6025-15WL	Women's Large, Black/Green	8 – 10.5	12.5 – 17"	9.5 – 11"	749852111372
6025-15WX	Women's X-Large, Black/Green	10.5 & over	13.5 – 17.5"	11 – 12.5"	749852111389
6025-18WS	Women's Small, Black/Pink	4 – 5	10.5 – 14.5"	6.5 – 8"	749852111396
6025-18WM	Women's Medium, Black/Pink	5.5 – 7.5	11.5 – 15.5"	8 – 9.5"	749852111402
6025-18WL	Women's Large, Black/Pink	8 – 10.5	12.5 – 17"	9.5 – 11"	749852111419
6025-18WX	Women's X-Large, Black/Pink	10.5 & over	13.5 – 17.5"	11 – 12.5"	749852111426

**APPLICATION INSTRUCTIONS:** Gather sock from top and pull over foot and calf. Smooth out wrinkles while pulling over the calf.

**CARE INSTRUCTIONS:** Hand or machine wash on gentle cycle in cold water using a mild detergent. Air-dry. Do not machine dry or use bleach, fabric softener or dryer sheets.

Not intended for persons with severe foot deformities or circulatory disorders resulting in infections or foot ulcerations. Always consult a medical professional when experiencing persistent pain, discomfort or swelling.



888.237.3668 • 513.825.1888  
 Email: info@powersteps.com  
 www.powerstep.com  
 Stable Step, LLC, West Chester, OH