

Tools

> Home > Healthcare > Laboratory Diagnostics > Point of Care Testing > Urinalysis > CLINITEK Status+ Analyzer

CLINITEK Status+ Analyzer

Need simple, enhanced clinical information for urinalysis results?

Overview > Features & Benefits > Technical Specifications



The CLINITEK Status[®]+ Analyzer features new automatic checks (Auto-Checks^{*}). The analyzer automatically checks each test strip for humidity exposure, common sample interferences^{**} and strip identification for Siemens test strips. Together, these provide improved clinical information.


- Enhances confidence in clinical decisions with (Auto-Checks) performed on every test
- Ensures consistent results interpretation with automatic timing and reading of every urine specimen
- Ready for connectivity to data management solutions (DMS), via a simple upgrade kit
- Offers a comprehensive point-of-care urinalysis menu for routine testing, detection of early kidney disease and hCG pregnancy testing

CLINITEK Status+ Analyzer
Need simple, enhanced clinical
information for urinalysis results?




Related Products, Services & Resources

>




>




Chemistry Series

>



>



Resources - Point of Care
Resources at Your Fingertips

*Only available when using Siemens test strips with IR or color bands

**Feature not available in the U.S.