

Medi-Vac® Suction and Fluid Collection Products Catalog

Guardian™ Hardware



CardinalHealth

Medi-Vac® Suction and Fluid Collection Products – Guardian™ Hardware

Our number one goal is to help you stay competitive in today's health-care environment. To that end, we offer families of products that allow you to meet the diverse needs of your facility and control the cost of doing business.

The Guardian™ suction canister system features single-use disposable hard locking lid attachments. Special brackets and stands are available for convenient mounting configurations and mobility. Designed for safety, economy and convenience, Guardian canister systems are available in a variety of sizes to meet all of your suction and collection needs.

For more information about the Guardian system or to place an order, please contact your Cardinal Health sales representative or call customer service at 1.800.964.5227.

Disposable Canisters

Guardian™ disposable hard canisters feature rigid blue lids with an automatic shutoff valve to help prevent cross contamination of regulators and wall vacuum outlets. Its 90° adapters allow tubing to connect at right angles to help prevent kinking and impeded fluid flow. The locking lid encourages the proper disposal of infectious liquid medical waste and enhances worker safety. Lid includes accessory and orthopedic ports.



1200cc Canister



2000cc Canister



3000cc Canister

Guardian Canisters

Cat. No.	Size	Qty.
65651-212	1200cc	40/cs
65651-220	2000cc	40/cs
65651-230	3000cc	40/cs

Canisters with PPV (Porous Plastic Valve)

The Guardian PPV suction canister functions in smoke-intense environments, reducing the likelihood of occlusion or premature shutoff.

Guardian PPV Canisters

Cat. No.	Size	Qty.
65651-412	1200cc	40/cs
65651-420	2000cc	40/cs
65651-430	3000cc	40/cs

Direct-to-Regulator Canister

Disposable hard canister with rigid blue lid mounts directly to regulator without the need for a ring bracket or additional tubing.

Guardian Direct-to-Regulator Canister

Cat. No.	Size	Qty.
65651-100	1200cc	36/cs



Direct-to-Regulator Canister

Ring Brackets

Support devices used to configure your Guardian™ suction system.

Ring Brackets

Reusable ring brackets are constructed of durable polycarbonate for years of maintenance-free use.

Brackets slide into a wall plate or roll stand. 1 each.



Guardian Ring Bracket

Guardian Ring Bracket

Reusable polycarbonate bracket used to support the Guardian canister in any slide mount such as the wall plate or mobile suction cart. Has preattached tubing for connection to canister.

- 65652-155 Used with 1200cc disposable canister
Wall to outer edge of ring: 7³/₄"
- 65652-135 Used with 2000 and 3000cc disposable canister
Wall to outer edge of ring: 9¹/₂"



Guardian Close-to-Wall Ring Bracket

Guardian Close-to-Wall Ring Bracket

Reusable polycarbonate bracket used to support the Guardian canister in any slide mount such as the wall plate or mobile suction cart. Designed to be close to the wall.

- 65652-148 Used with 1200cc disposable canister
Wall to outer edge of ring: 6³/₄"
- 65652-168 Used with 2000 and 3000cc disposable canister
Wall to outer edge of ring: 9¹/₂"



Guardian Ring Bracket with Built-in On/Off Valve

Guardian Ring Bracket with Built-in On/Off Valve

Reusable polycarbonate bracket used to support the Guardian canister in any slide mount such as the wall plate or mobile suction cart. Comes with on/off valve built into bracket. Has preattached tubing for connection to canister.

- 65652-159 Used with 1200cc disposable canister
Wall to outer edge of ring: 7³/₄"
- 65652-139 Used with 2000 and 3000cc disposable canister
Wall to outer edge of ring: 9¹/₂"



Guardian Ring Bracket Extended for Regulator

Guardian Ring Bracket Extended for Regulator

Reusable polycarbonate bracket used to support the Guardian canister in any slide mount such as the wall plate or mobile suction cart. Used to mount regulator on ring bracket.

- 65652-143 Used with 1200cc disposable canister
Wall to outer edge of ring: 9¹/₂"
- 65652-144 Used with 2000 and 3000cc disposable canister
Wall to outer edge of ring: 10"



Guardian Ring Bracket Extended for Regulator (with D.I.S.S. and Hose Assembly)

Guardian Ring Bracket Extended for Regulator (with D.I.S.S. and Hose Assembly)

Used to mount regulator on bracket assembly. Comes with the following preattached: 65652-131– Universal Female D.I.S.S. Preattached tubing for connection to canister.

- 65652-173 Used with 1200cc disposable canister
Wall to outer edge of ring: 9¹/₂"
- 65652-174 Used with 2000 and 3000cc disposable canister
Wall to outer edge of ring: 10¹/₄"

Large Volume Collection Canister (LVC)



Canister and 18" roll stand

Guardian™ Large Volume Collection Canister (LVC)

Tough, shatter-resistant hard plastic canister collects 12L. Sturdy, ergonomic handle and gradations allow for viewing and measuring contents. Available with 18" or 30" roll stand for easy mobility.

Large Volume Collection System

Cat. No.	Technical Information	Qty.
65651-120	Guardian LVC	4/cs
65652-583	Guardian LVC roll stand (30")	1
65652-598	Guardian LVC roll stand (18")	1

Stands

Stands

Holders for Guardian™ canisters.



Single-Canister Base Stand



Single-Canister Roll Stand



Double-Canister Roll Stand

Single-Canister Base Stand (65652-582, 65652-524, 65652-911)

Holds Guardian canister for mounting on various surfaces such as portable pumps, shelves and carts. Double-sided foam tape for mounting purposes is preattached.

Single-Canister Roll Stand (65652-161, 65652-163, 65652-924)

Holds one Guardian canister. Mounted on four casters for stability and mobility.

Double-Canister Roll Stand (65652-181, 65652-182)

Holds two Guardian canisters. Mounted on four casters for stability and mobility.

Stands

Cat. No.	Size	Qty.
65652-582	1200cc	1
65652-524	2000cc	1
65652-911	3000cc	1
65652-161	1200cc	1
65652-163	2000cc	1
65652-924	3000cc	1
65652-181	1200cc	1
65652-182	2000cc	1

Metal Roll Stands

Additional roll stands for Guardian canisters. 1 each.



Metal Roll Stand with Regulator Hardware, Single



Metal Roll Stand with Regulator Hardware, Double

Metal Roll Stand, Single (65652-164, not pictured)

Holds single 2000 or 3000cc canister.

Metal Roll Stand with Regulator Hardware, Single (65652-264, pictured)

Holds single 2000 or 3000cc canister.

Metal Roll Stand, Double (65652-183, not pictured)

Holds two 2000 or 3000cc canisters.

Metal Roll Stand with Regulator Hardware, Double (65652-283, pictured)

Holds two 2000 or 3000cc canisters.

Medi-Vac® Accessories

Medi-Vac® accessories make it easy to add auxiliary items to your system. Designed for use with CRD™, Flex Advantage® and Guardian™ systems, Medi-Vac accessories are interchangeable to help you select the most effective setup for your needs.

Specimen Socks

Used to collect and retain specimens for analysis.
10 per case.



Standard Specimen Sock

Standard Specimen Sock (65652-122)

Used to convert pour/accessory port into a male port specimen retainer. Accepts up to $\frac{3}{8}$ " I.D. tubing with male port. Length 9".



Short Specimen Sock

Short Specimen Sock (65652-123)

Used to convert pour/accessory port into a male port specimen retainer. Accepts up to $\frac{3}{8}$ " I.D. tubing with male port. Length 4 $\frac{1}{2}$ ".



Short Specimen Sock Without Stick

Short Specimen Sock Without Stick (65652-132)

Used to convert pour/accessory port into a male port specimen retainer. Accepts up to $\frac{3}{8}$ " I.D. tubing with male port. Length 4 $\frac{1}{2}$ ".



Large-Port Specimen Sock

Large-Port Specimen Sock (65652-124)

Used to convert pour/accessory port into a male port specimen retainer. Accepts $\frac{3}{8}$ " to $\frac{1}{2}$ " I.D. tubing with male port. Length 9".



Female Specimen Sock

Female Specimen Sock (65652-169)

Used to convert the pour/accessory port into a female port specimen retainer which will accept tubing sets with preattached hard plastic (or metal) male end fittings. Length 9".



Pour/Accessory Port Convertor

Pour/Accessory Port Convertor (65652-119)

Used to convert the pour/accessory port into an additional patient port. Accepts up to $\frac{3}{8}$ " I.D. tubing with male port.



Pour/Accessory Port Convertor/Large

Pour/Accessory Port Convertor/Large (65652-120)

Used to convert the pour/accessory port into an additional patient port. Accepts $\frac{3}{8}$ " to $\frac{1}{2}$ " I.D. tubing with male port.

Medi-Vac® Accessories *(continued)*

Tubing

Flexible, connecting tubing to adapt system for various setups.



Bulk Connecting Tubing

Bulk Connecting Tubing (516500)

Used to connect vacuum source to canisters and manifolds. 8mm I.D. 100' roll per box; 4 boxes per case.



CRD™ and Guardian™ Tandem Tube

CRD™ and Guardian™ Tandem Tube (65652-542A)

Used to connect canisters in a series for high-volume setups. Not compatible with Flex Advantage® liners. 18" long x 8mm I.D. 30 per case.

Critical Measure

Small-volume measuring device.



Universal Critical Measure

Universal Critical Measure (65651-700)

Used to closely monitor collected fluid volume when precise measurements are needed. Graduated in 5cc increments up to 300cc. Fits over pour/accessory port. Includes 12" tandem tube. 50 per case.

Right-Angle Adapters

Special-use connectors.



Right-Angle Adapter

Right-Angle Adapter (65652-107)

Used to connect patient tubing to CRD and Guardian patient ports at a 90° angle. 50 per case.



Right-Angle Adapter for Flex Advantage Port and CRD and Guardian Ortho Ports

Right-Angle Adapter for Flex Advantage Patient Port and CRD and Guardian Ortho Ports (65652-109)

Used in 90° right-angle connections. 50 per case.

Universal Hardware

Medi-Vac® universal hardware products are the common denominators for all Medi-Vac suction systems. The carts, valves, manifolds and connectors work in conjunction with CRD™, Flex Advantage® and Guardian™ suction canisters. These interchangeable parts make it easy to customize and/or change your suctioning systems.

Vacuum Manifolds

Used to connect one vacuum source to multiple canisters.



Two-Canister Manifold

Two-Canister Manifold (65652-002)

Used to connect one vacuum source to two canisters. 1 each.



Three-Canister Manifold

Three-Canister Manifold (65652-003)

Used to connect one vacuum source to three canisters. 1 each.



Four-Canister Manifold

Four-Canister Manifold (65652-004)

Used to connect one vacuum source to four canisters. 1 each.

Valve

Used to turn vacuum supply on and off as well as regulate the flow of the vacuum supply to the canister.



In-Line Vacuum On/Off Valve

In-Line Vacuum On/Off Valve (65652-580)

Used to turn vacuum supply on and off as well as regulate the flow of the vacuum supply to the canister. 1 each.

Connectors

Allow customization of suction system.



Pigtail Connector

Pigtail Connector (65652-104)

For connection to vacuum source on portable pumps or to connect large I.D. vacuum tubing to small tubing barbs. Contains latex. 3" latex tube attached to tubing barb. 1 each.



Tubing Barb

Tubing Barb (65652-108A)

Used to connect two pieces of tubing. Male to male. 50 per case.



Tee Connection

Tee Connection (65652-543)

Used to connect tubing for various applications. 1 each.

Roll Stands

Roll stands allow flexibility in the location and configuration of Medi-Vac® canister (CRD™, Flex Advantage® or Guardian™ system) installations.



Tall



Tall with Regulator Mount

Tall (65652-586),

Tall with Regulator Mount (65652-596)

30" from floor to canister top. Four slide mounting locations for canister placements. 1 each.



Medium



Medium with Regulator Mount

Medium (65652-500),

Medium with Regulator Mount (65652-510)

20" from floor to canister top. Four slide mounting locations for canister placements. 1 each.



Short



Short with Regulator Mount

Short (65652-587),

Short with Regulator Mount (65652-597)

11" from floor to canister top. Four slide mounting locations for canister placements. 1 each.

Not compatible with 3000cc Flex Advantage canisters.



Tall Eight-Canister
Dual-Head Roll Stand

Tall Eight-Canister Dual-Head Roll Stand (65652-686)

Each head has four slide mounting locations for canister placement (holds a total of eight canisters). 1 each.

Universal Hardware *(continued)*

D.I.S.S. Connection Accessories

Allow use with standard vacuum regulator Diameter Index Safety System (D.I.S.S.) connections.



Universal Female D.I.S.S.

Universal Female D.I.S.S. (65652-131)

Used to support a regulator on a ring bracket or roll stand or to suspend a regulator ring bracket from a regulator. $\frac{1}{8}$ " national pipe thread plastic adapter. 1 each.



Male D.I.S.S. (Metal)

Male D.I.S.S. (Metal) (65652-137)

Metal replacement for plastic adapter. $\frac{1}{8}$ " national pipe thread metal adapter. 1 each.

Wall Plates and Adapters

Attachment devices for use with ring brackets.



Wall Plate

Wall Plate (65652-145)

Used to attach ring brackets to a vertical surface. Predrilled with four $\frac{1}{4}$ " holes. 1 each.



Variable Height Adapter

Variable Height Adapter (Low Boy) (65652-134)

Used with existing wall plate to attach ring brackets at varying heights below the original wall plate height. 1 each.

Bracket

Support device for a regulator.



Regulator Bracket

Regulator Bracket (65652-149)

Used to attach regulator to wall plate or mobile suction cart. 1 each.

SAF-T Pump™ Alternate Disposal System

The SAF-T Pump™ system quickly and safely empties canisters containing infectious liquid medical waste into the sanitary sewer with no pouring required.



- Fast – drains (1) liter of fluid in just 4 seconds
- Powered entirely by tap water – no moving parts to wear out
- Easy to use – just insert a dip tube, connect the pump and turn it on
- Avoids the added expense of disposing of full canisters
- Accommodates all Medi-Vac® suction canister sizes – basket holds multiple canisters securely
- Provides a closed system when used with Medi-Vac canisters

SAF-T Pump System and Accessories

Cat. No.	Description	Qty. Cs
65652-000	SAF-T Pump	1
65652-008	Dip tube, 8"	50
65652-011	Dip tube, 11"	50
65652-016	Dip tube, 16"	50
65652-024	Pump tube, 24"	10
65652-044	Pump tube, 44"	10
65652-001	Dispenser for dip tubes	1

Dip Tube Compatibility Guide

	65652-008	65652-011	65652-016
1200cc Guardian™ canister	•		
2000cc Guardian canister	•		
3000cc Guardian canister		•	
Guardian large-volume collection canister (LVC)			•
All CRD™ canister sizes	•		
1000cc Flex Advantage® liners	•		
1500cc Flex Advantage liners	•		
3000cc Flex Advantage liners			•

Multiple Guardian™ Canister Assembly

This is a multiple Guardian™ canister setup with two suction sites available. It has been arranged with two 1200cc Guardian canisters in tandem and two 3000cc Guardian canisters for 8400cc total collection capacity. Other configurations are available to meet your specific suctioning requirements.



Multiple Guardian Canister Assembly



Guardian Ring Bracket
with Built-in On/Off Valve
65652-139



Bulk Tubing
516500



Tandem Tube
65652-542A



Four-Canister Manifold
65652-004



Roll Stand, Tall
65652-586

Assembly Options

The Guardian™ system is available in several assembly options, allowing for flexibility at the point of use.



Wall-Mount Single-Canister Assembly



Base-Mounted Assembly

Wall-Mount Single-Canister Assembly

This is a wall-mount single-canister installation typical of a patient room situation. The Guardian ring bracket is available with or without the on/off valve.

Base-Mounted Assembly

This is a typical base-mounted installation, commonly used on portable vacuum pumps, shelves, etc.


Guardian™ Hardware Catalog Information

Cat. No			Description (Unit of Measure)
1200cc	2000cc	3000cc	Guardian™ Hardware
<input type="checkbox"/> 65651-212	<input type="checkbox"/> 65651-220	<input type="checkbox"/> 65651-230	Disposable Canister (40/case)
<input type="checkbox"/> 65651-412	<input type="checkbox"/> 65651-420	<input type="checkbox"/> 65651-430	Guardian PPV Canister (40/case)
<input type="checkbox"/> 65651-100	—	—	Direct-to-Regulator Canister (36/case)
<input type="checkbox"/> 65652-155	<input type="checkbox"/> 65652-135	<input type="checkbox"/> 65652-135	Guardian Ring Bracket (1 each)
<input type="checkbox"/> 65652-148	<input type="checkbox"/> 65652-168	<input type="checkbox"/> 65652-168	Guardian Close-to-Wall Ring Bracket (1 each)
<input type="checkbox"/> 65652-159	<input type="checkbox"/> 65652-139	<input type="checkbox"/> 65652-139	Guardian Ring Bracket with Built-in On/Off Valve (1 each)
<input type="checkbox"/> 65652-143	<input type="checkbox"/> 65652-144	<input type="checkbox"/> 65652-144	Guardian Ring Bracket Extended for Regulator (1 each)
<input type="checkbox"/> 65652-173	<input type="checkbox"/> 65652-174	<input type="checkbox"/> 65652-174	Guardian Ring Bracket Extended for Regulator (with D.I.S.S. and Hose Assembly) (1 each)
<input type="checkbox"/> 65652-582	<input type="checkbox"/> 65652-524	<input type="checkbox"/> 65652-911	Single-Canister Base Stand (1 each)
<input type="checkbox"/> 65652-161	<input type="checkbox"/> 65652-163	<input type="checkbox"/> 65652-924	Single-Canister Roll Stand (1 each)
<input type="checkbox"/> 65652-181	<input type="checkbox"/> 65652-182	—	Double-Canister Roll Stand (1 each)
—	<input type="checkbox"/> 65652-164	<input type="checkbox"/> 65652-164	Metal Roll Stand, Single (1 each)
—	<input type="checkbox"/> 65652-264	<input type="checkbox"/> 65652-264	Metal Roll Stand with Regulator Hardware, Single (1 each)
—	<input type="checkbox"/> 65652-183	<input type="checkbox"/> 65652-183	Metal Roll Stand, Double (1 each)
—	<input type="checkbox"/> 65652-283	<input type="checkbox"/> 65652-283	Metal Roll Stand with Regulator Hardware, Double (1 each)
Guardian LVC			
		<input type="checkbox"/> 65651-120	65651-120 Guardian LVC (4/case)
		<input type="checkbox"/> 65652-583	65652-583 Guardian LVC Roll Stand (30") (1 each)
		<input type="checkbox"/> 65652-598	65652-598 Guardian LVC Roll Stand (18") (1 each)
Accessories			
		<input type="checkbox"/> 516500	Bulk Connecting Tubing (100' roll/box, 4 boxes/case)
		<input type="checkbox"/> 65652-542A	CRD™ and Guardian Tandem Tube (30/case)
		<input type="checkbox"/> 65652-122	Standard Specimen Sock (10/case)
		<input type="checkbox"/> 65652-123	Short Specimen Sock (10/case)
		<input type="checkbox"/> 65652-132	Short Specimen Sock Without Stick (10/case)
		<input type="checkbox"/> 65652-124	Large-Port Specimen Sock (10/case)
		<input type="checkbox"/> 65652-169	Female Specimen Sock (10/case)
		<input type="checkbox"/> 65652-119	Pour/Accessory Port Convertor (10/case)
		<input type="checkbox"/> 65652-120	Pour/Accessory Port Convertor/Large (10/case)
		<input type="checkbox"/> 65651-700	Universal Critical Measure (50/case)
		<input type="checkbox"/> 65652-107	Right-Angle Adapter for CRD and Guardian Patient Ports (50/case)
		<input type="checkbox"/> 65652-109	Right-Angle Adapter for Flex Advantage® Patient Port and CRD and Guardian Ortho Ports (50/case)
Universal Hardware			
		<input type="checkbox"/> 65652-002	Two-Canister Manifold (1 each)
		<input type="checkbox"/> 65652-003	Three-Canister Manifold (1 each)
		<input type="checkbox"/> 65652-004	Four-Canister Manifold (1 each)
		<input type="checkbox"/> 65652-580	In-Line Vacuum On/Off Valve (1 each)
		<input type="checkbox"/> 65652-587	Roll Stand, Short (1 each)
		<input type="checkbox"/> 65652-500	Roll Stand, Medium (1 each)
		<input type="checkbox"/> 65652-586	Roll Stand, Tall (1 each)
		<input type="checkbox"/> 65652-597	Roll Stand, Short with Regulator Mount (1 each)
		<input type="checkbox"/> 65652-510	Roll Stand, Medium with Regulator Mount (1 each)
		<input type="checkbox"/> 65652-596	Roll Stand, Tall with Regulator Mount (1 each)
		<input type="checkbox"/> 65652-686	Tall Eight-Canister Dual-Head Roll Stand (1 each)
		<input type="checkbox"/> 65652-104	Pigtail Connector (1 each)
		<input type="checkbox"/> 65652-108A	Tubing Barb (50/case)
		<input type="checkbox"/> 65652-543	Tee Connection (1 each)
		<input type="checkbox"/> 65652-131	Universal Female D.I.S.S. (1 each)
		<input type="checkbox"/> 65652-137	Male D.I.S.S. (Metal) (1 each)
		<input type="checkbox"/> 65652-145	Wall Plate (1 each)
		<input type="checkbox"/> 65652-134	Variable Height Adapter (Low Boy) (1 each)
		<input type="checkbox"/> 65652-149	Regulator Bracket (1 each)
SAF-T Pump™ System and Accessories			
		<input type="checkbox"/> 65652-000	SAF-T Pump (1 each)
		<input type="checkbox"/> 65652-008	Dip Tube, 8" (50/case)
		<input type="checkbox"/> 65652-011	Dip Tube, 11" (50/case)
		<input type="checkbox"/> 65652-016	Dip Tube, 16" (50/case)
		<input type="checkbox"/> 65652-024	Pump Tube, 24" (10/case)
		<input type="checkbox"/> 65652-044	Pump Tube, 44" (10/case)
		<input type="checkbox"/> 65652-001	Dispenser for Dip Tubes (1 each)

Medi-Vac® Products

To order, please contact your Cardinal Health sales representative or call customer service at 1.800.964.5227.

© Copyright 2004. Cardinal Health, Inc. or one of its subsidiaries.
All rights reserved.
Lit. No. MV00101 (0104/7.5M/MC1666)

 Printed on recycled paper with 10% post-consumer fiber.



Cardinal Health
Medi-Vac® Products
1500 Waukegan Road
McGaw Park, IL 60085

www.cardinal.com/mps