



ESSENTIAL. EFFECTIVE. EFFICIENT.

[Home](#) [Products](#) [Kidney & Stone Management](#) [Accessories](#) [Ellik Bladder Evacuator](#)

# Prescriptive Information

[read more](#)

# Training and Education

For additional training and education on a number of topics, visit our Training Center

[read more](#)

# CONTACT US

[read more](#)

# Ellik Bladder Evacuator



The BARD® Ellik Bladder Evacuator, available in disposable or reusable models, is designed for ease of maneuverability during the removal of tissue sections during transurethral prostatectomy, TURPs, bladder biopsy, etc.

# Features and Specifications

Description	Part Number	Size	Quantity	Natural Latex?
Wound Accessories Ellik Bladder Evacuator	000451	-	10/box	Yes

<b>Wound, Accessories, Ellik Bladder Evacuator</b>	004501	-	1 each	Yes
<b>Wound, Accessories, Ellik Bladder Evacuator, Glass Reservoir Only</b>	0450A	-	1 each	No
<b>Wound, Accessories, Ellik Bladder Evacuator, Medium Wall, Rubber Bulb</b>	0450M	-	1 each	Yes
<b>Wound, Accessories, Rubber Bulb, Medium Wall</b>	004501	-	1 each	Yes

Specifications are subject to change without notice.



## Bard Medical is excited to announce the acquisition of Rochester Medical

**Bard Medical and Rochester Medical are coming together, combining quality, innovative products and services to better serve Clinician and Patient needs.**



For information about Bard Medical Home Care products. [Click here](#)



For information about MAGIC 3® Hydrophilic Catheters. [Click here](#)

For Rochester Medical product orders and all other inquiries, please continue to contact Rochester Medical at 800.243.3315 or visit [www.rocm.com](http://www.rocm.com)

Copyright © 2013 C. R. Bard, Inc. All rights reserved. All products on this website may not be available for sale in all geographies. Please contact your local Bard Medical representative for information about products available in your area.

This site is intended for healthcare professionals. If you are a patient seeking more information, please consult your healthcare provider.