

The Power of Purple*



BAIRD
ACCESS SYSTEMS

*Advancing the Delivery of Health Care.**

Nursing Full Procedure Tray

Benefits Include:

- **Optimized component layout for clinician convenience during PICC placement procedure**
 - Tray is organized in user-friendly fashion designed to follow PICC insertion procedural steps
- **Standardized PICC line tray kits and components**
 - Provide familiarity with Bard Access Systems' PICC trays across all PICC lines
- **Needle-stick safety components**
 - Enhances clinician/patient safety and facilitates hospital compliance with federal safety regulations
- **Latex-free components**
 - Provides options for clinicians/patients with latex allergies




Multiple Configurations:

Single-Lumen	Dual-Lumen	Triple-Lumen
4 Fr.	5 Fr.	5 Fr.
5 Fr.	6 Fr.	6 Fr.


Features and Advantages:


- **Maximum injection rates** – allow injection of contrast media for CECT scans at a maximum rate of 5 ml/sec
- **Easy identification** – catheters are clearly labeled to differentiate them from all non-indicated power injection PICCs
- **Easy to read ID tags** – indicate maximum injection rates and prompt users to check catheter patency
- **CVP monitoring** – all PowerPICC* catheters are indicated for CVP monitoring
- **Patient satisfaction** – helps to improve patient comfort due to the potential for decreased needle sticks
- **Versatile vascular access device** – a single device that can be used for multiple therapies, power injection of contrast media, and CVP monitoring





PowerPICC* catheters combine the efficacy of PICC access with the efficiency of power injection for CECT scans





Ordering & Product Information

Product Specifications:

Type of Catheter	Product Code	Tray	Gravity Flow (ml/hr)	Priming Volume (ml)	Gauge
4 Fr. SL PowerPICC* Catheter	3174108 3174118	Full w/MI Basic w/MI	1,272	0.67	18
5 Fr. SL PowerPICC* Catheter	3175108 3175118	Full w/MI Basic w/MI	1,185	0.66	18
5 Fr. DL PowerPICC* Catheter	3275108 3275118	Full w/MI Basic w/MI	578/578	0.57/0.57	18/18
5 Fr. TL PowerPICC* Catheter	3385108 3385118	Full w/MI Basic w/MI	982/131/131	0.76/0.43/0.43	18/20/20
6 Fr. DL PowerPICC* Catheter	3276108 3276118	Full w/MI Basic w/MI	753/753	0.62/0.62	18/18
6 Fr. TL PowerPICC* Catheter	3386108 3386105	Full w/MI Basic w/MI	1,163/275/275	0.76/0.47/0.47	17/19/19

PowerPICC* Tray Components:

Description	Full Tray with MI	Basic Tray with MI
PICC	X	X
Microintroducer	X	X
Guidewire, 0.018 in.	X	X
Scalpel, safety	X	X
IV catheter, safety, 20 Ga. x 4.5 cm	X	X
Introducer needle, 21 Ga. x 7 cm	X	X
Needle hypodermic, safety, 25 Ga. x 0.63 in.	X	
Filter straw	X	
Wound closures (pkg of three)	X	
Transparent dressing	X	
Drape, fenestrated	X	
Towel, absorbent	X	
Poly underdrape	X	
Tourniquet, latex-free	X	
Face mask	X	
Surgical tape	X	
Scissors	X	
Measuring tape	X	
End cap (one per lumen)	X	X
Syringe, 12 ml luer-lock, latex-free	X	
Surgical gloves, sterile, latex-free	X	
Chloraprep* applicator	X	
Lidocaine	X	
Saline, 10 ml	X	
Alcohol wipe	X	
Gauze pad 4 in. x 4 in.	X	
Gauze pad 2 in. x 2 in.	X	
StatLock* catheter stabilization device	X	X

The Power of Purple*

Please consult product labels and inserts for any indications, contraindications, hazards, warnings, cautions and instructions for use.

Learn More - For more information on the PowerPICC* Catheter, please contact your Bard Access Systems representative or visit www.powerpicc.com. Call toll-free at 1-800-545-0890, or visit us on the web at www.bardaccess.com.

BARD

ACCESS SYSTEMS

Bard Access Systems, Inc.

Salt Lake City, Utah 84116 USA
801-522-5000
www.bardaccess.com

Customer Service: 800-545-0890
Clinical Information: 800-443-3385