

Orthopaedic Catalogue

Index	2
Order Form	6
Section 1 - Lower Limb	7
Section 2 - Spinal	31
Section 3 - Upper Limb	43
Section 4 - Paediatric	53
Section 5 - Compression Therapy	59
Section 6 - General Products	69



3D Lite Sheet	72	B	Circulation Socks	65		
3/4 Orthotic	10		Civic Collar	41		
3/4 E-Fit Orthotic	8		Clavicle Strap	50		
A		Bariatric Wheelchair	77	Collar & Cuff - OAPL	51	
	Abdominal Binder	32	Basic Camp Wrist with Thumb	46	Comfort Zone Orthotics	10
	Achilles Heel Pad	81	Basic Camp Wrist	46	Commode Chair Budget	77
	Acromioclavicular Splint	52	Basic Mouthguard	78	Commode Chair Folding	77
	A-Flex Protective Helmet	41	Bebax Boot	55	Component Orthotics	10
	AFO Socks	65	Becker Shoulder Holster	51	Compression Shorts	28
	AFTR DC	21	BioFoam Foam Impression Box	76	Compression Stockings	60 - 64
	AFTR GEL	21	BioSkin Ankle Bracing	20	Contender 16" knee	23
	AFTR	20	BioSkin Back Skin	33	Control Orthotics	8
	Air Stirrup Ankle Brace	20	BioSkin Elbow Bracing	49	Cool Fit Cinch	32
	All Gel Care Kit	83	BioSkin Knee Bracing	23	Corset Fronts for Kydex Braces	35
	Aluminum Crutches	76	BioSkin Thigh Bracing	28	CT Spine	34
	Aluminum Forearm Crutches	76	Blue Line Bandage	64	Cushion Ball of Foot Pad	81
	Amsterdam Footwear	12	Blue Rocker	16	Cushion Digital Cap	79
	Ankle Brace with Straps	19	Bobarth Sling	51	Cushion Orthotics	8
	Ankle Bracing	19	Body Disc	82	Cushion Toe Comb	80
	Ankle Foot Orthosis (AFO"S)	14	Boot Bumper	81	Cushion Toe Separators	80
	WalkAide	14	Braces - Mouthguard	78	Cushion Toe Spreader	80
	Dynamic Walk	15	Budget Airliner	18	Cut 4 Custom LSO	36
	Dictus Band	15	Budget Walking Stick	76	Cut 4 Custom TLSO	38
	Ypsilon	16	Budget Wheelchair	77	Cyberspine CO	40
	Blue Rocker	16	Bunion Gel Shield	84	Cyberspine TLSO	38
	Toe Off	16	Bunion Gel Sleeve	80	Cybertech 652 SI Belt	32
	Navigait	17	Bunion Sleeve	84	Cybertech CTO	39
	Supra-Lite	17	Bunionette Gel Shield	84	Cybertech Minerva	39
	Springleaf	17	C	Cybertech Mini Minerva	39	
	EFB Polyprop	17				
HeelPro	17	Calf Skin	25	D		
Heel Lift	17	Calibrated Midfoot				
Plantar Fasciitis Sock	17	Compression Wrap	21	Delta Dry	75	
OAPL Night Splint	17	California Ant Thoracic Extension	38	Deltacast Elite	74	
RCAI MPO Active	18	California CEO	38	Deluxe Rigid Post Op Shoe	11	
Ankle Walkers	18	California Compressor	34	Deluxe Wheelchair	77	
OAPL Standard	18	California Eco	37	Dennis Browne Bar	54	
Premium Walker	18	California Hip Orthosis	29	Diabetic Crew Socks	65	
Low Top Walker	18	California Mid Profile	35	Diabetic Insole	10	
Full Shell	18	California Soft Spinal	35	Dictus Band	15	
Anterior Closure Knee Wrap	23	Camp Combo	27	Digistrip	79	
Anterior Overlap Brace	36	Camp Corsets	34	Digital Cap with Toe Spreader	79	
Anti Embolism Stockings	60	Carpal Gel Pads	82	Digital Caps	79	
Arch Additions	9	Cast Boot	11	Digital Care Kit	83	
Arch Support Pair	80	Casting Products	72 - 76	Digital Dot Pads	84	
Arm Abduction Pillow	50	Cervical Collar Philli Style	40	Digital Gel Cap	84	
Arm Elevator	51	Cervical Roll	42	Digital Pads on Strip	79	
Arm Sling - Economy	51	Charnley Pillow	30	Digital Toe Prop	84	
Arm Sling Universal Immobilising	51	Check knee	25	Digitube Plus	79	
Arthritic Gel Sock	82	Child Dorsi Stretch AFO	54	Dorsal Carpal Tunnel Splint	47	
Asthma Products	78	Childrens Orthotics	8	Double Hernia Support	29	
Athletic Support	28					

DP2 Wrist Skin	47	Full Length Orthotic	10	Hinged Knee Support	23
DP3 Wrist Skin	47	Functional Resting Splint	47	Hip Flexion Assist Device	29
Dual Density Orthotics	8	G		Horton Click Bar	55
Dual Foam Orthotics	8	Gel Bun	80	Hot & Cold Therapy	78
Dynamic Sport Orthotic	10	Gel Care Advanced	83	I	
Dynamic Walk	15	Gel Cushions	81	Ice Pack	78
E		Gel E Roll	82	Indelible Pencils	76
EFB Polyprop AFO	17	Gel Fit Bun Sleeve	80	Industrial Lifting Belt	33
Elastic Ankle	19	Gel Fit Metatarsal Pad	81	Invisible Forefoot Cushion	85
Elastic Elbow	49	Gel Foot Cover	80	Invisible Gel Dots	85
Elastic Fig8 Ankle	19	Gel Max Mouthguard	78	Invisible Heel Cushion	85
Elastic Knee X Action	22	Gel Met Pad with Toe Spread	81	Invisible Heel Guard	85
Elastic Knee	22	Gel Metatarsal Pad	81	It Stays Glue	65
Elastic Maternity Belt	32	Gel Products	79 - 85	J	
Elastic Wrist Strap	45	Gel Squares	82	Jewitt Hyperextension Brace	38
Elastic Wrist Support	45	Gel Thong Protector	85	K	
Electric Plaster Cutting Saw	76	Gel Toe Prop	79	KAFO Socks	65
Essential JS Boots	54	Gel Toe Separators	80	Kassel Footwear	14
Even Up	19	Gel Toe Spreaders	80	KiddieGait	54
Exercise Pulley Set	42	Geriatric Hand	47	Knee Brace with Hyperext	22
Extended Length Humeral Fracture Brace	50	Gipsan Plaster Bandage	74	Knee Braces	22
EZ Arm Abduction Orthosis	50	Gladiator DT	25	Knee Cage Swedish	25
F		Gladiator XT	25	Knee Hinged	22
Fillauer Indelible Pencils	75	Gutter Cushion	27	Knee Skin with Straps	23
Finger Splint Ext with Wrist Support	44	Gypsona Extra 3	74	Knee Skin	23
Finger Splint Extension with Clockspring	44	H		Knee Spiral	22
Finger Splint Finger Knuckle Bender	44	H/S Mod Oppenheimer w/Thumb	48	Kydex Chairback Brace	35
Finger Splint Knuckle Bender Reverse	44	H/S Thomas Suspension Wen4	48	Kydex Knight Taylor Brace	37
Finger Splint Mallet	44	H/S Weniger 1 Cock-Up Splint	48	L	
Finger Splint Safety Pin Modified Wen 6A	44	Hallux Control Strap	21	Lace-Up Ankle Brace	20
Finger Splint Safety Pin Wen 6	44	Hallux Valgus Night Splint	21	Ladtofa Compression Stockings	63
Fleecy Web Extra	75	Hand Splint comp Spring & Elastic	44	Lastofa Cotton Compression Stockings	64
Flexifast	64	Hand Splint Knuckle Bender Flexion	44	Leipzig Footwear	13
Flexion Control Wrist	45	Hand Splint knuckle Bender Reverse	44	Lerman 3 Point Knee	25
Foam Cuff & Sling	51	Hans Exerciser	52	Long Shoe Horn	76
Foam Impression Box	76	Heel Cup	10	Low Profile SI Belt	32
Folding Walking Frame with Wheels	77	Heel Lift	17	Low Top Walker	18
Folding Walking Frame	77	Heel Pad	10	Lower Hernia Support	28
Folding Walking Stick	76	Heel Wedges	9	Lumbar Roll	42
Footbionic Orthotics	8 - 9	HeelPro	17	Lumbar Sacral Abdominal Bib	35
Forefoot Additions	9	Heidelberg Footwear	13	Lumbar Sacral Brace	35
Four Point Walking Stick	76	Helping Hand Euro Style	76	Lumbar Sacral Strap Kit	35
Fulflex Tubular Bandage	64	Hernia Support	28	Lympha Press Mini	68
		High Spinal Brace	37	Lympha Press Optimal	68
		Hinged Knee Skin Front Closure	25	Lymphedema Treatment	68 - 70
		Hinged Knee Skin	24		

M

Mackie Elbow	50
Mackie Knee	26
Maple Leaf Orthosis	56
Marseille Footwear	12
Maternity Belts	32
Max Compression Stockings	65
Mego Afek Products	67
Memory Compression Stockings	64
Metatarsal Pads	9
Minerva Collar	39
Mixture Adhesive Felt	75
Mouthguard Case	78
Mouthguards	78
MPO Active	18

N

Navigait	17
Neoprene Shorts	28
Neutral Position Resting	48
Newport 3 Hip System	29
Newport 4 Hip System	29
Newport Jnr Adjustable Abduction Bar	56
Newport Jnr Paediatric Orthosis	55
Night Resting Splint with Thumb	48
Night Splint - OAPL	17

O

OAPL Anti Embolism Stockings	60
OAPL Cast Boot	11
OAPL Collar & Cuff	51
OAPL Night Splint	17
OAPL Orthotics	10
OAPL Post Op Shoe	11
OAPL ROM Elbow	49
OAPL ROM Knee	26
OAPL Shoulder Immobiliser	51
OAPL Standard Walker	18
OAPL Travel Socks	60
OFA Bamberg Stockings	63 - 65
Order Form	6
Orthotykes Pavlik Harness	55
Overdoor Traction Kit	42

P

Paed Abduction Bar	55
Paediatric AFO	54
Paediatric Casting Sock STS	74
Paediatric leg & Elbow Wrap	57
Paediatric Walker	54
Patella Knee Stabiliser	23
Patella Stabiliser - Bioskin	24
Patella Tendon Strap	22
Pavlik Harness	55
Peak Flow Meter	78
Philli Collar	40
Phlebo Press DVT	67
PhysiFeet Invisible Gel	85
PhysiPod Gel Products	84
PhysiPod Orthotics	10
Plantar Fasciitis Sock	17
Plaster Shears Bruns	76
Plastic Collar with Chin	41
Plastic Collar	41
POP Arm Brace	50
Post Op Knee Immobiliser	26
Post Operative Carpal	82
Post Surgery Abdomen	83
Post Surgery Breast	83
Power Ultra Mouthguard	78
Prefabricated Orthotics	8 - 11
Premium Airliner	18
Premium Ankle Support	19
Premium Back	32
Premium Elbow Support	49
Premium Knee & Patella Support	22
Premium Knee Stabiliser	22
Premium Knee Support	22
Premium Thigh Support	28
Premium Walker	18
Premium Wrist	45
Professional Orthotic	9

Q

Q Baby	24
Q Brace Front Closure	24
Q Brace	24
Q Lok System	24
Quattromed	26

R



RCAI MPO Active	19
RCAI Paed Abduction Bar	55
RCAI POP Arm Brace	50
Rearfoot Additions	9
Regensburg Footwear	12
Replacement Crutch Tips	77
Rhino Brace	55
Rib/Abdominal Binder	32
Rollator Walker	77
ROM Elbow	49
ROM Knees	26 - 27

S

Sacro Cinch	33
SAE Brace Hinged	38
Sequential Compression Garments	67 - 70
Sequential Compression Pumps	67 - 70
Shoe Horn	76
Shoulder Abduction Sling	51
Shoulder Brace-Pipeer	50
Shoulder Immobiliser	51
Shower Chair for Baths	77
Shower Chair	77
Shower Stool	77
Silicone Masks	78
Silipos Products	79 - 83
Slimfit cushion Orthotic	9
Sling - Bobarth	51
Socks	65
Soffban Natural	75
Soffban Synthetic	75
Soft Collar	40
Soft Humeral Fracture Brace	50
Somi Neck Brace	39
Space Chamber	78
Specialist Plaster	74
Springleaf AFO	17
Sroufe Super Shoe II	12
Standard Ankle Skin	20
Standard Knee Skin	23
Stockholm Footwear	12
Stockings Compression	60 - 65
Strap Arm Sling	51
STS Paediatric Socks	75
STS Slipper Socks	74
STS Socks	75
Supercushion Orthotics	8

Supra-Lite AFO	17	Travel Socks	60	Weniger 13-Combined Coiled Oppenheimer	48
Swash Low Profile	57	Trilok Ankle	20	Weniger 1A Cock-Up Splint Outrigger	48
Swash Original	56	Tro Camp Belt	32	Weniger 2A Knuckle Bender	45
T		Tubular Gauze Bandage	64	Weniger 2B Knuckle Bender Splt	45
Tennis Elbow Skin	49	U		Weniger 9A-Rev Knuckle Bender	45
Tennis Elbow Strap	49	Ultima ROM Knee	27	Wheelchairs	77
Therafim Stockings	62	Universal Plantar Fasciitis Support	21	Wondercup	11
Therapeutic Footwear	11 - 14	University Collar	40	Wondersole	11
Thermic Ankle Support	19	UNO-WHO Thumb Spica	46	Wondersport	11
Thermic Elbow	49	UNO-WHO Wrist Brace	46	Wonderspur	11
Thermic Knee	22	V		Wrist and/or Ankle Restraints	52
Thermic Wrist Support	45	V Force Knee Brace	25	Wrist Braces	44 - 48
Thigh Skin	28	Vantage AP Control ROM Knee	27	Wrist/Hand Immobiliser	47
Three Point Walking Stick with Stool	76	Visco Knee Skin	24	Wrist/Hand/Thumb Immobiliser	47
Thumb Spica	46	W		Wrist/Thumb Spica	46
Toe Alignment Splint	21	WalkAide	14	X Y Z	
Toe Off	16	Webly Ankle	19	X-Lite Products	72 - 73
Toilet Seat Raiser	77	Weil Exerciser	21	Ypsilon	16
Tourniquet	52	Weil Osteomy Strap	21		
Traction Socks	65	Weniger 13B Comb Oppenh Composite Elastic	48		
Trading Terms	86				

OAPL have a range of catalogues available, if you don't see what you're looking for here please refer to our:



In addition to this OAPL are constantly sourcing the latest products to add to our range so contact your local representative or OAPL Customer Service today.



Lower Limb

**Prefabricated
Orthotics** 8

**Therapeutic
Footwear** 11

AFO's 15

Ankle Walkers 18

Ankle Bracing 19

Knee Support 22

Thigh Support 28

Hips 29

CONTROL

The baseline model using a balanced EVA and PE blend for excellent foot control and cushioning. The slimmest of this range is easy to fit and great in dress footwear. Available with the option of Black (polo-knit) or Tan (brushed nylon) top covers.



Code	Size	Colour
FB-CONB-XXS	XXSMALL	BLACK
FB-CONB-XS	XSMALL	BLACK
FB-CONB-S	SMALL	BLACK
FB-CONB-M	MEDIUM	BLACK
FB-CONB-L	LARGE	BLACK
FB-CONB-XL	XLARGE	BLACK
FB-CONB-XXL	XXLARGE	BLACK
FB-CONT-XXS	XXSMALL	TAN
FB-CONT-XS	XSMALL	TAN
FB-CONT-S	SMALL	TAN
FB-CONT-M	MEDIUM	TAN
FB-CONT-L	LARGE	TAN
FB-CONT-XL	XLARGE	TAN
FB-CONT-XXL	XXLARGE	TAN

DUAL DENSITY

The control model with a soft, comfortable, cushioning top layer of 2mm Ortholite. Perfect for most sporting activities and work boots, combining support with good cushioning. Immediately comfortable with no break-in period.



Code	Size
FB-DD-XXS	XXSMALL
FB-DD-XS	XSMALL
FB-DD-S	SMALL
FB-DD-M	MEDIUM
FB-DD-L	LARGE
FB-DD-XL	XLARGE
FB-DD-XXL	XXLARGE

3/4 EZ-FIT

The contoured underside allows for a nice fit into nearly all dress shoes you will encounter. This contour also means that the orthotic can be modified easily unlike a flat based device. The orthotic is also contoured to closely follow the angled shank of the shoe. The arch profile is not aggressive however still provides good overall foot support. A deep heel cup is retained despite the profile being low enough to fit into most dress shoes. Ez-Fit orthoses are made using discrete black foam and are covered in black brushed nylon for a fresh feel.



Code	Size
FB-EZF-XS	XSMALL
FB-EZF-S	SMALL
FB-EZF-M	MEDIUM
FB-EZF-L	LARGE
FB-EZF-XL	XLARGE

DUAL FOAM

EVA/PE base covered with a 2mm of soft, comfortable EVA. No fabric cover means this durable model can be used in wet or rugged conditions such as football, hockey, school shoes, work boots or gumboots.



Code	Size
FB-DF-XXS	XXSMALL
FB-DF-XS	XSMALL
FB-DF-S	SMALL
FB-DF-M	MEDIUM
FB-DF-L	LARGE
FB-DF-XL	XLARGE
FB-DF-XXL	XXLARGE

SUPERCUSHION

The superior cushioning and support at a value for money price, you can not look past this model. Supercushion is covered from tip to toe with 3mm medical grade Poron on a soft, durable base. Ideal for diabetic and sensitive foot types.



Code	Size
FB-SCUSH-XXS	XXSMALL
FB-SCUSH-XS	XSMALL
FB-SCUSH-S	SMALL
FB-SCUSH-M	MEDIUM
FB-SCUSH-L	LARGE
FB-SCUSH-XL	XLARGE
FB-SCUSH-XXL	XXLARGE

CUSHION

A softer accommodating foam blend including a soft 2mm Ortholite top cover. Designed where additional cushioning required for conditions such as arthritis, plantar pressure lesions and plantar fasciitis.



Code	Size
FB-CUSH-XXS	XXSMALL
FB-CUSH-XS	XSMALL
FB-CUSH-S	SMALL
FB-CUSH-M	MEDIUM
FB-CUSH-L	LARGE
FB-CUSH-XL	XLARGE
FB-CUSH-XXL	XXLARGE

FOOTBIONICS SIZING		
Size	Men (US)	Women (US)
XXS	1-2.5	3 - 4.5
XS	3 - 5.5	5 - 6.5
S	6 - 7.5	7 - 8.5
M	8 - 9.5	9 - 10.5
L	10 - 11.5	11 - 12.5
XL	12 - 13.5	
XXL	14 - 16	

*US sizing is 1/2 a size larger than AUS

SLIMFIT CUSHION

Using the exact same base foam without our 2mm Ortholite top-layer, Footbionics Cushion Slim is only 3.5mm thick at the forefoot and thinner overall in profile through the arch and the heel. Covered with soft brush nylon material for added comfort and a fresh feel, Footbionics Cushion Slim feel good and are surprisingly supportive. The Cushion Slim fits comfortably into virtually any shoe.



Code	Size
FB-SFC-XXS	XXSMALL
FB-SFC-XS	XSMALL
FB-SFC-S	SMALL
FB-SFC-M	MEDIUM
FB-SFC-L	LARGE
FB-SFC-XL	XLARGE
FB-SFC-XXL	XXLARGE

PROFESSIONAL

The most advanced model integrates three different foam densities including a firm stability cradle extending from heel to the forefoot and covered with 2mm Ortholite. A fantastic device for grossly pronated feet, or for heavier patients.



Code	Size
FB-PRO-XXS	XXSMALL
FB-PRO-XS	XSMALL
FB-PRO-S	SMALL
FB-PRO-M	MEDIUM
FB-PRO-L	LARGE
FB-PRO-XL	XLARGE
FB-PRO-XXL	XXLARGE

CHILDRENS

This is a new range of orthoses for children. They have fantastic shape that is well tolerated and are simple to fit and adjust. There are three sizes for children aged 3 to 11 years. Available in Red or Blue.



Code	Size	Shoe Size	Colour
FB-CHB-S	SMALL	CHILD 7 - 9	RED
FB-CHB-M	MEDIUM	CHILD 9 - 11	RED
FB-CHB-L	LARGE	CHILD 11 - 13	RED
FB-CHR-S	SMALL	CHILD 7 - 9	BLUE
FB-CHR-M	MEDIUM	CHILD 9 - 11	BLUE
FB-CHR-L	LARGE	CHILD 11 - 13	BLUE

ARCH ADDITIONS

To compliment the rear foot posts these 3D shaped arch additions attach neatly to Footbionics orthotics to provide additional midtarsal control and arch support. These are nearly seamless when glued in place to Footbionics, no need to finish with a grinder.



Code	Size	Unit
FB-AA-S	SMALL	Pk 5 Pairs
FB-AA-M	MEDIUM	Pk 5 Pairs
FB-AA-L	LARGE	Pk 5 Pairs

REARFOOT ADDITIONS

Constructed from firm density pressed foam that produces an accurate 3D shape that closely resembles a completed 6 degree rear foot post. Designed for use specifically for Footbionics, you'll love using these time saving devices.



Code	Size	Unit
FB-RFA-S	SMALL	Pk 5 Pairs
FB-RFA-M	MEDIUM	Pk 5 Pairs
FB-RFA-L	LARGE	Pk 5 Pairs

FOREFOOT ADDITIONS

Constructed from firm density pressed foam that produces an accurate 3D shape that closely resembles a completed 6 degree fore foot post. Designed for use specifically for Footbionics, you'll love using these time saving devices.



Code	Size	Unit
FB-FFA-S	SMALL	Pk 5 Pairs
FB-FFA-L	LARGE	Pk 5 Pairs

HEEL WEDGES

These additions compliment our range. They are constructed from 8mm medium density EVA. They are ideal for leg length discrepancies and reducing loading on the Achilles and calf. No self adhesion, must be glued.



Code	Size	Unit
FB-HR-S	SMALL	Pk 10
FB-HR-M	MEDIUM	Pk 10
FB-HR-L	LARGE	Pk 10

METATARSAL PADS

Self adhesive, available in both soft and firm densities, and in two convenient sizes. Our met pads are compression moulded and cleverly shaped to fit the metatarsal parabola.



Code	Size	Density	Unit
FB-MPS-S	SMALL	SOFT	Pk 10
FB-MPF-S	SMALL	FIRM	Pk 10
FB-MPS-L	LARGE	SOFT	Pk 10
FB-MPF-L	LARGE	FIRM	Pk 10

PHYSIPOD COMPONENT ORTHOTICS

The component PhysiPod Orthotic has been developed after many years of podiatric practitioner experience, to offer a fully customisable Orthotic to suit multiple foot and lower limb complaints. All health professionals can design the Orthotic to suit their patient needs with the components supplied. The kit provided offers multiple combinations of mechanical and palliative additions for tailoring the Orthotic.



Code	Size
AH-CCO-S	SMALL
AH-CCO-M	MEDIUM
AH-CCO-L	LARGE
AH-CCO-XL	XLARGE

PHYSIPOD FULL LENGTH ORTHOTICS

The PhysiPod Full Length Orthotic has been designed to biomechanically correct foot alignment and offer maximum support to the foot throughout the gait cycle. Using the latest technology and highest quality materials the orthotic offers optimum control and comfort. Ideal for the treatment of; flat feet and over pronation, plantar fasciitis, heel spur, metatarsalgia, ball of foot pain, Morton's neuroma, inner knee pain, shin splints and lower back pain. Comes with an unique arch filler for added arch support.



Code	Size
AH-FLO-S	SMALL
AH-FLO-M	MEDIUM
AH-FLO-L	LARGE
AH-FLO-XL	XLARGE

PHYSIPOD 3/4 LENGTH ORTHOTICS

The PhysiPod 3/4 Length Orthotic has been designed to biomechanically correct foot alignment and offer maximum support to the foot throughout the gait cycle. Using the latest technology and highest quality materials the orthotic offers optimum control and comfort. Ideal for the treatment of; flat feet and over pronation, plantar fasciitis, heel spur, metatarsalgia, Morton's neuroma, inner knee pain, shin splints and lower back pain. Comes with an unique arch filler for added arch support.



Code	Size
AH-TQO-S	SMALL
AH-TQO-M	MEDIUM
AH-TQO-L	LARGE
AH-TQO-XL	XLARGE

PHYSIPOD HEEL CUP


This range is made from the highest quality medical grade silicone and mineral oil gel for durability, comfort, support and lasting wear.



Code	Size
AH-SHC-M	MEDIUM
AH-SHC-L	LARGE
AH-SHC-XL	XLARGE

OAPL DYNAMIC SPORT ORTHOTIC

Designed and developed by a Podiatrist to biomechanically correct and stabilise the active foot. Using the latest combination of premium materials, the Dynamic Sport Orthotic offers optimum control and comfort for the most demanding sporting requirements. Featuring; Deep, contoured heel cup with Thermoplastic heel guard for added stabilisation. TRI ZONE SUPPORT SYSTEM, featuring EVA orthotic base, PU gel strike zones in the heel and forefoot and ultra soft latex arch cradle. This system controls foot motion through the gait cycle and further cushions the heel and forefoot by absorbing high impact forces during sports. The arch is cushioned and supported by the generous contoured arch filler. Advanced antimicrobial top cover controls moisture and friction.



Code	Size
AH-DSO-S	SMALL
AH-DSO-M	MEDIUM
AH-DSO-L	LARGE
AH-DSO-XL	XLARGE

OAPL COMFORT ZONE ORTHOTIC

Designed and developed by a Podiatrist to biomechanically correct and stabilise the active foot. Using the latest combination of premium materials, the Comfort Zone Orthotic offers optimum control and enhanced comfort and support. Featuring, Soft Polyurethane shell and base. Flexible TPU orthotic laminated in the Polyurethane shell, allowing increased support and dynamic flex for shock absorption. Premium foam layer covers the orthotic providing added cushioning and comfort. Advanced antimicrobial mesh top cover control moisture and friction.



Code	Size
AH-SZO-S	SMALL
AH-SZO-M	MEDIUM
AH-SZO-L	LARGE
AH-SZO-XL	XLARGE

OAPL ORTHOTICS SIZING

Size	UKMens	UKWomens	US Mens	US Womens
Small	5-6	6-7	5-7	7-8
Medium	7-8	8-9	8-9	9-10
Large	9-10	10-11	10-11	10-11
X Large	11-12	12-13	12-13	11-12

OAPL DIABETIC INSOLE

Podiatrist designed insole specifically for the Diabetic, Arthritic and Sensitive Foot. The OAPL Diabetic insole is made of a Tri-Laminate construction, using a medical grade Plastazote top cover, soft foam middle layer and a soft Polyurethane base. The medical grade Plastazote top cover provides superior self moulding properties to evenly contour and distribute weight across the foot, reducing friction and aiding in impact absorption. The soft foam central layer further adds to cushioning and a reduction of shearing and friction forces, needed in protecting the at risk foot. The soft Polyurethane base further offers longitudinal arch support, further shock absorption and cushioning.



Code	Size
AH-DBO	UNIVERSAL

CUSHION HEEL PAD

Helps relieve heel pain. Self adhesive backing for proper placement. Reduces discomfort caused by plantar fasciitis. Packaged in a polybag.

Code	Size
1592	SMALL
1593	MEDIUM
1594	LARGE



WONDERSPUR

To aid calcaneal spurs, plantar fasciitis, heel spurs, knee and back pain. Dual density orthotics. Maximum support and shock absorption using two densities of medical grade silicone. Clear silicone with a softer density blue dot, provides additional relief to sensitive areas.

Code	Size	Shoe Size
AH-4500	SMALL	M 2 - 4 W 4 - 5
AH-4501	MEDIUM	M 5 - 7 W 6 - 9
AH-4502	LARGE	M 8 - 10 W 10 - 12
AH-4503	XLARGE	M 11 & Up W 12 & Up



WONDERCUP

To relieve pain throughout the heel. Dual density orthotics. Maximum support and shock absorption using two densities of medical grade silicone. Clear silicone with a softer density of red dot, provides additional relief to sensitive areas.

Code	Size	Shoe Size
AH-4504	SMALL	M 2 - 4 W 4 - 5
AH-4505	MEDIUM	M 5 - 7 W 6 - 9
AH-4506	LARGE	M 8 - 10 W 10 - 12
AH-4507	XLARGE	M 11 & Up W 12 & Up



WONDERSPORT

To cushion both heel, arch and metatarsal. Dual density orthotics. Maximum support and shock absorption using two densities of medical grade silicone. Clear silicone with a softer density green dot, provides additional relief to sensitive areas.

Code	Shoe Size
AH-4508	M 5 - 7
AH-4509	M 8 - 9
AH-4510	M 10 - 11
AH-4511	M 12 - 13
AH-4512	M 14 & Up
AH-4513	W 4 - 5
AH-4514	W 6 - 9
AH-4515	W 10 - 12
AH-4516	W 12 & Up



WONDERSOLE

The full length insole designed to provide support to metatarsals, arch and heel. Helps to relieve symptoms related to metatarsalgia, plantar fasciitis, fat pad atrophy and arthritis. Packaged in a box.

Code	Size	Shoe Size	Unit
AH-4530	SMALL	M: 7 W: 5-7	1 PAIR/PKG
AH-4531	MEDIUM	M: 8-9 W: 8-9	1 PAIR/PKG
AH-4532	LARGE	M: 10 W: 10-11	1 PAIR/PKG
AH-4533	XLARGE	M: 11-13	1 PAIR/PKG



OAPL CAST BOOT

To protect plaster cast from soil and moisture and provide secure, non-slip footing. A canvas sandal with rubber rocker sole and velcro closure.

Code	Size	Length
11079B	PAEDIATRIC	
11079A	XSMALL	180mm
11079	SMALL	240mm
11080	MEDIUM	280mm
11081	LARGE	320mm
11082	XLARGE	360mm



OAPL POST OP SHOE

To protect foot and dressings post operatively and provide secure non-slip footing. A nylon shoe with an enclosed heel counter and velcro closures. A cushioned sole provides additional comfort.

Code	Size	Variation	Length
11082C	SMALL	WOMENS	240mm
11082D	MEDIUM	WOMENS	255mm
11082E	LARGE	WOMENS	270mm
11082F	SMALL	MENS	265mm
11082G	MEDIUM	MENS	285mm
11082H	LARGE	MENS	300mm
11082I	XLARGE	MENS	310mm



DELUXE RIGID POST OP SHOE WITH COLLAR

Features a rigid sole for total stabilisation. Heel counter provides firm support and the heel collar is padded to reduce slipping. Loop-lock and D-Ring closure adjust to accommodate dressings.

Code	Size	Variation	Length	Width
445-12	SMALL	WOMENS	241mm	84mm
445-13	MEDIUM	WOMENS	254mm	88mm
445-14	LARGE	WOMENS	266mm	92mm
445-22	SMALL	MENS	279mm	98mm
445-23	MEDIUM	MENS	292mm	101mm
445-24	LARGE	MENS	304mm	104mm



SROUFE SUPER SHOE II

Designed for wear after soft tissue procedures or post-operative care. Especially useful for geriatric and diabetic patients with insensate feet. Adjustable straps accommodate the bulkiest dressings. Removable trim-fit tongue provides a better fit. Semi-rigid moulded non-skid sole absorbs shock during ambulation.



Code	Size	Variation
442-12	SMALL	WOMENS
442-13	MEDIUM	WOMENS
442-14	LARGE	WOMENS
442-22	SMALL	MENS
442-23	MEDIUM	MENS
442-24	LARGE	MENS

MARSEILLE

Designed to completely offload the heel the Marseille is an ultra-light shoe, with a fully lined, nylon upper and a rigid EVA sole that allows for easy modification. The sole has a 5 degree wedge, without hind foot support and an exposed heel. It also has increased midfoot support to further offload the heel. Protective forefoot and heel caps are available.



Code	Size	Shoe Size	Length	Width
FO-2120-S	SMALL	35 - 37	245mm	88mm
FO-2120-M	MEDIUM	38 - 40	260mm	92mm
FO-2120-L	LARGE	41 - 42	275mm	98mm
FO-2120-XL	XLARGE	43 - 46	295mm	103mm

PROTECTIVE HEEL CAP

Suitable for the Marseille Shoe.

Code	Size
FO9320-S	SMALL
FO9320-M	MEDIUM
FO9320-L	LARGE
FO9320-XL	XLARGE



STOCKHOLM

An ultra light shoe, with a fully lined, nylon upper and a rigid EVA sole that allows for easy modification. The sole has an 8 degree wedge cutout to off load the fore foot.



Code	Size	Shoe Size	Length	Width
FO-1120-XS	XSMALL	30 - 32	210mm	83mm
FO-1120-S	SMALL	33 - 35	230mm	87mm
FO-1120-M	MEDIUM	36 - 38	250mm	91mm
FO-1120-L	LARGE	39 - 41	270mm	97mm
FO-1120-XL	XLARGE	42 - 44	290mm	102mm
FO-1120-XXL	XXLARGE	45 - 47	310mm	109mm

PROTECTIVE FOREFOOT CAP

Suitable for the Stockholm and Marseille.

Code	Size
FO9121-XS	XSMALL
FO9121-S	SMALL
FO9121-M	MEDIUM
FO9121-L	LARGE
FO9121-XL	XLARGE
FO9121-XXL	XXLARGE



AMSTERDAM

The Amsterdam has been designed as a functional post-op shoe. It features a strapping system high on the dorsum and forefoot for optimal suspension, rigid rocker sole to facilitate ambulation and removable insole.



Code	Size	Shoe Size
TS6100-XS	XSMALL	33 - 35
TS6100-S	SMALL	36 - 39
TS6100-M	MEDIUM	40 - 42
TS6100-L	LARGE	43 - 44
TS6100-XL	XLARGE	45 - 47
TS6100-XXL	XXLARGE	48 - 50

PROTECTIVE FOREFOOT CAP

Suitable for the Amsterdam shoe.

Code	Size
FO9621-XS	XSMALL
FO9621-S	SMALL
FO9621-M	MEDIUM
FO9621-L	LARGE
FO9621-XL	XLARGE
FO9621-XXL	XXLARGE



REGENSBURG

The Regensburg features an anti-microbial insole and rigid rocker sole. It has been designed to be worn with a external fixator, and can be trimmed to allow a customised fit as required. It can also be used as a regular post-op shoe.



Code	Size	Shoe Size	Length
TS6141-XS	XSMALL	33 - 35	220mm
TS6141-S	SMALL	36 - 39	240mm
TS6141-M	MEDIUM	40 - 42	260mm
TS6141-L	LARGE	43 - 44	280mm
TS6141-XL	XLARGE	45 - 47	300mm
TS6141-XXL	XXLARGE	48 - 50	330mm

LEIPZIG

Designed as a extra depth rehabilitation shoe the Leipzig has nylon fully lined upper, a flexible sole and 3 removable insoles. The upper insole is latex and both lower insoles are EVA. The patient can wear the shoe with wound dressings as the volume of the shoe can be adjusted via the touch tape fasteners. There is also a zipper in the Achilles tendon area for easy donning and doffing of the shoe.



Code	Shoe Size	Orientation
TS5512-35L	35	LEFT
TS5512-36L	36	LEFT
TS5512-37L	37	LEFT
TS5512-38L	38	LEFT
TS5512-39L	39	LEFT
TS5512-40L	40	LEFT
TS5512-41L	41	LEFT
TS5512-42L	42	LEFT
TS5512-43L	43	LEFT
TS5512-44L	44	LEFT
TS5512-45L	45	LEFT
TS5512-46L	46	LEFT
TS5512-47L	47	LEFT
TS5512-48L	48	LEFT

TS5512-35R	35	RIGHT
TS5512-36R	36	RIGHT
TS5512-37R	37	RIGHT
TS5512-38R	38	RIGHT
TS5512-39R	39	RIGHT
TS5512-40R	40	RIGHT
TS5512-41R	41	RIGHT
TS5512-42R	42	RIGHT
TS5512-43R	43	RIGHT
TS5512-44R	44	RIGHT
TS5512-45R	45	RIGHT
TS5512-46R	46	RIGHT
TS5512-47R	47	RIGHT
TS5512-48R	48	RIGHT

HEIDELBERG

Designed as a extra depth rehabilitation shoe the Heidelberg has nylon fully lined upper with neoprene top, a flexible sole and 3 removable insoles. The upper insole is latex and both lower insoles are EVA. Due to its nylon construction and the removable insoles, the shoe is ideal for use with patients with wounds or swollen painful feet.



Code	Shoe Size	Orientation
TS4552-35L	35	LEFT
TS4552-36L	36	LEFT
TS4552-37L	37	LEFT
TS4552-38L	38	LEFT
TS4552-39L	39	LEFT
TS4552-40L	40	LEFT
TS4552-41L	41	LEFT
TS4552-42L	42	LEFT
TS4552-43L	43	LEFT
TS4552-44L	44	LEFT
TS4552-45L	45	LEFT
TS4552-46L	46	LEFT
TS4552-47L	47	LEFT
TS4552-48L	48	LEFT

TS4552-35R	35	RIGHT
TS4552-36R	36	RIGHT
TS4552-37R	37	RIGHT
TS4552-38R	38	RIGHT
TS4552-39R	39	RIGHT
TS4552-40R	40	RIGHT
TS4552-41R	41	RIGHT
TS4552-42R	42	RIGHT
TS4552-43R	43	RIGHT
TS4552-44R	44	RIGHT
TS4552-45R	45	RIGHT
TS4552-46R	46	RIGHT
TS4552-47R	47	RIGHT
TS4552-48R	48	RIGHT

KASSEL

Designed as a definitive or wound care shoe the Kassel features enough room inside, enabling the patient to wear a larger dressing if necessary. It has a bonded upper, the firm outer material able to be separated from the soft inner material and an opening cut to relieve pressure. The extra rigid rocker sole can be re-ground, to customise the rocker, to better immobilise or relieve lesions. The EVA insert (40 shore A) has a thickness of 2cm in the forefoot area. Special pressure relief areas can be ground into the material. A Forefoot Protective cap is available as an optional accessory.



Code	Shoe Size	Orientation
TS5326-35L	35	LEFT
TS5326-36L	36	LEFT
TS5326-37L	37	LEFT
TS5326-38L	38	LEFT
TS5326-39L	39	LEFT
TS5326-40L	40	LEFT
TS5326-41L	41	LEFT
TS5326-42L	42	LEFT
TS5326-43L	43	LEFT
TS5326-44L	44	LEFT
TS5326-45L	45	LEFT
TS5326-46L	46	LEFT
TS5326-47L	47	LEFT
TS5326-48L	48	LEFT

TS5326-35R	35	RIGHT
TS5326-36R	36	RIGHT
TS5326-37R	37	RIGHT
TS5326-38R	38	RIGHT
TS5326-39R	39	RIGHT
TS5326-40R	40	RIGHT
TS5326-41R	41	RIGHT
TS5326-42R	42	RIGHT
TS5326-43R	43	RIGHT
TS5326-44R	44	RIGHT
TS5326-45R	45	RIGHT
TS5326-46R	46	RIGHT
TS5326-47R	47	RIGHT
TS5326-48R	48	RIGHT

PROTECTIVE FOREFOOT CAP

Code	Size
F09520-XXS	XXSMALL
F09520-XS	XSMALL
F09520-S	SMALL
F09520-M	MEDIUM
F09520-L	LARGE
F09520-XL	XLARGE
F09520-XXL	XXLARGE



WALKAIDE

The WalkAide is a battery-operated, single-channel electrical stimulator that can be used to address foot drop with functional electrical stimulation. WalkAide utilises a Tilt Sensor to control stimulation during normal gait which is programmed by the clinician using WalkAide software.
Indications: Cerebrovascular accident, incomplete spinal cord injury, traumatic brain injury, multiple sclerosis, cerebral palsy.
Contra-Indications: Heart disease, epilepsy, recent surgery. Patient Kit includes: WalkAide control unit, WalkAide Electrode Lead Cable, WalkAide Electrodes, pk 4.



Code	Size	Variation	Calf Circ.
20-1000		PATIENT KIT	
20-240B	XSMALL	PAEDIATRIC CUFF	190 - 254mm
20-0250	SMALL	BI FLEX CUFF	254 -381mm
20-0260	MEDIUM	BI FLEX CUFF	297 - 432mm
20-0270	LARGE	BI FLEX CUFF	305 - 483mm
20-0400	LEFT	PATIENT FOOT SENSOR	
20-0401	RIGHT	PATIENT FOOT SENSOR	
20-0290	Pk 10	DISPOSABLE LINER	

PREMIUM ELECTRODES

20-0520	2" SQUARE	Pk 4
20-0530	2" SQUARE	Ctn 10 Pk 4
20-0521	1.25" ROUND	Pk 4
20-0531	1.25" ROUND	Ctn 10 Pk 4
20-0540	1.875" ROUND	Pk 4
20-0550	1.875" ROUND	Ctn 10 Pk 4

DYNAMIC WALK

Made from a lightweight mouldable carbon fibre composite material with an open heel construction. This AFO is low profile and the high activity design allows for running, walking and squatting. Its dorsiflexed position enables the dynamic walk to maintain its position on the foot/leg dynamically, without the need for straps. Suitable for flaccid foot drop with minimal ankle instability. Wide ankle version available (Non stocked item).

Code	Size	Shoe Size	Variation	Length	Orientation
700111121	SMALL	M 3.5 - 6 W 4 - 6.5	DUAL UPRIGHT	230mm	LEFT
700111141	MEDIUM	M 6.5 - 8.5 W 7 - 9	DUAL UPRIGHT	250mm	LEFT
700111161	LARGE	M 10 - 13.5 W 10.5 - 14	DUAL UPRIGHT	270mm	LEFT
700113121	SMALL	M 3.5 - 6 W 4 - 6.5	MEDIAL UPRIGHT	230mm	LEFT
700113141	MEDIUM	M 6.5 - 8.5 W 7 - 9	MEDIAL UPRIGHT	250mm	LEFT
700113161	LARGE	M 10 - 13.5 W 10.5 - 14	MEDIAL UPRIGHT	270mm	LEFT
700114121	SMALL	M 3.5 - 6 W 4 - 6.5	LATERAL UPRIGHT	230mm	LEFT
700114141	MEDIUM	M 6.5 - 8.5 W 7 - 9	LATERAL UPRIGHT	250mm	LEFT
700114161	LARGE	M 10 - 13.5 W 10.5 - 14	LATERAL UPRIGHT	270mm	LEFT
700111122	SMALL	M 3.5 - 6 W 4 - 6.5	DUAL UPRIGHT	230mm	RIGHT
700111142	MEDIUM	M 6.5 - 8.5 W 7 - 9	DUAL UPRIGHT	250mm	RIGHT
700111162	LARGE	M 10 - 13.5 W 10.5 - 14	DUAL UPRIGHT	270mm	RIGHT
700113122	SMALL	M 3.5 - 6 W 4 - 6.5	MEDIAL UPRIGHT	230mm	RIGHT
700113142	MEDIUM	M 6.5 - 8.5 W 7 - 9	MEDIAL UPRIGHT	250mm	RIGHT
700113162	LARGE	M 10 - 13.5 W 10.5 - 14	MEDIAL UPRIGHT	270mm	RIGHT
700114122	SMALL	M 3.5 - 6 W 4 - 6.5	LATERAL UPRIGHT	230mm	RIGHT
700114142	MEDIUM	M 6.5 - 8.5 W 7 - 9	LATERAL UPRIGHT	250mm	RIGHT
700114162	LARGE	M 10 - 13.5 W 10.5 - 14	LATERAL UPRIGHT	270mm	RIGHT
700100220	SMALL	PADDING KIT			
700100240	MEDIUM	PADDING KIT			
700100260	LARGE	PADDING KIT			



Ankle Foot Orthosis

DICTUS BAND

Made from a sturdy leather construction, with latex bands (stronger and longer lasting than elastic) to provide a dynamic dorsiflexion moment and superior foot lift. The latex band can be set in various positions enabling directional pull to combat equino varus foot/ankle posture. Supplied with both the steel hook and Monohook providing patient's with greater footwear options. Suitable for flaccid footdrop with minimal ankle instability.

Code	Size	Ankle Circ.
530200-002	CHILD	120 - 170mm
530500-003	MEDIUM	150 - 200mm
530500-004	LARGE	200mm & Up

A mid-foot cuff compatible with the Dictus band for shoeless dorsiflexion assist. This item is an accessory - it can only be purchased separately and must be used in conjunction with the Dictus band.

530600-003	UNIVERSAL	Shoeless Upgrade	LEFT
530600-004	UNIVERSAL	Shoeless Upgrade	RIGHT

Additional or replacement latex bands can be purchased separately. The Monohook is an option for anchoring the band to the shoe where eyelets are not available.

530100-041	Latex Bands	EACH
530100-003	Steel Hooks	PAIR
530100-002	Monohook	



YPSILON

Carbon fibre and Kevlar design. Low profile single strap at proximal edge. Y shaped shell design to load either side of tibial crest. New reinforced padding for calf band. Designed to provide dynamic toe off assist, while allowing natural ankle movements for the stable ankle. Suitable for, foot drop resulting from neurological and neuromuscular deficit disorders, such as stroke, multiple sclerosis, post polio or post spinal cord injury, muscular dystrophy and post surgery complications. **A foot othosis MUST be adhered to the footplate of this device.**



Code	Size	Length	Orientation
28400-1011	SMALL	230mm	LEFT
28400-1012	MEDIUM	250mm	LEFT
28400-1013	LARGE	270mm	LEFT
28400-2011	SMALL	230mm	RIGHT
28400-2012	MEDIUM	250mm	RIGHT
28400-2013	LARGE	270mm	RIGHT

ACCESORIES

Code	Name	Size
28402L	YPSILON:EXTRA STRAP	LARGE
28402M	YPSILON:EXTRA STRAP	MEDIUM
28402S	YPSILON:EXTRA STRAP	SMALL
284120011	EXTRA STRAP/COIL	SML/MED
284120013	EXTRA STRAP/COIL	LRG/XL

BLUE ROCKER

The Blue Rocker is designed to be a stiffer version of the Toe Off. It's carbon fibre and kevlar open design allows to heel to perform its biomechanical role of shifting from eversion to inversion. The anterior plate assists with knee extension at the mid-stance phase of gait. Dynamic ankle foot orthoses designed to provide extra stability to improve both balance and posture. Developed for patients with more involved pathology in addition to footdrop. Suitable for, multiple involvement, foot drop resulting from neurological and neuromuscular deficit disorder, such as stroke, multiple sclerosis, post polio or post spinal cord injury, muscular dystrophy and patrial foot amputation. **A foot othosis MUST be adhered to the footplate of this device.**



Code	Size	Length	Orientation
284051011	SMALL	230mm	LEFT
284051012	MEDIUM	250mm	LEFT
284051013	LARGE	270mm	LEFT
284051014	XLARGE	285mm	LEFT
284052011	SMALL	230mm	RIGHT
284052012	MEDIUM	250mm	RIGHT
284052013	LARGE	270mm	RIGHT
284052014	XLARGE	285mm	RIGHT

TOE OFF

The ToeOff's carbon fibre and Kevlar open design allows to heel to perform its biomechanical role, shifting from eversion to inversion. The anterior plate assists with knee extension at the mid-stance phase of gait. Its dynamic design provides extra stability to improve both balance and posture. Developed for patients with more involved pathology in addition to footdrop. Suitable for, multiple involvement, foot drop resulting from neurological and neuromuscular deficit disorders, such as stroke, multiple sclerosis, post polio or post spinal cord injury, muscular dystrophy and patrial foot amputation. **A foot othosis MUST be adhered to the footplate of this device.**



Code	Size	Foot Plate	Orientation
28375-1010	XSMALL	210mm	LEFT
28375-1011	SMALL	230mm	LEFT
28375-1012	MEDIUM	245mm	LEFT
28375-1013	LARGE	270mm	LEFT
28375-2010	XSMALL	210mm	RIGHT
28375-2011	SMALL	230mm	RIGHT
28375-2012	MEDIUM	245mm	RIGHT
28375-2013	LARGE	270mm	RIGHT

TOE OFF COMFORT KIT

Code	Size
28398XS	XSMALL
28398S	SMALL
28398M	MEDIUM
28398L	LARGE

TOE OFF SOFTKIT LINER

A softkit liner for the Toeoff AFO. Available in small, medium and large.

Code	Size
28390S	SMALL
28390M	MEDIUM
28390L	LARGE

TOE OFF SOFTSHELL LINER

A softshell liner for the Toeoff AFO. Available in all sizes.

Code	Size	Colour
28393XS	XSMALL	BEIGE
28393S	SMALL	BEIGE
28393M	MEDIUM	BEIGE
28393L	LARGE	BEIGE

NAVIGAIT

Low profile design and padded shell make this a discrete and comfortable alternative for AFO users with mild foot drop. The amount of dorsiflexion is easily adjusted with an elastic coil to accommodate the users specific needs. The longer design provides better suspension for the device. Position device above the Malleolus.

Code	Size	Height
284100011	SML/MED	275mm
284100013	LRG/XL	315mm



SUPRA-LITE

The supralite's lightweight open heel design will limit plantar flexion, provide some heel stabilisation without inhibiting dorsiflexion. The cut-outs incorporated in this design provide a gentle spring action and allow the heel to rest comfortably in the shoe.

Code	Size	Length	Orient
13020	SML/MED	240mm	LEFT
13022	MED/LRG	280mm	LEFT
13021	SML/MED	240mm	RIGHT
13023	MED/LRG	280mm	RIGHT



SPRINGLEAF AFO

The Springleaf AFO limits plantar flexion whilst allowing dorsiflexion. A polypropylene posterior shell with velcro closure proximally. Often requires adjustment with heat gun at time of fitting.

Code	Size	Shoe Size	Orientation
8250	SMALL	M 3.5 & Below	LEFT
8251	MEDIUM	M 4 - 6	LEFT
8252	LARGE	M 6 - 8.5	LEFT
8253	XLARGE	M 9 & Up	LEFT
8255	SMALL	M 3.5 & Below	RIGHT
8256	MEDIUM	M 4 - 6	RIGHT
8257	LARGE	M 6 - 8.5	RIGHT
8258	XLARGE	M 9 & Up	RIGHT



EFB POLYPROP

To maintain the foot and ankle in a neutral position restricting both plantar flexion and dorsiflexion. Anterior trim lines provide mediolateral stability and some subtalar joint control. A polypropylene posterior shell with Velcro closure proximally. Often requires adjustment with heat gun at time of fitting.

Code	Size	Shoe Size	Orient
13001	SMALL	M 3.5 & Below	LEFT
13003	MEDIUM	M 4 - 6	LEFT
13005	LARGE	M 6 - 8.5	LEFT
13007	XLARGE	M 9 & Up	LEFT
13002	SMALL	M 3.5 & Below	RIGHT
13004	MEDIUM	M 4 - 6	RIGHT
13006	LARGE	M 6 - 8.5	RIGHT
13008	XLARGE	M 9 & Up	RIGHT



HEEL LIFT

Available in a smooth or convoluted (original) foam. The Heelift Suspension Boot is designed to provide an effective solution for prevention and assistive healing of pressure ulcers at the heel, while delivering positioning and alignment of the lower limb. It is made from latex free medical grade foam with a friction free tricot base to ensure it stays in place at all times as well as improving client mobility. Comes with an extra pad to control hip rotation, foot drop or provide added elevation. One size fits all.

Code	Size
11022	ORIGINAL
11023	CONVOLUTED



Ankle Foot Orthosis

HEELPRO

A foot and ankle positioning orthosis designed to unweigh the heel and protect the foot from potential ulcers and or loss of dorsiflexion range of motion. The unique design suspends the heel and reduces pressure on the toes from bed linen to facilitate comfort and protect skin. Designed primarily for patients spending long periods of time in recumbent position, the HeelPRO can improve clearance in swing phase and promote a neutral position during stance for limited ambulation.

Code	Size	Shoe Size	Foot Lgth	Calf
3530	SMALL	M 3 - 4.5	200 - 230mm	410mm
3531	MEDIUM	M 4.5 - 6	230 - 250mm	410mm
3532	LARGE	M 7.5 - 10	250 - 280mm	460mm
3533	XLARGE	M 10.5 - 13.5	280 - 310mm	510mm
3529.02	OPTIONAL DE-ROTATION BAR KIT			



PLANTAR FASCIITIS SOCK

The OAPL Plantar Fasciitis socks is an easy to use, inexpensive yet highly effective device. By keeping tension on the plantar fascia tissue at night it helps reduce the impact associated with initial stretch under body weight in the morning. The unique proximal tibial cuff reduces migration down the limb and broad strap attachment spreads load across the entire forefoot. Used alone or in conjunction with other treatment modalities, it can also help to decrease recovery time.

Code	Size	Shoe Size
WI-PFSS	SMALL	M Up to 6.5 W Up to 7.5
WI-PFSR	REGULAR	M 6.5 - 9 W 7.5 - 10.5
WI-PFSL	LARGE	M 9+ W 10.5+



OAPL NIGHT SPLINT

To immobilise ankle in plantar grade or slightly dorsiflexed position during rest. Lycra lined foam padded straps with centre release buckles. Foam lined kydex chassis with EVA walking sole and oblique dorsiflexion straps, suitable for transfers and short distance ambulation. Indicated for plantar fasciitis and patients at risk of losing dorsiflexion range.

Code	Size	Shoe Size
50023E	SMALL	M Up to 6.5 W Up to 7.5
50025E	MEDIUM	M 6.5 - 9 W 7.5 - 10.5
50027E	LARGE	M 9+ M 10.5+



RCAI MPO ACTIVE 2000

Dynamically corrects foot and ankle contractures, provides heel clearance for pressure relief. Features a moisture wicking liner, hook and loop closures. The rotator bar controls hip and leg rotation in a supine position. Bilateral Multi Pods Systems with a spreader bar provide controlled hip abduction. The low profile ambulatory attachment pad maintains proper hip alignment. Converts to a bed-positioning device with "twist and lock" removal of ambulatory attachment.



Code	Size	Shoe Size	Calf Circ.
O2MP-S	SMALL	178 - 203	318
O2MP-M	MEDIUM	203 - 229	432
O2MP-L	LARGE	229 - 254	432
O2MP-XL	XLARGE	254 - 279	483

The RCAI Spreader Bar provides 35° of controlled hip abduction, controls hip rotation and maintains proper alignment

12SB-P	PAEDIATRIC	SPREADER BAR
12SB-A	ADULT	SPREADER BAR

To increase the calf and foot circumference, especially useful in the case of edema.

O2EXT-S	SMALL	EXTENDER
O2EXT-M/L	MEDIUM/LARGE	EXTENDER
O2EXT-XL	XLARGE	EXTENDER

OAPL STANDARD WALKER

Fabricated using a combination of plastic shells with adjustable aluminium uprights, this walker is ideal where shape customisation is required.

Code	Size
OA0200	XSMALL
OA0400	SMALL
OA0600	MEDIUM
OA0800	LARGE
OA1000	XLARGE



PREMIUM WALKER

Comprising of a moulded plastic construction, this light-weight walker has been designed with a low heel height and smooth rocker sole profile for optimal forward progression.

Code	Size
EWK010A-17S	SMALL
EWK010A-17M	MEDIUM
EWK010A-17L	LARGE
EWK010A-17XL	XLARGE



LOW TOP WALKER

Identical in design to the OAPL Premium Walker, the OAPL Low Top features shorter uprights and is best suited to stable forefoot pathology.

Code	Size
EWK010A-11S	SMALL
EWK010A-11M	MEDIUM
EWK010A-11L	LARGE
EWK010A-11XL	XLARGE



OAPL FULL SHELL WALKER

The OAPL Full Shell Walker provides full foot and ankle protection with the added benefit of pneumatic support. This Walker features a rocker sole for added comfort and ease of ambulation as well as a wide foot base to aid comfort.

Code	Size
OAFSS	SMALL
OAFSM	MEDIUM
OAFSL	LARGE



STANDARD LINER

Replacement foam liner for OAPL Walkers.

Code	Size
OA0200LINER	XSMALL
OA0400LINER	SMALL
OA0600LINER	MEDIUM
OA0800LINER	LARGE
OA1000LINER	XLARGE



BUDGET AIRLINER

Inflatable liner with removable pump and valve.

Code	Size
OA0400BAIR	SMALL
OA0600BAIR	MEDIUM
OA0800BAIR	LARGE
OA1000BAIR	XLARGE



PREMIUM AIRLINER

Inflatable liner with in-built pump and valve.

Code	Size
OA0400AIR	SMALL
OA0600AIR	MEDIUM
OA0800AIR	LARGE
OA1000AIR	XLARGE



ACHILLES WEDGE

Heel raise consisting of 1cm layers laminated together. Can be used in walkers post surgery.

Code
45618

\$40.00 each for purchases of 10 or more, Standard Walkers, Premium Walkers & Low Top Walkers (can include multiple sizes and styles)

WALKER SIZING - ALL STYLES

Walker Size	Mens Shoe	Womens Shoe
X-Small	Up to 2.5	2 - 4
Small	3 - 5	4 - 6
Medium	5.5 - 8	6 - 9
Large	9 - 11	10 - 12
X-Large	11+	12+

WALKER SHELL

Anterior and posterior moulded plastic shells for additional control and support.

Code	Size
OA0400SH	SMALL
OA0600SH	MEDIUM
OA0800SH	LARGE



EVENUP

Evenup compensates for leg length changes while wearing an ankle walker or post op shoe. The Evenup solves the problem of back, hip, or knee stress by levelling the opposite shoe so that your patient will walk on an even platform, reducing body strain. Evenup has a removable lining to provide either 1.3cm or a 2.6cm shoe raise.

Code	Size	Length
EVEN-S	SMALL	254 - 286mm
EVEN-M	MEDIUM	292 - 324mm
EVEN-L	LARGE	330 - 368mm



ELASTIC ANKLE

To support and provide compression to the ankle. Pull-on elasticised anklet with heel cut-out.

Code	Size	Ankle Circ.
12007	SMALL	190 - 200mm
12008	MEDIUM	210 - 240mm
12009	LARGE	250 - 860mm
12010	XLARGE	270mm & Up



ELASTIC FIG8 ANKLE

To provide firm compression and support without restricting ankle motion. A Figure 8 design, pull on elasticised anklet with heel cut-out. Simulates bandaging without the bulk.

Code	Size	Ankle Circ.
12011	SMALL	190 - 200mm
12012	MEDIUM	210 - 240mm
12013	LARGE	250 - 860mm
12014	XLARGE	270mm & Up



THERMIC ANKLE SUPPORT

The OAPL Thermic Ankle Support, constructed from 3mm thick neoprene, provides support and warmth to promote rapid healing of soft tissue injuries. Maintaining body heat helps to increase blood flow promoting elasticity of muscles, tendons and ligaments. The mild compression actively reduces swelling as well as providing relief to arthritic and rheumatic pain.

Code	Size	Ankle Circ.
65015	SMALL	160 - 220mm
65016	MEDIUM	220 - 250mm
65017	LARGE	250 - 280mm
65018	XLARGE	280 - 320mm



PREMIUM ANKLE SUPPORT

The OAPL Premium Ankle Support provides compression and warmth to promote rapid healing of soft tissue injuries. Maintaining body heat helps to increase blood flow, promoting elasticity of muscles, tendons and ligaments. This ventilated neoprene support with moisture wicking technology and soft terry lining provides comfortable support which aids in reducing swelling and associated pain.

Code	Size	Ankle Circ.
65010	SMALL	230 - 260mm
65011	MEDIUM	270 - 300mm
65012	LARGE	310 - 340mm
65013	XLARGE	350 - 400mm



ANKLE BRACE WITH STRAPS

The OAPL Ankle Brace with Straps is constructed of heavy gauge ballistic nylon with non stretch, non fraying nylon straps. Featuring a neoprene tongue and Achilles relief padding to aid in comfort. Useful for the treatment and prevention of sprains and strains whilst providing compression to reduce swelling and restricting ankle inversion and eversion. A low profile support that allows good fit in shoes and boots.

Code	Size	Shoe Size
20611	XXSMALL	W 5 - 6
20612	XSMALL	M 3 - 4.5 W 5 - 6.5
20613	SMALL	M 5 - 6.5 W 7 - 8.5
20615	MEDIUM	M 7 - 8.5 W 9 - 10.5
20617	LARGE	M 9 - 10.5 W 11 - 12.5
20618	XLARGE	M 11 - 12.5 W 13 - 14.5



WEBLY ANKLE

The Webly ankle brace is designed to limit inversion/eversion without restricting plantar or dorsiflexion. Constructed of light-weight pack cloth with a NoSweat tongue for added flexibility. Featuring Achilles relief padding and are indicated either prophylactically or for rehabilitation of ankle injury.

Code	Size	Shoe Size
HW-306XXS	XXSMALL	W 5-6
HW-306XS	XSMALL	M 6-7 W 6-8
HW-306S	SMALL	M 7.5-9 W 8-10
HW-306M	MEDIUM	M 9.5-11 W 10-12
HW-306L	LARGE	M 11.5-13 W 12-14
HW-306XL	XLARGE	M 13.5-15



LACE-UP ANKLE BRACE

The OAPL Lace Up Ankle Brace is designed for the treatment of lateral ankle joint ligament instability and injury. This brace controls lateral movement without inhibiting normal ankle joint motion. It provides protection from re-injury during sporting and recreational activities. Viable and economical alternative to strapping. Can be worn comfortably inside most shoes.



Code	Size	Shoe Size
12018F	XXSMALL	W 3 - 4.5
12018G	XSMALL	M 6 - 7 W 5 - 6
12018A	SMALL	M7.5 - 9 W 6.5 - 8
12018B	MEDIUM	M 9.5 - 11 W 8.5 - 10
12018C	LARGE	M 11.5 - 13 W 10.5 - 12
12018D	XLARGE	M 13.5 - 15 W 12.5 - 14
12018H	XXLARGE	M 15+

AIR STIRRUP ANKLE BRACE

Unique stirrup ankle brace with air bladder. Features adjustable heel pad for accurate positioning and optimal control during physical activity. Provides secure mediolateral ankle support.



Code	Size	Orientation
12021	UNIVERSAL	SMALL
12022	UNIVERSAL	REGULAR

BIOSKIN PREMIUM BRACING

The result of combining innovative designs, based on a firm understanding of biomechanics, with materials that have been designed specifically for orthopaedic and sports medicine applications. The patented BioSkin material provides a high level of compression without limiting activity or creating any discomfort. BioSkin is thin and gives radial compression so it conforms to the curves of the body. BioSkin is also breathable, so perspiration evaporates right through it and is hypo-allergenic eliminating the chance of skin irritation.

STANDARD ANKLE SKIN

To support and provide compression, added stability and increased proprioception to the ankle. Hook and loop rear closure allows easy application, micro-fleece lining with anti-microbial finish for extra comfort, universal left and right.



Code	Size
B53500	XSMALL
B53501	SMALL
B53502	MEDIUM
B53503	LARGE
B53504	XLARGE
B53505	XXLARGE

STANDARD ANKLE SKIN WITH FIGURE 8 WRAP

To support and provide compression, added stability and increased proprioception. Hook and loop rear closure allows easy application, micro-fleece lining with anti-microbial finish for extra comfort, universal left and right. With the addition of a figure 8 wrap to provide added ankle support.



Code	Size
B53510	XSMALL
B53511	SMALL
B53512	MEDIUM
B53513	LARGE
B53514	XLARGE
B53515	XXLARGE

TRILOK ANKLE

A wrap-around brace that achieves forefoot control and ankle stabilization. An ideal alternate to strapping. The Footlok Strap protects the ankle in its most vulnerable position, plantar flexion. The unique Footlok Strap duplicates taping for protection and may be reversed for management of Tib Post dysfunction. This strap controls foot motion (inversion and eversion) and foot positioning (pronation and supination).



Code	Size	Shoe Size
B53601	SMALL	M 4.5-6 W 4.5-6
B53602	MEDIUM	M 6.5-8 W 6.5-9
B53603	LARGE	M 9-12 W 9.5+
B53604	XLARGE	M 12.5-15

AFTR

The AFTR provides support and stabilisation for the non-operative and post-operative treatment of many foot and ankle pathologies. It includes medial and lateral plastic stays as well as figure-8 stirrup straps to limit movement in the ankle and hind foot. Made from Bio Skin material so it perfectly contours to the anatomy during normal volume changes following injuries and surgery. Bio Skin is hypo-allergenic, breathable and thin.



Code	Size
B53801	XSMALL/SMALL
B53803	MEDIUM/LARGE
B53805	XLARGE/XXLARGE

BIOSKIN ANKLE SIZING



Measure ankle circumference from Heel to top of Foot and back to Heel

Size	Cirumference
XSmall	250 - 280mm
Small	280 - 310mm
Medium	310 - 330mm
Large	330 - 360mm
XLarge	360 - 380mm
XXLarge	380 - 410mm

BIOSKIN AFTR SIZING



Measure ankle circumference from Heel to top of Foot and back to Heel

Size	Cirumference
XSmall/Small	279 - 330mm
Medium/Large	318 - 356mm
XLarge/XXLarge	343 - 381mm

AFTR DC

The AFTR DC (dorsal closure) provide support and stabilisation for the non-operative and post-operative treatment of many foot and ankle pathologies. It includes medial and lateral plastic stays as well as figure-8 stirrup straps to limit movement in the ankle and hind foot. Made from Bio Skin material it perfectly contours to the anatomy during normal volume changes following injuries and surgery. Bio Skin is hypo-allergenic, breathable and thin.



Code	Size
B53861	XSMALL/SMALL
B53863	MEDIUM/LARGE
B53865	XLARGE/XXLARGE

AFTR GEL

The AFTR Gel provides support and stabilisation for the non-operative and post-operative treatment of many foot and ankle pathologies. Featuring gel pads behind the malleoli to reduce swelling. It include medial and lateral plastic stays as well as figure-8 stirrup straps to limit movement in the ankle and hind foot. Made from Bio Skin material it perfectly contours to the anatomy during normal volume changes following injuries and surgery. Bio Skin is hypo-allergenic, breathable and thin.



Code	Size
B53831	XSMALL/SMALL
B53833	MEDIUM/LARGE
B53835	XLARGE/XXLARGE

CALIBRATED MIDFOOT COMPRESSION WRAP

The Calibrated Midfoot Compression Wrap allows you to use any of the Weil Digit Straps without covering the entire foot. The material is treated with a hypoallergenic non-slip finish on the inside to effectively hold it in place. The low profile design is well suited for use inside normal footwear.



Code	Unit
B95361	SINGLE
B95362	PK 10

HALLUX CONTROL STRAP

The HCS may be used post-bunionectomy to maintain the hallux in a neutral position during critical soft tissue healing. After primary healing, the HCS may be used day or night to stretch connective tissue and increase range of motion. The HCS can be positioned to provide resistance for FHL strengthening.



*Shown with Midfoot Compression Wrap, also compatible with AFTR DC, AFTR and AFTR Gel. These must be purchased separately.

Code	Size	Unit
B95391	SMALL	SINGLE
B95392	SMALL	PK 10
B95393	LARGE	SINGLE
B95394	LARGE	PK 10

WEIL OSTEOTOMY STRAP

The Weil Osteotomy Strap helps prevent dorsal drift and floating toes by splinting the toe in a neutral or plantar-flexed position. The Weil Osteotomy Strap can be worn in a walking boot, post-op shoe, or normal footwear and my be also used for non-operatively to treat metatarsalgia or plantar plate tears.



*Shown with Midfoot Compression Wrap, also compatible with AFTR DC, AFTR and AFTR Gel. These must be purchased separately.

Code	Unit
B95371	SINGLE
B95372	PK 10

WEIL EXERCISER

The Weil Exerciser is applied to the lesser digits to provide a home therapy program to strengthen the plantarflexors. Each set comes with one high tension and one low-tension exerciser to help patients regain normal FDL strength.



Code	Variation	Unit
B95383	LOW TENSION	SINGLE
B95384	LOW TENSION	PK 10
B95385	HIGH TENSION	SINGLE
B95386	HIGH TENSION	PK 10

HALLUX VALGUS NIGHT SPLINT

A plastic foam lined orthosis indicated as a night splint that will rigidly hold the Hallux in medial deviation. Suitable for Hallux Valgus or as a post-operative orthosis to maintain correction following surgery.



Code	Shoe Size	Orientation
OB-150501L	M 3-4.5 W 4.5-6	LEFT
OB-150502L	M 5-7.5 W 6.5-9	LEFT
OB-150503L	M 8-12 W 9.5-14	LEFT
OB-150501R	M 3-4.5 W 4.5-6	RIGHT
OB-150502R	M 5-7.5 W 6.5-9	RIGHT
OB-150503R	M 8-12 W 9.5-14	RIGHT

TOE ALIGNMENT SPLINT

Designed to maintain proper toe alignment post surgery. Consists of an elastic cotton metatarsal band which will fit any foot securely, along with T-Strap and lock straps for the toes.

Code	Size
DA-TAS	UNIVERSAL



UNIVERSAL PLANTAR FASCIITIS SUPPORT

Provides pressure and continuous support to help reduce the pain caused by plantar fasciitis. The 11/2" elastic tension strap applies pressure over the medial calcaneal tubercle and plantar fascia. Reducing the pain associated with inflammation due to plantar fasciitis, allowing the patient to exercise and maintain proper muscle tone. Provide relief day or night and can be worn with shoes, slippers, sandals - even barefoot. The breathable body wicks away moisture for enhanced comfort.



Code	Size
61-1000	SMALL/MEDIUM
61-1001	LARGE/XLARGE

ELASTIC KNEE

A cotton elastic pull-on knee support.

Code	Size	Knee Circ.
14067	SMALL	290 - 330mm
14068	MEDIUM	330 - 380mm
14069	LARGE	380 - 420mm
14070	XLARGE	420mm & Up



X ACTION ELASTIC KNEE

To provide support to the knee without restricting motion. An elasticised panel anteriorly with two crossover posterior elastic strips.

Code	Size	Knee Circ.
15000	SMALL	290 - 330mm
15001	MEDIUM	330 - 380mm
15002	LARGE	380 - 420mm
15003	XLARGE	420mm & Up



KNEE HINGED

A cotton elastic support with stainless steel, single-axis hinges and upper and lower encircling straps.

Code	Size	Knee Circ.
14019	SMALL	290 - 330mm
14020	MEDIUM	330 - 380mm
14021	LARGE	380 - 420mm
14022	XLARGE	420mm & Up



KNEE BRACE HINGED HYPEREXT

A hinged knee brace with posterior cross over strapping behind the knee and a hinge with an extension stop. This brace is designed to prevent hyperextension.

Code	Size	Knee Circ.
14036	SMALL	290 - 330mm
14037	MEDIUM	330 - 380mm
14038	LARGE	380 - 420mm
14039	XLARGE	420mm & Up



THERMIC KNEE

3mm Neoprene support, 2 colour reversible support. Ideal for post injury sporting applications, provides support and warmth.

Code	Size	Knee Circ.
65065	SMALL	300 - 350mm
65066	MEDIUM	350 - 400mm
65067	LARGE	400 - 450mm
65068	XLARGE	450mm & Up



PREMIUM KNEE SUPPORT

The OAPL Premium Knee Support provides compression and warmth to promote rapid healing of soft tissue injuries. Maintaining body heat helps to increase blood flow, promoting elasticity of muscles, tendons and ligaments. This ventilated neoprene support with moisture wicking technology and soft terry lining provides comfortable support which aids in reducing swelling and associated pain.

Code	Size	Knee Circ.
65060	SMALL	290 - 330mm
65061	MEDIUM	330 - 380mm
65062	LARGE	380 - 420mm
65063	XLARGE	420mm & Up



PREMIUM KNEE & PATELLA SUPPORT

The OAPL Premium Knee and Patella Support provides compression and warmth to promote rapid healing of soft tissue injuries. Maintaining body heat helps to increase blood flow, promoting elasticity of muscles, tendons and ligaments. This ventilated neoprene support with moisture wicking technology and soft terry lining provides comfortable support which aids in reducing swelling and associated pain.

Code	Size	Knee Circ.
65080	SMALL	290 - 330mm
65081	MEDIUM	330 - 380mm
65082	LARGE	380 - 420mm
65083	XLARGE	420mm & Up



PREMIUM KNEE STABILISER

The OAPL Premium Knee Stabiliser provides compression and warmth to promote rapid healing of soft tissue injuries. Maintaining body heat helps to increase blood flow, promoting elasticity of muscles, tendons and ligaments. This ventilated neoprene support with moisture wicking technology and soft terry lining. Provides comfortable support which aids in reducing swelling and associated pain.

Code	Size	Knee Circ.
65006	SMALL	290 - 330mm
65007	MEDIUM	330 - 380mm
65008	LARGE	380 - 420mm
65009	XLARGE	420 - 460mm



PATELLA TENDON STRAP

To aid in the correction of patella tracking abnormalities. A soft material strap with firm rubber anterior insert and velcro cinch strap.

Code	Size	Knee Circ.
14075	SMALL	Up to 350mm
14076	MEDIUM	355 - 375mm
14077	LARGE	380mm & Up



PATELLA KNEE STABILIZER

The brace exerts a medially directed force to the lateral border of the patella, as well as pressure on the patella tendon throughout the normal range of knee flexion. Low profile neoprene wrap with firm buttress and counter force D-Ring/Touch Tape.

Code	Size	Knee Circ.
14119	XSMALL	305 - 330mm
14120	SMALL	330 - 355mm
14121	MEDIUM	355 - 380mm
14122	LARGE	380 - 405mm
141221	XLARGE	405 - 460mm
141222	XXLARGE	460 - 510mm
141223	3XLARGE	510 - 560mm
141224	4XLARGE	560 - 610mm



CONTENDER 16"

To provide medial and lateral support for the treatment of knee joint instability. Features a hinge for medial/lateral support and vel-stretch proximal strap for maximum flexibility and increased patient comfort. Manufactured from premium grade latex-free neoprene with suede hinge cover.

Code	Size	Mid Patella Circ
88-8453	SMALL	330 - 355mm
88-8454	MEDIUM	355 - 380mm
88-8455	LARGE	380 - 405mm
88-8456	XLARGE	405 - 460mm
88-8457	XXLARGE	460 - 510mm
88-8458	3XLARGE	510 - 560mm
88-8459	4XLARGE	560 - 610mm



Knee Support

HINGED KNEE SUPPORT

To provide medial and lateral support for the treatment of knee instability and/or degenerative joint disease. Features a sleeve design with heavy duty hinges for optimal medial/lateral support. Manufactured from premium grade latex-free neoprene.

Code	Size	Mid Patella
88-0322	XSMALL	305 - 330mm
88-0323	SMALL	330 - 355mm
88-0324	MEDIUM	355 - 380mm
88-0325	LARGE	380 - 405mm
88-0326	XLARGE	405 - 460mm
88-0327	XXLARGE	460 - 510mm
88-0328	3XLARGE	510 - 560mm
88-0329	4XLARGE	560 - 610mm



BIOSKIN PREMIUM BRACING

The result of combining innovative designs, based on a firm understanding of biomechanics, with materials that have been designed specifically for orthopaedic and sports medicine applications. The patented BioSkin material provides a high level of compression without limiting activity or creating any discomfort. BioSkin is thin and gives radial compression so it conforms to the curves of the body. BioSkin is also breathable, so perspiration evaporates right through it and is hypo-allergenic eliminating the chance of skin irritation.

BIOSKIN KNEE SIZING

Measure knee circumference at mid-patella

Size	Circumference
XSmall	310 - 330mm
Small	330 - 360mm
Medium	360 - 410mm
Large	410 - 430mm
XLarge	430 - 460mm
XXLarge	460 - 540mm

ANTERIOR CLOSURE KNEE WRAP

Features a wrap around design for ease of application with heavy duty hinges for optimal medial/lateral support. Manufactured from premium grade latex-free neoprene.

Code	Size	Mid Patella
88-1303	SMALL	305 - 355mm
88-1304	MEDIUM	355 - 405mm
88-1305	LARGE	405 - 460mm
88-1306	XLARGE	460 - 510mm
88-1307	XXLARGE	510 - 560mm
88-1308	3XLARGE	560 - 610mm
88-1309	4XLARGE	610 - 660mm



STANDARD KNEE SKIN

For mild compression and increased proprioception. Available in open or closed patella, along with a unique strapping system (optional) that minimises migration.

Code	Size
B41100	XSMALL
B41101	SMALL
B41102	MEDIUM
B41103	LARGE
B41104	XLARGE
B41105	XXLARGE



STANDARD KNEE SKIN WITH STRAPS

For mild compression and increased proprioception. Available in open or closed patella, along with a unique strapping system (optional) that minimises migration.

Code	Size
B41110	XSMALL
B41111	SMALL
B41112	MEDIUM
B41113	LARGE
B41114	XLARGE
B41115	XXLARGE



VISCO KNEE SKIN

Knee ring provides intermittent passive massage for soft tissue. Ideal after knee arthroscopy. Finger loops enable easy pull-on application, along with a unique strapping system that minimises migration.

Code	Size
B41700	XSMALL
B41701	SMALL
B41702	MEDIUM
B41703	LARGE
B41704	XLARGE
B41705	XXLARGE



VISCO KNEE SKIN WITH STRAPS

Knee ring provides intermittent passive massage for soft tissue. Ideal after knee arthroscopy. Finger loops enable easy pull-on application, along with a unique strapping system that minimises migration.

Code	Size
B41710	XSMALL
B41711	SMALL
B41712	MEDIUM
B41713	LARGE
B41714	XLARGE
B41715	XXLARGE



UNDER SKIN

Locks brace in position, prevents chafing or rubbing.

Code	Size
B41400	XSMALL
B41401	SMALL
B41402	MEDIUM
B41403	LARGE
B41404	XLARGE
B41405	XXLARGE



PATELLA STABILIZER

Effective for patellar stabilization and treatment anterior knee pain. Interchangeable horseshoe or donut buttresses, along with a unique strapping system that minimises migration.

Code	Size
B41250	XSMALL
B41251	SMALL
B41252	MEDIUM
B41253	LARGE
B41254	XLARGE
B41255	XXLARGE



Q BABY

Perfect for the treatment of Osgood-Schlatter and patellar tendonitis, Silicone pad may be placed over patellar tendon or tibial tuberosity.

Code	Size
B51652	XSMALL-MEDIUM
B51655	LARGE-XXLARGE



Q BRACE

To provide continuous therapy and patella tracking. Unique brace offers multi-directional tracking to the patella and achieves greater dynamic force as the knee is flexed. Strapping systems minimizes migration.

Code	Size
B41610	XSMALL
B41611	SMALL
B41612	MEDIUM
B41613	LARGE
B41614	XLARGE
B41615	XXLARGE



Q-BRACE FRONT CLOSURE

To provide continuous therapy and patella tracking. Unique brace offers multi-directional tracking to the patella and achieves greater dynamic force as the knee is flexed. Strapping systems minimizes migration.

Code	Size
B51630	XSMALL
B51631	SMALL
B51632	MEDIUM
B51633	LARGE
B51634	XLARGE
B51635	XXLARGE



Q LOK SYSTEM

Reduces patellofemoral pain by increasing surface contact area with the joint. Provides a traction force on lateral retinaculum, reduces patella tilt, seats patella in trochlear notch. Wrap around design with lightweight hinges contoured to the leg. Made from ultra comfortable and breathable Stratus™ material. Universal left or right.

Code	Size
B61640	XSMALL
B61641	SMALL
B61642	MEDIUM
B61643	LARGE
B61644	XLARGE
B61645	XXLARGE



HINGED KNEE SKIN

To provide medial and lateral support to the knee. The lightest compression knee support available. Finger loops enable easy pull-on application. Light weight hinges reduces pressure at the knee joint. Donut and horseshoe buttresses can be placed in a position to provide custom patella support.

Code	Size
B41300	XSMALL
B41301	SMALL
B41302	MEDIUM
B41303	LARGE
B41304	XLARGE
B41305	XXLARGE



HINGED KNEE SKIN WITH FRONT CLOSURE

To provide medial and lateral support to the knee. The lightest compression knee support available. Finger loops enable easy pull-on application. Light weight hinges reduces pressure at the knee joint. Donut and horseshoe buttresses can be placed in a position to provide custom patella support.

Code	Size
B51330	XSMALL
B51331	SMALL
B51332	MEDIUM
B51333	LARGE
B51334	XLARGE
B51335	XXLARGE

**CHECK KNEE**

To prevent genu recurvatum, provide mediolateral stability and permit knee flexion. Long uprights reduce pressure on the femur and tibia. Contours to the shape of the leg. Unique strapping mechanism, lightweight, slim uprights allow full flexion.

Code	Size	Variation
14004	UNIVERSAL	Aluminium
14004SS	UNIVERSAL	Stainless Steel
28651	UNIVERSAL	Padding Kit

**GLADIATOR XT**

To provide support for mild moderate ACL, PCL, MCL and LCL deficiencies. Pivot-point level lock design stabilises knee, easily adjustable ROM setting. Extension 0-90 degrees, flexion 0-90 degrees. Features custom moulded soft condylar padding.

Code	Size
B41340	XSMALL
B41341	SMALL
B41342	MEDIUM
B41343	LARGE
B41344	XLARGE
B41345	XXLARGE

**KNEE CAGE SWEDISH**

To prevent genu recurvatum, provide mediolateral stability and yet permit knee flexion. One piece aluminium frame with webbing straps. Recurvatum is controlled by providing an anteriorly directed force behind the knee and two posteriorly directed forces above and below the knee.

Code	Size
14005	CHILD
14006	ADULT

**LERMAN 3 POINT KNEE**

To apply a varus valgus corrective force to the knee. Single size is fully adjustable and can be used left or right. Thigh and calf cuffs have integral adhesion layers designed to prevent migration of the orthosis.

Code	Size
200444	UNIVERSAL

**GLADIATOR DT**

To provide support for mild moderate ACL, PCL, MCL and LCL deficiencies. Pivot-point level lock design stabilises knee, easily adjustable ROM setting. Extension 0-90 degrees, flexion 0-90 degrees. Features custom moulded soft condylar padding. Gladiator DT can be used for mid-level sports activities. Front Closure.

Code	Size
B51360	XSMALL
B51361	SMALL
B51362	MEDIUM
B51363	LARGE
B51364	XLARGE
B51365	XXLARGE

**V FORCE KNEE BRACE**

Streamlined limited motion hinge design. Custom or prefabricated models available. Easily adjusted unloading to cope with your patients changing requirements. Low maintenance with self lubricating hinge.

Code	Size	Orientation
46503	UNIVERSAL	LEFT
46505	UNIVERSAL	RIGHT

**CALF SKIN**

Contoured fit, controls muscle movement and reduces muscle vibration.

Code	Size
B46300	XSMALL
B46301	SMALL
B46302	MEDIUM
B46303	LARGE
B46304	XLARGE
B46305	XXLARGE



QUATTROMED

Stabilises knee joint to relieve pain. Corrects knee joint movement and prevents lateral movement of the patella. Lightweight, low profile design improves comfort and promotes compliance. Multiple flex points mimics the knee's natural movement and range of motion for maximum mobility. Comfortable thigh, knee and calf contours provide a custom fit for improved patient comfort and compliance.



Code	Size	Orientation
0416-2201	XSMALL	LEFT
0416-2202	SMALL	LEFT
0416-2203	MEDIUM	LEFT
0416-2204	LARGE	LEFT
0416-2205	XLARGE	LEFT
0416-2206	XSMALL	RIGHT
0416-2207	SMALL	RIGHT
0416-2208	MEDIUM	RIGHT
0416-2209	LARGE	RIGHT
0416-2210	XLARGE	RIGHT

QUATTROMED SIZING

Size	15cm above the Patella	Mid Patella	15cm below the Patella
XS	390 - 460 mm	330 - 390 mm	310 - 370 mm
S	460 - 530 mm	390 - 440 mm	370 - 420 mm
M	530 - 600 mm	440 - 500 mm	420 - 470 mm
L	600 - 670 mm	500 - 560 mm	470 - 530 mm
XL	670 - 750 mm	560 - 630 mm	530 - 590 mm

MACKIE KNEE

The Mackie Hinge Knee Brace uses a unique infinitely adjustable hinge mounted on a rigid, non-contact frame. This design offers simplicity of function, strong construction and superior patient comfort. Because human physiology varies greatly from individual to individual, the hinge adjustability allows the patient to control how much force is applied to achieve the best result. For this patient-directed approach to be successful, patient tolerance and compliance is crucial. Static-Progressive Stretch is a proven concept for effective and permanent gain in range of motion.



Code	Size	Knee Circ.
S4000	SMALL	298 - 369mm
M4001	MEDIUM	381 - 432mm
L4002	LARGE	432 - 508mm

OAPL ROM KNEE

A ROM knee brace featuring extendable arms, reducing inventory as there is no need for stocking different devices of varying length. Aluminium arms allows for a more intimate fit as they can be contoured to fit any patient's leg shape. Pinless hinge allows easy adjustment with no parts that can be lost.



Code	Variation	Length
E-KN098	STANDARD	500 - 620mm
E-KN098C	COOL FIT	500 - 620mm

UNIVERSAL POST OP KNEE

To maintain the knee in full extension whilst providing mediolateral support. A soft foam support with detachable panels and touch-tape closures. Removable medial, lateral and posterior aluminium stays can be contoured to the shape of the leg. Egg shell foam removable knee support/pad. Three Panels.



Code	Size	Length
14000E	UNIVERSAL	400mm
14001E	UNIVERSAL	500mm
14002E	UNIVERSAL	600mm
B14001E	BARIATRIC	500mm
B14002E	BARIATRIC	600mm

POST OP KNEE IMMOBILISER

Adjustable condyle pads; Increased patient comfort and control by providing snug padded M/L containment. Touch tape closure for stay pockets; Easy access to stays for adjustment during fitting. Secure pockets for stays preventing migration during use. "T" Section Aluminium stays; Increased strength providing greater immobilization. One Panel.



Code	Size	Length	Circumference mm
14011E	SMALL	480mm	Calf 300 Thigh 400
14012E	MEDIUM	480mm	Calf 400 Thigh 500
14013E	LARGE	480mm	Calf 500 Thigh 640
14014E	XLARGE	480mm	Calf 600 Thigh 770
14011AE	SMALL	580mm	Calf 300 Thigh 400
14012AE	MEDIUM	580mm	Calf 400 Thigh 500
14013AE	LARGE	580mm	Calf 500 Thigh 640
14014AE	XLARGE	580mm	Calf 600 Thigh 770

ROM POST OP 20' KNEE

Pinless ROM Hinge for ease of adjustment, flexion and extension stops at 10 degrees increments. Anti-tamper screw to prevent patients from accessing ROM adjustments. Composite uprights to allow contouring for the patient's limb and trimmable straps make this easy to customise.

Code	Size
ROM013BAS	UNIVERSAL



ROM 23 KNEE"

Pinless ROM Hinge for ease of adjustment, flexion and extension stops at 10 degrees increments. Anti-tamper screw to prevent patients from accessing ROM adjustments. Metal uprights to allow contouring to patient's limbs and snap of adjustments on both proximal and distal uprights. Adjustable strap carriers can be then place over uprights leaving no exposed rough edges. Automatic full extension locking feature and trimmable straps make this easy to customise.

Code	Size
ROM13000SB	UNIVERSAL



ULTIMA ROM & DL 23' KNEE

Pinless ROM Hinge for ease of adjustment, flexion and extension stops 10 degree increments. It features anti-tamper screw to prevent patient from accessing ROM adjustments. Metal uprights allow contouring to patient's limb and snap of adjustments on both proximal and distal uprights. Adjustable strap carriers can then be placed over uprights leaving no exposed rough edges. Automatic full extension locking feature and trimmable straps.



Code	Size
ROMC13000SB	UNIVERSAL

VANTAGE AP CONTROL ROM KNEE

Pinless ROM Hinge for ease of adjustment, flexion and extension stops at 10 degrees increments. Anti-tamper screw to prevent patients from accessing ROM adjustments. Metal uprights to allow contouring to patient's limbs and snap of adjustments on both proximal and distal uprights. Adjustable strap carriers can be then placed over uprights leaving no exposed rough edges and adjustable strap carriers allowing precise positioning of straps.



Code	Size
ROM130050SB	UNIVERSAL

CAMP COMBO

Combo is a lightweight, removable knee brace attachment for ToeOFF® (or BlueRocker™) that creates an innovative orthotic solution for management of footdrop combined with genu-recurvatum. Uprights easily shaped and adjusted to optimum height for patient, modular components allow you to accommodate variances in thigh and calf circumferences, low profile, lightweight, easy to don and off, easily removed when wearer is relaxing. Combo is fitted individually for each patient. The different components can be combined to accommodate patient's size and length. Size of thigh cuff is determined by length of thigh segment as well as circumference of the thigh. Size of Tibia connector is determined by the size of the ToeOFF or BlueRocker. Size of the Hyperextension adjustment is determined from knee width (M-L). All 5 product codes need to be ordered for a complete orthosis. Combo cannot be used for patients over 150 kg or 2 meters tall.



Code	Component	Size	Length	Orientation	Thigh Circ.
287800000	JOINTS & UPRIGHTS				
287811011	THIGH CUFF	SMALL	110 - 150mm	LEFT	300 - 400mm
287811012	THIGH CUFF	MEDIUM	130 - 150mm	LEFT	380 - 500mm
287811013	THIGH CUFF	LARGE	150 - 190mm	LEFT	480 - 600mm
287811014	THIGH CUFF	XLARGE	170 - 210mm	LEFT	580 - 700mm
287812011	THIGH CUFF	SMALL	110 - 150mm	RIGHT	300 - 400mm
287812012	THIGH CUFF	MEDIUM	130 - 150mm	RIGHT	380 - 500mm
287812013	THIGH CUFF	LARGE	150 - 190mm	RIGHT	480 - 600mm
287812014	THIGH CUFF	XLARGE	170 - 210mm	RIGHT	580 - 700mm
287820011	HYPEREXTENSION ADJUSTABLE	SMALL	70 - 105mm		
287820012	HYPEREXTENSION ADJUSTABLE	MEDIUM	95 - 130mm		
287820013	HYPEREXTENSION ADJUSTABLE	LARGE	125 - 160mm		
287830011	SOFTKIT	SMALL			
287830012	SOFTKIT	MEDIUM			
287830013	SOFTKIT	LARGE			
287830014	SOFTKIT	XLARGE			
287840011	TIBIA CONNECTOR	SMALL			
287840012	TIBIA CONNECTOR	MEDIUM			
287840013	TIBIA CONNECTOR	LARGE			

PREMIUM THIGH SUPPORT

To provide support and maintain warmth, assisting in repair of thigh injuries. Neoprene support covering the thigh.

Code	Size	Thigh Circ.
65050	SMALL	490 - 530mm
65051	MEDIUM	540 - 580mm
65052	LARGE	590 - 630mm
65053	XLARGE	640 - 680mm



NEOPRENE SHORTS

To provide support and maintain warmth, assisting in the repair of muscle and tendon injuries to the thigh and groin.

Code	Size	Thigh Circ.
70075	SMALL	760 - 860mm
70076	MEDIUM	860 - 960mm
70077	LARGE	960 - 1065mm
70078	XLARGE	1065 - 1170mm



BIOSKIN PREMIUM BRACING

The result of combining innovative designs, based on a firm understanding of biomechanics, with materials that have been designed specifically for orthopaedic and sports medicine applications. The patented BioSkin material provides a high level of compression without limiting activity or creating any discomfort. BioSkin is thin and gives radial compression so it conforms to the curves of the body. BioSkin is also breathable, so perspiration evaporates right through it and is hypo-allergenic eliminating the chance of skin irritation.

THIGH SKIN

To provide compression to quads and hamstrings. Pull-on thigh support available with our without leg cinch strap.

Code	Size	Thigh Circ.
B46100	XSMALL	410 - 460mm
B46101	SMALL	460 - 510mm
B46102	MEDIUM	510 - 560mm
B46103	LARGE	560 - 640mm
B46104	XLARGE	640 - 710mm
B46105	XXLARGE	710 - 810mm



THIGH SKIN W/CINCH STRAP

To provide compression to quads and hamstrings. Pull-on thigh support available with our without leg cinch strap.

Code	Size	Thigh Circ.
B46110	XSMALL	410 - 460mm
B46111	SMALL	460 - 510mm
B46112	MEDIUM	510 - 560mm
B46113	LARGE	560 - 640mm
B46114	XLARGE	640 - 710mm
B46115	XXLARGE	710 - 810mm



COMPRESSION SHORTS

The high-level compression provided by Bio Skin Compression Shorts has been shown to improve athletic performance and endurance by decreasing muscle vibration. The unique Bio Skin material allows moisture to evaporate, providing a cool, comfortable experience.

Code	Size	Waist Circ.
B46210	XSMALL	660 - 710mm
B46201	SMALL	710 - 810mm
B46202	MEDIUM	810 - 920mm
B46203	LARGE	920 - 102mm
B46204	XLARGE	1020 - 1120mm
B46205	XXLARGE	1120 - 1220mm



COMP SHORTS W CINCH STRAP

B46231	SMALL	710 - 810mm
B46232	MEDIUM	810 - 920mm
B46233	LARGE	920 - 102mm
B46234	XLARGE	1020 - 1120mm

COMP SHORTS W GROIN STRAP

B46211	SMALL	710 - 810mm
B46212	MEDIUM	810 - 920mm
B46213	LARGE	920 - 102mm
B46214	XLARGE	102 - 112mm

ATHLETIC SUPPORT

Offers sportsmen support and comfort throughout a full range of hip and trunk movement. A broad elastic wristband and cotton elastic leg straps connect to a soft open weave pouch.

Code	Size	Waist Circ.
44044	SMALL	700 - 800mm
44045	MEDIUM	810 - 910mm
44046	LARGE	920 - 102mm
44047	XLARGE	1020mm & Up



LOWER HERNIA SUPPORT

The Lower Hernia Support provides gentle and gradual pressure to reduce right or left sided inguinal hernias. It helps prevent recurrent bulge on exertion and strain involved in lifting and forward bending for example. The discreet support is fitted with adjustable and removable pads for user comfort and effectiveness. The support has reduced stitching to decrease likelihood of friction and irritation for comfort and diabetic use.



Code	Size	Orientation	Wasit Circ
601-1-L	SMALL	LEFT	650 - 800mm
601-2-L	MEDIUM	LEFT	800 - 1000mm
601-3-L	LARGE	LEFT	1000 - 1250mm
601-1-R	SMALL	RIGHT	650 - 800mm
601-2-R	MEDIUM	RIGHT	800 - 1000mm
601-3-R	LARGE	RIGHT	1000 - 1250mm

DOUBLE HERNIA SUPPORT

The Double Hernia Support provides gentle and gradual pressure to reduce bilateral lower/inguinal hernias. It helps prevent recurrent bulge on exertion and strain involved in lifting and forward bending for example. The discreet support is fitted with adjustable and removable pads for user comfort and effectiveness. The support has reduced stitching to decrease likelihood of friction and irritation for comfort and diabetic use.



Code	Size	Waist Circ
602S	SMALL	650 - 800mm
602M	MEDIUM	800 - 1000mm
602L	LARGE	1000 - 1250mm

NEWPORT 3 HIP SYSTEM

The Newport 3 Hip Orthosis has been the industry standard for over a decade. Innovations like the patented clamping disks, allowing precise contouring of the pelvic section and the Universal Virtual Joint, which provides unprecedented strength in a small component, have made the Newport Hip Orthosis the most popular choice for physicians, orthotists and patients worldwide.



Code	Size	Orientation
3640	SMALL	LEFT
3641	MEDIUM	LEFT
3642	LARGE	LEFT
3643	XLARGE	LEFT
3645	SMALL	RIGHT
3646	MEDIUM	RIGHT
3647	LARGE	RIGHT
3648	XLARGE	RIGHT

NEWPORT 3 HIP SYSTEM UNIVERSAL

3640.40	SMALL
3641.41	MEDIUM
3642.42	LARGE
3643.43	XLARGE

NEWPORT 4 HIP SYSTEM

The totally modular and patented Newport 4 Hip Orthosis has been designed and engineered to provide the most comfortable hip orthosis for your patients with innovative new options to increase function and control.



Code	Size	Orientation
3742	SMALL	LEFT
3743	MEDIUM	LEFT
3744	LARGE	LEFT
3745	XLARGE	LEFT
3746	SMALL	RIGHT
3747	MEDIUM	RIGHT
3748	LARGE	RIGHT
3749	XLARGE	RIGHT

CALIFORNIA HIP ORTHOSIS

The California Hip Orthosis is the latest innovation in Orthomerica's leading hip orthosis designs which helps to prevent dislocation. The California Hip Orthosis was created to facilitate the best pelvic and thigh control and minimise the bulk of conventional hip orthoses for patients at risk to dislocate.



Code	Size	Orientation
3960	SMALL	LEFT
3961	MEDIUM	LEFT
3962	LARGE	LEFT
3963	XLARGE	LEFT
3965	SMALL	RIGHT
3966	MEDIUM	RIGHT
3967	LARGE	RIGHT
3968	XLARGE	RIGHT

CALIFORNIA SIZING

Size	Waist Circ	Hip Circ	Mid Thigh
S	560 - 610 mm	660 - 810 mm	330 - 480 mm
M	710 - 860 mm	810 - 970 mm	380 - 530 mm
L	860 - 1020 mm	970 - 1120 mm	430 - 610 mm
XL	1020 - 1170 mm	1120 - 1270 mm	480 - 660 mm

NEWPORT SIZING

Size	Waist Circ	Hip Circ	Mid Thigh
S	560 - 710 mm	660 - 810 mm	330 - 480 mm
M	710 - 860 mm	810 - 960 mm	380 - 530 mm
L	860 - 1010 mm	960 - 1110 mm	430 - 610 mm
XL	1010 - 1160 mm	1110 - 1270 mm	480 - 660 mm

HIP FLEXION ASSIST DEVICE

Lightweight padded waist strap. Dual adjustable flexion assist cords. Improved patient comfort and compliance. Low profile design. Improved compliance as device can be worn under clothing. Universal shoe attachment. Increased flexibility for patient footwear selection.



Code	Size	Orientation
200HFAD-RS	SMALL	RIGHT
200HFAD-RM	MEDIUM	RIGHT
200HFAD-RL	LARGE	RIGHT
200HFAD-LS	SMALL	LEFT
200HFAD-LM	MEDIUM	LEFT
200HFAD-LL	LARGE	LEFT
200HFAD-BS	SMALL	BI-LATERAL
200HFAD-BM	MEDIUM	BI-LATERAL
200HFAD-BL	LARGE	BI-LATERAL

HIP FLEXION SIZING

Size	Waist Circ
Small	610 - 820 mm
Medium	830 - 1020 mm
Large	1030 - 1220 mm

CHARNLEY PILLOW

To maintain the hips in abduction whilst resting in bed. A foam cushion covered in vinyl with straps to hold legs in prison.

Code	Size
16015	UNIVERSAL





Spinal

Sacro Iliac 32

Lumbo Sacral 33

Thoraco Lumbo
Sacral 37

Cervical 40

Headgear 41

Other 42

ELASTIC MATERNITY BELT

All Elastic belt with adjustable touch tape closure to provides comfort and support by relieving the back and abdominal region from excessive pressure, during pregnancy.

Code	Hip Circ.
43031	810mm
43032	880mm
43033	990mm
43034	108mm
43035	1160mm
43036	1270mm



TROCHANTER BELT WEBBING WITH FOAM HIP PADS

Webbing strap belt featuring soft foam hip pads to cushion and provide extra comfort and support by relieving the back and abdominal region from excessive pressure during pregnancy. D ring and touch tape closure.

Code	Hip Circ.
43026	810mm
43027	880mm
43028	990mm
43029	108mm
43030	1160mm



TROCHANTER BELT CANVAS/ELASTIC

Canvas belt with an elastic back section, with D ring and touch tape closure. Provides comfort and support by relieving the back and abdominal region from excessive pressure during pregnancy.

Code	Hip Circ.
43021	810mm
43022	880mm
43023	990mm
43024	108mm
43025	1160mm



CYBERTECH 630 SI BELT

Significantly reduces motion about the sacroiliac joint. Instant pain relief in most cases. 5:1 mechanical advantage compression power provides unmatched pelvic sacral support. Flexible sacroiliac orthoses with breathable mesh and fabric closure for one size fits all. Features flexible pulley system which adjusts the brace to the body and adjustable pads which can be customised for the patient.

Code	Size	Circumference
G200	Universal	584 - 1280mm



LOW PROFILE SI BELT

The new low profile SI Belt provides stabilisation of the sacroiliac joints through direct circumferential compression utilising the 3 to 1 ratio mechanical advantage pulley system. The SI belt is prescribed for conservative treatment of sacroiliac strain or arthritis. Significantly reduces motion about the sacroiliac joint providing instant pain relief in most cases. Breathable mesh and adjustable pads create a very comfortable orthosis.

Code	Size	Circumference
CYB630-LB-2	SMALL	Up to 889mm
CYB630-LB-3	STANDARD	Over to 889mm



TRO CAMP BELT

Designed by camp, provides comfort and support by relieving the back and abdominal region from excessive pressure.

Code
210210-105



RIB/ABDOMINAL BINDER

The OAPL Rib / Abdominal Binder can be used to provide support to the ribs and / or abdominal area. Support of this area can reduce back pain and act as a reminder for correct posture and partially restrict motion. The OAPL Rib / Abdominal Binder features an all elastic belt with touch tape closure for correct fitting. Available in three widths to support a variety of applications.

Code	Size	Height	Circumference
OA-WA0506S	SMALL	150mm	800-1100mm
OA-WA0506M	MEDIUM	150mm	980-1350mm
OA-WA0506L	LARGE	150mm	1180-1580mm
OA-WA0506XL	XLARGE	150mm	1400-2200mm
OA-WA0509S	SMALL	230mm	800-1100mm
OA-WA0509M	MEDIUM	230mm	980-1350mm
OA-WA0509L	LARGE	230mm	1180-1580mm
OA-WA0509XL	XLARGE	230mm	1400-2200mm
OA-WA05012S	SMALL	300mm	800-1100mm
OA-WA05012M	MEDIUM	300mm	980-1350mm
OA-WA05012L	LARGE	300mm	1180-1580mm
OA-WA05012XL	XLARGE	300mm	1400-2200mm



PREMIUM BACK

The OAPL Premium Back Support provides compression and warmth to promote rapid healing of soft tissue injuries. Maintaining body heat helps to increase blood flow, promoting elasticity of muscles, tendons and ligaments. This ventilated neoprene support with moisture wicking technology and soft terry lining. Provides comfortable support which aids in reducing swelling and associated pain.

Code	Size	Waist Circ.
65090	SMALL	700 - 800mm
65091	MEDIUM	810 - 900mm
65092	LARGE	910 - 1000mm
65093	XLARGE	1000 - 1200mm



COOL FIT CINCH

To provide support to the lumbar spine. A well ventilated mesh support with elastic segment and touch tape adjustment.

Code	Size	Waist Circ.
21009	XSMALL	750 - 850mm
21010	SMALL	850 - 950mm
21011	MEDIUM	950 - 1050mm
21012	LARGE	1050 - 1150mm
21013	XLARGE	1150 - 1250mm
21014	XXLARGE	1250 - 1350mm



SACRO CINCH

Reduces intradiscal pressure by increasing intra-abdominal pressure. Acts as a reminder to partially restrict motion. An all elastic support with velcro closures and a wide encircling cinch pull strap. Metal stays, if indicated, fit into two pouches posteriorly which require individual shaping at time of fitting.

Code	Size	Waist Circ.
21005A	SMALL	710 - 810mm
21006A	MEDIUM	810 - 910mm
21007A	LARGE	910 - 1015mm
21008A	XLARGE	1015 - 1115mm
21009A	XXLARGE	1250 - 1350mm



INDUSTRIAL LIFTING BELT

To provide comfortable support to reduce back strain during lifting and work activity. A black elastic belt with flexible stays, adjustable shoulder straps, elastic compression pulls and a velcro closure anteriorly.

Code	Size	Waist Circ.
1060S	SMALL	860 - 960mm
1060M	MEDIUM	960 - 106mm
1060L	LARGE	106 - 116mm
1060XL	XLARGE	116-126mm



BIOSKIN BACK SKIN

To provide pain relief for chronic and acute back pain. Low profile brace that offers a cool and comfortable fit. The unique lumbar supports offer pain relieving support to the low spine. The vented elastic support offer pain relieving support to the low spine. The vented elastic side panels provide heat release and the removable double cinch belt gives flexibility and support. This brace can be customised with buttress options.



Code	Size	Waist Circ.
B47100	XSMALL	610 - 660mm
B47101	SMALL	660 - 760mm
B47102	MEDIUM	760 - 860mm
B47103	LARGE	860 - 970mm
B47104	XLARGE	970 - 1070mm
B47105	XXLARGE	1070 - 1170mm

WITH LUMBAR PAD

B47200	XSMALL	610 - 660mm
B47201	SMALL	660 - 760mm
B47202	MEDIUM	760 - 860mm
B47203	LARGE	860 - 970mm
B47204	XLARGE	970 - 1070mm
B47205	XXLARGE	1070 - 1170mm

WITH LUMBAR PAD & FLEXIBLE SUPPORT

B47240	XSMALL	610 - 660mm
B47241	SMALL	660 - 760mm
B47242	MEDIUM	760 - 860mm
B47243	LARGE	860 - 970mm
B47244	XLARGE	970 - 1070mm
B47245	XXLARGE	1070 - 1170mm

WITH POCKET FOR THERMAL INSERT

B47300	XSMALL	610 - 660mm
B47301	SMALL	660 - 760mm
B47302	MEDIUM	760 - 860mm
B47303	LARGE	860 - 970mm
B47304	XLARGE	970 - 1070mm
B47305	XXLARGE	1070 - 1170mm

WITH POCKET & THERMAL INSERT

B47390	XSMALL	610 - 660mm
B47391	SMALL	660 - 760mm
B47392	MEDIUM	760 - 860mm
B47393	LARGE	860 - 970mm
B47394	XLARGE	970 - 1070mm
B47395	XXLARGE	1070 - 1170mm

CALIFORNIA COMPRESSOR

The Compressor is a low-profile lumbar orthosis that provides maximum abdominal support with one easy pull. Its one-handed compression pull tab makes application simple and quick. Breathable comfortable mesh metal material minimises heat build-up common with other spinal orthoses. Indications include acute and chronic low back pain, spondylolisthesis and post-op low lumbar laminectomy. Its lightweight construction is optimal for patients requiring support for everyday activities such as lifting, long periods of standing, work related motions, golf and gardening.



Code	Size	Height	Waist Circ.
3219	XSMALL	200mm	630-760mm
3220	SMALL	200mm	760-890mm
3221	MEDIUM	200mm	890-1010mm
3222	LARGE	200mm	1010-1140mm
3223	XLARGE	200mm	1140-1270mm
3224	XXLARGE	200mm	1270-1390mm
3225	3XLARGE	200mm	1390-1520mm
3226	4XLARGE	200mm	1520-1640mm
3229	XSMALL	250mm	630-760mm
3230	SMALL	250mm	760-890mm
3231	MEDIUM	250mm	890-1010mm
3232	LARGE	250mm	1010-1140mm
3233	XLARGE	250mm	1140-1270mm
3234	XXLARGE	250mm	1270-1390mm
3235	3XLARGE	250mm	1390-1520mm
3236	4XLARGE	250mm	1520-1640mm

CT SPINE BRACE

First step intervention for the treatment of lower back pain, designed to give control immediately to the patient. Mechanical advantage pulley system is designed to enable adjustment throughout the day to provide real support when you needed and less support when not needed.



Code	Size	Waist Circ.
G100S	SMALL	630 - 760mm
G100M	MEDIUM	760 - 890mm
G100L	LARGE	890 - 1020mm
G100XL	XLARGE	1020 - 1140mm

CAMP CORSETS

Reduces intradiscal pressure by increasing intra-abdominal pressure. Acts as a reminder to partially restrict motion. Increased comfort and compliance due to the material's ability to conform and it's breathability. Cinch is made from eyeletted, non-stretch, nylon cinch straps allows the patient to easily adjust the degree of support applied to the spine.

CAMP IMPROVE BACK SACRO

Front Height 110mm, Back Height 185mm

Code	Size
23331S	SMALL
23331M	MEDIUM
23331L	LARGE
23331XL	XLARGE



CAMP IMPROVE BACK LUMBO

Front Height 185mm, Back Height 255mm.

Code	Size
23332S	SMALL
23332M	MEDIUM
23332L	LARGE
23332XL	XLARGE



CAMP IMPROVE LONG LUMBO

Front Height 230mm, Back Height 330mm.

Code	Size
23333S	SMALL
23333M	MEDIUM
23333L	LARGE
23333XL	XLARGE



CAMP XXI LUMBOSACRAL CORSET

Front Height 250mm, Back Height 320mm.

Code	Size
29433S	SMALL
29433M	MEDIUM
29433L	LARGE
29433XL	XLARGE



CAMP XXI THORACO LUMBAR CORSET

Front Height 240mm, Back Height 520mm.

Code	Size
29437S	SMALL
29437M	MEDIUM
29437L	LARGE
29437XL	XLARGE



CAMP CORSETS SIZING

Size	Waist Circ
S	800 - 950mm
M	900 - 1050mm
L	1000 - 1150mm
XL	1100 - 1250mm

OAPL LUMBAR SACRAL BRACE

To limit flexion and extension of the lumbosacral spine. The abdominal pad reduces intradiscal pressure by increasing intra-abdominal pressure. A lightweight aluminium posterior frame covered with leather and lined with felt. Attaches via webbing straps to the abdominal pad. Please note, strap kit and abdominal bib must be used with this brace and purchased separately. Length is 1" below scapula to 3" below PSIS>



Code	Variation	Length	Apx. Patient Hgt
22009		280mm	1630mm
22010		300mm	1630 - 1730mm
22011		330mm	1730 - 1830mm
22012		350mm	1830mm & Up
22009S	SHOWER	280mm	1630mm
22010S	SHOWER	300mm	1630 - 1730mm
22011S	SHOWER	330mm	1730 - 1830mm
22012S	SHOWER	350mm	1830mm & Up

OAPL LUMBAR SACRAL STRAP KIT

Is used in conjunction with the OAPL Lumbar Sacral Brace and OAPL High Spinal Brace to fit the brace to the patient.

Code	Size
84019AAA	XXSMALL
84019AA	XSMALL
84019A	SMALL
84019B	MEDIUM
84019C	LARGE
84019D	XLARGE
84019DD	XXLARGE



OAPL LUMBAR SACRAL ABDOMINAL BIB

Is used in conjunction with the OAPL Lumbar Sacral Brace and OAPL High Spinal Brace to fit the brace to the patient.

Code	Size	Variation
23001A	SMALL	
22015A	LARGE	
23001AS	SMALL	SHOWER
22015AS	LARGE	SHOWER



KYDEX CHAIRBACK BRACE

To limit flexion, extension and lateral movement of the lumbosacral spine. The corset front reduces intradiscal pressure by increasing intra-abdominal pressure. A polyfoam lined kydex plastic frame posteriorly, which attaches to a corset front anteriorly. Hygienic and lightweight. Use heatgun if reshaping is required.

Code	Size	Length	Hip Circ.
22001	SMALL	290mm	740 - 820mm
22002	MEDIUM	320mm	830 - 920mm
22003	LARGE	330mm	930 - 1020mm
22004	XLARGE	330mm	1030 - 1120mm



CORSET FRONTS FOR KYDEX BRACES

Open weave, mesh, full front corsets to be used in conjunction with lumbosacral and thoracolumbosacral Kydex braces. Corresponds with Kydex brace size. Xlarge corset can be used with the large Kydex braces when hip circumference exceeds 112cm.



Code	Size
22018	SMALL
22019	MEDIUM
22021	LARGE
22022	XLARGE

CALIFORNIA MID PROFILE

New laparoscopic surgical techniques and innovative rehabilitation programs are changing the orthopaedic management of spine disorders. Since these procedures have become less invasive and create less soft tissue disruption, there is a need for new orthotic management to complement the spinal treatment program. The Mid-Profile Orthosis provides stability and control only where it is indicated in order to reduce spinal muscle atrophy and deconditioning. This device can be upgraded to a TLSO using the California Ant Thoracic Extension.



Code	Size	Height	Width	Waist Circ.
3239	XSMALL	300mm	460mm	640-760mm
3240	SMALL	300mm	360mm	760-890mm
3241	MEDIUM	300mm	260mm	890-1020mm
3242	LARGE	340mm	560mm	1020-1140mm
3243	XLARGE	340mm	610mm	1140-1270mm
3244	XXLARGE	340mm	660mm	1270-1400mm
3245	3XLARGE	340mm	700mm	1400-1520mm

CALIFORNIA SOFT SPINAL

The pre-fabricated "Original" California Soft Spinal Orthosis has been designed to alleviate pain, limit undesirable motion and provide circumferential pressure to stabilise the spine. The basic system is modular, comprising rigid components inside a soft and comfortable interface. Indicated for spinal problems between L3 and S1. The orthosis can be effective for more superior spinal problems between T7 and L3, by adding either an anterior, posterior, or combined upper thoracic component. All components can be retrofitted and are modular and interchangeable. Velcro attachable lateral stabilising supports are also available when more control of lateral flexion and rotation is needed. The anterior and posterior ABS modular components can be heated and flared to provide a more intimate fit. This device can be upgraded to a TLSO using the California Ant Thoracic Extension.



Code	Size	Height	Waist Circ.
3360	XSMALL	270mm	710-810mm
3361	SMALL	280mm	760-860mm
3362	MEDIUM	280mm	830-940mm
3363	LARGE	280mm	910-106mm
3364	XLARGE	280mm	1040-1190mm
3365	XXLARGE	290mm	1160-1320mm
3366	3XLARGE	300mm	1290-1440mm

ANTERIOR OVERLAP BRACE

Pre and Post-Operative stabilisation of the spine. Prevents flexion, extension or lateral deviation. 1/8" Polyethylene cylinder, anterior opening, 15 Lordosis standard, closure kit included.



Code	Size	Hip Circ	Waist Circ	Chest Circ
210401	1	710mm	510mm	660mm
21041.5	1.5	735mm	560mm	685mm
210402	2	760mm	610mm	710mm
21042.5	2.5	785mm	660mm	735mm
210403	3	810mm	710mm	760mm
21043.5	3.5	837mm	735mm	800mm
210404	4	865mm	760mm	840mm
21044.5	4.5	890mm	785mm	865mm
210405	5	915mm	810mm	890mm
21045.5	5.5	940mm	837mm	915mm
210406	6	965mm	865mm	940mm
21046.5	6.5	990mm	890mm	965mm
210407	7	1015mm	915mm	990mm
21047.5	7.5	1052mm	940mm	1040mm
210408	8	1090mm	965mm	1090mm
21048.5	8.5	1122mm	990mm	1140mm

CUT 4 CUSTOM LSO

3D-LITE Cut 4 Custom, when you need a custom LSO or TLSO NOW! Comfortable, lightweight, breathable. Ability to contour to severe scoliosis or kyphotic spinal deformities. May be worn in the shower. Non-Toxic. No reports of skin irritation. Same day fabrication. Front and back Cut 4 Custom sections have straps pre-attached. Both sections are ready to be heated and custom moulded to patient. Backs, fronts and liners sold separately and interchangeable.



Code	Size	Variat.	Height	Waist Circ.
649190011	XS/S	FRONT	195mm	600-750mm
649190012	S/M	FRONT	195mm	700-900mm
649190013	M/L	FRONT	215mm	850-1050mm
649190014	L/XL	FRONT	245mm	950-1200mm
649190015	XL/XXL	FRONT	270mm	1150-1400mm
649210011	XS/S	BACK	300mm	600-750mm
649210012	S/M	BACK	315mm	700-900mm
649210013	M/L	BACK	335mm	850-1050mm
649210014	L/XL	BACK	350mm	950-1200mm
649210015	XL/XXL	BACK	365mm	1150-1400mm

LINNING

649290011	XS/S	LINING FRONT
649290012	S/M	LINING FRONT
649290013	M/L	LINING FRONT
649290014	L/XL	LINING FRONT
649290015	XL/XXL	LINING FRONT
649310011	XS/S	LINING BACK
649310012	S/M	LINING BACK
649310013	M/L	LINING BACK
649310014	L/XL	LINING BACK
649310015	XL/XXL	LINING BACK

OAPL HIGH SPINAL BRACE

To limit flexion and extension of the thoracolumbosacral spine. The abdominal pad reduces intradiscal pressure by increasing intra-abdominal pressure. Shoulder straps restrict upper thoracic flexion. A lightweight aluminium posterior framed covered with leather and lined with felt. Attaches via webbing straps to the abdominal pad. Please note, strap kit and abdominal bid must be used with this brace and purchased separately. Length measurement is spine of scapula to 3" below PSIS.



Code	Variation	Post Height
23001		460mm
23002		485mm
23003		510mm
23004		535mm
23005		560mm
23001S	SHOWER	460mm
23002S	SHOWER	485mm
23003S	SHOWER	510mm
23004S	SHOWER	535mm
23005S	SHOWER	560mm

OAPL LUMBAR SACRAL STRAP KIT

Is used in conjunction with the OAPL Lumbar Sacral Brace and OAPL High Spinal Brace to fit the brace to the patient.

Code	Size
84019AAA	XXSMALL
84019AA	XSMALL
84019A	SMALL
84019B	MEDIUM
84019C	LARGE
84019D	XLARGE
84019DD	XXLARGE



OAPL LUMBAR SACRAL ABDOMINAL BIB

Is used in conjunction with the OAPL Lumbar Sacral Brace and OAPL High Spinal Brace to fit the brace to the patient.

Code	Size	Variation
23001A	SMALL	
22015A	LARGE	
23001AS	SMALL	SHOWER
22015AS	LARGE	SHOWER



KYDEX KNIGHT TAYLOR BRACE

To limit flexion, extension and lateral movement of the thoracolumbosacral spine. Corset front reduces intradiscal pressure by increasing intra-abdominal pressure. A polyfoam lined Kydex plastic frame posteriorly, which attaches to a corset front anteriorly. Hygienic and lightweight. Use heat gun if reshaping is required.



Code	Size	Length	Hip Circ.
23014	SMALL	470mm	740 - 820mm
23015	MEDIUM	490mm	830 - 920mm
23016	LARGE	540mm	930 - 1020mm
23017	XLARGE	540mm	1030 - 1120mm

CORSET FRONTS FOR KYDEX BRACES

Open weave, mesh, full front corsets to be used in conjunction with lumbosacral and thoracolumbosacral Kydex braces.

Corresponds with Kydex brace size. Xlarge corset can be used with the large Kydex braces when hip circumference exceeds 112cm.



Code	Size
22018	SMALL
22019	MEDIUM
22021	LARGE
22022	XLARGE

CALIFORNIA ECO

Extension Compression Orthosis was developed to comfortably extend the thoracic spine and unweight the wedged vertebral body in order to promote healing and pain relief. The ECO stabilises the pelvis and lumbar spine to ease pain. Effective for individuals with Osteoporosis, vertebral compression fractures, thoracic strains and debilitating Kyphosis. The ECO is as simple to apply as a backpack or jacket, with easy-to-reach self-adjusting and padded shoulder straps. The metal spinal frame can be easily contoured to the patient's thoracic and lumbar spine.



Code	Size	Variation	Waist Circ.
3253	XSMALL	STANDARD	640-760mm
3254	SMALL	STANDARD	760-890mm
3255	MEDIUM	STANDARD	890-1020mm
3256	LARGE	STANDARD	1020-1140mm
3257	XLARGE	STANDARD	1140-1270mm
3273	XSMALL	TALL	640-760mm
3274	SMALL	TALL	760-890mm
3275	MEDIUM	TALL	890-1020mm
3276	LARGE	TALL	1020-1140mm
3277	XLARGE	TALL	1140-1270mm

CALIFORNIA ANT THORACIC EXTENSION

Orthomerica's new modular ATE interfaces seamlessly with the industry leading California Soft Spinal System, Cosmic, Torso, Airback and other pre-fabricated spinal modules. It also easily attaches to custom bi-valve LSOs or TLSOs.

Code	Variation	Length
2013	STANDARD	160-290mm
2014	LONG	220-350mm



CUT 4 CUSTOM TLSO

3D-LITE Cut 4 Custom, when you need a custom LSO or TLSO NOW! Comfortable, lightweight, breathable. Ability to contour to severe scoliosis or kyphotic spinal deformities. May be worn in the shower. Non-Toxic. No reports of skin irritation. Same day fabrication. Front and back Cut 4 Custom sections have straps pre-attached. Both sections are ready to be heated and custom moulded to patient. Backs, fronts and liners sold separately and interchangeable.



Code	Size	Variat.	Height	Wasit
649220010	S/M	FRONT	420mm	700-900mm
649220012	M\L	FRONT	480mm	850-1100mm
649220014	XL/XXL	FRONT	540mm	1050-1350mm
649230010	S/M	BACK	420mm	700-900mm
649230012	M\L	BACK	480mm	850-1100mm
649230014	XL/XXL	BACK	540mm	1050-1350mm

LINING

649320010	S/M	FRONT LINER
649320012	M\L	FRONT LINER
649320014	XL/XXL	FRONT LINER
649330010	S/M	BACK LINER
649330012	M\L	BACK LINER
649330014	XL/XXL	BACK LINER

CYBERSPINE TSLO

The Cyberspine Thoracal Lumbar Sacral Orthosis is designed to comfortably control and manage acute pain often associated with osteoporosis, compression fractures, spinal stenosis, strain, and excessive kyphosis. The CyberSpine TLSO provides superior triplanar motion control, extending from sacrococcygeal junction and terminating just inferior of the scapular spine. The four segmental rigid plastic shell components are ergonomically designed to restrict gross trunk motion in the sagittal, frontal and transverse planes. The anterior rigid plastic and thoracic extension extends from the symphysis pubis to the sternal notch. Overlapping plastic lateral panels provide lateral stabilizing strength and control.



Code	Size	Height	Circumference
CSTL-NB-2	SMALL	432mm	762 - 889mm
CSTL-NB-3	MEDIUM	432mm	889 - 1016mm
CSTL-NB-4	LARGE	483mm	1016 - 1143mm
CSTL-NB-5	XLARGE	483mm	1143 - 1270mm

CALIFORNIA CEO

The spinal hyperextension orthosis is specifically designed to extend the thoracic spine, prevent forward flexion and limit rotation. The patented compound closure system provides excellent low back support, additional comfort and adjustable relief.



Code	Size	Variation	Length	Waist mm
3380	XSMALL	STANDARD	310-490	640-760
3381	SML/MED	STANDARD	310-490	760-1020
3381.01	SML/MED	EXTENDED	410-580	760-1020
3382	LRG/XL	STANDARD	310-490	1020-1270
3382.01	LRG/XL	EXTENDED	410-580	1020-1270
3383	XXLARGE	STANDARD	310-490	1270-1400

SAE BRACE HINGED

A comfortable, lightweight, universally sized and anatomically shaped, hyperextension orthosis. Featuring adjustable vertical and horizontal anodized aluminium bars, a patient friendly anterior left or right velcro closure with kydex plastic pads lined in soft foam that can be easily modified with a heat gun.



Code	Height mm	Width mm
42215	340 - 460	290 - 330

JEWITT HYPEREXTENSION BRACE

Consists of a fully adjustable anterior aluminium frame with screw-on anterior pads. The posterior thoracolumbar pad is made of polyethylene and padded. Increased range of adjustment allows more accurate positioning of pads.



Code	Length mm	Chest & Waist Circ. mm
23009	340 - 410	C 630-740 W 410-560
23010	410 - 470	C 740-890 W 560-740
23011	470 - 530	C 890-1020 W 740-910
23012	530 - 570	C 1020-1140 W 910-1070
23013	530 - 570	C 1140-1270 W 1070-1270

CYBERTECH CTO

The Cervical Thoracic Orthosis "CTO" incorporates the best design features of past proven systems with today's materials for a functional comprehensive ALL adjustable non-invasive cervical thoracic immobilisation system. Utilised for the treatment and immobilisation disorders involving the cervical and upper thoracic spine.



Code	Size	Chin-Sternal	Mid Neck Circ
CTO-SG-2	SML/MED	63 - 101mm	280 - 380mm
CTO-SG-4	LRG/XLRG	83 - 134mm	380 - 510mm

MINERVA COLLAR

To control flexion, extension, lateral and rotary movement of the cervical spine. The anterior section consists of a sternal plate and a mandibular shell connected by a height adjustable aluminium upright with provision for a trachea tube. The posterior section comprises thoracic and occipital shells again connected by a height adjusted aluminium upright. The anterior and posterior sections are held firmly in place by padded shoulder and lateral chest straps. This orthosis has a removable washable foam lining within the polyethylene shells.



Code	Size
24049	PAEDIATRIC
24050	REGULAR

PADDING KIT

24050B

CYBERTECH MINERVA

Minerva Orthosis incorporates the best design features of past proven systems with today's materials for a functional comprehensive ALL adjustable non-invasive cervical thoracic immobilisation system. Utilised for the treatment and immobilisation disorders involving the cervical and upper thoracic spine.



Code	Size	Chin-Sternal	Mid Neck Circ
MCTO-SG-2	SML/MED	63 - 101mm	280 - 380mm
MCTO-SG-4	LRG/XLRG	83 - 134mm	380 - 510mm

CYBERTECH MINI MINERVA

The unique aluminum multi-post yoked anterior adjustment system allows complete linear adjustability for flexion; extension and neutral positioning and immobilization with unique adjustment pull tabs. The multi-post yoked design accommodates the passage of tracheotomy facilitation. The Y-Chin support members allow width and narrowing customization for rotational control and patient comfort.



The ergonomic molded posterior occipital component encompasses an aluminum adjustable spine element allowing the practitioner full customization and total contact conformity of the cervical arch contours

Code	Size	Chin-Sternal	Mid Neck Circ
MIMO-SG-2	SML/MED	63 - 101mm	280 - 380mm
MIMO-SG-4	LRG/XLRG	83 - 134mm	380 - 510mm

SOMI NECK BRACE

To control flexion, extension, lateral and rotary movements of the cervical spine. The foam lined aluminium sternal section is held firmly in place by shoulder bands with crossover straps attaching anteriorly. The occipital and mandibular pieces are height adjustable (to alter flexion/extension) and both are attached anteriorly on the sternal piece.



Code	Size	Variation
24001	SMALL	CHILD
24002	REGULAR	ADULT

SOFT COLLAR

Serves as a reminder to limit cervical flexion. Soft foam covered in skin tone stockinette with velcro closure.



Code	Size	Length	Height
24026	XSMALL	480mm	65mm
24027	SMALL	480mm	75mm
24028	SMALL/LONG	535mm	75mm
24029	MEDIUM	480mm	90mm
24030	MED/LONG	535mm	100mm
24031	LARGE	535mm	100mm
24032	XLARGE	585mm	100mm

CERVICAL COLLAR PHILADELPHIA STYLE

Flexion, extension and rotary control. Heat mouldable expanded polyethylene collar with Kydex anterior and posterior reinforcements and velcro closure. This collar maintains the head in the desired position of flexion or extension. Lightweight, hygienic, reduces pressure problems usually associated with more rigid types. May be fitted and adjusted while the patient is lying flat in bed. Preforated for ventilation. Trachea opening included for all sizes except 2 1/4" 57mm height.



Code	Size	Height	Neck Circ.
24005A	INFANT		152 - 203mm
24005B	CHILD		203 - 279mm
24006A	SMALL	2 1/4"	250 - 350mm
24007A	MEDIUM	2 1/4"	330 - 380mm
24008	LARGE	2 1/4"	405 - 455mm
24009A	SMALL	3 1/4"	250 - 350mm
24010A	MEDIUM	3 1/4"	330 - 380mm
24011A	LARGE	3 1/4"	405 - 455mm
24012A	SMALL	4 1/4"	250 - 350mm
24013A	MEDIUM	4 1/4"	330 - 380mm
24014A	LARGE	4 1/4"	405 - 455mm
24015A	SMALL	5 1/4"	250 - 350mm
24016A	MEDIUM	5 1/4"	330 - 380mm
24017A	LARGE	5 1/4"	405 - 455mm

PHILI LINER

Liner reduces chafing and helps keep cervical collars clean. Two-way stretch material, latex free.

Code	Size
24017B	SMALL
24017C	MEDIUM
24017D	LARGE



PHILI STABILIZER ADJUSTABLE

The Philadelphia Stabilizer is used with the Phili collar to greatly restrict flexion and extension of the cervical spine. As the patient recovers, the Stabilizer can be removed and the collar can remain to continue to provide the necessary immobilization.

Code
1056933



UNIVERSITY COLLAR

The University Cervical Collar is designed for multipurpose long term patient care, comfort and ease of adjustability. Universally sized with unique patent pending variable height positioning for fast and secure adjustment. Simply pull out the dual adjustment tabs, select height and reengage. The University Cervical Collar meets the specialized needs of each individual with unparalleled comfort, control and compliance.

Code
COLU-NG-U



CYBERSPINE CO

The ergonomic molded posterior occipital component encompasses an aluminum adjustable spine element allowing the practitioner full customization and total contact conformity of the cervical arch contours. Semi-rigid links connect the occipital and mandibular sections with full pivot locking joints improving position and helps eliminate AP/linear shearing migration of components and unique patent pending level locking slide buckles.

Code	Size	Chin-Sternal	Mid Neck Circ
CERV-SG-2	SML/MED	63 - 101mm	280 - 380mm
CERV-SG-4	LRG/XLRG	83 - 134mm	380 - 510mm



CIVIC COLLAR

To provide firm support against flexion, extension and rotation. The Civic Collar is shaped closely to the head. It is made from PVC and closed cell foam to that it is lightweight and comfortable for the patient. Adjustments can be made using the screws and the velcro straps can be constricted. The grey PVC may be heated gently and moulded to the desired contour.



Code	Size	Height	Neck Circ.
24045	SMALL	100mm	100 - 350mm
24046	MEDIUM	115mm	115 - 400mm
24047	LARGE	130mm	130 - 430mm

PLASTIC COLLAR 18M

Flexion control. An adjustable plastic collar, featuring padded edges and velcro closure.



Code	Size	Neck Circ.
24018	SMALL	380mm
24019	MEDIUM	405mm
24020	LARGE	430mm

PLASTIC COLLAR WITH CHIN

In addition to preventing cervical flexion this collar offers some resistance to lateral and rotational movements. An adjustable plastic collar, featuring padded edges and velcro closure. Front height is adjustable, ranging from 90mm to 110mm.



Code	Size	Neck Circ.
24022	SMALL	380mm
24023	MEDIUM	405mm
24024	LARGE	430mm

A-FLEX PROTECTIVE HEADGEAR

An adjustable alternative to other hard protective devices. The A-Flex Protective Headgear is fabricated from a flexible plastic that readily conforms to varying head shapes making it easy-to-fit upon demand. The protection is ideal for low impact forces that are distributed (like a shock absorber) across the entire surface without cracking or penetrating. The headgear is equipped with easy-to-adjust straps minimising the need for chin straps. The posterior strap "locks" the device below the mastoids and accommodates size variations while the lateral straps are used for easy donning and doffing. Lateral straps and side-slits are also optional. The A-Flex features ventilation holes to reduce trapped heat. It can be easily trimmed with a pair of scissors without the need for complex equipment and tools. Custom manufacture also available.




Code	Size	Colour	Circumf.
783	CUSTOM	BLUE	CUSTOM
783.01	PAEDIATRIC SMALL	BLUE	400-430mm
783.02	PAEDIATRIC MEDIUM	BLUE	430-460mm
783.03	PAEDIATRIC LARGE	BLUE	460-480mm
783.04	PAEDIATRIC XLARGE	BLUE	480-510mm
783.05	ADULT XSMALL	BLUE	510-530mm
783.06	ADULT SMALL	BLUE	530-560mm
783.07	ADULT MEDIUM	BLUE	560-580mm
783.08	ADULT LARGE	BLUE	580-610mm
783.09	ADULT XLARGE	BLUE	610-640mm
783.10	ADULT XXLARGE	BLUE	640-660mm
783.21	PAEDIATRIC SMALL	PINK	400-430mm
783.22	PAEDIATRIC MEDIUM	PINK	430-460mm
783.23	PAEDIATRIC LARGE	PINK	460-480mm
783.24	PAEDIATRIC XLARGE	PINK	480-510mm

BACKEZE CUSHION

Contoured pillow to relieve low back pain or fatigue while sitting. Foam towelling covered cushion contoured to fit the lower area of the spine.


Code	Size
46001	UNIVERSAL



CERVICAL ROLL

To provide support and maintain the neck in a position of rest when lying down.

Code	Size
24000	UNIVERSAL



LUMBAR ROLL

To provide support to the lower back and help in reducing lower back pain.

Code	Size
46013B	UNIVERSAL



OVERDOOR TRACTION KIT

Designed for effective traction which may be set up at home in conjunction with most swing-type doors. The kit consists of a 10lb plastic bag, a spreader bar, pulleys and cord, an overdoor traction pulley and disposable head halter.

Code	Size
43007	UNIVERSAL



EXERCISE PULLEY SET

Individual component for the Overdoor Traction Kit.

Code
24320



Other



Upper Limb

Hand/Finger 44

Hand/Wrist 45

Elbow Support 49

Shoulder/Arm 50

Slings 51

Misc 52

FINGER SPLINT Mallet

To passively maintain the distal interphalangeal joint in extension whilst permitting unrestricted movement of the proximal interphalangeal joint.

Code	Size	PIP Circ.
31001	1	40mm
31002	2	45mm
31003	3	50mm
31004	4	55mm
31005	5	60mm
31006	5.5	65mm
31007	6	70mm
31008	7	75mm



Finger Splint Safety Pin WEN 6

To passively extends interphalangeal joints.

Code	Size	PIP Circ.
31014	XSMALL	40 - 55mm
31015	SMALL	55 - 65mm
31016	MEDIUM	65 - 75mm
31017	LARGE	75 - 85mm
31018	XLARGE	85 - 100mm



Finger Splint Safety Pin Modified WEN 6A

To extend interphalangeal joints and allow active flexion of the proximal interphalangeal joint.

Code	Size	PIP Circ.
31019	XSMALL	40 - 55mm
31020	SMALL	55 - 65mm
31021	MEDIUM	65 - 75mm
31022	LARGE	75 - 85mm
31023	XLARGE	85 - 100mm



Finger Splint Knuckle Bender Reverse

To actively extend the proximal interphalangeal joint.

Code	Size	PIP Circ.
31036	XSMALL	40 - 55mm
31037	SMALL	55 - 65mm
31038	MEDIUM	65 - 75mm
31039	LARGE	75 - 85mm
31040	XLARGE	85 - 100mm



Finger Splint Extension With Clockspring

To actively extend finger to thumb.

Code	Size	PIP Circ.
31024	SMALL	55 - 65mm
31025	MEDIUM	65 - 75mm
31026	LARGE	75 - 85mm



Finger Splint Ext. With Wrist Support

Maintains interphalangeal joints in extension whilst allowing flexion of the proximal interphalangeal joints.

Code	Size	PIP Circ.
31027	XSMALL	40 - 55mm
31028	SMALL	55 - 65mm
31029	MEDIUM	65 - 75mm
31030	LARGE	75 - 85mm



Finger Splint Finger Knuckle Bender

To actively flex the proximal interphalangeal joint.

Code	Size	PIP Circ.
31031	XSMALL	40 - 55mm
31032	SMALL	55 - 65mm
31033	MEDIUM	65 - 75mm
31034	LARGE	75 - 85mm
31035	XLARGE	85 - 100mm



Hand Splint Knuckle Bender Flexion

To actively flex metacarpophalangeal joints.

Code	Size	PIP Circ.
31076	XSMALL	40 - 55mm
31077	SMALL	55 - 65mm
31078	MEDIUM	65 - 75mm
31079	LARGE	75 - 85mm
31080	XLARGE	85 - 100mm



Hand Splint Knuckle Bender Reverse

To actively extend the metacarpophalangeal joints.

Code	Size	PIP Circ.
31081	SMALL	55 - 65mm
31082	MEDIUM	65 - 75mm
31083	LARGE	75 - 85mm



Hand Splint Comp. Spring & Elastic

To actively extend the metacarpophalangeal and interphalangeal joints.

Code	Size	PIP Circ.
31084	SMALL	55 - 65mm
31085	MEDIUM	65 - 75mm
31086	LARGE	75 - 85mm



WENIGER 2A-KNUCKLE BENDER

To actively flex the metacarpophalangeal joints and extend the interphalangeal joints.

Code	Size	PIP Circ.
31114	XSMALL	40 - 55mm
31115	SMALL	55 - 65mm
31116	MEDIUM	65 - 75mm
31117	LARGE	75 - 85mm
31118	XLARGE	85 - 100mm

**WENIGER 2B-KNUCKLE BENDER SPLT**

To actively flex the metacarpophalangeal and interphalangeal joints. Can be extended to control the wrist.

Code	Size	MCP Width	PIP Circ.
31119	XSMALL	60mm	40 - 55mm
31120	SMALL	75mm	55 - 65mm
31121	MEDIUM	82mm	65 - 75mm
31122	LARGE	90mm	75 - 85mm
31123	XLARGE	100mm	85 - 100mm

**WENIGER 9A-REV KNUCKLE BENDER**

To actively extend the metacarpophalangeal and interphalangeal joints.

Code	Size	MCP Width	PIP Circ.
31124	XSMALL	60mm	40 - 55mm
31125	SMALL	75mm	55 - 65mm
31126	MEDIUM	82mm	65 - 75mm
31127	LARGE	90mm	75 - 85mm
31128	XLARGE	100mm	85 - 100mm

**ELASTIC WRIST STRAP**

To provide compression and support for mild sprains and strains without restricting motion. An elastic wrist sleeve with velcro adjustment.

Code	Size	Wrist Circ
32000	REGULAR	Up to 170mm
32000A	LARGE	170mm & Up

**ELASTIC WRIST SUPPORT**

To provide compression and support for mild sprains and strains without restricting motion.

Code	Size	Wrist Circ
31197	SMALL	140 - 165mm
31198	MEDIUM	165 - 190mm
31199	LARGE	190 - 215mm
31199A	XLARGE	215mm & Up

**THERMIC WRIST SUPPORT**

To immobilise the wrist but allow full range of motion for fingers and thumb, 3mm neoprene, 2 colour reversible support; blue/beige. Lightweight neoprene ideal for post injury sporting applications.

Code	Size	Wrist Circ
65035	SMALL	110 - 130mm
65036	MEDIUM	140 - 160mm
65037	LARGE	170 - 200mm
65038	XLARGE	210 - 250mm

**PREMIUM WRIST**

To provide support and maintain warmth, assisting in the repair of soft tissue injuries. Neoprene sleeve with thumb cutout.

Code	Size	Wrist Circ
65030	SMALL	110 - 130mm
65031	MEDIUM	140 - 160mm
65032	LARGE	170 - 200mm
65033	XLARGE	210 - 250mm

**WRIST-FLEXION CONTROL**

To provide support to the wrist and when indicated, limit movement of the wrist. An elasticised nylon brace with a removable aluminium palmar spoon, velcro closures dorsally and a unique non-elastic wrist strap. This together with the aluminium palmar spoon, provides a three point pressure system to effectively control wrist motion.

Code	Size	Palm Circ	Orien
32009	SMALL	Up to 190mm	LEFT
32010	MEDIUM	190 - 220mm	LEFT
32011	LARGE	220 - 240mm	LEFT
32012	XLARGE	240mm & Up	LEFT
32013	SMALL	Up to 190mm	RIGHT
32014	MEDIUM	190 - 220mm	RIGHT
32015	LARGE	220 - 240mm	RIGHT
32016	XLARGE	240mm & Up	RIGHT



BASIC CAMP WRIST

To provide support to the wrist and when indicated, limit movement of the wrist. Low profile elasticised support made from a durable and lightweight cotton/lycra mix. It features an aluminium palmar stay to allow for contouring for individual patients.

Code	Size	Orientation
27441XSL	XSMALL	LEFT
27441SL	SMALL	LEFT
27441ML	MEDIUM	LEFT
27441LL	LARGE	LEFT
27441XLL	XLARGE	LEFT
27441XSR	XSMALL	RIGHT
27441SR	SMALL	RIGHT
27441MR	MEDIUM	RIGHT
27441LR	LARGE	RIGHT
27441XLR	XLARGE	RIGHT

**UNO-WHO THUMB SPICA**

To provide support to the wrist and hand, removable stay pod for custom sizing and malleable for a custom fit. Constructed of polyester outer shell with tricot covered foam lining for comfort and durability. Stockinette tongue and dorsal hook and loop closure offer simple application. Useful in the treatment of de Quervain's tendonitis, gamekeeper injuries, scaphoid/navicular injuries and used for early cast removal.

Code	Orientation
HW-429L	LEFT
HW-429R	RIGHT

**BIOSKIN PREMIUM BRACING**

The result of combining innovative designs, based on a firm understanding of biomechanics, with materials that have been designed specifically for orthopaedic and sports medicine applications. The patented BioSkin material provides a high level of compression without limiting activity or creating any discomfort. BioSkin is thin and gives radial compression so it conforms to the curves of the body. BioSkin is also breathable, so perspiration evaporates right through it and is hypo-allergenic eliminating the chance of skin irritation.

BOOMERANG

Restricts motion of the wrist. Universal size, easy attachment with open thumb. Universal left and right.

Code	Size
B54200	UNIVERSAL

**THUMB SPICA**

Restricts MP Joint. Low profile, comfortable thumb spica. Universal left and right. To size measure length of thumb from the base knuckle to the nail, and consult size below.

Code	Size	Thumb Lgth
B54502	XSML-MED	80 - 110mm
B54505	LRG-XXL	110 - 140mm

**WRIST/THUMB SPICA**

Restricts MP Joint. Low profile, comfortable thumb spica.

Code	Size	Orientation
B54521-L	XS -SML	LEFT
B54523-L	MED-LRG	LEFT
B54525-L	XL-XXL	LEFT
B54521-R	XS -SML	RIGHT
B54523-R	MED-LRG	RIGHT
B54525-R	XL-XXL	RIGHT

**BASIC CAMP WRIST WITH THUMB**

Provide wrist and 1st metacarpal control. Low profile elasticised support made from a durable and lightweight cotton/lycra mix. It features an aluminium palmar stay to allow for contouring for individual patients.

Code	Size	Orientation
27442XSL	XSMALL	LEFT
27442SL	SMALL	LEFT
27442ML	MEDIUM	LEFT
27442LL	LARGE	LEFT
27442XLL	XLARGE	LEFT
27442XSR	XSMALL	RIGHT
27442SR	SMALL	RIGHT
27442MR	MEDIUM	RIGHT
27442LR	LARGE	RIGHT
27442XLR	XLARGE	RIGHT

**UNO-WHO WRIST BRACE**

To provide support to the wrist and hand, removable stay pod for custom sizing and malleable for a custom fit. Constructed of polyester outer shell with tricot covered foam lining for comfort and durability. Stockinette tongue and dorsal hook and loop closure offer simple application.

Code	Orientation
HW-428L	LEFT
HW-428R	RIGHT

**BIOSKIN WRIST SIZING**

Size	Wrist Circ
XSmall - Small	100 -150mm
Medium - Large	150 - 200mm
XLarge - XXLarge	200 -250mm

DP2 WRIST SKIN

Provides support to the wrist for carpal tunnel syndrome. Cock-up wrist support, soft circumferential wrist strap from stabilization. Removable palmar and dorsal stays and ease of application with hook and loop closures. DP2 provides 6" length, DP3 provides 8" length.



Code	Size	Orientation
B54231-L	XS - SML	LEFT
B54233-L	MED - LRG	LEFT
B54235-L	XL - XXL	LEFT
B54231-R	XS - SML	RIGHT
B54233-R	MED - LRG	RIGHT
B54235-R	XL - XXL	RIGHT

DP3 WRIST SKIN

B54241-L	XS - SML	LEFT
B54243-L	MED - LRG	LEFT
B54245-L	XL - XXL	LEFT
B54241-R	XS - SML	RIGHT
B54243-R	MED - LRG	RIGHT
B54245-R	XL - XXL	RIGHT

WRIST/HAND IMMOBILISER

To immobilise the wrist but allow full range of motion for fingers and thumb. Heat mouldable polyethylene with soft plastazote lining. Radiolucent, easy application, waterproof and washable.



Code	Size	Orientation	Circumference
31180	SMALL	LEFT	Up to 200
31182	MEDIUM	LEFT	200 - 220
31184	LARGE	LEFT	220 & Up
31181	SMALL	RIGHT	Up to 200
31183	MEDIUM	RIGHT	200 - 220
31185	LARGE	RIGHT	220 & Up

WRIST/HAND/THUMB IMMOBILISER

To immobilise the wrist and the first metacarpal phalangeal joint of the thumb, allowing full range of finger motion. Heat mouldable polyethylene with soft plastasote lining. Radiolucent, easy application, waterproof and washable.



Code	Size	Orientation	Circumference
31188	SMALL	LEFT	Up to 200
31190	MEDIUM	LEFT	200 - 220
31192	LARGE	LEFT	220 & Up
31189	SMALL	RIGHT	Up to 200
31191	MEDIUM	RIGHT	200 - 220
31193	LARGE	RIGHT	220 & Up

DORSAL CARPEL TUNNEL SPLINT

Design to limit wrist flexion and extension during repetitive hand motion. Made of lightweight, yet durable Kydex plastic, the splint can be moulded for patient customising and positioning. Can also be used as a wrist extension splint to support a weak or flaccid wrist.



Code	Size	Orientation	Palm Width
44DCT-ES-L	XSMALL	LEFT	51 - 64
44DCT-S-L	SMALL	LEFT	64 - 76
44DCT-M-L	MEDIUM	LEFT	76 - 89
44DCT-L-L	LARGE	LEFT	89 - 102
44DCT-XL-L	XLARGE	LEFT	102 - 114
44DCT-ES-R	XSMALL	RIGHT	51 - 64
44DCT-S-R	SMALL	RIGHT	64 - 76
44DCT-M-R	MEDIUM	RIGHT	76 - 89
44DCT-L-R	LARGE	RIGHT	89 - 102
44DCT-XL-R	XLARGE	RIGHT	102 - 114

GERIATRIC HAND

To provide stable wrist and finger positioning. Lightweight, durable Kydex construction can be heat moulded to customise. It features an adjustable palmar stay, adjustable web space abduction with foam roll.



Code	Size	Orientation	Palm Width
32GHK-S-L	SMALL	LEFT	64 - 76
32GHK-M-L	MEDIUM	LEFT	76 - 89
32GHK-L-L	LARGE	LEFT	89 - 102
32GHK-XL-L	XLARGE	LEFT	102 - 114
32GHK-S-R	SMALL	RIGHT	64 - 76
32GHK-M-R	MEDIUM	RIGHT	76 - 89
32GHK-L-R	LARGE	RIGHT	89 - 102
32GHK-XL-R	XLARGE	RIGHT	102 - 114

FUNCTIONAL RESTING SPLINT

To maintain a resting position of the wrist, hand and thumb. A plastic palmar support with dorsal velcro closures. The wrist/hand is held in resting position. The fingers and thumb are held in a comfortable, slightly flexed position.



Code	Size	Orientation	Palm Circ
35360	SMALL	LEFT	Up to 170
35361	MEDIUM	LEFT	170 - 190
35362	LARGE	LEFT	195 - 220
35365	SMALL	RIGHT	Up to 170
35366	MEDIUM	RIGHT	170 - 190
35367	LARGE	RIGHT	195 - 220

NIGHT RESTING SPLINT WITH THUMB

To maintain a resting position of the wrist, hand and thumb. To assist in the prevention of flexion contractures and radial or ulnar deviation. A heat mouldable plastic palmar support with a soft polyfoam lining, dorsal velcro closures and thumb rest.



Code	Size	Orientation	Circumference
32025	XSMALL	LEFT	Up to 170mm
32026	SMALL	LEFT	170 - 190mm
32027	MEDIUM	LEFT	195 - 220mm
32028	LARGE	LEFT	220mm & Up
32029	XSMALL	RIGHT	Up to 170mm
32030	SMALL	RIGHT	170 - 190mm
32031	MEDIUM	RIGHT	195 - 220mm
32032	LARGE	RIGHT	220mm & Up

NEUTRAL POSITION RESTING

To maintain a resting position of the wrist, hand and thumb. A plastic palmar support with dorsal velcro closures. The wrist/hand is held in a mid pronation/supination position. The fingers and thumb are held in a comfortable, slightly flexed position.



Code	Size	Orientation	Circumference
32032A	SMALL	LEFT	Up to 170mm
32032B	MEDIUM	LEFT	170 - 190mm
32032C	LARGE	LEFT	195 - 220mm
32032D	XLARGE	LEFT	220mm & Up
32032E	SMALL	RIGHT	Up to 170mm
32032F	MEDIUM	RIGHT	170 - 190mm
32032G	LARGE	RIGHT	195 - 220mm
32032H	XLARGE	RIGHT	220mm & Up

H/S WENIGER 1 COCK-UP SPLINT

To actively dorsiflex the wrist.



Code	Size	MP Width
32033	SMALL	75mm
32034	MEDIUM	82mm
32035	LARGE	90mm

H/S MOD.OPPENHEIMER W/THUMB

Oppenheimer spring wire splint to actively dorsiflex the wrist and thumb attachment to actively abduct the thumb.



Code	Size	MP Width
32036	SMALL	75mm
32037	MEDIUM	82mm
32038	LARGE	90mm

H/S THOMAS SUSPENSION WEN 4

To provide active wrist dorsiflex and thumb abduction for radial palsy.



Code	Size	MP Width
32039	SMALL	75mm
32040	LARGE	90mm

WENIGER 13-COMBINED COILED OPPENHEIMER

An active orthosis consisting of a spring Oppenheimer to dorsiflex the wrist, knuckle bender to flex the metacarpophalangeal joints and outrigger to extend the interphalangeal joints.



Code	Size	MP Width
32055	XSMALL	60mm
32056	SMALL	75mm
32057	MEDIUM	82mm
32058	LARGE	90mm
32059	XLARGE	100mm

WENIGER 1A-COCK-UP SPLINT OUTRIGGER

To actively dorsiflex the wrist and flex the metacarpophalangeal joints.



Code	Size	MP Width
32060	SMALL	75mm
32061	MEDIUM	82mm
32062	LARGE	90mm

WENIGER 1A-COCK-UP SPLINT OUTRIGGER

An active orthosis consisting of a spring Oppenheimer to dorsiflex the wrist, reverse knuckle bender to extend the metacarpophalangeal joints and outrigger to extend the interphalangeal joints.



Code	Size	MP Width
32063	XSMALL	60mm
32064	SMALL	75mm
32065	MEDIUM	82mm
32066	LARGE	90mm
32067	XLARGE	100mm

WENIGER 13B-COMB OPPENH COMPOSITE ELASTIC

An active orthosis consisting of a spring Oppenheimer to dorsiflex the wrist, combined with a spring and elastic splint to extend the metacarpophalangeal and interphalangeal joints.



Code	Size	MP Width
32068	SMALL	75mm
32069	MEDIUM	82mm
32070	LARGE	90mm

TENNIS ELBOW STRAP

This strap prevents the forearm muscles from expanding to their maximum during contraction, thus relieving the tension on the tendon. A nylon support with slight elastic properties that fits immediately distal to the elbow joint. Tension is adjustable by a velcro strap.

Code	Size	Forearm Circ
46007	SMALL	200 - 240mm
46008	MEDIUM	240 - 280mm
46009	LARGE	280 - 330mm



ELASTIC ELBOW

To provide compression and support to the elbow without restricting movement. An elastic pull on elbow sleeve.

Code	Size	Elbow Circ
34001	SMALL	200 - 260mm
34002	MEDIUM	260 - 320mm
34003	LARGE	320 - 380mm
34004	XLARGE	380mm & Up



THERMIC ELBOW

To support and maintain warmth assisting in repair of soft tissue injuries. 3mm neoprene support. 2 colour reversible support. Lightweight neoprene ideal for post injury sporting applications.

Code	Size	Elbow Circ
65025	SMALL	220 - 260mm
65026	MEDIUM	260 - 300mm
65027	LARGE	300 - 320mm
65028	XLARGE	330 - 360mm



PREMIUM ELBOW SUPPORT

The OAPL Premium Elbow Support provides compression and warmth to promote rapid healing of soft tissue injuries. Maintaining body heat helps to increase blood flow, promoting elasticity of muscles, tendons and ligaments. This ventilated neoprene support with moisture wicking technology and soft terry lining. Provides comfortable support which aids in reducing swelling and associated pain.

Code	Size	Elbow Circ
65020	SMALL	230 - 260mm
65021	MEDIUM	270 - 300mm
65022	LARGE	310 - 350mm
65023	XLARGE	360mm & Up



BIOSKIN PREMIUM BRACING

The result of combining innovative designs, based on a firm understanding of biomechanics, with materials that have been designed specifically for orthopaedic and sports medicine applications. The patented BioSkin material provides a high level of compression without limiting activity or creating any discomfort. BioSkin is thin and gives radial compression so it conforms to the curves of the body. BioSkin is also breathable, so perspiration evaporates right through it and is hypo-allergenic eliminating the chance of skin irritation.

STANDARD ELBOW SKIN

Full elbow coverage, comfortable, low-profile fit.

Code	Size
B42100	XSMALL
B42101	SMALL
B42102	MEDIUM
B42103	LARGE
B42104	XLARGE



TENNIS ELBOW SKIN

Silicone pad fits over painful area, full elbow coverage.

Code	Size
B42200	XSMALL
B42201	SMALL
B42202	MEDIUM
B42203	LARGE
B42204	XLARGE



TENNIS ELBOW SKIN

Silicone pad fits over painful area, fits comfortably below elbow.

Code	Size
B42500	XSMALL
B42501	SMALL
B42502	MEDIUM
B42503	LARGE
B42504	XLARGE



BIOSKIN ELBOW SIZING

Size	Forearm Circ
XSmall	200 - 230mm
Small	230 - 260mm
Medium	260 - 280mm
Large	280 - 310mm
XLarge	310 - 330mm

OAPL ROM ELBOW

Features extendable Arms, reducing inventory as there is no need for different length devices. Pinless Hinge, allowing easy adjustment with no parts to be lost. Flexible Cuffs to accommodate large size range providing intimate fit. Also features alligator tabs touch tape closure and easily adjustable strap length to accommodate sizing variations.

Code	Orientation
E-L092L	LEFT
E-L092R	RIGHT



RCAI POP ARM BRACE

Post surgery and/or related injuries where limited ROM or locked control of the elbow is required. Open, lightweight design for patient comfort. This brace is construction from one piece for easy application and removal. The washable “wicking” liner absorbs moistures. Lateral hinge arm prevents rubbing. Flexion settings 0-120 degrees, extension from 0-90 degrees. Plastic humeral and forearm cuffs provide support.



Code	Size	Orien	Mid Humeral
74APOP-S-L	SMALL	LEFT	279 - 356mm
74APOP-S-R	SMALL	RIGHT	279 - 356mm
74APOP-M-L	MEDIUM	LEFT	305 - 406mm
74APOP-M-R	MEDIUM	RIGHT	305 - 406mm
74APOP-L-L	LARGE	LEFT	330 - 457mm
74APOP-L-R	LARGE	RIGHT	330 - 457mm

MACKIE ELBOW

Proven and reliable system for providing incremental stretch to reduce soft tissue contractures of the elbow. Allows patients to control their own progress and maximise treatment effectiveness. Pre-fabricated sizes designed to provide faster patient delivery. Achieve optimal fit and function with minimal inventory. Developed in conjunction with the Mayo Clinic, Rochester, Minnesota.



Code	Size	Orientation	Circumference
S2000	SMALL	LEFT	200 - 300mm
M2001	MEDIUM	LEFT	260 - 360mm
L2002	LARGE	LEFT	320 - 420mm
S2003	SMALL	RIGHT	200 - 300mm
M2004	MEDIUM	RIGHT	260 - 360mm
L2005	LARGE	RIGHT	320 - 420mm

SOFT HUMERAL FRACTURE BRACE

A prefabricated polyethylene, foam lined shell with velcro closures. Of universal design, fitting the right or left arm, this design allows unrestricted range of motion at the elbow and shoulder.



Code	Size	Length	Distal Circ	Proximal Circ
1229X	XSMALL	280mm	150 - 200mm	200 - 250mm
1230	SMALL	300mm	180 - 230mm	230 - 280mm
1231	MEDIUM	330mm	230 - 280mm	280 - 330mm
1232	LARGE	330mm	280 - 330mm	330 - 380mm
1233	XLARGE	330mm	330 - 380mm	380 - 430mm

EXTENDED LENGTH HUMERAL FRACTURE BRACE

A prefabricated foam lined polyethylene, two part shell with proximal shoulder cap, anterior and posterior openings with polyethylene tongues and velcro closures. Adjustable to maintain good containment.



Code	Size	Distal Circ	Proximal Circ
1089	XSMALL	130-180mm	180-230mm
1090	SMALL	180-230mm	230-280mm
1091	MEDIUM	230-330mm	280-330mm
1092	LARGE	280-330mm	330-380mm
1093	XLARGE	330-380mm	380-430mm

SHOULDER BRACE-PIPEER

To retract the shoulders, encouraging an improved posture. The shoulder is made from a ventilated cushioned material featuring posterior stays, padded shoulder straps, webbing chest strap and anterior buckle adjustment.



Code	Size	Chest Circ
33006	XSMALL	610 - 710mm
33007	SMALL	760 - 810mm
33008	MEDIUM	860 - 910mm
33009	MED/LRG	960 - 1010mm
33010	LARGE	1070 - 1120mm
33010A	XLARGE	1170 - 1220mm

EZ ARM ABDUCTION ORTHOSIS

To provide a variable range of abduction, flexion and rotation at the shoulder joint. Lightweight metal orthosis with soft padded lining and webbing straps for suspension. Shoulder abduction and degree of rotation is adjustable as is the angle of elbow flexion.



Code	Size	Wholesale
109300	ADULT	\$700.00

CLAVICLE STRAP

To retract the shoulders. A padded figure-eight straps with posterior velcro fasteners.

Code	Size	Chest Circ
33002	SMALL	750 - 850mm
33003	MEDIUM	850 - 950mm
33004	LARGE	950 - 1050mm
33005	XLARGE	1050mm & Up



ARM ABDUCTION PILLOW

To comfortably maintain the shoulder in 45 degree or 30 degree abduction. Foam pillow contoured to body shape with comfortable velofoam lining. Maintained in position by a circumferential body strap and opposing shoulder strap.



Code	Size	Variation
33015A	UNIVERSAL	30 DEGREE
33015B	UNIVERSAL	45 DEGREE

WILMER CARRYING ORTHOSIS

To prevent glenohumeral subluxation. Lightweight metal, leather and cloth design. Weighs less than 200 grams. Allows passive exo/endo rotation of the humerus. Suppresses oedema of the forearm and hand (postural drainage).



Code	Size	Shoe Size	Forearm Circ
WCO-S-N-L	NORMAL	LEFT	230 - 290
WCO-S-L-L	LARGE	LEFT	280 - 340
WCO-S-N-R	NORMAL	RIGHT	230 - 290
WCO-S-L-R	LARGE	RIGHT	280 - 340

ARM SLING - ECONOMY

To maintain forearm in position of rest, limiting movement of the arm. Mesh sleeve with adjustable webbing strap.

Code	Size
33014	UNIVERSAL



ARM SLING UNIVERSAL IMMOBILISING

To immobilise the shoulder by securing the humerus against the trunk. Maintains the forearm in a position of rest. Can be used to elevate the upper limb.

Code	Size
33019	UNIVERSAL



OAPL SHOULDER IMMOBILISER

The low cost, multipurpose sling featuring, touch tape compatible straps with alligator closure providing ease of adjustment. Padded straps increasing patient comfort, leading to increased compliance. This lightweight cotton sling provides increased comfort due to breathability and conformability of material. The thumb loop reduces slippage of sling and improved wrist support. Sling can be removed without undoing straps to allow easier use. The high attachment point of waist strap provides improved immobilization of shoulder unobstructed wrist movement.



Code	Size	Variation
OAPLSI-P	Paediatric	
OAPLSI-S	SMALL	
OAPLSI-M	MEDIUM	
OAPLSI-L	LARGE	
OAPLSI-PD	Paediatric	Pack 5
OAPLSI-SD	SMALL	Pack 5
OAPLSI-MD	MEDIUM	Pack 5
OAPLSI-LD	LARGE	Pack 5

SHOULDER ABDUCTION SLING

To maintain shoulder abduction, with abduction pillow. Universal left and right sling reduces inventory while ensuring that you have a sling available. Soft foam straps and alligator closure provide patient comfort and ease of adjustment.



Code	Size
OASSS	SMALL
OASSM	MEDIUM
OASSL	LARGE
OASSXL	XLARGE

OAPL COLLAR & CUFF

Features, alligator closure and touch tape compatible straps that provide ease of adjustment, straps can be cut to length and alligator closures applied. No need for sharps such as pins or able ties for closure. Enabling changes in positions due to volume changes or stretching of sling easily adjusted. Padded straps provide increased patient compliance. Dispenser box creates easy storage and access with alligator touch tape housed in its one section.



Coming Soon!

FOAM CUFF & SLING

A continuous length of stockinette covered foam used for arm slings, etc.



Code	Variation	Size	Width
72030	STANDARD	10m Roll	50mm
720305	2 ROLLS	2 x 5m Roll	50mm
72030B	WITH 20 TIES	2 x 6m Roll	50mm
72030GB	PAEDIATRIC	20m Roll	30mm

STRAP ARM SLING

Designed to maintain forearm in position of rest across the chest, limiting motion of the shoulder, elbow and forearm.

Code	Size	Wholesale
33020	UNIVERSAL	\$12.05



ARM ELEVATOR

Designed to elevate the upper limb. Canvas sleeve with velcro closures and webbing strap for elevation.

Code	Size
33040	UNIVERSAL



BECKER SHOULDER HOLSTER

The Becker Shoulder Holster is fully adjustable in all planes and is designed to enhance functional recovery, while effectively controlling shoulder movement. The orthosis may be used postoperatively, or in conditions where the immobilization of the shoulder is the desired goal.

Code	Size
200185	UNIVERSAL



SLING - BOBARTH

To immobilise the shoulder by firmly securing the humerus against the trunk. Comfortable velofoam sling and binder with simple adjustment straps.

Code	Size
33021	UNIVERSAL



ACROMIOCLAVICULAR SPLINT

To align and maintain the acromioclavicular joint.

Code	Size
33001	UNIVERSAL



TOURNIQUET

Unique buckle design providing simple one handed use, plus full stretch elastic webbing band.

Code	Size
55160	UNIVERSAL



WRIST AND/OR ANKLE RESTRAINTS

To restrain patient to bed/ wheelchair etc.

Code	Size
32075	UNIVERSAL



HAND EXERCISER

A hand exerciser that is ergonomically shaped. Designed to fit the contours of the hand, it has a shape and feel that appeals to patients. Available in four strength.

Code	Size	Variation
31200	EXTRA SOFT	PINK
31200	SOFT	BLUE
31200	MEDIUM	GREEN
31200	FIRM	ORANGE
31200	EXTRA FIRM	BLACK



Misc



+oapl.TM

Paediatric

Lower Limb 54

Upper Limb 57

KIDDIEGAIT

Carbon composite with anterior shell, lateral strut and an engineered footplate with open calcaneus. It is designed to allow integration of your custom foot orthotic UCB or SMO to help control ground-up forces. It also features an anterior shell to assist in the management of proximal deficits and improve proprioceptive response. To provide ground reaction control and dynamic assist for toddlers.

Code	Size	Orientation
283721011	SMALL	LEFT
283721012	MEDIUM	LEFT
283721013	LARGE	LEFT
283722011	SMALL	RIGHT
283722012	MEDIUM	RIGHT
283722013	LARGE	RIGHT

28397011	SMALL	SOFTKIT FOR KIDDIEGAIT
28397012	MEDIUM	SOFTKIT FOR KIDDIEGAIT
28397013	LARGE	SOFTKIT FOR KIDDIEGAIT



PAEDIATRIC AFO

This specially designed range of polypropylene orthoses for children maintain the foot and ankle in a neutral position restricting both plantar flexion and dorsiflexion. Medial lateral stability is also achieved. Trim lines can be adjusted to permit varying range of motion.

Code	Size	Length mm	Orien	Calf Circ.
13030	XSMALL	165 - 100	LEFT	150 - 190
13032	SMALL	190 - 115	LEFT	190 - 230
13034	MEDIUM	230 - 120	LEFT	220 - 240
13036	LARGE	250 - 130	LEFT	230 - 265
13038	XLARGE	280 - 150	LEFT	250 - 290
13031	XSMALL	165 - 100	RIGHT	150 - 190
13033	SMALL	280 - 115	RIGHT	190 - 230
13035	MEDIUM	230 - 120	RIGHT	220 - 240
13037	LARGE	250 - 130	RIGHT	230 - 265
13039	XLARGE	280 - 150	RIGHT	250 - 290



CHILD DORSI STRETCH AFO

To immobilise ankle in plantar grade or slightly dorsiflexed position during rest. Formfit padding (allows strategic increases in padding thickness). Lycra lined foam padded straps with centre release buckles. Kydex chassis with EVA walking sole.

Code	Size	Length
27321	SMALL	140mm
27322	MEDIUM	160mm
27323	LARGE	180mm



PAEDIATRIC WALKER

To provide support to the ankle, foot and lower leg during ambulation. Aluminium alloy chassis and uprights. Breathable foam liner. Rocker sole.

Code	Size	Length
27326	SML-MED	152-228mm
27327	LRG-XLRG	191-305mm
OA0200	XSMALL	W 2 - 5



ESSENTIAL JS BOOTS

Flexible inner liner of boot allows them to contour to the foot, providing a more intimate fit. It also provides improved comfort and decreased slippage for smaller sizes. The clear liner enables both parent and clinician to easily access any pressure areas and foot position. The shoes incorporate heel relief and reduce pressure over the heel area. The sizing allows interchanging between JS and Talipes Boot. The JS Boots have identical attachment components for the Denis Browne Bar as the Markell Boots and Talipes Boots.

Code	Shoe Size
11069JS	00000
11070JS	0000
11071JS	000
11072JS	00
11073JS	0
11074JS	1
11075JS	2
11076JS	3
11077JS	4
11078JS	5
11079JS	6
11080JS	7
11081JS	8



DENIS BROWNE BAR

To correct infants rotational problems in the feet, tibia or hips and in the management of Talipes Equino Varus. It may also be used to hold hips in slight abduction. Attaches directly to straight last or equino varus boots or may be fitted to patients own shoes.

Code	Length
11049	100mm
11050	125mm
11051	150mm
11052	200mm
11053	250mm
11054	100 - 280mm



HORTON CLICK BAR

The Horton Click is easier to fit on the child because you fit the shoes on independently then click them on the foot rotation bar. This allows for better positioning of the foot in the shoe which reduces the child from being able to kick the shoes off.



Code	Length
6-8ADJ	6" - 8"
8-12ADJ	8" - 12"

BEBAX BOOT

The Bebax Boot is a prewalker shoe that allows progressive correction of congenital forefoot deformities in all three space planes. The Bebax bootie is made of leather and composed of two soft, pliable plastic inner shells, foamed and connected by a multi-directional hinge.



Code	Length	Orientation
652385	75mm	LEFT
652407	80mm	LEFT
652423	85mm	LEFT
652440	90mm	LEFT
652466	95mm	LEFT
652482	100mm	LEFT
648183	105mm	LEFT
648205	115mm	LEFT
648221	125mm	LEFT
648248	135mm	LEFT
652393	75mm	RIGHT
652415	80mm	RIGHT
652431	85mm	RIGHT
652458	90mm	RIGHT
652474	95mm	RIGHT
652490	100mm	RIGHT
648191	105mm	RIGHT
648213	115mm	RIGHT
648230	125mm	RIGHT
648256	135mm	RIGHT

PAVLIK HARNESS

To maintain hips in desired degree of flexion and abduction, whilst encouraging active motion. The harness fits over the torso and straps embrace the legs. Pelite booties and chest strap are perforated for ventilation. Leg straps are adjustable. Buckle closure.

Code	Size	Chest Circ
16016	NEW BORN	250 - 330mm
16017	SMALL	330 - 405mm
16018	MEDIUM	405 - 500mm
16019	LARGE	500 - 605mm



ORTHOTYKES PAVLIK HARNESS

To maintain hips in desired degree of flexion and abduction, whilst encouraging active motion. Velcro closure. Webbing and felt lined. Open and closed toe arrangement. Markings for straps to prevent confusion.



Code	Size	Chest Circ
16016O	NEW BORN	305 - 356mm
16017O	SMALL	356 - 406mm
16018O	MEDIUM	406 - 457mm
16019O	LARGE	457 - 534mm

RHINO BRACE

To hold hips in the correct abduction position, but flexible enough to allow walking. A moulded plastic orthosis with foam lining. Anterior velcro closures strap around the abdomen and thigh shells.

Code	Size	Waist Circ
16016R	XSMALL	356 - 406mm
16017R	SMALL	406 - 457mm
16018R	MEDIUM	457 - 508mm
16019R	LARGE	508 - 559mm



RCAI PAED ABDUCTION BAR

Designed to provide adjustable abduction and adduction control. Nylon ball and socket joint duplicate the anatomical motion of the hip. Breathable neoprene liners discourage limb migration during gait training or limited ambulation.

Code	Size	Thigh Circ.
26AB-A-N	ADULT	432 - 635mm
26AB-P-N	PAEDIATRIC	305 - 457mm
26AB-P-W	PAEDIATRIC	305 - 457mm
26AB-P-XN	PAEDIATRIC	216 - 330mm



NEWPORT JR PAEDIATRIC ORTHOSIS

The Newport jr Hip System with the new pj, Virtual Paediatric Joint enhances hip stability by supporting and positioning the hips in optimal alignment while allowing controlled ranged of motion. The function during activities of daily living. Improved alignment through the hip encourages appropriate acetabular and femoral modelling and reduces the tendency for subluxation and dislocation.

Code	Size	Variation	Length	Hip Circ.
3919	XSMALL	EX-SHORT	180-220mm	450-480mm
3920	SMALL	SHORT	210-260mm	480-510mm
3921	MEDIUM	SHORT	240-300mm	510-560mm
3922	LARGE	LONG	290-360mm	550-610mm
3923	XLARGE	LONG	350-420mm	610-650mm



ADJUSTABLE ABDUCTION BAR

The abduction bar self-adjusts to the angle of the thigh cuffs and attaches easily. Designed to provide rigidity similar to that provided by a hip spica cast, the Newport Jr with abduction bar is the perfect post-operative orthotic choice. Features of the new abduction bar include: Dual sliding bar design: One bar is slotted for maximum adjustability, and the other has pre-drilled holes for stable fixation. Thumb screw stabilizes the desired setting and can be released quickly to facilitate entry into cars, doorways, etc. Adjustable for up to 30 degrees of abduction on each side.



Code	Size	Variation	Length
3950	XS-SML	SHORT ATTACHED	150-250mm
3950.01	XS-SML	LONG KIT ONLY	150-250mm
3951	MED-LRG-XL	SHORT ATTACHED	230-430mm
3951.01	MED-LRG-XL	LONG KIT ONLY	230-430mm

MAPLE LEAF ORTHOSIS

The re-designed, anatomically contoured proximal section provides greater control of the lumbar spine and pelvic region, improved comfort and patient compliance. Improved design allows the pelvic section to be adjusted by simply loosening the screws and sliding the pelvic section to the desired width. Hip joints have been offset to provide greater clearance. Rotational adjustments can easily be made to accurately position the hip in the transverse plane. The thigh cuff design with extended medial contour provides greater leverage for abducting the patient's hip.



Code	Size	Pelvic Diameter
200147-XS	XSMALL	105 - 190mm
200147-S	SMALL	200 - 240mm
200147-M	MEDIUM	247 - 292mm
200147-L	LARGE	300 - 345mm
200147-C	CUSTOM	CUSTOM
200147-SP	SMALL	PADDING KIT
200147-MP	MEDIUM	PADDING KIT
200147-LP	LARGE	PADDING KIT

SWASH ORIGINAL

The SWASH Original is a hip stabilisation and sitting orthosis, offering significant benefits for the non-ambulatory child. The SWASH works not just for CP but other neuromuscular disorders and for children with mild to severe involvement. The Original version features: Plastic pelvic section extends from L2/L3 to the coccyx, providing maximum torso support for candidates with flaccidity or very limited trunk control. Single strap fixed front closure to secure torso support. Increased posterior and side support. "Childproof" thigh cuff closure. Proven results in several different independent studies. Available with 115° or 123° uprights.



Code	Size	Variation	Wasit Circ.
288000001	SIZE 1	WAIST PART	400 - 465mm
288000002	SIZE 2	WAIST PART	460 - 550mm
288000003	SIZE 3	WAIST PART	540 - 630mm
288000004	SIZE 4	WAIST PART	620 - 720mm
288000005	SIZE 5	WAIST PART	710 - 810mm
288000006	SIZE 6	WAIST PART	800 - 900mm

Code	Size	Variation	Distal Thigh
288030001	SIZE 1	THIGH CUFF	250 - 290mm
288030002	SIZE 2	THIGH CUFF	250 - 290mm
288030003	SIZE 3	THIGH CUFF	290 - 330mm
288030004	SIZE 4	THIGH CUFF	330 - 380mm
288030005	SIZE 5	THIGH CUFF	380 - 440mm
288030006	SIZE 6	THIGH CUFF	440 - 500mm

Code	Size	Variation	Length	Dia.
288010001	SIZE 1	LEG BARS 115°	310mm	6
288010002	SIZE 2	LEG BARS 115°	375mm	7
288010003	SIZE 3	LEG BARS 115°	440mm	7
288010004	SIZE 4	LEG BARS 115°	480mm	7
288010005	SIZE 5	LEG BARS 115°	540mm	7
288010006	SIZE 6	LEG BARS 115°	600mm	7
288020001	SIZE 1	LEG BARS 123°	310mm	6
288020002	SIZE 2	LEG BARS 123°	375mm	7
288020003	SIZE 3	LEG BARS 123°	440mm	7
288020004	SIZE 4	LEG BARS 123°	480mm	7
288020005	SIZE 5	LEG BARS 123°	540mm	7
288020006	SIZE 6	LEG BARS 123°	600mm	7

SWASH LOW PROFILE

The SWASH Low Profile is a hip stabilisation and sitting orthosis, offering significant benefits for the non-ambulatory child. The SWASH works not just for CP but other neuromuscular disorders and for children with mild to severe involvement. The Low Profile version features: Malleable spring steel "iliac roll" that rests over the pelvis and is secured in position by attaching to a small plastic reinforced abdominal pad. Double strap front closure allows symmetrical application. 8 mm diameter uprights in all sizes, extra strong for extra high tone. Interchangeable thigh cuff sizes. Simplified thigh cuff closure for easy application and removal. Contoured posterior joints reduce interference with posterior walkers and other standing/walking/sitting aids. Available with 115° or 123° uprights.



Code	Size	Variation	Wasit Circ.
------	------	-----------	-------------

288350001	SIZE 1	WAIST PART	400 - 465mm
-----------	--------	------------	-------------

288350002	SIZE 2	WAIST PART	460 - 550mm
-----------	--------	------------	-------------

288350003	SIZE 3	WAIST PART	549 - 630mm
-----------	--------	------------	-------------

288350004	SIZE 4	WAIST PART	620 - 720mm
-----------	--------	------------	-------------

288350005	SIZE 5	WAIST PART	710 - 810mm
-----------	--------	------------	-------------

288350006	SIZE 6	WAIST PART	800 - 900mm
-----------	--------	------------	-------------

Code	Size	Variation	Distal Thigh
------	------	-----------	--------------

288420001	SIZE 1	THIGH CUFF	210 - 250mm
-----------	--------	------------	-------------

288420002	SIZE 2	THIGH CUFF	250 - 290mm
-----------	--------	------------	-------------

288420003	SIZE 3	THIGH CUFF	290 - 330mm
-----------	--------	------------	-------------

288420004	SIZE 4	THIGH CUFF	330 - 380mm
-----------	--------	------------	-------------

288420005	SIZE 5	THIGH CUFF	380 - 440mm
-----------	--------	------------	-------------

288420006	SIZE 6	THIGH CUFF	440 - 500mm
-----------	--------	------------	-------------

Code	Size	Variation	Length	Dia.
------	------	-----------	--------	------

288400001	SIZE 1	LEG BARS 115°	310mm	6
-----------	--------	---------------	-------	---

288400002	SIZE 2	LEG BARS 115°	375mm	7
-----------	--------	---------------	-------	---

288400003	SIZE 3	LEG BARS 115°	440mm	7
-----------	--------	---------------	-------	---

288400004	SIZE 4	LEG BARS 115°	480mm	7
-----------	--------	---------------	-------	---

288400005	SIZE 5	LEG BARS 115°	540mm	7
-----------	--------	---------------	-------	---

288400006	SIZE 6	LEG BARS 115°	600mm	7
-----------	--------	---------------	-------	---

288410001	SIZE 1	LEG BARS 123°	310mm	8
-----------	--------	---------------	-------	---

288410002	SIZE 2	LEG BARS 123°	375mm	8
-----------	--------	---------------	-------	---

288410003	SIZE 3	LEG BARS 123°	440mm	8
-----------	--------	---------------	-------	---

288410004	SIZE 4	LEG BARS 123°	480mm	8
-----------	--------	---------------	-------	---

288410005	SIZE 5	LEG BARS 123°	540mm	8
-----------	--------	---------------	-------	---

288410006	SIZE 6	LEG BARS 123°	600mm	8
-----------	--------	---------------	-------	---

PAEDIATRIC LEG & ELBOW WRAP

Velfoam lined wraps with metal stays and Velcro/D ring closure. Designed for paediatric use only. Item comes in either blue or pink, please specify colour choice when ordering.

Code	Size	Length
14018F	6"	150mm
14018G	9"	225mm
14018B	12"	300mm
14018C	14"	350mm
14018D	16"	400mm
14018E	18"	450mm



OAPL SHOULDER IMMOBILISER

The low cost, multipurpose sling featuring, touch tape compatible straps with alligator closure providing ease of adjustment. Padded straps increasing patient comfort, leading to increased compliance. This lightweight cotton sling provides increased comfort due to breathability and conformability of material. The thumb loop reduces slippage of sling and improved wrist support. Sling can be removed without undoing straps to allow easier use. The high attachment point of waist strap provides improved immobilization of shoulder unobstructed wrist movement.



Code	Variation
OAPLSI-P	
OAPLSI-PD	Pack 5

LERMAN SARO

Designed to provide pre and post-operative management of the paediatric shoulder joint. Easily adjustable and lightweight. The adjustment at the shoulder joint and elbow joint allow the clinician to precisely align the child's shoulder joint and upper extremity. The SARO may be used for the right or left.

Code	Size	Variation
200190-S	SMALL	Under 4 Years
200190-M	MEDIUM	4 - 8 Years







Compression Therapy

Compression Stockings 60

Socks 65

Bandages 66

Electro Stimulators 67

Sequential Compression Pumps 67

OAPL TRAVEL SOCKS

OAPL Travel Socks minimise the risk of DVT forming in the legs during travel by increasing circulation. The latex free, hospital grade stockings feature graduated compression, 14-18mmHg. Comprising of 76% Nylon, 24% Spandex these Australian made stockings are available in Black or Beige, knee length. They are suitable for tired, achy legs (heaviness), ankle swelling and moderate to pronounced spider and varicose veins. Hand wash in warm water with mild soap.

Code	Size	Colour	Height	Calf
41032	SMALL	BLACK	Up to 400mm	250 - 300mm
41033	MEDIUM	BLACK	Up to 420mm	300 - 370mm
41034	LARGE	BLACK	Up to 450mm	370 - 440mm
41035	XLARGE	BLACK	Up to 480mm	440 - 500mm
41036	SMALL	BEIGE	Up to 400mm	250 - 300mm
41037	MEDIUM	BEIGE	Up to 420mm	300 - 370mm
41038	LARGE	BEIGE	Up to 450mm	370 - 440mm
41039	XLARGE	BEIGE	Up to 480mm	440 - 500mm



OAPL ANTI-EMBOLISM STOCKING BELOW KNEE

OAPL Anti Embolism stockings are designed to support and aid in venous return, recommended for use by every recumbent patient to prevent DVT Prophylaxis. OAPL Anti Embolism stockings are latex free graduated compression, 18mmHg at the ankle, comprising of 76% Nylon, 24% Lycra Spandex. OAPL Anti Embolism stockings should be used with patients with a history of phlebitis or thrombosis or immobile, aged patients suffering venous insufficiency. Also used post-operatively for oedema management. Stocking can be washed or sterilised (Autoclave at 121 degrees) or hand wash in warm water with mild soap and drip dried.

BELOW KNEE

Code	Size	Variation	Height	Calf
41031J	XSMALL	REGULAR	Up to 400mm	Under 250mm
41031A	SMALL	REGULAR	Up to 400mm	250 - 300mm
41031C	MEDIUM	REGULAR	Up to 420mm	300 - 370mm
41031E	LARGE	REGULAR	Up to 450mm	370 - 440mm
41031G	XLARGE	REGULAR	Up to 480mm	440 - 500mm
41031B	SMALL	LONG	Over 400mm	250 - 300mm
41031D	MEDIUM	LONG	Over 420mm	300 - 370mm
41031F	LARGE	LONG	Over 450mm	370 - 440mm
41031H	XLARGE	LONG	Over 480mm	440 - 500mm



OAPL ANTI-EMBOLISM STOCKING THIGH LENGTH

Code	Size	Variation	Height - Gluteal Fold	Calf
41028A	SMALL	REGULAR	Up to 820mm	250 - 300mm
41028C	MEDIUM	REGULAR	Up to 820mm	300 - 370mm
41028E	LARGE	REGULAR	Up to 820mm	370 - 440mm
41028G	XLARGE	REGULAR	Up to 820mm	440 - 500mm
41028B	SMALL	LONG	Over 820mm	250 - 300mm
41028D	MEDIUM	LONG	Over 820mm	300 - 370mm
41028F	LARGE	LONG	Over 820mm	370 - 440mm
41028H	XLARGE	LONG	Over 820mm	440 - 500mm



OAPL ANTI-EMBOLISM STOCKING with WAIST BELT

Code	Size	Variation	Height - Gluteal Fold	Calf
41027A	SMALL	REGULAR	Up to 820mm	250 - 300mm
41027C	MEDIUM	REGULAR	Up to 820mm	300 - 370mm
41027E	LARGE	REGULAR	Up to 820mm	370 - 440mm
41027G	XLARGE	REGULAR	Up to 820mm	440 - 500mm
41027B	SMALL	LONG	Over 820mm	250 - 300mm
41027D	MEDIUM	LONG	Over 820mm	300 - 370mm
41027F	LARGE	LONG	Over 820mm	370 - 440mm
41027H	XLARGE	LONG	Over 820mm	440 - 500mm



OAPL STOCKING MEASURING GUIDE


Sizing the Stockings

Care should be taken in determining the correct size of the stocking.

Measuring

Knee Length Stockings


- Start by measuring the circumference of the leg at the mid calf 'A'
- Next measure the length of the leg 'H' to 'K' measuring from the heel to the back of the knee. (popliteal fossa)
- **Please Note:** Knee Length stockings should finish 3cm below the popliteal fossa



Knee Length

Thigh Length & Thigh with Waist Belt Stockings

- Measure circumference of the leg at the mid calf 'A'
- Next measure the length of the leg from the heel, 'H' to the gluteal fold 'G'
- **Please Note:** Thigh Length stockings should end at Gluteal fold.



Thigh Length

Important

For best results, instructions for applying the stockings should be carefully followed. Failure to do so could result in restricted circulation or reduced effect.

OAPL stock the full range of Jobst Stockings. For a comprehensive Jobst guide please contact customer service or your local representative.

THERAFIRM WOMEN'S KNEE HIGH 15-20mmHg

Designed to help relieve minor ankle or leg swelling, minor varicosities and leg fatigue. Soft comfort band is non-restrictive while staying in place all day. Reinforced toe and heel for durability. Sheer, two-way stretch nylon and spandex. Available in closed toe style. Colours available are Sand and Black.

Code	Size	Variation
68128	SMALL	BLACK
68124	MEDIUM	BLACK
68125	LARGE	BLACK
68132	XLARGE	BLACK
68200	SMALL	SAND
68201	MEDIUM	SAND
68202	LARGE	SAND
68203	XLARGE	SAND



THERAFIRM 15-20mmHg SIZING

Size	Ankle	Calf	Length
Small	160-180mm	280-340mm	310-460mm
Medium	190-220mm	330-390mm	330-480mm
Large	240-270mm	380-440mm	360-510mm
X-Large	280-290mm	450-470mm	380-560mm
XX-Large	300-320mm	460-500mm	380-560mm

THERAFIRM THIGH HIGH WITH LACE 15-20mmHg

Designed to help relieve minor ankle or leg swelling, minor varicosities and leg fatigue. Over-the-knee style stockings feature non-restrictive lace band with silicone border to keep the stockings in place. reinforced toe and heel for durability. Made from sheer, two way stretch nylon and spandex. Available in closed toe style.

Code	Size	Variation
68450	SMALL	SAND
68451	MEDIUM	SAND
68452	LARGE	SAND
68453	XLARGE	SAND



THERAFIRM 15-20mmHg THIGH SIZING

Size	Ankle	Calf	Thigh	Length
Small	160-180mm	280-340mm	280-340mm	580-640mm
Medium	190-230mm	330-390mm	410-440mm	640-710mm
Large	230-270mm	380-440mm	440-480mm	710-790mm
X-Large	280-290mm	450-470mm	480-530mm	790-860mm
XX-Large	300-320mm	460-500mm	530-580mm	860-910mm

THERAFIRM MEN SOCKS 15-20mmHg

Designed to help relieve minor ankle or leg swelling, minor varicosities and leg fatigue. Stylish men's dress socks feature over-the-calf ribbed style graduated compression with a non-binding, knit-in band that keeps socks in place all day. Features comfort heel and toe. Two-way stretch nylon and spandex. Available in closed toe style.

Code	Size	Variation
68310	SMALL	BLACK
68311	MEDIUM	BLACK
68312	LARGE	BLACK
68340	XLARGE	BLACK
68353	SMALL	BROWN
68313	MEDIUM	BROWN
68314	LARGE	BROWN
68342	XLARGE	BROWN



THERAFIRM MENS SIZING

Size	Ankle	Calf	Length
Small	180-230mm	270-380mm	280-360mm
Medium	230-280mm	280-420mm	310-380mm
Large	250-360mm	310-430mm	330-410mm
X-Large	280-380mm	330-480mm	360-460mm

THERAFIRM KNEE HIGH 20-30mmHG CLASS 1

Available in Closed & Open toe. Class 1 stockings are designed to offer moderate to firm support. These stocking may be prescribed for relief from chronic ankle swelling, prophylaxis or clotting after uncomplicated vein stripping, mild lymphedema and more.

Code	Size	Variation
79320	SMALL	OPEN - BEIGE
79321	MEDIUM	OPEN - BEIGE
79322	LARGE	OPEN - BEIGE
79323	XLARGE	OPEN - BEIGE
79324	XXLARGE	OPEN - BEIGE
79325	XXXLARGE	OPEN - BEIGE
79326	XXXXLARGE	OPEN - BEIGE
75112	SMALL	CLOSED - BEIGE
75113	MEDIUM	CLOSED - BEIGE
75114	LARGE	CLOSED - BEIGE
75115	XLARGE	CLOSED - BEIGE
75116	XXLARGE	CLOSED - BEIGE



THERAFIRM KNEEHIGH SIZING

Size	Ankle	Calf	Length
Small	160-180mm	270-330mm	280-380mm
Medium	190-220mm	290-360mm	310-430mm
Large	220-250mm	320-380mm	330-460mm
X-Large	250-290mm	340-410mm	360-480mm
XX-Large	230-320mm	370-430mm	410-530mm
XXX-Large	320-350mm	420-520mm	410-530mm
XXXX-Large	350-380mm	510-610mm	410-530mm

THERAFIRM KNEE HIGH 30-40mmHG CLASS 2

Available in Closed & Open toe. Class 2 stocking are designed to offer extra firm support. These stocking may be prescribed in cases of severe or chronic venous insufficiency related to acute leg swelling, lymphedema etc.

Code	Size	Variation
76505	SMALL	OPEN - BEIGE
76506	MEDIUM	OPEN - BEIGE
76507	LARGE	OPEN - BEIGE
76508	XLARGE	OPEN - BEIGE
76509	XXLARGE	OPEN - BEIGE
76600K	SMALL	CLOSED - BEIGE
76601K	MEDIUM	CLOSED - BEIGE
76602K	LARGE	CLOSED - BEIGE
76603K	XLARGE	CLOSED - BEIGE
76604K	XXLARGE	CLOSED - BEIGE



OFA BAMBERG LASTOFA BELOW KNEE CLASS 2

Code	Size	Variation
OB-5401-1	1	OPEN
OB-5401-2	2	OPEN
OB-5401-3	3	OPEN
OB-5401-4	4	OPEN
OB-5401-5	5	OPEN
OB-5401-6	6	OPEN
OB-5401-2-BLK	2	OPEN-BLACK
OB-5401-3-BLK	3	OPEN-BLACK
OB-5401-4-BLK	4	OPEN-BLACK
OB-5401-5-BLK	5	OPEN-BLACK
OB-5401-6-BLK	6	OPEN-BLACK
OB-5401-1-CT	1	CLOSED
OB-5401-2-CT	2	CLOSED
OB-5401-3-CT	3	CLOSED
OB-5401-4-CT	4	CLOSED
OB-5401-5-CT	5	CLOSED
OB-5401-6-CT	6	CLOSED



OFA BAMBERG LASTOFA

Lastofa stands for compression stockings that are effective and a pleasure to wear. An all-round talent that is most popular because of its optimal fit. Modern materials ensure enhanced breathability, which is above all very pleasant in summer. Lastofa is also especially recommended for allergy sufferers, as the fine polyamide threads do not allow mould fungus, pollen etc to deposit themselves in the fibre spaces. All stockings are particularly kind to the skin and have therefore been accredited the strict European eco-label "Öko-Tex-Standard 100".

OFA BAMBERG LASTOFA BELOW KNEE CLASS 1

Code	Size	Variation
OB-5201-1	1	OPEN
OB-5201-2	2	OPEN
OB-5201-3	3	OPEN
OB-5201-4	4	OPEN
OB-5201-5	5	OPEN
OB-5201-6	6	OPEN
OB-5201-1-CT	1	CLOSED
OB-5201-2-CT	2	CLOSED
OB-5201-3-CT	3	CLOSED
OB-5201-4-CT	4	CLOSED
OB-5201-5-CT	5	CLOSED
OB-5201-6-CT	6	CLOSED



OFA BAMBERG LASTOFA THIGH HIGH CLASS 1

Code	Size	Variation	Wholesale
OB-5204-1	1	OPEN	\$98.14
OB-5204-2	2	OPEN	\$98.14
OB-5204-3	3	OPEN	\$98.14
OB-5204-4	4	OPEN	\$98.14
OB-5204-5	5	OPEN	\$98.14
OB-5204-6	6	OPEN	\$98.14



OFA BAMBERG LASTOFA THIGH HIGH CLASS 2

Code	Size	Variation
OB-5404-1	1	OPEN
OB-5404-2	2	OPEN
OB-5404-3	3	OPEN
OB-5404-4	4	OPEN
OB-5404-5	5	OPEN
OB-5404-6	6	OPEN
OB-5404-2-BLK	2	OPEN - BLACK
OB-5404-3-BLK	3	OPEN - BLACK
OB-5404-4-BLK	4	OPEN - BLACK
OB-5404-5-BLK	5	OPEN - BLACK
OB-5404-6-BLK	6	OPEN - BLACK



LASTOFA SIZING

Size	Ankle	Calf	Length
1	19.5 - 20cm	28 - 33cm	38.1 - 43cm
2	20 - 22cm	31 - 36cm	38.1 - 43cm
3	22 - 24.5cm	32 - 39cm	38.1 - 43cm
4	24.5 - 27cm	35 - 43cm	38.1 - 43cm
5	27 - 30cm	37 - 46cm	38.1 - 43cm
6	30 - 33cm	39 - 49cm	38.1 - 43cm

OFA BAMBERG LASTOFA COTTON

The robust compression stocking with a high cotton content. The cotton has been processed in such a way that it is on the inside of the stocking and therefore directly on the skin. Particularly suitable for allergy sufferers and patients who prefer a firmer quality. The optimal fit stands for full medical efficacy. Lastofa with cotton is especially kind to the skin and has been accredited the eco-label "Öko-Tex-Standard 100".

OFA BAMBERG LASTOFA COTTON BELOW KNEE CLASS 2

Code	Size	Variation
OB-3401-1	1	OPEN
OB-3401-2	2	OPEN
OB-3401-3	3	OPEN
OB-3401-4	4	OPEN
OB-3401-5	5	OPEN
OB-3401-6	6	OPEN
OB-3401-1-CT	1	CLOSED
OB-3401-2-CT	2	CLOSED
OB-3401-3-CT	3	CLOSED
OB-3401-4-CT	4	CLOSED
OB-3401-5-CT	5	CLOSED
OB-3401-6-CT	6	CLOSED



OFA BAMBERG LASTOFA COTTON BELOW KNEE CLASS 3

Code	Size	Variation
OB-3501-1	1	OPEN
OB-3501-2	2	OPEN
OB-3501-3	3	OPEN
OB-3501-4	4	OPEN
OB-3501-5	5	OPEN
OB-3501-6	6	OPEN
OB-3501-1-CT	1	CLOSED
OB-3501-2-CT	2	CLOSED
OB-3501-3-CT	3	CLOSED
OB-3501-4-CT	4	CLOSED
OB-3501-5-CT	5	CLOSED
OB-3501-6-CT	6	CLOSED



OFA BAMBERG LASTOFA COTTON THIGH HIGH CLASS 2

Code	Size	Variation
OB-3404-1	1	OPEN
OB-3404-2	2	OPEN
OB-3404-3	3	OPEN
OB-3404-4	4	OPEN
OB-3404-5	5	OPEN
OB-3404-6	6	OPEN



OFA BAMBERG LASTOFA COTTON THIGH HIGH CLASS 3

Code	Size	Variation
OB-3504-1	1	OPEN
OB-3504-2	2	OPEN
OB-3504-3	3	OPEN
OB-3504-4	4	OPEN
OB-3504-5	5	OPEN
OB-3504-6	6	OPEN



OFA BAMBERG MEMORY

Memory spoils you with its wonderfully soft wearing comfort. High quality functional yarn makes it feel most friendly to the skin. The most modern knitting technology ensures a perfect fit. You will sense the compression as less strong, while it is still reliably effective.

The soft toe end allows the toes much freedom to move. The anatomically moulded heel adapts perfectly to the shape of your foot. Hardly any constriction in the ankle area nor in the hollow of the knee even if the leg is positioned at 90 degrees.

OFA BAMBERG MEMORY BELOW KNEE CLASS 1

Code	Size	Variation
OB-8101-2	2	OPEN
OB-8101-3	3	OPEN
OB-8101-4	4	OPEN
OB-8201-2	2	OPEN
OB-8201-3	3	OPEN
OB-8201-4	4	OPEN
OB-8101-2-CT	2	CLOSED
OB-8101-3-CT	3	CLOSED
OB-8101-4-CT	4	CLOSED
OB-8201-2-CT	2	CLOSED
OB-8201-3-CT	3	CLOSED
OB-8201-4-CT	4	CLOSED



OFA BAMBERG MEMORY THIGH CLASS 1

Code	Size	Variation
OB-7104-2	2	OPEN
OB-7104-3	3	OPEN
OB-7104-4	4	OPEN
OB-7204-2	2	OPEN
OB-7204-3	3	OPEN
OB-7204-4	4	OPEN



MEMORY SIZING

Size	Ankle	Calf	Length
Small	19.5 - 21.5cm	29 - 35cm	38 - 43cm
Medium	21.5 - 25.5cm	34 - 41cm	38 - 43cm
Large	25.5 - 30cm	37 - 48cm	38 - 43cm

OFA BAMBERG MAX BELOW KNEE

The noble compression stocking for the modern, health-conscious man. The compression stocking that can be hardly distinguished from a normal knee stocking; 50% cotton make Max such a pleasure to wear. On account of their special skin compatibility this line in stockings has been accredited the eco-label "Öko-Tex-Standard 100". Class 1 & 2. Available in black only.



Code	Size	Variation
OB-6501-S	SMALL	CLASS 1
OB-6501-M	MEDIUM	CLASS 1
OB-6501-L	LARGE	CLASS 1
OB-6601-S	SMALL	CLASS 2
OB-6601-M	MEDIUM	CLASS 2
OB-6601-L	LARGE	CLASS 2

MAX SIZING

Size	Ankle	Calf	Length
Small	19.5 - 21.5cm	29 - 35cm	38 - 43cm
Medium	21.5 - 25.5cm	34 - 41cm	38 - 43cm
Large	25.5 - 28cm	37 - 46cm	38 - 43cm

IT STAYS GLUE

Gentle roll-on body adhesive that washes off with water leaving no stain. It is pliable and moves with the skin. Used to support hosiery, socks straps and orthopaedic devices.

Code	Size
112013	60ml



CIRCULATION SOCK

Featuring an unrestrictive top band, hand-linked seams and thickened sole, the PhysiPod Circulation Socks have been designed with ultimate comfort in mind. They are ideal for every day wear as well as for feet prone to swelling and circulatory issues.



Code	Size	Shoe Size	Variation
AH-DS-S-W	SMALL	M 5-7 W 6-8	WHITE
AH-DS-M-W	MEDIUM	M 8-10 W 9-11	WHITE
AH-DS-L-W	LARGE	M 11-3 W 12-14	WHITE
AH-DS-S-N	SMALL	M 5-7 W 6-8	NAVY
AH-DS-M-N	MEDIUM	M 8-10 W 9-11	NAVY
AH-DS-L-N	LARGE	M 11-3 W 12-14	NAVY
AH-DS-S-BR	SMALL	M 5-7 W 6-8	BROWN
AH-DS-M-BR	MEDIUM	M 8-10 W 9-11	BROWN
AH-DS-L-BR	LARGE	M 11-3 W 12-14	BROWN
AH-DS-S-B	SMALL	M 5-7 W 6-8	BLACK
AH-DS-M-B	MEDIUM	M 8-10 W 9-11	BLACK
AH-DS-L-B	LARGE	M 11-3 W 12-14	BLACK

TRACTION SOCKS

OAPL Traction Socks have been designed to increase traction with the floor whilst walking. The increase in traction may help prevent you from falling as a result of slipping when you are unable to wear your normal footwear. The socks feature a slip resistant sole on a lightweight breathable sock. They are a latex free product constructed from 80% Acrylic and 20% Nylon material.



Code	Size	Shoe Size	Colour
CR7	SMALL	PAEDIATRIC	RED
CR8	MEDIUM	M 3-4.5 W 3-6	GREEN
CR11	LARGE	M 5-6.5 W 6.5-7.5	LIGHT BLUE
CR14	XLARGE	M 7-8.5 W 8-9.5	TAN
CR16	XXLARGE	M 9-11.5 W 10-11.5	BLUE
CR17	XXXLARGE	M 12-13.5 W 11.5+	GREY

DIABETIC CREW SOCKS

Ideal for diabetic, arthritic and sensitive feet. Promotes dry, healthy feet, with a non-binding Halo top. Helps prevent skin irritations; patented Smartknit seamless socks eliminate pressure points found in other diabetic or traditional socks.



Code	Size	Shoe Size	Variation
71122	SMALL	M 5-7 W 5-8	WHITE
71132	MEDIUM	M 8-10 W 9-11	WHITE
71142	LARGE	M 11-13 W 12-14	WHITE
71152	XLARGE	M 14 & Up	WHITE
71722	SMALL	M 5-7 W 5-8	BLACK
71732	MEDIUM	M 8-10 W 9-11	BLACK
71742	LARGE	M 11-13 W 12-14	BLACK
71752	XLARGE	M 14 & Up	BLACK

AFO & KAFO SOCKS

Form fitting interface specifically sized for AFO & KAFO use. Protects skin health and increases success of orthosis use. Moisture-wicking fibres keep feet cool, dry and healthy and assist in inhibiting odour. Core-spun, high-stretch fibres hug feet comfortably like a protective second skin. Helps prevent skin irritations and pressure over toes. Patented SmartKnit seamless socks eliminate pressure points found in other traditional socks. Non-binding Halo Top, fits comfortably to minimise indentations. White Only.



AFO SOCKS

Code	Size
1SP1-AF-NF	INFANT
1SP1-AF-CS	CHILD SMALL
1SP1-AF-CR	CHILD REGULAR
1SP1-AF-AS	ADULT SMALL
1SP1-AF-AR	ADULT REGULAR

KAFO SOCKS

1SP1-AF-CL	CHILD
1SP1-NA-XL	NARROW
1SP1-RG-XL	REGULAR

BLUE LINE BANDAGE

High quality elastic bandage indicated for the treatment of varicose veins.

Code	Length
72008	25000mm
72009	3000mm



FULFLEX TUBULAR BANDAGE 10 METRE

A tubular bandage designed to provide compression and support to body segments.



Code	Size	Colour	Width
72010	SIZE A	BEIGE	40mm
72011	SIZE B	BEIGE	65mm
72012	SIZE C	BEIGE	70mm
72013	SIZE D	BEIGE	75mm
72014	SIZE E	BEIGE	85mm
72015	SIZE F	BEIGE	100mm
72016	SIZE G	BEIGE	125mm
72016H	SIZE H	BEIGE	150mm
72017	SIZE J	BEIGE	175mm
72018	SIZE K	BEIGE	215mm
72019	SIZE L	BEIGE	250mm

FULFLEX TUBULAR BANDAGE 1 METRE

A tubular bandage designed to provide compression and support to body segments.



Code	Size	Colour	Width
720101	SIZE A	BEIGE	40mm
720111	SIZE B	BEIGE	65mm
720121	SIZE C	BEIGE	70mm
720131	SIZE D	BEIGE	75mm
720141	SIZE E	BEIGE	85mm
720151	SIZE F	BEIGE	100mm
720161	SIZE G	BEIGE	125mm

SHAPED FULFLEX - ABOVE KNEE

Code	Size	Colour	Length	Width
72011B/D	SIZE B-D	BEIGE	1200mm	650-750mm
72012C/E	SIZE C-E	BEIGE	1200mm	700-850mm
72013D/F	SIZE D-F	BEIGE	1200mm	750-100mm
72014E/G	SIZE E-G	BEIGE	1200mm	850-125mm

SHAPED FULFLEX - BELOW KNEE

Code	Size	Colour	Length	Width
72011B/C	SIZE B-C	BEIGE	700mm	65 - 70mm
72012C/D	SIZE C-D	BEIGE	700mm	70 - 75mm
72013D/E	SIZE D-E	BEIGE	700mm	75 - 85mm
72014E/F	SIZE E-F	BEIGE	700mm	85 - 100mm

FULFLEX TUBULAR BANDAGE 10 METRE NATURAL

Code	Size	Colour	Width
72011N	SIZE B	NATURAL	65mm
72012N	SIZE C	NATURAL	70mm
72013N	SIZE D	NATURAL	75mm
72014N	SIZE E	NATURAL	85mm
72015N	SIZE F	NATURAL	100mm
72016N	SIZE G	NATURAL	125mm
72016HN	SIZE H	NATURAL	150mm

TUBULAR GAUZE BANDAGE 20 METRE



Code	Size	Variation	Width
SL-SM600	00	Small fingers and toes.	8mm
SL-SM601	01	Fingers and toes.	18mm
SL-SM602	12	Bulk dress fingers.	25mm
SL-SM603	34	Child limb.	40mm
SL-SM604	56	Adult limb.	65mm
SL-SM605	78	Large limb.	85mm

FLEXIFAST

Tubular retention bandage designed to provide light compression to limb or torso. 10 metre roll.



Code	Size	Variation	Width	Circ
FFR-10	Child ankles, wrists arms.	RED	35mm	100 - 160mm
FFG-10	Small hands, wrists or limbs.	GREEN	50mm	140 - 240mm
FFB-10	Adult hands, arms, legs.	BLUE	75mm	210 - 360mm
FFY-10	Adult elbows, arms, legs.	YELLOW	110mm	300 - 500mm

VEINOPLUS

To increase blood flow in the legs, reducing venous stasis and minimising the risk of DVT. Conveniently sized stim unit with electrodes included. Suitable for foot and leg oedemas, prevention of DVT, pregnancy, obesity, air and ground travel, prolonged periods of sitting or standing.



Code	Variation	Unit
VEINO	Veinoplus Kit	Each
VEINOE	Replacement Electrodes	P2

MEGO AFEK COMPRESSION THERAPY

A global leader in effective non-drug and non-surgical treatment of blood and lymph flow disorders. Mego Afek specializes in manufacturing non-invasive medical devices, pressure control systems, and high-frequency welded products including inflatable structures and components.

Their medical division develops and produces medical devices for the support and treatment of lymph disorders, circulatory disorders, Angio vascular diseases and prevention of Deep Vain Thrombosis.

Numerous clinical studies are available on the efficacy of these sequential compression devices. Please contact OAPL marketing to request a copy of the studies today. marketing@oapl.com.au

PHLEBO PRESS DVT 601

Phlebo Press® DVT mimics the action of the muscle pump and the foot pump on the veins of the leg. Its unique design provides maximum blood flow enhancement. Featuring; easy to follow colour scheme, making the application process simple and faultless. Minimal noise generated from machine. Optimal compression cycle for maximum blood transition in deep veins. Comfortable and easy to wear, pleasant and smooth massage action. Audio and visual alarm, easy operation and monitoring.



Code
L10000D3

DISPOSABLE SLEEVE CALF

This 4 cell compression garment covers the calf and is connected via the colour coded hose valves to the Phlebo Press DVT pump, providing a compression sequence to secure full drainage of venous blood from the deep veins of the calf. Single use, disposable device constructed from non-allergenic material.



Code	Size	Length	Calf Circ
L4ZD0110	SMALL	260mm	300 - 420mm
L4ZD0120	MEDIUM	330mm	420 - 480mm

REUSABLE SLEEVE CALF & FOOT

This 4 cell compression garment covers the foot and calf and is connected via the colour coded hose valves to the Phlebo Press DVT pump, providing a compression sequence to secure full drainage of venous blood from the deep veins of the foot and calf.



Code	Size	Length	Calf Circ
L4ZD01310	SMALL	380mm	380 - 420mm
L4ZD01320	MEDIUM	450mm	420 - 480mm
L4ZD01330	LARGE	480mm	500 - 540mm

DISPOSABLE SLEEVE CALF & FOOT

This 4 cell compression garment covers the foot and calf and is connected via the colour coded hose valves to the Phlebo Press DVT pump, providing a compression sequence to secure full drainage of venous blood from the deep veins of the foot and calf. Single use, disposable device constructed from non-allergenic material.



Code	Size	Length	Calf Circ
L4ZD0310	SMALL	380mm	300 - 420mm
L4ZD0320	MEDIUM	420mm	420 - 480mm
L4ZD0330	LARGE	480mm	500 - 540mm

DISPOSABLE SLEEVE CALF & THIGH

This 4 cell compression garment covers the calf and thigh and is connected via the colour coded hose valves to the Phlebo Press DVT pump, providing a compression sequence to secure full drainage of venous blood from the deep veins of the calf and thigh. Single use, disposable device constructed from non-allergenic material.



Code	Size	Length	Ankle Circ	Knee Circ	Thigh Circ
L4ZD0430	SMALL	400mm	270 - 330mm	360 - 420mm	400 - 450mm
L4ZD0420	MEDIUM	544mm	320 - 400mm	420 - 500mm	480 - 560mm
L4ZD0410	LARGE	544mm	410 - 460mm	510 - 570mm	580 - 630mm

LYMPHA PRESS MINI

Treating lymphedema effectively means lifelong commitment to a regimen of care. Home therapy with the Lympha Press Mini® is easy and effective, saves clinic visits, and increases patient compliance for better treatment results. Features : Up to 24 overlapping pressure chambers, Fully variable pressure range 20 to 80 mmHg for safe use at home, User friendly - easy to use, Quiet, Compact, easy to store, Sturdy and reliable, Portable, Can treat up to two limbs simultaneously, Compatible with comprehensive range of quality garments, Easy reading controls (Pressure and Treatment duration). Indications: Lymphedema of legs and arms, Post mastectomy arm lymphedema, Healing of venous stasis ulcers, Promotion of wound recovery, Relief from muscular and nervous pain, Relief of discomfort from leg fatigue.



Code
L1000024A

LYMPHA PRESS OPTIMAL

The Lympha Press Optimal® is a pressure controlled device featuring Pretherapy treatment for home care patients as well as for clinics. The Pretherapy Treatment is a gentle massage of the proximal lymphatics, with the purpose of priming the lymph vessels prior to the full treatment, based on MLD (manual lymph drainage) techniques. Pretherapy starts the treatment at the proximal end of the treated limb. Lympha Press® Optimal® then moves on to the main treatment, either sequential build-up pressure, or Wave peristaltic massage. Features: Pressure controlled device, Pretherapy cycle preceding all treatments, Pretherapy cycle override option, Informative LCD screen, Remote control pause switch, Quick-connect connector for hose bundles, Up to 24 overlapping air cells per compression garment, Fully variable pressure range 20 to 90 mmHg, User friendly - easy to use, Quiet, Compact, easy to store, Sturdy and reliable, Portable, Can treat up to two limbs simultaneously, Compatible with most Lympha Press garments, Easy reading controls.



Code	Wholesale
L10002C6	\$2,093.00

COMFY LEG SLEEVE

Comfy Leg Sleeve is a multiple air-cell compression garment for the treatment of full legs - foot, calf & thigh. Use comfy sleeve with Lympha Press Optimal and Lympha Press Mini.



Code	Groin	Mid Thigh	Calf	Ankle
L41L2285	700mm	600mm	470mm	420mm
L41L2295	780mm	680mm	510mm	430mm
L41L2365	810mm	730mm	600mm	560mm
L41L2375	850mm	760mm	610mm	560mm
L41L2385	900mm	800mm	630mm	560mm

COMFY LEG SLEEVE PARALLEL EXPANDER

Comfy Leg Sleeve Parallel Expander allows an increase in the circumference of the Comfy leg Sleeve. Each expander adds 100mm from the top to the bottom of the sleeve.

Code	Use with Sleeve
L415A097	L41L2385
L415A087	L41L2285, L41L2375
L415A101	L41L2295
L415A078	L41L2265

COMFY LEG SLEEVE TAPERED EXPANDER

Comfy Leg Sleeve Tapered Expander allows an increase in the circumference of the Comfy leg Sleeve. Each expander adds 100mm above the knee tapered to 40mm at the ankle.

Code	Size
L415T087	L41L2285, L41L2375
L415T094	L41L2295
L415T078	L41L2365
L415T097	L41L2385

COMFY ARM SLEEVE

The Comfy Arm Sleeve 1 - 75 offer a new and greatly improved way to treat lymphedema in the upper limbs. The treated area does not stop at the upper arm but now includes the armpit as well as the upper shoulder, one side of the back, chest and abdomen.



Code	Size	Orien	Torso Circ
L41LS175L	SMALL	LEFT	1050 - 1150mm
L41LM175L	MEDIUM	LEFT	1110 - 1250mm
L41LL175L	LARGE	LEFT	1220 - 1400mm
L41LS175R	SMALL	RIGHT	1050 - 1150mm
L41LM175R	MEDIUM	RIGHT	1110 - 1250mm
L41LL175R	LARGE	RIGHT	1220 - 1400mm

LYMPHA PANTS

Designed to give best lymph drainage for both legs and lower abdomen. Very usefull and efficient where both or only one leg suffers from lymphedema. One size only expanders available to increse fit.


Code
L4MP2121



LYMPHA PANTS EXPANDER

Allows the expansion of Lympha Pants, each expander offer 100mm increase in circumference. Up to three expanders can be used on each side of pants and are sold in a pair.

Code	Unit
L4Z5LP00	Pair




LYMPHA JACKET

Designed to give best lymph drainage for both arms and upper abdomen. Very usefull and efficient where both or only one leg suffers from lymphedema.



Code	Style	Arm 1	Arm 2	Chest Circ
L4MJ2126	SHORT	800mm	570mm	1050 - 1230mm
L4MJ2136	LONG	960mm	720mm	1050 - 1230mm



VELCRO LEG SLEEVE

Velcro Sleeve FL is a multiple air-cell compression garment for the treatment of full legs - foot, calf & thigh.




Code	Ankle Circ	Thigh Circ
L4131060	300 - 380mm	430 - 560mm
L4131085	300 - 380mm	540 - 640mm
L4132050	320 - 420mm	450 - 570mm
L4132060	320 - 420mm	480 - 600mm
L4132070	300 - 440mm	480 - 620mm
L4132075	340 - 440mm	560 - 640mm
L4132085	340 - 440mm	570 - 670mm
L4132090	340 - 440mm	580 - 700mm
L4132100	340 - 440mm	590 - 730mm
L4133065	430 - 550mm	620 - 740mm
L4133080	420 - 560mm	660 - 780mm
L4133085	420 - 560mm	660 - 800mm
L4133095	420 - 560mm	700 - 840mm
L4133100	420 - 560mm	740 - 880mm
L4133110	420 - 560mm	780 - 910mm
L4134070	650 - 740mm	800 - 940mm
L4134075	600 - 740mm	840 - 980mm
L4134080	640 - 740mm	860 - 1000mm

VELCRO LEG SLEEVE EXPANDER

Allows Extension of the circumference of the Velcro Leg Sleeve.

Code	Length	Width
L4Z5A030	300mm	300mm
L4Z5A035	350mm	300mm

ZIPPER LEG SLEEVE



Code	Length	Ankle Circ	Thigh Circ
4119290	760mm	300 - 400mm	510 - 620mm
L41192A0	770mm	290 - 390mm	530 - 640mm
L41192B0	850mm	270 - 370mm	540 - 650mm
4119370	620mm	400 - 530mm	620 - 740mm
4119380	720mm	400 - 530mm	650 - 770mm
4119390	770mm	400 - 530mm	660 - 790mm
L41193A0	830mm	400 - 530mm	690 - 820mm
L41113B0	910mm	360 - 490mm	680 - 820mm
L41114A0	830mm	700 - 840mm	990 - 1120mm
4112140	320mm	230 - 310mm	300 - 400mm
4112250	370mm	240 - 300mm	370 - 460mm
4112370	490mm	240 - 330mm	370 - 460mm
4112480	560mm	240 - 330mm	390 - 490mm

BOOT SLEEVE

Zipper leg sleeve attachment used to treat the foot area.


Code	Size	Length	Toe Circ	Ankle Circ
L4122111	PAED 1	210mm	180mm	120 - 260mm
L4122110	PAED 2	270mm	220mm	150 - 260mm
L4121110	REG	300mm	280mm	180 - 340mm
L4121210	MEDIUM	280mm	230mm	300 - 460mm
L4121310	LARGE	370mm	280mm	360 - 600mm
L4121410	XLARGE	400mm	320mm	420 - 700mm

ZIPPER SLEEVE EXPANDER

Zipper leg and arm sleeve attachment, with a zipper on each side to increase the overall circumference of the garment by 65mm.

Code	Use with Sleeve	Length
L4159370	4119370, 4112480	500mm
L4159290	4119380	620mm
L4150680	All arm sleeves; 4119290, L41192A0, 4119390	680mm
L4150750	L41193A0, L41114A0	750mm
L4150850	L41113B0	850mm

ZIPPER ARM SLEEVE



Code	Length	Wrist Circ	Bicep Circ
L41111A6	750mm	240 - 320mm	350 - 450mm
L41111A5	750mm	250 - 330mm	400 - 500mm
L41111A4	750mm	260 - 350mm	440 - 550mm
L41111B6	850mm	260 - 350mm	360 - 460mm
L41111B5	850mm	250 - 330mm	410 - 520mm
L41111B4	850mm	250 - 330mm	460 - 560mm

ZIPPER ARM SLEEVE SHOULDER ATTACHMENT

Secures the Zipper Arm Sleeve to the upper torso and compresses the shoulder area.

Code	Size	Orientation
L4171L2S	SMALL	LEFT
L4171L2M	MEDIUM	LEFT
L4171L2L	LARGE	LEFT
L4171R2S	SMALL	RIGHT
L4171R2M	MEDIUM	RIGHT
L4171R2L	LARGE	RIGHT



General Products

Casting 72

Mobility & Rehabilitation 76

Mouthguards 78

Asthma Management 78

Hot/Cold Packs 78

Gel Products 79

X-LITE®

X-LITE is an Airy and Lightweight fabrication material engineered to meet your splinting, bracing, and casting demands. The material is made from a 100% cotton mesh that has been impregnated with a thermoplastic. All reactive substances in the thermoplastic have been cured which gives us a NON-Toxic and Environmental friendly material.

BONDING - No solvents or "roughing " the surface required.

EXCELLENT CONFORMABILITY - Forms easily to the anatomy - no need to stretch or tug.

VARIABLE RIGIDITY - Add layers for additional rigidity or add strips just to areas where additional reinforcement or support is desired.

100% REVERSIBLE - Just return to heated water for a few seconds and the splint will convert to the flat X-LITE material.

NON-TOXIC - No gloves required. No reports of skin irritation.

VENTILATED & COMFORTABLE - The open structure in the material gives great ventilation and comfort to the user.

X-LITE CLASSIC SPLINT

X-LITE sheets come in many different widths and lengths to fit all the bracing and splinting needs that can come across in a day's work. Sheets should be immersed into water heated to 71°C for 1 minute. Working time is about 1 - 1.5 minutes. Cooling time about 2 minutes depending on number of layers.



Code	Colour	Length	Width	Unit
780W1	WHITE	500mm	25mm	50 Sheets
793W	WHITE	375mm	75mm	15 Sheets
794W	WHITE	375mm	100mm	15 Sheets
795W	WHITE	375mm	125mm	15 Sheets
796W	WHITE	375mm	150mm	15 Sheets
790W	WHITE	375mm	200mm	15 Sheets
797W	WHITE	900mm	150mm	15 Sheets
792W	WHITE	900mm	200mm	15 Sheets
798W	WHITE	500mm	450mm	15 Sheets
799W	WHITE	500mm	500mm	15 Sheets
791W	WHITE	1200mm	200mm	15 Sheets
793G	GREEN	375mm	75mm	15 Sheets
794G	GREEN	375mm	100mm	15 Sheets
795G	GREEN	375mm	125mm	15 Sheets
796G	GREEN	375mm	150mm	15 Sheets
790G	GREEN	375mm	200mm	15 Sheets
797G	GREEN	900mm	150mm	15 Sheets
792G	GREEN	900mm	200mm	15 Sheets
798G	GREEN	500mm	450mm	15 Sheets
799G	GREEN	500mm	500mm	15 Sheets
791G	GREEN	1200mm	200mm	15 Sheets

X-LITE CLASSIC SPLINT DISPENSER

Code	Length	Width
XSC0710DW	1000mm	75mm
XSC1010DW	1000mm	100mm
XSC1510DW	1000mm	150mm
XSC2010DW	1000mm	200mm

X-LITE PREMIUM SPLINT SHEETS

Premium very suitable for small adult anatomy, particularly phalanges and paediatrics. Sheets should be immersed into water heated to 71°C for 1 minute. Working time is about 1 - 1.5 minutes. Cooling time about - 2 minutes depending on number of layers.



Code	Length	Width	Unit
XSP0737SS	375mm	75mm	10 Sheets
XSP10372S	375mm	100mm	10 Sheets
XSP15372S	375mm	150mm	10 Sheets
XSP20372S	375mm	200mm	10 Sheets
XSP45600S	375mm	450mm	4 Sheets

X-LITE PREMIUM SPLINT DISPENSER

Code	Length	Width
XSP0710DS	1000mm	75mm
XSP1010DS	1000mm	100mm
XSP1510DS	1000mm	150mm
XSP2010DS	1000mm	200mm

3D-LITE SHEET


3D-LITE is an open weave polyester material impregnated with thermoplastic. The open structure gives an airy and comfortable material for the user that is Non-Toxic and user friendly. 3D-LITE is woven three-dimensionally to provide: stretch in only one direction, rigidity in all directions, optimum draping and mouldability, a smooth surface, less weight and breathability for patient comfort, no shrinkage. 3D-LITE comes in three different thicknesses, each one suitable for different applications:

3mm - Small joint and small anatomy. Excellent for rheumatoid patients.

5mm - When more stability is needed, humeral fracture bracing, tibia fracture bracing, elbow splints, paediatric TLSOs.

7mm - An extremely rigid material. Very suitable for adult TLSO's and other applications.

Where maximum rigidity is required.



Code	Variation	Length	Widt	Unit
3DLITE3	3mm	600mm	600mm	Sheet
3DLITE5	5mm	600mm	600mm	Sheet
3DLITE7	7mm	600mm	600mm	Sheet

X-LITE FINISHING SYSTEM

Same proprietary manufacturing process as X-LITE Classic and X-LITE Premium but cotton gauze replaces the cotton mesh. Since both the splinting material and gauze are cotton, they bond almost instantly to each other. Edging material should be immersed into water heated to 71°C for 10 - 15 seconds. Working time is about 15 - 20 seconds. Cooling time about 5 - 10 seconds.



Code	Colour	Length	Width	Unit
770DG	GREEN	7500mm	18mm	1 Dispenser
771DG	GREEN	7500mm	18mm	1 Dispenser
772DG	GREEN	7500mm	50mm	1 Dispenser
770G	GREEN	500mm	18mm	50 Sheets
771G	GREEN	500mm	25mm	50 Sheets
772G	GREEN	500mm	50mm	50 Sheets
770DW	WHITE	7500mm	18mm	1 Dispenser
771DW	WHITE	7500mm	25mm	1 Dispenser
772DW	WHITE	7500mm	50mm	1 Dispenser
700W	WHITE	500mm	18mm	50 Sheets
771W	WHITE	500mm	25mm	50 Sheets
772W	WHITE	500mm	50mm	50 Sheets

X-LITE THERMOPLASTIC HOOK

An easy solution to attach your pile closure straps to your splint or brace. Bonded to the back of the "hook" is the same cotton gauze used in the Finishing/Edging Tape with the same excellent bonding properties. Thermoplastic hook should be immersed into water heated to 71°C for 15 - 20 seconds. Working time is about 20 - 30 seconds. Cooling time about 15 - 20 seconds.



Code	Length	Width	Unit
AXH0505DW	7500mm	50mm	1 Dispenser
AXH05501W	7500mm	50mm	5 Sheets

X-LITE ADHESIVE FOAM WHITE

3D-SAFE is a self-adhesive soft foam material that adheres easily to X-LITE or 3D-LITE to crease a soft edge around your orthosis.



Code	Length	Width	Unit
495	2000mm	500mm	1 Sheet

X-LITE CLASSIC CAST ROLLS

The only low temperature thermoplastic offered in casting rolls for easy application of cylindrical casts - a practitioner and patient friendly alternative to fiberglass. The rolls have a separating plastic in between layers to prevent sticking together when heated. Available in blue, green, and white. Cast rolls should be immersed into water heated to 71°C for 2 minutes. Working time is about 2 - 3 minutes. Cooling time between 4 - 15 minutes depending on number of layers.



Code	Colour	Length	Width	Unit
103B	BLUE	1800mm	50mm	20 Rolls
133B	BLUE	1800mm	75mm	20 Rolls
143B	BLUE	1800mm	100mm	10 Rolls
153B	BLUE	1800mm	125mm	10 Rolls
163B	BLUE	1800mm	150mm	10 Rolls
103G	GREEN	1800mm	50mm	20 Rolls
133G	GREEN	1800mm	75mm	20 Rolls
143G	GREEN	1800mm	100mm	10 Rolls
153G	GREEN	1800mm	125mm	10 Rolls
163G	GREEN	1800mm	150mm	10 Rolls
103W	WHITE	1800mm	50mm	20 Rolls
133W	WHITE	1800mm	75mm	20 Rolls
143W	WHITE	1800mm	100mm	10 Rolls
153W	WHITE	1800mm	125mm	10 Rolls
163W	WHITE	1800mm	150mm	10 Rolls

X-LITE PREMIUM CAST ROLLS

Premium very suitable for small adult anatomy, particularly phalanges and paediatrics. Sheets should be immersed into water heated to 160°F for 1 minute. Working time is about 1 - 1.5 minutes. Cooling time about - 2 minutes depending on number of layers.



Code	Length	Width	Unit
XCP071W	180mm	75mm	5 Rolls
XCP101W	180mm	100mm	5 Rolls
XCP121W	180mm	125mm	5 Rolls
XCP151W	180mm	150mm	5 Rolls

GYPSONA EXTRA 3

Gypsona Extra 3 bandages combine fast setting, high quality plaster of Paris with the conformable and stable interlocked weave of leno cloth. The bandage demonstrates minimum plaster loss whilst the extra rich creamy texture allows easy handling.



Code	Length	Width	Unit
SN-87071035	3500mm	75mm	Box of 24
SN-87071535	3500mm	100mm	Box of 24
SN-87072035	3500mm	150mm	Box of 12
SN-87077535	3500mm	200mm	Box of 12

GIPSAN PLASTER BANDAGE

Gipsan Bandages are made of high quality plaster and a supporting gauze. The Gipsan products combine the traditional application technique with durability and lightness in weight, hardness and fast drying. It is recommended for all kinds of immobilization of single parts of the body. The modern forms and different sizes of Gipsan products enable quick and qualitative treatment of injuries.



Code	Length	Width	Unit
GI-200-15	3500mm	50mm	Box of 12
GI-200-10	3500mm	100mm	Box of 24

SPECIALIST PLASTER

Specialist Plaster of Paris uses a formula and manufacturing process that produces a consistent plaster with high dry/wet plaster retention and an early green strength. It also offers minimal plaster loss to ensure maximum strength and a smooth cast finish.



Code	Length	Width	Unit
SN-007367	4600mm	100mm	Box of 12
SN-007368	4600mm	150mm	Box of 12

DELTACAST ELITE

Delta-Cast Elite has been designed with ease of application and patient comfort in mind. Our patented fibreglass casting material combined strength and rigidity with soft edges for enhanced patient comfort. Easy to apply, clear radiograph viewing, very comfortable, rigid support.



Code	Variation	Width
72279-00	WHITE	50mm
72279-01	WHITE	75mm
72279-02	WHITE	100mm
72279-03	WHITE	125mm
72275-00	RED	50mm
72275-01	RED	75mm
72275-02	RED	100mm
72276-00	BLUE	50mm
72276-01	BLUE	75mm
72276-02	BLUE	100mm
72277-00	PURPLE	50mm
72277-01	PURPLE	75mm
72277-02	PURPLE	100mm
72278-00	BLACK	50mm
72278-01	BLACK	75mm
72278-02	BLACK	100mm

STS SLIPPER SOCK

The Slipper Sock is an exciting new innovative casting product that will benefit both you and your patients when casting for foot orthotics. With the Slipper Sock and its ultra-fast setting resin, you can now obtain a superior foot mould without the plaster mess in less than 2 minutes! Among its benefits: Reduces casting time for foot orthotics by 50%. No more hassle with plaster dust from separating messy splinting material. No more concerns about plaster clogging drains or getting on your clothing, furniture or floor. No more time wasted holding the foot in "neutral position" waiting for the plaster to harden. No more staff time wasted cleaning your patients' feet and legs or cleaning up the room. No more need for special "casting rooms". No more waiting until the next day to ship the casts to your lab. No more special packing of plaster casts to prevent damage during shipping.



Code	Size	Shoe Size	Unit
903-S	SMALL	W 4-5.5	Box 10
903-M	MEDIUM	W 7-10.5, M 6-9.5	Box 10
903-L	LARGE	W 11-12, M 10-11.5	Box 10
903-XL	XLARGE	M 12-15	Box 10

STS PAEDIATRIC CASTING SOCK

The Paediatric AFO Casting Sock is designed for use in obtaining lower extremity moulds of children. The casting socks are available in four sizes. Reduces casting time. Cure time in less than 4 minutes. More accurate than fiberglass rolls. Virtually eliminates all wrinkling and provides excellent conformity to the foot and ankle contours. Glass-smooth finish is ideal for CAD/CAM. Cutting strip and Channel tube included for safe and innovative removal of toe contractures. No more concerns about plaster clogging drains or getting on your clothing, furniture or floor. No more cast cutter. No more staff time wasted cleaning room. No more need for special "casting rooms". No more waiting until the next day to ship the cast to your central fabrication laboratory. No more special packing of plaster casts to prevent damage during shipping.



Code	Size	Length	Height	Unit
904-S	SMALL	85mm	100mm	Box 10
904-M	MEDIUM	100mm	125mm	Box 10
904-L	LARGE	115mm	160mm	Box 10
904-XL	XLARGE	140mm	190mm	Box 10

STS SOCK

The Fitted Polyester Casting Socks are used for obtaining moulds for custom shoes, boots and ankle foot orthotic devices. The various sizes cover the broad range of applications needed by the practitioner to achieve the best possible fit and finish for their respective moulds. Reduces casting time - Cure time in less than 5 minutes. More accurate than fiberglass rolls. Virtually eliminates all wrinkling and provides excellent conformity to the foot and ankle contours. Glass-smooth finish is ideal for CAD/CAM. Cutting strip and Channel tube included for safe and innovative removal. No toe contractures. No more concerns about plaster clogging drains or getting on your clothing, furniture or floor. No more cast cutter. No more staff time wasted cleaning room. No more need for special "casting rooms". No more waiting until the next day to ship the cast to your central fabrication laboratory. No more special packing of plaster casts to prevent damage during shipping. Other sizes available upon request.



Code	Size	Shoe Size	Variation	Unit
900-M	MEDIUM	W 7-10, M 7-9	ANKLE	Box 10
901-S	SMALL	W 4-6.5	MID LEG	Box 10
902-M	MEDIUM	W 6-10, M 6-9	BERMUDA	Box 10
902-L	LARGE	W 11-12, M 10-12	BERMUDA	Box 10
902-XL	XLARGE	M 12-15	BERMUDA	Box 10

SOFFBAN NATURAL

Soffban Natural tears and tears apart easily, aiding rapid application and reducing the risk of stressing the fracture. It is also autoclavable, making it suitable for use in theatre or sterile procedures.



Code	Length	Width	Unit
SN-71472-00	2700mm	50mm	Box of 12
SN-71473-00	2700mm	75mm	Box of 12
SN-71469-00	2700mm	100mm	Box of 12
SN-71470-00	2700mm	150mm	Box of 12

SOFFBAN SYNTHETIC

Soffban is a comfortable, soft and smooth orthopaedic padding. It features a high loft for maximum cushioning and protection and feathers out easily, giving the padding a smooth finish. Soffban is highly conformable for easy shaping around awkward body contours and joints. The padding is resistant to moisture absorbency, minimising the risk of skin maceration. Soffban is easily tearable for quick, trouble-free application and is also autoclavable, making it suitable for use in theatre or sterile procedures.



Code	Length	Width	Unit
SN-36000025	2700mm	50mm	Box of 24
SN-36000026	2700mm	75mm	Box of 12
SN-36000027	2700mm	100mm	Box of 12
SN-36000028	2700mm	150mm	Box of 6

DELTA DRY

The unique patented technology in Delta-Dry has created a water resistant cast padding which allows easy application by the clinician and affords the patient all the lifestyle benefits of a cast they can get wet.



Code	Length	Width	Unit
SN-73443-00	2400mm	50mm	Box of 12
SN-73443-01	2400mm	75mm	Box of 12
SN-73443-02	2400mm	100mm	Box of 12

MIXTURE ADHESIVE FELT

Semi compressed weave, wool and synthetic fibres.



Code	Thickness	Length	Width
DLT-PD01	5mm	225mm	450mm
DLT-PD02	7mm	225mm	450mm

PURE WOOL ADHESIVE FELT

Semi compressed weave. Pack of 4.



Code	Thickness	Length	Width
DLT-PD00	2mm	225mm	450mm
DLT-PD29	3mm	225mm	450mm
DLT-PD30	5mm	225mm	450mm
DLT-PD31	7mm	225mm	450mm
DLT-PD32	10mm	225mm	450mm

FLEECY WEB EXTRA

Brushed fleece cotton, 12.5% larger sheet size.



Code	Length	Width
DLT-PD05X	225mm	450mm

MOLESKIN

Roll

Code	Size
DLT-PD26	280 x 580mm
DLT-PD27	280 x 2500mm



FILLAUER INDELIBLE PENCILS

Code	Unit
EBE705	Box of 12



PLASTER SHEARS BRUNS

Code	Length	Unit
LR-LPSB23	2300mm	Pair



ELECTRIC PLASTER CUTTING SAW

Features: Swing, Frequency: >12,500ft/min; Range: >5 degrees Power Supply: 220V AC, 50Hz, Input Power: 125W, Blade Diameter: 62MM, Power Cord Length: 2M. Low noise total optimized mechanical structure and parts design. Alloy blade is suitable for constant use. Function of parenchyma protection is safe and reliable. Suitable for cutting macromolecule plastic material.



Code	Blade
LR-LSEL	62mm

BIOFOAM IMPRESSION FOAM

BIO-FOAM® Impression Foam is an easily compressible foam, engineered for accurate orthotic, podiatric and prosthetic modelling. BIO-FOAM Impression Foam accurately captures and reproduces the plantar surface of your patient's foot - in only minutes! A clinically accepted alternative to plaster casting, BIO-FOAM Impression Foam is praised for its ease of handling, durability, cleanliness and unmatched quality. Because of its closed-cell structure, BIO-FOAM Impression Foam is moisture-proof to prevent distortion and contamination. Our strict manufacturing and quality control procedures ensure the highest quality impression foam - consistency and reliability with each and every use. Discount applies if purchasing by the Carton, 6 units per carton.



Code
199329

HELPING HAND EURO STYLE

A lightweight, extension clasper with a squeeze trigger. This product is designed to aid in picking up products and lifting them of high shelves. It is excellent for people suffering from arthritis in the hands, or people who find it difficult to pick things up from the floor.



Code	Length
10209	670mm

LONG SHOE HORN

Long handled plastic shoehorn with easy grip top to assist with putting shoes on for the mobility impaired.



Code	Length
10921/24ABS	600mm

FOLDING WALKING STICK

A lightweight aluminium walking stick that conveniently folds into four for easy storage.

Code	Height
KY-927L	820 -920mm



BUDGET WALKING STICK

A lightweight aluminium walking stick.

Code	Height
KY-939L	820 -920mm



THREE POINT WALKING STICK WITH STOOL

A light weight aluminium walking stick with three points for added stability and a fold down stool seat for resting on the go.

Code	Height
KY911L	860mm



FOUR POINT WALKING STICK

A sturdy base of four points adds stability to this aluminium walking stick.

Code	Height
KY-932L	740-970mm



ALUMINIUM CRUTCHES

Adjustable underarm crutches in lightweight aluminium. Weight capacity 75kg.

Code	Size	Height
KY-925S	SMALL	95 -116mm
KY-925M	MEDIUM	114 -135mm
KY-925L	LARGE	135 - 156mm



ALUMINIUM FOREARM CRUTCHES

Adjustable forearm crutches in lightweight aluminium. Max weight capacity 75kg.

Code	Size
KY-9331L	LARGE ONLY



REPLACEMENT CRUTCH TIPS

Replacement rubber tips for our aluminium crutches.

Code	Size
CTRG19	19mm
CTRG22	22mm
CT225	25mm

**FOLDING WALKING FRAME**

A light weight, folding aluminium walking frame that easily folds together for transportation. Max weight capacity 110kg.

Code	Height	Width
KY-913L	760-850mm	580mm

**FOLDING WALKING FRAME WITH WHEELS**

A light weight, folding aluminium walking frame with two 5inch rubber wheels for added ease of movement. Max weight capacity 110kg.

Code	Height	Width
KY-912L-5	780-960mm	570mm

**SHOWER STOOL**

Stool style aluminium and plastic shower chair with adjustable legs. Max weight capacity 90kg.

Code	Height
KY-797L	400 - 500mm

**SHOWER CHAIR**

A full chair with back and armrests, featuring adjustable height legs, constructed from light weight aluminium and plastic. Max weight capacity 110kg.

Code	Height	Width
KY-792	350-450mm	520mm

**SHOWER CHAIR FOR BATHS**

A full chair with back and armrests, featuring adjustable length legs to sit on a bath tub, constructed from light weight aluminium and plastic. Max weight capacity 110kg.

Code	Width
KY-793S	550mm

**ROLLATOR WALKER**

A lightweight aluminium Rollator walker featuring; Easy, one hand folding. Special lever lock that prevents a persons fingers from "snapback" contact when releasing the brake. Manufactured with wrap-around welds providing additional strength and support. Comes with universal basket that can be mounted under seat or in front of rollator. Max weight capacity 90kg. 30 day warranty.

Code	Height	Variation
KY-9144LB	900 -960mm	BLUE
KY-9144LR	900 -960mm	RED

**BUDGET WHEELCHAIR**

Chromed steel budget wheelchair featuring a seat width of 41cm, fixed arm height and 59cm wheel. Max weight capacity 100kg.

Code
KY-809

**BARIATRIC WHEELCHAIR**

This wheelchair caters for larger patients with a powder coated steel heavy duty frame. Generous 700mm seat width and 500mm seat depth. 590mm PU wheels and heavy duty PU Castors. With a fixed armrest and footrest for added stability. Also features a removable seat section and commode tray. Max weight capacity 150 kgs.

Code	Height	Width
KY-607-70	700mm	900mm

**DELUXE WHEELCHAIR**

A comfortable padded wheelchair, constructed from a lightweight aluminium frame. Featuring 590mm, quick release PU wheel and 190mm heavy duty PU castors. Flip up armrest and detachable footrest aid in comfort and usability. Max weight capacity 110 kgs.

Code	Height	Width
KY-953LQ	500mm	410mm

**TOILET SEAT RAISER**

A plastic over toilet seat, that provides added height for ease of access for the mobility impaired.

Code	Height	Width
LC-1664	100mm	400mm

**FOLDING COMMODE BUDGET**

Lightweight aluminium folding commode chair with plastic armrest and removable bucket with lid. Max weight capacity 110kg.

Code	Height	Width
KY-899	780mm	440mm

**COMMODE CHAIR FOLDING**

Lightweight aluminium folding commode chair with padded armrest and removable bucket with lid. Height adjustable legs. Max weight capacity 110kg.

Code	Height	Width
KY-894	780-880mm	450mm



MOUTHGUARD PRO

Proven Dual Composite Construction with internal shock-absorbing frame offers protection and fit that outperforms conventional mouthguards. Shock absorbing jaw pads help protect against concussions and create a channel for improved breathing. Shock absorbing bumper absorbs and deflects impact away from front teeth. Works well with braces. Guaranteed for your protection.



Code	Size
MG5100Y	YOUTH
MG5100A	ADULT

GEL MAX

Heavy duty rubber exoskeletal shock frame and gel-fit provide maximum protection, fit and comfortable in an easy to fit triple layer design. Breathing channel allows for easy breathing while clenched. Guaranteed for your protection.



Code	Size	Variation
MG6100Y	YOUTH	BLUE
MG6200Y	YOUTH	PINK
MG6100A	ADULT	BLUE
MG6200A	ADULT	PINK

POWER ULTRA

Patented semi-rigid shock transfer core design, provides, ultimate protection, fit and lower jaw stability in a compact size for easier breathing and stability. Compression fit system promotes easy and precise moulding for a super tight fit.



Code	Size
MG7100Y	YOUTH
MG7100A	ADULT

BRACES

With a high percentage of young athletes wearing braces, Shock Doctor's introduction of the silicone Instafit Braces mouthguard has been a welcome innovation to parents and kids everywhere. Constructed of 100% medical-grade silicone for ultimate comfort, this extremely durable mouthguard utilizes an integrated ortho-channel that fits securely over braces, protecting the mouth from laceration. Guaranteed for your protection.



Code	Size
MG4100Y	YOUTH
MG4100A	ADULT

BASIC MOUTHGUARD

A basic clear mouthguard which easily conforms to the users teeth with boil and bite technology.



Code	Size
MG600A	YOUTH
MG600Y	YOUTH

MOUTHGUARD CASE

Keep you mouthguard clean and secure with this sturdy mouthguard case.



Code
MG102C

SILICONE MASKS

Designed for attachment to the Space Chamber aerosol spacer as an alternative to inhaling directly from the mouthpiece. Flexible design to fit the contours of the face. Creates an effective seal to optimise the delivery of medication. Exceptional patient comfort and non-collapsible.



Code	Size	Variation
11717	ADULT	ADULT
11718	CHILD	CHILD
11719	INFANT	INFANT

SPACE CHAMBER

The Space Chamber is designed to deliver more medication into the lungs with significantly less resistance to help patients take control of their asthma and Chronic Obstructive Pulmonary Disease (COPD).



Code	Variation
11711	STANDARD

PEAK FLOW METER

Expertly designed to measure asthma symptoms to enable early treatment and to help asthma patients take control of their asthma and COPD.



Code	Variation
11701	STANDARD
11702	LOW-RANGE

INSTANT ICE COMPRESS

Instant cold relief provided for one hour. Easy use, lightly ouch to activate cold action, shake to distribute and apply to affected area.



Code	Length	Width
SP-1005	210mm	125mm

REUSABLE HOT & COLD PACK

Hot & Cold pack that can be kept proven for cold relief or heated in hot water or the microwave, for 60 seconds on a low setting.



Code	Size	Length	Width
SP-7200	MEDIUM	260mm	105mm
SP-7503P	LARGE	280mm	160mm

SILIPOS DIGITAL PAD

Protects toes and fingers from friction and abrasion. Helps to alleviate footwear pressure on toes. Use where corns, calluses, or hammertoes cause discomfort, packed in a polybag.

Code	Size	Unit
10055	MINI	Pk 12
10065	SML/MED	Pk 12
10075	LRG/XL	Pk 12

**SILIPOS CORN PADS**

Helps to relieve painful corns and calluses. Each pad is 3cm in length and contains 2mm thick mineral oil gel pad. Gel pad allows for custom placement. Packaged in a polybag.

Code	Size	Unit
90055	NARROW	Pk 100
90065	MEDIUM	Pk 100
90075	WIDE	Pk 100

**SILIPOS DIGITAL CAP**

Fully lined with mineral oil gel for all around digital protection. Helps protect against pain caused by pressure and friction. Can be used to help treat nail or skin problems due to dryness. Helps reduce discomfort caused by corns, blisters and ingrown nails size. Packaged in a polybag.

RIBBED KNIT FABRIC

Code	Size	Unit
10195	MINI	Pk 6
10205	SML/MED	Pk 6
10206	LRG/XL	Pk 6

**MESH FABRIC - GREATER COMPRESSION**

10235	SML/MED	Pk 8
-------	---------	------

SILIPOS ALL GEL DIGITAL CAP

Helps protect against pain caused by pressure and friction. Can be used to help treat nail or skin problems due to dryness. Helps reduce discomfort caused by corns, blisters and ingrown nails. Wear with sandals or open-toe shoes. Packaged in a polybag.

Code	Size	Unit
10495	SML/MED	Pk 2
10505	LRG/XL	Pk 2

**DIGITAL CAP WITH TOE SPREADER**

Traditional Digital Cap with added benefit of built-in Gel Toe Spreader. Fully lined with mineral oil gel for all around digital protection. Helps protect against pain caused by pressure and friction. Comfortable and simple solution for patients suffering from multiple ailments. Packaged in a polybag.

Code	Size	Unit
1596	LRG/XL	Pk 2

**CUSHION DIGITAL CAP**

Soft, breathable foam cap provides lightweight comfort and protection. Helps reduce discomfort caused by corns, blisters and ingrown nails. Packed in a polybag.

Code	Size	Unit
1589	SML/MED	Pk 4
1590	LRG/XL	Pk 4

**SILIPOS DIGITAL PADS ON STRIP**

Mineral oil gel pads on a continuous strip. Elastic tubing is easy to cut to preferred size. Once cut, product is turned inside out so gel can be applied to skin. Provides a more cost effective alternative to Digital Pads. Packed in a polybag.

Code	Size	Unit
1105	MINI	17 Pads per Strip, 2 Strips Pk
1106	SML/MED	16 Pads per Strip, 2 Strips Pk
1107	LRG/XL	14 Pads per Strip, 2 Strips Pk

**SILIPOS DIGISTRIP**

Half coated with mineral oil gel strip. Elastic tubing is easy to cut to preferred size. Once cut, product is turned inside out so gel can be applied to skin. Provides a more cost effective alternative to our Corn Pads. Approximately 30 pads per strip (when cut to 3m lengths). Packed in a polybag.

Code	Size	Length	Unit
1100	NARROW	600mm	2 Strips Pkg
1101	MEDIUM	600mm	2 Strips Pkg
1102	WIDE	600mm	2 Strips Pkg

**SILIPOS DIGITUBE PLUS**

Fully lined with mineral oil gel. Provides all around digital protection. Not as bulky as conventional foam tubing. Durable and will not bottom out as foam materials do. Can aid in reducing scar tissue. Cut to desired size. Packaged in a polybag.

RIBBED KNIT FABRIC

Code	Size	Diameter	Length	Unit
10865	NARROW	20mm	150mm	Pk 2
10875	WIDE	25mm	150mm	Pk 2

**MESH FABRIC - GREATER COMPRESSION**

10855	NARROW	20mm	150mm	Pk 3
10845	WIDE	76mm	254mm	Pk 1

SILIPOS GEL TOE PROP

Helps to relieve pressure and pain caused by hammer, claw and mallet toes. Silicone gel designed to absorb shock and provide even pressure distribution. Reinforced toe loop to ensure proper fit. Will not bottom out like foam products. Available in polybag or box packaging.

Code	Size	Orientation	Unit
10415	SMALL	LEFT	Pk 3
10435	LARGE	LEFT	Pk 3
10425	SMALL	RIGHT	Pk 3
10445	LARGE	RIGHT	Pk 3



GEL TOE SEPARATORS

Anatomically designed to fit comfortably between toes. Relieves pressure and friction between toes. Helps correct overlapping toes. Packaged in a polybag.

Code	Size	Unit
11605	SMALL	Pk 15
11615	LARGE	Pk 15



CUSHION TOE SEPARATORS

Relieves pressure and friction between toes. Helps correct overlapping toes. Soft breathable foam provides lightweight comfort and protection. Packaged at a polybag.

Code	Size	Unit
1584	SMALL	Pk 6
1585	MEDIUM	Pk 6
1586	LARGE	Pk 6



GEL TOE SEPARATORS

Conforms and comfortably fits between toes. Larger in size compared to the Gel Toe Separators. Used as an inter-digital wedge for bunion deformity, overlapping toes, or toe drift. Packaged in a polybag.

Code	Size	Unit
11505	SMALL	Pk 4
11515	MEDIUM	Pk 4
11525	LARGE	Pk 4



CUSHION TOE COMB

Separates, cushions and protects toes. Helps relieve discomfort caused by crooked or overlapping toes. Soft breathable foam provides lightweight comfort and protection. Packaged in a polybag.

Code	Size	Unit	Wholesale
1591	UNIVERSAL	Pk 6	\$7.50



CUSHION TOE SPREADER

Conforms and comfortably fits between toes. Helps relieve pain caused by bunions. Soft breathable foam provides lightweight comfort and protection. Packed in a polybag.

Code	Size	Unit
1587	MEDIUM	Pk 4
1588	LARGE	Pk 4



SILIPOS BUNION GEL SLEEVE

Anatomically designed to comfortably conform to the shape of the bunion. Mineral oil gel pad cushions, protects and reduces pressure on affected area. Helps to reduce scarring after surgery. Packaged in a polybag.

Code	Size	Unit
10325	SML/MED	Pk 1
10335	LRG/XL	Pk 1



SILIPOS GEL FIT BUN SLEEVE

Elastic bunion sleeve with gel pad incorporated inside two layers of a comfortable and conforming fabric. The covered gel pad shields the painful and sensitive bony protrusion of the forefoot. Product helps relieve shoe pressure, friction and shear forces on the bunion. May also be worn on the hand to protect the thumb (thenar eminence) from discomfort during gripping activities. Product's design is thin and compact making it perfect for everyday use and fit with most dress, casual or athletic shoes and underneath hosiery. Packaged in a polybag.

Code	Size	Unit
1110	SMALL	Pk 1
1111	MEDIUM	Pk 1



SILIPOS GEL BUNION

Anatomically designed to comfortably conform to the shape of the bunion. Mineral oil gel pad cushions protects and reduces pressure on affected area. Helps to reduce scarring after surgery. Packaged in a small clamshell.

Code	Size	Unit
10085	UNIVERSAL	Pk 5



ARCH SUPPORT PAIR

Elastic sleeve provides gentle compression and comfort to arch of foot. Comfortably lifts and supports the plantar fascia. Helps to relieve heel and arch pain. Packaged in a polybag.

Code	Size	Unit
1595	UNIVERSAL	Pk 2



SILIPOS GEL FOOT COVER

Thin layer of gel covered with fabric offers complete toe protection. Extends to provide protection and cushioning to the metatarsal heads. Designed to be worn comfortably with everyday footwear and underneath hosiery. Ideal for activities that place stress on toes and metatarsals. Packaged in a polybag.

Code	Size	Size	Unit
10275	SMALL	95 x 125mm	Pk 1
10285	MEDIUM	100 x 135mm	Pk 1
10295	LARGE	110 x 150mm	Pk 1



SILIPOS UNI GEL STRAP

Comfortable elastic sleeve with anatomically designed gel pad to cover all five metatarsal heads. Helps reduce pain and discomfort in the forefoot. Serves as a soft tissue supplement for fat pad atrophy. Multi-functional use; cut stitching to fit to instep, bunions, elbow hand or wrist.

Code	Size	Variation	Unit
10155	SML/MED	UNCOVERED	Pk 1
10165	LRG/XL	UNCOVERED	Pk 1
10175	SML/MED	COVERED	Pk 1
10185	LRG/XL	COVERED	Pk 1



SLIM GEL-FIT METATARSAL PAD

Elastic gel metatarsal sleeve with gel pad incorporated inside two layers of a comfortable and conforming fabric. The covered gel pad eases pain under the ball of the foot caused by metatarsalgia, calluses, rheumatism, neuromas or from long hours on your feet. Product design is thin and compact making it perfect for everyday use and fit with most dress, casual or athletic shoes underneath hosiery. Great for use with high heels that put extra pressure on metatarsal heads. Packaged in a polybag.

Code	Size	Unit
1113SIL	SMALL	Pk 1
1114	MEDIUM	Pk 1

**SILIPOS BALL OF FOOT PAD**

Thin fabric covered gel pad (2mm) to protect against abrasion, friction and pressure. Excellent for relieving the burning sensation under the metatarsal heads. May be used as a soft tissue supplement for metatarsalgia. Ideal for calluses and small neuromas. Thin enough to fit in any shoe type. Packaged in a polybag.

Code	Size	Unit
10025	UNIVERSAL	Pk 2

**SILIPOS GEL METATARSAL PAD**

Absorbs shock, impact and provides even distribution of pressure in the forefoot area. Anatomically designed for comfort and support. Helps reduce the burning sensation in the metatarsal heads. Ideal for conditions such as fat pad atrophy, metatarsalgia and Morton's neuroma. Silicone construction with reinforced toe loop to ensure proper fit. Packed in a polybag.

Code	Size	Unit
10465	UNIVERSAL	1 Pair

**GEL MET PAD W/TOE SPREAD**

Anatomically designed for comfort and protection metatarsal heads. Helps reduce pain and burning sensations under metatarsal heads. Ideal for fat pad atrophy, metatarsalgia and Morton's neuroma. Silicone construction with built-toe spreader helps reduce discomfort between toes while securing product in place. Ideal for activities that place stress on toes and metatarsals. Packaged in a box.

Code	Size	Unit
4020	SML/MED	1 Pair
4021	LRG/XL	1 Pair

**OAPL GEL CUSHIONS**

Absorbs shock and pressure on the ball of foot. Gently massages the bottom of the foot. Helps prevent calluses by reducing friction and abrasion. Packed in a polybag.

Code	Size	Unit
4247	UNIVERSAL	1 Pair

**CUSHION BALL OF FOOT PAD**

Absorbs shock and pressure on the ball of foot. Gently massages the bottom of the foot. Helps prevent calluses by reducing friction and abrasion. Packaged in a polybag.

Code	Size	Unit
1583	UNIVERSAL	1 Pair

**SILIPOS ACHILLES HEEL PAD**

Soft conforming gel pad protects the posterior of the heel from pressure, impact and friction. Provides gentle compression and dissipation of pressure around inflamed area. Effective in protecting against discomfort caused by ski boots, skates, cleats, other athletic footwear and shoes or boots. Reducing pain caused by footwear pressure on irregular heel bumps. Packaged in a polybag.

Code	Size	Unit
10385	SML/MED	Pk 1
10395	LRG/XL	Pk 1

**SILIPOS MALLEOLAR GEL SLEEVE**

Designed with two discs of gel moulded on both the medial and lateral malleolae. Cushions and protects the ankle against impact, pressure and abrasion. Provides gentle compression to the ankle. Effective in protecting against the discomfort caused by ski boots, skates, cleats and other athletic footwear. Packaged in a polybag.

Code	Size	Unit
15005SIL	UNIVERSAL	Pk 1

**SILIPOS SOFT SKIN GEL HEEL**

Thin gel pad anatomically moulded to fit the heel. Designed to moisturize and soften dry, cracked skin and help to reduce heel fissures. Beneficial for dermatitis and eczema. Can be used in conjunction with skin moisturizing creams and lotions. Recommended as a non-ambulatory therapy product not to be worn with shoes. Packaged in a polybag.

Code	Size	Unit
15265	UNIVERSAL	Pk 1

**SILIPOS HEEL/ELBOW SLIPOVER**

Versatile mineral oil gel can be positioned almost anywhere around foot and heel. Aids in the healing and skin conditioning of dry crackled heels and elbows. Effective in helping to prevent heel and elbow ulcerations. Packaged in a polybag.

Code	Size	Unit
15225	MED/SML	Pk 1
15235	LRG/XL	Pk 1

**BOOT BUMPER**

Helps protect the foot and ankle from abrasions without the added bulk. Effective in protecting against discomfort caused by ski boots, skates, cleats and other athletic footwear. Provides light compression to the foot and ankle. Packaged in a plastic clamshell.

Code	Size	Unit
1770	SML/MED	Pk 1
1772	LRG/XL	Pk 1



SILIPOS ARTHRITIC GEL SOCK

Single layer sock with 2mm of gel throughout the sole. Helps prevent foot ulcerations. Non-restrictive cuff. Can be used with or without diabetic shoes. Provides optimum protection for the neuropathic foot. Helps reduce friction, abrasion and shear forces preventing callus formation. Serves as a soft tissue supplement for fat pad atrophy. Packaged in a box.

Code	Shoe Size	Colour	Unit
1701	M 7-9	White	1 Pair
1702	M 9-11	White	1 Pair
1703	M 11-13	White	1 Pair
1706	M 7-9	Black	1 Pair
1707	M 9-11	Black	1 Pair
1708	M 11-13	Black	1 Pair



SILIPOS SOFT SOCK WHITE

Double layer sock with 6mm of gel throughout the sole. Excellent for any condition that requires a soft, friction absorbing foot bed. Indicated for neuropathic & rheumatoid arthritis sufferers. Fluid like properties of the gel help massage and comfort the foot while walking. Effecting in reducing shock and impact on the heel. Serves as a soft tissue supplement for fat pad atrophy. Requires extra depth footwear. Packaged in a polybag.

Code	Size	Shoe Size	Unit
10605	SMALL	M 7-9	1 Pair
10615	MEDIUM	M 9-11	1 Pair
10625	LARGE	M 10-13	1 Pair
10635	XLARGE	M 13-16	1 Pair



SILIPOS PARTIAL FOOT SOCK

For use with transtmetatarsal amputees. CoolMax Material provides elasticity, comfort and breathability. May be used with foot orthoses or shoes with fillers. Available in 2 sizes. Mineral oil gel covers site of amputation. Packaged in a polybag.


Code	Size
10805	REGULAR
10815	EXTRA LONG



SILIPOS HALF COATED BODY SLEEVE

New and improved fabric for better fit and flexion. Knit sleeving with moisturising mineral oil gel coating on one side of the sleeve only. Available in foot rolls allowing sleeve to be dispensed and cut to desired length. Excellent for use on upper and lower extremity allowing for man universal uses. Miniaturisation and mild compression attributes make this product excellent for: skin grafts, burns, scar treatment and Lymphedema applications. Packaged in a polybag.

Code	Length	Width	Unit
10915	3100mm	83mm	1 Roll
10925	3100mm	100mm	1 Roll
10935	3100mm	135mm	1 Roll



SILIPOS GEL E ROLL

Elastic bandage with one side coated in a thin layer of mineral oil gel. Conforms to difficult or irregular surfaces of the body and will not slide. Gel moves with the skin and does not promote skin irritation. Packaged in a polybag. May be used; to reduce scar tissue, for cold therapy, to provide compression, as a residual limb shrinker, to aid in muscle strains or arthritic pain, to moisturise areas of dry skin, for burn care and to control oedema.

Code	Length	Width	Unit
10355	1200mm	100mm	Pk 1



POST OPERATIVE CARPAL

Anatomically designed gel pad moulded to elastic sleeve contours to shape of palm. Helps to prevent and reduce shock, impact and vibration. Reduces discomfort caused by hypersensitivity. For placement directly over carpal tunnel post-operative scar area. Mineral oil gel pad helps prevent and reduce scar tissue. Can be worn under gloves and/or splints. Packaged in a polybag.

Code	Size	Orientation	Unit
14115	SMALL	LEFT	Pk 1
14135	LARGE	LEFT	Pk 1
14125SIL	SMALL	RIGHT	Pk 1
14145	LARGE	RIGHT	Pk 1



SILOPAD BODY DISC

Multi-purpose discs of mineral oil gel. Apply to parts of the body suffering from discomfort caused by pressure and friction. Can be trimmed to customised size with scissors. Mineral oil gel moisturises and conditions skin. Packaged in a polybag.

Code	Diameter	Unit
10305	63mm	Pk 2
10315	100mm	Pk 2



CARPAL GEL PADS

Anatomically designed gel pad (same as in Carpal Gel Sleeve) with self-adhesive backing. Can be incorporated into custom made or pre-fabricated hand wrist braces. Packaged in a polybag.

Code	Lengt	Width	Unit
14205	152.4mm	50mm	Pk 3



SILIPOS GEL SQUARES


Versatile solution for discomfort related to pressure and friction. Gel square is 3mm thick and 10x10cm in size. Can be trimmed to customised size with scissors. Mineral oil gel moisturises and conditions skin. Can be used as an interface between the skin and a brace. Packaged in a polybag.

Code	Size	Variation	Unit
15505	100x100mm	SELF-ADHESIVE BACKING	Pk 2
15515	100x100mm	NO BACKING	Pk 2



SILIPOS PRESSURE RELIEF PADDING


2mm thick gel pad with self-adhesive backing. Used to reduce shock, impact and vibration. Can be used as a replacement for moleskin or felt products. Excellent use with custom ankle foot orthosis devices. Perfect cushion protection in hand splints. Packaged in a cardboard tube.



Code	Length	Width	Unit
11205	910mm	100mm	Pk 1
11215	910mm	200mm	Pk 1

GEL CARE ADVANCED

Thick layer of self-adhesive mineral oil gel with conforming material backing. Available in discs, sheets or strips. Moisturises the skin and helps to reduce scar tissue. Hypoallergenic, washable and reusable. Cost effective solution for scar management. Excellent for use in the prevention and treatment of hypertrophic scars, keloid scars, skin barrier repair, injury, surgery, burns, acne and stretch marks.



Code	Variation	Length	Width	Unit
GEL-610	SHEET	150mm	120mm	Pk 1
GEL-611	SHEET	150mm	120mm	Pk 10
GEL-613	STRIPS	250mm	40mm	Pk 6
GEL-615	DISCS	40mm	40mm	Pk 6

GEL CARE ADVANCED

Thick layer of self-adhesive mineral oil gel with conforming material backing. Available in tow sizes. Moisturises the skin and helps to reduce scar tissue. Hypoallergenic, washable and reusable. Flexibility of application on hard to reach and bendable areas of the body. Excellent for use in the prevention and treatment of hypertrophic scars, keloid scars, skin barrier repair, injury, surgery, burns, acne and stretch marks.



Code	Length	Width	Unit
GEL-620	910mm	75mm	Pk 1
GEL-621	910mm	150mm	Pk 1

SILIPOS GEL SHEET

2mm thick layer of self-adhesive mineral oil gel with conforming material backing. Moisturises the skin and helps to reduce scar tissue. Hypoallergenic, washable and reusable. Cost effective solution for scar management.



Code	Length	Width	Unit
760	100mm	100mm	Pk 1
761	200mm	140mm	Pk 1
762	250mm	40mm	Pk 2

SILIPOS POST SURGERY BREAST

0.5mm thick layer of self-adhesive mineral oil gel with conforming nylon backing. For use after breast augmentation procedures. Allows patient to wear clothing next to hyper sensitive skin areas. Packaged in a reusable tray.



Code	Size	Length	Width	Unit
740	SMALL	914mm	812mm	Pk 1
741	LARGE	1067mm	914mm	Pk 1

SILIPOS POST SURGERY ABDOMEN

0.5 mm thick layer of self-adhesive mineral oil gel with conforming nylon backing. For use after abdominoplasty procedures. Allows patient to wear clothing next to hyper sensitive skin areas. Packaged in a reusable tray.



Code	Size	Length	Width	Unit
744	SMALL	60mm	44mm	Pk 1
745	LARGE	90mm	45mm	Pk1

SILIPOS DIGITAL CARE KIT

A convenient and economical way to store, display and dispense our most popular digital care products. Special compartments for each product allows you to easily select the best solution for the condition you are treating.



Code	Unit
30365	1 Kit

ALL GEL CARE KIT

Offers an assortment of our most popular gel foot care products in a compact kit for easy in-office dispensing. An excellent compliment to the original Digital Care Kit.



Code	Unit
31744	1 Kit

PHYSIPOD GEL DIGITAL DOT PADS

PhysiPod Gel Digital Dot Pads Are easy to wear fabric covered gel pads, which can be applied to clawed and hammered toes. The elastic fabric cover conforms and stretches to most toes and easily holds the soft gel pad in place. Gel Digital Dot Pads absorb pressure and friction, to cushion and protect from painful corns, callouses and blisters.



Code	Unit
AH-GDDP	Pk 2

PHYSIPOD DIGITAL GEL CAP

PhysiPod Digital Gel Caps are designed to protect the apices of fingers and toes from direct pressure. The elastic fabric cover conforms and stretches to most fingers and toes and envelops the digit with a thick soft gel button on the apex. Ideal for use in clawed and hammered toes and fingers, by cushioning the tip of the digit.



Code	Unit
AH-DGC	Pk 2

PHYSIPOD GEL TOE SEPARATORS

PhysiPod Gel Toe Separators are designed to separate small toes and relieve pressure and friction from clawed and hammered toes, or tight footwear. They gently separate the toes and absorb pressure, reducing irritation to the toes and the interdigital lesions. They are conveniently made and supplied in two sizes.



Code	Size	Unit
AH-GTS-S	SMALL	Pk 2
AH-GTS-L	LARGE	Pk 2

PHYSIPOD GEL TOE WEDGE

The PhysiPod Gel Toe Wedges are designed to fit between the first and second toe, to straighten and re-align overriding toes and relieve pressure from bunions. The soft gel palliates the toes and also interdigital lesions, due to the moulded wedge shape. Painful lesions are also softened by the mineral oil properties released to the skin.



Code	Unit
AH-GTW	Pk 1

PHYSIPOD DIGITAL TOE PROP

PhysiPod Toe Props are designed to lift and palliate the apices of the lesser toes. The tapered design follows the web space of the toes and is ideal for the clawed and hammered toe deformities. Made of a medium durometer Silicone, the toe prop eases pain and discomfort to corns and callous by lifting and supporting the toes. Sold separately as left and right.



Code	Orientation	Unit
AH-SDTP-L	LEFT	Pk 1
AH-SDTP-R	RIGHT	Pk 1

PHYSIPOD BUNIONETTE GEL SHIELD

The PhysiPod Bunionette Gel Shield is tailored to protect and palliate the last toe area from friction and pressure from footwear. The thin and contoured shell conforms to the Bunionette and fits most footwear. Easy application and durable, the PhysiPod Bunionette shield can be worn on either foot.



Code	Unit
AH-BGS-5	Pk 1

PHYSIPOD BUNION GEL SHIELD

The PhysiPod Bunion Gel Shield is tailored for the painful bunion, by offering a shaped and contoured shell to protect and palliate the joint from the shoe. Conveniently thin and tapering, the bunion shield can be worn in most fashion footwear. Easy application and durable, the PhysiPod Bunion Gel Shield relieves friction and pressure from footwear pressure, and may be worn on either foot.



Code	Unit
AH-BGS-1	Pk 1

PHYSIPOD BUNION SLEEVE

The PhysiPod Silicone Bunion Sleeve offers a soft gel disc within a flesh coloured sleeve to cushion and protect painful bunions. The gel disc is held in place by the soft and conforming elasticised sleeve. Ideal for most footwear due to a new thin design.



Code	Unit
AH-SBS	Pk 1

PHYSIPOD FOREFOOT CUSHION SLEEVE

The Silicone Forefoot Cushion Sleeve offers a soft get disc across the whole forefoot, to protect against friction and pressure. The gel disc serves as a soft tissue supplement where there is atrophy of the plantar fat pad or bony prominences. The gel disc is held in place by a soft elasticised flesh coloured sleeve covering the forefoot. Ideal for most footwear due to a new thin design.



Code	Unit
AH-SFCS	Pk 2

**PHYSIFEET INVISIBLE
GEL DOTS**


As the name suggests, PhysiFeet Invisible Gel Dots are clear gel dots, which can be placed over pressure areas where footwear rubs or irritates the foot. Easily placed, with sticky side to the shoe, the gel dots protect and cushion the foot. Thin and small, PhysiFeet Invisible Gel Dots are ideal for difficult and tight spots.



Code	Unit
AH-IGD	Pk 6

**PHYSIFEET INVISIBLE
HEEL GUARD**

PhysiFeet Invisible Heel Guard discreetly sits on the counter of the heel to cushion and protect from footwear pressure and friction. The Heel Guard prevents slipping and rubbing as the contoured shape holds the heel comfortable in place. Easily placed, with the sticky side to the heel counter, the Invisible Heel Guard offers a lasting comfortable wear.



Code	Unit
AH-IHG	PAIR

**PHYSIFEET INVISIBLE
HEEL CUSHION**

The PhysiFeet Invisible Heel Cushion provides relief to tired and aching heels, especially for long hours of standing (and dancing) in fashion footwear. The thin clear gel cushion offers a non slip textured base, which fits easily into all shoes. Easily placed in the base of the heel area, the Heel Cushion is recommended for high heel shoes.



Code	Unit
AH-IHC	PAIR

**PHYSIFEET INVISIBLE
FOREFOOT CUSHION**

The PhysiFeet Invisible Forefoot Cushion is for the tired and burning forefoot. When the pressure from high heels on the ball of the foot becomes too much, the thin cushion fits easily into place to offer relief. A non slip textured base prevents friction and shearing and fits easily sticky side down to cushion the forefoot.



Code	Unit
AH-IFC	PAIR

**PHYSIFEET GEL THONG
PROTECTOR**

The PhysiFeet Gel Thong Protector is a perfect accessory for all straps found on thongs and sandals. The soft gel is designed and shaped to fit over the toe strap, to prevent friction and rubbing to the toes. In addition PhysiFeet Gel Thong Protectors will alleviate rubbing, blistering and add comfort to the wearing of thongs and sandals.



Code	Unit
AH-GTP	PAIR