

PRO PLUS 603HS

One of our most popular acoustic scopes, the 603 combines outstanding acoustic performance with rugged durability.



The Pro Plus 603HS features:

- Combination adult size chestpiece precisely machined to exacting tolerances from surgical stainless steel for outstanding performance and rugged durability
- Extra large bell (1³/₈" diameter) for unsurpassed low frequency response
- Ultra-sensitive diaphragm (1⁷/₈" diameter) for greater amplification and crisper high frequency response
- Color coordinated non-chill bell and diaphragm retaining rim for patient comfort
- Stainless steel binaurals with double leaf internal spring fixed at 15° angle for greater comfort
- Reinforcing yoke molded into flexible 22" PVC "Y" tubing
- Deluxe PVC eartips (in small and large) for the ultimate in wearing comfort and acoustic seal
- Spare diaphragm included
- Available in 3 striking colors for departmental coding or personal identification
- Scope ID Tag included
- Weighs 6.5 oz.
- Overall length 31"
- Lifetime Warranty

PRIVATE LABELING	
COMPONENT	PACKAGING
N/A	2 Pc. Box/Label



Feature	Benefit
Exclusive source has been producing premium scopes longer than any foreign mfg	Unsurpassed quality and consistency
CNC machined chestpieces	Manufacturing and acoustic consistency
Stainless steel chestpiece AND binaurals	The ultimate in durability and acoustic performance
Non chill bell and diaphragm rims color coordinated to the tubing	For patient comfort and striking style
Binaurals fixed at 15°	Ergonomically correct position for better acoustic seal.
Dual leaf internal binaural spring	Tested to 10,000 cycles for durability
Reinforcing yoke at tube Y	Reinforces stress point for greater durability
Combination diaphragm and bell chestpiece	Allows attenuating of desired frequencies. Diaphragm for high frequency sounds, bell for lower frequency sounds
Threaded eartips	Stay on securely – for the ultimate in wearing comfort and acoustic seal
100% inspected in the USA	To ensure quality, and performance of EVERY scope
Lifetime warranty	To protect the consumer's investment

ORDERING INFORMATION					
Schein Item#	ADC Item#	Product Description	Latex Content	U/M	Cs Pk
9004816	603BDHS	PRO PLUS 603, Adult, Burgundy	Latex-Free	ea.	10
9004815	603BKHS	PRO PLUS 603, Adult, Black	Latex-Free	ea.	10
9004817	603NHS	PRO PLUS 603, Adult, Navy	Latex-Free	ea.	10