

Tools

> Home > Healthcare > Laboratory Diagnostics > Immunoassay > Immunoassay Systems > IMMULITE 1000 Immunoassay System

IMMULITE 1000 Immunoassay System

Reliability that keeps your lab running

> Overview | **Features & Benefits** > Assays > Technical Specifications

Key features of the IMMULITE® 1000 Immunoassay System:

Dependable operation provides laboratories with confidence

- The reliability allows laboratories confidence in assay turnaround times
- Cost of ownership is low due to minimal service, operator interventions, and operator time

Improving laboratory efficiency through the Ease of Use

- Few consumables
- Daily maintenance in less than 10 minutes
- Sample and reagent level sensors
- Easy entry of assay parameters, expiration dates, and calibration adjustment intervals
- Integration of the computer, keyboard, and monitor provide a consolidated footprint

Flexible Solutions to adapt to a laboratory's needs

- Throughput up to 120 tests per hour
- Turbo assays (Troponin I, CK-MB, Myoglobin, Intact PTH, and HCG) allow results in 15 minutes for critical tests
- 12 position reagent carousel

Related Products, Services & Resources



1000 *Turbo*
Assays



Management
System



Remote Service



Management