

# Magnitude

## ACCESSORY GUIDE

REV A



Model 3150M  
REF 865128  
April 2011









## Frequently Ordered Items

Part Number	Description	Page	QTY
<b>ECG Accessories</b>			
989803152291	Prep Gel	4	
989803152321	Standard Quadrode Electrodes, Box of 50	4	
989803152341	Neonatal Quadrode Electrodes, Box of 10	4	
989803152361	CV Quadrode Electrodes, Box of 10	4	
<b>SpO<sub>2</sub> Accessories</b>			
989803162501	Adjustable Grip Sensor, Adult/Neonatal	5	
989803162511	Adjustable Grip Sensor, Infant/Pediatric	5	
989803166551	SpO <sub>2</sub> Quick Connect, Adult, Box of 20	6	
989803166561	SpO <sub>2</sub> Quick Connect, Pediatric, Box of 20	6	
989803166571	SpO <sub>2</sub> Quick Connect, Infant, Box of 20	6	
989803166581	SpO <sub>2</sub> Quick Connect, Neo, Box of 20	6	
989803167111	SpO <sub>2</sub> Quick Connect, Starter Kit, (5 sensors)	6	
<b>NIBP Accessories</b>			
989803170431	Disposable Infant Cuff, Box of 10	7	
989803170441	Disposable Pediatric Cuff, Box of 10	7	
989803170451	Disposable Small Adult Cuff, Box of 10	7	
989803170461	Disposable Adult Cuff, Box of 10	7	
989803170471	Disposable Adult Extra-Long Cuff, Box of 10	7	
989803170481	Disposable Large Adult Cuff, Box of 10	7	
989803170491	Disposable Large Adult Extra-Long Cuff, Box of 10	7	
989803169421	Neonatal NIBP Hose (Single Lumen)	8	
989803170401	Neonatal Disposable Cuff, #3	8	
989803170421	Neonatal Disposable Cuff, #5	8	
<b>Accessories for EtCO<sub>2</sub>/Agent System</b>			
989803154101	Sample Line Kit	9	
989803152561	Adult Cannula	9	
989803152571	Pediatric Cannula	9	
989803152581	Small Infant Cannula	9	
989803152591	Infant Cannula	9	
989803152601	Divided Adult Cannula	9	
989803152611	Divided Small Infant Cannula	9	
989803152621	Divided Infant Cannula	9	
989803152631	Divided Pediatric Cannula	9	
<b>Temperature Accessories</b>			
989803152801	FiberTEMP Surface Temperature Sensor	11	






## Accessory Guide Contents by Parameter Type









<b>Accessory Type</b>	<b>Page</b>
<b>ECG Accessories</b>	<b>4</b>
<b>SpO<sub>2</sub> Accessories</b>	<b>5</b>
<b>NIBP Accessories</b>	<b>7</b>
<b>Respiration &amp; Gas Accessories</b>	<b>9</b>
<b>Temperature Accessories</b>	<b>11</b>
<b>Power and Batteries</b>	<b>12</b>
<b>Documentation</b>	<b>13</b>
<b>Gating Accessories</b>	<b>13</b>
<b>Miscellaneous</b>	<b>13</b>

## ECG Accessories

Part Number	Description	Image	QTY
989803152291 (9009)	<b>Prep Gel</b> Package of 3 tubes (453564120521). Must be used when applying MRI electrodes		
989803152311 (9303A)	<b>Standard Quadrode (Box of 10)</b> MRI Electrodes compatible with the 9224 Standard Lead Cable		
989803152321 (9303N)	<b>Standard Quadrode (Box of 50)</b> MRI Electrodes compatible with the 9224 Standard Lead Cable		
989803152341 (9372A)	<b>Neonatal Quadrodes (Box of 10)</b> MRI Electrodes compatible with the 9222 Neonatal Lead Cable		
989803152361 (9371A)	<b>CV Quadrode (Box of 10)</b> MRI Electrodes compatible with the CV Lead Cable 9223, Advanced Filter Cable 989803170121 and Special Application Cable 452213305301		
989803154751 (9218)	<b>CV ECG Leadwire Set for 3150</b> For use with 931A Quadtrodes: Cable Set Includes: ECG/Leadwires (9218B) and Cable (9240B)		
989803154761 (9219)	<b>Neo Quadrode ECG Leadwire Set for 3150</b> For use with 9272A Quadtrods Only: Cable Set Contains: ECG/Leadwires (9219A) and Cable (9240B)		
989803154811 (9240)	<b>Standard Quadrode ECG Leadwire Set for 3150</b> For use with 9303A Quadtrodes only. Cable set Includes: ECG/Leadwires (9240A) and Cable (9240B)		











## SpO<sub>2</sub> Accessories







Part Number	Description	Image	List Price
<b>989803155081</b>  (9367)	<b>Reusable Clip Sensor for 3150</b> <i>Adult sized clip sensor for 3150. Not compatible with wireless SpO<sub>2</sub></i>		
<b>989803155201</b>  (9399A)	<b>Grip Sensor Cable for 3150</b> <i>Includes one each of Adult, Pediatric, Infant, and Neonatal Grips (9467) - Is not compatible with wireless SpO<sub>2</sub></i>		
<b>989803162491</b>	<b>Adjustable SpO<sub>2</sub> Starter Pack (Box of 2)</b> Contains one legacy reusable adjustable adult/neonatal grip and one legacy reusable adjustable infant/pediatric grip. Designed for adult or neonatal patients. Compatible with the 9399B Sensor		
<b>989803162501</b>	<b>Adult/Neonatal SpO<sub>2</sub> Grip (Box of 2)</b> Contains two legacy reusable adjustable adult/neonatal grips. Designed for adult or neonatal patients. Compatible with the 9399B Sensor		
<b>989803162511</b>	<b>Infant/Pediatric SpO<sub>2</sub> Grip (Box of 2)</b> Contains two legacy reusable adjustable infant/pediatric grips. Designed for infant or pediatric patients. Compatible with the 9399B Sensor		

Part Number	Description	Image	List Price
989803174181	<b>Quick Connect Adapter for 9399A</b> <i>Adapts the 9399A cable to accept the Quick Connect SpO<sub>2</sub> accessories</i>		
989803166551	<b>Adult Quick Connect Disposable SpO<sub>2</sub> Applicator (Box of 20)</b> <i>For use with 989803174181 Quick Connect SpO<sub>2</sub> Adapter Only</i>		
989803166561	<b>Pediatric Quick Connect Disposable SpO<sub>2</sub> Applicator (Box of 20)</b> <i>For use with 989803174181 Quick Connect SpO<sub>2</sub> Adapter Only</i>		
989803166571	<b>Infant Quick Connect Disposable SpO<sub>2</sub> Applicator (Box of 20)</b> <i>For use with 989803174181 Quick Connect SpO<sub>2</sub> Adapter Only</i>		
989803166581	<b>Neonatal Quick Connect Disposable SpO<sub>2</sub> Applicator (Box of 20)</b> <i>For use with 989803174181 Quick Connect SpO<sub>2</sub> Adapter Only</i>		
989803167111	<b>Quick Connect SpO<sub>2</sub> Starter Kit</b> <i>For use with 989803174181 Quick Connect SpO<sub>2</sub> Adapter Only. Includes 5 each of Quick Connect Disposable Grips (Adult, Pediatric, Infant, and Neonatal)</i>		
989803166531	<b>Adult Quick Connect Reusable SpO<sub>2</sub> Clip</b> <i>For use with 989803174181 Quick Connect SpO<sub>2</sub> Adapter Only</i>		
989803166541	<b>Pediatric Quick Connect SpO<sub>2</sub> Clip</b> <i>For use with 989803174181 Quick Connect SpO<sub>2</sub> Adapter Only</i>		














## NIBP Accessories







Part Number	Description	Image	QTY
989803152451 (94014)	<b>“Quick Connect” NIBP Starter Kit</b> <i>Includes: MRI NIBP Adult Hose (989803169411), Adult Cuff (989803169451) and Pediatric Cuff (989803169441)</i>		
989803169411	<b>“Quick Connect” NIBP Hose (Twin Lumen) 18’ MRI Hose</b> <i>For use with Quick Connect NIBP Cuffs</i>		
989803170431	<b>“Quick Connect” Infant NIBP Disposable Cuff (Box of 10)</b> <i>Circumference range 10 to 15 cm compatible with 989803169411 MRI NIBP Hose</i>		
989803170441	<b>“Quick Connect” Pediatric NIBP Disposable Cuff (Box of 10)</b> <i>Circumference range 14 to 21 cm compatible with 989803169411 MRI NIBP Hose</i>		
989803170451	<b>“Quick Connect” Small Adult NIBP Disposable Cuff (Bx of 10)</b> <i>Circumference range 20 to 28 cm compatible with 989803169411 MRI NIBP Hose</i>		
989803170461	<b>“Quick Connect” Adult NIBP Disposable Cuff (Box of 10)</b> <i>Circumference range 27 to 36 cm compatible with 989803169411 MRI NIBP Hose</i>		
989803170471	<b>“Quick Connect” Adult Extra Long NIBP Disposable Cuff (Box of 10)</b> <i>Circumference range 27 to 36 cm compatible with 989803169411 MRI NIBP Hose</i>		
989803170481	<b>“Quick Connect” Large Adult NIBP Disposable Cuff (Bx of 10)</b> <i>Circumference range 35 to 45 cm compatible with 989803169411 MRI NIBP Hose</i>		
989803170491	<b>“Quick Connect” Large Adult Extra-Long NIBP Disposable Cuff (Box of 10)</b> <i>Circumference range 35 to 45 cm compatible with 989803169411 MRI NIBP Hose</i>		
989803169431	<b>“Quick Connect” Infant NIBP Cuff (Twin Lumen)</b> <i>Circumference range 5 to 15 cm compatible with 989803169411 MRI NIBP Hose</i>		

989803169441	<p><b>“Quick Connect” Pediatric NIBP Cuff (Twin Lumen)</b> <i>Circumference range 9 to 25 cm compatible with 989803169411 MRI NIBP Hose</i></p>		
989803169451	<p><b>“Quick Connect” Adult NIBP Cuff (Twin Lumen)</b> <i>Circumference range 17 to 41 cm compatible with 989803169411 MRI NIBP Hose</i></p>		
989803169461	<p><b>“Quick Connect” Large Adult NIBP Cuff (Twin Lumen)</b> <i>Circumference range 39 to 62 cm compatible with 989803169411 MRI NIBP Hose</i></p>		
989803169421	<p><b>“Quick Connect” Neonatal NIBP Hose (Single Lumen)</b> <i>15' MRI Neonatal Hose for use with 989803170401 and 989803170421 NIBP Cuffs</i></p>		
989803170401	<p><b>“Quick Connect” Disposable Neonatal NIBP Cuff #3 (Single Lumen) - QTY 1 -</b> <i>Circumference range 6 to 11 cm compatible with 989803169421 MRI Neonatal NIBP Hose</i></p>		
989803170421	<p><b>“Quick Connect” Disposable Neonatal NIBP Cuff #5 (Single Lumen) - QTY 1 -</b> <i>Circumference range 8 to 15 cm compatible with 989803169421 MRI Neonatal NIBP Hose</i></p>		





## Respiration & Gas Accessories





Part Number	Description	Image	QTY
989803154101 (9010D)	<b>ETCO2 SAMPLING KIT,20 FT.</b> <i>EtCo2 sample line kit for 3150. Includes sample line and nafion tubing.</i>		
989803152561 (9012)	<b>Adult Disposable Cannula</b>		
989803152571 (9013)	<b>Pediatric Disposable Cannula</b>		
989803152581 (9014)	<b>Small Infant Disposable Cannula</b>		
989803152591 (9015)	<b>Infant Disposable Cannula</b>		
989803152601 (9016)	<b>Divided Adult Disposable Cannula</b> <i>Disposable Adult Divided Cannula for simultaneous delivery of O<sub>2</sub> and EtCO<sub>2</sub> sampling</i>		
989803152611 (9016A)	<b>Divided Small Infant Disposable Cannula</b> <i>For simultaneous delivery of O<sub>2</sub> and EtCO<sub>2</sub> sampling</i>		
989803152621 (9016B)	<b>Divided Infant Disposable Cannula</b> <i>For simultaneous delivery of O<sub>2</sub> and EtCO<sub>2</sub> sampling</i>		
989803152631 (9016C)	<b>Divided Pediatric Disposable Cannula</b> <i>For simultaneous delivery of O<sub>2</sub> and EtCO<sub>2</sub> sampling</i>		
989803152681 (9471)	<b>Scavenger Hose</b>		
989803152691 (9025)	<b>Endotracheal Airway Adapter (Box of 50)</b>		

<p><b>989803152861</b> (HP176)</p>	<p><b>Zero Port Filter</b> <i>Replacement filter</i></p>		
<p><b>989803153461</b> (9445)</p>	<p><b>O2 Cell</b> <i>Replacement O2 Cell for Precess &amp; 3155A systems</i></p>		
<p><b>989803152781</b> (9461)</p>	<p><b>O2 Sensor Kit - Disposable</b> <i>Replacement O2 Cell and Assembly for 3155A</i></p>		
<p><b>989803154121</b> (9010H)</p>	<p><b>NAFION DRYER TUBE</b> <i>Dryer tube used with 3150 and 3155A systems</i></p>		
<p><b>989803154281</b> (9026)</p>	<p><b>EtCO2 Disc Filter (50 PK)</b> <i>protective hydrophobic disc filter used with 3150 and 3155A systems</i></p>		
<p><b>989803155491</b> (9456)</p>	<p><b>Disposable Water Trap</b> <i>Includes water trap, connection tubing, nafion tube dryer and 0.08 micron hydrophobic disk filter</i></p>		

## Temperature Accessories

Part Number	Description	Image	QTY
<b>989803152801</b>  (9320)	<b>FiberTEMP Sensor</b> <i>For surface temperature only</i>		
<b>989803152811</b>  (9321)	<b>FiberTemp Sensor Applicators (Box of 30)</b> <i>Only for surface application</i>		





## Power Accessories

Part Number	Description	Image	QTY
<b>989803152221</b>  (AC517A)	<b>5-Foot MRI Power Cable</b> Connects MR Fusion Power Supply (989803171721) or White Power Supply (989803169201) to Monitor Base. For use with Magnitude system		
<b>989803152231</b>  (AC517A)	<b>25-Foot MRI Power Cable</b> Connects MR Fusion Power Supply (989803171721) or White Power Supply (989803169201) to Monitor Base. For use with Magnitude system		
<b>989803169201</b>	<b>Traditional White MRI Power Supply</b> Must be used with a line cord and AC517 cables		
<b>989803168211</b>	<b>North American Line Cord</b> Hospital grade 125v cored with US plug type for use with all Invivo monitor power supplies.		



## Documentation

Part Number	Description	Image	QTY
989803155611 (9546CD)	3150 Service Manual	Image currently unavailable	
989803155721 (9572)	3150 and 3155 Operators Manual	Image currently unavailable	

## Gating Options/Accessories

Part Number	Description	Image	QTY
989803154871 (9270)	<b>Siemens Gating Cable for 3150</b> Interfaces Magnitude with Siemens Espree/Avanto/Trio/Symphony/Sonata/Harmony/Allegra MRI systems. - Provides cardiac and peripheral gating.		
989803154881 (9278)	<b>GE Excite Gating Interface</b> Connects Magnitude cardiac and peripheral gating output to GE Excite, MR 450, MR750 and Horizon MRI systems		
989803154871 (9277)	<b>Hitachi/Toshiba Gating Interface</b> Connects Magnitude cardiac and peripheral gating output to Hitachi and Toshiba MRI systems		
989803162671 (9297)	<b>Vector Gating Interface</b> Connects Magnitude cardiac and peripheral gating option to MRI systems equipped with vector gating such as the Philips Achieva MRI systems		
989803162831 (AC429)	<b>Philips Intera Gating Interface</b> Connects Magnitude cardiac and peripheral gating output to Philips Intera MRI systems	Image currently unavailable	

## Miscellaneous Accessories

Part Number	Description	Image	QTY
989803154301 (9029)	POLE,IV,NON-MAGNETIC	Image currently unavailable	
989803155521 (9456)	KIT,BASKET,POLE MOUNT,MRI	Image currently unavailable	
989803156161 (HM66T)	HANDLEBAR,3150/M	Image currently unavailable	
989803152871 (9059)	<b>2.4GHz Passive Antenna Kit</b> <i>Cable lengths (AC485) will need to be known before ordering. Maximum distance is 100 ft Additional installation charge may be required</i>		
989803152911 (9180T)	<b>Printer Paper (Box of 10)</b>		



inspired by a better way

For more information about **Invivo MRI Patient Monitoring Systems** or any of the Complete Solution products from Invivo, please visit our website at:

[www.invivocorp.com](http://www.invivocorp.com) or  
[www.invivoexpression.com](http://www.invivoexpression.com)

#### INVIVO CORPORATION

ORLANDO, FL USA

T: 800-331-3220

E: [INFO@INVIVOCORP.COM](mailto:INFO@INVIVOCORP.COM)

EUROPE

T: +31 (0)499 378624

E: [CONTACT@INVIVOCORP.COM](mailto:CONTACT@INVIVOCORP.COM)

[WWW.INVIVOCORP.COM](http://WWW.INVIVOCORP.COM)

[WWW.INVIVOEXPRESSION.COM](http://WWW.INVIVOEXPRESSION.COM)

# PHILIPS

Invivo is a Philips Company

Products, services and technical specifications are subject to change without notice.

© 2011 Invivo Corporation All Rights Reserved.