

WALKERS

2014 Product Catalog

Drive Medical's broad line of walkers exceeds the market demand for quality, selection, reliability & value.



for the road ahead™
drive



About Drive

The Company markets and distributes its products to customers located throughout North America, Europe, Central America, South America, the Middle East and Asia.

Drive Medical Design and Manufacturing is one of the fastest growing global manufacturers and distributors of durable medical equipment. The Company markets and distributes its products to customers located throughout North America, Europe, Central America, South America, the Middle East and Asia.

Philosophy/Mission

Our mission is to manufacture and market world class products which focus on innovation, functionality and value to improve the quality of life and to promote independence for the individuals that use our products. The Company is dedicated to its customers by consistently offering the best quality and most functional products across a broad spectrum of health care providers' needs; better value and aggressive pricing; world-class customer service and technical support; friendly people to work with; and ongoing innovation.

Facilities

Drive has corporate offices and distribution facilities located throughout the United States, Canada, the United Kingdom, Germany, Romania, China and Taiwan. The Company's North America Technical Service Center is located in Bohemia, New York.

Office Locations

Corporate Headquarters

99 Seaview Boulevard
Port Washington, NY 11050
t: 877.224.0946 | 516.998.4600

CHAD® / Inovo

401 Leonard Blvd N
Lehigh Acres, FL 33971

Distribution Centers

Southeast and Central US

Distribution Center

4131 Southmeadow Parkway West
Suite 100
Atlanta, GA 30349

West Coast Distribution Center

3550-B Jurupa Street
Ontario, CA 91761

North East Distribution Center

200 Docks Corner Road
Suite 241
South Brunswick, NJ 08810

International Distribution Center

Drive Medical Canada Inc.
50 Vogell Road, Unit 7
Richmond Hill, Ontario L4B 3K6

Technical Service Center

1101 Lakeland Avenue
Suite 100
Bohemia, NY 11716
t: 866.694.5085

For international operations visit
www.drivemedical.com

Icon Library



New
Product



1-Pack
Available



Allergy
Free



Environment
Friendly



Fashion
Colors



Latex
Free



Retail
Packaged





Tool-free
Assembly



Weight
Capacity



Bariatric Weight
Capacity



LiftWalker™ with Retractable Stand Assist Bars

| Item # | Product Description | UOM |
|---------|--|------|
| 10277LW | LiftWalker™, Two Button, Adult, with Retractable Assist Bars | 1/cs |

- LiftWalker™ will provide convenient, safe, assisted independence anywhere, anytime it's needed.
- For use by individuals who need assistance with sitting/standing on their own.
- Safer than being pulled by an attendant or caregiver.
- Reduces risk of injury to an attendant or caregiver.
- The retractable stand assist bars are easy to deploy. (Patented)
- The LiftWalker™ makes it easy to stand because the individual uses biceps, chest, back, stomach and leg muscles all at once. (Figure A - C)
- Easy to lock/unlock foot brake.
- Easy push-button mechanism may be operated by fingers, palms, or side of hand.
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing.
- Comfortable, soft foam-padded hand grips.
- Comes with 5" wheels.
- New retail packaging.

specifications

| | |
|------------------------|---|
| Width | 17" inside hand grip 21" inside back legs |
| Overall Width | 26" |
| Depth (Opened at Base) | 21" |
| Height | 32.5" - 38.5" |
| Weight | 11 lbs. |
| Carton Shipping Weight | 13 lbs. |
| Weight Limit | 300 lbs. |
| Warranty | Limited Lifetime |



A



A

Deluxe Folding Walker, Two Button

| Item # | Product Description | Product Height | UOM |
|---------|---------------------|----------------|------|
| 10200-4 | Adult | 32" - 39" | 4/cs |
| 10201-4 | Junior | 25" - 32" | 4/cs |

- Easy push-button mechanisms may be operated by fingers, palms or side of hand.
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing.
- Unique U-shaped frame design creates greater clearance.
- The patient is able to bring the walker closer to assist in standing because of the improved crossframe. (Figure A)
- Sturdy 1" diameter, aluminum construction ensures maximum strength while remaining lightweight.

specifications

| | |
|--|---|
| Width | 17" inside hand grip 19" inside back legs |
| Overall Width | 22" |
| Depth (Opened at Base) | 15.5" |
| Weight | 6 lbs. |
| Carton Shipping Weight (10200 10201) | 26 24 lbs. |
| Weight Limit | 350 lbs. |
| Warranty | Limited Lifetime |

Universal (Adult/Junior) Deluxe Folding Walker, Two Button

| Item # | Product Description | UOM |
|----------|--------------------------|------|
| 10248N-1 | Universal (Adult/Junior) | 1/cs |
| 10248N-4 | Universal (Adult/Junior) | 4/cs |

- Walker height can be adjusted from 28.25" - 38.25".
- Two walkers in one.
- Side brace adjusts with height of walker to provide additional stability.
- New easy-release closure. (Figure A)
- Easy push-button mechanisms may be operated by fingers, palms or side of hand. Patent Pending.
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing.
- Sturdy 1" diameter aluminum construction ensures maximum strength while remaining lightweight.

specifications

| | |
|--|---|
| Width | 17" inside hand grip 19" inside back legs |
| Overall Width | 21.75" |
| Depth (Opened at Base) | 15.75" |
| Height | 28.25" - 38.25" |
| Weight | 6 lbs. |
| Weight Limit | 300 lbs. |
| Carton Shipping Weight (10248N-1 10248N-4) | 8 28.5 lbs. |



A

B



A

Universal (Adult/Junior) Folding Walker, Two Button with 5" Wheels

| Item # | Product Description | UOM |
|---------|--------------------------|------|
| 10253-1 | Universal (Adult/Junior) | 1/cs |
| 10253-4 | Universal (Adult/Junior) | 4/cs |

- Two walkers in one. Height adjusts from 28.75" - 38.5".
- Side brace adjusts with height of walker to provide additional stability.
- Easy push-button mechanisms may be operated by fingers, palms or side of hand. Patent Pending.
- New easy-release closure. (Figure A)
- Each side operates independently to allow easy movement through narrow spaces and greater stability.
- Sturdy 1" diameter, aluminum construction ensures maximum strength while remaining lightweight.
- Newly designed rear glide cap allows walker to glide easily and smoothly over most surfaces. (Figure B)

specifications

| | |
|--|--|
| Width | 17" inside hand grip 18.75" inside back legs |
| Overall Width | 24" |
| Depth (Opened at Base) | 17.75" |
| Height | 28.75" - 38.5" |
| Weight | 7 lbs. |
| Carton Shipping Weight (10253-1 10253-4) | 8 30.5 lbs. |
| Weight Limit | 300 lbs. |
| Warranty | Limited Lifetime |


Deluxe Folding Walker, Two Button with 3" Wheels

| Item # | Product Description | UOM |
|---------|---------------------|------|
| 10206-4 | Adult | 4/cs |

- Easy push-button mechanisms may be operated by fingers, palms or side of hand.
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing.
- Sturdy 1" diameter, aluminum construction ensures maximum strength while remaining lightweight.
- Newly designed rear glide cap allows for walker to glide easily and smoothly over most surfaces. (Figure A)
- Vinyl-contoured hand grip.
- Comes with wheels. (Replacement wheels item # 10109).

specifications

| | |
|------------------------|---|
| Width | 17" inside hand grip 20" inside back legs |
| Overall Width | 24" |
| Depth (Opened at Base) | 17.5" |
| Height | 32" - 39" |
| Weight | 7 lbs. |
| Carton Shipping Weight | 31.5 lbs. |
| Weight Limit | 350 lbs. |
| Warranty | Limited Lifetime |



Deluxe Folding Walker, Two Button with 5" Wheels

| Item # | Product Description | Product Height | UOM |
|---------|---------------------|----------------|------|
| 10210-4 | Adult | 32" - 39" | 4/cs |
| 10211-4 | Junior | 25" - 32.25" | 4/cs |

- Easy push-button mechanisms may be operated by fingers, palms or side of hand.
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing.
- Sturdy 1" diameter aluminum construction ensures maximum strength while remaining lightweight.
- Newly designed rear glide cap allows for walker to slide easily and smoothly over most surfaces. (Figure A)
- Vinyl-contoured hand grip.
- Comes with wheels. (Replacement wheels item # 10109).

specifications

| | |
|--|---|
| Width | 17" inside hand grip 20" inside back legs |
| Overall Width | 24" |
| Depth (Opened at Base) (10210 10211) | 17.5" 16.25" |
| Weight (10210 10211) | 7.5 7 lbs. |
| Carton Shipping Weight (10210 10211) | 30 29 lbs. |
| Weight Limit | 350 lbs. |
| Warranty | Limited Lifetime |

Bariatric Aluminum Folding Walker, Two Button

(Wider & Deeper Frame Design)

| Item # | Product Description | Product Height | UOM |
|---------|---------------------|----------------|------|
| 10220-1 | Adult | 32" - 39" | 1/cs |
| 10220-2 | Adult | 32" - 39" | 2/cs |
| 10223-2 | Junior | 25" - 32" | 2/cs |

- For Bariatric use.
- Steel legs and side braces.
- Precision design provides additional strength while adding minimal weight.
- Wider and deeper frame design accommodates individuals up to 500 lbs.
- Easy push-button mechanisms may be operated by fingers, palms or side of hand.
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing.
- Vinyl-contoured hand grip.

specifications

| | |
|--|---|
| Width | 20.5" inside hand grip 23" inside back legs |
| Overall Width (10220 10223) | 25.25" 24.5" |
| Depth (Opened at Base) | 20" |
| Weight (10220 10220) | 9 8 lbs. |
| Carton Shipping Weight (10220-1 10220-2 10223-1) | 12 21 18 lbs. |
| Weight Limit | 500 lbs. |
| Warranty | Limited Lifetime |



drive™ GUARDIAN®
10220HD G30754B

Bariatric Aluminum Folding Walker, Two Button

(Wider & Deeper Frame Design) Heavy Duty Wheels

| Item # | Product Description | UOM |
|-----------|---------------------|------|
| 10220-2WW | Adult | 2/cs |

- For Bariatric use.
- Steel legs and side braces.
- Precision design provides additional strength while adding minimal weight.
- Wider and deeper frame design accommodates individuals up to 500 lbs.
- Easy push-button mechanisms may be operated by fingers, palms or side of hand.
- Each side operates independently to allow easy movement through narrow spaces and greater stability.
- Vinyl-contoured hand grip.

specifications

| | |
|------------------------|---|
| Width | 20.5" inside hand grip 23" inside back legs |
| Overall Width | 27" |
| Depth (Opened at Base) | 23" |
| Height | 32" - 39" |
| Weight | 11 lbs. |
| Carton Shipping Weight | 24 lbs. |
| Weight Limit | 500 lbs. |
| Warranty | Limited Lifetime |

Bariatric, Extra Wide and Deeper Frame, Folding Walker, Two Button

| Item # | Product Description | UOM |
|---------|---------------------|------|
| 10220HD | Heavy-Duty, Adult | 2/cs |

- For Bariatric use.
- Precision design provides additional strength while adding minimal weight.
- Wider and deeper frame design accommodates individuals up to 700 lbs.
- Easy push-button mechanisms may be operated by fingers, palms or side of hand.
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing.
- Heavy-duty vinyl-contoured hand grip.

specifications

| | |
|------------------------|---|
| Width | 19.5" inside hand grip 21.5" inside back legs |
| Overall Width | 24.5" |
| Depth (Opened at Base) | 21.75" |
| Height | 32.5" - 39.5" |
| Weight | 9 lbs. |
| Carton Shipping Weight | 20 lbs. |
| Weight Limit | 700 lbs. |
| Warranty | Limited Lifetime |



**Lightest
Walker on
the Market!**

Push Palm Release Folding Walker

| Item # | Product Description | Product Height | UOM |
|--------|---------------------|----------------|------|
| 10291 | Adult | 30" - 37" | 4/cs |
| 10292 | Junior | 26" - 29" | 4/cs |


- The unique Push Palm Release makes it easy to open and close the walker from a seated or standing position.
- (Figure A)
- The release mechanism pops up with an audible “click” that ensures the side frame is completely opened.
- Weighs less than 5 lbs.
- Manufactured utilizing a unique process that eliminates the use of fasteners to prevent any component from becoming loose over time. (Figure B)
- Sturdy 1" diameter aluminum is lightweight and ensures maximum strength. (Figure C)
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing.
- Deluxe soft-ribbed contoured hand grips.
- Patent Pending.

specifications

| | |
|--|------------------|
| Width inside hand grip (10291 10292) | 17.5" 17.25" |
| Width inside back legs | 22" 20.5" |
| Depth (Opened at Base) (10291 10292) | 17.5" 16.25" |
| Weight (10291 10292) | 4.8 4.3 lbs. |
| Carton Shipping Weight (10291 10292) | 23 21 lbs. |
| Weight Limit | 300 lbs. |
| Warranty | Limited Lifetime |

accessories

| | | |
|-------|-----------------------|------|
| 10126 | Walker Wheels, Adult | 1/pr |
| 10127 | Walker Wheels, Junior | 1/pr |



**Lightest
Walker on
the Market!**

Push Palm Release Folding Walker with 5" Wheels

| Item # | Product Description | Product Height | UOM |
|--------|---------------------|----------------|------|
| 10293 | Adult | 33" - 39" | 4/cs |
| 10294 | Junior | 28" - 30" | 4/cs |

- The unique Push Palm Release makes it easy to open and close the walker from a seated or standing position. (Figure A)
- The release mechanism pops up with an audible "click" that ensures the side frame is completely opened.
- Weighs only 4.9 lbs.
- Manufactured utilizing a unique process that eliminates the use of fasteners to prevent any component from becoming loose over time. (Figure B)
- Sturdy 1" diameter aluminum is lightweight and ensures maximum strength. (Figure C)
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing.
- Deluxe soft-ribbed contoured hand grips.
- Patent Pending.

specifications

| | |
|--|------------------|
| Width inside hand grip (10293 10294) | 17.25" 17.25" |
| Width inside back legs (10293 10294) | 21.75" 20.25" |
| Depth (Opened at Base) (10293 10294) | 18.5" 17.25" |
| Weight (10293/10294) | 5.3 4.9 lbs. |
| Weight Limit | 300 lbs. |
| Carton Shipping Weight (10293 10294) | 26 24 lbs. |
| Warranty | Limited Lifetime |

accessories

| | | |
|-------|-----------------------|------|
| 10126 | Walker Wheels, Adult | 1/pr |
| 10127 | Walker Wheels, Junior | 1/pr |



Flame Blue



A



Flame Blue



A

Deluxe, Trigger Release Folding Walker

| Item # | Product Description | Color | Product Height | UOM |
|---------|---------------------|-----------------------|----------------|------|
| 10224-4 | Adult | Flame Blue Crossbrace | 32" - 39" | 4/cs |
| 10225-4 | Junior | Flame Blue Crossbrace | 25.5" - 32" | 4/cs |

- Trigger release feature for opening and closing allows users to close the walker without removing their hands from the handle and is ideal for individuals with limited finger dexterity. (Figure A) Trigger release can be reached with hands of all sizes.
- Durable, Composite Trigger Release resists cracking.
- Sturdy 1" diameter anodized, extruded aluminum construction ensures maximum strength while remaining lightweight.
- Patented.

specifications

| | |
|--|--|
| Width | 18" inside hand grip 20.75" inside back legs |
| Overall Width | 23.5" |
| Depth (Opened at Base) (10224 10225) | 16.75" 16.25" |
| Weight (10224 10225) | 6.5 lbs. 6 lbs. |
| Carton Shipping Weight (10224 10225) | 27 26 lbs. |
| Weight Limit | 350 lbs. |
| Warranty | Limited Lifetime |

Deluxe, Trigger Release Folding Walker with 5" Wheels

| Item # | Product Description | Color | Product Height | UOM |
|---------|---------------------|-----------------------|----------------|------|
| 10226-4 | Adult | Flame Blue Crossbrace | 32"-39" | 4/cs |
| 10227-4 | Junior | Flame Blue Crossbrace | 26"-33.5" | 4/cs |

- Trigger release feature for opening and closing allows users to close the walker without removing their hands from the handle and is ideal for individuals with limited finger dexterity. Trigger release can be reached with hands of all sizes. Patented.
- Composite Trigger Release resists cracking.
- Sturdy 1" diameter anodized, extruded aluminum construction ensures maximum strength. Newly designed rear glide cap allows for walker to glide easily and smoothly over most surfaces. (Figure A)

specifications

| | |
|--|--|
| Width | 18" inside hand grip 20.75" inside back legs |
| Overall Width | 25" |
| Depth (Opened at Base) (10226 10227) | 16.75" 16.25" |
| Weight (10226 10227) | 7/6 lbs. |
| Carton Shipping Weight (10226 10227) | 29 27 lbs. |
| Weight Limit | 350 lbs. |
| Warranty | Limited Lifetime |



COMPETITIVE EDGE Line

Deluxe, Trigger Release Folding Walker with 3" Wheels

| Item # | Product Description | UOM |
|--|---------------------|------|
| 10228G | Adult | 4/cs |
| <ul style="list-style-type: none"> • Trigger release feature for opening and closing allows users to close the walker without removing their hands from the handle and is ideal for individuals with limited finger dexterity. (Figure A) Trigger release can be reached with hands of all sizes. • Durable, Composite Trigger Release resists cracking and breaking. • Sturdy 1" diameter anodized, extruded aluminum construction ensures maximum strength while remaining lightweight. • Comes with 3" wheels (Replacement item # 10111) • Patented. | | |

specifications

| | |
|------------------------|--|
| Width | 18" inside hand grip 20.75" inside back legs |
| Overall Width | 24.5" |
| Depth (Opened at Base) | 16.75" |
| Height | 32" - 39" |
| Weight | 7 lbs. |
| Carton Shipping Weight | 32 lbs. |
| Weight Limit | 350 lbs. |
| Warranty | Limited Lifetime |

Deluxe Folding Walker, One-Button

| Item # | Product Description | Product Height | UOM |
|---------|---------------------|----------------|------|
| 10202-2 | Adult | 32" - 39" | 2/cs |
| 10203-2 | Junior | 25" - 32" | 2/cs |

- Single-button release allows quick, one-hand folding.
- Cross brace is located near hand grips to allow for full stride.
- Sturdy 1" diameter anodized, extruded aluminum construction ensures maximum.
- Vinyl-contoured hand grip.

specifications

| | |
|--------------------------------------|---|
| Width | 17" inside hand grip 20.5" inside back legs |
| Overall Width | 24.5" |
| Depth (Opened at Base) (10202/10203) | 18.5 17" |
| Weight (10202/10203) | 6 5 lbs. |
| Carton Shipping Weight (10202/10203) | 13 12 lbs. |
| Weight Limit | 300 lbs. |
| Warranty | Limited Lifetime |



Deluxe Folding Walker, One Button with 5" Wheels

| Item # | Product Description | UOM |
|---------|---------------------|------|
| 10212-2 | Adult | 2/cs |

- Single-button release allows quick, one-hand folding.
- Cross brace is located near hand grips to allow for full stride.
- Sturdy 1" diameter anodized, extruded aluminum construction ensures maximum strength.
- Newly designed rear glide cap allows for walker to glide easily and smoothly over most surfaces. (Figure A)
- Vinyl-contoured hand grip.
- Comes with wheels. (Replacement item #10109)

specifications

| | |
|------------------------|---|
| Width | 17" inside hand grip 20.5" inside back legs |
| Overall Width | 25.5" |
| Depth (Opened at Base) | 18.5" |
| Height | 32" - 39" |
| Weight | 7 lbs. |
| Carton Shipping Weight | 15 lbs. |
| Weight Limit | 300 lbs. |
| Warranty | Limited Lifetime |





10248NBL-1



10248NP-1



Flame Blue



Flame Red



Universal (Adult/Junior) Deluxe Folding Walker, Two Button Black Grips and Tips

| Item # | Product Description | Color | UOM |
|------------|---------------------|-------------------------------------|------|
| 10248NBL-1 | Universal, Adult | Flame Blue Crossbar with Blue Frame | 1/cs |
| 10248NP-1 | Universal, Adult | Flame Red Crossbar with Pink Frame | 1/cs |

- Attractive Blue fashion or Pink fashion color. (Figure A)
- Two walkers in one. Height adjusts from 28.25" - 38.25".
- Side brace adjusts with height of walker to provide additional stability.
- New easy-release closure.
- Easy push-button mechanisms may be operated by fingers, palms or side of hand.
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing.
- Sturdy 1" diameter aluminum construction ensures maximum strength while remaining lightweight.
- Vinyl-contoured hand grip.
- Patent Pending.

specifications

| | |
|------------------------|---|
| Width | 17" inside hand grip 19" inside back legs |
| Overall Width | 21.75" |
| Depth (Opened at Base) | 15.75" |
| Height | 28.25" - 38.25" |
| Weight | 6 lbs. |
| Carton Shipping Weight | 8 lbs. |
| Weight Limit | 300 lbs. |
| Warranty | Limited Lifetime |



Flame Red





Deluxe Folding Travel Walker with 5" Wheels and Fold-Up Legs

| Item # | Product Description | Color | UOM |
|-------------|---------------------|-----------------------------------|------|
| RTL10263KDR | Adult | Flame Red Crossbar with Red Frame | 1/cs |

- Easily folds tool-free to $\frac{1}{2}$ the size of a traditional walker and can be stored in a "Ready Set Go" pouch provided with each walker. (Figure A - C)
- Trigger release feature allows the user to close the walker without removing their hands from the handle and is ideal for individuals with limited finger dexterity. Trigger release can be reached with hands of all sizes.
- Legs are attached with strong elastic cords so legs can be folded. (Figure D - E)
- Comes with 5" front wheels and rear glide skis that are ideal for most surfaces.
- Two side-storage pouches with large pockets for personal items are held in place by Velcro® tabs to both side frames.
- Sturdy 1" diameter anodized, extruded aluminum construction ensures maximum strength while remaining lightweight.
- Deluxe soft-ribbed contoured hand grips.
- Retail Packaged.

specifications

| | |
|------------------------|---|
| Width | 17" inside hand grip 22" inside back legs |
| Overall Width | 26.5" |
| Depth (Opened at Base) | 19.75" |
| Height | 29.5 - 37" |
| Weight | 8 lbs. |
| Carton Shipping Weight | 10 lbs. |
| Weight Limit | 300 lbs. |
| Warranty | Limited Lifetime |



Economy Knee Walker

| Item # | Product Description | UOM |
|--------|---------------------|------|
| 780 | Adult | 1/cs |

- Ideal for individuals recovering from foot surgery, breaks, sprains, amputation and ulcers of the foot.
- Deluxe lever brake ensures safety. Can be adjusted for either right, or left-handed use. (Figure A)
- Tool-free leg and handle-height adjustments. (Figure B - C)
- 8" casters are ideal for indoor/outdoor use.
- 300 lb. Weight Capacity.

specifications

| | |
|------------------------|--|
| Handle Height | 32.5" - 41" |
| Unit Dimensions | 14" (W) x 32" (D) |
| Knee Pad Height | 18" - 22" |
| Knee Pad Dimensions | 7.5" (W) x 14.25" (D) |
| Weight | 17 lbs. |
| Carton Shipping Weight | 21 lbs. |
| Weight Limit | 300 lbs. |
| Warranty | Limited Lifetime, Lifetime Warranty on Brake Cable |

accessories

| | | |
|------------|-----------------------------|------|
| 780 BASKET | Basket for 780 Knee Walkers | 1/ea |
|------------|-----------------------------|------|



DV8 Steerable Aluminum Knee Walker

| Item # | Product Description | UOM |
|--------|---------------------|------|
| 790 | Adult | 1/cs |

- Knee walker can be steered for increased maneuverability.
- Ideal for individuals recovering from foot surgery, breaks, sprains, amputation and ulcers of the foot.
- Deluxe dual-braking system. Sturdy dual frame adds stability and durability. Leg pad with "channel" provides added stability and comfort.
- Tiller can be easily folded for storage and transportation.
- New foldable steering column allows for easy storage.
- Tool-free height adjustment of seat and handle. (Figure A - B)
- 8" casters are ideal for indoor/outdoor use.
- Silver vein finish (Figure C) is easy to maintain.

specifications

| | |
|------------------------|--|
| Handle Height | 33" - 37.5" |
| Unit Dimensions | 17" (W) x 32" (D) |
| Knee Pad Height | 17.5" - 21.5" |
| Pad Dimensions | 7.25" (W) x 14.25" (D) |
| Weight | 19 lbs. |
| Carton Shipping Weight | 26 lbs. |
| Weight Limit | 300 lbs. |
| Warranty | Limited Lifetime, Lifetime Warranty on Brake Cable |



Clever-Lite Walker, Adult, with 5" Casters

with Seat and Loop Locks

| Item # | Product Description | Color | UOM |
|--------|---------------------|-------|------|
| 10230 | Adult | Blue | 1/cs |

- Front wheels can be set in either "swivel" (Figure A) or "fixed" (Figure B) position.
- Special loop lock made of internal aluminum casting operates easily and ensures safety.
- Sturdy 1" diameter anodized aluminum construction.
- Comes with flip-up seat with built-in carry handle.
- Allows individual to be seated or use as traditional walker by raising seat to step inside the frame.
- Patented dual-lever side-paddle release. (Figure C)

specifications

| | |
|------------------------|---|
| Width | 18.5" inside hand grip 19.5" inside back legs Overall 26" |
| Depth | 23" Opened at Base |
| Seat Dimensions | 18.75" (W) x 11" (D) |
| Seat Height | 20.5" - 24.5" |
| Handle Height | 32.5" - 36.5" |
| Weight | 14 lbs. |
| Carton Shipping Weight | 17 lbs. |
| Weight Limit | 300 lbs. |
| Warranty | Limited Lifetime |

accessories

| | | |
|--------|--------|------|
| 10230B | Basket | 1/ea |
|--------|--------|------|

Clever-Lite Walker, Junior, with 5" Casters

with Seat and Loop Locks

| Item # | Product Description | Color | UOM |
|--------|---------------------|-------|------|
| 10230J | Junior | Blue | 1/cs |


- Front wheels can be set in either "swivel" (Figure A) or "fixed" (Figure B) position.
- Special loop lock made of internal aluminum casting operates easily and ensures safety.
- Sturdy 1" diameter anodized aluminum construction.
- Comes with flip-up seat with built-in carry handle.
- Allows individual to be seated or use as traditional walker by raising seat to step inside the frame.
- Patented dual-lever side-paddle release. (Figure C)

specifications

| | |
|------------------------|---|
| Width | 18.5" inside hand grip 19.5" inside back legs 26.5" Overall |
| Depth | 18" opened at base |
| Seat Dimensions | 18.75" (W) x 11" (D) |
| Seat Height | 17.5" - 21.5" |
| Handle Height | 29.5" - 33.5" |
| Weight | 14 lbs. |
| Carton Shipping Weight | 17 lbs. |
| Weight Limit | 300 lbs. |
| Warranty | Limited Lifetime |

accessories

| | | |
|--------|--------|------|
| 10230B | Basket | 1/ea |
|--------|--------|------|



Clever-Lite Walker, Adult, with 8" Casters with Seat and Loop Locks

| Item # | Product Description | Color | UOM |
|--------|---------------------|-------|------|
| 10243 | Adult | Blue | 1/cs |

- Front wheels can be set in either “swivel” (Figure A) or “fixed” (Figure B) position.
- Special loop lock made of internal aluminum casting operates easily and ensures safety.
- Sturdy 1" diameter anodized aluminum construction.
- Comes with flip-up seat with built-in carry handle.
- Allows individual to be seated or use as traditional walker.
- Folds with patented side-paddle release. (Figure C)

specifications

| | |
|------------------------|---|
| Width | 19" inside hand grip 21" inside back legs 27" Overall |
| Depth | 20.5" opened at base |
| Seat Dimensions | 19" (W) x 11" (D) |
| Seat Height | 23" - 27.5" |
| Handle Height | 35" - 39" |
| Weight | 18 lbs. |
| Carton Shipping Weight | 21 lbs. |
| Weight Limit | 300 lbs. |
| Warranty | Limited Lifetime |

accessories

| | | |
|--------|--------|------|
| 10230B | Basket | 1/ea |
|--------|--------|------|

Clever-Lite LS, Adult Walker with Seat and Push-Down Brakes

| Item # | Product Description | UOM |
|--------|---------------------|------|
| 10233 | Adult | 1/cs |

- Rear glide push-down brakes. (Figure A)
- Sturdy 1" diameter anodized, extruded aluminum construction.
- Comes with flip-up seat with built-in carry handle.
- Allows individual to be seated or use as traditional walker by raising seat to step inside the frame.
- Easily folds with patented dual-lever side-paddle release. (Figure B)
- Soft, flexible backrest provides comfort and security while seated.
- 5" fixed front wheels.

specifications

| | |
|------------------------|---|
| Width | 18.5" inside hand grip 21" inside back legs 27" Overall |
| Depth | 21" opened at base |
| Seat Dimensions | 18.75" (W) x 11" (D) |
| Seat Height | 19" - 25.5" |
| Handle Height | 29.5" - 39" |
| Weight | 11 lbs. |
| Carton Shipping Weight | 15 lbs. |
| Weight Limit | 300 lbs. |
| Warranty | Limited Lifetime |



KA 1200N



KA 2200N



KA 3200N



KA 4200N



KA 5200N

drive
Nimbo1P
ACK75
lb85
lb150
lb190
lb250
lbFlame Red
HCPCS
E1399

Nimbo Lightweight Gait Trainer

| Item # | Color | Product Size | Product Width Inside Handlebar | Product Handlebar Height | Product Dimensions (open at base) | Product Weight | Product Weight Capacity | UOM |
|----------|------------------|--------------|--------------------------------|--------------------------|-----------------------------------|----------------|-------------------------|------|
| KA 1200N | Goldenrod Yellow | Tyke | 13.5" | 15.5" - 19.5" | 24" (W) x 20" (D) | 9 lbs. | 75 lbs. | 1/ea |
| KA 2200N | Lime Green | Junior | 13.5" | 19" - 25" | 24" (W) x 27" (D) | 10 lbs. | 85 lbs. | 1/ea |
| KA 3200N | Cornflower Blue | Youth | 15" | 23" - 30.5" | 26" (W) x 31" (D) | 13 lbs. | 150 lbs. | 1/ea |
| KA 4200N | Midnight Blue | Young Adult | 16.5" | 28" - 36" | 26.5" (W) x 34" (D) | 14 lbs. | 190 lbs. | 1/ea |
| KA 5200N | Flame Red | Adult | 17.75" | 36" - 41" | 29" (W) x 37" (D) | 16 lbs. | 250 lbs. | 1/ea |

- Height-adjustable aluminum frame in bright colors.
- Soft rubber wheels adhere to any surface and will not allow the Nimbo to slide backwards.
- One-directional 5" rear wheels with aluminum ratchet and pin mechanism.
- Anti-reverse override bracket disengages the one-directional rear wheels and allows forward and reverse mobility.

- 5" front swivel wheels' locking mechanism easily switches from swivel to non-swivel.
- The KA 2200N - KA 5200N can be used in both anterior and posterior positions.
- Latex-free handgrips with flanged ends.
- Folds easily for convenient transporting.

specifications

| | |
|-----------------------------------|------------------|
| Carton Shipping Weight (KA 1200N) | 11 lbs. |
| Carton Shipping Weight (KA 2200N) | 13 lbs. |
| Carton Shipping Weight (KA 3200N) | 14 lbs. |
| Carton Shipping Weight (KA 4200N) | 18 lbs. |
| Carton Shipping Weight (KA 5200N) | 20 lbs. |
| Warranty | Limited Lifetime |

Nimbo Lightweight Gait Trainers are available with seats. Please see our Wenzelite section of the catalog for details.

Pelvic Stabilizer



Nimbo



HCPCS
E1399

| Item # | For Use With | UOM |
|---------|-------------------------------------|------|
| KA 1050 | All Nimbo Lightweight Gait Trainers | 1/bx |

- The pelvic stabilizer is adjustable in width and features back and lateral pads and a support belt. It is recommended for users who have an asymmetrical posture when standing or walking and, therefore, cannot stay centered in the gait trainer.

Dial-A-Speed Tabs



Nimbo



HCPCS
E1399

| Item # | For Use With | UOM |
|---------|-------------------------------------|---------|
| KA 1400 | All Nimbo Lightweight Gait Trainers | 1 pr/bx |

- Dial-A-Speed tabs provide variable resistance on one or both rear wheels. They allow user to control gait trainer speed and aligns asymmetrical gait patterns.

Forearm Platforms



Nimbo



HCPCS
E1399

| Item # | Product Description | Product Height (from handlebar to armrest) | Product Armrest Length | UOM |
|-------------|---------------------|--|------------------------|---------|
| KA 1035 FPS | Small | 3.5" - 4.5" | 5" | 1 pr/bx |
| KA 1035 FPL | Large | 5" - 7.5" | 8" | 1 pr/bx |

- Constructed of lightweight aluminum.
- Height and depth adjustable.
- Handgrips extend forward or back, supinate, pronate and are angle-adjustable.
- Contoured armrest with a flexible edge.
- Velcro® straps to prevent arms from slipping.
- Choice of mounting positions—inside or outside of the handlebar for lateral adjustment. Spacers are included for additional lateral range.

Bumper Guards



Nimbo



HCPCS
E1399

| Item # | For Use With | UOM |
|---------|-------------------------------------|---------|
| KA 8000 | All Nimbo Lightweight Gait Trainers | 1 pr/bx |

- Bumper guards fit over front casters and protect walls.

Soft Seat Harness



Nimbo

HCPCS
E1399

| Item # | Product Description | UOM |
|----------|---------------------|------|
| CE 1070S | Small | 1/bx |
| CE 1070L | Large | 1/bx |

- The Seat Harness provides cushioned comfort. It moves with user's gait pattern allowing the user to shift weight from one leg to another.
- Straps can be lengthened, shortened and positioned towards the front or back of the gait trainer's frame
- It assists in weight bearing, properly positions user's pelvis and can be used as a safety sling.

Seat



Nimbo

HCPCS
E1399

| Item # | For Use With | Product Seat-to-Floor Height Range | UOM |
|---------|--------------|------------------------------------|------|
| KA 1285 | KA 1200N | 9" - 13" | 1/bx |
| KA 2285 | KA 2200N | 11.5" - 17" | 1/bx |
| KA 3285 | KA 3200N | 12" - 19.5" | 1/bx |
| KA 4285 | KA 4200N | 14.5" - 22" | 1/bx |
| KA 5285 | KA 5200N | 21.5" - 27.5" | 1/bx |

- Solid seats fold up for standing and walking and flip down for convenient seating.

Pouch



Nimbo

HCPCS
E1399

| Item # | Product Size | UOM |
|----------|--------------|------|
| 10258-1 | Small | 1/bx |
| 10258L-1 | Large | 1/bx |



- Nylon carry pouch features two pockets.

9" Front Wheels



Nimbo

HCPCS
E1399

| Item # | For Use With | UOM |
|---------|-------------------|---------|
| KA 8100 | KA 2200N-KA 5200N | 1 pr/bx |

- Non-swivel front wheels with inner legs.
- Height adjustable.
- Easily used outdoors and over rough grounds.
- Adds 1" to overall height of walker.



10219FRD-1



10222FRD-1



10221FRD-1

Glider Walker

| Item # | Product Size | Product Handle Height Range | Product Depth (opened at base) | Product Width (inside back legs) | Product Weight | Product Weight Capacity | UOM |
|------------|--------------|-----------------------------|--------------------------------|----------------------------------|----------------|-------------------------|------|
| 10219FRD-1 | Standard | 29" - 37" | 25" | 24" | 8 lbs. | 175 lbs. | 1/bx |
| 10222FRD-1 | Junior | 22" - 33" | 25" | 19" | 7 lbs. | 175 lbs. | 1/bx |
| 10221FRD-1 | Pediatric | 12" - 21" | 20" | 17.5" | 5.5 lbs. | 175 lbs. | 1/bx |

- Aluminum frame with flame-red crossbar.
- Height and angle-adjustable handles.
- Handles rotate inward and outward for adjustable width.
- 90° handle on standard and junior model; 45° handle on pediatric model.
- Rubber wheel allows walker to roll easily and smoothly over irregular surfaces.
- Tool-free adjustments.
- Wide base provides stability and balance.

specifications

| | |
|-------------------------------------|------------------|
| Carton shipping Weight (10219FRD-1) | 11 lbs. |
| Carton Shipping Weight (10222FRD-1) | 10 lbs. |
| Carton Shipping Weight (10221FRD-1) | 7 lbs. |
| Warranty | Limited Lifetime |



HCPCS
E0141



A

Platform Attachment

| Item # | Product Description | UOM |
|---------|---------------------|------|
| 10219PA | Standard | 1/pr |
| 10222PA | Junior | 1/pr |
| 10221PA | Pediatric | 1/pr |

- Height adjustable.
- Can be mounted anywhere along the side bars.
- The handgrips extend, pronate and supinate.
- Can be mounted with an anterior or posterior tilt.
- Brackets can be mounted facing the inside or the outside of the side bars allowing for lateral adjustability.

Side Walker

| Item # | Product Description | UOM |
|---------|---------------------|------|
| 10240-1 | Adult (Hemi-Walker) | 1/ea |

- Designed for individuals with the use of only one hand or arm.
- Lighter than a walker and more stable than a cane.
- Folds easily with one hand. (Figure A)
- Ideal for users who do not require a walker but need a wider base for support.

specifications

| | |
|------------------------|---------------------|
| Overall Dimensions | 18.5" (W) x 16" (D) |
| Height | 28.5" - 35.5" |
| Weight | 4 lbs. |
| Carton Shipping Weight | 5 lbs. |
| Weight Limit | 300 lbs. |
| Warranty | Limited Lifetime |



A



A

Bariatric Side Walker

| Item # | Product Description | UOM |
|--------|---------------------|------|
| 10242 | Adult (Hemi-Walker) | 1/ea |

- Designed for individuals with the use of only one hand or arm.
- Lighter than a walker and more stable than a cane.
- Folds easily with one hand. (Figure A)
- Ideal for users who do not require a walker but need a wider base for support.
- 500 lb. Weight Limit.

specifications

| | |
|------------------------|---------------------|
| Overall Dimensions | 18.5" (W) x 16" (D) |
| Height | 29.5" - 36" |
| Weight | 4 lbs. |
| Carton Shipping Weight | 5 lbs. |
| Weight Limit | 500 lbs. |
| Warranty | Limited Lifetime |

Platform Walker/Crutch Attachment

| Item # | Product Description | For Use With | Product Height | UOM |
|---------|---|--|----------------|------|
| 10105-1 | | Adult and Junior Walkers and Aluminum Crutches | 37" - 53" | 1/ea |
| 10105-2 | | Adult and Junior Walkers and Aluminum Crutches | 37" - 53" | 2/cs |
| 10116 | Universal Platform Walker/Crutch Attachment | All Drive walkers and Invacare® Paddle Release Walkers | 36" - 55" | 1/ea |

- Designed for those who are unable to grip the walker.
- Cradles forearm with soft vinyl padding for comfort.
- Large top bracket allows the platform attachment to be placed over the hand grip of the walker.
- Adjustable strap securely holds arm in place.

specifications

| | |
|--|------------------|
| Overall Width (10105 10116) | 5.75" 6" |
| Overall Depth | 18" |
| Weight (10105 10116) | 2 7 lbs. |
| Carton Shipping Weight (1010501 10105-2 10116) | 3 6 10 lbs. |
| Weight Limit | 300 lbs. |
| Warranty | Limited Lifetime |



Bariatric Platform Walker/Crutch Attachment

| Item # | Product Description | For Use With | UOM |
|-----------|--|--|------|
| 10105HD-2 | Universal Platform Walker and Aluminum Crutch Attachment | Adult and Junior Walkers and Aluminum Crutches | 2/cs |

- Designed for those who are unable to grip the walker.
- Cradles forearm with soft vinyl padding for comfortable weight bearing.
- Large top bracket allows the platform attachment to be placed over the hand grip of the walker.
- Adjustable strap securely holds arm in place.
- Attachment for Underarm Crutch ([Figure A](#)) for item #10105 must be ordered separately. (Cannot be used with Bariatric model)
- Bariatric 500 lb. Weight Limit.

specifications

| | |
|------------------------|----------------------|
| Overall Dimensions | 6.5" (W) x 18.5" (D) |
| Height | 40.5" - 55" |
| Weight | 4 lbs. |
| Carton Shipping Weight | 9 lbs. |
| Weight Limit | 500 lbs. |
| Warranty | Limited Lifetime |

Universal 5" Walker Wheels

Dual-Hole Adjustment Column and Two Sets of Rear Glides

| Item # | Color | Product Size | UOM |
|----------|---------------------------------------|--------------|------|
| 10109 | Silver tubing, gray tire and silencer | 5" Wheels | 1/pr |
| 10109-BL | Blue tubing, black tire and silencer | 5" Wheels | 1/pr |
| 10109-G | Green tubing, black tire and silencer | 5" Wheels | 1/pr |
| 10109-P | Pink tubing, black tire and silencer | 5" Wheels | 1/pr |

- Converts folding walker into wheeled walker.
- Fits most manufacturers' walkers.
- Allows for 8 height adjustments.
- Rubber wheel allows walker to roll easily and smoothly over irregular surfaces.
- Comes with rear glide caps, (item #10107) and glide covers (item #10107C) allowing use on all surfaces.
- 350 lb. Weight Limit.

specifications

| | |
|--------------|----------|
| Weight Limit | 350 lbs. |
|--------------|----------|



5" Walker Wheels For Adult Walker

5" Walker Wheels For Junior Walker

| Item # | For Use With | UOM |
|--------|------------------------|------|
| 10126 | Models 10291 and 10293 | 1/pr |

- Converts folding walker into wheeled walker.
- Allows for 8 height adjustments.
- Rubber wheel allows walker to roll easily and smoothly over irregular surfaces.
- 300 lb. Weight Limit.

| Item # | For Use With | UOM |
|--------|------------------------|------|
| 10127 | Models 10292 and 10294 | 1/pr |

- Converts folding walker into wheeled walker.
- Allows for 4 height adjustments.
- Rubber wheel allows walker to roll easily and smoothly over irregular surfaces.
- 300 lb. Weight Limit.

specifications

| | |
|--------------|----------|
| Weight Limit | 300 lbs. |
|--------------|----------|

specifications

| | |
|--------------|----------|
| Weight Limit | 300 lbs. |
|--------------|----------|



3" Universal Walker Wheels with Two Sets of Rear Glides

| Item # | Color | UOM |
|-----------|---------------|------|
| 10111 | Silver Tubing | 1/pr |
| 10100NBL | Blue Tubing | 1/pr |
| 10100NG | Green Tubing | 1/pr |
| 10100NP | Pink Tubing | 1/pr |
| 10100NBLK | Black Tubing | 1/pr |

- Converts folding walker onto wheeled walker.
- Allows for 8 height adjustments.
- Rubber wheel allows walker to roll easily and smoothly over irregular surfaces.
- Comes with rear glide caps, (item #10107) (Figure A) and glide covers (item #10107C) (Figure B), allowing use on all surfaces.
- For use with most leading manufacturers' 1" walkers.
- 350 lb. Weight Limit.

specifications

| | |
|--------------|----------|
| Weight Limit | 350 lbs. |
|--------------|----------|

5" Walker Wheels with Two Sets of Rear Glides

| Item # | For Use With | Color | UOM |
|---------|---------------------------------|---------------------------------------|------|
| 10128 | 10248N and 10253 Series Walkers | Silver tubing, gray tire and silencer | 1/pr |
| 10128BL | 10248N and 10253 Series Walkers | Blue tubing, black tire and silencer | 1/pr |
| 10128P | 10248N and 10253 Series Walkers | Pink tubing, black tire and silencer | 1/pr |

- Converts folding walker into wheeled walker.
- Allows for 11 height adjustments.
- Rubber wheel allows walker to roll easily and smoothly over irregular surfaces.
- Comes with rear glide caps (item #10107) (Figure A) and glide covers (item #10107C) (Figure B), allowing use on all surfaces.
- 350 lb. Weight Limit.

specifications

| | |
|--------------|----------|
| Weight Limit | 350 lbs. |
|--------------|----------|



5" Bariatric Walker Wheels with Two Sets of Rear Glides

| Item # | Product Description | Color | UOM |
|---------|--------------------------|------------------------------|------|
| 10118SV | 5" Bariatric Wheels | Silver Vein with Gray Wheels | 1/pr |
| 10118T | 5" Tall Bariatric Wheels | Silver Vein with Gray Wheels | 1/pr |

- Converts heavy-duty walker into wheeled walker.
- Allows for 8 height adjustments.
- 10118T adds 4" to height of walker.
- Dual rubber wheels allow walker to roll easily and smoothly over irregular surfaces.
- Comes with rear glide caps, (item #10107) (Figure A) and glide covers (item #10107C) (Figure B), allowing use on all surfaces.
- For use with items # 10220-2 and 10223-2.
- 500 lb. Weight Limit.

specifications

| | |
|--------------|----------|
| Weight Limit | 500 lbs. |
|--------------|----------|

5" Bariatric Walker Wheels, Combo Pack

| Item # | Product Description | UOM |
|----------|---|---------|
| 10118CSV | Combo Pack Includes 1 pair Bariatric Walker Wheels (item # 10118TSV) and 1 Pair Extension Legs | 1Set/bx |

- Converts heavy-duty walker into wheeled walker.
- Allows for 8 height adjustments.
- Dual rubber wheels allow walker to roll easily and smoothly over irregular surfaces.
- Comes with rear glide caps, (item #10107) (Figure A) and glide covers (item #10107C) (Figure B), allowing use on all surfaces.
- For use with items #10220-2 and 10223-2.
- 500 lb. Weight Limit.

specifications

| | |
|--------------|----------|
| Weight Limit | 500 lbs. |
|--------------|----------|



Swivel Wheel with Lock & Two Sets of Rear Glides

| Item # | Product Description | UOM |
|--------|---------------------------|--------|
| 10115 | 5" Swivel Wheel with Lock | 1pr/bx |

- Can be locked for fixed wheel or unlocked for a swivel wheel. (Figure A - B)
- Allows for 5 adjustments.
- Rubber wheels allow walker to roll easily and smoothly over irregular surfaces.
- Comes with rear glide caps, (item #10107) and glide covers (item #10107C), allowing use on all surfaces.
- 350 lb. Weight Limit.

specifications

| | |
|--------------|----------|
| Weight Limit | 350 lbs. |
|--------------|----------|



Swivel Wheel with Lock & Two Sets of Rear Glides

| Item # | Product Description | UOM |
|--------|---------------------------|--------|
| 10113 | 3" Swivel Wheel with lock | 1pr/bx |

- Can be locked for fixed wheel or unlocked for a swivel wheel. (Figure A - B)
- Allows for 8 height adjustments.
- Rubber wheels allow walker to roll easily and smoothly over irregular surfaces.
- Comes with rear glide caps, (item #10107) and glide covers (item #10107C), allowing use on all surfaces.
- 350 lb. Weight Limit.

specifications

| | |
|--------------|----------|
| Weight Limit | 350 lbs. |
|--------------|----------|



10108W

10108



Tall Extension Legs with Wheels, Combo Pack (Adds 4")

| Item # | Product Description | UOM |
|---------|--|----------|
| 10108WC | Combo Pack 1 pair Extension Legs with 5" Wheels and 1 pair Extension Legs | 1 set/bx |

- Allows for height adjustments from 36" to 43" (walker handgrip to floor). (Adds 4" to height of walker).
- 5" rubber wheels allow walker to roll easily and smoothly over irregular surfaces.
- For use with all Drive and most leading manufacturers' walkers.
- Includes rear glide caps, (item #10107) and glide covers (item #10107C), allowing use on all surfaces.
- 350 lb. Weight Limit.

specifications

| | |
|--------------|----------|
| Weight Limit | 350 lbs. |
|--------------|----------|

Tall Extension Legs (Adds 4")

| Item # | UOM |
|--------|------|
| 10108 | 4/bx |

- Makes height adjustments from 36" to 43" on a standard adult walker. (Adds 4" to height of walker)
- Allows for 8 height adjustments.
- Sturdy anodized, extruded aluminum.
- For use with all Drive and most leading manufacturers' walkers.
- 350 lb. Weight Limit.

specifications

| | |
|--------------|----------|
| Weight Limit | 350 lbs. |
|--------------|----------|



Rear Glide Brakes with 3" Wheels

| Item # | UOM |
|--------|------|
| 10106 | 1/pr |

- For use on back of wheeled walker.
- 3" wheel glides along surface until weight is applied to hand grip, which activates rubber-tip brake.
- Allows for 8 height adjustments.
- 350 lb. Weight Limit.

specifications

| | |
|--------------|----------|
| Weight Limit | 350 lbs. |
|--------------|----------|

Rear Glide Brakes

| Item # | UOM |
|--------|------|
| 10104 | 1/pr |

- For use on back of wheeled walker.
- Nylon tip glides along surface until weight is applied to hand grip, which activates rubber-tip brake.
- Allows for 8 height adjustments.
- 350 lb. Weight Limit.

specifications

| | |
|--------------|----------|
| Weight Limit | 350 lbs. |
|--------------|----------|



RTL6062



RTL6063



Fashion Walker Tote Bags

| Item # | Product Description | UOM |
|---------|-----------------------|------|
| RTL6062 | Solid, Quilt Fabric | 3/cs |
| RTL6063 | Pattern, Quilt Fabric | 3/cs |

- Double sided with two pockets 14" (Wide) by 14" (Deep), held in place with hook and loop.
- Easier to carry lighter items such as magazines, newspapers or wallet.

specifications

| | |
|----------|------------------|
| Warranty | Limited 6 Months |
|----------|------------------|

Carry Pouch for Walkers and Rollators, Large

| Item # | UOM |
|----------|------|
| RTL10250 | 1/ea |

- For use with most folding walkers.
- 2 large pockets with Velcro®-type closure.
- Detachable mesh pocket.
- Made of durable, easy-to-clean nylon.
- Size: 17.5" (L) x 18.5" (W) x 1" (D).
- Easily attaches to walker.

specifications

| | |
|----------|-----------------|
| Warranty | Limited 90 Days |
|----------|-----------------|



Walker Basket Carry Liner

| Item # | UOM |
|---------|------|
| 10259-1 | 1/ea |

- Easily attaches to walker basket.
- Inside zippered pouch.
- 2 inside mesh pouches.
- Made of durable, easy-to-clean nylon.
- Size: 6" (L) x 16" (W) x 7" (D).

specifications

| | |
|----------|-----------------|
| Warranty | Limited 90 Days |
|----------|-----------------|

Carry Pouch for Walker

| Item # | Product Size | UOM |
|----------|--------------|------|
| 10258L-1 | Large | 1/bx |
| 10258-1 | Small | 1/bx |

- Easily attaches to walker.
- For use with most folding walkers.
- Large zippered pouch.
- 3 mesh pouches.
- Made of durable, easy-to-clean nylon.
- Large size: 14" (W) x 10.5 (L) x 2" (D). Small size: 9" (W) x 8.5" (L) x 1" (D).

specifications

| | |
|----------|---------|
| Warranty | 90 Days |
|----------|---------|



Side Walker Carry Pouch

| Item # | UOM |
|---------|------|
| 10255-1 | 1/ea |

- Easily attaches to side frame of walker.
- For use with most folding walkers.
- 2 large pockets with Velcro®-type closure.
- Made of durable, easy-to-clean nylon.
- Detachable mesh pocket.
- Size: 14.5" (W) x 2.5" (D) x 15" (L).

specifications

| | |
|----------|-----------------|
| Warranty | Limited 90 Days |
|----------|-----------------|



Tennis Ball Glides with Replaceable Glide Pads

| Item # | Product Description | UOM |
|--------|------------------------------------|---------|
| 10119 | Tennis Ball with Glides in Can | 1pr/can |
| 10121 | Tennis Ball Box | 1/pr |
| 10123 | Replacement Tennis Ball Glide Pads | 4/pkg |

- Provides a quiet, smooth and durable glide experience when used with a walker.
- Easy and safe to install.
- Lasts longer than plastic glide cap.
- Patent Pending.

specifications

| | |
|----------|----------------|
| Warranty | Limited 1 Year |
|----------|----------------|



RTL100014



Sneaker Walker Glides

| Item # | UOM |
|-----------|---|
| RTL100014 | 1 pr of Sneaker Walker Glides Per Package, 4 pks/cs |
| 100014 | 1/pr |

- Stylish tennis shoe glides personalize the walker.
- Allows walker to glide easily and smoothly over most surfaces.
- Fits all size tubing. Knob turns on bottom of glide for different tubing diameter. (Figure A - B)
- Easy, tool-free installation.
- Replaces rubber tip.

Glide Cover



10107C

- Tip style glide cover for use over existing rubber tip.
- Allows walker to "glide" across most surfaces.
- Easily added to walker. Makes less noise.
- Very quiet on outdoor surfaces.

| Item # | Color | UOM |
|---------|-------|------|
| 10107C | Black | 1/pr |
| 10107CG | Gray | 1/pr |

Glide Cap



- Allows walker to glide easily and smoothly over most surfaces.
- Replaces rubber tip.

| Item # | For Use With | UOM |
|--------|--|------|
| 10117C | Model #s 10291 10292 10293 10294 | 1/pr |

Glide Cap



- Allows walker to glide easily and smoothly over most surfaces.
- Replaces rubber tip.

| Item # | UOM |
|--------|------|
| 10107 | 1/pr |

Mushroom Glide Cap



| Item # | Product Description | For Use With | Color | UOM |
|-----------|---------------------|--|-------|------|
| 10107MC-G | Mushroom Caps | 10253, 10206, 10210, 10211, 10226, 10227 | Gray | 1/pr |
| 10107MC-B | Mushroom Caps | 10253, 10206, 10210, 10211, 10226, 10227 | Black | 1/pr |

Glide Ski



| Item # | Color | UOM |
|--------|-------|------|
| 10110 | White | 1/pr |
| 10110B | Black | 1/pr |

- For use with all 1" folding walkers.
- Easily added to walker.
- Durable composite construction ensures long life.
- Ideal on most surfaces.
- For use with all Drive and most leading manufacturers' walkers.


Walker Tips



| COMPARE TO | |
|------------|----------|
| drive™ | SUNRISE® |
| RTL10389G | 1202 |

- Provides stability to walker.
- For use with all Drive and most leading manufacturers' with 1" diameter walkers.

| Item # | Color | Product Size | UOM |
|-----------|-------|--------------|---------------------------|
| RTL10389G | Gray | 1" | 1pr/bx, 8 retail boxes/cs |



Walker Basket

| Item # | Product Description | UOM |
|-----------|------------------------------|--|
| 10200B | Insert Included with Basket | 1/ea |
| 10200BR-6 | Inserts Included with Basket | 6 Flat Retail Boxes Included in Carton, 6/cs |

- For use with all 1" folding walkers.
- 10200B comes with plastic insert tray with cup holder (10200B-1). (Figure A)
- Includes retail boxes in shipping carton. (Figure B).

accessories

| | | |
|----------|----------------------------------|------|
| 10200B-1 | Walker Basket Insert Replacement | 1/ea |
|----------|----------------------------------|------|


Folding Walker Basket

| Item # | UOM |
|------------|------|
| RTL10200FB | 1/bx |

- For use with most manufacturers' 1" walkers.
- Folds flat for storage.
- Comes with plastic insert tray with cup holder (item #10200B-1). (Figure A)

accessories

| | | |
|----------|----------------------------------|------|
| 10200B-1 | Walker Basket Insert Replacement | 1/ea |
|----------|----------------------------------|------|



Clever Lite Rollator Basket

Universal Cup Holder

| Item # | For Use With | UOM |
|--------|----------------------|------|
| 10230B | 10230, 10230J, 10243 | 1/ea |

- For use with Clever Lite rollators.
- Easy to attach and remove.

| Item # | UOM |
|-----------|------|
| STDS1040S | 1/ea |

- New clamp allows cup holder to swivel to keep contents level.
- Easily attaches to wheeled walkers, rollators, transport chairs/wheelchairs with tubing diameter from $\frac{5}{8}$ " to 1".
- Special clamp insert can be added or removed as necessary to accommodate different size tubing.
- Made with durable, lightweight plastic.

specifications

| | |
|-------------|------------------|
| Dimensions: | 3" (W) x 3" (L). |
| Warranty | Limited 1 Year |



COMPARE TO
drive™
JOBAR®
RTL10131 JB5770

Light & Go Mobility Light

| Item # | UOM |
|---------|-------|
| RTL1100 | 10/cs |

- Patented motion and light-activated technology.
- For hands-free use with wheelchairs (Figure A), Walkers (Figure B) and Rollators. (Figure C)
- Comes with protective rubber boot.
- 3 bright, long-life LEDs for enhanced visibility.
- Runs on 2 standard AA batteries (included).
- Easy-to-use clamp with foam padding for secure grip.
- Retail Packaged.

E-Z Walker Caddy

| Item # | UOM |
|----------|------|
| RTL10131 | 3/cs |

- Allows personal items to be carried from room to room.
- Includes a removable top that conceals personal items.
- Fits most manufacturers' 2-button walkers.
- Easily attaches with side clips.
- Made of easy-to-clean durable plastic.
- Easy to install.
- Designed by Chris Clarke.
- Size: 15½" (W) 5" (H) x 12 ½" (D).

specifications

Warranty

Limited 1 Year



10124



COMPARE TO
drive ROSE HEALTHCARE®
10124 1212



10125



Universal Walker Tray with Cup Holder

| Item # | UOM |
|--------|------|
| 10124 | 1/ea |

- Allows personal items to be carried from room to room.
- Contains one cup holder.
- Fits most manufacturers' walkers.
- Made of easy-to-clean durable plastic.
- Easy to install.
- Size: 23" (W) x 1.5" (H) x 17" (D).

specifications

| | |
|----------|----------------|
| Warranty | Limited 1 Year |
|----------|----------------|

Walker Tray with Cup Holders

| Item # | UOM |
|--------|------|
| 10125 | 1/ea |

- Allows personal items to be carried from room to room.
- Contains two cup holders.
- Fits most manufacturers' 2-button walkers.
- Easily folds down when not in use.
- Made of easy-to-clean, durable plastic.
- Easy to install.
- Size: 16" (W) x 2" (H) x 12" (D).

specifications

| | |
|----------|----------------|
| Warranty | Limited 1 Year |
|----------|----------------|



2014 WALKERS PRODUCT CATALOG

Drive Medical Design & Manufacturing

99 Seaview Boulevard
Port Washington, NY 11050

t: 877.224.0946 | f: 516.998.4601

www.drivemedical.com



PRINTED IN CHINA

Copyright© 2014 Medical Depot, Inc. dba Drive Medical Design and Manufacturing®, 99 Seaview Boulevard, Port Washington, NY 11050. All trademarks used in association with the sale of products of Drive Medical Design and Manufacturing are trademarks owned by Medical Depot, Inc. All other trademarks, trade names, service marks and logos referenced herein belong to their respective companies.