

AIRCAST®

AirSport+™
Simple. Stable. Smarter.



The Aircast® AirSport+™ is a clinically proven* stirrup design that aids in ankle stabilization for moderate to severe ankle sprains and chronic instabilities.

The anatomically correct design and durable, breathable materials offer a stable, comfortable fit.

Along with the Aircast technology, simple step-in design and intuitive, unified closure system, the AirSport+ provides a functional management option to help address ankle injuries and rehabilitation.

Recommended for those with:

Chronic ankle instability - Weak ankles due to ligament damage and over stressed ligaments

Acute ankle sprains & strains - Over stretched and/or torn ankle ligaments

Prophylaxis - Helping to prevent reoccurring ankle injury and returning to daily activities

POSTERIOR HIGH ANKLE STRAP

- + Helps provide compression to the tibia/fibula to aid in improved ankle stabilization

EASY, COMFORTABLE FIT

- + Simple, easy application with posterior step-in design
- + Breathable materials

ATF CROSS STRAP

- + Anterior TaloFibular (ATF) cross strap & forefoot wrap provides additional compression & stability to the ankle joint, helping to minimize inversion



STIRRUP DESIGN

- + Bilateral anatomically-shaped shells help minimize the risk of ankle rollover (inversion/eversion)
- + Helps to provide medial/lateral support for unstable ankles
- + Air/foam combination enhances comfort

BOA® FIT SYSTEM FOR CUSTOMIZABLE COMPRESSION

- + Allows patients to customize compression levels to their needs
- + Criss-cross laces generate uniform compression designed to minimize pressure points





AirSport+™

Simple. Stable. Smarter.

Part Number		Shoe Size		
Right	Left	Size	Male	Female
02RXSR	02RXSL	XS	< 3	< 5
02RSR	02RSL	S	3 – 5	5 – 6 ½
02RMR	02RML	M	5 – 8 ½	6 ½ – 10
02RLR	02RLL	L	8 ½ – 12	10 – 13 ½
02RXLR	02RXLL	XL	>12	>13 ½



T 800.526.8785 | D 760.727.1280 | F 800.457.4221
 2900 Lake Vista Drive | Dallas, TX 75067 | U.S.A.
aircast.com

Copyright © 2018 by DJO, LLC
 MKT00-7796 Rev A

* Data on file.
 Boa® is a registered trademark of Boa® Technology Inc.
 Individual results may vary. Neither DJO Global, Inc. nor any of its subsidiaries dispense medical advice. The contents of this brochure do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.