

# Polyisoprene Surgical Gloves

## From A Name You Can Trust.

*Designed and manufactured to provide strength, protection, and value.*



### Features & Benefits:

- ✓ *Similar structure to natural rubber latex* - supreme comfort
- ✓ *Fully anatomical shape* - reduced hand fatigue and great fit
- ✓ *Skin friendly accelerators* - Low allergenic potential
- ✓ *Micro-rough texture* - superior instrument handling and grip
- ✓ *Multi-layer coating* - better donning capability
- ✓ *Extended & Beaded Cuff* - prevents roll down effect

### Ordering Information

Size	Item No.
5.5	570-1119
6.0	570-1120
6.5	570-1121
7.0	570-1122
7.5	570-1123
8.0	570-1124
8.5	570-1125
9.0	570-1126

50 pairs per box,  
4 boxes per case

### Specifications:

**Application:** Surgical Glove

**Chemo Rated:** N/A

**Color:** Cream

**Cuff Style:** Extended

**Glove Exterior:** Fully Micro Textured

**Glove Interior:** Polymer Coated

**Glove Mil Thickness:** 7.6 mils

**Material:** Polyisoprene

**Powder Content:** Powder Free

**Sterility:** Sterile

**Usage:** Single Use

**Latex Free Indicator:** Not made with Natural Rubber Latex



**RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!**

# Criterion<sup>®</sup> CR

## Powder-Free Polychloroprene Surgical Gloves



Item Code	Description
(112-7079)	Criterion <sup>®</sup> CR (Surgical) 5½ Box, 50 pair
(112-7080)	Criterion <sup>®</sup> CR (Surgical) 6 Box, 50 pair
(112-7081)	Criterion <sup>®</sup> CR (Surgical) 6½ Box, 50 pair
(112-7082)	Criterion <sup>®</sup> CR (Surgical) 7 Box, 50 pair
(112-7083)	Criterion <sup>®</sup> CR (Surgical) 7½ Box, 50 pair
(112-7085)	Criterion <sup>®</sup> CR (Surgical) 8 Box, 50 pair
(112-7086)	Criterion <sup>®</sup> CR (Surgical) 8½ Box, 50 pair
(112-7088)	Criterion <sup>®</sup> CR (Surgical) 9 Box, 50 pair

<b>Application:</b>	Surgical Glove
<b>Chemo Rated:</b>	N/A
<b>Color:</b>	White
<b>Cuff Style:</b>	Beaded Cuff
<b>Glove Exterior:</b>	Fully Textured
<b>Glove Interior:</b>	Polymer Coated
<b>Glove Mil Thickness:</b>	Min. 5.9 mil
<b>Material:</b>	Polychloroprene
<b>Powder Content:</b>	Powder-Free
<b>Sterility:</b>	Sterile
<b>Usage:</b>	Single Use
<b>Latex-Free Indicator:</b>	Not Made with Natural Rubber Latex

### Glove Box Holders



**Acrylic, Vertical**  
**(102-7201)** Four boxes  
**(102-8344)** Three boxes  
**(102-6468)** Two boxes  
**(102-6467)** One box

**Wire, Vertical**  
**(112-5474)** Three boxes  
**(112-5472)** Two boxes

**Wire, Horizontal**  
**(112-5475)** Three boxes  
**(112-5473)** Two boxes  
**(112-5471)** One box



**RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!**

# Criterion® Powder-Free Latex Surgical Gloves



Item Code	Description
(570-0635)	Criterion® 5½ Box, 50 pair
(570-0636)	Criterion® 6 Box, 50 pair
(570-0637)	Criterion® 6½ Box, 50 pair
(570-0638)	Criterion® 7 Box, 50 pair
(570-0639)	Criterion® 7½ Box, 50 pair
(570-0640)	Criterion® 8 Box, 50 pair
(570-0641)	Criterion® 8½ Box, 50 pair
(570-0642)	Criterion® 9 Box, 50 pair

<b>Application:</b>	Surgical Glove
<b>Chemo Rated:</b>	N/A
<b>Color:</b>	Natural
<b>Cuff Style:</b>	Beaded Cuff
<b>Glove Exterior:</b>	Fully Finger Textured
<b>Glove Interior:</b>	Chlorinated
<b>Glove Mil Thickness:</b>	Min. 7.1 mil
<b>Material:</b>	Latex
<b>Powder Content:</b>	Powder-Free
<b>Sterility:</b>	Sterile
<b>Usage:</b>	Single Use
<b>Latex-Free Indicator:</b>	Not Made with Natural Rubber Latex

### Glove Box Holders



**Acrylic, Vertical**  
**(102-7201)** Four boxes  
**(102-8344)** Three boxes  
**(102-6468)** Two boxes  
**(102-6467)** One box

**Wire, Vertical**  
**(112-5474)** Three boxes  
**(112-5472)** Two boxes

**Wire, Horizontal**  
**(112-5475)** Three boxes  
**(112-5473)** Two boxes  
**(112-5471)** One box



**RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!**