



In alignment with the FDA Guidance Document, Unique Device Identification: Policy Regarding Compliance Dates for Class I and Unclassified Devices, Direct Marking, and Global Unique Device Identification Database Requirements for Certain Devices (September 24, 2022), End-Users may follow this information to assemble the UDI.

1. For devices without direct markings, users may assemble the UDI information of the product by referencing the serial numbers engraved on the devices itself, in addition to referencing the table listed below. The serial numbers are located on the devices themselves as an engraving and/or name plates. Below are several examples of serial numbers depicted on the products for reference.

Example 1. – Essentials E-Type Lowspeed Contra Angle Sheath



HS REF: 112-7304

Item Description: Essentials E-Type Lowspeed Contra Angle Sheath

Serial Number: 21120441

Engraved Serial Number

Example 2. – Maxima XTEND AX Maintenance System



HS REF: 570-1174

Item Description: Maxima XTEND AX Maintenance System

Serial Number: 22070019

Serial Number



Example 3. – Acclean Ultrasonic Scaler



HS REF: 570-2771

Item Description: Acclean Ultrasonic Scaler

Serial Number: 22016704



Serial Number

2. Users may then assemble the UDI using the serial number obtained from step 1 and referencing the corresponding table below with associated part number and GTIN.

UDI may be obtained as follows: '(01)'GTIN + '(21)'SN = UDI.

- In example #1, the human-readable UDI Information for 112-7304 will be:
(01)00304040058268(21)21120441
- In example #2, the human-readable UDI Information for 570-1174 will be:
(01)00304040099520(21)22070019
- In example #3, the human-readable UDI Information for 570-2771 will be:
(01)00304040155837(21)22016704

3. This information can be reprinted from the product web page on the Henry Schein website at www.HenrySchein.com.



UDI Codes Cross-Matching Summary Table

HS REF	Item Description	GTIN Number
1127297	Essentials C-Type Push-Button Handpiece	00304040058190
1127298	Essentials C-Type Handpiece	00304040058206
1127299	Essentials C-Type Mini Push-Button Handpiece	00304040058213
1127300	Essentials C-Type Mini Handpiece	00304040058220
1127301	Essentials E-Type Lowspeed System	00304040058237
1127302	Essentials E-Type Lowspeed 20,000 RPM Motor	00304040058244
1127303	Essentials E-Type Lowspeed Straight Attachment	00304040058251
1127304	Essentials E-Type Lowspeed Contra Angle Sheath	00304040058268
5700072	Maxima XTEND Maintenance	00304040058282
5700105	Maxima ELITE Handpiece	00304040058305
5700106	Maxima ELITE Handpiece Mini	00304040058312
5700741	Precision-Torque 2Cr SS	00304040064375
5701072	Maxima PRO 2 Highspeed Handpiece	00304040090046
5701073	Maxima PRO 2L Highspeed Optic Handpiece	00304040090060
5701074	Maxima PRO 2S Highspeed Swivel Handpiece	00304040090053
5701174	Maxima XTEND AX Maintenance	00304040099520
5701431	Essentials Pro 2 Handpiece	00304040110089
5701432	Essentials Pro 2M Handpiece	00304040110096
5701433	Essentials Pro 2S Handpiece	00304040110102
5701434	Essentials Pro 2SM Handpiece	00304040110119
5701449	Maxima Pro 2L HP Mini F/O	00304040110898
5701450	Maxima Pro 2S HP Mini N/O	00304040110881
5701451	Maxima Pro 2L HP Mini F/O	00304040110874
5701531	MaxPro Personalized	00304040117385
5701532	MaxPro ML Personalized	00304040117415
5700107	Maxima ELITE L Standard Head	00304040058329
5700401	Precision Torque 2 Push Button	00304040058374
5700402	Precision Torque 2L PB FO	00304040058381
5700486	Master 3 Standard Push Button	00304040058398
5700487	Master 3 Standard HighSpeed HP	00304040058404
5700488	Master 3 Fiber-Optic PB	00304040058411
5700489	Master 3 Mini Push Button	00304040058428
5700490	Master 3S PB Swivel Type	00304040058435
5700491	Master Torque 3 E-Type LS HP	00304040058442
5700492	Master Torque 3 E-Type LS HP	00304040058459
5700494	Master Torque 3 E-Type Motor	00304040058466
5700495	Master Torque 3 E-Type Motor	00304040058473



HS REF	Item Description	GTIN Number
5700496	Master 3 Mini High-Speed HP	00304040058480
5700497	Master Torque 3 Straight	00304040058497
5700498	Master Torque 3 Contra Angle	00304040058503
5700502	Master Torque 3 4:1 Straght f/	00304040058541
5700684	Maxima Electric System	00304040116104
5700685	Maxima Electric Micromotor	00304040116111
5700687	PDS MaxPro Highspeed Handpiece	00304040062067
5700688	PDS MaxPro M Highspeed Handpiece	00304040062074
5700691	MaxPro Highspeed Handpiece	00304040062012
5700692	MaxPro M Highspeed Handpiece	00304040062029
5700697	MaxPro L Highspeed Handpiece	00304040062098
5700698	MaxPro ML Highspeed Handpiece	00304040062081
5701533	MaxPro M Personalized	00304040117392
5701534	MaxPro L Personalized	00304040117408
5702094	Maxima PRO 2L AR Highspeed Handpiece	00304040133682
5702095	Maxima PRO 2S AR Handpiece	00304040133699
5702771	Acclean Ultrasonic Scaler	00304040155837
5702855	Acclean Hygiene Handpiece 2S	00304040160091
5702861	Acclean Hygiene Handpiece 2	00304040160107
5702863	PRO 2-K Lowspeed Attachment Contra Angle Sheath	00304040159613
5702864	PRO 2-M Lowspeed Attachment Straight Attachment	00304040159637
5702865	PRO 2-M Lowspeed Attachment Contra Angle Sheath	00304040159644
5702866	PRO 2-N Lowspeed Attachment Straight Attachment	00304040159651
5702867	PRO 2-N Lowspeed Attachment Contra Angle Sheath	00304040159668
5702868	PRO 2-S Lowspeed Attachment Straight Attachment	00304040159675
5702869	PRO 2-S Lowspeed Attachment Contra Angle	00304040159682
5702871	T-Style Lowspeed Attachment Contra Angle	00304040159705