



PowerSpin™ Centrifuge Fixed Speed

- Henry Schein Brand
- Fixed speed at 3,400 RPM
- 6 place fixed angle rotor
- 30 minute timer
- Stainless steel rotor chamber
- Lid latch safety switch

Order model # 112-6905



PortaSpin™ Mobile Centrifuge

- Henry Schein Brand
- Plugs into a car 12V Power Outlet
- Fixed speed at 3,500 RPM
- 6 place fixed angle rotor
- 30 minute timer
- Stainless steel carrying handle
- Optional AC power adapter

Order model # 112-6906

PowerSpin™ Centrifuge Variable

Variable speed centrifuge with color-coded panel for blood, urine and fecal specimens.

- Variable speed 500-4,000 RPM
- 6 or 8 place fixed angle rotor
- 30 minute timer
- Stainless steel rotor chamber
- Lid latch safety switch

Order model # 112-6907 6 place

Order model # 112-6909 8 place

Henry Schein Brand



Microhematocrit Centrifuge

- Henry Schein Brand
- 24 place rotor for 75mm capillary tubes
- Fixed speed at 12,000 RPM
- 5 minute timer
- All metal housing
- Safety switch
- Microhematocrit Reader included

Order model # 112-6910





Laboratory Rotator

- 11x11" platform
 - Variable speed 25-250 RPM
 - 30 minute timer
 - Humidity cover included
- Order model # 112-6912

Vortex Tube Mixer

- Variable speed up to 3,000 RPM
 - Switch from Auto to Touch mode
 - Power "on" LED indicator
- Order model # 112-6913



Blood Tube Rocker


- 15 place platform for up to 120mm tubes
 - 25 RPM fixed speed
 - Compact size for minimal footprint
- Order model # 112-6911



Phlebotomy Tray

- Large carrying capacity
 - Includes racks (sharps and alcohol dispenser not included)
 - Handy pull out drawer
- Order model # 112-6916



<i>Distributed by (in US only):</i>	
 <small>(manufacturer)</small>	HENRY SCHEIN INC. 135 Duryea Road Melville, NY 11747 USA EC REP
	Henry Schein U.K. Holdings Ltd. Medcare House, Centurion Close Gillingham Business Park Gillingham ME8 0SB U.K.

www.henryschein.com




Binocular Microscope

- Henry Schein Brand
- Siedentopf Binocular head
- Dual diopter adjustments
- 10X WF Highpoint Eyepieces
- 4X, 10X, 40X, 100X Oil Objectives
- Variable LED illumination
- 1.25 NA Abbe Condenser

Order model # 112-6914

Laboratory Microscope

- Henry Schein Brand
- Siedentopf Binocular head
- 10X WF Highpoint Eyepieces
- 4X, 10X, 40X, 100X Oil Objectives
- 1.25 NA Abbe Condenser
- Koehler Illumination
- Variable Halogen Illumination

Order model # 112-6915



