

# Electrodes

From A Name You Can Trust.




## ECG Monitoring Electrodes

Foam Electrode w/ Metal Snap  
570-1165...50/Pk

### Features & Benefits:

- Not made with natural rubber latex
- Solid gel
- Stainless steel snap stops corrosion from occurring, eliminating artifact in the tracing
- 40 x 35 mm Rectangle
- Hypoallergenic
- 24 month shelf life



## Electrode Resting Tabs

570-1173...100/Pk

### Features & Benefits:

- Conductive Adhesive Gel
- Unique liner card provides ease of removal and application
- Latex-Free
- Adhesive solid gel requires minimal clean-up
- Pliable substrate conforms to the body
- 24 month shelf life



**RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!**



## ECG Monitoring Electrodes

Foam Electrode w/ Metal Snap, Diaphoretic  
570-1166...50/Pk

### Features & Benefits:

- Not made with natural rubber latex
- Solid gel
- Hypoallergenic
- 40 x 36 mm Tear Drop
- Stainless steel snap stops corrosion from occurring, eliminating artifact in the tracing
- 24 month shelf life



## ECG Monitoring Electrodes

Radiolucent Foam w/ Carbon Snap  
570-1167...50/Pk

### Features & Benefits:

- Can be worn in X-Ray/MRI environment
- Not made with natural rubber latex
- Radiolucent Ag/AgCl sensor
- Hypoallergenic
- Solid gel
- 40 x 36 mm Tear Drop
- 24 month shelf life



## ECG Monitoring Electrodes

Clear Tape w/ Metal Snap  
570-1169...50/Pk

### Features & Benefits:

- Not made with natural rubber latex
- Hypoallergenic
- Solid gel
- 40 x 36 mm Tear Drop
- Stainless steel snap stops corrosion from occurring, eliminating artifact in the tracing
- 24 month shelf life



**RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!**