


Emesis Bags

From A Name You Can Trust.

Features & Benefits

- Eliminate exposure to emesis
- Minimize spillage
- Opening supported by plastic notched ring to provide secure grip and unobstructed access into the bag.
- Twist bag and secure into notched ring after use
- Measured up to 40 fluid ounces or 1000 CC/mL
- Compatible with following dispensers:
 - 116-0999 - Xodus Medical
 - 582-3154 - Cardinal Health



**Emesis Bags
570-1178**

Specifications

- Quantity: 24 per pack
- Dimensions: 16.5 cm x 38.1 cm
- Thickness: 1.35 mil thickness
- Material: High-density Polyethylene



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!