

High Frequency Desiccator



Item Number: 570-0583

40 watt high frequency desiccator

Less to do and less to buy makes this the most cost effective and user friendly unit on the market

Adjustable in 1/10th watt increments in the first 10 watts, in monopolar and bipolar

More accurate power with Universal Power Control 100-240 VAC

The Henry Schein Desiccator is a ground reference unit which means no grounding pad is required.



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!



STANDARD EQUIPMENT

- 570-0583 (1) Henry Schein High Frequency Generator
- 895-0124 (1) 3-Button Autoclavable Handpiece
- 112-7171 (5) non-sterile sharp dermal tips
- 112-7172 (2) sterile sharp dermal tips
- 112-7173 (5) non-sterile blunt dermal tips
- 112-7175 (2) sterile blunt dermal tips
- 895-0127 (2) non-sterile DERM-Elite™ Premium blunt dermal tips
- 895-0126 (2) sterile DERM-Elite Premium blunt dermal tips
- 112-7176 (2) sterile disposable handpiece sheath
- 570-0331 (2) non-sterile disposable handpiece sheath
- 926-3686 (1) wall mount kit
- 110 VAC Hospital-grade power cord*
- User's guide*



895-0124



112-7171, 112-7172, 112-7173,
112-7175, 895-0127, 895-0126



926-3686



112-7176



570-0331

COMPLEMENTARY PRODUCTS †

† Sold Separately



102-7058
Table Top Stand †



600-5943
Footswitch †

NEW



5700888, 5700890
Sharp H-Type Tips †

NEW



5700891, 5700889
Blunt H-Type Tips †



118-3875
Force Cable †



147-2819
Reusable Bipolar Forceps †



682-5062
Ball Electrode †



650-3112
Needle Electrode †



329-9379
Top Tray for Mobile Stand †



330-8942
Mobile Stand †

PRODUCT SPECIFICATIONS

OUTPUT POWER
TO 40 WATTS (0-10 WATTS IN
1/10TH WATT INCREMENTS
ON BOTH MONOPOLAR
AND BIPOLAR)

DIMENSIONS
HEIGHT: 188MM (7.4 IN.)
WIDTH: 228MM (8.98 IN.)
DEPTH: 105MM (4.13 IN.)
WEIGHT: < 2.26 KG (< 5 LBS.)

OUTPUT FREQUENCY
328 – 492KHZ

LINE VOLTAGE
570-0583: 100-240 VAC

LINE FREQUENCY
50 - 60 HZ

WARRANTY
FOUR YEARS



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!