

Surgical Sutures

From A Name You Can Trust.

Maxima[®] Sutures by Henry Schein[®]

- Silicone-coated 300 series stainless steel needles give finer point geometry for superior penetration
- Precision needles are strong enough to resist bending, yet ductile enough to resist breaking
- Large breadth of suture and needle combination for a range of surgeries
- Significant cost savings for your clinic compared to the leading brand


RACETRACK FEATURES AND BENEFITS

POLYLITH CARD

Helps reduce moisture content, an advantage with absorbable suture materials that degrade by hydrolysis.

LARGER DIAMETER

Larger overall diameter to ensure less suture memory and minimize the number of suture rotations around track.



NEEDLE PARK

Two needle parks to accommodate a wide range of needle types, sizes, and double-arm sutures. Needle park also protects needles from damage.




PRESS-PUSH NEEDLE LIFT

The tab, when pressed from behind, lifts needle for easy removal.







RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!




PGA Absorbable Suture

Precision Reverse Cutting Needles			PGA ABSORBABLE SUTURE								
			Material	Color	Length	6-0	5-0	4-0	3-0	2-0	0
3/8 circle 13mm	NP-3		Braided	Undyed	18" 45cm		570-1048 12/Box	570-1049 12/Box			
3/8 circle 19mm	NPS-2		Braided	Undyed	18" 45cm			570-1045 12/Box	570-1046 12/Box		
3/8 circle 19mm	NFS-2		Braided	Violet	18" 45cm			570-1399 12/Box			
				Undyed	30" 75cm		570-1397 12/Box	570-1041 36/Box	570-1047 36/Box		

Nylon Non-absorbable Suture

Precision Reverse Cutting Needles			NYLON NON-ABSORBABLE SUTURE								
			Material	Color	Length	6-0	5-0	4-0	3-0	2-0	0
3/8 circle 13mm	NP-3		Monofilament	Black	18" 45cm		570-1039 12/Box	570-1040 12/Box			
3/8 circle 19mm	NPS-2		Monofilament	Black	18" 45cm			570-1042 12/Box	570-1044 36/Box		
3/8 circle 19mm	NFS-2		Monofilament	Black	18" 45cm		570-1398 12/Box	570-1043 12/Box			
3/8 circle 24mm	NFS-1		Monofilament	Black	18" 45cm				570-1050 12/Box		

Polypropylene Non-absorbable Suture

Precision Reverse Cutting Needles			POLYPROPYLENE NON-ABSORBABLE SUTURE								
			Material	Color	Length	6-0	5-0	4-0	3-0	2-0	0
3/8 circle 13mm	NP-3		Monofilament	Blue	18" 45cm		570-1591 12/Box	570-1593 12/Box			
3/8 circle 19mm	NPS-2		Monofilament	Blue	18" 45cm			570-1592 12/Box	570-1594 36/Box		
3/8 circle 19mm	NFS-2		Monofilament	Blue	18" 45cm			570-1595 12/Box			

HSI-PUB1001 R02 Effective Date: XX/XX/XXXX



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!