

Surgical Dressings

From A Name You Can Trust.

Laparotomy Sponge

Features & Benefits:

- 18" x 18"
- Pre-washed, 100% high quality cotton with tight weave
- Cross stitched with super locked edges to minimize linting and loose threads
- X-ray detectable
- Banded 5's
- Not made with natural rubber latex



570-1115 - Laparotomy Sponge, 5/Pk

OR Towels

Features & Benefits:

- 17" x 27"
- Pre-washed, 100% high quality cotton
- Highly absorbent and virtually lint free
- All sterile towels are CSR wrapped and have a special corner fold that allows for aseptic presentation into a sterile field
- Not made with natural rubber latex



570-1113 - OR Towel Non-Sterile, 100/Ca

570-1117 - OR Towel Sterile, Blue 2s, 2/Pk

570-1116 - OR Towel Sterile, Blue 4s, 4/Pk

570-1114 - OR Towel Sterile, Blue 6s, 6/Pk

USP Type VII Gauze

Features & Benefits:

- 10/Pack
- Latex-Free
- Sterile
- Item number 570-1098 & 570-1099 are X-Ray detectable, the rest are not X-Ray detectable



570-1096 - USP Type VII Gauze, 4x4, 12 Ply

570-1097 - USP Type VII Gauze, 4x4, 8 Ply

570-1098 - USP Type VII Gauze, 4x4 16 Ply (X-Ray)

570-1099 - USP Type VII Gauze, 4x4 12 Ply (X-Ray)

570-1100 - USP Type VII Gauze, 4x4 16 Ply



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!