



Pressure-Activated Safety Lancets

From A Name You Can Trust.

Henry Schein Safety Lancets can be used by all trained healthcare professionals to draw capillary blood samples. These lancets are considered to be an integral component in preventing sharps-related injuries at all medical facilities.



1. Carefully pull off the protective cap.



2. Push the lancet firmly onto the chosen sample site to fire.



3. Dispose of the used lancet in a suitable sharps container.

Pressure-Activated Safety Lancets				
Item Number	Color	Size	Specification	Package
(900-7993)	Green	18g x 1.8mm	High Flow	100/Bx
(900-7981)	Gray	18g x 2.3mm	High Flow	100/Bx
(900-7983)	Yellow	21g x 1.0mm	Medium Flow	100/Bx
(900-7992)	Light Blue	21g x 1.8mm	Medium Flow	100/Bx
(900-7980)	Orange	21g x 2.2mm	Medium Flow	100/Bx
(900-7982)	Rose	21g x 2.8mm	High Flow	100/Bx
(570-0674)	White	23g x 1.8mm	Medium Flow	100/Bx
(570-0675)	Yellow	26g x 1.8mm	Low Flow	100/Bx
(900-7991)	Purple	28g x 1.8mm	Low Flow	100/Bx
(570-0676)	Turquoise	30g x 1.8mm	Low Flow	100/Bx

Features & Benefits:

Designed and manufactured to provide strength, protection, and value.

- Applicable to glucose and other in vitro tests.
- Needle is shielded before and after use.
- Convenient pre-loaded, pressure-activated design.
- Comfortable ultra-sharp beveled needle.
- Various sizes available for accurate sampling.



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!