

CRITERION[®] GLOVES PORTFOLIO

High-quality, value-priced exam and surgical gloves



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!

Nothing should come between your skilled hands and your patients...



...except Criterion® gloves from Henry Schein®

The health of your patients is in your hands. Protect those hands and your patients with gloves from Henry Schein's Criterion® portfolio. With a wide range of options to suit all of your medical and surgical needs, Criterion® gloves deliver the fit, texture, and durability of national brands at a competitive price. Criterion® gloves are backed by the Henry Schein brand you trust.

Choose the gloves that best suit your needs

	Latex	Vinyl	Nitrile	Poly-chloro-prene	Polyiso-prene	Powder-Free	Chemo Rated	Sterile	Exam	Surgical	Notes
Criterion® (Latex)	✓					✓		✓		✓	Excellent grip in wet or dry conditions.
Criterion® (Nitrile)			✓			✓	✓		✓		
Criterion® (Vinyl)		✓				✓			✓		Supple touch for maximum tactile sensitivity.
Criterion® CL	✓					✓			✓		Provides comfortable fit and low protein levels.
Criterion® COATS			✓			✓			✓		Features an all-natural colloidal oatmeal USP coating.
Criterion® CR (Exam)				✓		✓			✓		More natural feel and fit to hand.
Criterion® CR (Surgical)				✓		✓		✓		✓	Can be worn for extended periods without hand fatigue.
Criterion® EC			✓			✓	✓		✓		Extended cuff.
Criterion® N200			✓			✓	✓		✓		
Criterion® N250			✓			✓	✓		✓		
Criterion® N300 ULTRA			✓			✓	✓		✓		Eco-friendly packaging reduces waste and saves space.
Criterion® PC	✓					✓			✓		
Criterion® PI					✓	✓		✓		✓	Low allergenic potential due to skin-friendly accelerators.
Criterion® Pure Freedom			✓			✓	✓		✓		Reduced potential for dermatitis for glove wearers.

Criterion® (Latex) Powder-Free Latex Surgeons' Gloves



Item Code	Description
(570-0635)	Criterion® (Latex) 5½ Box, 50 pair
(570-0636)	Criterion® (Latex) 6 Box, 50 pair
(570-0637)	Criterion® (Latex) 6½ Box, 50 pair
(570-0638)	Criterion® (Latex) 7 Box, 50 pair
(570-0639)	Criterion® (Latex) 7½ Box, 50 pair
(570-0640)	Criterion® (Latex) 8 Box, 50 pair
(570-0641)	Criterion® (Latex) 8½ Box, 50 pair
(570-0642)	Criterion® (Latex) 9 Box, 50 pair

Application:	Surgical Glove
Chemo Rated:	N/A
Color:	Natural
Cuff Style:	Beaded Cuff
Glove Exterior:	Full Finger Textured
Glove Interior:	Chlorinated
Glove Mil Thickness (at Palm):	Min. 7.1 mil
Material:	Latex
Powder Content:	Powder-Free
Sterility:	Sterile
Usage:	Single Use
Latex-Free Indicator:	Made with Natural Rubber Latex



Acrylic, Vertical
(102-7201) Four boxes
(102-8344) Three boxes
(102-6468) Two boxes
(102-6467) One box

Glove Box Holders

Wire, Vertical
(112-5474) Three boxes
(112-5472) Two boxes



Wire, Horizontal
(112-5475) Three boxes
(112-5473) Two boxes
(112-5471) One box



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!

Criterion® (Nitrile) Powder-Free Nitrile Examination Gloves



Item Code	Description
(111-8533)	Criterion® (Nitrile) XS Box, 100 ct
(111-8535)	Criterion® (Nitrile) S Box, 100 ct
(111-8536)	Criterion® (Nitrile) M Box, 100 ct
(111-8537)	Criterion® (Nitrile) L Box, 100 ct
(111-8539)	Criterion® (Nitrile) XL Box, 90 ct

Application:	Exam Glove
Chemo Rated:	Chemo Rated
Color:	Blue
Cuff Style:	Beaded Cuff
Glove Exterior:	Fingertip Textured
Glove Interior:	Chlorinated
Glove Mil Thickness (at Palm):	Min. 1.97 mil
Material:	Nitrile
Powder Content:	Powder-Free
Sterility:	Non-Sterile
Usage:	Single Use
Latex-Free Indicator:	Not Made with Natural Rubber Latex

Chemotherapy Drugs and Concentration:

- Cisplatin 1.0 mg/mL (1,000 ppm)
- Cyclophosphamide (Cytoxan), 20.0 mg/mL (20,000 ppm)
- Dacarbazine (DTIC) 10.0 mg/mL (10,000 ppm)
- Doxorubicin Hydrochloride 2.0 mg/mL (2,000 ppm)
- Etoposide (Toposar) 20.0 mg/mL (20,000 ppm)
- Fluorouracil 50.0 mg/mL (50,000 ppm)
- Methotrexate 25.0 mg/mL (25,000 ppm)
- Mitomycin C 0.5 mg/mL (500 ppm)
- Paclitaxel (Taxol) 6.0 mg/mL (6,000 ppm)
- Vincristine Sulfate 1.0 mg/mL (1,000 ppm)

Minimum Breakthrough Detection Time

- > 240 min.
- > 240 min.
- > 240 min.
- > 240 min.
- > 240 min.
- > 240 min.
- > 240 min.
- > 240 min.
- > 240 min.
- > 240 min.



Acrylic, Vertical
(102-7201) Four boxes
(102-8344) Three boxes
(102-6468) Two boxes
(102-6467) One box

Glove Box Holders

Wire, Vertical
(112-5474) Three boxes
(112-5472) Two boxes

Wire, Horizontal
(112-5475) Three boxes
(112-5473) Two boxes
(112-5471) One box



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!

Criterion® (Vinyl) Powder-Free Vinyl Examination Gloves



Item Code	Description
(570-0087)	Criterion® (Vinyl) XS Box, 100 ct
(570-0088)	Criterion® (Vinyl) S Box, 100 ct
(570-0089)	Criterion® (Vinyl) M Box, 100 ct
(570-0090)	Criterion® (Vinyl) L Box, 100 ct
(570-0091)	Criterion® (Vinyl) XL Box, 90 ct

Application:	Exam Glove
Chemo Rated:	N/A
Color:	Clear
Cuff Style:	Beaded Cuff
Glove Exterior:	Smooth
Glove Interior:	N/A
Glove Mil Thickness (at Palm):	Min. 3.15 mil
Material:	Vinyl
Powder Content:	Powder-Free
Sterility:	Non-Sterile
Usage:	Single Use
Latex-Free Indicator:	Not Made with Natural Rubber Latex



Acrylic, Vertical
(102-7201) Four boxes
(102-8344) Three boxes
(102-6468) Two boxes
(102-6467) One box



Wire, Horizontal
(112-5475) Three boxes
(112-5473) Two boxes
(112-5471) One box

Glove Box Holders

Wire, Vertical
(112-5474) Three boxes
(112-5472) Two boxes



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!

Criterion[®] CL Powder-Free Latex Examination Gloves



Item Code	Description
(102-5354)	Criterion [®] CL XS Box, 100 ct
(102-8296)	Criterion [®] CL S Box, 100 ct
(102-5355)	Criterion [®] CL M Box, 100 ct
(102-7763)	Criterion [®] CL L Box, 100 ct
(102-5357)	Criterion [®] CL XL Box, 100 ct

Application:	Exam Glove
Chemo Rated:	N/A
Color:	Natural
Cuff Style:	Beaded Cuff
Glove Exterior:	Fully Textured
Glove Interior:	Chlorinated
Glove Mil Thickness (at Palm):	Min. 4.3 mil
Material:	Latex
Powder Content:	Powder-Free
Sterility:	Non-Sterile
Usage:	Single Use
Latex-Free Indicator:	Made with Natural Rubber Latex



Acrylic, Vertical
(102-7201) Four boxes
(102-8344) Three boxes
(102-6468) Two boxes
(102-6467) One box

Glove Box Holders

Wire, Vertical
(112-5474) Three boxes
(112-5472) Two boxes



Wire, Horizontal
(112-5475) Three boxes
(112-5473) Two boxes
(112-5471) One box



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!

Criterion® COATS

Powder-Free Nitrile Examination Gloves



Item Code	Description
(112-6220)	Criterion® COATS XS Box, 200 ct
(112-6221)	Criterion® COATS S Box, 200 ct
(112-6222)	Criterion® COATS M Box, 200 ct
(112-6223)	Criterion® COATS L Box, 200 ct
(112-6224)	Criterion® COATS XL Box, 200 ct

Application:	Exam Glove
Chemo Rated:	N/A
Color:	White
Cuff Style:	Beaded Cuff
Glove Exterior:	Fingertip Textured
Glove Interior:	Chlorinated
Glove Mil Thickness (at Palm):	Min. 2.0 mil
Material:	Nitrile
Powder Content:	Powder-Free
Sterility:	Non-Sterile
Usage:	Single Use
Latex-Free Indicator:	Not Made with Natural Rubber Latex



Acrylic, Vertical
(102-7201) Four boxes
(102-8344) Three boxes
(102-6468) Two boxes
(102-6467) One box

Glove Box Holders

Wire, Vertical
(112-5474) Three boxes
(112-5472) Two boxes



Wire, Horizontal
(112-5475) Three boxes
(112-5473) Two boxes
(112-5471) One box



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!

Criterion® CR (Exam)

Powder-Free Polychloroprene Examination Gloves



Item Code	Description
(112-6857)	Criterion® CR (Exam) XS Box, 100 ct
(112-6858)	Criterion® CR (Exam) S Box, 100 ct
(112-6859)	Criterion® CR (Exam) M Box, 100 ct
(112-6860)	Criterion® CR (Exam) L Box, 100 ct
(112-6861)	Criterion® CR (Exam) XL Box, 90 ct

Application:	Exam Glove
Chemo Rated:	N/A
Color:	Green
Cuff Style:	Beaded Cuff
Glove Exterior:	Fingertip Textured
Glove Interior:	Polymer Coated
Glove Mil Thickness (at Palm):	Min. 3.15 mil
Material:	Polychloroprene
Powder Content:	Powder-Free
Sterility:	Non-Sterile
Usage:	Single Use
Latex-Free Indicator:	Not Made with Natural Rubber Latex



Acrylic, Vertical
(102-7201) Four boxes
(102-8344) Three boxes
(102-6468) Two boxes
(102-6467) One box

Glove Box Holders

Wire, Vertical
(112-5474) Three boxes
(112-5472) Two boxes



Wire, Horizontal
(112-5475) Three boxes
(112-5473) Two boxes
(112-5471) One box



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!

Criterion® CR (Surgical) Powder-Free Polychloroprene Surgical Gloves



Item Code	Description
(112-7079)	Criterion® CR (Surgical) 5½ Box, 50 pair
(112-7080)	Criterion® CR (Surgical) 6 Box, 50 pair
(112-7081)	Criterion® CR (Surgical) 6½ Box, 50 pair
(112-7082)	Criterion® CR (Surgical) 7 Box, 50 pair
(112-7083)	Criterion® CR (Surgical) 7½ Box, 50 pair
(112-7085)	Criterion® CR (Surgical) 8 Box, 50 pair
(112-7086)	Criterion® CR (Surgical) 8½ Box, 50 pair
(112-7088)	Criterion® CR (Surgical) 9 Box, 50 pair

Application:	Surgical Glove
Chemo Rated:	N/A
Color:	White
Cuff Style:	Beaded Cuff
Glove Exterior:	Fully Textured
Glove Interior:	Polymer Coated
Glove Mil Thickness (at Palm):	Min. 5.9 mil
Material:	Polychloroprene
Powder Content:	Powder-Free
Sterility:	Sterile
Usage:	Single Use
Latex-Free Indicator:	Not Made with Natural Rubber Latex



Acrylic, Vertical
(102-7201) Four boxes
(102-8344) Three boxes
(102-6468) Two boxes
(102-6467) One box

Glove Box Holders

Wire, Vertical
(112-5474) Three boxes
(112-5472) Two boxes



Wire, Horizontal
(112-5475) Three boxes
(112-5473) Two boxes
(112-5471) One box



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!

Criterion® EC Powder-Free Nitrile Examination Gloves



Item Code	Description
(112-6785)	Criterion® EC XS Box, 50 ct
(112-6786)	Criterion® EC S Box, 50 ct
(112-6787)	Criterion® EC M Box, 50 ct
(112-6788)	Criterion® EC L Box, 50 ct
(112-6790)	Criterion® EC XL Box, 50 ct

Application:	Exam Glove
Chemo Rated:	Chemo Rated
Color:	Aqua Blue
Cuff Style:	Beaded Cuff
Glove Exterior:	Fingertip Textured
Glove Interior:	Chlorinated
Glove Mil Thickness (at Palm):	Min. 2.8 mil
Material:	Nitrile
Powder Content:	Powder-Free
Sterility:	Non-Sterile
Usage:	Single Use
Latex-Free Indicator:	Not Made with Natural Rubber Latex

Chemotherapy Drugs and Concentration

- Cisplatin 1.0 mg/mL (1,000 ppm)
- Cyclophosphamide (Cytoxan), 20.0 mg/mL (20,000 ppm)
- Dacarbazine (DTIC) 10.0 mg/mL (10,000 ppm)
- Doxorubicin Hydrochloride 2.0 mg/mL (2,000 ppm)
- Etoposide 20.0 mg/mL (20,000 ppm)
- Fluorouracil 50.0 mg/mL (50,000 ppm)
- Methotrexate 25.0 mg/mL (25,000 ppm)
- Mitomycin C 0.5 mg/mL (500 ppm)
- Paclitaxel (Taxol) 6.0 mg/mL (6,000 ppm)
- Vincristine Sulfate 1.0 mg/mL (1,000 ppm)

Minimum Breakthrough Detection Time

- > 240 min.
- > 240 min.
- > 240 min.
- > 240 min.
- > 240 min.
- > 240 min.
- > 240 min.
- > 240 min.
- > 240 min.
- > 240 min.



- Acrylic, Vertical**
- (102-7201) Four boxes
 - (102-8344) Three boxes
 - (102-6468) Two boxes
 - (102-6467) One box

Glove Box Holders

- Wire, Vertical**
- (112-5474) Three boxes
 - (112-5472) Two boxes

- Wire, Horizontal**
- (112-5475) Three boxes
 - (112-5473) Two boxes
 - (112-5471) One box



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!

Criterion N200

Powder-Free Nitrile Examination Gloves



Application:	Exam Glove
Chemo Rated:	Chemo Rated
Color:	Blue
Cuff Style:	Beaded Cuff
Glove Exterior:	Fingertip Textured
Glove Interior:	Chlorinated
Glove Mil Thickness (at Palm):	Min. 1.97 mil
Material:	Nitrile
Powder Content:	Powder-Free
Sterility:	Non-Sterile
Usage:	Single Use
Latex Free Indicator:	Not Made with Natural Rubber Latex

Item Code	Description
(900-7437)	Criterion N200 XS Box, 200 ct
(900-7438)	Criterion N200 S Box, 200 ct
(900-7439)	Criterion N200 M Box, 200 ct
(900-7440)	Criterion N200 L Box, 200 ct
(900-7441)	Criterion N200 XL Box, 180 ct

Chemotherapy Drugs and Concentration

- Cisplatin 1.0 mg/mL (1,000 ppm)
- Cyclophosphamide (Cytoxan), 20.0 mg/mL (20,000 ppm)
- Dacarbazine (DTIC) 10.0 mg/mL (10,000 ppm)
- Doxorubicin Hydrochloride 2.0 mg/mL (2,000 ppm)
- Etoposide (Toposar) 20.0 mg/mL (20,000 ppm)
- Fluorouracil 50.0 mg/mL (50,000 ppm)
- Methotrexate 25.0 mg/mL (25,000 ppm)
- Mitomycin C 0.5 mg/mL (500 ppm)
- Paclitaxel (Taxol) 6.0 mg/mL (6,000 ppm)
- Vincristine Sulfate 1.0 mg/mL (1,000 ppm)

Minimum Breakthrough Detection Time

- > 240 min.
- > 240 min.
- > 240 min.
- > 240 min.
- > 240 min.
- > 240 min.
- > 240 min.
- > 240 min.
- > 240 min.
- > 240 min.



- Acrylic, Horizontal**
- (102-7201)** Four boxes
 - (102-8344)** Three boxes
 - (102-6468)** Two boxes
 - (102-6467)** One box

Glove Box Holders

- Wire, Vertical**
- (112-5474)** Three boxes
 - (112-5472)** Two boxes



- Wire, Horizontal**
- (112-5475)** Three boxes
 - (112-5473)** Two boxes
 - (112-5471)** One box



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!

Criterion® N250 Powder-Free Nitrile Examination Gloves



Item Code	Description
(570-0099)	Criterion® N250 XS Box, 250 ct
(570-0100)	Criterion® N250 S Box, 250 ct
(570-0101)	Criterion® N250 M Box, 250 ct
(570-0102)	Criterion® N250 L Box, 250 ct
(570-0103)	Criterion® N250 XL Box, 250 ct

Application:	Exam Glove
Chemo Rated:	Chemo Rated
Color:	Cornflower Blue
Cuff Style:	Beaded Cuff
Glove Exterior:	Fingertip Textured
Glove Interior:	Chlorinated
Glove Mil Thickness (at Palm):	Min. 1.97 mil
Material:	Nitrile
Powder Content:	Powder-Free
Sterility:	Non-Sterile
Usage:	Single Use
Latex-Free Indicator:	Not Made with Natural Rubber Latex

Chemotherapy Drugs and Concentration

- Amethopterin Hydrate 25.0 mg/mL (25,000 ppm)
- Carmustine (BCNU) 3.3 mg/mL (3,300 ppm)
- Cisplatin 1.0 mg/mL (1,000 ppm)
- Cyclophosphamide (Cytoxan) 20.0 mg/mL (20,000 ppm)
- Dacarbazine (DTIC) 10.0 mg/mL (10,000 ppm)
- Doxorubicin HCl (Adriamycin) 2.0 mg/mL (2,000 ppm)
- Etoposide (Toposar) 20.0 mg/mL (20,000 ppm)
- 5-Fluorouracil (Adrucil) 50.0 mg/mL (50,000 ppm)
- Methotrexate 25.0 mg/mL (25,000 ppm)
- Mitomycin C 0.5 mg/mL (500 ppm)
- Paclitaxel (Taxol) 6.0 mg/mL (6,000 ppm)
- Thiotepa 10.0 mg/mL (10,000 ppm)
- Vincristine Sulfate 1.0 mg/mL (1,000 ppm)

Minimum Breakthrough Detection Time

- > 240 min.
- > 10.1 min.
- > 240 min.
- > 240 min.
- > 240 min.
- > 240 min.
- > 240 min.
- > 240 min.
- > 240 min.
- > 240 min.
- > 240 min.
- > 20.2 min.
- > 240 min.



- Acrylic, Vertical**
- (102-7201) Four boxes
 - (102-8344) Three boxes
 - (102-6468) Two boxes
 - (102-6467) One box

Glove Box Holders

- Wire, Vertical**
- (112-5474) Three boxes
 - (112-5472) Two boxes

- Wire, Horizontal**
- (112-5475) Three boxes
 - (112-5473) Two boxes
 - (112-5471) One box



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!

Criterion® N300 ULTRA Powder-Free Nitrile Examination Gloves



Item Code	Description
(112-6829)	Criterion® N300 ULTRA XS Box, 300 ct
(112-6830)	Criterion® N300 ULTRA S Box, 300 ct
(112-6831)	Criterion® N300 ULTRA M Box, 300 ct
(112-6832)	Criterion® N300 ULTRA L Box, 300 ct
(112-6833)	Criterion® N300 ULTRA XL Box, 250 ct

Application:	Exam Glove
Chemo Rated:	Chemo Rated
Color:	Tiffany Blue
Cuff Style:	Beaded Cuff
Glove Exterior:	Fingertip Textured
Glove Interior:	Chlorinated
Glove Mil Thickness (at Palm):	Min. 1.97 mil
Material:	Nitrile
Powder Content:	Powder-Free
Sterility:	Non-Sterile
Usage:	Single Use
Latex-Free Indicator:	Not Made with Natural Rubber Latex

Chemotherapy Drugs and Concentration

- Cisplatin 1.0 mg/mL (1,000 ppm)
- Cyclophosphamide (Cytoxan), 20.0 mg/mL (20,000 ppm)
- Dacarbazine (DTIC) 10.0 mg/mL (10,000 ppm)
- Doxorubicin Hydrochloride 2.0 mg/mL (2,000 ppm)
- Etoposide 20.0 mg/mL (20,000 ppm)
- Fluorouracil 50.0 mg/mL (50,000 ppm)
- Methotrexate 25.0 mg/mL (25,000 ppm)
- Mitomycin C 0.5 mg/mL (500 ppm)
- Paclitaxel (Taxol) 6.0 mg/mL (6,000 ppm)
- Vincristine Sulfate 1.0 mg/mL (1,000 ppm)

Minimum Breakthrough Detection Time

- > 240 min.
- > 240 min.
- > 240 min.
- > 240 min.
- > 240 min.
- > 240 min.
- > 240 min.
- > 240 min.
- > 240 min.
- > 240 min.



Acrylic, Vertical
(102-7201) Four boxes
(102-8344) Three boxes
(102-6468) Two boxes
(102-6467) One box

Glove Box Holders

Wire, Vertical
(112-5474) Three boxes
(112-5472) Two boxes



Wire, Horizontal
(112-5475) Three boxes
(112-5473) Two boxes
(112-5471) One box



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!

Criterion® PC Powder-Free Latex Examination Gloves



Item Code	Description
(102-5418)	Criterion® PC XS Box, 100 ct
(102-5419)	Criterion® PC S Box, 100 ct
(102-5421)	Criterion® PC M Box, 100 ct
(102-6730)	Criterion® PC L Box, 100 ct
(102-5422)	Criterion® PC XL Box, 90 ct

Application:	Exam Glove
Chemo Rated:	N/A
Color:	White
Cuff Style:	Beaded Cuff
Glove Exterior:	Fully Textured
Glove Interior:	Polymer Coated
Glove Mil Thickness (at Palm):	Min. 3.15 mil
Material:	Latex
Powder Content:	Powder-Free
Sterility:	Non-Sterile
Usage:	Single Use
Latex-Free Indicator:	Made with Natural Rubber Latex



Acrylic, Vertical
(102-7201) Four boxes
(102-8344) Three boxes
(102-6468) Two boxes
(102-6467) One box

Glove Box Holders

Wire, Vertical
(112-5474) Three boxes
(112-5472) Two boxes



Wire, Horizontal
(112-5475) Three boxes
(112-5473) Two boxes
(112-5471) One box



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!

Criterion® PI

Polyisoprene (PI) Surgical Glove



Application:	Surgical Glove
Chemo Rated:	N/A
Color:	Cream
Cuff Style:	Extended Cuff
Glove Exterior:	Fully Microtextured
Glove Interior:	Polymer Coated
Glove Mil Thickness (at Palm):	Min. 7.6 mil
Material:	Polyisoprene
Powder Content:	Powder-Free
Sterility:	Sterile
Usage:	Single Use
Latex-Free Indicator:	Not Made with Natural Rubber Latex

Item Code	Description
(570-1119)	Criterion® PI Surgical Glove, 5.5 Box, 50 pair
(570-1120)	Criterion® PI Surgical Glove, 6.0 Box, 50 pair
(570-1121)	Criterion® PI Surgical Glove, 6.5 Box, 50 pair
(570-1122)	Criterion® PI Surgical Glove, 7.0 Box, 50 pair
(570-1123)	Criterion® PI Surgical Glove, 7.5 Box, 50 pair
(570-1124)	Criterion® PI Surgical Glove, 8.0 Box, 50 pair
(570-1125)	Criterion® PI Surgical Glove, 8.5 Box, 50 pair
(570-1126)	Criterion® PI Surgical Glove, 9.0 Box, 50 pair



Acrylic, Vertical
(102-7201) Four boxes
(102-8344) Three boxes
(102-6468) Two boxes
(102-6467) One box

Glove Box Holders

Wire, Vertical
(112-5474) Three boxes
(112-5472) Two boxes



Wire, Horizontal
(112-5475) Three boxes
(112-5473) Two boxes
(112-5471) One box



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!

