

# Hinged Knee Brace

## From A Name You Can Trust.

*Our Knee Brace Hinged contain all the features necessary to ensure your patients heal with a proper alignment, for a quick and complete recovery, while having the best patient experience available.*



Hinged Knee Brace			
Item Number	Brace Size	Thigh Width	Calf Width
(570-1029)	X-Small	16.5"	14"
(570-1030)	Small	17.5"	15"
(570-1031)	Medium	18.5"	16"
(570-1032)	Large	19.5"	17"
(570-1033)	X-Large	21"	18"
(570-1034)	2X-Large	23"	20"
(570-1035)	3X-Large	25"	22"

### Features & Benefits:

*Designed and manufactured to provide strength, protection, and value.*

- Removable dual axis polycentric hinges mimic the true anatomy of the knee and provide both medial and lateral support
- Front closure provides easy application
- Condyle pads for additional comfort
- Wrap around configuration and strap provide warmth and equal compression
- Brushed neoprene and seamed edges increase patient comfort and ease of washing



**RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!**