

Laboratory Equipment From A Name You Can Trust.



PowerSpin™ Centrifuge Variable
(112-6907) 6-Place
(112-6909) 8-Place

- Variable speed centrifuge with color-coded panel for blood, urine, and fecal specimens.
- Variable speed 500–4,000 RPM
 - 6- or 8-place fixed angle rotor
 - 30-minute timer
 - Stainless steel rotor chamber
 - Lid latch safety switch



PowerSpin™ Centrifuge Fixed Speed
(112-6905)

- Fixed speed at 3,400 RPM
- 6-place fixed angle rotor
- 30-minute timer
- Stainless steel rotor chamber
- Lid latch safety switch



PortaSpin™ Mobile Centrifuge
(112-6906)

- Plugs into a car 12V power outlet
- Fixed speed at 3,500 RPM
- 6-place fixed angle rotor
- 30-minute timer
- Stainless steel carrying handle
- Optional AC power adapter



Microhematocrit Centrifuge
(112-6910)

- 24-place rotor for 75mm capillary tubes
- Fixed speed at 12,000 RPM
- 5-minute timer
- All-metal housing
- Safety switch
- Microhematocrit reader included



PowerSpin™ DXB Centrifuge
(570-0483)

- Variable speed 500–5,500 RPM
- 6-place horizontal rotor
- Digital timer with recall
- Lid latch safety switch
- Compatible with BD Vacutainer™ Barricor™ tubes



Phlebotomy Tray
(112-6916)

- Large carrying capacity
- Includes racks (sharps and alcohol dispenser not included)
- Handy pull out drawer



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!



Laboratory Microscopes
(112-6914) LED Illumination
(112-6915) Halogen Illumination

- Binocular head
- 10x WF Highpoint Eyepieces
- 4X, 10X, 40X, 100X Oil Achromat Objective
- 1.25 NA Abbe Condenser



Advanced Microscope
(570-0671)

- Binocular head
- 10x WF Highpoint Eyepieces
- 4X, 10X, 40X, 100X Oil Plan Objectives
- LED Illumination



Clinical Microscopes
(570-0437) Binocular Head
(570-0438) Trinocular Head

- 10X WF Highpoint Eyepieces
- 2X, 4X, 10X, 40X, 100X Infinity Plan Objectives
- Flip-out-Abbe Condenser
- Halogen Illumination



Blood Tube Rocker
(112-6911)

- 15-place platform for up to 120mm tubes
- 25 RPM fixed speed
- Compact size for minimal footprint



Laboratory Rotator
(112-6912)

- 11" x 11" platform
- Variable speed 25-250 RPM
- 30-minute timer
- Humidity cover included



Vortex Tube Mixer
(112-6913)

- Variable speed up to 3,000 RPM
- Switch from Auto to Touch mode
- Power "on" LED indicator



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!