

Put your practice **ONE STEP** ahead
of today's frantic pace with

 HENRY SCHEIN®

One Step+ hCG Pregnancy Tests



Count on it for...

EASE OF USE

- Just 3 drops of sample
- CLIA-waived (urine)
- Test with urine or serum

SPEED

- True 1-step urine test
- Positive result in as little as 1 minute

RELIABILITY

- Proven patented technology
- Built-in procedural controls for accuracy



Item # 900-4072 OneStep+ hCG Urine Strip Test - 25 Strips
Item # 900-4073 OneStep+ hCG Urine Cassette Test - 25 Cassettes
Item # 900-4074 OneStep+ hCG Urine Cassette Test - 100 Cassettes
Item # 900-4075 OneStep+ hCG Combo Test - 25 Cassettes
Contact Your Sales Representative for More Information.
www.henryschein.com

 HENRY SCHEIN® BRAND
Quality Products From a Company You Trust!

hCG Pregnancy Tests



hCG COMBO TEST COMPARISON CHART

Product Name	Henry Schein® One Step+ hCG Combo Test	Quidel QuickVue® Combo	Genzyme OSOM® hCG Combo
Format	25 Cassettes	50 Cassettes	25 Cassettes
CLIA Status	Waived: Urine Moderate: Serum	Waived: Urine Moderate: Serum	Waived: Urine Moderate: Serum
Sample	Urine or serum	Urine or serum	Urine or serum
Sensitivity	25 mIU/mL	25 mIU/mL	Urine: 20 mIU/mL Serum: 10 mIU/mL
Specificity	>99%	>99%	>99%**
Procedural Steps	1	1	1
Time to Result	(+) = 1 minute* (-) Urine = 3 minutes (-) Serum = 5 minutes	Urine = 3 minutes Serum = 5 minutes	Urine = 3 minutes Serum = 5 minutes
Shelf Life	24 months	18 months	24 months

*Positive results may appear in as little as 1 minute. **Agreement on Negative Samples



hCG PREGNANCY TEST PROCEDURES

hCG Cassette



Read result in 3 minutes for urine or 5 minutes for serum.
A positive result may be observed in as little as 1 minute.

hCG Strip



Read result in 3 minutes.
A positive result may be observed in as little as 1 minute.

Item # 900-4072 OneStep+ hCG Urine Strip Test - 25 Strips
 Item # 900-4073 OneStep+ hCG Urine Cassette Test - 25 Cassettes
 Item # 900-4074 OneStep+ hCG Urine Cassette Test - 100 Cassettes
 Item # 900-4075 OneStep+ hCG Combo Test - 25 Cassettes
 Contact Your Sales Representative for More Information.
www.henryschein.com

HENRY SCHEIN® BRAND
 Quality Products From a Company You Trust!