




Topical Skin Adhesive



Flexible Skin Adhesive with Applicator Tip

Blue
(570-0351).....ea.
(112-7158)..... 10/box.

- Beveled Tip allows for a smoother, broader layer of skin adhesive
- The flexible skin adhesive provides BOTH the tensile strength of a Butyl formula and the flexibility of an Octyl formula
- Ampule can be used with or without the tip

0.5mL Ampule Topical Skin Adhesive

Blue
(112-6111).....ea.

Clear
(112-6106).....ea.

All versions of our Topical Skin Adhesives are

- FDA-approved liquid adhesive consisting of n-butyl-2 cyanoacrylate
- Highest bond strength of all topical tissue adhesives on the market
- Inhibits growth of bacteria
- Sterile, single-patient-use plastic ampules
- Provide a precise, controlled application
- Will not clog or dry when opened



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!