

SBA 258

Product: Henry Schein Waterline Maintenance Tablets		
CAS No.	Name	Function
77-92-9	Citric Acid	pH Adjuster
151-21-3	Sodium Lauryl Sulfate	Surfactant
144-558	Sodium Bicarbonate	Builder
9003-39-8	PVP	Binder
7440-22-4	Silver	Active Ingredient
15792-67-3	Acid Blue Aluminum Salt	Colorant

LINKS TO DESIGNATED LISTS

1. California Proposition 65	2. EU CMRs
3. EU Endocrine Disruptors	4. EPA IRIS Neurotoxicants
5. EPA IRIS Carcinogens	6. EU PBTs and vPvBs
7. Canadian DSL PBiTs	8. EU Respiratory Sensitizers
9. IARC Carcinogens	10. ATSDR Neurotoxicants
11. EPA NWMP PBTs	12. NTP OHAT Reproductive or Developmental Toxicants
13. EPA TRI PBTs	14. Washington PBTs
15. NTP Carcinogens	16. California NLs
17. California MCLs	18. California TACs
19. California 303(c) and (d) Priority Pollutants	20. OEHHA RELs
21. California Priority Chemicals	22. OSPAR Priority Action Part A
23. Non-Functional Constituents	24. EU Fragrance Allergens