



MOLECULAR. IN MINUTES.™

ID NOW™ INFLUENZA A & B 2


The fastest molecular flu test
just got faster and
even better.

- Now with early detection — positive results in as few as 5 minutes¹
- Detects up to 42% more true positives than RADTs²
- Room temperature storage
- CLIA waived




ID NOW™ Influenza A & B 2

SPEED – THE FASTEST MOLECULAR INFLUENZA A & B TEST FOR EXPEDITING WORKFLOW



- Positive results detected in as few as 5 minutes, with 95% detected within 7 minutes¹
- Negative results in 13 minutes²

PERFORMANCE – SUPERIOR SENSITIVITY ENHANCES CLINICAL CONFIDENCE



- Detects up to 42% more true positives than RADTs²
- Reduces unnecessary antibiotic usage

ID NOW INFLUENZA A & B 2 VS. PCR³

Direct Swab	Sensitivity	Specificity
Flu A	96.3%	97.4%
Flu B	100%	97.1%

VTM	Sensitivity	Specificity
Flu A	92.8%	98.5%
Flu B	100%	97.7%

EASE OF USE - SIMPLIFYING THE TEST PROCESS

- Room temperature storage eliminates need for refrigeration
- CLIA waived nasal or nasopharyngeal swabs, for both Direct and VTM
- Integrated pipette eliminates manual pipetting for direct samples
- Standardized respiratory trio (Influenza A & B, RSV, Strep A) across all sites of care
- Bi-directional connectivity to industry leading RALS[®] platform facilitates satellite testing and remote management

ORDER INFORMATION

PRODUCT NAME	PRODUCT CODE
ID NOW™ INFLUENZA A & B 2 - 24 TEST KIT	427-000
ID NOW™ INFLUENZA A & B - CONTROL SWAB KIT	425-080

CPT® CODE: 87502-QW

FOR MORE INFORMATION, CONTACT YOUR LOCAL ABBOTT REPRESENTATIVE OR VISIT ALERE.COM.

1. Moore N, et al. Evaluation of the ID NOW™ Influenza A & B 2 Assay. 2018 ASM Clinical Virology Symposium, West Palm Beach, FL. Poster.
 2. Chartrand et al, Accuracy of Rapid Influenza Diagnostic Tests A Meta-analysis; Ann-Intern Med. 2012; 156:500-511
 3. ID NOW Influenza A & B 2 Product Insert
 4. Veritor System for Rapid Detection of Flu A + B Product Insert, Table 2 (U.S. Sites)

VTM, viral transport medium

© 2019 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners. Any photos displayed are for illustrative purposes only. Any person depicted in such photos is a model. IDDM187625 120004452-04 01/19

