

**CRITERION™ CHLOROPRENE
EXAM GLOVES**

**EXPERIENCE
THE SOFTEST
CRITERION™
CHLOROPRENE
EXAM GLOVE.**



For more information about the full line of Criterion Personal Protective Equipment,
scan here or visit HENRYSCHEINDENTAL.COM/CRITERION

CRITERION™ CHLOROPRENE EXAM GLOVES

- Criterion's softest chloroprene exam glove
- Chloroprene material provides the comfort and fit of latex but with the synthetic strength of nitrile
- Precision performance provides dexterity and extended-wear comfort
- Puncture and tear resistance for high-risk procedures
- Microtextured fingertips for a secure grip and improved control during procedures
- Interior polymer coating allows for easy donning when hands are wet, damp, or dry
- Beaded cuff helps prevent rollback. Cuff length at 240mm (9.5")
- 3.54 mil/ 0.09 mm thickness at palm
- Ambidextrous
- Conforms to ASTM D3577-19 and ISO 10282-2014



AQUA BLUE COLOR

SIZE	QUANTITY	ITEM CODE
XS	100/box	570-5400
S	100/box	570-5401
M	100/box	570-5402
L	100/box	570-5403
XL	100/box	570-5404



CRITERION™ provides a smart selection of PPE products that let you outfit your practice from head-to-toe with quality you can rely on.

Order today at HENRYSCHEINDENTAL.COM/CRITERION