

Ordering Information

Part Number Description

Triton® 6M

4769M1XUS Triton® 6M traction table

Now chose your options:

(X) Upholstery Color Options

Blue
(x = 02)Black
(x = 08)

When ordering, replace 'x' in Part Number with desired 2-digit colour code above. 4769M108 = Triton 6M, **black**
NOTICE: Color charts are produced within the limitations of lithographic color printing and are intended as a general guide only.

Technical Specifications

Mains Power: 120 V~ 50/60 Hz

Current: Max 2A with 120V

Dimensions:	Total*	90.5 x 37 x 31 inches
	Head / back section	28 inches
	Middle section	8.6 inches
	Pelvic section	10 inches
	Leg section	30 inches
	Width	27.5 inches

Adjustment:	Height	19 inches
		37 inches
	Angle Inclination	Head Section -15° to +40°
		Pelvic Section 0° to 15°
Displacement	Leg Section 0° to 45°	
	Leg Section 4 inches	

Dynamic lifting capacity: 440 LBS

Maximum workload: 605 LBS

Patient Weight: 363 LBS

Electrical Safety Class: Class I, Type B

Safety Tests: CE, Class I Rule 12

Table Weight: 428 LBS

Average Shipping Weight: 484 LBS

Shipment Dimensions: 92.5 x 37.4 x 31.5 inches

* Dimensions are listed as L x W x H

Triton® 6M
Traction Table

Features

Comprised of 6 sections, the Triton 6M traction table uses a combination of gas springs and electric power to aid manual adjustment, allowing for easy control of table positioning. Each section is adjustable, providing a wealth of achievable patient positioning options. Together with a Triton Traction Unit (see pages 9-10), the Triton 6M is a solid choice for the basis of any therapist's practice.

- 6-section treatment table provides maximum patient comfort in prone and supine positions to aid manual and traction therapy
- Pelvic tilt section can be adjusted from 0° to 15° to aid flexion of lumbar spine and to help address pelvic asymmetry
- Head piece can be repositioned from -15° to +40° to allow easy access to thoracic and cervical spine
- Tuck-away grab bars for hanging traction without thoracic restraint
- Leg section is adjustable from 0° to 45°
- Gas springs assist in friction-free movement of head, pelvic, chest and leg sections which are easily accessible from either side of the table
- Foot/hand controller allows easy adjustment of table height and includes a lockable safety mechanism
- Complete range of accessories to suit every patient need

