

# Maxiwipe® LA Disinfecting Towelettes

## From A Name You Can Trust.

*Your first line of defense for surface disinfection*



### Features & Benefits:

- Low-alcohol disinfectant towelettes
- Clean and disinfect hard, nonporous surfaces and patient-care equipment
- Contains 17.2% alcohol
- For use as a virucide, including HIV-1, HBV and HCV applications
- Also effective against methicillin-resistant *Staphylococcus aureus* (MRSA), vancomycin-resistant *Enterococcus faecalis* (VRE), and *Staphylococcus aureus* with reduced susceptibility to vancomycin
- Use one towelette to completely preclean surface of all gross debris
- Use a second towelette to thoroughly wet the surface
- Surface should remain visibly wet for 2 minutes at room temperature (68°F/20°C) for efficiency
- Repeated use of the product may be required to ensure the surface remains visibly wet for 3 minutes at room temperature for tuberculocidal activity and effectiveness against *S. aureus*, *Pseudomonas aeruginosa*, *Trichophyton mentagrophytes*, and *Salmonella enterica*

Large-6" x 6-3/4"  
**(112-6843)**.....160 per Canister

X-Large-9" x 12"  
**(112-6844)**.....65 per Canister



**RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!**