



VersaSure™ Cleaner Disinfectant Wipes



PRODUCT SOLUTIONS



Kills viruses* and bacteria in 2 minutes



Versatility to use everywhere with the **assurance** of broad spectrum disinfection

Efficacy

- Kills bacteria and fungus, including MRSA and TB, in 2 min
- Kills viruses, including mumps and measles, in 30 sec†

Aesthetics

- Alcohol-free
- Low odor
- Low residue

Compatibility

- Good compatibility on surfaces commonly found in the healthcare setting

	VersaSure™	PDI® Sani-Cloth® AF3	PDI® Super Sani-Cloth®	CaviWipes
Bacteria and virus* kill times	30 sec—2 min	3 min	2 min	3 min
Total kill claims	44	44	29	14
Toxicity Signal Word	Caution	Caution	Warning	Caution

*See virus claims table on product
†Rotavirus is 2 minutes

Clorox Healthcare® VersaSure™ Cleaner Disinfectant Wipes


Stays wet for contact time compared to other quat/quat-alcohol wipes



Still wet after 2 minutes

VersaSure™

Contact time 2 min



Almost dry after 2 minutes

PDI® Sani-Cloth®
AF3 Wipe

Contact time 3 min



Almost dry after 90 seconds

PDI® Super-Sani Cloth®
Wipes (Purple)

Contact time 2 min



Almost dry after 90 seconds

CaviWipe

Contact time 3 min

Features patented technology that enhances quat activity on surface to provide faster kill times without alcohol


Designed for Use on Common Healthcare Surfaces

- ✓ Common plastics
- ✓ Stainless steel
- ✓ Glass
- ✓ High-end finishes (Corian™, marble)


Ideal for Use in Patient Areas

- ✓ Little residue
- ✓ Excellent cleaning power
- ✓ Stays wet for full 2 minutes
- ✓ Minimal odor

Low residue profile compared to other quat/quat-alcohol wipes after 30 wipes




VersaSure™



PDI® Sani-Cloth®
AF3 Wipe




PDI® Super-Sani Cloth®
Wipes (Purple)
(PVC)



CaviWipe


Best in class cleaning power (after 10 scrubs)




VersaSure™




PDI® Sani-Cloth®
AF3 Wipe



PDI® Super-Sani
Cloth® Wipes (Purple)



CaviWipe



Control:
dried porcine blood

VersaSure™ Cleaner Disinfectant Wipes

FOR HARD SURFACES



Multipurpose Wipes

Designed for: Equipment and hard-surface disinfecting

Item No. 31757, 85 ct. canister, 6/case (6.75" x 8")

Item No. 31758, 150 ct. canister, 6/case (6" x 5")

Item No. 31760, 30 ct. pouch, 24/case (6.75" x 8")

Terminal Wipes

Designed for: Everyday cleaning and disinfecting, daily patient room cleaning

Item No. 31759, 110 ct. bucket, 2/case (12" x 12")

Item No. 31761, 110 ct. refill pouch, 2/case (12" x 12")



For more information, contact your Clorox sales representative or call **800-234-7700**.
Email: healthcare@clorox.com
Visit us: www.cloroxhealthcare.com

© 2018 Clorox Professional Products Company, 1221 Broadway, Oakland, CA 94612