



POINT OF CARE



LAB-ACCURATE TESTING WITHOUT THE WAITING.



piccolo *xpress*
chemistry analyzer

TRANSFORM YOUR EFFICIENCY WHILE MAINTAINING QUALITY AND COMPLIANCE.

- **Test results in minutes**—Clinicians have results during a patient visit, accelerating treatment decisions and avoiding call-backs and follow-up visits
- **Accuracy you can trust**—Lab-accurate results enable clinicians to make decisions with confidence
- **Maintain control and compliance**—Automated quality control checks on every test
- **Easy to use**—Fully automated system makes testing fast and simple for staff
- **Flexible connectivity**—Capable of bi-directional interfacing



ACCELERATE TESTING. ASSURE COMPLIANCE.

- Accuracy and precision comparable to larger laboratory analyzers
- iQC® system performs continuous quality control checks on every run, helping maintain compliance while eliminating daily control testing requirements
- Fully automated system—minimal testing experience required
- Small blood sample—just 100 µL of whole blood, serum, or plasma

A 2014 study in the American Journal of Clinical Pathology¹ found that point-of-care testing significantly improved clinical operations, resulting in:

- 21% reduction in tests ordered
- 89% reduction in follow-up phone calls
- 85% reduction in follow-up letters
- 61% reduction in follow-up visits
- Practice efficiency cost savings totaling \$24.64/patient²

FROM SAMPLE TO RESULT IN 3 SIMPLE STEPS



Step 1: Add sample



Step 2: Insert rotor



Step 3: Read results

MONITOR A RANGE OF PARAMETERS WITH ONE, EASY-TO-USE PLATFORM



Comprehensive Metabolic Panel

ALB
ALP
ALT
AST
BUN
Ca
Cl⁻
CRE

GLU
K⁺
Na⁺
TBIL
tCO₂
TP
eGFR*



07P02-08



Basic Metabolic Panel

BUN
Ca
Cl⁻
CRE
GLU

K⁺
Na⁺
tCO₂
eGFR*



07P02-04



Moderately Complex Panels



Hepatic Function Panel

ALB
ALP
ALT
AST

DBIL
TBIL
TP



07P02-06

Basic Metabolic Panel Plus

BUN
Ca
Cl⁻
CRE
GLU

K⁺
Na⁺
tCO₂
Mg
eGFR*



07P02-13

Lipid Panel

CHOL
CHOL/HDL*
nHDLc*
HDL
LDL*
TRIG
VLDL*



07P02-05



Lipid Panel Plus

ALT
AST
CHOL
HDL*

nHDLc*
GLU
HDL
LDL*
TRIG
VLDL*



07P02-12



Liver Panel Plus

ALB
ALP
ALT
AMY

AST
GGT
TBIL
TP



07P02-10



MetLyte Plus CRP

BUN
CK
Cl⁻
CRE
GLU

K⁺
Na⁺
tCO₂
CRP
eGFR*



07P02-14

Biochemistry Panel Plus

ALB
ALP
ALT
AMY
AST
BUN
Ca

CRE
CRP
GGT
GLU
TP
UA
eGFR*



07P02-15

MetLac 12 Panel

BUN
Ca
Cl⁻
CRE
GLU

LAC
K⁺
Mg
Na⁺
PHOS



07P02-16

General Chemistry 6

ALT
AST
BUN
CRE



07P02-01



General Chemistry 13

ALB
ALP
ALT
AMY
AST
BUN
Ca

CRE
GGT
GLU
TBIL
TP
UA
eGFR*



07P02-11



Electrolyte Panel

Cl⁻
K⁺
Na⁺
tCO₂



07P02-02



Kidney Check

BUN
CRE
eGFR*



07P02-09



Renal Function Panel

BUN
Ca
Cl⁻
CRE
GLU

K⁺
Na⁺
PHOS
tCO₂
eGFR*



07P02-07



MetLyte 8 Panel

BUN
CK
Cl⁻
CRE
GLU

K⁺
Na⁺
tCO₂
eGFR*



07P02-03



*Calculated. Note: Each panel comes in a package of 10.

1. Crocker, J et al, Implementation of Point-of-Care Testing in an Ambulatory Practice of an Academic Medical Center, Am J Clin Pathol November 2014;142:640-646 2. Actual savings can vary depending on state-specific factors.



REAL-TIME CHEMISTRY. REAL-WORLD BENEFITS.

- Optimize efficiency and productivity with fast, accurate, on-site testing
- Make decisions with confidence with lab-accurate results
- Maintain compliance with automated quality control features
- Increase patient satisfaction by reducing the need for call-backs
- Meet a range of testing needs with a single, versatile platform

BROAD CLINICAL APPLICATIONS.

Piccolo Xpress® meets a wide range of clinical testing needs across a variety of healthcare settings and specialties.

- Urgent Care
- Internal Medicine
- Cardiology
- Hospital Clinics
- Primary/Family Care
- Long-Term Care
- Oncology
- Pediatrics
- Ambulatory Surgery
- Dermatology

piccolo *xpress*
chemistry analyzer



Scan this code to see the Piccolo Xpress in action.

To learn more about how the Piccolo Xpress can transform your practice, contact your Abbott Point of Care Representative, or visit www.abbottpointofcare.com

For in vitro diagnostic use only © Abbott Point of Care Inc. 046406 Rev. A 07/16. Piccolo Xpress®, Piccolo® and iQC® are registered trademarks of Abaxis, Inc. Abaxis Part Number 888-3293 Rev. A.

CHOOSE TRANSFORMATION™

See Inspiration in Action at abbottpointofcare.com

