

P.B.S.® Diamond Burs

Nail Reduction & Tissue Debridement

Professional treatment system for hyperkeratotic tissue and hypertrophic nails. The exclusive P.B.S.® bonding system permanently bonds diamonds to the bur surface, providing unsurpassed product life and performance. P.B.S.® burs prove more efficient than steel or carbide instruments and generate less heat for greater patient comfort.

All diamond burs fitted with $\frac{3}{32}$ " shaft size.



Barrel 2026005
(To reduce thick, deformed nails.)



Short 2026020
(For better control in treating smaller digits.)



Umbrella 2026035
(For the reduction of hyperkeratoses, tylomas and heel callus.)



Pear 2026010
(To smooth the rough edges of nails.)



Olive 2026025
(Similar to the PEAR and BUD.)




Small Olive 2026040
(For working smaller areas.)



Bud 2026015
(The most popular pattern for nail reduction.)



Taper 2026030
(For the treatment of incurvated nails or the nail groove.)



Small Cylinder 2026045
(For gross reduction of large nail areas.)

Diamond Bur System Kit 2026050

A convenient kit of four nail and callus skin burs plus technique information.

Includes one each:
Bud, Short, Taper, Umbrella

* P.B.S. is a registered trademark of the product manufacturer, Abrasive Technology, Inc.