



PRO 760HS SERIES

Our most popular aneroid, the 760HS is the model of choice for demanding institutions.

The 760HS series features:

- Contemporary black enamel 300mmHg CE no-pin stop manometer, provides day-in, day-out dependability
- Nylon cuff with proprietary Size Guide™ marking system, prevents mis-cuffing
- PVC latex free inflation bladder withstands heavy use
- Chrome-plated brass air release valve
- Leatherette carrying case with durable nylon zipper
- Assembled, inspected, and packaged in the U.S.A. from foreign components



PRIVATE LABELING

COMPONENT	PACKAGING
Gauge & Cuff	1 Pc. Box/Label

Feature

Benefit

Nissei Made Aneroid Gauge EXCLUSIVELY - The MOST experienced and respected Asian Gauge manufacturer	Unsurpassed accuracy, quality, and reliability. Our gauges are randomly tested to 30,000 cycles, 3x MORE than required under the National Standard (ANSI SP 10)
Gauge is 100% visually inspected for 0 point, 5% re-inspected for accuracy.	To ensure conformance with National Standard ANSI SP 10
Accurate to +/-2mmHg	Exceeds industry requirements by 33%
Gauge features extended graduation marks to prevent parallax error	Makes it easier to read from oblique angles
Cuff has a properly positioned artery mark at the bladder midpoint	Ensures proper compression of the brachial artery. Most low cost imports fail to meet this requirement from ANSI SP 10
Size Guide marking system (points out when to choose larger or smaller cuffs, properly positioned artery mark)	Prevents mis-cuffing which can lead to measurement error
Leak tested Valve	Valves are tested to ensure conformance with ANSI SP 10 for dump rate and leakage.
Premium filter equipped valve and end valve	Large valve is more ergonomic, filters prevent damaging dust build up
100% inspected, assembled and packaged in the USA	Ensures unmatched consistency, quality, and performance in a value priced BP instrument
20 year calibration, 3 year system warranty	To protect the consumer's investment

ORDERING INFORMATION

Schein Item#	ADC Item#	Product Description	Latex Content	U/M	Cs Pk
1126071	760-10SABKHS	PRO 760, Small Adult, Black	Latex-Free	ea.	20
1126072	760-11ABKHS	PRO 760, Adult, Black	Latex-Free	ea.	20
1126073	760-9CBKHS	PRO 760, Child, Black	Latex-Free	ea.	20
1126074	760-13TBKHS	PRO 760, Thigh, Black	Latex-Free	ea.	12
1126075	760-12XBKHS	PRO 760, Lrg. Adult, Black	Latex-Free	ea.	20