



PRO PLUS 700HS

Our finest pocket aneroid, and it is the instrument of choice for the discriminating professional.

The PRO PLUS 700HS series features:

- Chrome-plated precision crafted 300mmHg CE no-pin stop manometer
- Luminescent dial features bold numerals and extended graduation marks for easy reading in all light conditions and from virtually any angle
- Nylon cuff with proprietary Size Guide™ marking system, prevents mis-cuffing
- PVC Latex free inflation bladder withstands heavy use
- Latex free inflation bulb fitted with filter screen end valve to reduce dust build-up
- Chrome-plated brass air release valve with filter screen
- Leatherette carrying case with durable nylon zipper
- Assembled, inspected, and packaged in the U.S.A. from foreign components

PRIVATE LABELING	
COMPONENT	PACKAGING
Gauge & Cuff	2 Pc. Box/Label



Feature	Benefit
Nissei Made Aneroid Gauge EXCLUSIVELY - The MOST experienced and respected Asian Gauge manufacturer	Unsurpassed accuracy, quality, and reliability. Our gauges are randomly tested to 30,000 cycles, 3x MORE than required under the National Standard (ANSI SP 10)
Gauge is 100% visually inspected for 0 point, 5% re-inspected for accuracy.	To ensure conformance with National Standard ANSI SP 10
Accurate to +/-2mmHg	Exceeds industry requirements by 33%
Gauge features extended graduation marks to prevent parallax error	Makes it easier to read from oblique angles
Gauge is chrome plated	More durable finish
Gauge features Luminescent dial	Easier to read in low light conditions
Cuff has a properly positioned artery mark at the bladder midpoint	Ensures proper compression of the brachial artery. Most low cost imports fail to meet this requirement from ANSI SP 10
Size Guide marking system (points out when to choose larger or smaller cuffs, properly positioned artery mark)	Prevents mis-cuffing which can lead to measurement error
Unit features Latex Free bulb and bladder	Eliminates risk of latex reaction
Leak tested Valve	Valves are tested to ensure conformance with ANSI SP 10 for dump rate and leakage.
Premium filter equipped valve and end valve	Large valve is more ergonomic, filters prevent damaging dust build up
100% inspected, assembled and packaged in the USA	Ensures unmatched consistency, quality, and performance in a value priced BP instrument
Lifetime calibration, 3 year system warranty	To protect the consumer's investment

ORDERING INFORMATION					
Schein Item#	ADC Item#	Product Description	Latex Content	U/M	Cs Pk
1126064	700-10SABKHS	PRO PLUS 700, Sm. Adult, Black	Latex-Free	ea.	20
1126065	700-11ABKHS	PRO PLUS 700, Adult, Black	Latex-Free	ea.	20
1126066	700-12XABKHS	PRO PLUS 700, Lg. Adult, Black	Latex-Free	ea.	20