

Upgrade to the TRUE METRIX[®] PRO

Professional Monitoring Blood Glucose System

The TRUE METRIX[®] portfolio* achieves a level of performance driven by science, research and technological advancements.

Featuring  **TRIPLE SENSE TECHNOLOGY**[®]

Each system automatically detects, analyzes and corrects variables in each blood sample to ensure proven accuracy and confidence in results.

Proven Accuracy

The TRUE METRIX[®] Systems exceed the new, more stringent minimum ISO 15197:2013 performance requirements.¹

With TRUE METRIX[®]—
Accuracy is Confidence



Products	MFR Item #	HSI Item #
Professional Monitoring Blood Glucose System (Intended for multiple-patient use in professional healthcare settings)		
TRUE METRIX[®] PRO-Professional Monitoring Blood Glucose Meter Kit (Includes: Meter, Installed 3-Volt Battery, 10 ct. Test Strips, Instructions for Use, Log Book, Compact Carrying Case, Control Solution)	RE4099P-43	570-0329
TRUE METRIX [®] PRO Professional Monitoring Blood Glucose Test Strips, 50 ct.	R3099P-450	570-0327
Self-Monitoring Blood Glucose System		
TRUE METRIX[®] Self-Monitoring Blood Glucose Meter Kit, Not For Retail Sale (NFRS) (Includes: Meter, Installed 3-Volt Battery, 10 ct. Test Strips, 10 Lancets and Lancing Device, Instructions for Use, Log Book, Compact Carrying Case, Control Solution)	RE4099-45	570-0330
TRUE METRIX [®] Self-Monitoring Blood Glucose Test Strips, 50 ct.	R3099-450	570-0326
TRUE METRIX [®] Self-Monitoring Blood Glucose, 50 ct. Medicare/Medicaid	R3H01-350	125-4037
Control Solution		
TRUE METRIX [®] Control Solution Level 1	R5H01-1	125-4034
TRUE METRIX [®] Control Solution Level 2	R5H01-2	125-4035
TRUE METRIX [®] Control Solution Level 3	R5H01-3	125-4036
Lancets		
TRUEplus [®] Lancets 28g, 100 ct. (Single-use)	743500	890-0175

* TRUE METRIX[®] is intended for self-monitoring blood glucose only and not for multiple patient use. Only TRUE METRIX[®] PRO is intended for multiple patient use.

Features:

- No coding
- Tiny 0.5 microliter sample size
- As fast as 4 seconds
- Stores 500 results
- Data management
- Download capabilities



All test strips are made in the U.S.A.

Extensive education, training and support for you and your patients

Your facility can be the leader in advanced education, service and products with the TRUE METRIX® PRO Support Program



24/7 Customer Care

- Product experts available 24/7 to answer all your patients' product questions – in both English and Spanish – Call 1-800-803-6025



Field Sales Specialists

- Field Sales Specialists in your area are ready to provide in-person support and education on the advanced technology this system offers
- The team can also provide access to all patient and product support tools



Resource Guides

- Resource Guides are available for the system
- Overview of the meter and test strip, testing guide, troubleshooting guide, training programs and more



TRUEmanager® Diabetes Management Software


Track a Healthier Course to Better Diabetes Management

- Patients and healthcare professionals can access and analyze results from their meters
- View Record Summary, Log Book, Glucose Trend, and Conformance Reports



Online Support

- Find information on our advanced diabetes products, professional resources, and programs at trivdiahealth.com
- Easy, three-step training and certification program for you and your staff who test point-of-care blood glucose in multiple-patient facilities



Patient Education & Training Materials

- Educational brochures and product instruction guides are available for you and your patients
- Accessible online at trivdiahealth.com

NICO-3544 © 08/16 Trividia Health, Inc. TRUE METRIX, TRIPLE SENSE TECHNOLOGY, TRUEinsight, TRUEplus and TRUEmanager are trademarks of Trividia Health, Inc. All other products names and brand names are trademarks of their respective owners. 1. Data on file